



WORLD  
PRECISION  
INSTRUMENTS

ANIMA  LAB  
animal facility and laboratory equipment • animal research models

# LABORATORY EQUIPMENT AND SUPPLIES

catalog 2019

QUALITY COST-EFFECTIVE  
**EQUIPMENT**  
FOR WHEREVER  
SCIENCE TAKES  
**YOU**





WORLD  
PRECISION  
INSTRUMENTS



## Let us Be Your Partner in Discovery

We are at the dawn of an exciting new day of science—improving the lives of all humanity with breakthroughs like pheno-genomic testing and therapy. Daily, researchers like you are discovering new models, pathways and techniques which impact our health and our environment. These discoveries result in life-changing—and in many cases life-prolonging—therapies and revolutionary protocols. We are also working with researchers in monitoring environmental influences in our oceans and streams. New surgical instruments are being created to assist in specific surgeries and diagnoses. This year, for example, we added a new line of tools specifically for work under a microscope. They are non-reflective and are also modified with special tips to enable surgeries in arteries and veins.

At WPI, our solutions enable pioneers on the cusp of discovery to impact the quality of life for all of us for generations to come. WPI provides solutions for a wide range of disciplines, including electrophysiology, molecular biology, cellular biology, surgical procedures and physiology, to advance the life science knowledge base. As researchers continue the quest for answers to life's mysteries, we provide the instruments to make laboratory research effective, reliable and reproducible. We aim to be your partner in the scientific adventure and instill confidence in your journey of discovery. Your research methods and ideas, combined with our scientists and engineers, help us develop novel solutions for your work.

Our customers' journey of discovery, and the promise it brings of life-changing impact, inspires our team. Let us be your research partner with quality, cost-effective instruments for wherever science takes you.



# Table of Contents

<b>New Product Summary</b>				<b>2</b>
<b>Cell &amp; Tissue</b>				<b>3</b>
Epithelial Physiology Muscle Physiology Biofluorometer	Micromanipulators Mounting Devices Magnetic Holding Devices	Live Cell Imaging Vibration Isolation		
<b>Microinjection</b>				<b>39</b>
Microinjection Systems Microinjection Accessories	Pinpoint Injector Picoliter Injectors			
<b>Fluid Handling</b>				<b>49</b>
Microperfusion Peristaltic Pumps UltraMicroPumps	Syringes UMP3 Application Nanoliter Injectors	Syringe Pumps Microfluidics Pipettors	Pipetter Tips Vacuum Pumps Luer Fittings	
<b>Electrophysiology</b>				<b>73</b>
Amplifiers Stimulators Stimulus Isolators	Electrophysiology Accessories Microelectrodes Wire, Cables & Connectors	Pipette Holders Capillary Glass Specialty Wire & Half Cells	Micropipettes Pullers, Microforges Micropressure Systems	
<b>Animal Physiology</b>				<b>117</b>
Stereotaxic Anesthesia Blood Pressure Monitoring	Thermometers Temperature Control Neuroscience Cannulas	Analgesia Adhesives		
<b>Biosensing</b>				<b>143</b>
Biosensors Free Radical Analyzers Oxygen Detection Systems	Data Acquisition pH Detection Biosensing Accessories			
<b>Optical Measurement Systems</b>				<b>161</b>
Qualified Optic Cables Optical Detection Flow Cells	Sample Holders Dipping Probes Cuvettes	Spectrometers Light Sources Biofluorometer	Microscopy & Cameras Illumination	
<b>Microdissection Instruments</b>				<b>192</b>
Forceps Scissors Needle Holders	Black Instruments Punches Clamps, Scalpels	Rodent Accessories Microdissection Kits Surgical Equipment	Care and Cleaning	
<b>Lab Supplies</b>				<b>212</b>
Adhesives Glass Capillaries Micropipettes	Luer Fittings BNC Cables & Connectors Wire	Bases & Mounts Calipers & Micrometers Pipettors & Tips	Illumination Culture Dishes	
<b>Index</b>				<b>226</b>

All specifications are subject to change.  
For research only.

Copyright © 2019, World Precision Instruments, Inc. All rights reserved.



# Products at WPI

## PINPOINT INJECTOR



Preliminary Design

The new **MICRO-ePORE™** cell penetrator is a simple and versatile pin-point controlled electroporation system that can be used for efficient microinjection of reagents into oocytes and pre-implantation stage mammalian embryos. Patent pending Flutter Electrode Technology assists in small, clean, precise membrane penetration without tearing or damaging the membrane. See page 42.

## QUALIFIED OPTICAL FIBERS

The use of qualified optical fibers ensures better fibers, better assemblies, better probes. WPI excels at making optical fibers with high batch-to-batch performance consistency. WPI UV fiber has industry leading degradation resistance to deep UV light. WPI fiber composition and manufacturing is unique so that we can offer you qualified and guaranteed UV performance. Our optical fiber assemblies offer excellent optical performance and low lot-to-lot variation for scientific and process applications. Our manufacturing know-how and optical designs give you the ability to order custom assemblies designed to meet the specifications of your application. See page 162.



## FLUORESCENCE PROBE

WPI offers fluorescence probes to match your application, both single and double emission. These probes may be used for the detection of the transient response of free ion concentrations, like calcium, potassium, sodium or magnesium. Fluorescence probes can also be designed to detect pH or membrane potential. Auto-fluorescence like the detection of ATPase activity via NADH or FAD is another application. For details, talk with a WPI application specialist. See page 178.



## OPTICAL OXYGEN SENSOR



Now you can detect oxygen in small sample in living tissues using the phase shift between a reference signal and a measured signal. **BioOxy** is a new and innovative technology for measuring oxygen in gaseous and aqueous phase. **BioOxy** is an optical oxygen sensor with important advantages over using common Clark type electrodes. These sensors have a rapid response time with the  $t_{90} \leq 3$  seconds (in gas phase). No oxygen is consumed making measurements. These sensors are perfect for bioprocess control and making oxygen measurement in small samples. See page 152.

## TEER MEASUREMENT CHAMBERS



ENDOHM-6G

ENDOHM-12G

ENDOHM-24GSNAP

Our new EndOhm chambers are now made of glass, and the new insert holder with 120° tri-supports can accommodate three leg inserts! The crystal clear glass chamber allows visualization of apical electrode positioning. See page 5.

## STERILIZATION TRAYS AND BASKETS

With top quality structural integrity, these trays are ideal for the handling and storage of all standard microdissection and surgical instruments. The unique "grid" design of the base makes it easy to install the silicone finger mats required to protect the delicate instruments. They have burr-free edges, and we offer a large selection to fit your needs. See [www.wpiinc.com/baskets](http://www.wpiinc.com/baskets). See page 211.



## OPHTHALMIC KITS AND INSTRUMENTS

We offer kits for specific research and veterinary procedures. And we have over 100 ophthalmic instruments each sold separately. See [www.wpiinc.com/eyes](http://www.wpiinc.com/eyes). See page 198.



503364

## CERAMIC COATED INSTRUMENTS

Our black instruments are coated with titanium nitride (TiN), an extremely hard ceramic material. The TiN coating hardens and protects the cutting edge. Ceramic coated anti-reflective instruments are perfect for microscopy and microsurgical applications. Coating surgical instruments with a black ceramic adds a thin layer to the metal instrument, making the instrument harder and giving you greater precision. This anti-glare surface minimizes reflections off the surface of your instruments. The incredibly smooth coating improves the instruments' resistance to corrosion and minimizes friction. The ceramic coating is virtually impenetrable, because the raw material is bonded to the instrument both physically and chemically. These instruments are much more resilient to the pressure of daily use and chemical processing. Coated instruments last considerably longer. See page 203.



504487

## STERI-LITE MICROBEAD STERILIZER

Ideal for the sterilization of small research instruments, the chamber is filled with glass beads and heated up to 300°C. Eliminate bacteria, spores and other microorganisms. Chamber temperature is displayed on the large LED control panel and may be adjusted with the control knob. See page 210.



ST5191

### WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • [wpiuk@wpi-europe.com](mailto:wpiuk@wpi-europe.com) • [www.wpi-europe.com](http://www.wpi-europe.com)  
Brazil: 011 55 13 40629703 • [info@brazil.wpiinc.com](mailto:info@brazil.wpiinc.com) • [www.wpiinc.com](http://www.wpiinc.com)

Germany: +49 (0)30 61 88845 • [wpide@wpi-europe.com](mailto:wpide@wpi-europe.com) • [www.wpi-europe.com](http://www.wpi-europe.com)  
China: +86 21 6888 5517 • [chinasales@china.wpiinc.com](mailto:chinasales@china.wpiinc.com) • [www.wpiinc.net](http://www.wpiinc.net)

# Cell & Tissue



## Comprehensive Products for Cell & Tissue Applications

WPI's cell and tissue products cover a broad range of applications. Our products are tested and optimized to provide you with instruments that are reliable in supporting your research requirements. We offer a range of products for TEER measurement, live cell imaging, microscopy studies and muscle physiology applications. WPI's products have been cited in 100s of reference papers and in many cases are the only products available of their kind. WPI's **EVOM2™** is the original instrument designed specifically to perform non-destructive TEER measurement on epithelial monolayer cell cultures.

# Epithelial Volt/Ohm (TEER) Meter

## Non-destructively test for epithelial monolayer confluence in 2D cell cultures

### Features

- Measures trans-epithelial electrical resistance or trans-epithelial voltage
  - Compatible with 12 and 24 well culture plate systems out of the box
- Includes industry standard **STX2** hand held "chopstick" electrodes

- Analog output for recording resistance or voltage measurements
- Auto ranging from 0-10 K $\Omega$
- Battery powered

### Benefits

- You can verify performance and calibrate the meter for TEER function using provided test resistor
- Battery powered meter is portable

A variety of accessory electrodes are available for measuring TEER in 6- and 96-well fixed (HTS) and removable well culture systems (See **STX100** series (page 6) and **Endohm** (page 5) electrodes)

### Applications

- TEER and trans-epithelial voltage measurements in 2D cell cultures

The **EVOM™** was the first instrument designed specifically to perform routine Trans Epithelial Electrical Resistance (TEER) measurement in tissue culture research. **EVOM2™** is the next generation, redesigned for ease of use. The **EVOM2™** not only qualitatively measures cell monolayer health, but also quantitatively measures cellular confluence. The unique electronic circuit of the **EVOM2™** and the included **STX2** electrode detect the confluence of the cellular monolayer. When combined with WPI's Endohm chamber, the **EVOM2™** can also be used to perform more accurate quantitative measurements or lower resistance measurements like transendothelial electrical resistance measurements.

#### Isolated battery power for 10 hours of use


The isolated power source of the **EVOM2™** was specifically designed to avoid adverse effects on tissue and the formation of electrode metal deposits, even when it is plugged into a standard wall outlet. Now, the **EVOM2™** is always on when you need it. In addition, its rechargeable battery allows up to 10 hours of mobile use.

#### Accurate reading every time

The four-and-a-half digit readout provides a range of 1-10,000  $\Omega$ . The included test electrode lets you calibrate the resistance measurements for an accurate reading every time, and the voltage meter never needs calibration. An analog BNC output is standard with the **EVOM2™**, providing an output port for recording data or remote display of the **EVOM2™** output.

#### Electrode pair to measure voltage, pass current

**EVOM2™** comes complete with the popular **STX2** "chopstick" electrodes, 4 mm wide and 1 mm thick. Each stick of the electrode pair contains a silver/silver-chloride pellet for measuring voltage and a silver electrode for passing current. The small size of each electrode is designed to facilitate placement of the electrodes into a variety of standard cell culture wells.



### References

- Di, S., Gujie, M., & Thomas, W. (2016). Magnetic ferri-liposomes for triggered drug release across the blood-brain barrier. *Frontiers in Bioengineering and Biotechnology*, 4. <http://doi.org/10.3389/conf.FBIOE.2016.01.00061>
- Paradis, A., Leblanc, D., & Dumais, N. (2016). Optimization of an *in vitro* human blood-brain barrier model: Application to blood monocyte transmigration assays. *MethodsX*, 3, 25–34. <http://doi.org/10.1016/j.mex.2015.11.009>



EVOM2™

Meenach, S. A., Tsoras, A. N., McGarry, R. C., Mansour, H. M., Hilt, J. Z., & Anderson, K. W. (2016). Development of three-dimensional lung multicellular spheroids in air- and liquid-interface culture for the evaluation of anticancer therapeutics. *International Journal of Oncology*, 48(4), 1701–9. <http://doi.org/10.3892/ijo.2016.3376>

Turdaliev, A., Solandt, J., Shambetova, N., Xu, H., Blom, H., Brismar, H., ... Fu, Y. (2016). Bioelectric and Morphological Response of Liquid-Covered Human Airway Epithelial Calu-3 Cell Monolayer to Periodic Deposition of Colloidal 3-Mercaptopropionic-Acid Coated CdSe-CdS/ZnS Core-Multishell Quantum Dots. *PLoS One*, 11(2), e0149915. <http://doi.org/10.1371/journal.pone.0149915>

Hirano, M., & Hirano, K. (2016). Myosin di-phosphorylation and peripheral actin bundle formation as initial events during endothelial barrier disruption. *Scientific Reports*, 6(1), 20989. <http://doi.org/10.1038/srep20989>

Marvin, S. A., et al. (2016). Type I Interferon Response Limits Astrovirus Replication and Protects against Increased Barrier Permeability In Vitro and In Vivo. *Journal of Virology*, 90(4), 1988–1996. <http://doi.org/10.1128/JVI.02367-15>

Jacobs, M. E., et al. (2016). SGK1 regulation by miR-466g in cortical collecting duct cells. *American Journal of Physiology-Renal Physiology*, 310(11), F1251–F1257. <http://doi.org/10.1152/ajprenal.00024.2016>

Yu, C., et al. (2015). The Effects of Glucagon-like Peptide-2 on the Tight Junction and Barrier Function in IPEC-J2 Cells through Phosphatidylinositol 3-kinase-Protein Kinase B-Mammalian Target of Rapamycin Signaling Pathway. *Asian-Australasian Journal of Animal Sciences*, 29(5), 731–738. <http://doi.org/10.5713/ajas.15.0415>

### SPECIFICATIONS

	EVOM2™	300523
MEMBRANE VOLTAGE RANGE	±200 mV	±200 mV
RESOLUTION	0.1 mV	0.1 mV
RESISTANCE RANGE	0 to 9999 $\Omega$	0 to 100 K $\Omega$
RESISTANCE RESOLUTION	1 $\Omega$	10 $\Omega$
AC SQUARE WAVE CURRENT	±10 $\mu$ A nominal at 12.5 Hz	1 $\mu$ A nominal at 12.5 Hz
POWER	Internal rechargeable 6V NiMH 2700 mAh battery with external 12 VDC supply for recharging	
NOMINAL BATTERY RUN TIME	10 hours	
BNC OUTPUT	1-10 V (1 mV/ $\Omega$ )	1-10 V (1 mV/10 $\Omega$ )
DIMENSIONS	19 x 11 x 6 cm (7.25" x 4.25" x 2.30")	
WEIGHT	1.4 kg (3 lb.)	
ELECTRODE CONNECTION	RJ-11 connector (telephone style)	
TEST RESISTOR	External, 1000 $\Omega$	
ENVIRONMENTAL RANGE	10-38°C (50-100°F)	
	0-90% non-condensing relative humidity	

### ORDERING INFORMATION

<b>EVOM2</b>	Epithelial Tissue Volt/ohmmeter 9999 $\Omega$ range (includes STX2)
<b>300523</b>	EVOM2™ 100 K $\Omega$ Range (includes STX2)
<b>91799</b>	Epithelial Tissue Volt/ohmmeter 9999 $\Omega$ range (includes STX3)

### OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>STX2</b>	Replacement "Chopstick" Electrode Set
<b>STX3</b>	Adjustable Tip Spread "Chopstick" Electrode Set
<b>3993</b>	Electrode Adapter (for electrodes with 2 mm pins)
<b>91736</b>	Replacement Battery, Rechargeable NiMH
<b>91750</b>	EVOM2™ Test Resistor

### WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • [wpiuk@wpi-europe.com](mailto:wpiuk@wpi-europe.com) • [www.wpi-europe.com](http://www.wpi-europe.com)  
Brazil: (013) 406-29703 • [info@brazil.wpiinc.com](mailto:info@brazil.wpiinc.com) • [www.wpiinc.com](http://www.wpiinc.com)

Germany: +49 (0)6031 1602171 • [wptide@wpi-europe.com](mailto:wptide@wpi-europe.com) • [www.wpi-europe.com](http://www.wpi-europe.com)  
China: +86 21 6888 5517 • [chinasales@china.wpiinc.com](mailto:chinasales@china.wpiinc.com) • [www.wpiinc.net](http://www.wpiinc.net)

# Reproducible Resistance of Endothelial Tissue

## TEER measurement of epithelial and endothelial cell cultures

### Features

- Compatible with original **EVOM™** and **EVOM2™** meters
- Adjustable apical electrode height
- Crystal clear glass chamber allows visualization of apical electrode positioning
- New insert holder with 120° tri-supports for three leg inserts
- Three sizes cover a range of well cup sizes from a variety of manufacturers

### Benefits

- Stability and reproducibility superior to the **STX2** electrodes to 1% tolerance
- Can be used with 6, 12 or 24 well plates with removable inserts
- Symmetrical electrode pattern disperses test current uniformly
- Simple test procedure to verify electrode performance

### Applications

- TEER measurement for removable culture cup systems using **EVOM2™** meters for endothelial and epithelial cell cultures

Using WPI's **EVOM2™** resistance meter, Endohm chambers provide reproducible resistance measurements of endothelial and epithelial monolayers in culture cups. Transfer cups from their culture wells to the Endohm chamber for measurement rather than using hand-held electrodes. The chamber and the cap each contain a pair of concentric electrodes: a voltage-sensing silver/silver chloride pellet in the center plus an annular current electrode. The height of the top electrode can be adjusted to fit cell culture cups of different manufacturers.

### Make more precise measurements with Endohms

Endohm's symmetrically opposing circular disc electrodes, situated above and beneath the membrane, allow a more uniform current density to flow across the membrane than with **STX2** electrodes. The background resistance of a blank insert is reduced from 150 Ω (when using WPI's hand-held **STX2** electrodes) to less than 5 Ω. With Endohm's fixed electrode geometry, variation of readings on a given sample is reduced from 10-30 Ω with **STX2** electrodes (depending on your experience) to 1-2 Ω. Compared with other resistance measurement methods, Endohm with **EVOM2™** offers a much more convenient and economic solution to "leaky tissue" measurement. Because of the uniform density of the AC square wave current from **EVOM2™**, errors caused by electrode polarization or membrane capacitance are largely eliminated. Endohm together with **EVOM2™** offers the most accurate and economical endothelial ohm meter now available. To date, cups from Costar, Millipore, and Falcon have been tested. Endohm chambers may be sterilized with ETO, alcohol or a bactericide. Do not autoclave.

### References

Srinivasan, B., Kolli, A. R., Esch, M. B., Abaci, H. E., Shuler, M. L., & Hickman, J. J. (2015). TEER measurement techniques for in vitro barrier model systems. *Journal of Laboratory Automation*, 20(2), 107-26. <http://doi.org/10.1177/2211068214561025>

ENDOHM-12 COMPATIBILITY CHART				
Corning	Millipore	Membrane Diameter	Growth Surface Area	Membrane Pore Size
3401		12 mm	1.12 cm <sup>2</sup>	0.4 μm
3402	PITP01250	12 mm	1.12 cm <sup>2</sup>	3.0 μm
3403	PITP01250	12 mm	1.12 cm <sup>2</sup>	3.0 μm
3493		12 mm	1.12 cm <sup>2</sup>	0.4 μm
3494		12 mm	1.12 cm <sup>2</sup>	3 μm
3460	PIHT15R48* PET Insert	12 mm	1.12 cm <sup>2</sup>	0.4 μm
	PIRP15R48* PET Insert	12 mm	1.12 cm <sup>2</sup>	1 μm
3462	PISP15R48* PET Insert	12 mm	1.12 cm <sup>2</sup>	3 μm
	PIMP15R48* PET Insert	12 mm	1.12 cm <sup>2</sup>	5 μm
	PIEP30R48* PIEP15R48* PIEP15R48* PET Insert	12 mm	1.12 cm <sup>2</sup>	8 μm



ENDOHM-24 COMPATIBILITY CHART			
Corning	Millipore	Material	Pore Size (μm)
3407		Polycarbonate	0.4
3801		Polycarbonate	0.4
3802		Polycarbonate	3.0
3412	PIHT30R48*	Polycarbonate	0.4
3414		Polycarbonate	3.0
	PITT03050	Polycarbonate	3.0
3428		Polycarbonate	8.0
3450		Polyester	0.4
3452		Polyester	3.0
3491		Collagen	0.4
3492		Collagen	3.0
	PICMORG50	Organotypic Insert	0.4
	PIHA03050	HA Insert	0.45
	PIHP03050	PCF Insert	0.4
	PICM03050	HA Mixed Cellulose Esters	0.4
	PIHT30R48*	PET Insert	0.4
	PIRP30R48*	PET Insert	1.0
	PISP30R48*	PET Insert	3.0
	PIMP30R48*	PET Insert	5.0
	PIEP30R48*	PET Insert	8

ENDOHM-6 COMPATIBILITY CHART				
Corning	Millipore	Membrane Diameter	Growth Surface Area	Membrane Pore Size
3470		6.5 mm	0.33 cm <sup>2</sup>	0.4 μm
3472	PITP01250	6.5 mm	0.33 cm <sup>2</sup>	3 μm
3413	PCF Insert	6.5 mm	0.33 cm <sup>2</sup>	0.4 μm
3415	PITP 01250 PCF Insert	6.5 mm	0.33 cm <sup>2</sup>	3 μm
3421		6.5 mm	0.33 cm <sup>2</sup>	5 μm
3422	PIEP 01250 PCF Insert	6.5 mm	0.33 cm <sup>2</sup>	8 μm
3495	PISP12R48* PIHT12R48* PET Insert	6.5 mm	0.33 cm <sup>2</sup>	0.4 μm
	PIHA012 50 (HA Insert)	6.5 mm	0.33 cm <sup>2</sup>	0.45 μm
	PICM012 50 (CM Insert)	6.5 mm	0.33 cm <sup>2</sup>	0.4 μm
3496	PISP12R48* PET Insert	6.5 mm	0.33 cm <sup>2</sup>	3 μm
	PIRP12R48* PET Insert	6.5 mm	0.33 cm <sup>2</sup>	1 μm
	PIMP12R48* PET Insert	6.5 mm	0.33 cm <sup>2</sup>	5 μm
	PIEP12R48* PET Insert	6.5 mm	0.33 cm <sup>2</sup>	8 μm
	PIXP01250 PCF Insert	6.5 mm	0.33 cm <sup>2</sup>	12 μm
	PITP01250 PIHP 01250	6.5 mm	0.33 cm <sup>2</sup>	1.0 μm 3.0 μm

\* New Chambers Support Tri-Leg Inserts

ORDERING INFORMATION	
<b>ENDOHM-6G</b>	Endohm for 6 mm Culture Cup (24 wells per plate)
<b>ENDOHM-12G</b>	Endohm for 12 mm Culture Cup (12 wells per plate)
<b>ENDOHM-24GSNAP</b>	Endohm for 24 mm & Costar Snapwell™ Cup (6 wells per plate)
	Requires EVOM2™, EVOM™, EVOMX™ or Millicell ERS-2
<b>53330-01</b>	Replacement Endohm Cable

# HTS Electrodes for Use with EVOM™ & EVOM2™

For High Throughput Screening (HTS) cell culture filter plates

## Features

- Designed for use with 24-well HTS plate (Corning Costar and BD Falcon) and with 96-well plates (Millipore)
- Improved accuracy down to 5 Ω
- Sterilize with EtO, alcohol or bactericide

## Benefits

- Smaller tip size than the **STX2** electrode constructed for durability fits neatly into the keyhole shaped filter well
- Electrode design reduces chance of contamination

## Applications

- **STX-100** Electrodes are designed for TEER measurement in HTS culture plates using the **EVOM2™**

Semi-permeable HTS (high throughput screening) culture plates have become a standard tool for pharmaceutical and institutional research in epithelial transport. HTS culture plates have well cups which are bonded together into a single assembly and are not removable. This makes the plates ideal for automated applications, but imposes a significant inconvenience when TEER must be measured in the absence of an automated system. WPI's **STX100** series electrodes provide a cost-effective alternative to automation, allowing HTS well plates to be measured manually using a hand-held electrode.

The spatial orientation of an electrode during TEER measurement can have a significant effect on the resistance reading. When compared to the **STX-2** "chopsticks" electrode, the **STX100** series electrodes offer a technical advantage. The design of the **STX100** series electrodes guarantees spatial repeatability. Each **STX100** electrode is designed to self-align based on the form factor of the apical and basal access ports of the HTS plate. Each manufacturer of HTS culture systems has a unique form factor, so the **STX100** electrodes are manufacturer specific.

WPI developed **STX100** electrodes for 24 and 96 well plates specific to Corning, BD Falcon and Millipore. Refer to the charts for compatibility information or contact your WPI sales representative for assistance. If an automated system is preferable for your application, WPI's automated TEER measurement system (REMS) on page 7 is an affordable solution.



**STX100M**

### STX100 M COMPATIBILITY CHART

Millipore	Description	Pore Size	Membrane
PSRP004R1	96-Well Plate	1.0 μm	PET
PSHT004R5	96-Well Plate	0.4 μm	PCF
PSRP004R5	96-Well Plate	1.0 μm	PET
PSHT004S5	96-Well Plate	0.4 μm	PCF
PSHT004R1	96-Well Plate	1.0 μm	PCF

### STX100 C96 COMPATIBILITY CHART

Corning	Description	Pore Size	Membrane
3380	HTS Transwell-96 System	1.0 μm	PC
3392	HTS Transwell-96 System	1.0 μm	PET
3381	HTS Transwell-96 System	0.4 μm	PC
3391	HTS Transwell-96 System	0.4 μm	PC
3385	HTS Transwell-96 Well Plate	3.0 μm	PC
3386	HTS Transwell-96 Well Plate	3.0 μm	PC
3387	HTS Transwell-96 Well Plate	5.0 μm	PC
3388	HTS Transwell-96 Well Plate	5.0 μm	PC
3374	HTS Transwell-96 Well Plate	8.0 μm	PET
3384	HTS Transwell-96 Well Plate	8.0 μm	PET

### STX100 C COMPATIBILITY CHART

Corning	Description	Pore Size	Membrane
3379	HTS Transwell-24	0.4 μm	PET
3378	HTS Transwell-24	0.4 μm	PET
3396	HTS Transwell-24	0.4 μm	
3397	HTS Transwell-24	0.4 μm	
3398	HTS Transwell-24	3.0 μm	
3399	HTS Transwell-24	3.0 μm	

### STX100F COMPATIBILITY CHART

Falcon	Description	Pore Size	Membrane
351180	BD Falcon (24 well) HTS Multiwell Insert System	1.0 μm	PET
351181	BD Falcon (24 well) HTS Multiwell Insert System	1.0 μm	PET
351182	BD Falcon (24 well) HTS Multiwell Insert System	3.0 μm	PET
351183	BD Falcon HTS (24 well) Multiwell Insert System	3.0 μm	PET
351184	BD Falcon HTS (24 well) Multiwell Insert System	8.0 μm	PET
351185	BD Falcon (24 well) HTS Multiwell Insert System	8.0 μm	PET
354803	BD BioCoat (24 well) HTS Fibrillar Collagen Multiwell Insert System	1.0 μm	PET
354804	BD BioCoat (24 well) HTS Fibrillar Collagen Multiwell Insert System	1.0 μm	PC



### ORDERING INFORMATION

<b>STX100C</b>	STX100 for Corning Costar HTS Transwell-24
<b>STX100F</b>	STX100 for BD Falcon 24 well HTS Multiwell Insert System
<b>STX100M</b>	STX100 for Millipore Multiscreen™ HTS 96-Well Plate
<b>STX100C96</b>	STX100 for Corning HTS 96-Well Plate

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
Brazil: (013) 406-29703 • info@brazil.wpiinc.com • www.wpiinc.com

Germany: +49 (0)6031 1602171 • wvide@wpi-europe.com • www.wpi-europe.com  
China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net



# Automated TEER Measurement System

PC-controlled high throughput TEER measurement for epithelial monolayer

## Features

- Automates TEER measurement and data logging for use with HTS well plates
- PC controlled positioning
- Data acquisition in LabView
- Manufacturer specific electrodes available for 24- and 96-well HTS plates
- Plate configuration files and sample sequences are user-definable
- Two user-defined rinse locations
- Manual mode

## Benefits

- Speed—capable of acquiring TEER data on a 96-well plate in less than five minutes
- Automation reduces the possibility for human error

## Applications

- Automated measurement and data logging of TEER for 24 and 96 well HTS culture plates



SYS-REMS

The **REMS** AutoSampler automates measurements of TEER epithelial or endothelial monolayers cultured on HTS well plates. It is a PC controlled tissue resistance measurement system that offers reproducibility, accuracy, flexibility and ease-of-operation. Automated measurement of tissue resistance in cell culture microplates provides the advantages of speed, precision, decreased opportunity for contamination and the rapid availability of measured resistance data.

The main components of the **REMS** AutoSampler include: the robotic sampler that moves the electrode over each well of the microplate, the electrode which is located on the robotic arm, a base plate for the 24- and 96-well tray, a Windows-based data acquisition card, the **REMS** electrode interface unit and the **REMS** software to operate the system on a Windows-based computer.

### Automate TEER measurements

The **REMS** AutoSampler automates TEER measurements that would otherwise be performed manually with WPI's **EVOM2™** Epithelial VoltOhmmeter. Automated tissue resistance measurements up to 20 kΩ can be performed on 24- or 96-well HTS microplates. See [www.wpiinc.com/REMS](http://www.wpiinc.com/REMS) for manufacturer plate compatibility.

### Precisely and reproducibly positions electrode

The **REMS** AutoSampler will automatically measure and record tissue resistance from a user-specified matrix of culture wells on the microplate. According to the specified sequence, the robotic arm moves over the identified wells taking TEER measurements. By means of an x-y-z locating system, the electrode is positioned precisely into the well. The ability of the **REMS** AutoSampler to reproducibly position the electrode contributes to consistent TEER measurements. TEER data is incrementally stored as the electrode moves from one well to the next.

### Compact electrode pair

The use of AC current to measure resistance provides several advantages over DC current, including:

- Absence of offset voltages on measurements
- There is a zero net current being passed through the membrane and, therefore, it is not adversely affected by a current charge
- No electrochemical deposition of electrode metal.

### Rinse and calibration check stations

The **REMS** AutoSampler also features two rinse stations. If occasional rinsing of the **REMS** electrode is required, it may be sent to a rinse station by pressing the rinse station button on the menu bar.

## References

- Gallego, M., Grootaert, C., Mora, L., Aristoy, M. C., Van Camp, J., & Toldrá, F. (2016). Transepithelial transport of dry-cured ham peptides with ACE inhibitory activity through a Caco-2 cell monolayer. *Journal of Functional Foods*, 21, 388–395. <http://doi.org/10.1016/j.jff.2015.11.046>
- Kamiloglu, S., Grootaert, C., Capanoglu, E., Ozkan, C., Smagghe, G., Raes, K., & Van Camp, J. (2016). Anti-inflammatory potential of black carrot (*Daucus carota* L.) polyphenols in a co-culture model of intestinal Caco-2 and endothelial EA.hy926 cells. *Molecular Nutrition & Food Research*. <http://doi.org/10.1002/mnfr.201600455>
- Cain, M. D., Salimi, H., Gong, Y., Yang, L., Hamilton, S. L., Heffernan, J. R., ... Klein, R. S. (2017). Virus entry and replication in the brain precedes blood-brain barrier disruption during intranasal alphavirus infection. *Journal of Neuroimmunology*, 308, 118–130. <http://doi.org/10.1016/j.jneuroim.2017.04.008>

## REMS AUTOSAMPLER SPECIFICATIONS

MEMBRANE RESISTANCE RANGE	Auto-ranging or 0–2000 Ω and 0–20 kΩ
AC SQUARE WAVE CURRENT	± 20 μA @ 12.5 Hz
ELECTRODE POSITIONING	Resolution in X, Y and Z: ± 1 mm
ELECTRODE ARM SPEED	X- and Y-axis: 250 mm/sec Z-axis: 247.3 mm/sec
TYPICAL MEASUREMENT TIME	
24-WELL	1 min., 10 sec.
SCAN PATTERN	Preset or user-defined
LINE VOLTAGE (User Specified)	100/120 V or 220/240 V
DIMENSIONS	53.5 × 43.7 × 37.1 cm (21½ × 17¼ × 14¾ in.)
WEIGHT	24 kg (52 lb.)

## ORDERING INFORMATION

**SYS-REMS** Automated Tissue Resistance Measuring System

*Includes robot sampler, data acquisition board; computer, display, keyboard, mouse; software for Windows 7 or 10; and electrode (SPECIFY WHEN ORDERING).*

## OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>REMS-24</b>	Replacement Electrode for 24-well HTS Plate
<b>REMS-96</b>	Replacement Electrode for Millipore™ 96-well Plate
<b>REMS-96C</b>	Replacement Electrode for Corning 96-well HTS plate
<b>REMS-24M</b>	Replacement Electrode for Millipore 24-well HTS plate

*See [www.wpiinc.com/REMS](http://www.wpiinc.com/REMS) detailed information.*

# Ussing System for Epithelial Research

## Non-destructive TEER measurement for epithelial tissue

### Features

- Direct connect low-resistance electrodes
- Simple operation, easy to control temperature and clean after use
- Luer type leak-free attachment of tubing and electrodes
- Recessed electrode ports avoid bubble formation
- Secure membrane holding by sharp stainless steel pins or O-ring
- Specialized chamber adapts cell culture insert (Costar Snapwell) for monolayer cell culture
- Chambers with rectangular openings for tubular tissues from small animals

### Benefits

- Leak free design of Ussing chambers
- Can be used with monolayer cell culture inserts
- Optional drains for quick evacuation of radioactive or toxic substances
- Circulation reservoirs available in two sizes
- Control temperature with a circulating water bath (option) available

### Applications

- Ion transport studies
- Nutrient transport studies

WPI's Ussing System offers researchers a quick, effective means of making low-resistance electrical connections to the Ussing chamber without need of long agar bridges or Calomel half-cells. Ag/AgCl half-cells screw into short tubes which plug firmly into place in the chamber's Luer ports. These direct-connect electrodes eliminate the inconvenience and expense of Calomel half-cells in open liquids. The system includes one Ussing Chamber (eight sizes available), support stand, electrode kit, glass circulation reservoir (two sizes available), and a tubing start-up kit (25 feet of 0.375-in. tubing, 10 feet of 0.156-in. tubing, plus four male Luer fittings, two compressor clamps, one Y-connector, and one clip). Sixteen possible system configurations are available. Components are also available separately. (Preamplifier in photo not included.)

### Leak free design of Ussing chambers


WPI's classical Ussing Chambers are well established perfusion chambers that are easy to operate, easy to control temperature, and easy to clean after use. Ussing Chambers are machined from solid acrylic with eight entry ports for fluid lines, electrodes, or agar bridges. For easy, leak-free attachment of tubing and electrodes, all eight ports are Luer type. The four ports for voltage and current electrodes are recessed to prevent formation of air bubbles in the chamber. The fluid compartments in each side of the chamber are separated by the epithelial membrane being studied. Sharp stainless steel pins on one side of the chamber hold the membrane in position and mate with holes in the opposite chamber interface. (In the **CHM4**, tissue is held by an O-ring instead of pins.)

### Can be used with monolayer cell culture inserts

The **CHM5** chamber adapts the Costar 12 mm Snapwell, a cell culture insert for monolayer cell culture, into WPI's "classical" epithelial voltage clamp system. Classical Ussing Chambers have not been widely used for monolayer cell culture inserts, because most inserts have a very deep profile, limiting good fluid perfusion at the surface of the membrane and limiting voltage electrodes from measuring the potential close to the surface of the membrane. **CHM5** solves these problems: Perfusion fluid is introduced into the chamber at an angle so that it flows directly to the surface of the membrane. The voltage electrode is also inserted into the chamber at an angle to reduce the distance between the surface of the membrane and the electrode.

### Rectangular openings for tubular tissue

Two small chambers with rectangular openings are designed for tubular tissue from small animals such as the mouse intestinal tract membrane



Complete Ussing System includes stand, glass reservoir, electrodes, Ussing chamber and tubing.

(**CHM6**) and rat intestinal tract membrane (**CHM7**). The rectangular opening more closely matches the shape of the tissue than would a circular opening, significantly increasing the membrane area available for testing. The larger membrane area increases the transport rate of low permeability chemicals. It also reduces the electrical resistance of the system for easier current clamping.

### Optional drains

Drains may be added to Ussing chambers to allow quick and complete evacuation of radioactive or toxic substances. To have drains added at the time of order, add a "D" to the part number (such as "USS1LD"). The cost of the drain will be added to the cost of the chamber or system ordered.

### Cartridge electrodes

The Electrode Kit contains four voltage/current electrodes, plus four Luer-tipped cartridges. Electrodes are threaded and screw securely into the end of each cartridge. The Luer tip then plugs securely into the Luer openings of the chamber. The cable from each electrode terminates with a 2 mm pin which may be plugged into mos voltage/current clamps.

The miniature electrode-gel cartridge is a small plastic tube with a male Luer tip identical to those at the tip of hypodermic syringes. The tube may be filled with different gel materials. Agar is commonly used, but other gel materials may also be satisfactory.




EKV and EKC Cartridge Electrodes

US Patent No. 4,912,060

### WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
Brazil: (013) 406-29703 • info@brazil.wpiinc.com • www.wpiinc.com

Germany: +49 (0)6031 1602171 • wpide@wpi-europe.com • www.wpi-europe.com  
China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net



Assembled chambers are 101.6 mm (4 in.) long.

## References

M.Khera, G.T. Somogyi, S. Kiss, T.B. Boone, C.P. Smith "Botulinum toxin A inhibits ATP release from bladder urothelium after chronic spinal cord injury" *Neurochemistry International* 45, 2004: 987-993

McLamb, B. L., Gibson, A. J., Overman, E. L., Stahl, C., & Moeser, A. J. (2013). Early Weaning Stress in Pigs Impairs Innate Mucosal Immune Responses to Enterotoxigenic E. coli Challenge and Exacerbates Intestinal Injury and Clinical Disease. *PLoS ONE*, 8(4), e59838. <http://doi.org/10.1371/journal.pone.0059838>

Khera, M., Somogyi, G. T., Kiss, S., Boone, T. B., & Smith, C. P. (2004). Botulinum toxin A inhibits ATP release from bladder urothelium after chronic spinal cord injury. *Neurochemistry International*, 45(7), 987-993. <http://doi.org/10.1016/j.neuint.2004.06.001>

## ORDERING INFORMATION

### USSING SYSTEMS, LARGE RESERVOIR

USS1L	Medium Chamber, Stand, Reservoir, Electrodes, Tubing
USS2L	Small Chamber, Stand, Reservoir, Electrodes, Tubing
USS3L	Large Chamber, Stand, Reservoir, Electrodes, Tubing
USS4L	Extra Small Chamber, Stand, Reservoir, Electrodes, Tubing
USS5L	Snap Chamber, Stand, Reservoir, Electrodes, Tubing
USS6L	Small Rectangular Chamber, Stand, Reservoir, Electrodes, Tubing
USS7L	Large Rectangular Chamber, Stand, Reservoir, Electrodes, Tubing
USS8L	Extra Small Chamber, Stand, Reservoir, Electrodes, Tubing

### USSING SYSTEMS, SMALL RESERVOIR

USS1S	Medium Chamber, Stand, Reservoir, Electrodes, Tubing
USS2S	Small Chamber, Stand, Reservoir, Electrodes, Tubing
USS3S	Large Chamber, Stand, Reservoir, Electrodes, Tubing
USS4S	Extra Small Chamber, Stand, Reservoir, Electrodes, Tubing
USS5S	Snap Chamber, Stand, Reservoir, Electrodes, Tubing
USS6S	Small Rectangular Chamber, Stand, Reservoir, Electrodes, Tubing
USS7S	Large Rectangular Chamber, Stand, Reservoir, Electrodes, Tubing
USS8S	Extra Small Chamber, Stand, Reservoir, Electrodes, Tubing

### OPTIONAL ACCESSORIES/REPLACEMENT PARTS

System components also available separately:

xxxxD	Drain option (add "D" to part number of chamber or system)
CHM1	Medium Chamber
CHM2	Small Chamber
CHM3	Large Chamber
CHM4	Extra Small Chamber with O-Ring Seal
CHM5	Snap Chamber (fits Costar Snapwell cups)
CHM6	Small Rectangular Chamber
CHM7	Large Rectangular Chamber
CHM8	Extra Small Chamber with Mounting Pins
EK1	Ussing Electrode Kit (2 voltage, 2 current)
EKC	Extra Ussing Current Electrode (red) (each)
EKV	Extra Ussing Voltage Electrode (blue) (each)
5210	Large Glass Circulation Reservoir, (20 mL per side)
5233	Replacement Condenser for 5210
5362	Small Glass Circulation Reservoir, (10 mL per side)
5361	Replacement Condenser for 5362
3955	EKV Cartridges, 35 mm (pkg. of 12)
3960	EKC Cartridges, 58 mm (pkg. of 12)
3669	Tubing Kit (flexible hose and Luer fittings)
3579-20	Replacement Luer fittings for tubing connections (pkg. of 20)
5153	Support Stand
3485	Post Mounting Kit for Preamp
505060	Julabo Circulating Bath, 5L volume, 115V
505061	Julabo Circulating Bath, 13L volume, 115V
505062	Julabo Circulating Bath, 5L volume, 230V
505063	Julabo Circulating Bath, 13L volume, 230V



### Circulation reservoirs available in two sizes


Hand-blown borosilicate glass reservoirs with jacketed chambers for temperature control are available in two sizes — 5210 holds 20 mL per side, and 5362 (at left) holds 10 mL per side (useful when expensive chemicals are involved). Reservoir condenser caps prevent air bubbles and turbulence in fluid reservoirs.



### Control temperature with a circulating water bath

The Julabo circulating bath is ideal for controlling temperatures of external systems. With a powerful 15 L/min flow rate, the pump provides optimum heat exchange. The tap water cooling feature is standard with a range of 20-100°C.

### Reservoir openings and pins



## SPECIFICATIONS

	CHM1 (Medium)	CHM2 (Small)	CHM3 (Large)	CHM4 (Extra Small)	CHM5 (Snap)	CHM6 (Rect., Small)	CHM7 (Rect., Large)	CHM8 (Extra Small)
RESERVOIR OPENING	12 mm	9 mm	13.5 mm	4 mm	12 mm	5 x 14.5 mm	7 x 30 mm	4 mm
HALF-CHAMBER VOLUME	1.0 mL	0.75 mL	1.2 mL	0.5 mL	1.7 mL	0.8 mL	5.5 mL	0.5 mL
PIN CIRCLE DIAMETER	17 mm	12 mm	18.5 mm	6 mm*	N/A	7 x 16.5 mm	9 x 32 mm	5.5 mm
SURFACE AREA	113 mm <sup>2</sup>	63.5 mm <sup>2</sup>	143 mm <sup>2</sup>	12.6 mm <sup>2</sup>	113 mm <sup>2</sup>	67.1 mm <sup>2</sup>	199.5 mm <sup>2</sup>	12.6 mm <sup>2</sup>

\*O-ring diam.

# Muscle Testing Platform

Customizable system for muscle physiology research

## Features

- Heated cuvette and micrometer
- Study intact muscle fiber bundles\*
- Data recording and analysis included
- Modular design
- Corrosion-free

## Benefits

- System is completely customizable with a host of accessories

## Applications

- Electrical stimulation, including twitch and tetany, to analyze force response in amplitude; and kinetics, like contraction and relaxation times and velocities
- Intracellular calcium concentration/ distribution and muscle force in intact muscle fiber bundles.\* (Requires a Biofluorometer.)

These applications may require specific components or other electronics that are compatible with this system.

The **SI-MT-L** is the standard muscle research platform, which can be configured to study intact muscle fiber bundles, muscle strips and small whole muscles. Modular design allows the system to be configured for turnkey solutions for specific applications. The system is built on a solid platform making precise mechanical and optical measurements easy. Like the **Cell Tester** and the **SI-HTB** systems, **SI-MT-L** and **SI-MTM-L** use **SI-KG** optical force transducers and are constructed with corrosion-free materials (stainless steel, anodized aluminum or plastic).

The system is supplied with a **LabTrax-MDAC** data acquisition system to:

- Record signals from the force transducer, the motor position monitor or any further relevant physiological signal using up to eight analog inputs.
- Control the stimulator, or an external stimulus isolator, through an analog output
- Control the position of the linear motor through a second analog output
- Control other devices through two digital outputs
- Program numerous customized protocols


## Customize the system to suit your needs

The flexibility of the Muscle Testing Platform allows for complete customization. For example, you can add:

- Linear motor for muscle length perturbation studies
- Isolated constant current stimulator for twitch and tetanus stimulation
- Biofluorometer for intracellular calcium studies
- Additional accessories like force transducers, data acquisition system



SI-MT-L



MDAC Protocol segment, showing the combination of basic signal patterns to generate an experimental protocol.

## ORDERING INFORMATION

### SYSTEMS AVAILABLE FOR INTACT MUSCLE SAMPLES INCLUDE:

- SI-MTM-L** System with a Univette heated cuvette & linear motor
- SI-MT-L** System with a Univette heated cuvette & digital micrometer

Systems include: Base plate; force transducer and stand, digital micrometer or linear motor (optional), micrometer or motor stand, cuvette and table, oxygenation system, evacuation system, signal conditioning system with transducer amplifier and temperature control module, anti-oscillation module and linear motor amplifier (if applicable), LabTrax-MDAC data acquisition system with MDAC software for recording force and position signals and for controlling stimulation and the movement of the linear motor.

\*Contact WPI for help evaluating your specific experimental requirements.

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
Brazil: (013) 406-29703 • info@brazil.wpiinc.com • www.wpiinc.com

Germany: +49 (0)6031 1602171 • wpide@wpi-europe.com • www.wpi-europe.com  
China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net

# Horizontal Tissue Bath

Research system for higher throughput of complex pharmacological/physiological assays

## Features

- Two channel system for increased productivity, easily expanded to add channels
- Fully independent heating and fluid control for each channel
- Low profile/small footprint
- Variable volume, chemically inert Teflon bath with shape configurations from variable to fixed 500  $\mu$ L–10 mL\*
- Modular, space-saving, blade-style electronics
- Large variety of force transducers covering mN–N forces
- Can be combined with automated fluid control systems
- Add an electrometer like the WPI Duo773 or Electro705

## Benefits

- Tissue bath for volumes as small as 500  $\mu$ L
- Low profile, small footprint to reduce bench space requirements
- Versatile system is easy to upgrade with added channels and options

## Applications

- Tissue fluorescence (calcium, ATPase via NADH, FAD). (Requires the multi-purpose Biofluorometer\* (SI-BF-100). See page 20.)

These applications may require specific components or other electronics that are compatible with this system.

The SI-H Horizontal Tissue/Organ Bath system (SI-HTB) combines the ease of use and productivity of a traditional vertical organ bath with the more advantageous features of single tissue physiology platforms.

### Bath options 500 $\mu$ L up

The SI-HTB system breaks through the (large) volume limitations of the traditional organ bath, allowing volumes as low as 500  $\mu$ L in an inert, Teflon-based bath. The bath design allows multiple shape options for thick, long, flat and thin tissue. When pharmaceuticals are available in precious, small amounts, you will appreciate this standard feature. A wide range of transducers and tissue mounting supports allow a great variety of tissue shape, volume and size.


### Low profile

The low profile and small footprint of the bath system, combined with the modular, space-saving, chassis-mounted design of the electronics, reduces the bench space requirement up to 4-fold when compared with standard 4-channel organ baths.

### Versatile system, easy to upgrade

The SI-HTB combines advanced physiological techniques with the throughput needed in pharmacological assays in one flexible platform. Upgrades to four or more channels are easy and economical.

The motor option (SI-MOT) turns your system into a tissue work-out station with isotonic and eccentric force measurement capabilities. Nearly all established myo-mechanical tests from stretch-release to work-loops and muscle fatigue are now possible in a single organ bath system. Some of these procedures require automated length changing methods, which can be conveniently programmed using WPI's LabTrax-MDAC data acquisition system.



## Options

The solid horizontal tissue bath design is ideal for combination with electrophysiology on the same platform. Intracellular measurements can share the stable solid base of the bath system.

WPI's new fiber-optic based, multi-channel Biofluorometer allows for tissue fluorescence measurements (calcium, NADH, FAD) on the SI-HTB platform.

**Now, you can design a system to meet your needs and budget. And, it is fully upgradeable in the future.**

WPI's 16-bit, full speed, Labview-based Muscle Data Acquisition system LABTRAX-MDAC is perfect for this platform.

\*Contact WPI for help evaluating your specific experimental requirements.

## ORDERING INFORMATION

**SI-HTB2** Horizontal Tissue Bath, 2-Channel System  
2-Channel SI-HTB platform for isometric force (1), SI-KG Force Transducers (2), SI-BAM21-LCB Optical Force Transducer Amplifiers (2), SI-TCM2B 2-Channel Temperature Controller (1), SI-BMFA power frame enclosure (1), LabTrax 8/16 with MDAC Data Acquisition software

**SI-HTB4** Horizontal Tissue Bath, 4-Channel System  
2-Channel SI-HTB platform for isometric force (2), SI-KG Force Transducers (4), SI-BAM21-LCB Optical Force Transducer Amplifiers, (4) SI-TCM2B 2-Channel Temperature Controller (2), SI-BMFA power frame enclosure (1), LabTrax 8/16 with MDAC Data Acquisition software (1)

## OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>LABTRAX-MDAC</b>	LabTrax 8/16 8-Channel Data Acquisition with MDAC Software
<b>SI-AOSUB</b>	Anti-Oscillation Module
<b>SI-PF100</b>	Programmable Filter
<b>MINISTAR</b>	Miniature Peristaltic Pump, 1-channel
<b>PERIPRO-2HS</b>	Peri-Star™ Pro, 2-channel, High Rate, Large Tubing (110-220V)
<b>PERIPRO-4HS</b>	Peri-Star™ Pro, 4-channel, High Rate, Large Tubing (110-220V)
<b>PERIPRO-4LS</b>	Peri-Star™ Pro, 4-channel, Low Rate, Small Tubing (110-220V)
<b>PERIPRO-8LS</b>	Peri-Star™ Pro, 8-channel, Low Rate, Small Tubing (110-220V)
<b>801566</b>	Mini Vacuum Pump (110V)
<b>801963</b>	Mini Vacuum Pump (220V)
<b>SI-COLUB</b>	Constant Load Module
<b>97204</b>	Pinger for the Anti-Oscillation Module

# Cell Tester System

Perform muscle physiology tests on single skeletal muscle fibers

## Features

- Fits ANY inverted microscope
- Unique rotational stage—improves the experimental throughput
- Two integral piezo manipulators included

## Benefits

- Simultaneous use of the multi-purpose biofluorometer\* (**SI-BF-100**) for tissue fluorescence (calcium, ATPase via NADH, FAD)
- For use with small multi-cellular preparations or skinned muscle fibers\*
- Stimulate and perform perturbations of muscle fibers
- Options available for customization of the system

## Applications

- Stretching and relaxation of cells with nanometer resolution
- Classic cross bridge cycling studies in cardiac tissue
- Skeletal muscle fiber axial stretch and isometric force
- Intracellular calcium concentration/distribution in muscle fibers as muscle force is measured
- ATPase activity as muscle force is measured

The **SI-CTS200** system, which is the result of blending the latest technologies in electronics, mechanics and optics, permits researchers to investigate living systems at a new level of observation. The Cell Tester provides researchers with the comprehensive ability to investigate and characterize the physiological, biomechanical and biophysical properties of single isolated living cells.

### Study multi-cellular or skinned fibers

The **SI-CTS200** systems can be used on small multi-cellular preparations or skinned muscle fibers or strips. It is designed to sit on the stage of any standard, research-level, inverted microscope while maintaining the optical path of the microscope for simultaneous fluorescence or confocal imaging.

### Stimulate and perform perturbations of cells

The Cell Tester is an integrated system of components needed to maintain and handle muscle fibers, stimulate and perform perturbations of the cells, and detect, amplify, and record signals, like contractile force.

### Customize your system with options

- Optical force transducer for measuring with nanonewton sensitivity and integrated cell attachment system.
- Equipped with a nanomotor for stretching and relaxing cells with nanometer resolution.
- A rotating cuvette system for easy alignment of cells increases productivity. It is designed to orient cells in the XY plane so that no physical manipulation of the position of the tissue itself is required.
- Interchangeable bath inserts provide a range of options for the handling of live fibers.
- **LabTrax-MDAC** data acquisition system records and control your experiments, using up to eight analog inputs.
- Biofluorometer (See page 20.)



SI-CTS200

## References

**Prosser BL, Ward CW, Lederer WJ.** X-ROS signaling: rapid mechano-chemo transduction in heart. *Science*. 2011 Sep 9;333(6048):1440-5. PMID: 21903813

### FORCE TRANSDUCER SPECIFICATIONS

RANGE	0–5 mN (0–0.5 g)
FORCE RESOLUTION	20 nN at 10× gain
COMPLIANCE	10 μm/mN
NOISE	0.3 μN
RESONANCE FREQUENCY	250 Hz (This is eliminated from the measurement by the AOSUB)
TIME RESOLUTION	7 ms

*Resolutions were determined while using the SI-AOSUB anti-oscillation filter.*

### NANOMOTOR SPECIFICATIONS

TOTAL TRAVEL	±90 μm
RESOLUTION	20 nm
SMALLEST STEP	60 nm
INPUT	±10V (calibrated at 10 μm/V)

### CELL TESTER PLATFORM SPECIFICATIONS

DISTANCE (bottom of the base plate to the bottom of glass) 0.508mm

### ORDERING INFORMATION

**SI-CTS200** Complete Cell Tester System

*System Includes: Rotating Cuvette System; Micromanipulator System; Signal Conditioning Amplifier with four modules: Optical Transducer Amplifier; Temperature Controller; Anti-Oscillation Unit; Position Controller; data acquisition system with MDAC Software for recording, controlling stimulation and nanomotor position; Force Transducer of choice; Nanomotor; Glass Fiber Cell Mounts (1 set)*

**SI-CTS200B** Cell Tester, Non-Rotating, no Micromanipulator System

**SI-CTS200A** Cell Tester, Manual Platform, no Micromanipulator System

For additional information, including application notes, go to [wpiinc.com/celltester](http://wpiinc.com/celltester)

\*Contact WPI for help evaluating your specific experimental requirements.

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • [wpiuk@wpi-europe.com](mailto:wpiuk@wpi-europe.com) • [www.wpi-europe.com](http://www.wpi-europe.com)  
Brazil: (013) 406-29703 • [info@brazil.wpiinc.com](mailto:info@brazil.wpiinc.com) • [www.wpiinc.com](http://www.wpiinc.com)

Germany: +49 (0)6031 1602171 • [wptide@wpi-europe.com](mailto:wptide@wpi-europe.com) • [www.wpi-europe.com](http://www.wpi-europe.com)  
China: +86 21 6888 5517 • [chinasales@china.wpiinc.com](mailto:chinasales@china.wpiinc.com) • [www.wpiinc.net](http://www.wpiinc.net)

# Signal Conditioning Amplifier System

Choose the amplifier modules you need to measure nearly anything!



SI-BMFA power frame enclosure plus optional modules (pages 14-17).

## Features

- Ergonomic design
- 8-Channel
- Small footprint
- Backplane design includes provision for configurable communication between modules

## Benefits

- Up to eight modules connected through the backplane of the chassis allows researchers to assemble the set of electronics they need for their own custom applications
- Standard configuration options are available for Muscle Tester platforms

## Applications

- Process transduction of physical signals, displacement and optical force transducer outputs

The transduction of physical signals in the last decade has increasingly moved in the direction of the computer with an electronic amplifier. Further signal conditioning and analysis of sampled raw data is then handled efficiently in software, like WPI's MDAC package. The software preserves the raw data and is highly reliable. Operations such as integration, differentiation, filtration and even waveform generation are now efficiently handled in software. On the other hand, however, the transduction of physical signals such as bio-potentials, force, temperature, pressure or ionic concentrations must be measured with an electronic amplifier.

SIH/WPI's physiology amplifier system focuses on this idea and provides a flexible electronic platform intended to process the transduction of physical signals, displacement transducer outputs from force transducer signals. This platform simply focuses on the reliable transduction of the electronic signal and provides a convenient passage for the translation of real world signals to a computer for analysis.

### Eight modules connected through backplane

The system consists of a chassis with eight slots on the backplane and includes an ultra quiet, shielded power supply. All of the module outputs

are routed to rear panel connectors. If you prefer, outputs may be routed internally to the inputs of other modules. The system has a small footprint and may be stacked to provide as many channels as you need.

The **SI-BMFA** Power Frame is the foundation of the SI modular physiology suite. It incorporates a robust power supply that can accommodate up to eight physiology modules, which can be mixed or matched in any combination. Modules are quick and easy to install, thanks to an innovative and mechanically solid card rail system.

The system is flexible and configurable. A variety of modules are available for the Signal Conditioning Amplifier System, and you can mix and match the modules to suit your requirements.

**Optical Transducer Amplifier—SI-BAM21-LCB** (page 15)

**Programmable Filter Module—SI-PF100B** (page 17)

**Linear Motor Control Module—SI-MOTDB** (page 14)

**Temperature Control Module—SI-TCM2B** (page 16)

**Anti-Oscillation Module—SI-AOSUB** (page 17)

**Constant Load Module—SI-COLUB** (page 16)

### Standard options available for muscle testers

When the system is ordered with **SI-MT** (Muscle Tester) system, the Signal Conditioning Amplifier System (chassis) is configured with an **SI-BAM21-LCB**. Optional modules include an **SI-TCM2B** Temperature Control Module, an **SI-MOTDB** Linear Motor Controller, an **SI-PF100B** Programmable Filter Module, the **SI-AOSUB** Anti-Oscillation Module and the **SI-COLUB** Constant Load Unit. The Temperature Control Module and Linear Motor Controller require two slots each on the chassis backplane.

## ORDERING INFORMATION

<b>SI-BMFA</b>	Power Frame Enclosure
<b>SI-MOT-MT</b>	Linear Motor for SI-HTBM system
<b>SI-MOTDB</b>	Linear Motor Controller
<b>SI-BAM21-LCB</b>	Optical Transducer Amplifier
<b>SI-TCM2B</b>	2-Channel Temperature Control Module
<b>SI-COLUB</b>	Constant Load Module
<b>SI-AOSUB</b>	Anti-Oscillation Module
<b>SI-PF100B</b>	Programmable Filter Module

# Signal Conditioning Amplifier System

Choose the amplifier modules you need to measure nearly anything!

## Linear Motor Control Module

### Features

- Powers the motor and provides an output indicating the actual motor position
- Connects to Analog to Digital Converter output of the computer or data acquisition system (like LabTrax-MDAC) to allow the control of the programmed waveform and timing of the motor control
- Input range of  $\pm 10$  VDC
- Linear motor position is determined by a DC value applied from the Position In port
- 2-Slot control module for a linear motor

### Benefits

- Over current protection that automatically shuts down when the supply voltage dips below the reference value
- External  $\pm 10$  VDC position command input for control by a constant load module or data acquisition system

### Applications

- Slack test
- Isotonic release
- Constant velocity release
- Stretch release
- Eccentric/concentric contractions (intact muscle)



SI-MOTDB

The SI-H Linear Motor Controller is designed for use with the SI-H line of muscle physiology research platforms. For systems that require a linear motor, this unit provides the precision control of the motor. A linear motor is required to perform mechanical muscle testing such as slack-tests, isotonic releases, constant velocity releases, stretch releases, after-loaded contractions and eccentric contractions (intact muscle). The position of the linear motor is determined by a combination of the data from the controller indicating the current position and the DC value applied to the front panel at the Position In port. The applied Position In signal can be provided by a data acquisition system (LabTrax-MDAC to use standard or customized protocols). The data acquisition analog output signal is set to define the waveform and timing pattern of force to be applied to the sample.



The SI-MOT motor can be used with any SI-MTM-L Muscle Tester.

### SI-MOTDB SPECIFICATIONS

POWER REQUIREMENTS	12 V DC provided by the chassis
INPUT	$\pm 10$ V DC
TRAVEL	1 mm/2 VDC current
MAXIMUM TRAVEL	7 mm ( $\pm 3.5$ mm from center of travel)

### ORDERING INFORMATION

**SI-MOT-MT** Linear Motor for SI-MTM-L system  
**SI-MOTDB** Linear Motor Controller

## Biofluorometer

Reliable, simplified and affordable LED based fluorometer

Perfect for:  
Ratiometric Calcium & ATPase

Recent advancements in optics and LED technology simplify ratiometric calcium imaging, making this equipment more affordable. A breakthrough in WPI patented technology allows the SI-BF-100 to use wavelengths below 380 nm and produce more light in those spectra. This technology significantly cuts the cost of photometric calcium imaging without sacrificing resolution or quality.

### Streamlined system to reduce errors

Up till now, calcium imaging systems have been required to compensate for errors and noise introduced by the complexity of their design. The systems require mechanical filters and use expensive xenon or mercury light sources. The beauty of the SI-BF-100 is its simplicity. The elegance of its design reduces the noise introduced into the system and the errors inherent in traditional systems.

This single wavelength spectrophotometer can be customized for your specific application. See page 20 for details.



SI-BF-100



### WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
Brazil: (013) 406-29703 • info@brazil.wpiinc.com • www.wpiinc.com  
Germany: +49 (0)6031 1602171 • wptide@wpi-europe.com • www.wpi-europe.com  
China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net



## Optical Force Transducer Amplifier Module

### Features

- Designed for use with SI-H Muscle Tester Platforms
- Rapid auto zeroing function with fine offset adjustment
- Offset indicator LED's
- Multiple gain ranges with adjustable fine tuning for precise calibration
- 1X, 2X, 5X, 10X gains. Optional factory setting allows for 10X, 20X, 50X and 100X gains

### Benefits

- Single amplifier that spans the entire SI-KG optical force transducer line
- Manual calibration to the  $\pm 10$  V measurement range, for covering maximum forces ranging from 50 mN to 2N. When using the **SI-KG7A**, can calibrate as small as 5 mN.

### Applications

- Muscle physiology studies, including isotonic studies

The **SI-BAM21-LCB** amplifier for SI-KG Optical Force Transducers is used in conjunction with the SI-H muscle physiology systems. The **SI-BAM21-LCB** powers the force transducer and converts the output of the transducer to an amplified analog voltage that is proportional to the force applied to the force transducer. The output signal can be multiplied by a factor of 1, 2, 5 or 10 to provide better resolution for a minimal change in applied force.

NOTE: An optional factory setting increases the multiplier by a factor of 10, allowing the signal to be multiplied by 10, 20, 50 and 100.

The **SI-BAM21-LCB** amplifier works with SI-KG optical force transducers to:

- Generate an analog output (-10 VDC to +10 VDC) that is proportional to the force applied to the tissue sample.
- Supply a DC voltage that powers the SI-KG force transducer to which it is connected.

Also available in a single stand-alone enclosure, either version provides an incredibly quiet, linear and stable transducer signal to your data recording system.



**SI-BAM21-LCB**

## How the SI-BAM21-LCB amplifier works

In a typical setup, a muscle is held by a force transducer. The force transducer is connected to the **SI-BAM21-LCB**. As the muscle contracts or releases, the transducer converts the force into an electrical current signal which is proportional to the force applied to the transducer. The **SI-BAM21-LCB** converts the current signal into a voltage signal that can be displayed on the screen of the recording device.

Before initiating an experiment, the **SI-BAM21-LCB** must first be zeroed. This sets the baseline for measurements to follow.

The output signal is buffered and multiplied by 1, 2, 5 or 10, depending on the Gain switch setting on the front panel of the amplifier module. The  $\times 10$  setting is useful when output signals are extremely small. Finally, the force proportional signal is sent through the output amplifier circuit.

The analog output has a range of  $-10$  V to  $+10$  V that drives the **LABTRAX-MDAC** data acquisition system, multimeter or oscilloscope.

NOTE: When the Signal Conditioning Amplifier System is configured at the factory for an **SI-HTB** or **SI-MT** Muscle Tester system, the signal is routed internally from the **SI-BAM21-LCB** module to the **SI-AOSUB** module.

### SI-BAM21-LCB SPECIFICATIONS

INPUT CONFIGURATION	Current to voltage converter
GAIN	1 $\times$ , 2 $\times$ , 5 $\times$ , 10 $\times$ Optional factory setting: 10 $\times$ , 20 $\times$ , 50 $\times$ , 100 $\times$
INPUT OFFSET ADJUSTMENT	$\pm 2.0$ VDC
OUTPUT IMPEDANCE	470 $\Omega$
POWER	12 VDC provided by chassis
OUTPUT RANGE	$\pm 10$ VDC

### ORDERING INFORMATION

<b>SI-BAM21-LCB</b>	Optical Transducer Amplifier
<b>SI-BAM21-LC</b>	Stand alone Optical Force Transducer Amplifier

### OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>LABTRAX-MDAC</b>	LabTrax 8/16 with MDAC software
<b>2851</b>	BNC Cable
<b>SI-KG2</b>	0-2 N Force Transducer
<b>SI-KG2B</b>	0-200 mN Force Transducer
<b>SI-KG4</b>	0-50 mN Force Transducer
<b>SI-KG7a</b>	0-5 mN Force Transducer
<b>SI-KG7b</b>	0-10 mN Force Transducer
<b>SI-KGxx</b>	contact WPI for specialty transducers with different range

See *Optical Force Transducers* on page 22.



**SI-BAM21-LC**

The stand-alone **SI-BAM21-LC** Force Transducer Amplifier has all the capabilities of the Signal Conditioning Amplifier module (**SI-BAM21-LCB**).

# Signal Conditioning Amplifier System

Choose the amplifier modules you need to measure nearly anything!

## Temperature Control Module

### Features

- Uses PID control to maintain a constant temperature with  $\pm 0.1^\circ\text{C}$  tolerance
- Easy to control with simple interface
- Also available as a stand-alone device
- 2-Slot control module for maintaining temperature control

### Benefits

- Controls two cuvettes simultaneously
- User defined high and low alarm warnings

### Applications

- Muscle physiology applications using SI-H muscle tester platforms



SI-TCM2B

The SI-H Temperature Control Unit is designed for use with the SI-H line of muscle physiology research platforms. It maintains the temperature of an SI-H cuvette up to  $45^\circ\text{C}$ . It is accurate to  $0.1^\circ\text{C}$ . The circuit is appropriate to RTD (resistive temperature device) applications. It linearly converts a temperature reading to a voltage that is displayed as a temperature on the SI-TCM2 and can be recorded. This unit is available in a stand-alone model and as a module for the Signal Conditioning Amplifier System backplane.

Call for details and pricing information.

### SI-TCM2 SPECIFICATIONS

INPUT CONFIGURATION	Current to voltage converter
POWER REQUIREMENTS	12 V DC at 2.5A 50/60 Hz wall adaptor, 2.5 mm ID/5.5 mm OD with positive center DC barrel (included-WPI #801513)
OPERATING TEMPERATURE RANGE	Room temperature
MAXIMUM TEMPERATURE	$45^\circ\text{C}$
DISPLAY PRECISION	$0.1^\circ\text{C}$
CONTROLLER RESOLUTION	$0.1^\circ\text{C}$
CUVETTE TEMPERATURE SENSOR	1000 $\Omega$ RTD (1000 $\Omega$ at $0^\circ\text{C}$ )

### ORDERING INFORMATION

SI-TCM2B	2-Channel Temperature Control Module
SI-TCM2	2-Channel Temperature Control, Stand-Alone

## Constant Load Module

### Features

- Offers three modes including Constant Load, External Loop and Bypass
- Can be configured using LabTrax-MDAC data acquisition system

### Benefits

- Precise time resolution for feedback control in isotonic muscle testing

### Applications

- Muscle physiology applications using SI-H muscle tester platforms


The SI-COLUB Constant Load Module for performing constant load experiments, has augmented flexibility. In its primary mode (Constant Load) the unit takes an external trigger command from the force transducer to perform a constant load cycle. In addition, the module allows for a different external trigger or you can completely bypass the module without having to switch cabling.

The Constant Load Module lets you maintain a constant force, muscle length or sarcomere length rather than keep the total length of the preparation constant during an isotonic test. This is accomplished using a feedback loop.

The SI-COLUB monitors a designated parameter to determine how much force is necessary. It also monitors a feedback signal. The motor position command signal driving the motor is constantly adjusted to drive the feedback signal to the commanded setpoint.



SI-COLUB



### SI-COLUB SPECIFICATIONS

COMMAND REQUEST	$\pm 10$ V
FEEDBACK	$\pm 10$ V
MOTOR OUTPUT	$\pm 10$ V
POWER REQUIREMENTS	12 V DC provided by the chassis

### ORDERING INFORMATION

SI-COLUB	Constant Load Module
----------	----------------------

## Anti-Oscillation Module

### Features

- Neutralizes nearly 100% of the unavoidable transducer oscillation in optical force transducer
- High time-resolution of the corrected force signal in ms range
- Additional powerful signal smoothing with high time fidelity for all higher harmonics

### Benefits


- Goes beyond traditional low-pass filtering of an oscillated signal
- No phase-shift in the force signal
- High time fidelity in following the applied step signal

### Applications

- All types of optical-based force measurements

Within a force transducer, the applied force displaces an elastic component of the transducer. This displacement is transformed into an electrical signal. If the displacement of the elastic component is small and if the displacement is linearly transformed into an electrical signal, the external force and electrical signal are linearly related to each other.

However, shortly after the force changes rapidly the electrical signal does not follow the applied force linearly, as the elastic element and the force transducer pin have a mass. So, the system acts not only as an elastic device but also as a harmonic oscillator, notably when the force changes rapidly, resulting in oscillation of the system around the new displacement level (ringing phenomenon). To address this problem, WPI's Anti-Oscillation Unit (**SI-AOSUB**) uses a genuine electronic approach to remove the unavoidable transducer oscillation, based on electronic differentiation of the force signal and subsequent low-pass filtering of the resulting high-frequency signal.



The upper trace is a force transient obtained directly from the bridge amplifier output, and the lower trace shows the signal after it passes through the anti-oscillation module.



To adjust the anti-oscillation filter properly, the transducer is excited at its resonance frequency using a magnetic driver or pulser (WPI #97204). The Pinger (WPI #97204) is included with the anti-oscillation module.



**SI-AOSUB**

## Programmable Filter Module

### Features

- Low Pass Filter
- Bessel and Butterworth filters
- Programmable cutoff frequency
- Signal may be routed through the backplane or through the front panel BNCs

### Benefits

- General purpose filter module for nearly all type of noisy signal

### Applications

- Muscle physiology applications using SI-H muscle tester platforms

When you use a motor, an **SI-PF100B** Programmable Filter is necessary to minimize the natural vibration. It is designed so you can eliminate the resonance frequency without affecting the signal of interest. It is a low pass filter set to pass signals of interest below the specified frequency. It can be calibrated from 5 to 1,000 Hz.

You may select either a Bessel or a Butterworth filter. Then, you must carefully select the cutoff frequency based on the typical resonance frequency of your force transducer and your own experimental setup.

When the Signal Conditioning Amplifier System electronics are configured at the factory with an **SI-PF100B** Programmable Filter, the signal is routed internally from the amplifier module (**SI-BAM21-LCB**) to the **SI-PF100B**. If you prefer, the signal may be routed from the amplifier through the ports on the front panel of the Programmable Filter using a standard BNC cable.



### SI-PF100 SPECIFICATIONS

POWER	12 VDC provided by chassis
INPUT	± 10 VDC
CUTOFF FREQUENCY RANGE	5-1,000 Hz
FILTER TYPES	BESSEL, BUTTERWORTH

### ORDERING INFORMATION

<b>SI-PF100B</b>	Programmable Filter Module
<b>SI-PF100</b>	Programmable Filter in Stand-Alone Enclosure



The Programmable Filter Module as stand-alone unit (**SI-PF100**) can be used for filtering two independent input signals, with the installation of a second **SI-PF100B** module (not included).

### SI-AOSUB SPECIFICATIONS

POWER	12 VDC provided by chassis
INPUT	± 10 VDC
PULSER OUTPUT	0 – 10 VDC adjustable 0.1 Hz – 4.0 KHz
DAMPING FREQUENCY RANGE	0.1 Hz – 2.0 KHz
OUTPUT RANGE	± 10 VDC

### ORDERING INFORMATION

<b>SI-AOSUB</b>	Anti-Oscillation Module
<b>97204</b>	Pinger for Anti-Oscillation Module
<b>2851</b>	BNC Cable

# Optical Sarcomere Spacing System

Measure sarcomere spacing with nanometer accuracy

## Features

- Monitors high-speed dynamic sarcomere length changes
- Sarcomere spacing measurement with nanometer accuracy
- NIH open-source  $\mu$ Manager\* software interface between USB 3.0 camera and microscope
- Synchronous use of Sarcomere Spacing plug-in and  $\mu$ Manager software
- User defined Region of Interest (ROI) for image capturing from whole cell down to a minimum of 6 consecutive sarcomeres
- Compatible with Windows 10

## Benefits



- Camera image capturing using the standard  $\mu$ Manager software familiar to researchers
- Use Image toolboxes for complementary image processing
- Live image processing and display of length changes in sarcomere spacing of contracting muscle cells and fibers
- Synchronous storage of live images and processed data in user-friendly MS Excel format, for further offline analysis
- User defined ROI setting allows focusing on specific cell regions and changing image capturing speed
- High-speed USB 3.0 camera fits to all standard inverted microscopes with C-Mount adapter

## Applications

- Cardiac muscle physiology
- Skeletal muscle physiology
- Quality control of muscle slices in meat industry

WPI's Optical Sarcomere Spacing System quantifies live images to detect length changes in sarcomere spacing in nanometer accuracy during twitch or tetanus contractions of muscle cells/fibers. The NIH  $\mu$ Manager software interface and WPI's image processing plug-in for sarcomere spacing detection monitor and measure sarcomere length changes of the contracting muscle cells in real time with an optimal frame rate of 500 FPS.

The Optical Sarcomere Spacing System fits perfectly with WPI's Cell Tester system (SI-CTS200) to capture twitch contraction of cardio myocytes, synchronous with the force signal. Extend the quantification by adding measurements of  $Ca^{2+}$  or ATPase concentration of contracting muscle cells/fibers using the WPI's Biofluorometer (SI-BF-100) simultaneously.



*$\mu$ Manager software interface and Optisarc image processing plug-in for sarcomere spacing detection run simultaneously during live image capturing of contracting cardio myocyte on selected ROI (Image with courtesy of Dr. Michael Kohlhaas, Universitätsklinikum des Saarlandes, Homburg).*

The Optical Sarcomere Spacing System can also be extended for capturing sarcomere length changes in slow-twitch and fast-twitch skeletal muscle cells/fibers.

Sarcomere spacing detection of any muscle slice is also possible using genuine image processing methods for high reproducible measurements of selected sarcomere regions in static images.

## OPTISARC CAMERA SPECIFICATIONS

IMAGE SENSOR	1/2.5" 5 megapixel monochrome CMOS
ACTIVE PICTURE ELEMENTS	2,592 (H) x 1,944 (V).
MAXIMUM FRAME RATE	2,106 FPS @ 32 x 32.
OPTIMAL FRAME RATE	500 FPS @ 296 X 148
CONNECTION	USB 3.0 Micro-B

## ORDERING INFORMATION

SI-OSARC	USB 3.0 CMOS Camera with OptiSarc software plug-in for live image capture to detect sarcomere spacing
INV-101	Inverted Microscope with C-mount Adapter
SI-BF-100**	Biofluorometer
SI-CTS200**	Cell Tester

$\mu$ Manager\* was developed at Vale laboratory at UCSF, funded by an NIH grant R01-EB007187 from the National Institute of Biomedical Imaging and Bioengineering (NIBB). <https://micro-manager.org/> and <https://imagej.nih.gov/ij/>

\*Arthur Edelstein, Nenad Amodaj, Karl Hoover, Ron Vale and Nico Stuurman (2010): *Computer Control of Microscopes using  $\mu$ Manager*. Current Protocols in Molecular Biology, Chapter 14, Unit 14.20, 22 pages.

\*\* For simultaneously use with SI-BF-100 and SI-CTS200, ask WPI about selecting an appropriate microscope.

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • [wpiuk@wpi-europe.com](mailto:wpiuk@wpi-europe.com) • [www.wpi-europe.com](http://www.wpi-europe.com)  
Brazil: (013) 406-29703 • [info@brazil.wpiinc.com](mailto:info@brazil.wpiinc.com) • [www.wpiinc.com](http://www.wpiinc.com)

Germany: +49 (0)6031 1602171 • [wpi@wpi-europe.com](mailto:wpi@wpi-europe.com) • [www.wpi-europe.com](http://www.wpi-europe.com)  
China: +86 21 6888 5517 • [chinasales@china.wpiinc.com](mailto:chinasales@china.wpiinc.com) • [www.wpiinc.net](http://www.wpiinc.net)

# Ca<sup>2+</sup> Detection in Muscle Tissue using Fluorescence Spectroscopy

The use of fluorescent probes in cell physiology has emerged as indispensable tool in the analysis of cell functioning over recent years. The physics underlying fluorescence is illustrated by the electronic-state diagram (so-called Jablonski diagram, see Fig. 1), showing the three-stage process to create the fluorescent signal (Excitation - Excited/State Lifetime - Fluorescence Emission) in a fluorophore/indicator and simplified described below.





Fig. 1- Jablonski diagram illustrating the processes of fluorescence by absorption of higher photon energy by a fluorophore and subsequent emission of lower photon energy, resulting in fluorescence during the fluorescence-lifetime.

Fluorescence is obtained when an excitation photon ( $h\nu_{EX}$ ) from an external source, such as a high-power LED, is absorbed by a fluorophore that elevates its energy ( $S1'$ ). During the fluorescence-lifetime, the elevated energy ( $S1'$ ) decays to a lower energy state  $S1$ . Then, fluorescence results in the emission of a photon with lower energy ( $h\nu_{EM}$ ) and at a lower wavelength. Fundamental in spectroscopy is the difference in energy or wavelength represented by ( $h\nu_{EX}-h\nu_{EM}$ ), which is called the Stokes shift. The Stokes shift allows efficient discrimination of the excitation, making fluorescence a very sensitive technique and able to be detected against a low background, isolated from excitation photons.

Fig. 2- Typical excitation and emission spectra, showing the lower wavelength of the excitation source and higher wavelength of the fluorescence.



Four essential elements of fluorescence signaling can be then identified to build up a detection system:

- Excitation light source adapted to the absorption bandwidth of the fluorophore (e.g. high-power LED of specific wavelength)
- A fluorophore/indicator (e.g. Fura-8 for free Ca<sup>2+</sup> detection in muscle tissue)
- Emission wavelength filters to limit the bandwidth of the emission photons or overlapping bands
- A detector system that registers the fluorescence light and produces a recordable output as an electrical signal (e.g. Photomultiplier tubes).

Regardless of the application, compatibility of these four elements is essential for optimizing fluorescence detection.

## Example of free Ca<sup>2+</sup> detection in muscle tissue

Typically, a fluorescent dye is introduced into tissue or single cells to obtain a fluorescent response of the labeled molecule. A typical example is the detection of the transient increase in the cytoplasmic/myoplasmic



free calcium concentration ( $\Delta[Ca^{2+}]$ ) as the intermediate signaling event of the excitation-contraction coupling. The quantification of  $\Delta[Ca^{2+}]$  is done using a monochromatic light to excite the dye labeled Ca<sup>2+</sup> molecule in a tissue/cell sample either in a tissue bath or microscopic experimental set-up. The emitted fluorescence signal from the indicator dye can be then used to monitor the amplitude and time-course of the  $\Delta[Ca^{2+}]$  detected by sensitive detectors, such as highly sensitive photomultiplier tubes (PMT module) or cameras.

The ratiometric indicator dye Fura-8 was selected for the detection of free calcium concentration ( $\Delta[Ca^{2+}]$ ) in heart muscle tissue slices. Fura-8 was excited at 365nm and 410 nm wavelengths and the emission recorded at 535 nm wavelength in dual excitation/ single emission mode. The advantages of choosing this ratiometric measurement technique with dual excitation/single emission using Fura-8 were minimization of movement artifact, cancellation of possible effects of uneven loading, inhomogeneous distribution of fluorescence indicator in the cells or indicator bleaching in the detection of free calcium concentration ( $\Delta[Ca^{2+}]$ ) in the muscle tissue.

This allowed quantification and comparison between:

- High spatial versus high time resolution techniques on the human left ventricular slices
- The possibility to measure free calcium concentration ( $\Delta[Ca^{2+}]$ ) transients in a horizontal tissue bath on human left ventricular slices or murine slices.




Fig. 3- Qualitative representation some results of free Ca<sup>2+</sup> detection in human heart slices using the SI-BF-100 system.

Average fluorescence intensities of Fura-8 loaded human left ventricular slices detected at 525 nm, when excited at 340 nm and 410 nm, respectively, and ratios calculated (lower trace) from the imaging data of a Rolera EM-C2 camera (left). Right, the response of the SI-BF-100 detected at 525 nm, when excited at 365 nm and 410 nm wavelength using two aperture settings and calculated ratios (lower traces). Furthermore, fluorescent data collected with the small aperture setting and low-pass filtered at 50 Hz is shown (from: Belz et al., Proc. SPIE 9702, 2016).

## References

- Belz M., et al. Fiber optic Biofluorometer for physiological research on muscle slices. Proc. SPIE 9702, Optical Fibers and Sensors for Medical Diagnostics and Treatment Applications XVI, 2016.
- Spectrophotometry. Wikipedia, the free encyclopedia, 2017 (Cross-references).

# Biofluorometer

## Reliable, 2-channel, LED-based fluorometer

### Features

- Two photomultiplier inputs
- Light excitation with high power LEDs
- Modes: single excitation & single emission, single excitation & dual emission, dual excitation & single emission
- 2 channel mode (two single excitation & single emission)
- Optical connections: Liquid Light Guides (LLGs) and SMA terminated fibers
- Sampling rates up to 1 kHz (1000 ratios/second)
- Automatic LED light drift correction for long term measurements
- Automatic room light correction
- Optional fiber optic probes for horizontal tissue bath applications
- Optional imaging probes for Langendorff systems
- Optional attachments for direct connection to fluorescence microscopes via epifluorescence port (excitation) and C-Mount (Emission) via liquid light guides

### Benefits

- Versatile instrument for horizontal tissue bath, Langendorff and microscope applications
- Warmup time < 1 minute
- Low bleaching mode (5% LED On time)
- Can be combined with imaging based Sarcomere Detection System (OptiSarc)
- SMA and Liquid Light Guide connections
- Single and dual emission or excitation detection methods
- 2-channel instrument for single excitation/single emission dyes
- Customized analysis techniques in WPI's MDAC data acquisition software

### Applications

- Ratiometric calcium measurement and ATPase
- Fluorometric applications in neuroscience and cell biology
- The **SI-BF100** is an LED-based fluorometer for life science applications. It is ideally suited for ratiometric calcium detection (FURA-8) and ATPase detection (via NADH fluorescence). With up to three LED modules (wavelengths), the **SI-BF100** covers many fluorometric applications in neuroscience and cell biology.

The **SI-BF-100** enables the detection and analysis of fluorescence signals in four different modes:

- **Single excitation/single emission**—In this classical mode, a fluorophore is excited at one wavelength and the fluorescence signal is detected at a single higher wavelength using one photomultiplier. The concentration of the analyte is directly proportional to the intensity of the detected signal.
- **Dual excitation/single emission**—A fluorophore is excited at two wavelengths and the fluorescence signal is detected at one wavelength using one photomultiplier. The concentration of the analyte is proportional to the ratio of the two detected fluorescence signals. This ratiometric concept minimizes the effect of indicator dye bleaching and motion artifact in experiments. A typical example is the detection of free calcium in muscle tissue using the indicator dye Fura-8™.
- **Single excitation/dual emission**—A fluorophore is excited at one wavelength and the fluorescence signal is detected at two wavelengths using two photomultipliers.
- **Dual excitation/dual emission**—Two separate fluorophores are excited at different wavelengths and the fluorescence signal of each fluorophore is detected at two separate wavelengths using two photomultipliers.




SI-BF-100


### Experimental Setups

Possible experimental setups for the **SI-BF-100** include:


**A:** This **SI-BF-100LLG** setup depicts the attachment to a fluorescent microscope.



**B:** This **SI-BF-100LLG** setup shows a single fiber based detection system with a dichroic mirror and an imaging system at the distal end of a liquid light guide. Dichroic mirrors are used in A and B to separate excitation and emission light.



**C:** The **SI-BF-100SMA** setup has a round shaped fiber bundle that is used to deliver excitation light to and pick up fluorescent light from the sample (adapted from Belz et al., 2016).




### Single and dual emission probes

All probes use fibers with a 300 µm core diameter. Excitation fibers have 1000 µm SMA connectors for Excitation and Double Emission Probes and 1500 µm SMA connectors for Single Emission Probes.



This SMA Single Emission Probe has a for muscle strips, where the middle fibers are used for excitation and the rectangular fiber pattern, optimized surrounding fibers are used for detection.

The Biofluorometer can be connected to the epifluorescence port of an inverted microscope and its high intensity LED light source is used for excitation of the corresponding dye.



Epifluorescence Port

Emission Fiber

Excitation Fiber

### WORLD PRECISION INSTRUMENTS


UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
Brazil: (013) 406-29703 • info@brazil.wpiinc.com • www.wpiinc.com

Germany: +49 (0)6031 1602171 • wptide@wpi-europe.com • www.wpi-europe.com  
China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net

# Biofluorescence of calcium use in heart muscle slices with Fura-8

## Calcium sensing

A specific target molecule in muscle contraction is  $Ca^{2+}$  as key intermediate signaling event between excitation and contraction of muscle fibers, and thus essential for the analysis of the force development in muscles. On the cellular level, force production is therefore directly related to the transient increase in the myoplasmic free calcium concentration. More specifically, the assessment of both parameters simultaneously is therefore critical in the evaluation and interpretation of the force development characteristics. Fluorescence techniques used in conjunction with muscle research systems (like WPI's **SI-MT-L** or **SI-HTB2**) to record muscle force is an **innovative technique** in cardiac muscle and skeletal muscle physiology. WPI's Biofluorometer (**SI-BF-100**) was specifically developed to monitor rapid changes in  $Ca^{2+}$  transients, i.e.  $\Delta[Ca^{2+}]_i$ .



Organ bath setup: Average fluorescence intensities of Fura-8 loaded murine myocardium slices excited at 365/410 nm wavelength and detected at 525 nm. Left, using two apertures for fluorescence detection via the microscope. Right, same experimental condition using optical fiber fluorescence detection with the SI-BF-100SMA in the organ bath, sampled at 1 kHz. The organ bath data are low-pass filtered at 200 Hz or 50 Hz cut-off frequency. Adapted and reprinted with permission from Belz et al., Proc of SPIE Vol 9702, 97020Q-1 - 97020Q-11, 2016.

## Advantages using Fura-8 dye:

- More sensitive to calcium than Fura-2
- Higher signal-to-noise ratio than Fura-2
- Simplified dye loading by incubation for 1 hour at room temperature
- Emission peak response is shifted to longer wavelength (peak at 525 nm)
- Red-shift dual excitation wavelength (354 nm and 415 nm)

## Application in heart muscle slices

The Biofluorometer opens a wide field in functional fluorescence research, by studying the fundamental and/or applied aspects of the underlying energetics and signaling aspects of muscle contraction. This is notably useful in:

- Pre-clinical & Toxicological Studies:
  - screening of potential drugs
  - evaluating the side effects of drugs
  - evaluating models of cardiac disease
- Sports & Rehabilitation:
  - disuse vs. overuse
  - muscle damage
  - function for heart transplantation



The SI-BF-100 Organ/tissue bath setup with an SMA single emission probe for direct sensing of calcium fluorescence signal. Photo courtesy of Professor Andreas Dendorfer, Walter Brendel Zentrum, München, Germany.

# System configuration

The **SI-BF-100** setup for detecting  $Ca^{2+}$  transients via Fura-8 depends on the excitation and emission wavelengths of the indicator dye Fura-8.

The **SI-BF-100** setup depends on your experimental paradigm and the sample size of interest:

- **BF-100SMA-C-Fura 8**: direct sensing of calcium via direct measurement in tissue sample
- **BF-1002XLG-C-Fura 8**: direct sensing of calcium via microscope setup in muscle cell or in small tissue sample
- **BF-100CAM-C-Fura 8**: calcium imaging via microscope setup with camera in muscle cell or in small tissue sample

BF-100SMA-C-Fura 8	Includes
SI-BF-100 Main Unit with SMA connection	(1) SI-BF-100SMA
HP LED 365 nm	(1) 99209-1
HP LED 420 nm	(1) 99209-4
Optical Filter 535 nm, 43 nm BP	(1) 802238
Single Emission Probe	(1) 94650

BF-100CAM-C-Fura 8	Includes
SI-BF-100 Main Unit with LLG connection	(1) SI-BF-100LLG
HP LED 365 nm	(1) 99209-1
HP LED 420 nm	(1) 99209-4
Optical Filter 535 nm, 43 nm BP	(1) 802238
Liquid Light Guide Ø 3mm	(3) 802407
CM-CAM-2XLG Packaging	(1) 99259

BF-1002XLG-C-Fura 8	Includes
SI-BF-100 Main Unit with LLG connection	(1) SI-BF-100LLG
HP LED 365 nm	(1) 99209-1
HP LED 420 nm	(1) 99209-4
Optical Filter 535 nm, 43 nm BP	(1) 802238
Liquide Light Guide Ø 3mm	(2) 802407
CM-2XLG Packaging	(1) 99261

## SI-BF-100 SPECIFICATIONS

MEASUREMENT PRINCIPLE	Fiber optic fluorometer with 2 inputs and 1 output
DETECTOR INPUTS	2 PMTs
EXCITATION	High Power LED Modules: 365 nm, 420 nm, 470 nm, 530 nm (select any 3 modules, when ordering)
ANALOG OUTPUT RANGE	0-10V
OPTICAL CONNECTIONS	Choice of Liquid Light Guide (LLG) or SMA connections
POWER	12 V/2 A (includes external 100 – 240 V / 50 – 60 Hz power supply)
DIMENSIONS (h x w x d)	3.5 x 17 x 13 in. (88 x 431 x 330 mm)

## ORDERING INFORMATION

<b>SI-BF-100LLG</b>	Biofluorometer with LLG Optical Connections
<b>SI-BF-100SMA</b>	Biofluorometer with SMA Fiber Optic Connections
<b>99261</b>	C-Mount Microscope Attachment for 2x PMTs
<b>99259</b>	C-Mount microscope attachment for 1x camera & 2x PMTs <i>Includes 1x camera C-Mount adapter with adjustable aperture</i>

## OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>SI-BF-SMA-UPGRADE</b>	Biofluorometer Upgrade Kit for SMA optical probes
<b>SI-BF-LLG-UPGRADE</b>	Biofluorometer Upgrade Kit for LLG connections
<b>802407</b>	Liquid Light Guide (LLG), 3mm diameter, 6 ft. long
<b>M3301</b>	Manual Manipulator for securing the probe
<b>M10</b>	Magnetic Base
<b>94650</b>	Single Emission, Small Tissue Probe
<b>94689</b>	Dual Emission, Small Tissue Probe

Check our website for new LED modules, emission filters and dichroic mirrors for specific applications.

# KG Optical Force Transducers

- Simple calibration
- Different models to accommodate a wide range of forces and sensitivities
- Nearly insensitive to changes in temperature/ambient light
- Extremely high level of linearity
- Virtually indestructible with normal use
- KG transducers are required for use with **SI-BAM21** amplifiers



SI-KG4

## TRANSDUCER SPECIFICATIONS UNLOADED TRANSDUCER WITHOUT TISSUE MOUNTING SUPPORT

	Force Range	Range	Noise	Compliance	Resonance Frequency
<b>SI-KG7</b>	0-5 mN	0-0.5 g	0.2 $\mu$ N	10 nm/mN	250 Hz
<b>SI-KG7a</b>	0-5 mN	0-0.5 g	0.4 $\mu$ N	5 nm/mN	500 Hz
<b>SI-KG7b</b>	0-10 mN	0-1 g	1 $\mu$ N	1.5 nm/mN	550 Hz
<b>SI-KG2</b>	0-2 N	0-200 g	250 $\mu$ N	150 nm/mN	1.3 kHz
<b>SI-KG2B</b>	0-0.2 N	0-20 g	80 $\mu$ N		590 Hz
<b>SI-KG4</b>	0-50 mN	0-5 g	15 $\mu$ N	0.5 nm/mN	1.2 kHz
<b>SI-KGxx</b>	Contact WPI for specialty transducers with different ranges.				

## Tissue Mounting Hooks

Mounting hooks can be used in a variety of combinations, depending on the type of tissue to be examined.

Mounting hooks are sold in kits. Currently, there are 11 kit configurations, each available in four different sizes. The mounting hook size that is required depends on the force transducer used.

Vascular hooks are available for mounting blood vessels (rings). They are normally used with a pair of blunt hooks (**SI-TM8**).




For larger muscles, screw clamps (**SI-TM11**) and spring clips (**SI-TM9**, **SI-TM10**) are available.

The micrometer and motor receive a large (**SI-KG4** size) tissue mount. If a smaller tissue mount is used, the **97909** adapter is required. This adapter is included with every **SI-MT** or **SI-HTB** system.


### Ordering

When ordering tissue mounts, specify the tissue mount configuration and force transducer to be used.



Force Transducer mount pictured on the left; Motor/micrometer mount pictured on the right.

**Tweezer    Basket    Pointed Hook    Blunt Hook    Spring Clip    Screw Clip    Tendon Hook**



\*The **97909** tissue mount adapter tube (OD:0.096", ID:0.035") allows you to use **SI-KG4** size mounting hooks with **SI-KG2** size force transducers.

## ORDERING INFORMATION

	Force Transducer Mount	Micrometer/Motor Mount	Force Transducers
<b>SI-TM1</b> Papillary Muscle	Basket	Pointed Hook	Available for all force transducers
<b>SI-TM2</b> General Purpose	Pointed Hook	Pointed Hook	Available for all force transducers
<b>SI-TM3</b> Small Skeletal Muscle	Tweezer	Tendon Hook	Available for all force transducers
<b>SI-TM4</b> Trabeculae	Basket	Basket	Available for all force transducers
<b>SI-TM5</b> Papillary Muscle	Basket	Tweezer	Available for all force transducers
<b>SI-TM6</b> Muscle Rings	Blunt Hook	Blunt Hook	Available for all force transducers
<b>SI-TM7</b> General Purpose	Tweezer	Tweezer	Available for all force transducers
<b>SI-TM8</b> Muscle Rings	Blunt Hook/Vascular Hook	Blunt Hook/Vascular Hook	Available for all force transducers
<b>SI-TM9</b> Strong Skeletal Muscle	Tweezer	Spring Clip	<b>SI-KG2, SI-KG2A, SI-KG2B</b> Only
<b>SI-TM10</b> Strong Skeletal Muscle	Pointed Hook	Spring Clip	<b>SI-KG2, SI-KG2A, SI-KG2B</b> Only
<b>SI-TM11</b> Very Strong Skeletal Muscle	Pointed Hook	Screw Clamp	<b>SI-KG2, SI-KG2A, SI-KG2B</b> Only

### WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
Brazil: (013) 406-29703 • info@brazil.wpiinc.com • www.wpiinc.com

Germany: +49 (0)6031 1602171 • wpide@wpi-europe.com • www.wpi-europe.com  
China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net



# 4-Channel Data Acquisition System with Software

Low noise, high resolution system with 8 analog input and 3 analog output channels

## Features

- Powerful low-noise (<1 mV RMS) and high-resolution (16 bits) data acquisition system for sampling up-to 8 analog input channels and 3 analog output channels simultaneously, using standard BNC connections
- MDAC software provides easy to use interface controlling, with extensible standard and customized Data Processing and Analysis Tools



## Benefits

- Online Channel Math operations, general purpose Fast Fourier analysis (FFT) and digital filtering of Analog In channels
- Numerous basic signal forms can be combined to design experimental protocols, for most physiological applications
  - Factory designed standard or customized protocols
  - Semi-automated data analysis toolbox
- Protocol repeat function to avoid time consuming protocol programming of extended experiments

## Applications

- Muscle physiology (Can be used with **SI-MTM** Muscle Testing Platform, **SI-CTS200** Cell Tester System, **SI-HTB2** Horizontal Tissue Bath and **SI-BF-100** Biofluorometer)
- Stand alone general data recorder for Spectroscopy, Neuroscience and Electrophysiology (Can be used with **TBR4100** Free Radical Analyzer, Extracellular Bioamplifiers like **SYS-DAM50**, **SYS-DAM80**, **SYS-900A**, **ISO-80**, **EVOM2™** Volt Ohm Meter, **ATC2000** Animal Temperature Controller, **BP-1** Blood Pressure Monitor or the **BAT-12** Microprobe Thermometer)
- Instrument control for software triggered devices like **A365/A385/A395** Constant Current Stimulators, **MPS-2** Perfusion System, **SYS-PV820/SYS-PV830** Pneumatic PicoPumps, **Duo 773** Intracellular amplifiers, and the **SYS-TBM4M** Transbridge Transducer Amplifier (e.g. for FORT force transducers)

Knowledge of the physiological characteristics of muscle tissue can be useful to quantify beneficial or adverse effects of drug supply on muscle function in pre-clinical and toxicological studies, evaluating muscle dystrophies, training effects in sports and rehabilitation (disuse vs. overuse) and advanced physiology and biomedical research.

This is usually achieved by quantifying the contractile and/or the elastic properties of muscle tissue. This needs the programming of different and specific experimental protocols (isometric, concentric and eccentric, isokinetic or isotonic), so that the physiological structure of interest can be quantified. **LabTrax-MDAC** data acquisition software was designed for use with WPI's Muscle Physiology line to test physiological characteristics of muscle tissues in various conditions, using factory designed standard or customized protocols. The semi-automated Data Analysis Toolbox of standard protocols gives quick access to user-friendly, readable and interpretable results of the experiments.

### Variety of muscle physiology applications

The physiological response of muscle tissue to training, disuse, nutrition, drug supply and others factors may be studied by adding accessories to the system, like:

- Study of the muscle's force production capacities in combination with the Ca<sup>2+</sup> release from the sarcoplasmic reticulum (SR) and ATPase consumption. The perfect instrument for this is WPI's Biofluorometer (**SI-BF-100**) in combination with any system of WPI's Muscle Physiology line, controlled via **LABTRAX-MDAC**.

- Study of the muscle's force production capacities from direct muscle or peripheral nerve stimulation. For this experiment, use WPI's programmable isolated current stimulators (**A365**, **A380** or **A395**), controlled via **LABTRAX-MDAC**.

**LABTRAX-MDAC** provides easily used continuous stimulation protocols, so that especially cardiac cells/tissue remain intact during experimental resting periods.

**LABTRAX-MDAC** is also well suited for other software triggered instruments or as a stand-alone general data recorder for selected WPI Instruments.



The back panel of the Lab-Trax-8/16 has four analog outputs, digital inputs or outputs, a USB port, power socket and power switch.

## LAB-TRAX-8/16 SPECIFICATIONS

ANALOG INPUTS	8 BNC connections
INPUT RANGE	± 10V
SYSTEM NOISE	< 1 mV RMS
ISOLATION	1,500V
OPERATING CURRENT	800 mA maximum
ANALOG OUTPUTS	4 BNC connections
OUTPUT RANGE	± 10V
IMPLEMENTED FILTER	5 <sup>th</sup> order low-pass Bessel filter with 3dB cut-off frequency
OUTPUT IMPEDANCE	100Ω
OUTPUT CURRENT	15 mA
DIGITAL I/O	16/16 TTL (BNC or DB-9 Connector)
LOGIC HIGH VOLTAGE	3.3 V minimum
LOGIC LOW VOLTAGE	1.0 V maximum
ANALOG & DIGITAL INPUTS	Operating voltage protected to ±30V
PC INTERFACE	USB 2.0
RESOLUTION	16 bits
POWER SOURCE	12V DC

## ORDERING INFORMATION

<b>LABTRAX-MDAC</b>	Lab-Trax-8/16 with MDAC software
<b>LAB-TRAX-4</b>	4-Channel General Data Acquisition System
<b>2851</b>	BNC to BNC Cable

# Micromanipulators

Micromanipulators are used when precision work is conducted under a microscope. A micropipette, electrode or probe can be mounted on a micromanipulator and move as little as a micron at a time. This tool can be used for *in vitro* fertilization, patch clamp experimentation, extracellular recording, microinjection and any application requiring fine mechanical placement (resolution).

In addition to micromanipulators, WPI offers tilt bases, piezo translators and a variety of stands.

Micromanipulators can be broken out into three broad categories

- Manual
- Manual/Motorized
- Motorized



## MICROMANIPULATOR COMPARISON

Micromanipulator	Manual or Motorized	Resolution	Travel	Stands	Tilt Base	Piezo Translator	Notes
<b>SU-TRIO235</b>	Motorized	< 100 nm	X: 25 mm Y: 25 mm Z: 25 mm Diagonal: 50 mm	Fixed platform stage with imperial/standard holes, chamber insert and gantry supports	N/A	N/A	Carries up to a kilogram
<b>SU-TRIO245</b>	Motorized	< 100 nm	X: 25 mm Y: 25 mm Z: 25 mm	Fixed platform stage with imperial/standard holes, chamber insert and gantry supports	N/A	N/A	Carries up to a kilogram
<b>SU-QUAD</b>	Motorized	< 100 nm	X: 25 mm Y: 25 mm Z: 25 mm Diagonal: 30 mm	Fixed platform stage with imperial/standard holes, chamber insert and gantry supports	N/A	N/A	Quiet mode eliminates electrical noise
<b>M3301</b>	Manual	0.01 mm (X fine) 0.1 mm (X, Y, Z)	X (fine): 10 mm X: 37 mm Y: 20 mm Z: 25 mm	M9, M10, M10L, 501622, 501623	TBS, M-3	MPM20	
<b>KITE</b>	Manual	0.1 mm	X (fine): 10 mm X: 35 mm Y,Z: 20 mm	M9, M10, M10L, 501622, 501623	TBS, M-3		
<b>M325</b>	Manual	10 µm	X: 25 mm Y,Z: 10 mm	M9, M10, M10L, 501622, 501623	TBS, M-3		
<b>MMJ</b>	Manual	0.1 mm	X: 37 mm Y: 20 mm Z: 25 mm	M9, M10, M10L, 501622, 501623	TBS, M-3		Joystick control
<b>MD4</b>	Manual	10 µm (X fine) 100 µm (X, Y, Z)	X (fine): 10 mm X: 37 mm Y: 20 mm Z: 25 mm	M9, M10, M10L, 501622, 501623	TBS, M-3		Holds two electrodes
<b>MM1</b>	Manual (mini)	1.0 µm	3 mm				225 g load
<b>MM3</b>	Manual (mini)	1.5 µm	13 mm				340 g load
<b>MM1-3</b>	Manual (mini)	1.0 µm	3 mm	MB-2			225 g load
<b>MM3-3</b>	Manual (mini)	1.5 µm	13 mm	MB-2 with MM1A (adapter)			340 g load

### WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
Brazil: (013) 406-29703 • info@brazil.wpiinc.com • www.wpiinc.com

Germany: +49 (0)6031 1602171 • wpide@wpi-europe.com • www.wpi-europe.com  
China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net

# Manual Micromanipulators

Popular manual micromanipulator

Economy manual micromanipulator



**M3301L**

**M-10**  
magnetic  
stand sold  
separately.

Left or right-handed versions of the **M3301** are supplied with a standard 12mm clamp (**M2**) and one microelectrode holder (**M3301EH**).

## Features

- The most widely used micromanipulator
- Lightweight 550 g
- Sure, repeatable movement without drift

## Benefits

- Control knobs clustered in 8 cm area in a single vertical plane for quick resolution
- Right and left hand orientation options available

## Applications

- Microinjection
- Electrophysiology recording

Weighing just 550 g and employing a slim space-saving design, this well-built micromanipulator outsells all others worldwide for high precision experiments where magnification is in the range of up to 250x. Its design allows units to stand tightly grouped, since all control knobs project to the rear. And because control knobs are clustered within an 8 cm area in a single vertical plane, resolution is quick. The hand works blindly while the eye monitors the microscopic image. Vernier scales allow readings to 0.1 mm; X-axis fine control allows readings to 10 μm.

The instrument employs rack-and-pinion drive, V-shaped guide ways and cross roller bearings, so all movement is sure and repeatable, without drift, side play, backlash or sticking. Contact parts are milled of hardened steel for high performance and long life.

### M3301 SPECIFICATIONS

	TRAVEL RANGE	RESOLUTION
X-axis Fine	10 mm	0.01 mm
X-axis	37 mm	0.1 mm
Y-axis	20 mm	0.1 mm
Z-axis	25 mm	0.1 mm

### ORDERING INFORMATION

<b>M3301R</b>	Manual Manipulator (right-handed)
<b>M3301L</b>	Manual Manipulator (left-handed)
<b>M3301-M3-R</b>	Manual Manipulator (right-handed) & Tilting Base*
<b>M3301-M3-L</b>	Manual Manipulator (left-handed) & Tilting Base*
<b>502105</b>	Axis Adjustment Tool

\*Requires 5464 weight



**KITE-R**

## Features

- Vernier scales allow readings to 0.1 mm
- X-axis fine control allows readings to 10 μm
- Choice of optional **M3** Tilting base which can be mounted to a table with M6 screws
- Left or right-handed versions of the **KITE** micromanipulator are supplied with a standard 12 mm clamp and electrode holder **M3301EH**

## Benefits

- Control knobs clustered in the same plane for quick resolution
- Right and left hand orientation options available

## Applications

- Microinjection
- Electrophysiology recording

The **KITE** is an economical manual micromanipulator that is perfect for student use. It is available in right or left hand versions and is also sold bundled with an optional **M3** Tilting Base. Our 5 lb. weight (WPI #5464) is frequently sold with this unit to provide stability.

### KITE SPECIFICATIONS

	TRAVEL RANGE	RESOLUTION
X-axis Fine	10 mm	0.01 mm
X-axis	35 mm	0.1 mm
Y-axis	20 mm	0.1 mm
Z-axis	20 mm	0.1 mm

### ORDERING INFORMATION

<b>KITE-R</b>	Kite Manual Manipulator (right-handed)
<b>KITE-L</b>	Kite Manual Manipulator (left-handed)
<b>KITE-M3-R</b>	Kite (right-handed) + Tilting Base Combo*
<b>KITE-M3-L</b>	Kite (left-handed) + Tilting Base Combo*

\*Requires 5464 weight



### OPTIONAL ACCESSORIES

<b>M3301EH</b>	Replacement Electrode Holder (14 cm long)
<b>15873</b>	Optional Angled Electrode Holder (13 cm long)
<b>M-3</b>	80° Tilting Base M6 x 1 mm screw
<b>5464</b>	5-lb Weight for Tilting Base (Shipping weight: 7 lb (3 kg))
<b>500475</b>	Ball Joint, 7 cm long, for Ø 8 mm Holder
<b>500476</b>	Ball Joint, 4 cm long, for Ø 4 mm Holder
<b>M4C</b>	Microscope Stage Adapter

Also see magnetic stands (page 28)

# Miniature Micropositioners

One and three axis positioners for small jobs



## MM1 and MM1-3

### Features

- Precise and smooth motion
- Less than 1  $\mu\text{m}$  maximum wobble

### Benefits

- Compact
- Provides precise and smooth motion with no backlash

### Applications

- Oocyte injection with Nanoliter Injector

Single stage measures only 5 × 11 × 26 mm with 3 mm travel. Provides precise and smooth motion with no backlash, positive spring loaded carriage, straight within 1  $\mu\text{m}$  and less than 1  $\mu\text{m}$  maximum wobble. Features fine 80 TPI screw adjustment. 10 mm square mounting surface has a 3.9 mm tapped center hole for transmission and/or mounting. Available in single X (**MM1**) and X-Y-Z (**MM1-3**) axis configurations.

## MM3 and MM3-3

### Features

- Opens 0–4mm
- Precise and smooth motion
- Less than 1  $\mu\text{m}$  maximum wobble

### Benefits

- Compact
- Provides precise and smooth motion with no backlash

### Applications

- Oocyte injection with Nanoliter Injector

Single stage measures only 7 × 17 × 44 mm with 13 mm travel. Offers precise and smooth motion with no backlash, positive spring-loaded carriage, straight within 1.5  $\mu\text{m}$ , and less than 1.5  $\mu\text{m}$  maximum wobble. Features fine 80 TPI screw adjustment. 13 mm square mounting surface has a 7mm tapped center hole for transmission and/or mounting. Available in single X (**MM3**) and X-Y-Z (**MM3-3**) axis configurations.

### MINI-MICROPOSITIONER SPECIFICATIONS

	MM1	MM1-3	MM3	MM3-3
AXIS	X	X-Y-Z	X	X-Y-Z
STRAIGHT LINE ACCURACY	Within 1 $\mu\text{m}$ over 3 mm travel	Within 1 $\mu\text{m}$ over 3 mm travel	Within 1.5 $\mu\text{m}$ over 13 mm travel	Within 1.5 $\mu\text{m}$ over 13 mm travel
CLEAR APERTURE	3.9 mm tapped hole, 8-32 thread	3.9 mm tapped hole, 8-32 thread	7 mm tapped hole, 5/16-16 thread	7 mm tapped hole, 5/16-16 thread
LOAD CAPACITY	255 g Normal	255 g Normal	340 g Normal	340 g Normal
FINISH	Black Anodized	Black Anodized	Black Anodized	Black Anodized
WEIGHT	3 g/axis	12 g/axis	14 g/axis	48 g/axis
TYPE	Fine Screw	Fine Screw	Fine Screw	Fine Screw
TRAVEL	3 mm	3 mm	13 mm	13 mm

### ORDERING INFORMATION

<b>MM1</b>	Mini Micropositioner, one axis, 3 mm travel
<b>MM1-3</b>	Mini Micropositioner, three axes, 3 mm travel
<b>MM1-A</b>	Mounting Adapter for MM1 and MM1-3
<b>MM1-C</b>	Clamp for MM1 and MM1-3
<b>MM3</b>	Micropositioner, one axis, 13 mm travel
<b>MM3-3</b>	Micropositioner, three axes, 13 mm travel
<b>MM3-A</b>	Mounting Adapter for MM3 and MM3-3
<b>MM3-C</b>	Clamp for MM3 and MM3-3
<b>MM3-ALL</b>	Complete 3-Axis Micropositioner & Magnetic Stand
<b>MM1-ALL</b>	Complete 3-Axis Mini Micropositioner & Magnetic Stand

### WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
Brazil: (013) 406-29703 • info@brazil.wpiinc.com • www.wpiinc.com

Germany: +49 (0)6031 1602171 • wpide@wpi-europe.com • www.wpi-europe.com  
China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net

# TRIO™ 3-Axis Manipulators

Highly stable, motorized, 3-axis manipulator with 25 mm of movement on each axis



SU-TRIO245



SU-TRIO235

## Features

- Three independent axes - 25 mm orthogonal travel in X, Y and Z
- Sub-micron (less than 100 nm) resolution
- Carries up to a kilogram
- DIP switches on the Control box select direction of movement produced by turn of the Control box knob
- USB interface for computer control
- The TRIO™ comes standard with a universal mounting system suitable for the most popular headstages or pipette holders.

## Benefits

- User selectable angle from 0 - 90° via Control box input
- Fast movement with a top speed of 3mm/sec (while homing)
- Mechanically robust construction for high stability
- Compact, fanless, user-friendly Control box controller preserves bench and rack space
- Push button control of multiple functions - WORK, HOME, LOCK, PULSE and RELATIVE
- Mechanical designs and software algorithms allow the motors to power down completely during recording

## Applications

- Suited for *in vivo* and *in vitro* electrophysiological recording
- Applications with noise-sensitive recordings

The new **SU-TRIO245** is a highly stable 3-axis manipulator with 25 mm of travel on each axis. The TRIO's synthetic 4th axis can be set in software as any angle between 0 and 90° for diagonal movement. Based on a lead-screw design with a smaller overall size and footprint than most manipulators, the SU-TRIO245 is ideal for applications that require 2 pipettes in one setup or for setups where space is limited.

The compact design of the integrated control box requires minimal bench space. It provides quiet, fan-free operation, and is easy to use. No rack mounted controller is required. Position coordinates, in relative or absolute values, are displayed directly on the control box. The TRIO manipulators use a logarithmic acceleration algorithm that eliminates the need for speed selection. As the knobs on the control box are turned faster, acceleration ramps up. This allows for smooth and intuitive motion control of electrode position without the need to stop and change speeds or lift your hand from the knobs. A Y-axis lockout function (accessible by DIP switch) is also available, allowing X/Z-only axial movement during HOME and WORK repositioning.

An alternate configuration, the SU-TRIO235, removes the Z axis and replaces it with an adjustable diagonal axis. This configuration pairs a traditional X and Y axis with a diagonal axis that can be adjusted in angle with a set screw. The diagonal features 50 mm of travel with X and Y axes that have 25 mm of travel. On this model, we have added software to create a synthetic "Z-axis" by combining the diagonal and X axes in reverse of how we create a synthetic diagonal axis on our other manipulators. The SU-TRIO235 and SU-TRIO245 both feature the same precision bearing and lead-screw design.

### TRIO SPECIFICATIONS

TRAVEL	SU-TRIO235:	25 mm on X and Y, 50 mm on diagonal
	SU-TRIO245:	25 mm on X, Y and Z axes
CONTROL BOX		5.5 x 7.5 x 4" (14 x 19 x 10.2 cm)
WEIGHT		2.2 lbs. (1.0 kg)
POWER		115/230 V, 50/60 Hz

### ORDERING INFORMATION

<b>SU-TRIO235</b>	3-Axis Motorized Micromanipulator with X, Y, D Control
<b>SU-TRIO245-L</b>	3-Axis Motorized Micromanipulator with X, Y, Z Control (left-hand)
<b>SU-TRIO245-R</b>	3-Axis Motorized Micromanipulator with X, Y, Z Control (right-hand)
<b>505069</b>	Z-Axis Vertical Extension for Quad/Trio
<b>505071</b>	Rotating Base for MP225 Mechanical
<b>505073</b>	4-in. Dovetail Extension
<b>505074</b>	Mounting Adapter Plate
<b>505085</b>	Rod Holder for Motorized Manipulator

*Basic systems include the manipulator, controller, rod holder, 4-inch dovetail extension, mounting adapter plate, Z-axis vertical extension, cables, and power supply*

# QUAD® 4-Axis Motorized Manipulator

4-axis motorized manipulator with 25 mm of movement on each axis



SU-QUAD-R

## Features

- Four independent axes - 30 mm travel in diagonal for coaxial pipette movement, 25 mm travel in X, Y and Z
- Sub-micron 100 nm resolution
- Quiet mode eliminates electrical noise
- True diagonal assures coaxial movement
- USB interface for computer control

## Benefits

- Compact, fanless, user-friendly control box preserves bench and rack space
- Push button control of multiple functions - WORK, HOME, LOCK, PULSE and RELATIVE
- Robotic HOME and WORK position moves for easy automated pipette exchange
- Display indicates coordinates in relative or absolute
- Five conveniently located control buttons
- When you are ready to begin collecting data, the motor drive electronics can be suppressed

## Applications

- Suited for *in vivo* and *in vitro* electrophysiological recording
- Cell penetration
- Applications with noise-sensitive recordings

The QUAD® motorized micromanipulator is easy to use and features four independent axes. Each axis has a 25 mm range of motion, a digital display of position and a control box. The compact, intuitive controller takes up minimal bench space, is fan free and is easy to use. Three axes provide the X, Y, Z-orthogonal motion typical of most motorized micromanipulators. In the QUAD® a true fourth axis moves the electrode coaxially at exactly the desired angle of approach. The fourth axis has 30 mm of travel that significantly extends the range of travel for the system.

The QUAD® is very easy to operate. It includes a control box with digital display of position. The control box inputs for each axis allow facile manual control of electrode position. Five conveniently located buttons control all of the basic functions that you will need in normal operation. Press and hold the **WORK** button to quickly store a work position; press **WORK** again and the manipulator will return to the same location. **HOME** sends the manipulator to an initial location or any user defined position, which is useful for changing electrodes rapidly. When you are ready to begin collecting data the motor drive electronics can be suppressed by pressing **LOCK**. Display coordinates can be zeroed at any location by pressing **RELATIVE**; go back to absolute coordinates by pressing **RELATIVE** again. Finally, **PULSE** activates a pulse movement mode that produces small, rapid bursts of motion that can be advantageous for sharp electrode cell penetration.

## QUAD SPECIFICATIONS

TRAVEL	30 mm on the diagonal 4th axis 25 mm on the X, Y and Z axes
CONTROL BOX	5.5 x 7.5 x 4 in. (14 x 19 x 10.2 cm)
WEIGHT	2.2 lb. (1.0 kg)
POWER	115/230 V, 50/60 Hz

## ORDERING INFORMATION

<b>SU-QUAD-L</b>	4-Axis Motorized Micromanipulator (left-hand)
<b>SU-QUAD-R</b>	4-Axis Motorized Micromanipulator (right-hand)
<b>505069</b>	Z-Axis Vertical Extension for Quad/Trio
<b>505070</b>	Mounting Adapter Plate for Quad
<b>505071</b>	Rotating Base for MP225 Mechanical
<b>505073</b>	4-inch Dovetail Extension
<b>505085</b>	Rod Holder for Motorized Manipulator


## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
Brazil: (013) 406-29703 • info@brazil.wpiinc.com • www.wpiinc.com

Germany: +49 (0)6031 1602171 • wpide@wpi-europe.com • www.wpi-europe.com  
China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net

# Specialty Micromanipulators

## Joystick Controlled Manual Micromanipulator



M3301EH  
electrode holder  
included

MMJR (right-handed controls)

### Features

- X axis can be tilted 90°
- Easy steering motion that translates normal hand movement into smooth sub-millimeter shifts

### Benefits

- Fine adjustment for the X and Y axes can be controlled by the joystick
- Joystick manipulation

### Applications


- Oocyte injection with Nanoliter Injector

### ORDERING INFORMATION

<b>MMJR</b>	Joystick Micromanipulator (right-handed)
<b>MMJL</b>	Joystick Micromanipulator (left-handed)

Options and replacement parts listed below

## Micrometer Slide Micromanipulator



The M325 can be configured from left- to right-handed version in 60 seconds!

Toolholder rotates around two axes, can hold tools from Ø3.0 mm to Ø12.7 mm.

M3301EH  
electrode holder  
included

Mounts onto horizontal or vertical rods — any diameter from 10 mm to 12.7 mm.

### Features

- Built of precision micrometer-actuated linear slides
- Can be configured for right or left hand use
- Resolution is 10 µm

### Benefits

- Each slide is comprised of a large micrometer head and a spring-return linear slide
- Micrometer head is graduated in 10 µm steps which enables repeatable positioning to an accuracy of 2 µm

### Applications

- Electrophysiology recording and injection

### ORDERING INFORMATION

<b>M325</b>	3-Axis Fine Controlled Manual Micromanipulator
-------------	--

Options and replacement parts listed below

## Dual Tool-Holder Micromanipulator



M3301EH  
electrode holders  
included

MD4L with left-  
hand controls

M-3 Tilting Base sold separately.

### Features

- Scales allow coarse adjustment readings with an accuracy of 100 µm
- X-axis fine control is achieved with a micrometer screw

### Benefits

- Dual electrode holders

### Applications

- Electrophysiology differential amplification

### ORDERING INFORMATION

<b>MD4R</b>	Double-Holder Micromanipulator (right)
<b>MD4L</b>	Double-Holder Micromanipulator (left)
<b>MD4-M3-R</b>	Double-Holder Micromanipulator (right) + Tilting Base
<b>MD4-M3-L</b>	Double-Holder Micromanipulator (left) + Tilting Base

### OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>M3301EH</b>	Replacement Electrode Holder (14 cm 5 Ø 7.2 mm)
<b>15873</b>	Angled Electrode Holder (13 cm long)
<b>M4C</b>	Microscope Stage Adapter
<b>500475</b>	Ball Joint, 7 cm long, for Ø 8 mm Holder
<b>500476</b>	Ball Joint, 4 cm long, for Ø 4 mm Holder

# Magnetic Holding Devices



M9

## M9

Mechanical clamp tightens three rotatable joints simultaneously with one locking knob. Arm adjusts without distortion. Base exerts 100kg magnetic force for great stability. Fine adjustment for precise operations included.

### Magnetic Base:

50 (w) x 60 (l) x 55 (h) mm  
(2.2 x 2.4 x 2.2 in.)

### Vertical Holding Power:

100 kgf (220 lb force)

### Arms:

L1: 119 mm (4.7 in.)  
L2: 106 mm (4.2 in.)  
L3: 25 mm (0.98 in.)  
Ø 12 mm (0.472 in.)

### Clamp Hole:

none

### Weight:

1.8 kg (4 lb.)

### ORDERING

M9 Magnetic Stand



M10

## M10

Similar to M1 with a 12mm diameter sub pole (fits 12mm clamp supplied with M3301, DC3001, MD4 and MMJ manipulators).

### Magnetic Base:

50 (w) x 58 (l) x 55 (h) mm  
(2.0 x 2.3 x 2.2 in.)

### Vertical Holding Power:

80 kgf (176 lb force)

### Main Pole:

diameter: 14 mm (0.55 in.)  
length: 178 mm (7 in.)

### Sub Pole:

diameter: 12 mm (0.47 in.)  
length: 165 mm (6.5 in.)

### Clamp Hole:


Adjustable from 4.5 mm to 6.5 mm

### Weight:

1.8 kg (4 lb.)

### ORDERING

M10 Magnetic Stand



M10L

## M10L

Same as M10, but equipped with a taller (14-inch) vertical main pole.

### Magnetic Base:

50 (w) x 58 (l) x 55 (h) mm  
(2.0 x 2.3 x 2.2 in.)

### Vertical Holding Power:

80 kgf (176 lb force)

### Main Pole:

diameter: 14 mm (0.55 in.)  
length: 356 mm (14 in.)

### Sub Pole:

diameter: 12 mm (0.47 in.)  
length: 165 mm (6.5 in.)

### Clamp Hole:

Adjustable from 4.5 mm to 6.5 mm

### Weight:

1.8 kg (4 lb.)

### ORDERING

M10L Magnetic Stand



M11

## M11

Useful as a probe holder, the M11 is not suitable for heavier items such as micromanipulators.

Bends freely for maximum flexibility. The connecting arm twists and bends like a snake. Lock the arm in position with a flick of the controlling lever.

### Magnetic Base:

50 (w) x 58 (l) x 55 (h) mm  
(2.0 x 2.3 x 2.2 in.)

### Vertical Holding Power:

80 kgf (176 lb force)

### Main Pole:

diameter: 16 mm (0.63 in.)  
length: 315 mm (12.4 in.)

### Sub Pole:

none

### Clamp Hole:

Adjustable from 6 mm to 8 mm

### Weight:

1.4 kg (3 lb.)

### ORDERING


M11 Magnetic Stand

# Powerful Ball Joint Rare Earth Magnet

## Construct holding devices for small parts/equipment

Small but very powerful. Holds 2 kilograms (~5 pounds)!  
Steel ball rotates freely 360° on a 180° axis  
M3 mounting screw on ball for attachment to equipment  
Magnet base threaded (M3) for mounting onto a base or equipment

This novel magnetic ball joint has phenomenal holding power for up to 2kg of attached weight while permitting the ball a full 360° rotation on a 180° axis. You can freely orient your equipment to an infinite number of positions within this rotation. This is made possible by the combination of a steel ball (10



mm diameter) and a powerful rare earth magnet contained in the magnet cylinder (Ø10 x 20 mm). Convenient M3 attachment sites are provided on both the ball (male) and the magnet base (female). For use with micromanipulators for the positioning and holding of optical instruments including various lighting sources and lasers, pipettes and any small parts that would benefit from the flexibility offered by this new magnetic ball joint.



500871

### ORDERING INFORMATION

500871 Magnetic Ball Joint


## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
Brazil: (013) 406-29703 • info@brazil.wpiinc.com • www.wpiinc.com

Germany: +49 (0)6031 1602171 • wptide@wpi-europe.com • www.wpi-europe.com  
China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net



The base of each stand exerts a powerful magnetic force that holds it solidly on ferrous metal surfaces — even vertically or upside down



## MB2

Mechanical clamping type tightens three joints simultaneously. Arm is freely adjustable without distortion. Equipped with fine adjuster and medium size magnet for stabilizing the base. Suitable for performing precision operation. Fits clamp for **MM1-ALL** and **MM3-ALL**.

**Magnetic Base:**

30 (w) x 35 (l) x 35 (h) mm  
(1.2 x 1.4 x 1.4 in.)

**Vertical Holding Power:**

17 kgf (37 lb. force)

**Arm:**

L1: 46 mm (1.8 in.)  
L2: 46 mm (1.8 in.)  
L3: 39 mm (1.5 in.)

**Clamp Hole:**

Adjustable from 5 to 6.5 mm

**Weight:**

0.38 kg (0.83 lb)

### ORDERING

**MB2** Compact Magnetic Stand

## M8

A ball joint at the base of the main post allows 360° rotation, offering considerable versatility. The second arm adopts angles up to 75°.

**Magnetic Base:**

50 (w) x 58 (l) x 55 (h) mm  
(2.0 x 2.3 x 2.2 in.)

**Vertical Holding Power:**

80 kgf (176 lb. force)

**Main Pole:**

diameter: 12 mm (0.47 in.)  
length: 194 mm (7.6 in.)

**Sub Pole:**

diameter: 10 mm (0.39 in.)  
length: 165 mm (6.5 in.)

**Clamp Hole:**

Adjustable from 4.5 mm to 6.5 mm

**Weight:**

1.8 kg (4 lb.)

### ORDERING

**M8** Magnetic Stand

## M1

A precision base providing stable support for such devices as electrodes and manipulators. Adjustable second arm adopts a variety of angles.

**Base:**

50 (w) x 58 (l) x 55 (h) mm  
(2.0 x 2.3 x 2.2 in.)

**Vertical Holding Power:**

80 kgf (176 lb. force)

**Main Pole:**

diameter: 12 mm (0.47 in.)  
length: 176 mm (6.9 in.)

**Sub Pole:**

diameter: 10 mm (0.39 in.)  
length: 165 mm (6.5 in.)

**Clamp Hole:**

diameter: 4.5 mm and 6.5 mm

**Weight:**

1.8 kg (4 lb.)

### ORDERING

**M1** Magnetic Stand

## M1L

Same base and support arm as M1, but equipped with a longer (14-inch) vertical post.

**Base:**

50 (w) x 58 (l) x 55 (h) mm  
(2.0 x 2.3 x 2.2 in.)

**Vertical Holding Power:**

80 kgf (176 lb. force)

**Main Pole:**

diameter: 12 mm (0.47 in.)  
length: 356 mm (14 in.)

**Sub Pole:**

diameter: 10 mm (0.39 in.)  
length: 165 mm (6.5 in.)

**Clamp Hole:**

diameter: 4.5 mm and 6.5 mm

**Weight:**


1.8 kg (4 lb.)

### ORDERING

**M1L** Magnetic Stand

# Versatile Manipulator Mounting Clamp

Mount your Micromanipulator



M5 — Ø10 mm clamp

Three of the stands above — **M1**, **M1L** and **M8** — have 10 mm diameter mounting rods. The standard mount on several WPI manipulators (**DC3001**, **KITE**, **M3301**, **MMJ**, and **MD4**) accommodates a 12 mm rod. In order to use one of these three stands, you will need to replace the manipulator's standard 12mm mounting clamp with the optional M5 clamp.

### ORDERING INFORMATION

**M5** Ø10 mm Clamp

# Magnetic Holding Devices

*A solid platform for mounting your manipulators*



501651

## Round Base

An ideal accessory for optical tables and vibration-free platform. Reduces experimental set-up time by allowing free positioning and instant clamp down of optical components. Switchable ON/OFF magnetic circuit permits fine adjustment and precise positioning.

- Easy ON/OFF operation using lever
- Thin and powerful magnetic force
- Generous array of tap holes

### Holding Power:

20 kgf (44 lb. force)

### Dimension:

75 (OD) x 20 (h) mm  
2.9 (OD) x 0.8 (h) in.

### Mounting Hole:

4-M4 x 0.7, depth 6 mm\*  
M8 x 1, depth 6 mm  
Span 35 mm

### Weight:

0.7 kg (1.5 lb.)

### ORDERING

**501651** Magnetic Base,  
75 mm diameter

**503568** Magnetic Base,  
50 mm diameter

\* Posts with M4-threads not available from WPI.



503570

## Square Base

An ideal accessory for optical tables and vibration-free platform. Reduces experimental set-up time by allowing free positioning and instant clamp down of optical components. Switchable ON/OFF magnetic circuit permits fine adjustment and precise positioning.

- Easy ON/OFF operation using lever
- Thin and powerful magnetic force
- Generous array of tap holes

### Holding Power:

20 kgf (44 lb. force)

### Dimension:

65 (w) x 65 (l) x 20 (h) mm  
2.6 (w) x 2.6 (l) x 0.8 (h) in.

### Mounting Hole:

8-M4 x 0.7, depth 6 mm\*  
M8 x 1, depth 6 mm  
Span 25 mm

### Weight:

0.6 kg (1.3 lb.)

### ORDERING

**501653** Magnetic Base,  
65x65 mm

**503569** Magnetic Base,  
45x45 mm

**503570** Magnetic Base,  
90x90 mm

**503571** Magnetic Base,  
120x120 mm

\* Posts with M4-threads not available from WPI.



501652

## MOBITY

MOBITY™ is a new magnetic clamping system. With its ease of use, only one hand is needed to operate the attractive power. The MOBITY™ has a strong 88lbf pull, yet weighs only 1.5 lbs. MOBITY™ meets various applications with 4 tapped holes on the top surface. Requires (1) 9V alkaline battery (included).

### Holding Power:

40 kgf (88 lb. force)

### Dimension:

55 (w) x 73 (l) x 50 (h) mm  
2.2 (w) x 2.9 (l) x 2.0 (h) in.

### Mounting Hole:

3-M4, depth 20 mm\*  
M8, depth 15 mm

### Weight:

0.7 kg (1.5 lb.)

### ORDERING

**501652** MOBITY Magnetic  
Clamping System

\* Posts with M4-threads not available from WPI.



M7

## M7

A small holder ideal for use where space is limited. Main post unscrews from base which may then be used alone as a switchable magnetic holder.

### Magnetic Base:

30 (w) x 35 (l) x 35 (h) mm  
1.2 (w) x 1.4 (l) x 1.4 (h) in.

### Vertical Holding Power:

20 kgf (44 lb. force)

### Main Pole:

Diameter: 7 mm (0.28 in.)  
Length: 52 mm (2 in.)

### Clamp Hole:

Diameter: 6 mm

### Weight:


0.36 kg (0.8 lb.)

### ORDERING

**M7** Compact Magnetic  
Stand

## Base Plates

A magnetic stand requires a steel mounting surface. WPI's steel base plates have plenty of mass to give stability to your experimental setup. Beveled edges make them easy to handle. Rubber feet hold them off the benchtop, making them easier to grasp when moving. The special black coating provides a durable, protective, non-reflective finish.



**5479**  
30.5x61x1 cm (12x24x38-in.)



**5052**  
21.6x30.5x1.3 cm (8½x12x½ in.)

### ORDERING INFORMATION

**5052** Steel Base Plate, 4.5 kg (10 lb)

**5479** Steel Base Plate, 14.5 kg (32 lb)

## WORLD PRECISION INSTRUMENTS


UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
Brazil: (013) 406-29703 • info@brazil.wpiinc.com • www.wpiinc.com

Germany: +49 (0)6031 1602171 • wpide@wpi-europe.com • www.wpi-europe.com  
China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net

# Mounting a Manipulator

## Tilting Base with Screw Adjustment

At last! A tilt base you can operate with one hand!



### Features

- High precision and stability
- Tilt angle up to 35°
- 1.6 kg (3.5 lb) steel base provides stability for your micromanipulator
- Holes also allow permanent mounting to your bench top
- Adjust tilt angle by turning knob
- Manipulator mounting bracket included

### Benefits

- 13x15 cm (5x6 in.) footprint saves space in your work area
- Two sets of mounting holes are pre-drilled for WPI manipulators (M3301R shown) but steel platform may be drilled for mounting other devices

### Applications


- Microinjection
- Electrophysiology amplification
- Stimulation

#### ORDERING INFORMATION

<b>TBS</b>	Tilt Base with Screw Adjustment
	Shipping Weight 7 lb.

## Ball-joint holder attachment

Here's a new angle for mounting




#### ORDERING INFORMATION

<b>500475</b>	Ball joint, 7 cm long, for O.D. 5-9 mm Electrode Holder (shown)
<b>500476</b>	Ball joint, 4 cm long, for O.D. 2.8-4.5 mm Electrode Holder

## Manipulator Clamps

Interchangeable clamps allow manipulators to be mounted on a variety of supports.



#### ORDERING INFORMATION

<b>M2</b>	Ø 12 mm Clamp
<b>M5</b>	Ø 10 mm Clamp
<b>M6</b>	Ø 12.7 mm (1/2-in.) Clamp

## Mount a manipulator on your microscope stage



**M4C Microscope Stage Adapter**  
— clamps onto microscope stage to provide a stable support for any manipulator with a Ø12 mm clamp.



The M4C is shown mounted on a microscope stage.

#### ORDERING INFORMATION

<b>M4C</b>	Microscope Stage Adapter
------------	--------------------------

# Environmental Control for Live Cell Microscopy

*Incubation and environmental control for long term imaging of cellular growth*

## Features

- Unique diffusion grid, combined with air input and return vents, provide an air flow pattern for consistent, even heating, with no hot or cold spots in the chamber
- External heater that can be placed far enough from the system to eliminate electrical and vibrational interference from the heater
- High degree of temperature precision and stability
- Minimal focal drift after equilibrium is achieved—accuracy  $\pm 0.1^{\circ}\text{C}$  at the sample itself, and  $0.2^{\circ}\text{C}$  across the microscope stage (allowing for uniform heating of multiwell dishes)
- Airflow pattern and temperature uniformity eliminate dramatic changes in environmental temperature when the incubator door opens
- Ergonomic design for ease of use—The focus and x/y stage controls are outside of the incubator itself. Large doors allow easy access to the specimen and small ones for cords, tubing, etc.
- Precision, shielded temperature probe
- Simple, one person setup of the system

## Benefits


- Choice of controller to manage air flow, heat, carbon dioxide and oxygen
- Control and monitor an external heater or lens warmer with **ECU or control and monitor temperature and monitor humidity with ECU-HOC**
- Diffuser grid and proper air venting insures consistent air flow inside the acrylic microscope chamber

## Applications


- Live cell imaging under a microscope
- Time lapse video research

This Live Cell Microscope Incubator was extensively tested in laboratories. When compared with other systems, it offers dramatic advantages. For example, other incubators for live cell microscopy rely on passive, random

diffusion of heated air from a single source to maintain the desired temperature setpoint. With no hot air return vent, the heated air escapes from the system through cracks at the microscope/incubator junction in an uncontrolled, random fashion. These systems offer no temperature uniformity, suffer from focus drift and often experience electrical and vibrational interference from the heater. You will also notice dramatic temperature drifts when the imaging environment is disturbed.



*Incubator on a Nikon Eclipse TE-2000 Microscope. Incubators for Zeiss and Olympus microscopes are also available, as well as versions that accommodate confocal modules. All incubators are compatible with all commercially available cameras, light sources, filter wheels, motorized stages and motorized nose pieces.*



*This diagram shows how the Stagetop environment connects with the ECU-HOC to control heat, air flow, carbon dioxide and oxygen.*

\*Aux. heater could be a lens warmer or an AirTherm Satellite heater (for a full microscope enclosure).

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
Brazil: (013) 406-29703 • info@brazil.wpiinc.com • www.wpiinc.com

Germany: +49 (0)6031 1602171 • wpide@wpi-europe.com • www.wpi-europe.com  
China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net

## Choice of controller to manage heat, CO<sub>2</sub> and O<sub>2</sub>

This unique, acrylic Live Cell Imaging chamber, combined with an Environmental Control Unit (ECU) and an AirTherm controller, ensure precision control of your incubator environment.

The ECU comes in different varieties so you have all the control you require.

- With the **ECU-H5**, you can control air flow and heating.
- In addition to air flow and temperature control, the **ECU-HC** lets you control the carbon dioxide level. It has an internal sensor.
- The **ECU-HOC** adds control of the oxygen level, which is accomplished by displacing the oxygen with nitrogen.

## Control external heater and monitor humidity

The first four ECU units are capable of controlling a simple, external heater, like the AirTherm Satellite (AirTherm-SAT) or a microscope lens warmer. The AirTherm-SMT can monitor and control temperature and monitor humidity level inside the microscope chamber.

## Diffuser grid and venting for consistent air flow

Air flow affects the temperature uniformity of incubators. The red arrows on Fig. 1 and Fig. 3 indicate air flow. The Live Cell Microscope Incubator uses a diffuser grid and proper venting to insure consistent air flow. Traditional incubators with poor air flow suffer with hot and cold spots in the incubator, as seen in thermal images (Fig. 2 and Fig. 4). Warmer temperatures are indicated by red and cooler temperatures by blue.




Fig. 1—Single air input and no venting causes random air flow in a traditional incubator.




Fig. 2—Hot and cold spots result from inconsistent flow.




Fig. 3—A diffusion grid with air input and exhaust vents yields consistent air flow.



Fig. 4—Consistent air flow means uniform heating.



CE

AIR-THERM-SMT



AIRTHERM-SAT

CE



This humidifier may be used inside a microscope chamber for controlling the humidity inside the stagetop environment chamber.

## ORDERING AN ENVIRONMENTAL CHAMBER

Acrylic enclosures are essentially custom-built. When ordering a system, you will need to provide the following information:

Microscope	Stage	Stage-Up
Perfect Focus	Camera	Left Port
Right Port	Analyzer	Fluor Attachment
Tirf	White Light Tirf	Binocular D Head
Tilting Head	Filter Wheels	Excitation
Emission	Dual Lamphouse	Transmitted Light Shutter
Cells	35 mm/60 mm	Wells
	Coverslips	

## ORDERING INFORMATION

### FULL CHAMBER SYSTEM WITH ECU CONTROLLER & NEW AIRTHERM-SAT/SMT

*All Systems Include: Proprietary Humidification Module, Stage Adapter, Stage Dish with Optical Grade Glass*

<b>FC-HC</b>	Microscope Environmental Chamber, Heat & CO <sub>2</sub> Controller <i>Requires 100% CO<sub>2</sub> &amp; ambient air supplies.</i>
<b>FC-HOC</b>	Microscope Environmental Chamber, Heat Controller, CO <sub>2</sub> and O <sub>2</sub> Controller
<b>FC-H5</b>	Microscope Environmental Chamber, Heat Controller, Digital Control for Bottle Gas
<b>FC-SMT</b>	Microscope Environmental Chamber, Heat Controller (AIRTHERM-SMT instead of ECU) <i>Requires pre-mixed 5% CO<sub>2</sub> gas supply</i>
<b>300626</b>	Optional Humidifier for AirTherm-SMT

# Stagetop Environmental Control

## Control temperature and CO<sub>2</sub> in a microscope stagetop environment

### Features

- Four programmable digital control loops:
  - Independent incubator base temperature PID control with  $\pm 0.1^\circ\text{C}$  precision
  - Independent incubator lid temperature PID control with  $\pm 0.1^\circ\text{C}$  precision
  - CO<sub>2</sub> digital PID control with  $\pm 0.1\%$  precision
  - Airflow digital PID control from 0–900 SCCM
- USB-based remote control and data logging
- Electronic flow meter
- Programmable alarm for out of tolerance condition on all four channels
- Compact and lightweight

### Benefits

- Compact housing that fits most inverted microscope stages and holds standard culture well plates
- Control temperature and CO<sub>2</sub>, O<sub>2</sub> using the Environmental Control Unit (ECU)
- Control and monitor system parameters using the ECU

### Applications

- For short term or long term studies of living cell cultures under a microscope (Live Cell Imaging)
- Time lapse video research

Time lapse video research requires a microscope stagetop incubator.

#### Compact housing fits most inverted microscope stages

Perfect for live cell imaging, **STEV2** (the stagetop environmental control platform) is a compact environmental case that houses your culture wells and fits on a microscope stage inside the live cell microscope incubator.

#### Control temperature and CO<sub>2</sub> and O<sub>2</sub>

This system offers precision control of both temperature and carbon dioxide, as well as remote control and data logging via a USB connection. The system is flexible and easy to configure for a variety of experimental conditions.

#### Control and monitor system parameters

The system includes the Environmental Control Unit electronics which use four programmable loops to control the temperatures of the case and the lid, CO<sub>2</sub> within the environmental case and airflow within the incubator.



The Environmental Control Unit home screen shows the real time readings for the parameters you want to see.

### ECU SPECIFICATIONS

Operating Temperature (ambient)	10 – 50°C (50 – 122°F)
Operating Humidity (ambient)	15 – 70% RH, non-condensing
Warm up Time	20 minutes
Computer Interface	USB via external USB/RS232 converter
Sensor	Non-dispersive infrared (NDIR), dual beam, 20s response time
Power	110/240V, 50/60Hz
CO <sub>2</sub> Sensor (ST-HC only)	
Sensor Range	0 – 20% CO <sub>2</sub>
Control Range	0 – 20%
Control Precision	0.1% CO <sub>2</sub>
Control Accuracy	0.1 – 3% of reading
Drift	<2.5% reading/year
O <sub>2</sub> Sensor (ST-HOC only)	
Sensor Type	Zirconium Dioxide, diffusion, 4 sec response time
Sensor Range	0–25%
Control Range	0–25%*
Control Precision	0.1% O <sub>2</sub>
Control Accuracy	$\pm 0.5\%$ (2% of the full scale)

### ORDERING INFORMATION

#### STAGETOP ENVIRONMENTAL CONTROL SYSTEMS

<b>ST-H5</b>	Chamber, Controller, Heat, Digital Flow Control for Bottle Gas
<b>ST-HC</b>	Chamber, Controller with CO <sub>2</sub> & Heat (using CO <sub>2</sub> internal sensor)
<b>ST-HOC</b>	Chamber, Controller with CO <sub>2</sub> , O <sub>2</sub> & Heat (using CO <sub>2</sub> internal sensor)

#### ACCESSORIES

<b>AIRTHERM-SAT-1W</b>	AirTherm Satellite Heater (110V)
<b>AIRTHERM-SAT-2W</b>	AirTherm Satellite Heater (230V)

### WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
Brazil: (013) 406-29703 • info@brazil.wpiinc.com • www.wpiinc.com

Germany: +49 (0)6031 1602171 • wptide@wpi-europe.com • www.wpi-europe.com  
China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net

# Heat Controller/Humidity Monitor

Smart, electrically quiet air heater for live cell imaging & custom incubators

## Features

- Precision heat controller for use in live cell imaging and custom incubators
- Control heat and monitor humidity (optional) with a single controller
- Electrically and acoustically quiet
- Quick, precise response to thermal change

## Benefits

- PID control algorithms allow for tight control of temperature in the environmental chamber
- System includes hoses and probes so that it is ready to connect with your microscope chamber, incubator or controlled environment

## Applications

- Live cell imaging systems fitted with full microscope environmental chamber enclosures
- Control environment in custom incubators

The **AirTherm™ SMT** is a new generation of heat control system from WPI designed to be used in Live Cell Imaging applications with microscopes fitted with a full microscope environmental chamber enclosure. The standard **AirTherm™ SMT** controls temperature and, as an option, monitors humidity.

### PID algorithm for maintaining tight control

The **AirTherm™ SMT** uses a PID control algorithm to tightly control the temperature and monitor humidity of the controlled environment


With **AirTherm™ SMT**, the temperature of the sample and microscope optics can be controlled to within 0.2°C. During operation, air is drawn out of the chamber through a flexible hose, heated by the **AirTherm™ SMT** heater and re-circulated to the chamber by the return hose.

### Ready to connect to your microscope chamber

The system is typically used in a closed loop configuration.

The **AirTherm™ SMT** system includes:

- Two coil-reinforced heater hose pieces and hose clamps.
- Temperature sensor for remote placement in environmental chamber
- Humidity probe for monitoring chamber humidity available as an option.



A typical **AirTherm™ SMT** installation places the heated air inflow at the bottom and the cold air return at the top of the microscope chamber.



AIR-THERM-SMT

CE

## References

X. Yin, D.A. Knecht, M.A. Lynes "Meallothionein mediates leukocyte chemotaxis" *BMC Immunology* 6. 2005: 6--21

## AIRTHERM SPECIFICATIONS

AIR FLOW RATE	20–50 CFM (0.55–1.4 m <sup>3</sup> /minute)
CONTROL TEMPERATURE RANGE	Ambient to 45°C
TEMPERATURE RESOLUTION	0.1°C
TEMPERATURE ACCURACY	0.2°C
ANALOG OUTPUT FOR CHART RE-CORDER	0.5°C resolution; 0–10 V represents 0–100°C
HEATING VOLUME	Less than 50 CF (1400 L), re-circulating
TEMPERATURE SENSOR TYPE	Platinum RTD 1000 Ω
HUMIDIFIER TYPE	Ultrasonic
HUMIDIFIER TANK CAPACITY	0.5 gallons
HUMIDIFIER DAILY OUTPUT	2 gallons
FUSE	For 120 VAC, 8A 250 V 5x20 mm metric For 230 VAC, 4A 250 V 5x20 mm metric
POWER	450 W, 95–135 V or 220–240 V, 50/60 Hz
DIMENSIONS	6½ x 8 x 7½ in. (15.5 x 21 x 19 cm)

## ORDERING INFORMATION

**AIRTHERM-SMT-2W** AirTherm™ SMT Heater, 230V  
**AIRTHERM-SAT-2W** AirTherm Satellite Heater, 230V

## OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>15590</b>	Replacement Hoses, 2.5" diam., 4.5 ft
<b>98727</b>	Replacement Temperature Probe
<b>98728</b>	Humidity Probe
<b>300626</b>	AirTherm SMT Humidifier Kit with Sensor & Hose (110 V)
<b>300411</b>	Replacement Humidifier for AirTherm SMT (110 V)
<b>5389</b>	Inlet and Outlet Hoses
<b>96824</b>	3-Way Hose Assembly (3 hoses and T-connector)

# Bench Top Vibration Isolation Platforms

*Simple setup and adjustment*

## Features

- Only 4.6" thin
- Portable
- Available in a variety of payload ranges

**MK-BM-8100**, 50-105 lb.  
payload weight range  
(23-48 kg)



## Benefits

- Better performance than an air table
- Simple setup and adjustment
- Requires no air or electricity

## Applications

- Microscopy and other research applications that require the isolation of the preparation from vibrations of any kind

These bench top platforms offers 10-100 times better performance than a full size air table in a package only 4.6 inches tall, and without air or electricity! These vibration isolation platforms are extremely easy to use and offer extreme performance — 1.5 Hz horizontal natural frequency and 0.5 Hz vertical natural frequency. There are only two adjustments. This is the thinnest, most portable and most user-friendly isolator ever offered that is capable of delivering this level of performance.

Weight: Approximately 40 lb. (16 kg)  
Dimensions: 18" W x 20" D x 4.6" H (457 x 508 x 117 mm)

## Performance

- Horizontal frequencies are weight dependent
- Horizontal frequency of 1.5 Hz is achieved at or near the upper limits of the payload range
- At the lower limits of the payload range the horizontal frequency is approximately 2.5 Hz
- Vertical frequency is tunable to 0.5 Hz throughout the payload range

## ORDERING INFORMATION

MODEL	PAYLOAD RANGE
<b>MK-BM-825</b>	Vibration Platform, 10 - 30 lb. (4.5 - 14 kg)
<b>MK-BM-850</b>	Vibration Platform, 25 - 55 lb. (11 - 25 kg)
<b>MK-BM-8100</b>	Vibration Platform, 50 - 105 lb. (23 - 48 kg)
<b>MK-BM-8125</b>	Vibration Platform, 90 - 130 lb. (40 - 59 kg)
<b>MK-BM-8150</b>	Vibration Platform, 125 - 155 lb. (57 - 70 kg)
<b>MK-BM-8175*</b>	Vibration Platform, 150 - 180 lb. (68 - 81.5 kg)
<b>MK-BM-8200*</b>	Vibration Platform, 175 - 205 lb. (79.5 - 93 kg)
<b>MK-BM-8225*</b>	Vibration Platform, 200 - 230 lb. (90.5 - 104 kg)
<b>MK-BM-8250*</b>	Vibration Platform, 225 - 255 lb. (102 - 115.5 kg)

\* Weight: Approximately 47 lb. (21 kg)/(same dimensions)

# Vibration-Free Tables

## Vibration-Free Workstation

### Features

- Vertical and horizontal vibration isolation
- High performance active-air suspension
- Automatic leveling
- VibraDamped steel
- Class 100 clean room compatible
- Leveling feet



Custom order.  
Call for pricing and  
part numbers.

### Benefits

- Eliminate inconsistent and unreliable performance
- Reduce excessive wear, maintenance and fatigue failures
- Protect sensitive instruments and equipment from faulty operation or failure using Vibration-Free Platforms and Workstations

### Applications

- Patch clamping, cell injection, analytical balances and optical microscopes

All buildings vibrate. Activities of people, machinery, heating and ventilation systems, and nearby truck or rail traffic cause all types of vibrations. These vibrations, though acceptable to occupants, cannot be tolerated by equipment used in patch clamping, cell injection, analytical balances and optical microscopes. *Additional tabletop sizes and finishes are available, as well as optional accessories such as side rails and casters. See WPI website for available configurations.*

# Universal Manipulator Stand

*Mount manipulators at angles and heights with infinite flexibility*



501622

## UMS SPECIFICATIONS

### DIMENSIONS

Base Plate	10.0 x 12.5 x 1.5 cm (LxWxH)
Stand	4.0 x 4.0 x 30 cm (LxWxH) (501622) 4.0 x 4.0 x 45 cm (LxWxH) (501623)

### Mounting Holes

- English: 1/4 20 x 1" (2 bolts supplied)
- Metric: M6 x 25 mm grid (2 bolts supplied)

### SHIPPING WEIGHT

<b>501622</b>	9 lb. (4 kg)
<b>501623</b>	11 lb. (5 kg)

Vibration-free Platform (VFP)  
not included

## ORDERING INFORMATION

<b>501622</b>	Universal Micromanipulator Stand 30 cm (includes one clamp)
<b>501623</b>	Universal Micromanipulator Stand 45 cm (includes one clamp)
<b>501624</b>	Additional Rotation Clamp

The **M3301**, **KITE**, **DC3001** and **SM325** (but not **M325**) can be used with these mounts.


## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
Brazil: (013) 406-29703 • info@brazil.wpiinc.com • www.wpiinc.com

Germany: +49 (0)6031 1602171 • wptide@wpi-europe.com • www.wpi-europe.com  
China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net



# Microinjection



MICROINJECTION

## We offer a complete microinjection system

WPI has been in the microinjection business since the very beginning. We were first on the market with pressure injection equipment, and we have been engineering and delivering innovative instruments ever since.

For essential fluidic control in neurophysiology, we offer solutions across pumping platforms including syringe, peristaltic, microfluidic and injection pumps. For picoliter to microliter ranges, our tools are designed for applications like patch clamp, oocyte work, infusion, perfusion, cellular injection, continuous flow and non-pulsatile flow.

With zebrafish (*Danio Rerio*) rapidly gaining in popularity as bio-medical research subjects, WPI meets the need for lab budgets with efficient, cost-effective stations for microinjection and transfection.



# WPI Microinjection System

The system depicted at left includes components often favored by researchers, indicated with \* in the list below: **PV820** Pneumatic PicoPump, **PUL-1000** Micropipette Puller, **M4C** stand, **M3301R** micromanipulator, **5430-XX** PicoNozzle Kit with a  $\mu$ Tip, **PZMTIII** microscope, lighted base, articulating mirror and optional **PRO-300 HDS** camera and view screen, **E2XX** micropipette storage jar, **Z-MOLDS** Microinjection and Transplantation Molds, **14003-G** Vannas spring scissors, glass capillaries, **77020** glass tweezers and Fluorodish optical grade glass bottom culture dishes. Whatever your needs, WPI offers a range of equipment to fill your requirements.

## Options for Customizing Your System

### INJECTORS FOR GLASS PIPETTES

- \* **PV820** Pneumatic PicoPump with Hold Pressure
- **PV830** Pneumatic PicoPump with Hold Pressure and Vacuum
- **NL2010MC2T** Nanoliter Injector with SMARTouch Controller



Designed to simplify intracellular injection and a variety of other microinjection tasks, WPI's **PicoPumps** use accurate timing and carefully regulated air pressures for securing cells and injecting. Injected volumes range from picoliters to nanoliters.



Microprocessor-controlled Nanoliter Injector **NL2010MC2T** with SMARTouch controller uses positive displacement injection, eliminating the need for pipette calibration. System uses glass micropipettes.

### PINPOINT INJECTOR

- \* **MICRO-ePORE™** Cell Penetrator



A simple and versatile pinpoint controlled electroporation system for clean efficient injection of CRISPER reagents into 2-stage embryos. The cell membrane is easily penetrated, without cell disruption or tear drop formation. The technique has been used in animal embryos, including mice, rats and bovine, as well as injecting morpholino oligomers into Zebrafish tails.




### INJECTORS FOR GLASS PIPETTES OR METAL NEEDLES

The versatile UltraMicroPump **UMP3** injector is a Syringe-based micropump using metal needles as small as 36 gauge and microsyringes to deliver picoliter volumes.

A glass pipette injection kit is also available

Use with the NANOFIL syringe and 1mm OD glass for your microinjection application."



#### WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com Germany: +49 (0)6031 1602171 • wpide@wpi-europe.com • www.wpi-europe.com  
Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • www.wpiinc.com China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net

## MICROSCOPES

- \* **PZMTIII** Precision Stereo Zoom Trinocular Microscope
- **PZMIII** Precision Stereo Zoom Binocular Microscope
- **PZMIII-MI** Binocular Microscope Head with LED illuminated base (articulating mirror)

**PZMTIII-MI** Trinocular Microscope Head with LED illuminated base (articulating mirror) and variable light intensity. Dual reflection lens/mirror system provides transmitted brightfield/pseudo-darkfield illumination.



The Trinocular Microscope may be used with the **PRO-300 HDS**.

Photo Courtesy of Benjamin Dubansky, PhD, University of North Texas

*NOTE: #Base can be used with other manufacturer's microscope optics*

## PULLERS

- \* **PUL-1000** Microprocessor-Controlled 4-Step Micropipette Puller
- **SU-P97** Flaming/Brown Pipette Puller
- **SU-P1000** Next Generation Pipette Puller

**PUL-1000** is a microprocessor controlled horizontal puller for making glass micropipettes or microelectrodes used in intracellular recording, microperfusion or microinjection. It offers programmable sequences of up to four steps with heating, force, movement and cooling time. Perfect for long taper pipettes.



*WPI is an authorized distributor for Sutter Instruments*

## MANIPULATORS

- \* **M3301** Manual Micromanipulator
- **KITE** Manual Micromanipulator
- **M4C** Microscope Stage Adapter
- **SU-TRIO** 3-Axis Micromanipulator
- **SU-QUAD** 4-Axis Micromanipulator



Weighing just 550 grams, the **M3301** is a well-built micromanipulator that outsells all others worldwide for high precision experiments where magnification is in the range of up to 250x.




**M4C** Microscope Stage Adapter

## GLASS MICROPIPETTES

- Pre-Pulled Long Taper Glass Pipettes (**TIP15FLT**, **TIP12FLT**)
- Glass Capillaries
- **µTip** Pre-Pulled Pipettes with Calibrated Tip ID



Calibrated Tip ID. Get pre-pulled, pipettes with or without a filament. Available in a variety of sizes. Luer connect option is available.




Eliminate the cost and trouble of pulling your own pipettes. Get new low cost, long taper pipettes with or without a filament. Taper length: 12-14 mm.

## ACCESSORIES

- **PRO-300 HDS** Camera and Monitoring Screen
- **5430-XX** PicoNozzle Kit
- **NanoFil** Microliter Syringes - 10 or 100 µL
- **NFINHLD-G10** 1.0 mm Glass Pipette Holder for NanoFil Syringe
- **MicroFil** for backfilling glass needles
- Pipettors
- **E2XX** Micropipette Storage Jar
- **801566/801963** Vacuum Pump for use with the **PV830**
- **FluoroDish** Optical Glass Bottom Dishes
- **M10** or **M-3** Manipulator Base
- **Z-MOLDS** Microinjection and Transplantation Molds
- Many Surgical Instruments

# Pinpoint Injector

The next generation injection method for increased knock-in efficiency



## Features

- Touch-screen display—resistive touch panel for use with gloves
- Injection control with foot switch or manually with touch screen
- Intuitive user interface
- User adjustable frequency and voltage with touch screen or knobs
- Small footprint
- Four user-programmable protocols
- Adjustable audio continuity tone indicates active probe
- Injection counter indicates total number of injections
- Electrodes include reference electrode option
- Electrode holder compatible with most common manipulators

## Benefits

- High yield (80-85%) of viable, 2-cell stage mouse embryos after injection
- Integrates with WPI's PV820/PV830 or the Eppendorf Femtojet® (Femtojet® 4i is a registered Trademark of Eppendorf AG)
- Hands-free operation with foot switch control
- Works with most inverted microscopes

## MICRO-ePore successes in university labs




Photo courtesy of Dr Janet Rossant, The Hospital for Sick Children in Toronto, Canada. An example of the MICRO-ePORE™ electrode setup on an inverted microscope. MICRO-ePORE™ electrode design ensures compatibility with the leading microscope manufacturers.

## Applications

- Microinjection into oocytes and pre-implantation stage mammalian embryos, including microinjection of CRISPR-Cas9 reagents into the cytoplasm of two-cell stage embryos
- Pronuclear rodent zygote microinjection
- Gene silencing in zebrafish



Preliminary Design

98991 - Auxiliary Knob Controller sold separately

## References

Gu, B., Posfai, E., & Rossant, J. (2018). Efficient generation of targeted large insertions by microinjection into two-cell-stage mouse embryos. *Nature Biotechnology*, 36(7), 632–637. <http://doi.org/10.1038/nbt.4166>

## MICRO-ePORE SPECIFICATIONS

VOLTAGE PARAMETERS	0–3.0 V, at 1 mV increments
FREQUENCY PARAMETERS	50–3000 Hz, at 1 Hz increments
PIPETTE RESISTANCE MAXIMUM	40 MΩ
DIMENSIONS	19.7 × 12.7 × 7.6 cm (7.75 × 5 × 3 in.)
WEIGHT	0.9 kg (2 lb)
CERTIFICATIONS	CE, RoHS

## ORDERING INFORMATION

### MICRO-ePORE MICRO-ePore System™

Includes Electrode Holder • Microinjection Headstage, including adapter for Femtojet™ • Holding Pipette • Platinum Reference Electrode • Vacuum Tubing Set • Foot Switch (Auxiliary Knob Controller sold separately)

<b>98991</b>	Auxiliary Knob Controller
<b>SYS-PV820</b>	Pneumatic PicoPump Microinjector
<b>SYS-PV830</b>	Pneumatic PicoPump Microinjector with Vacuum Pressure
<b>1B100-3*</b>	Single-Barrel Standard Borosilicate Glass Tubing, 3 in. (76mm) long, 1.0 mm OD, 0.58 mm ID, package of 500
<b>TIP1TW1*</b>	Pre-Pulled Glass Pipettes, 1 μm tip ID, 1.0 mm glass OD, Thin wall, fire polished

\*Many varieties and choices are available for glass capillaries and pre-pulled μTip micropipettes.

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • [wpiuk@wpi-europe.com](mailto:wpiuk@wpi-europe.com) • [www.wpi-europe.com](http://www.wpi-europe.com)  
Brazil: 011 55 13 40629703 • [info@brazil.wpiinc.com](mailto:info@brazil.wpiinc.com) • [www.wpiinc.com](http://www.wpiinc.com)

Germany: +49 (0)6031 1602171 • [wptide@wpi-europe.com](mailto:wptide@wpi-europe.com) • [www.wpi-europe.com](http://www.wpi-europe.com)  
China: +86 21 6888 5517 • [chinasales@china.wpiinc.com](mailto:chinasales@china.wpiinc.com) • [www.wpiinc.net](http://www.wpiinc.net)

# Reliable Pneumatic PicoPump

## *Pneumatic microinjector with vacuum pressure*



PV830

Microinjection of sub-nanoliter volumes through a very small opening in a glass micropipette is often accomplished best with a controlled burst of gas pressure. Designed to simplify intracellular injection and a variety of other microinjection tasks, WPI's PicoPumps use precisely regulated pressures for securing cells and injecting them with fluid. Injected volumes range from picoliters to nanoliters.

The **PV830** model provides separate regulation and pressure gauges for the eject and hold pressures. In addition to the features of the **PV820**, the **PV830** model also includes a vacuum pressure regulator and vacuum gauge.

See page 44 for more information.

# Microelectrode Puller

## *A compact, versatile and reliable workhorse*



CE

PUL-1000


**PUL-1000** is a microprocessor controlled horizontal puller for making glass micropipettes or microelectrodes used in intracellular recording, patch clamp studies, microperfusion or microinjection. The puller was designed with tight mechanical specifications and precision electronics for complete control of the pulling process and accurate reproducibility. It offers programmable sequences of up to four steps with heating, force, movement and cooling time. This allows graduated cycles for applications like patch clamp recording.

This puller is a reasonably priced, compact, versatile and reliable workhorse. The microprocessor, combined with the LCD display, makes the **PUL-1000** easy to use.

See page 108 for more information.

# Glass Capillaries

## *Quality glass, superior prices for microinjection/microelectrodes*



Standard, no Filament

Standard Filament

Thin Wall, no Filament

Thin Wall Filament

Our quality borosilicate glass capillaries are available with a large variety of options, including fire polished, filaments, thin wall, specialty glass and multi-barrel.

Borosilicate glass capillaries: Close dimensional tolerances assure microelectrode uniformity and reproducibility. Capillaries are available in 1, 2, 3, 5 and 7-barrel configurations, complete range of single barrel thin-wall sizes and a variety of special configurations. Capillaries with filaments contain a solid filament fused to the inner wall, which speeds filling of electrodes. Capillaries with or without inner filaments are available for making microelectrodes in a wide range of diameters.

For details, see page 47.

# Reliable Pneumatic PicoPump

## Pneumatic microinjector with vacuum pressure



PV830



### Features

- Inject into a single cell with picoliter volumes
- Regulated hold, ejection and vacuum pressure
- Carefully regulated air pressures for securing cells and injecting them with fluid
- Optional independent vacuum regulation
- Pressure Input: 0-120 psi
- Pressure Output: 0.3-90 psi
- Choice of two models

### Benefits

- Hold pressure prevents backfilling of the pipette by capillary action
- PV830 includes vacuum pressure for filling pipettes from the tip or securing floating cells during injection

### Applications

- Intracellular injection in the picoliter to nanoliter range

Although syringe and piston techniques are popular for microinjection in the nanoliter range, microinjection of sub-nanoliter volumes through a very small opening in a glass micropipette is often better accomplished with a controlled burst of gas pressure. Designed to simplify intracellular injection and a variety of other microinjection tasks, WPI's PicoPumps use precisely regulated pressures for securing cells and injecting them with fluid. Injected volumes range from picoliters to nanoliters. Separate ports supply positive and negative pressure—positive pressure for high-pressure ejection, and suction for supporting the cell or for filling the pipette from the tip. In the time between injections, a secondary default "hold" pressure is applied to the pipette to prevent fluid uptake through capillary action or diffusion. The pressure pulses are typically controlled by a precision timer. Timing, ejection pressure, holding pressure and suction are adjusted independently by control knobs and monitored with indicator gauges on the front panel. Injection pressure is controlled by a 20-turn regulator on the front panel. The built-in timing function controls the amount of time that the injection pressure is applied with millisecond resolution. Time intervals can range from 10 s down to 10 ms or less. The pressure burst can also be controlled manually or triggered by an external source.

The PV830 model provides separate regulation and pressure gauges for the eject and hold pressures. In addition to the features of the PV820, the PV830 model also includes a vacuum pressure regulator and vacuum gauge.

#### Hold Pressure Prevents Backfilling

Like the PV830, the PV820 was designed to simplify intracellular injection. The hold pressure prevents backfilling of the pipette by capillary action. Simple to use. Reliable every time. Each PicoPump is



PicoNozzle Kit **5430-ALL** (included) allows micropipettes to be securely mounted in micro-positioners for stable axial air delivery. (Handle diameter 7 mm × 100 mm long.)



The optional **5430-XX** PicoNozzle Kit has a Luer-fitted pipette holder for quick and easy changes. Choose your kit based on the pipette OD (XX indicates OD: 1.0, 1.2, 1.5 or 2.0 mm). The red handle fits easily into most manipulators. This kit contains:

- (1) MPH6S microelectrode holder
- (1) Handle for the MPH6S (4" hollow tube with male Luer fitting at both ends—diameter 6.25 mm × 100 mm long)
- (1) 5' tubing (0.060" ID, 0.120" OD, male locking Luer fitting on one end and a female locking Luer fitting at the other end, rated for 200 PSI and 86 durometer shore A)

supplied with a **5430-ALL** kit that includes two PicoNozzles and tubing to connect the holders to the pressure and vacuum ports.

#### Vacuum Pressure for Additional Applications

Eject pressure supplies a high-pressure pulse for injecting fluid. Hold pressure, which is not sufficient to cause fluid ejection, is used to prevent back filling of the pipette by capillary action or diffusion when the solenoid is inactive. Vacuum can be applied through the ejection port, but vacuum in the PV820 must be externally regulated (see Mini Vacuum Pump, page 71). To fill pipettes from the tip or through a secondary port to secure a floating cell during microinjection, vacuum may be switched from regulated vacuum to atmosphere by using the switch on the front panel.

### References

Di Cara F, Sheshachalam A, Braverman NE, Rachubinski RA, Simmonds AJ (2017). Peroxisome-Mediated Metabolism Is Required for Immune Response to Microbial Infection. *Immunity*. 2017 Jul 18;47(1):93-106.e7. <https://www.ncbi.nlm.nih.gov/pubmed/28723556>

Carrington, B., Varshney, G. K., Burgess, S. M., & Sood, R. (2015). CRISPR-STAT: an easy and reliable PCR-based method to evaluate target-specific sgRNA activity. *Nucleic Acids Research*, 43(22), e157. <http://doi.org/10.1093/nar/gkv802>

# Pneumatic picoliter injector

WPI's **PV820** Pneumatic PicoPump has been providing scientists with precision and repeatable microinjection in volumes ranging from picoliters to nanoliters for decades. The trusted Pneumatic PicoPump is the preferred Zebrafish microinjection pump. The **PV820** offers the same functionality as the **PV830**, minus the independent vacuum regulator.



## PICOPUMP SPECIFICATIONS

	PV820	PV830
<b>PRESSURE</b>		
PRESSURE INPUT	0 to 120 psi	0 to 120 psi
PRESSURE OUTPUT	0.3 to 90 psi *	0.3 to 90 psi
PRESSURE BURST TIMER (10-turn dial)	10 ms to 10 s in Timed Mode	10 ms to 10 s in Timed Mode
REGULATOR ACCURACY	0.1% (20-turn dial) *	0.1% (20-turn dial) *
REGULATOR REPEATABILITY	0.05 psi *	0.05 psi *
GAUGE ACCURACY	3% at full scale *	3% at full scale *
INPUT CONNECTOR	Quick Connect (¼ in. OD Tubing)	Quick Connect (¼ in. OD Tubing)
OUTPUT CONNECTOR	Barbed (⅙-in. ID Tubing)	Barbed (⅙-in. ID Tubing)
<b>VACUUM</b>		
VACUUM INPUT	0 to 30.0 in. Hg	0 to 30.0 in. Hg
VACUUM OUTPUT	Unregulated	0.2 to 29.9 in. Hg
LOWEST REGULATED VACUUM	Unregulated	3 in. water
REGULATOR ACCURACY	Unregulated	0.1% (20-turn dial)
REGULATOR REPEATABILITY	Unregulated	0.03 in. Hg
GAUGE ACCURACY	None	3% at full scale
INPUT CONNECTOR	Quick Connect (¼ in. OD Tubing)	Quick Connect (¼ in. OD Tubing)
OUTPUT CONNECTOR	Barbed (⅙ in. ID Tubing)	Barbed (⅙ in. ID Tubing)
CONTROL	Manual	Manual
VENT	Atmosphere	Atmosphere
<b>CONNECTIONS INCLUDED</b>		
INPUT KIT	10-ft nylon tubing (0.25-in. OD, 1000 psi), one ½-inch female NPT adapter	
OUTPUT KIT	Two PicoNozzle Kits	
<b>PHYSICAL SPECIFICATIONS</b>		
POWER	95-135 V or 220-240 V, 50/60 Hz	95-135 V or 220-240 V, 50/60 Hz
DIMENSIONS	17 x 3.5 x 9.5 in. (43 x 9 x 24 cm)	17 x 5.25 x 9.5 in. (43 x 13 x 24 cm)
SHIPPING WEIGHT	11 lb (5 kg)	14 lb (6.3 kg)

\*Both Hold and Eject Pressure

- Ramezani, T., Laux, D. W., Bravo, I. R., Tada, M., & Feng, Y. (2015). Live Imaging of Innate Immune and Preneoplastic Cell Interactions Using an Inducible Gal4/UAS Expression System in Larval Zebrafish Skin. *Journal of Visualized Experiments*, (96), e52107–e52107. <http://doi.org/10.3791/52107>
- Chokshi, P. (2015). Analysis of scn5Laa and scn5Lab Gene Function in Danio rerio (Zebrafish) Heart Development through TALENs/CRISPR-CAS9-mediated Gene Knockout. Retrieved March 14, 2016, from [http://thescholarship.ecu.edu/bitstream/handle/10342/4923/Chokshi\\_ecu\\_06000\\_11444.pdf?sequence=1](http://thescholarship.ecu.edu/bitstream/handle/10342/4923/Chokshi_ecu_06000_11444.pdf?sequence=1)
- Ewart, M.-A., Kennedy, S., MacMillan, D., Raja, A. L., Watt, I. M., & Currie, S. (2014). Altered vascular smooth muscle function in the ApoE knockout mouse during the progression of atherosclerosis. <http://doi.org/10.1016/j.atherosclerosis.2014.02.014>
- Henson, H. E., Parupalli, C., Ju, B., & Taylor, M. R. (2014). Functional and genetic analysis of choroid plexus development in zebrafish. *Frontiers in Neuroscience*, 8, 364. <http://doi.org/10.3389/fnins.2014.00364>
- Huang, L., Xiao, A., Wecker, A., McBride, D. A., Choi, S. Y., Zhou, W., & Lipschutz, J. H. (2014). A Possible Zebrafish Model of Polycystic Kidney Disease: Knockdown of *wnt5a* Causes Cysts in Zebrafish Kidneys. *Journal of Visualized Experiments*, (94), e52156–e52156. <http://doi.org/10.3791/52156>
- Konantz, J., & Antos, C. L. (2014). Reverse Genetic Morpholino Approach Using Cardiac Ventricular Injection to Transfect Multiple Difficult-to-target Tissues in the Zebrafish Larva. *Journal of Visualized Experiments*, (88), e51595–e51595. <http://doi.org/10.3791/51595>

## ORDERING INFORMATION


- SYS-PV820** PicoPump w/ hold pressure  
**SYS-PV830** PicoPump w/ hold pressure and vacuum  
 Specify line voltage All PicoPumps require external vacuum source.

## OPTIONAL ACCESSORIES/REPLACEMENT PARTS

- 3260** Foot Switch  
**2932** Rack Mount Kit, 3.5-in. high (PV820)  
**2933** Rack Mount Kit, 5.25-in. high (PV830)  
**5430-10** PicoNozzle Kit (MPH6S for 1.0 mm pipette & 5' tubing)  
**5430-12** PicoNozzle Kit (MPH6S for 1.2 mm pipette & 5' tubing)  
**5430-15** PicoNozzle Kit (MPH6S for 1.5 mm pipette & 5' tubing)  
**5430-20** PicoNozzle Kit (MPH6S for 2.0 mm pipette & 5' tubing)  
**5430-ALL** PicoNozzle Kit (for 1.0, 1.2, 1.5 and 1.65 mm pipettes & 5-ft. tubing)  
**75122-110** Gaskets for PicoNozzle , 1.0 mm, green, package of 10  
**75122-210** Gaskets for PicoNozzle — 1.2 mm, black, package of 10  
**75122-310** Gaskets for PicoNozzle — 1.5 mm, blue, package of 10  
**75122-410** Gaskets for PicoNozzle — 1.65 mm, red, package of 10  
**MPH6S** Micropipette Holder (specify 1.0, 1.2, 1.5 or 2.0 mm)  
**MPH6R** Micropipette Holder (specify 1.0, 1.2, 1.5 or 2.0 mm)  
**3316** Replacement Input Kit

# Manual Microsyringe Pumps

*Precise manual injection of fluid using glass pipettes*



## Features

- Accepts syringe sizes: 10  $\mu$ L–1 mL gas tight Luer tip
- Micrometer head has resolution of 10  $\mu$ m per division
- Solid, stainless steel frame

## Benefits

- Choice of models, with or without digital display
- Cost-effective solution

## Applications

- Use for perfusion or withdrawal of liquids

**MMP** and **DMP** are convenient tools for precise manual injection of fluid using glass pipettes or similar injection devices. The design allows visual feedback of flow at the pipette tip. They can also be used as a manual micro syringe pump for perfusion or withdrawal of liquids. The resolution of the injection volume can be continuously varied from 10 nL to the microliter range, depending on the syringe used.

### With or without digital display

DMP comes with a digital micrometer that will allow the reading of piston advancement easily with a 0.001 mm resolution. MMP has the traditional mechanical micrometer head with a resolution of 10  $\mu$ m per division and advances 500  $\mu$ m per revolution.

### Constructed of stainless steel

The entire frame body of the injector is constructed of stainless steel for excellent stability and durability. The piston of the micrometer can be adjusted to the syringe's plunger position. Small diameter PTFE tubing is used for the accuracy and solution compatibility. The unique design of the pipette holder can securely hold any pipette with an outer diameter of between 1.0 mm and 1.5 mm. All necessary accessories for removing air and filling the syringe and tubing with liquid are included. The system comes complete with a 100  $\mu$ L gas tight syringe. Other syringe sizes may be purchased separately.

## DMP & MMP SPECIFICATIONS

TRAVEL DISTANCE	25 mm
ADVANCES RESOLUTION	0.001 mm for DMP and 0.01 mm for MMP
SYRINGE SIZE	10 $\mu$ L to 1 mL gas tight Luer tip syringe
TUBING	1.5 m of PTFE tubing with 0.5 mm ID
PIPETTE HOLDER	0.24" x 5.2"
PIPETTE HOLDER FITS	1.0 to 1.5 mm OD pipette

## ORDERING INFORMATION

<b>MMP</b>	Manual Microsyringe Pump
<b>DMP</b>	Manual Microsyringe Pump with Digital Display

## OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>MMP-KIT</b>	Injection Assembly Parts Kit Not including valve—see <b>14057-10</b> , page 220)
----------------	---

# Microinjection and Transplantation Molds

*A simpler method for handling zebrafish embryos*

## Features

- For high throughput microinjection research
- Make impressions in agarose gel to facilitate embryo alignment
- Four molds per kit
- Reusable

## Benefits

- Organizes and immobilizes embryos for microinjection

## Applications


- Zebrafish research

Mold your agarose, pipette in your embryos, and watch them auto aligning in the grooves. **Z-MOLDS** Microinjection and Transplantation Molds (4 per kit) are designed for zebrafish research.

**Proteomics & Large Screening**—Inject many embryos, up to 1000. The grooves made by the mold in the agarose gel enable the embryos to self align.



**Xenograft and Larval Injection**—The sloped ridges make perfect angles in the agarose gel, which then makes it easier to do microinjections in the larvae.



**Transplantation**—Increase the speed of doing microinjections. Simply turn the petri dish as you are injecting.



**Standard Microinjection**—Use for blastomere transplantation.



## ORDERING INFORMATION

**Z-MOLDS** Microinjection & Transplantation Molds

## References

Kitambi SS, Toledo EM, Usoskin D, Wee S, Harisankar A, Svensson R, Sigmundsson, K, Kalderén C, Niklasson M, Kundu S, Aranda S, Westermark B, Uhrbom L, Andäng M, Damberg P, Nelander S, Arenas E, Artursson P, Walfridsson J, Forsberg Nilsson K, Hammarström LG, Ernfors P. Vulnerability of glioblastoma cells to catastrophic vacuolization and death induced by a small molecule. *Cell*. 2014 Apr 10;157(2):313-28.



# Glass Capillaries

Quality glass, superior prices for microinjection/microelectrodes

## Features

- Quality borosilicate glass capillaries
- Large variety available, including fire polished, filaments, thin wall, specialty glass and multi-barrel

## Benefits


- Superior pricing
- Most glass orders ship within 48 hours

## Applications

- Microinjection
- Electrophysiology
- Patch clamp
- Fluid Handling

## Fire Polishing

Fire-polished glass capillaries are easier to insert into microelectrode holders without damaging the gasket. More importantly, fire-polished glass won't scratch the chloridized wire used in a recording electrode. Fire-polishing does not affect the glass's mechanical or electrical properties.



## Making Uniform, Reproducible Microelectrodes



Borosilicate glass capillaries: Close dimensional tolerances assure microelectrode uniformity and reproducibility. Capillaries are available in 1, 2, 3, 5 and 7-barrel configurations, complete range of single barrel thin-wall sizes and a variety of special configurations. Capillaries with filaments contain a solid filament fused to the inner wall, which speeds filling of electrodes. Capillaries with or without inner filaments are available for making microelectrodes in a wide range of diameters.

## Filament Glass Capillaries


Single barrel standard wall thickness capillaries are offered either with or without inner filaments for quick filling in a variety of lengths and diameters.

## Thin Wall Glass Capillaries

Thin wall single barrel capillaries are offered both with or without inner filaments.



Specialty glass is also available. See page 217.



## ORDERING INFORMATION

	Length	OD(mm)	ID(mm)	Filament	Fire-Polished	Quantity	Item
Single Barrel Standard Borosilicate Glass	3 in. (76 mm)	1.0	0.58	✓		500	<b>1B100F-3</b>
	3 in. (76 mm)	1.0	0.58			500	<b>1B100-3</b>
	3 in. (76 mm)	1.2	0.68	✓		350	<b>1B120F-3</b>
	3 in. (76 mm)	1.2	0.68			350	<b>1B120-3</b>
	3 in. (76 mm)	1.5	0.84	✓		225	<b>1B150F-3</b>
	3 in. (76 mm)	1.5	0.84		✓	300	<b>1B150-3</b>
	4 in. (100 mm)	1.0	0.58	✓	✓	500	<b>1B100F-4</b>
	4 in. (100 mm)	1.0	0.58		✓	500	<b>1B100-4</b>
	4 in. (100 mm)	1.2	0.68	✓	✓	400	<b>1B120F-4</b>
	4 in. (100 mm)	1.2	0.68		✓	350	<b>1B120-4</b>
	4 in. (100 mm)	1.5	0.84	✓	✓	300	<b>1B150F-4</b>
	4 in. (100 mm)	1.5	0.84		✓	300	<b>1B150-4</b>
	4 in. (100 mm)	2.0	1.12	✓		125	<b>1B200F-4</b>
	4 in. (100 mm)	2.0	1.12		✓	200	<b>1B200-4</b>
	6 in. (152 mm)	1.0	0.58	✓		500	<b>1B100F-6</b>
	6 in. (152 mm)	1.0	0.58			500	<b>1B100-6</b>
6 in. (152 mm)	1.2	0.68	✓		350	<b>1B120F-6</b>	
6 in. (152 mm)	1.2	0.68			350	<b>1B120-6</b>	
6 in. (152 mm)	1.5	0.84	✓		225	<b>1B150F-6</b>	
6 in. (152 mm)	1.5	0.84			225	<b>1B150-6</b>	
6 in. (152 mm)	2.0	1.12	✓		125	<b>1B200F-6</b>	
6 in. (152 mm)	2.0	1.12			125	<b>1B200-6</b>	
Thin-Wall Single-Barrel Standard	3 in. (76 mm)	1.0	0.75	✓		500	<b>TW100F-3</b>
	3 in. (76 mm)	1.0	0.75			500	<b>TW100-3</b>
	3 in. (76 mm)	1.2	0.90	✓	✓	400	<b>TW120F-3</b>
	3 in. (76 mm)	1.2	0.90			350	<b>TW120-3</b>
	3 in. (76 mm)	1.5	1.12	✓		225	<b>TW150F-3</b>
	3 in. (76 mm)	1.5	1.12		✓	300	<b>TW150-3</b>
	4 in. (100 mm)	1.0	0.75	✓		500	<b>TW100F-4</b>
	4 in. (100 mm)	1.0	0.75		✓	500	<b>TW100-4</b>
	4 in. (100 mm)	1.2	0.90	✓		350	<b>TW120F-4</b>
	4 in. (100 mm)	1.2	0.90			350	<b>TW120-4</b>
	4 in. (100 mm)	1.5	1.12	✓		225	<b>TW150F-4</b>
	4 in. (100 mm)	1.5	1.12		✓	300	<b>TW150-4</b>
	6 in. (152 mm)	1.0	0.75	✓		500	<b>TW100F-6</b>
	6 in. (152 mm)	1.0	0.75		✓	500	<b>TW100-6</b>
	6 in. (152 mm)	1.2	0.90	✓	✓	400	<b>TW120F-6</b>
	6 in. (152 mm)	1.2	0.90			350	<b>TW120-6</b>
6 in. (152 mm)	1.5	1.12	✓		225	<b>TW150F-6</b>	
6 in. (152 mm)	1.5	1.12		✓	300	<b>TW150-6</b>	

Single barrel glass is Kimble N51A. All thin wall glass is Schott Duran 8330.

# Microelectrode Puller

## Cost-effective and compact micropipette puller

### Features

- Microprocessor controlled
- Program sequences up to four steps
- Store up to 95 programs in memory for easy recall
- Multiple factory programs installed

### Benefits

- Tempered glass cover to reduce the effects of humidity on puller reproducibility
- Switchable power supply ensures that line voltage fluctuations don't affect reproducibility

### Applications

- Pull your own microelectrodes and micropipettes



PUL-1000

CE

**PUL-1000** is a microprocessor controlled horizontal puller for making glass micropipettes or microelectrodes used in microperfusion or microinjection. The puller was designed with tight mechanical specifications and precision electronics for complete control of the pulling process and accurate reproducibility. It offers programmable sequences of up to four steps with heat, pulling force, pull length and cooling time.

This puller is a reasonably priced, compact, versatile and a reliable workhorse. The microprocessor, combined with the LCD display, makes the **PUL-1000** easy to use.

#### Tempered Glass Cover

The cover of the pulling chamber is made with tempered glass to minimize the temperature effect on the reproducibility of pulled pipettes.

#### Switchable Power Supply

**PUL-1000** has a high quality switching power supply for use anywhere in the world without worry about the line voltage differences. Pulling reproducibility is unaffected by line voltage fluctuation. Heating voltage can be controlled to within 0.1% accuracy even when line voltage fluctuates from 90 to 240 VAC

### References

Plautz, C. Z., Williams, H. C., & Grainger, R. M. (2016). Functional Cloning Using a Xenopus Oocyte Expression System. *Journal of Visualized Experiments* (107), e53518–e53518. <http://doi.org/10.3791/53518>

Komarova, Y., Peloquin, J., & Borisy, G. (2011). Components of a microinjection system. *Cold Spring Harbor Protocols*, 2011(8), 935–9. <http://doi.org/10.1101/pdb.ip27>

### PUL-1000 SPECIFICATIONS

HEATER ELEMENT	Platinum/Iridium
PULLING FORCE	Solenoid, adjustable
TAPER LENGTH	1–10 mm
CAPILLARY OD RANGE	1.0–2.0 mm
MAXIMUM CAPILLARY LENGTH	170 mm
MINIMUM CAPILLARY LENGTH	55 mm
MEMORY SETS	95
AUTO SHUT-OFF TIME	90 s
POWER	90–240VAC, 50/60 Hz
DIMENSIONS	34 x 24 x 12 cm
SHIPPING WEIGHT	16 lb

### ORDERING INFORMATION

#### PUL-1000 Micropipette Puller

System includes PUL-1000 Puller; TW100-4 Thin Wall Glass Capillaries (package of 500); MF34G-5 package of MicroFil

<b>13834</b>	Filament, Platinum/Iridium, 2.5mm wide, 2.5mm Square Box
<b>504951</b>	Filament, Platinum/Iridium, 4.5mm wide, Trough

## Also see these other pullers available from WPI

**SU-P1000**  
See page 110.



**SU-P97**  
See page 111.



### WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • [wpiuk@wpi-europe.com](mailto:wpiuk@wpi-europe.com) • [www.wpi-europe.com](http://www.wpi-europe.com) Germany: +49 (0)6031 1602171 • [wptide@wpi-europe.com](mailto:wptide@wpi-europe.com) • [www.wpi-europe.com](http://www.wpi-europe.com)  
Brazil: 011 55 13 40629703 • [info@brazil.wpiinc.com](mailto:info@brazil.wpiinc.com) • [www.wpiinc.com](http://www.wpiinc.com) China: +86 21 6888 5517 • [chinasales@china.wpiinc.com](mailto:chinasales@china.wpiinc.com) • [www.wpiinc.net](http://www.wpiinc.net)

# Fluid Handling



## A large variety of pumps, from picoliter to microliter

We offer a large variety of pumps for everything from microfluidics to general fluid handling applications. Our peristaltic pumps are easy to setup and clean, offer continuous flow with virtually “infinite” volume (depending only on the capacity of your source), require no contact with metal or the pump and are good for large volume pumping. Our syringe pumps provide accurate volume control and are an excellent choice for lower volume applications. Our pneumatic PicoPumps are designed for delivery of very small volumes.

# Which Pump is Right for Me?


## PUMPS & FLUID HANDLING COMPARISON CHART

	Fluid Range	Channels	Special features	Page
<b>PERISTALTIC PUMPS</b>				
<b>MiniStar™</b>	0.006- 37 mL/min	1	Compact design, remote control	53
<b>Peripro-2HS</b>	0.8 - 300 mL/min	2	Calibrated output, replaceable tubing cartridges	52
<b>Peripro-4HS</b>	0.8 - 300 mL/min	4	Calibrated output, replaceable tubing cartridges	52
<b>Peripro-4LS</b>	0.01-80 mL/min	4	Calibrated output, replaceable tubing cartridges	52
<b>Peripro-8LS</b>	0.01-80 mL/min	8	Calibrated output, replaceable tubing cartridges	52
<b>SINGLE-CHANNEL AND WHOLE-CELL SOLUTION EXCHANGE</b>				
<b>MPS2</b>	Up to 250 µL/min.	8	Programmable control; low dead volume	51
<b>LABORATORY SYRINGE PUMPS</b>				
<b>AL-1000</b>	0.73 µL/hr to 1699 mL/hr	1	Push/pull	65
<b>AL-2000</b>	0.73 µL/hr to 1699 mL/hr	2	Push/pull (2 networked pumps)	65
<b>AL-1000 HP</b>	1.459 µL/hr - 6120 mL/hr	1	Infuse/withdraw, High Pressure	65
<b>AL-4000</b>	1.459 µL/hr to 6120 µL/hr	2	Infuse/withdraw	65
<b>AL-6000</b>	0.73 µL/hr to 1699 mL/hr	6	Infuse/withdraw	65
<b>AL-8000</b>	0.454 µL/hr to 163 mL/hr	8	Infuse/withdraw	65
<b>SPLG100</b>	1.26 pL/min to 88.32 mL/min	1	Infuse only	72
<b>SPLG101</b>	1.26 pL/min to 25.99 mL/min	2	Infuse only	72
<b>SPLG110</b>	1.26 pL/min to 88.28 mL/min	1	Infuse/withdraw	73
<b>SPLG200</b>	0.5 pL/min to 220.97 mL/min	2	Infuse only	74
<b>SPLG210</b>	0.5 pL/min to 220.97 mL/min	2	Infuse/withdraw	73
<b>SPLG212</b>	0.5 pL/min to 220.97 mL/min	2	Infuse/withdraw programmable	73
<b>SPLG270</b>	0.5 pL/min to 220.97 mL/min	2+2	Push-pull	73
<b>SPLG272</b>	0.5 pL/min to 220.97 mL/min	2+2	Push-pull programmable	73
<b>SP100i</b>	0.0001-519 mL/hr	1	Basic single channel	-
<b>SP101i</b>	0.001 µL/hr - 35 mL/min	2	Micro dialysis application	-
<b>SP120p</b>	0.1 µL/hr - 127 mL/hr	1+1	Push pull, single cycle	-
<b>SP200i</b>	0.001 µL/hr - 145 mL/min	2	RS232 TTL/foot switch	-
<b>SP210c</b>	0.001 µL/hr - 86 mL/min	2+2	RS232 push pull, continuous	-
<b>SP210iw</b>	0.001 µL/hr - 145 mL/min	2	RS232 infuse/withdraw	-
<b>SP220i</b>	0.001 µL/hr - 21 mL/min	10	RS232 infuse only	-
<b>SP230iw</b>	0.001 µL/hr - 21 mL/min	10	RS232 infuse/withdraw	-
<b>SP250i</b>	0.001 µL/hr - 21 mL/min	4	RS232 infuse only	-
<b>SP260p</b>	0.001 µL/hr - 86 mL/min	2+2	RS232 push pull, single cycle	-
<b>MICRO SYRINGE PUMP / STEREOTAXIC INJECTION</b>				
<b>UMP3T-1</b>	0.03 nL/min - 10 µL/sec	1	Ultra micro infuse/withdraw RS232	56
<b>MMP</b>	Manual 100 µL-1 mL syringe	1	Manual	64
<b>DMP</b>	Manual 100 µL-1 mL syringe	1	Digital readout micrometer	64
<b>MICROINJECTION</b>				
<b>PV820</b>	Injected volumes from pL to nL	1	Injection pressure and holding pressure	63
<b>PV830</b>	Injected volumes from pL to nL	1	Injection pressure and holding pressure and vacuum	62
<b>NANOLITER2010</b>	Bolus, 2.3-69 mL/Injection	1	Oocyte injector, infuse only	61
<i>Microinjection Systems (Zebrafish, C. Elegans, Drosophila, Xenopus oocytes)</i>				
<b>MICROFLUIDICS</b>				
<b>ExiGo</b>	50 nL/min - 10 mL/min	1	Infuse only, feedback via integrated flow sensor, includes iPad mini which can control up to four pumps	77
<b>Mirus</b>	100nL/min - 10 mL/min ±1%	8	Microchip perfusion. Infuse only, reversible flow, ~600 µL dead volume, PC control	78
<b>Kima</b>	15 - 35 mL/hr ±4%	1	Microchip perfusion. Infuse only, recirculating pump controlled by iPod Touch, wi-fi communication, <300 µL dead volume	77
<b>Unigo</b>	1 µL/min-1 mL/min; unidirectional (push)	1	A precision, microfluidic, single-channel, pressure pump. Can expand up to 4 channels.	78
<b>SUPPLIES &amp; TOOLS</b>				
Pipetters — standard and autoclavable				81
Miniature Vacuum Pump				82

Available online: [wpiinc.com](http://wpiinc.com)

# Multichannel Perfusion System

For single ion channel and whole-cell solution exchange



## Features

- Fast solenoid valve (response time 2 ms)
- Manual or Programmable PC control with user-friendly GUI interface (software included)
- 8-Channels-10 mL syringes
- Gravity feed
- Quick cleaning and priming valves included

## Benefits

- Color-coded polyurethane ribbon style tubing for easy identification
- Super low dead volume (<100 nL) micromanifold
- Economically priced
- Organized, color-coded tubing system without contaminants of PVC tubing
- Low flow resistance and dead volume with the unique micromanifold

## Applications

- Single channel perfusion applications
- Whole cell patch perfusion applications
- Drug delivery under a microscope (in patch clamp application)

Micromanipulator, tilt base and petri dish not included.

**MPS-2** is a programmable 8-channel perfusion system designed for single channel and whole-cell patch preparations. Offering the best combination of performance and value, the **MPS-2** incorporates the same high quality solenoid valves found on similar but much more expensive systems. Unlike other perfusion systems on the market, which often compromise performance to fit every possible application, the **MPS-2** is the only perfusion system designed and optimized specifically for single-channel and whole-cell patch perfusion applications.

### Options for Manual and Automatic Control

The system can be controlled manually via membrane switches on the front panel or through a PC. Two different manual control modes are offered. One controls each channel independently and the other mode allows you to assign a master channel that will keep the system flow when all other channels are switched off. User-friendly graphic timing software is included, and the programmed perfusion sequence can be started by computer, a patch clamp amplifier or other external trigger, or manually.

### Organized, Color-coded, Ribbon Tubing

The perfusion fluid flows through specially designed color-coded polyurethane ribbon style tubing. The color-coding allows you to easily trace each channel for diagnostic checks or set up, and the ribbon style of tubing keeps the system neat and organized. Unlike PVC based tubing, polyurethane tubing contains no plasticizer, which can cause contamination. The tubing ribbon is designed as an economical disposable item, which is often critical when cleanliness is needed.

### Low Flow Resistance and Dead Volume

The most unique feature of the **MPS-2** is its perfusion micromanifold. Using the latest microfluidic techniques, the injection molded micromanifold provides the least flow resistance and dead volume of any product on the market. The flow channel inner diameter is approximately 1.0 mm, except for the last 5 mm before the junction point. This design allows a fast flow rate without using a pressured system. Small channels and a unique design at the merging point further reduce the chance of cross contamination. Dead volume is less than 100 nL at junction.

### MPS-2 SPECIFICATIONS

CHANNELS	8
VALVE RESPONSE TIME	2 ms
VALVE CONTROL	USB, TTL, external start via software
SYRINGE RESERVOIR VOLUME	10 mL
<b>MANIFOLD</b>	8 to 1
OUTPUT TUBING ID	250 $\mu$ m and 100 $\mu$ m.
	100 $\mu$ m ID tip, 8 $\mu$ L/min.
	Tubing cap: 7.85 nL/mm
MAXIMUM FLOW RATES (with 50 mm output tubing, gravity fed)	250 $\mu$ m ID tip, 250-500 $\mu$ L/min.
	Tubing cap: 49.1 nL/mm
DEAD VOLUME	< 100 nL excluding the single outlet tubing

### ORDERING INFORMATION


**MPS-2** Multichannel Perfusion System & Control Software

System includes: stand, PS-2 controller, valve console, syringe holder, power cord, USB cable, DB9-to-BNC 8-cable assembly, (2) 1A fuses, (10) 10 mL syringes, (10) 3-way stopcocks, (10) Luer fitting with barb for 1/16-in. ID tubing, color coded polyurethane tubing ribbon (5 ft), tubing for making micromanifold cleaning adaptor (3 in.), micromanifold holding rod, (2) micromanifolds with 100  $\mu$ m id tip, (2) micromanifolds with 250  $\mu$ m id tip and installation software

### OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>502109-15</b>	Color-coded PU Tubing, 1/16" ID x 8 Channels, 15 ft
<b>502110</b>	Micromanifold, 100 $\mu$ m ID tip, 2 pcs/pk
<b>502125</b>	Micromanifold, 250 $\mu$ m ID tip, 2 pcs/pk

Specify line voltage and Micromanifold tip OD when ordering.



Micromanifold closeup: Fluid-filled passages are shown in magenta.

# Affordable High Performance Peristaltic Pump

*Digital Peri-Star™ Pro peristaltic pump is affordably priced!*

## Features

- Available in 2-, 4- and 8-channel versions
- Display either rotation speed (RPM) or flow rate (mL/min)
- Wide flow range: 0.01 - 280 mL/min
- Accuracy of flow rate: 0.5% using self calibration function
- Accuracy of speed: 0.1 rpm
- Large backlit digital LCD display
- Programmable for all tubing sizes between 0.8 mm and 6.4 mm ID
- Easy and fast tubing replacement using snap-on cartridges
- Membrane keypad allows easy programming while protecting controls from fluid entry
- Actively driven rollers by planetary gears for long lasting tubing life

## Benefits

- Backlit display with water resistant cover
- Easy setup and calibration with simple instructions on screen
- Both high and low flow pumps employ planetary gears for minimal pulsations and greater accuracy

## Applications

- Chemical transfer
- Pharmaceutical processing

**Peri-Star™ Pro** peristaltic pumps provide accurate and precise pumping with convenience and versatility. Peri-Star Pro can be run in either flow rate mode (mL/min) or rotation speed mode (rpm). For good laboratory practice, pumps must be calibrated after changing the tubing and solution. You can easily calibrate Peri-Star Pro to deliver flow as accurate as 0.5% in a wide flow range from 0.01 mL/min to 280 mL/min. Under rotation speed mode, the digitally controlled stepping motor provides accurate and reproducible operation with 0.1% rpm both forward and in reverse.

### Backlit display with water resistant cover

Large backlit digital LCD display provides readouts of rotation direction, flow rate or rotation speed, tubing ID, drive status and remote control mode simultaneously. Water resistant membrane keypad allows easy programming while protecting LCD display and controls from fluid entry.

### Easy set up and calibrate the pump

Built-in Human Machine Interface (HMI) with screen instructions in plain English steps you through initial setup, calibration and operating procedures. The user-friendly interface reduces the need to frequently check the printed manual for instruction and reference.

### Unique planetary gears for accuracy

Peri-Star Pro is available in two versions: a 4-roller version for high flow and an 8-roller version for lower volumes which provides high pressure with minimal pulsations.

A unique planetary gear design with eight actively driven rollers (four rollers for higher flow rate model), together with independent tubing compression fine adjustment, enables high flow accuracy and prolongs tubing life. Snap-on cartridges allow tubing to be changed quickly without cross contamination of solutions.



CE

PERIPRO-4LS

## PERI-STAR PRO SPECIFICATIONS

	Peri-Star Pro 2H / 4H (High Rate)	Peri-Star Pro 4L / 8L (Low Rate)
NUMBER OF ROLLERS	4	8
NUMBER OF CHANNELS	2-4	4-8
ROTOR SPEED RANGE	1-100 rpm	1-100 rpm
FLUID FLOW RANGE	0.8-280 mL/min #17 Tubing: 3.5-280 mL/min	0.01-80 mL/min #14 Tubing: 0.2-18 mL/min
TUBING RANGE	3.1-6.4 mm ID	0.5-2.4 mm ID
SELF-CALIBRATION	Yes	Yes
WORKING ENVIRONMENT	0-45°C, Humidity < 80%	0-45°C, Humidity < 80%
POWER	110 V or 220 V AC, 50 - 60 Hz	110 V or 220 V AC, 50 - 60 Hz
DIMENSIONS	190 x 162 x 275 mm	190 x 162 x 275 mm
SHIPPING WEIGHT	11 lb. (5 kg)	11 lb. (5 kg)

## ORDERING INFORMATION

### PERI-STAR PRO PUMPS

<b>PERIPRO-2HS</b>	Peri-Star™ Pro, 2-channel, High Rate, Large Tubing (110-220V)
<b>PERIPRO-4HS</b>	Peri-Star™ Pro, 4-channel, High Rate, Large Tubing (110-220V)
<b>PERIPRO-4LS</b>	Peri-Star™ Pro, 4-channel, Low Rate, Small Tubing (110-220V)
<b>PERIPRO-8LS</b>	Peri-Star™ Pro, 8-channel, Low Rate, Small Tubing (110-220V)

### OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>503049</b>	Replacement Tubing Cartridge, Large
<b>503050</b>	Replacement Tubing Cartridge, Small
<b>503022</b>	Replacement Silicone Tubing, 1 m, 1.6 mm I.D., #14, with stops
<b>503023</b>	Replacement Silicone Tubing, 1 m, 6.4 mm I.D., #17
<b>503120</b>	TTL Control Module


## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
Brazil: (013) 406-29703 • info@brazil.wpiinc.com • www.wpiinc.com


Germany: +49 (0)6031 1602171 • wpide@wpi-europe.com • www.wpi-europe.com  
China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net

# Mini Peristaltic Pump

*Lightweight, portable Mini★Star™*




Pump head for  
**PERIPRO-4HS**  
(4 channels, 4 rollers)



Pump head for  
**PERIPRO-8LS**  
(8 channels, 8 rollers)



Pump head for  
**PERIPRO-2HS**  
(2 channels, 4 rollers)



Pump head for  
**PERIPRO-4LS**  
(4 channels, 8 rollers)

## Features

- Flow Range: 0.06-14.0 mL/min
- Speed: 1-50.0 rpm, forward/reverse
- Stand and clamp shown in the image are included
- Includes tubing

## Benefits

- Remote control (wired)
- Small footprint
- Quiet operation
- Compact and easily transported

## Applications

- Small volume chemical transfer
- Pharmaceutical processing
- Sample perfusion



MINISTAR

This compact, lightweight peristaltic pump fits just about anywhere. It can be mounted directly on the bench, in a regular rack or to a post. The speed can be adjusted from 1 to 50 rpm. With recommended silicone tubing, the volume can be set from 0.06 to 14.0 mL/min. The **MiniStar™** also features a hand held remote control that allows you to start and stop the pump, purge or adjust its speed and direction.

## MINISTAR™ SPECIFICATIONS

CHANNEL	1
SPEED	1-50.0 rpm, forward/reverse
FLOW RANGE	0.06~14.0 mL/min
RESOLUTION	1 rpm (0.1 rpm computer control)
SPEED CONTROL	Remote control
DISPLAY	Indicators for status and speed
POWER	12 V DC (110/220 VAC adapter incl.)
WORKING CONDITION	Temperature 0-40°C, humidity < 80%
TUBING (Two-stop Silicone)	
Wall Thickness	0.8~1.0 mm
Outer Diameter	≤ 4.8 mm
DIMENSION OF DRIVER	135×72×72 mm (L×W×H)
DIMENSION OF REMOTE CONTROL	105×50×16 mm (L×W×H)
WEIGHT OF DRIVER	0.5 kg

## ORDERING INFORMATION

<b>MINISTAR</b>	Miniature Peristaltic Pump, 1-channel
<b>504011</b>	MiniStar™ and Stand (as pictured above)
<b>503120</b>	TTL Control Module
<b>503121</b>	Silicone Tubing with stops, 2.4mm ID x 0.8mm wall x 1 m (5-pk)
<b>503122</b>	Silicone Tubing with stops, 1 mm ID x 1 mm wall x 1 m (5-pk)

# Microinjection Syringe Pump

*Delivering picoliter volumes precisely*

## Features

- Graphic display with SMARTouch touch screen controller for "intelligent", easy to use interface controlling up to four syringe pumps
- Splash proof touch screen
- User configurable mounting bar
- Dual mode motor drive
- Compatible with all UMP, UMP2 and UMP3 pumps
- Optional foot switch available
- 5-Digit display

## Benefits

- Accepts a wide variety of microinjection syringes
- Manual or automated injections
- Quiet operation for electrophysiology recordings
- Mounts directly on micromanipulator or stereotaxic frame
- Nominal injections down to 1.0 nL
- Rapid setup with intuitive touchscreen controller

## Applications

- Microinjection
- Neuroscience
- Microfluidics
- Micro delivery of biochemical agents or dyes

The UltraMicroPump 3 (**UMP3**) is a versatile pump which uses micro syringes to deliver picoliter to milliliter volumes. The pump is optimum for applications that require injections of precise and small amounts of liquid. With its touchscreen controller, UMP3 can dispense as little as 0.53  $\mu\text{L}/\text{step}$  (using 10  $\mu\text{L}$  syringe with 60 mm scale length).

The new SMARTouch™ controller for the UltraMicroPump features Patent Pending technology which includes:

- Total system calibration – Calibrate the syringe and the controller together as a system. This feature eliminates the variability of the syringes and delivers the calibrated volume.
- Smart smoothness – The controller can be set to automatically adjust microstepping according to the injection rate to deliver the smoothest flow.
- User defined travel limits – Set the limits for a specific syringe in the software. This prevents the pump from over-driving the plunger into the syringe, potentially causing syringe breakage.

The **MICRO2T** SMARTouch™ controller is feature rich. All operations are controlled through interactive touch screen. It has a graphical indication of the flow and the volume remaining in the syringe. It offers automatic end stop detection that is dependent on the syringe volume. You can control two pumps independently from one controller with its dual display. It also has automatic pump detection and a Pause/Resume feature that allows dosing during infusion/withdrawal. The volume accumulated is displayed on screen, as well as the percentage of volume left in the syringe. The SMARTouch controller is fully compatible with all earlier versions of the UltraMicroPump.

### Low Fluid Dead Volume


Syringes may be filled externally and then inserted into the pump or filled while mounted in the pump. Fluids injected or withdrawn are held entirely within the micro syringe to maintain a low fluid dead volume.

### Flexibility in Mounting

For positioning, the UMP3 may be attached to any of several WPI micro-positioners such as the **M3301** (manual), **SU-QUAD** or **SU-TRIO** (motorized) or any manual stereotaxic manipulator.

### Rapid Setup with Intuitive Touchscreen Controller

An integral component in the UMP3 system is the SMARTouch™ touchscreen controller, which provides an "intelligent" and easy-to-use interface to up to four (or two) syringe pumps. Operating parameters are



**UMP3** shown mounted to the **MTM-3** Motorized Stereotaxic Frame (not included).

set with the touchscreen panel. You can save your parameters for instant recall. An optional foot switch offers "hands free" start/stop operation.

**Computer Control**—A USB port on the rear of the controller can be used to connect it to a computer for scripted protocols.

**NOTE:** UMP3 accepts glass syringes with barrel diameters from 5.5 to 9 mm.

## References

Augustad, I. L., Nyman, A. K. G., Costa, A. I., Barnett, S. C., Sandvig, A., Håberg, A. K., & Sandvig, I. (2017). Effects of Neural Stem Cell and Olfactory Ensheathing Cell Co-transplants on Tissue Remodelling After Transient Focal Cerebral Ischemia in the Adult Rat. *Neurochemical Research*, 1–11. <http://doi.org/10.1007/s11064-016-2098-3>

Eleftheriadou, I., Dieringer, M., & Poh, X. (2017). Selective transduction of astrocytic and neuronal CNS subpopulations by lentiviral vectors pseudotyped with Chikungunya virus envelope. *Biomaterials*. Retrieved from <http://www.sciencedirect.com/science/article/pii/S0142961217300376>

Job, M. O., & Kuhar, M. J. (2017). Cart Peptide in the Nucleus Accumbens Regulates Psychostimulants: Correlations Between Psychostimulant and Cart Peptide Effects. *Neuroscience*. <http://doi.org/10.1016/j.neuroscience.2017.02.012>

Mosberger, A. C., et al. (2017). Axotomized Corticospinal Neurons Increase Supra-Lesional Innervation and Remain Crucial for Skilled Reaching after Bilateral Pyramidotomy. *Cerebral Cortex* (New York, 1991). <http://doi.org/10.1093/cercor/bhw405>



The SMARTouch™ controller can control up to two UMP3 pumps.

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • [wpiuk@wpi-europe.com](mailto:wpiuk@wpi-europe.com) • [www.wpi-europe.com](http://www.wpi-europe.com)  
Brazil: (013) 406-29703 • [info@brazil.wpiinc.com](mailto:info@brazil.wpiinc.com) • [www.wpiinc.com](http://www.wpiinc.com)

Germany: +49 (0)6031 1602171 • [wpide@wpi-europe.com](mailto:wpide@wpi-europe.com) • [www.wpi-europe.com](http://www.wpi-europe.com)  
China: +86 21 6888 5517 • [chinasales@china.wpiinc.com](mailto:chinasales@china.wpiinc.com) • [www.wpiinc.net](http://www.wpiinc.net)





UMP3 shown mounted on a standard micromanipulator (not included) and using the new SMARTouch controller.

## ULTRAMICROPUMP SPECIFICATIONS

(based on 10 µL syringe)

### NORMAL MODE

TRAVEL	62 mm
MINIMUM DISPENSING VOLUME	0.58 nL / step (10 µL syringe)
LINEAR MOTION PER STEP	3.175 µm/half step
WEIGHT	325 g (11.5 oz)
MOUNTING ROD DIAMETERS	7.9 mm (0.31 in.)
MAINS POWER SUPPLY	90-264VAC @ 47-63Hz
DIMENSIONS	∅ 32 mm x 190 mm (∅ 1.3 in. x 7.5 in.)
MAXIMUM SYRINGE SIZE	1 mL

### MICROSTEPPING MODE

Precision is increased eight-fold

## ORDERING INFORMATION

<b>UMP3T-1</b>	UltraMicroPump III (one) and <b>Micro2T</b> Controller
<b>UMP3T-2</b>	UltraMicroPump III (two) and <b>Micro2T</b> Controller
<b>UMP3-3</b>	UltraMicroPump III (three) and <b>Micro4</b> Controller
<b>UMP3-4</b>	UltraMicroPump III (four) and <b>Micro4</b> Controller
<b>UMP3</b>	UltraMicroPump III (without controller)
<b>MICRO2T</b>	SMARTouch Controller, Two-Channel

## OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>13142</b>	Foot Switch for MICRO2T*
<b>502201</b>	V-clamp for Stereotaxic Frame
<b>503301</b>	Extension Cable, miniDIN (male-female) 10 ft
<b>503207</b>	Small Base Stand and Clamps

\*13142 is the foot pedal for use with the MICRO2T ONLY. It is not compatible with the old MICRO4 controller, which uses the 15867 foot pedal. The two foot pedals are not cross-compatible.

## Microvolume Syringes

### ORDERING INFORMATION

#### Syringes with Luer Fitting (no needle)

Order No.	Volume	Description	O.D.	SCALE LENGTH	UMP3	UMP2
<b>ILS005LT</b>	5 µL	ILS 5 µL Gas-tight Luer tip	6.5 mm	54.1 mm	Y	Y
<b>ILS010LT</b>	10 µL	ILS 10 µL Gas-tight Luer tip	6.5 mm	54.1 mm	Y	Y
<b>ILS025LT</b>	25 µL	ILS 25 µL Gas-tight Luer tip	8.0 mm	60 mm	Y	Y
<b>SGE050TLL</b>	50 µL	SGE 50 µL Gas-tight Teflon Luer Lock	8.0 mm	60 mm	Y	Y
<b>SGE100TLL</b>	100 µL	SGE 100 µL Gas-tight Teflon Luer Lock	8.0 mm	60 mm	Y	Y
<b>SGE250TLL</b>	250 µL	SGE 250 µL Gas-tight Teflon Luer Lock	8.0 mm	60 mm	Y	N

#### Syringes with Replaceable Beveled Needles

Order No.	Volume	Description	O.D.	SCALE LENGTH	UMP3	UMP2
<b>SGE0005RN*</b>	0.5 µL	SGE 0.5 µL 23 ga (0.63 mm), 70 mm long needle	8.0 mm	54.1 mm	Y	Y
<b>SGE001RN*</b>	1.0 µL	SGE 1.0 µL 26 ga (0.47 mm), 70 mm long needle	8.0 mm	54.1 mm	Y	Y
<b>SGE005RN</b>	5 µL	SGE 5 µL 23 ga (0.63 mm), 70 mm long needle	8.0 mm	54.1 mm	Y	Y
<b>SGE010RNS</b>	10 µL	SGE 10 µL 26 ga (0.47 mm), 50 mm long needle	8.0 mm	54.1 mm	Y	Y
<b>SGE025RN</b>	25 µL	SGE 25 µL 25 ga (0.50 mm), 50 mm long needle	8.0 mm	60 mm	Y	Y
<b>SGE050RN</b>	50 µL	SGE 50 µL 25 ga (0.50 mm), 50 mm long needle	8.0 mm	60 mm	Y	Y
<b>SGE100RN</b>	100 µL	SGE 100 µL 25 ga (0.50 mm), 50 mm long needle	8.0 mm	60 mm	Y	Y

\* The capacity of this syringe is so small that the entire sample is contained within the needle. The plunger extends to the tip of the needle, displacing the full sample during injection — which gives the syringe zero dead volume.

SGE and ILS are respective trademarks of Scientific Glass Engineering and Innovative Labor Systeme.

#### Replacement Needles

<b>RN0005</b>	For Syringe SGE0005RN, 23 ga (0.63 mm) 70 mm long
<b>RN001</b>	For Syringe SGE001RN, 26 ga (0.47 mm) 70 mm long
<b>RN005</b>	For Syringe SGE005RN, 23 ga (0.63 mm) 50 mm long
<b>RN010</b>	For Syringe SGE010RN(S), 26 ga (0.47 mm) 50 mm long, 5-pack
<b>RN025</b>	For Syringes SGE025RN, SGE050RN, SGE0100RN, 26 ga (0.47 mm) 50 mm long, 5-pack

# Sub-Microliter Injection System

*Includes the smallest dead volume injection when the 10 µL syringe is used with WPI needles 34-36g*

## Features

- The smallest dead volume injection syringe
- Various needle sizes available from 26 ga. and 33-36 ga.
- Blunt or sharp needles
- Compatible with WPI's UMP3 microinjection system

## Benefits

- Low dead volume (0.5 µL or less)
- Switching the syringe tip during an experiment is easy
- Variety of tips.

## Applications

- Animal research
- Capillary electrophoresis
- Versatile research applications — RPE and IO Kits

**NanoFil™** is a specially designed 10 µL syringe developed in response to customer requests for improved microinjection in mice and other small animals. It makes quantitative nanoliter injection much easier and more accurate than any other method currently in use.

### Low Dead Volume

NanoFil's low dead volume eliminates the need for oil backfilling, a messy process which risks contamination of the injected sample. Injection is now simpler, and less messy, and there is no possibility of oil contamination in critical applications such as ophthalmology research (see the Retinal Pigment Epithelial (RPE) and Intra Ocular (IO) injection kits listed below).

### Easily Switch Syringe Tip

When the inner tip diameter of a conventional syringe is reduced to less than 100 µm, it is very difficult to front fill the solution at a reasonable speed. NanoFil solves this problem by using a tip coupling mechanism that makes it possible to change the syringe tip during the experiment. Simply load the sample using a larger tip, such as the 26 gauge needle provided with the syringe, and then replace it with a micro tip for sample injection. On a conventional 10 µL syringe, a solid ring or bushing is permanently bonded to the tubing. Replacing the tip in the middle of the experiment is not practical. With NanoFil, tips can be exchanged by a simple twist of the brass lock, gently pulling out the tip, and replacing with the desired new tip.

### Holds Metal Tips and Quartz Tubing

To secure the tip, NanoFil uses an olive-shaped silicon gasket that is similar to, but much sturdier than, some of the microelectrode holders used for electrophysiology recording. The silicone gasket makes it possible to hold not only metal needles but also **Silflex** tubing. Many types of tubing can be easily connected to the syringe as long as the outer diameter (OD) is close to, but not more than, the barrel inner diameter (ID) of 460 µm. Flexible quartz capillaries used in Gas Chromatography (GC) and Capillary Electrophoresis (CE) can also be easily coupled to the syringe.

### Variety of Tips

Specially designed needles as small as 36 gauge (110 µm OD) are offered in both blunt and beveled styles. Our studies have shown that these needles will cause less trauma to the tissue. NanoFil has a unique coupling mechanism that allows many different forms of small tubing and tips to be coupled with the syringe barrel.

### Selecting the correct tip for your application

The replaceable needles used with NanoFil are available with either blunt or beveled tips. The blunt tip is used for injection into soft tissue and when a uniform solution distribution is needed. The beveled style is



used for applications that involve the penetration of a tough tissue.

One of the main factors that can affect the resolution and accuracy of a microinjection in the low nanoliter range is diffusion from the needle opening. When the tip ID is equal to or larger than 100 µm, the error caused by tip diffusion can be in the nanoliter range level [(100 micron)<sup>3</sup> = 1 nanoliter]. With a 36 gauge needle installed on the NanoFil, the error caused by diffusion will be reduced, making accurate injection of a nanoliter possible.

All of WPI's beveled needles have a unique 25° tri-surface bevel that is optimized for microinjection. A 10 degree single-surface beveled tip penetrates better than one with a 25° angle, however the distance between the upper opening to the tip (the dimension "F" in the table on page 59) is longer. As a result, it requires deeper penetration of the tip. This can be an issue for microinjection into very small areas where the dimensions of the anatomy can't accommodate the required depth of insertion. WPI's unique 25° beveled tip reduces the required insertion depth by incorporating two extra beveled surfaces. The edge of a single surface beveled tip is actually a blade instead of a point. It dulls very quickly. In contrast, the tri-surfaced tip has a sharp point. It penetrates more easily and is more durable. Our tests show that our 33 gauge, 25 degree beveled tip penetrates easier and lasts longer than other manufacturers' 33 gauge, 10 degree single beveled tips. With a 35 gauge tri-surface beveled tip, the resistance to the penetration becomes even less. Each of our needles undergo a penetration test before leaving the factory to guarantee the best results for our customers.

### Available Tips

**33 gauge:** This tip is similar to Hamilton's 7762 and 7803 series removable needles in both tip length and outer diameter. However, our beveled tip version is shorter, more durable, and penetrates better due to the special tri-surface grinding technique. In the past, 33 gauge tips were the smallest size sold by other manufacturers and were frequently cited in literature. However, our new 35 gauge tip is much better for injections involving small animals, especially mice. Compared with Hamilton's 33 gauge, 10 degree beveled tip, our 35 gauge 25 degree beveled tip can reduce the depth of penetration by almost 80%. The distance between the tip and the upper rim of the opening (dimension F on the drawing) is 348 µm for the 33 gauge tip. The distance for our 35 gauge tip is only 230 µm. In addition, the smaller tip size significantly reduces the required penetration force. In nearly all applications, a 33 gauge tip can be replaced with our 35 gauge tip and produce better results.

### 34 gauge

This is a transitional size between the 33 gauge and 35 gauge. If the 35 gauge is too weak and the 33 gauge is too large, this makes a good alternative.

### 35 gauge

This was the most popular and preferred tip of most scientists during our field trial. The combination of its strength, length, durability and clogging resistance creates a balance with very little compromising of the individual properties. It is much smaller than the 33 gauge tip offered by other manufacturers. It is only slightly larger than the 36 gauge tip but is much stronger and less likely to be clogged. Samples can be directly loaded with this tip. Its 5 mm length is sufficient enough for almost all injection applications in mice.

### 36 gauge

This is the smallest tip that is commercially available. The tip is so small that it can be inserted into the opening of the 33 gauge needle tip. Because this is pushing the limits of what current technology can

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
Brazil: (013) 406-29703 • info@brazil.wpiinc.com • www.wpiinc.com

Germany: +49 (0)6031 1602171 • wptide@wpi-europe.com • www.wpi-europe.com  
China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.com



glass barrel:  $\varnothing$  7 mm

**NANOFIL**

produce, there are some limitations to consider before using it. Its thin diameter makes it necessary to limit its length to 2.5 to 3 mm in order to maintain a usable strength. Since the tip ID is in the 25 to 50  $\mu$ m range, it is very easily clogged. Therefore, only well filtered solutions can be used. Depending on the viscosity of the sample, you might also need to pre-load the syringe with a regular tip before switching to this tip for injection. We recommend using the 35 gauge tip instead of the 36 gauge unless it is absolutely necessary.

**Flexible Quartz Tubing**

The flexible quartz tubing tip is made of 160  $\mu$ m OD polyimide coated quartz tubing with a special adapter sleeve mounted at the end. It is designed for filling glass capillary electrodes or pipettes, just like WPI's traditional MF34G MicroFil. However, unlike the traditional MicroFil, which has about 50  $\mu$ L of dead volume in its Luer hub, the dead volume of this tip is less than 0.589  $\mu$ L. It is useful for loading electrodes with solutions that have a limited volume or are too expensive to waste.

**References**

Arriaga, G., Macopson, J. J., & Jarvis, E. D. (2015). Transsynaptic Tracing from Peripheral Targets with Pseudorabies Virus Followed by Cholera Toxin and Biotinylated Dextran Amines Double Labeling. *Journal of Visualized Experiments*, (103), e50672-e50672. <http://doi.org/10.3791/50672>

Delotterie, D. F., Mathis, C., Cassel, J.-C., Rosenbrock, H., Dorner-Ciossek, C., & Marti, A. (2015). Touchscreen tasks in mice to demonstrate differences between hippocampal and striatal functions. *Neurobiology of Learning and Memory*, 120, 16-27. <http://doi.org/10.1016/j.nlm.2015.02.007>


Han, Z., Banworth, M. J., Makkia, R., Conley, S. M., Al-Ubaidi, M. R., Cooper, M. J., & Naash, M. I. (2015). Genomic DNA nanoparticles rescue rhodopsin-associated retinitis pigmentosa phenotype. *FASEB Journal : Official Publication of the Federation of American Societies for Experimental Biology*, 29(6), 2535-44. <http://doi.org/10.1096/fj.15-270363>

**Using NanoFil™ in different configurations**

**Direct injection by hand:** This is the simplest and most economical way to inject. Any of our tips can be inserted directly into the NanoFil™ syringe. Even the SilFlex tubing can be inserted to switch from hand injection to the other methods listed below. This method is limited by the accuracy of plunger movement that is achievable with a human hand.

**Installed on WPI's UMP3 microsyringe pump:** This will allow the user to achieve nanoliter resolution and reproducibility. For neural system injection, mount the UMP3 on a stereotaxic frame.

**SilFlex tubing and holder:** The needle is mounted on a small plastic holder that is connected to the NanoFil by a 35 cm length of flexible tubing. The NanoFil syringe is mounted on the UMP3 pump. This configuration allows the user to hold the animal in one hand and insert the needle with the other. When the needle reaches the desired location, activate the pump using the foot switch and the pre-programmed injection volume will be delivered. This configuration gives a nanoliter level of accuracy and reproducibility. It is best suited for applications such as the RPE and IO injection.



SPECIFICATIONS										
Tip Order Number	Gauge	Tip O.D. "A"	Tip I.D. "B" nominal	Tip Length "C"	Total Length "D"	Shank O.D. "E"	Bevel Length "F"	Total Dead Volume	Tip Material	
NF33BV-2	33	210 $\mu$ m	115 $\mu$ m	10 mm	40 mm	460 $\mu$ m	$\approx$ 348 $\mu$ m	0.416 $\mu$ L	Stainless Steel	
NF34BV-2	34	185 $\mu$ m	85 $\mu$ m	5 mm	35 mm	460 $\mu$ m	$\approx$ 290 $\mu$ m	0.199 $\mu$ L	Stainless Steel	
NF35BV-2	35	135 $\mu$ m	55 $\mu$ m	5 mm	35 mm	460 $\mu$ m	$\approx$ 204 $\mu$ m	0.435 $\mu$ L	Stainless Steel	
NF36BV-2	36	120 $\mu$ m	35 $\mu$ m	3 mm	33 mm	460 $\mu$ m	$\approx$ 156 $\mu$ m	0.340 $\mu$ L	Stainless Steel	
NFQ34-5	34	160 $\mu$ m	100 $\mu$ m	55 mm	75 mm	460 $\mu$ m	n/a	0.589 $\mu$ L	Quartz	
NF33BL-2	33	210 $\mu$ m	115 $\mu$ m	10 mm	34 mm	460 $\mu$ m	$\approx$ 0	0.416 $\mu$ L	Stainless Steel	
NF34BL-2	34	185 $\mu$ m	85 $\mu$ m	5 mm	29 mm	460 $\mu$ m	$\approx$ 0	0.199 $\mu$ L	Stainless Steel	
NF35BL-2	35	135 $\mu$ m	55 $\mu$ m	5 mm	29 mm	460 $\mu$ m	$\approx$ 0	0.435 $\mu$ L	Stainless Steel	
NF36BL-2	36	120 $\mu$ m	35 $\mu$ m	3 mm	27 mm	460 $\mu$ m	$\approx$ 0	0.340 $\mu$ L	Stainless Steel	
<b>Silflex</b>			100 $\mu$ m		35 cm			2.749 $\mu$ L		
NF26BV-2	26	460 $\mu$ m	110 $\mu$ m	3 mm	40 mm	460 $\mu$ m		0.380 $\mu$ L		

**ORDERING INFORMATION**

**NANOFIL** NanoFil™ Syringe, 10 microliter  
Includes two 28-gauge MicroFil needles MF28G; one 1cc plastic syringe; and one 26-gauge beveled needle NF26BV.

**NANOFIL-100** NanoFil™ Syringe, 100 microliter  
26-gauge beveled needle NF26BV included.

**REPLACEMENT NEEDLE**

**NF26BV-2** 26G Beveled Needle, 460  $\mu$ m nominal diam. (package of 2)

**NanoFil™ NEEDLES**

- NF33BL-2 33 G blunt NanoFil™ Needle (pkg. of 2)
- NF34BL-2 34 G blunt NanoFil™ Needle (pkg. of 2)
- NF35BL-2 35 G blunt NanoFil™ Needle (pkg. of 2)
- NF36BL-2 36 G blunt NanoFil™ Needle (pkg. of 2)
- NF33BV-2 33 G beveled NanoFil™ Needle (pkg. of 2)
- NF34BV-2 34 G beveled NanoFil™ Needle (pkg. of 2)
- NF35BV-2 35 G beveled NanoFil™ Needle (pkg. of 2)
- NF36BV-2 36 G beveled NanoFil™ Needle (pkg. of 2)
- NF33-36BL Assortment of 4 blunt NanoFil™ Needles
- NF33-36BV Assortment of 4 beveled NanoFil™ Needles

**OPTIONAL ACCESSORIES/REPLACEMENT PARTS**

- NFINHLD NanoFil™ Injection Holder
- SILFLEX-2 SilFlex Tubing 35 cm long (pkg. of 2) (dead volume = 2.74  $\mu$ L)
- NFGSK-5 Spare Silicone Gasket for NanoFil™ & Holder (pkg. of 5)
- NFQ34-5 34 Gauge Flexible Quartz Tubing for filling (pkg. of 5)

# NanoFil™ Application Kits

*Designed for eye research*

## Features

- Includes SilFlex tubing, gasket, holder and tip assortment
- NanoFil™ syringe sold separately
- Less than 3 µL dead volume
- Achieve accurate, repetitive, oil-free injections

## Benefits

- No oil back filling necessary
- Comes with four needle sizes included in the kit

## Applications


- Intraocular injection
- Retinal pigment epithelium injection
- Mouse brain injection

These kits are specially designed for eye research for injecting retinal pigment epithelium (RPE) and intraocular (IO) in addition to brain injection in mice. They are used exclusively with a **NanoFil™** syringe and **UMP3** to achieve accurate and repeatable oil free injections down to submicroliter ranges. Each kit includes two pieces of Silflex tubing, a holder assembly, spare gaskets, and an assortment of four needles — blunt for the RPE kit and beveled tips for the IO kit. Each kit comes with one each of 33, 34, 35 and 36 gauge needles so that first time users can find the best size for their application.

The Silflex tubing is a very important component of the kit. This 35 cm long, flexible tubing has a precise outer diameter for airtight fitting with the syringe. It also has a small inner diameter to minimize dead volume. The SilFlex is coupled to the injection tip with a seal system similar to that of the NanoFil. The dead volume of the entire kit (including the tubing) is less than 3 microliters. All of the components in the kit are constructed of inert, solvent resistant materials for easy cleaning after viral injection.



RPE-KIT



RPE-KIT Holder assembly

## ORDERING INFORMATION

<b>RPE-KIT</b>	Retinal Pigment Epithelium (RPE) Injection Kit <i>SilFlex tubing, gasket, holder, and blunt tip mix</i>
<b>IO-KIT</b>	Intraocular (IO) Injection Kit <i>SilFlex tubing, holder, gasket, 4 beveled tips (33g, 34g, 35 g, 36g)</i>
<b>503207</b>	Stand & Clamps
<b>NFINHLD</b>	NanoFil Injection Holder
<b>NFINHLD-G10</b>	1.0 mm Glass Pipette Holder for NanoFil Syringe

## References

Park, S. W., Kim, J. H., Park, W. J., & Kim, J. H. (2015). Limbal Approach-Subretinal Injection of Viral Vectors for Gene Therapy in Mice Retinal Pigment Epithelium. *Journal of Visualized Experiments*, (102), e53030–e53030. <http://doi.org/10.3791/53030>



NFINHLD-G10

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • [wpiuk@wpi-europe.com](mailto:wpiuk@wpi-europe.com) • [www.wpi-europe.com](http://www.wpi-europe.com) Germany: +49 (0)6031 1602171 • [wptide@wpi-europe.com](mailto:wptide@wpi-europe.com) • [www.wpi-europe.com](http://www.wpi-europe.com)  
 Brazil: (013) 406-29703 • [info@brazil.wpiinc.com](mailto:info@brazil.wpiinc.com) • [www.wpiinc.com](http://www.wpiinc.com) China: +86 21 6888 5517 • [chinasales@china.wpiinc.com](mailto:chinasales@china.wpiinc.com) • [www.wpiinc.net](http://www.wpiinc.net)

# Nanoliter Injector

For oocyte injection and applications in the 2 to 70nL range

## Features

- Microprocessor-controlled injector
- Direct piston displacement in capillary glass
- Filling and injection speeds: 23nL/sec and 46 nL/sec
- Optional standard or **Micro2T** controller
- Optional foot switch

## Benefits

- No syringes required
- Graphical representation of volume status
- Automatically calculated injections based on syringe volume and scale length
- Intuitive touchscreen interface

## Applications

- Oocyte injection

WPI's microprocessor-controlled **Nanoliter 2010** uses direct piston displacement. By either pushing the injection button on the control box or pressing on the optional foot switch, a discrete volume will be injected. Choice of capillary filling and injection speeds are 23nL/sec or 46nL/sec (emptying speed is 92nL/sec). Maximum fluid ejection is 5 µL. Each unit comes with sufficient glass to pull at least 300 micropipettes (see **PUL-1000**, page 48). Glass is 1.14mm O.D. (nominal) and 0.5 mm I.D.

By setting the DIP switch, the injection volume can be changed from 2.3 to 69.0nL in 16 steps. Up to 100 injections may be performed from a single filling of the micropipette. Since the volume of a normal *Xenopus* oocyte is about 500nL, the instrument has the capability to inject from less than 1% to over 10% of the total volume of the oocyte in one preset step increment.

Included: 1 vial 3.5 in. capillaries (300), replacement "O" rings, Allen wrench, MicroFil **MF34G** backfilling needle and two sample **µTip™** pre-pulled micropipettes.

### New SMARTouch controller

Already own a unit with the standard controller (shown below) or a **Micro4** digital controller? Now you can upgrade to touchscreen control with the new **SMARTouch**.

## References

**Chouchane, M., et al.** (2017). Lineage Reprogramming of Astroglial Cells from Different Origins into Distinct Neuronal Subtypes. *Stem Cell Reports*, 9(1), 162–176. <http://doi.org/10.1016/j.stemcr.2017.05.009>

**Moreno-Juan, V., et al.** (2017). Prenatal thalamic waves regulate cortical area size prior to sensory processing. *Nature Communications*, 8, 14172. <http://doi.org/10.1038/ncomms14172>

**Tappert, L., et al.** (2017). Sublethal doses of imidacloprid disrupt sexual communication and host finding in a parasitoid wasp. *Scientific Reports*, 7(1), 42756. <http://doi.org/10.1038/srep42756>

**Lynagh, T., et al.** (2017). Unique Contributions of an Arginine Side Chain to Ligand Recognition in a Glutamate-gated Chloride Channel. *The Journal of Biological Chemistry*, 292(9), 3940–3946. <http://doi.org/10.1074/jbc.M116.772939>

**Li, S., et al.** (2017). Small RNA-Seq analysis reveals microRNA-regulation of the lmd pathway during *Escherichia coli* infection in *Drosophila*. *Developmental & Comparative Immunology*, 70, 80–87. <http://doi.org/10.1016/j.dci.2017.01.008>

The standard configuration of the Nanoliter 2010 includes the small controller which is simple to setup and operate.

NANOLITER2010



Micromanipulator sold separately.

## NANOLITER 2010 SPECIFICATIONS

REMOTE CONTROL	Yes
GLASS OD	1.14 mm
GLASS ID	0.5 mm
STEP	12.7 µm/step
INJECTION SPEED	
Slow	23 nL/sec
Fast	46 nL/sec
FILL SPEED	
Slow	23 nL/sec
Fast	46 nL/sec
EMPTY SPEED	92 nL/sec
SINGLE STEP VOLUME RANGE	2.3 - 69.0 nL
SMALLEST VOLUME	2.3 nL
LARGEST SINGLE STEP VOLUME	69 nL
TO CHANGE VOLUME	Set switch
INJECTIONS PER FILLING, MAX.	100 injections
SHIPPING WEIGHT	3 lb. (1.1 kg)

## ORDERING INFORMATION

**NANOLITER2010** Nanoliter 2010 and Standard Controller  
System includes Nanoliter 2010 injector, small controller, 1 vial 3.5 in. capillaries (300), replacement "O" rings, Allen wrench, MicroFil™ MF34G backfilling needle and two sample µTip™ pre-pulled micropipettes.

**NL2010MC2T** Nanoliter Injector & SMARTouch® Controller  
(small controller not included)  
System includes Nanoliter 2010 injector, MICRO2T controller, 1 vial 3.5 in. capillaries (300), replacement "O" rings, Allen wrench, MicroFil™ MF34G backfilling needle and two sample µTip™ pre-pulled micropipettes.

## OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>13142</b>	Foot Switch for Nanoliter 2010
<b>504127</b>	Replacement Nanoliter 2010 Injector Head
<b>504949</b>	Replacement 3.5-in. Glass Capillaries, 1.14 mm OD (300)
<b>504950</b>	Replacement 7-in. Glass Capillaries, 1.14 mm OD (300)
<b>TIP10XV119</b>	Micropipettes (10) Pre-pulled 10 µm
<b>40239</b>	Adapter for SMARTouch
<b>300521</b>	Replacement O-rings (five)
<b>500778</b>	Replacement Nanoliter Injector Mounting Adapter
<b>300033</b>	Adapter for Micro4

# Touchscreen Programmable Pump Series

*Legato syringe pumps for reliable delivery with the ease of a touch screen*

## Features

- High resolution color touch screen
- Real time clock
- Touch pad "lock" feature
- Full metal chassis
- Built in syringe table with up to 75 lb. linear force
- Built in RS-485 interface to link multiple pumps
- USB port & RS232 Interface, TTL interface
- Continuous mode of operation
- Protection with a spill dam
- Analog control option
- CE, UL, CSA, CB Scheme, EU RoHS compliance

## Benefits

- Automatic dispensing of small volumes
- Very precise flow rate control
- Hands free operation with foot pedal
- Better flow performance with accuracy  $\pm 0.35\%$

## Applications

- Drug administration
- Chemical applications with slow incorporation of fixed volumes of fluids

The large touch screen color display lets you see all of the pump's operating parameters to ensure proper operation during the experiments. Syringe size and flow rate are easily displayed, as well as the volume delivered and elapsed time. Set up is easy using the icon-driven software. An LED on the front panel makes it easy to see if the pump is running. Advanced micro stepping techniques are employed to further reduce the step angle to eliminate flow pulsation. Accuracy is  $\pm 0.5\%$ . A wide dynamic flow range from picoliters per minute to milliliters per minute can be programmed into the pump. These versatile pumps can be connected through an RS485 interface. Add the new Adagio software to maximize the use of the pump's functions and features. Adagio allows you to configure the pump through the software, as well as operate one or multiple pumps. LabVIEW drivers are available on the National Instruments website.

### SPLG100 SERIES FLOW RATES

Syringe	Diameter	Minimum	Maximum
0.5 $\mu$ L	0.103 mm	1.260 pL/min	1.325 $\mu$ L/min
1 $\mu$ L	0.146 mm	2.520 pL/min	2.651 $\mu$ L/min
2 $\mu$ L	0.206 mm	5.100 pL/min	5.299 $\mu$ L/min
5 $\mu$ L	0.343 mm	14.100 pL/min	14.690 $\mu$ L/min
10 $\mu$ L	0.485 mm	28.260 pL/min	29.380 $\mu$ L/min
25 $\mu$ L	0.729 mm	63.900 pL/min	66.370 $\mu$ L/min
50 $\mu$ L	1.03 mm	127.600 pL/min	132.500 $\mu$ L/min
100 $\mu$ L	1.457 mm	255.20 pL/min	265.100 $\mu$ L/min
250 $\mu$ L	2.304 mm	638.300 nL/min	662.900 $\mu$ L/min
500 $\mu$ L	3.256 mm	1.275 nL/min	1.324 mL/min
1000 $\mu$ L	4.608 mm	2.553 nL/min	2.652 mL/min
1 mL	4.699 mm	2.655 nL/min	2.757 mL/min
3 mL	8.585 mm	8.863 nL/min	9.204 mL/min
5 mL	11.989 mm	17.290 nL/min	17.950 mL/min
10 mL	14.427 mm	25.030 nL/min	25.990 mL/min
20 mL	19.05 mm	43.640 nL/min	45.320 mL/min
30 mL	21.59 mm	56.050 nL/min	58.210 mL/min
60 mL	26.594 mm	85.050 nL/min	88.320 mL/min

## SPLG100 Infuse-Only Syringe Pump

The **SPLG100** was the first single-syringe infusion-only pump with a touchscreen interface. The **SPLG100** has a wide flow rate range from 1.26 pL/min to 88.32 mL/min, depending on syringe size. It accommodates a single syringe from 0.5  $\mu$ L to 60 mL. Any type of syringe, including glass, plastic or stainless steel are held securely in place.



CE


SPLG100

### SPLG100 SPECIFICATIONS

SYRINGE SIZE	0.5 $\mu$ L to 60 mL
POWER	100-240 VAC: 50/60 Hz, 50W. 0.5 A fuse
MOTOR DRIVE CONTROL	Microprocessor with 1/16 microstepping
LINEAR FORCE (MAXIMUM)	13.6 kg (30 lb.) @ 100% force selection
NUMBER OF MICROSTEPS PER REVOLUTION OF LEAD SCREW	15,360
STEP RATE (MIN.)	27.5 sec/ $\mu$ step
STEP RATE (MAX.)	26 $\mu$ sec/ $\mu$ step
DRIVE MOTOR	0.9 degree Stepper Motor
PUSHER TRAVEL RATE (MIN.)	0.15 $\mu$ m/min
PUSHER TRAVEL RATE (MAX.)	159 mm/min
FLOW RATE (MIN.)	1.26 pL/min (0.5 $\mu$ L syringe)
FLOW RATE (MAX.)	88.32 mL/min (60 mL syringe)
DIMENSIONS	22.6 x 19.05 x 15 cm (9 x 7.5 x 5 in)
WEIGHT	2.66 kg (5.9 lb.)
CONNECTORS	RS485 - IEEE-1394 6 pos, USB Type B

## SPLG101 Dual Infuse-Only Syringe Pump

The **SPLG101** is ideal for applications where dual syringes are required with small volumes up to 10 mL. It accommodates two syringes from 0.5  $\mu$ L to 10 mL. The **SPLG101** has a wide flow rate range from 1.26pL/min to 25.99 mL/min, depending on syringe size.



CE

SPLG101

### SPLG101 SPECIFICATIONS

SYRINGE SIZE	0.5 $\mu$ L to 10 mL
POWER	100-240 VAC: 50/60 Hz, 50W. 0.5 A fuse
MOTOR DRIVE CONTROL	Microprocessor with 1/16 microstepping
LINEAR FORCE (MAX.)	13.6 kg (30 lbs.) @ 100% force selection
NUMBER OF MICROSTEPS PER REVOLUTION OF LEAD SCREW	15,360
STEP RATE (MIN.)	27.5 sec/ $\mu$ step
STEP RATE (MAX.)	26 $\mu$ sec/ $\mu$ step
DRIVE MOTOR	0.9 degree Stepper Motor
PUSHER TRAVEL RATE (MIN.)	0.15 $\mu$ m/min
PUSHER TRAVEL RATE (MAX.)	159 mm/min
FLOW RATE (MIN.)	1.26 pL/min (0.5 $\mu$ L syringe)
FLOW RATE (MAX.)	25.99 mL/min (10 mL syringe)
DIMENSIONS	22.6 x 19.05 x 15 cm (9 x 7.5 x 5 in)
WEIGHT	2.66 kg (5.9 lbs.)
CONNECTORS	RS-232 - 9 Pin D-Sub connector, RS485 - IEEE-1394 6 pos, USB - Type B


## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
 Brazil: (013) 406-29703 • info@brazil.wpiinc.com • www.wpiinc.com

Germany: +49 (0)6031 1602171 • wvide@wpi-europe.com • www.wpi-europe.com  
 China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net

## SPLG110 Infuse/Withdraw Syringe Pump

The **SPLG110** offers infuse/withdraw flow control and programmability for up to two multi-step programs of 50 steps each. The **SPLG110** has a wide flow rate range from 1.26pL/min to 88.28mL/min, depending on syringe size. The **SPLG110**



SPLG110



accommodates a single syringe from 0.5 µL to 60 mL. Any type of syringe can be used in the unit including glass, plastic or stainless steel. The pump is ideal for more complex multi-step dosing and has multi-mode operation including infusion only, withdrawal only, infusion and withdrawal and withdrawal/infusion modes.

### SPLG110 SPECIFICATIONS

SYRINGE SIZE	0.5 µL to 60 mL
POWER	100-240 VAC: 50/60 Hz, 50W. 0.5 A fuse
MOTOR DRIVE CONTROL	Microprocessor with 1/16 microstepping
LINEAR FORCE (MAX.)	13.6 kg (30 lb.) @ 100% Force Selection
NUMBER OF MICROSTEPS PER REVOLUTION OF LEAD SCREW	15,360
STEP RATE (MIN.)	27.5 sec/µstep
STEP RATE (MAX.)	26 µsec/µstep
DRIVE MOTOR	0.9 degree Stepper Motor
PUSHER TRAVEL RATE (MIN.)	0.15 µm/min
PUSHER TRAVEL RATE (MAX.)	159 mm/min
FLOW RATE (MIN.)	1.26 pL/min (0.5 µl syringe)
FLOW RATE (MAX.)	88.28 mL/min (60 ml syringe)
DIMENSIONS	22.6 x 19.05 x 15 cm (9 x 7.5 x 5 in)
WEIGHT	2.66 kg (5.9 lb.)
CONNECTORS	RS-232 - 9 Pin D-Sub connector, RS485 - IEEE-1394 6 pos, USB - Type B

## SPLG210 Infuse/Withdraw Syringe Pump

The **SPLG210** Infuse/Withdraw syringe pump is easy to use with the high resolution touch screen. The basic model works with one syringe or two (from 0.5 µL to 140 mL) and can be reconfigured in the field to be used with multiple syringes. A protective cover over the display prevents leakage into the display. To optimize your bench space,



SPLG210



the **SPLG210** can be placed on its side to reduce the footprint to only 3.5 x 9.75 in. The display also rotates with the change to allow you to operate the pump vertically. The programmable model offers maximum flexibility for configuring and running different programs. Up to 40 programs of 20 steps each can be configured and stored for quick recall with the touch of a button.

### SPLG210 SPECIFICATIONS

SYRINGE SIZE	0.5 µL to 140 mL
POWER	100-240 VAC: 50/60 Hz, 50W. 0.5 A fuse
MOTOR DRIVE CONTROL	Microprocessor with 1/16 microstepping
LINEAR FORCE (MAX.)	34 kg (75 lb.) @ 100% force selection
NUMBER OF MICROSTEPS PER REVOLUTION OF LEAD SCREW	6400
STEP RATE (MIN.)	27.5 sec/µstep
STEP RATE (MAX.)	26 µsec/µstep
DRIVE MOTOR	1.8 degree stepper motor
PUSHER TRAVEL RATE (MIN.)	0.36 µm/min
PUSHER TRAVEL RATE (MAX.)	190.80 mm/min
FLOW RATE (MIN.)	5 pL/min (0.5 µL syringe)
FLOW RATE (MAX.)	215.803 mL/min (140 mL syringe)
DIMENSIONS	8.89 x 25.4 x 27.94 cm (3.5 x 10 x 11 in))
WEIGHT	4.9 kg (10.75 lb.)
CONNECTORS	RS-232 - 9 Pin D-Sub connector, RS-485 - IEEE-1394 6 pos, USB - Type B, I/O & TTL - 15 Pin D-Sub connector

## SPLG270 Push/Pull Continuous

The **SPLG270** is a Push-Pull syringe pump. It accommodates two syringes from 0.5 µL to 140 mL for infusion and two syringes for withdrawal. This model supports infusion and withdrawal simultaneously at user-defined flow rates and with selectable target volumes to control the total volume pumped. It also supports infuse only, withdraw only, infuse/withdraw, withdraw/infuse and continuous mode. The touch screen interface lets you quickly create configurations and recall them for easy use. The 4.3" TFT color display with touch pad interface presents all the pump operating parameters on one easy-to-view run screen. Protective cover over the display prevents leakage into the display.



SPLG270



### SPLG270 SPECIFICATIONS

SYRINGE SIZE	0.5 µL to 140 mL
POWER	100-240 VAC: 50/60 Hz, 50W. 0.5 A fuse
MOTOR DRIVE CONTROL	Microprocessor with 1/16 microstepping
LINEAR FORCE (MAX.)	34 kg (75 lb.) @ 100% force selection
NUMBER OF MICROSTEPS PER REVOLUTION OF LEAD SCREW	6400
STEP RATE (MIN.)	27.5 sec/µstep
STEP RATE (MAX.)	26 µsec/µstep
DRIVE MOTOR	1.8° Stepper Motor
PUSHER TRAVEL RATE (MIN.)	0.36 µm/min
PUSHER TRAVEL RATE (MAX.)	190.80 mm/min
FLOW RATE (MIN.)	5 pL/min (0.5 µL syringe)
FLOW RATE (MAX.)	215.803 mL/min (140 mL syringe)
DIMENSIONS	8.89 x 25.4 x 27.94 cm (3.5 x 10 x 11 in)
WEIGHT	4.9 kg (10.75 lb.)
CONNECTORS	RS-232 - 9 Pin D-Sub connector, RS-485 - IEEE-1394 6 pos, USB - Type B, I/O & TTL - 15 Pin D-Sub connector

## Adagio Software

- Low cost, simple installation
- Flow evolution graph
- Import and export programs
- Quick and easy manual pump control
- Monitor one or more pumps
- Program data logging

The manual pump control tool allows easy direct control of the pump. Pump commands can be entered directly into the log. Multiple programs can be opened at the same time. The program's progression is tracked, and can be stored in a file for later access.

**Intuitive Run Screen** — Combining multiple parameters simultaneously with internationally recognizable icons allow the Legato™ Series to provide a new level of intuitive syringe pump operation.

### ORDERING INFORMATION

<b>SPLG100</b>	Legato 100 Syringe Pump, Infuse-Only
<b>SPLG101</b>	Legato 101 Syringe Pump, Dual Infuse-Only
<b>SPLG110</b>	Legato 110 Syringe Pump, Infuse/Withdraw
<b>SPLG200</b>	SPL Syringe Pump, Infuse Only
<b>SPLG210</b>	SPL Syringe Pump, Infuse/Withdraw
<b>SPLG212</b>	SPL Syringe Pump, Infuse/Withdraw Programmable
<b>SPLG270</b>	SPL Syringe Pump, Push-Pull
<b>SPLG272</b>	SPL Syringe Pump, Push-Pull Programmable
<b>504576</b>	Small Syringe Multi Rack (for six 30-60 mL syringes or ten 0.5 µL-20 mL syringes)
<b>504577</b>	Large Syringe Multi Rack (for up to four 60-140 mL plastic syringes)
<b>504578</b>	Software Adagio/USB Key

# Touchscreen Programmable Pump Series

*Legato syringe pumps for reliable delivery with the ease of a touch screen*

## Features

- High resolution color touch screen
- Real time clock
- Touch pad "lock" feature
- Full metal chassis
- Built in syringe table with up to 75 lb. linear force
- Built in RS-485 interface to link multiple pumps
- USB port & RS232 Interface, TTL interface
- Continuous mode of operation
- Protection with a spill dam
- Analog control option
- CE, UL, CSA, CB Scheme, EU RoHS compliance

## Benefits

- Automatic dispensing of small volumes
- Very precise flow rate control
- Hands free operation with foot pedal
- Better flow performance with accuracy  $\pm 0.35\%$



SPLG200

CE

## SPLG200 Infuse-Only Syringe Pump

The **SPLG200** Infuse only syringe pump is easy to use with the high resolution touch screen. The basic model works with one syringe or two (from 0.5  $\mu$ L to 140 mL) and can be reconfigured in the field to be used with multiple syringes. To optimize your bench space, the **SPLG200** can be placed on its side to reduce the footprint to only 3.5 x 9.75 inches. The display automatically reorients itself with the change to allow the user to operate the pump vertically. The pump features user definable flow rates with selectable target volumes and time values to control the infusion rate and the total volume. Up to 40 programs of 20 steps each can be configured and stored in the unit for quick recall with the touch of a button.

## Applications

- Drug administration
- Chemical applications with slow incorporation of fixed volumes of fluids



*When mounted vertically, the display screen of the SPLG series pumps automatically reorients for ease of use.*

## SPLG200 SERIES FLOW RATES

Syringe	Diameter	Minimum	Maximum
0.5 $\mu$ L	0.103 mm	3.12 pL/min	1.589 $\mu$ L/min
1 $\mu$ L	0.146 mm	6.18 pL/min	3.180 $\mu$ L/min
2 $\mu$ L	0.206 mm	12.301 pL/min	6.358 $\mu$ L/min
5 $\mu$ L	0.343 mm	33.96 pL/min	17.630 $\mu$ L/min
10 $\mu$ L	0.485 mm	67.72 pL/min	35.249 $\mu$ L/min
25 $\mu$ L	0.729 mm	153.42 pL/min	79.640 $\mu$ L/min
50 $\mu$ L	1.03 mm	306.24 pL/min	158.984 $\mu$ L/min
100 $\mu$ L	1.457 mm	612.72 pL/min	318.126 $\mu$ L/min
250 $\mu$ L	2.304 mm	1.533 nL/min	795.51 $\mu$ L/min
500 $\mu$ L	3.256 mm	3.06 nL/min	1.588 mL/min
1000 $\mu$ L	4.608 mm	6.129 nL/min	3.181 mL/min
1 mL	4.699 mm	6.373 nL/min	3.308 mL/min
3 mL	8.585 mm	21.272 nL/min	11.044 mL/min
5 mL	11.989 mm	41.485 nL/min	21.539 mL/min
10 mL	14.427 mm	60.073 nL/min	31.19 mL/min
20 mL	19.05 mm	104.74 nL/min	54.383 mL/min
30 mL	21.59 mm	134.533 nL/min	69.852 mL/min
50 mL	26.594 mm	204.122 nL/min	105.985 mL/min
100 mL	35.7 mm	367.839 nL/min	190.992 mL/min
140 mL	38.4 mm	415.623 nL/min	215.803 mL/min

## SPLG200 SPECIFICATIONS

SYRINGE SIZE	0.5 $\mu$ L to 140 mL
POWER	100-240 VAC: 50/60 Hz, 50W. 0.5 A fuse
MOTOR DRIVE CONTROL	Microprocessor with 1/16 microstepping
LINEAR FORCE (MAX.)	34 kg (75 lbs.) @ 100% force selection
NUMBER OF MICROSTEPS PER REVOLUTION OF LEAD SCREW	6400
STEP RATE (MIN.)	27.5 sec/ $\mu$ step
STEP RATE (MAX.)	26 $\mu$ sec/ $\mu$ step
DRIVE MOTOR	1.8° stepper motor
PUSHER TRAVEL RATE (MIN.)	0.36 $\mu$ m/min
PUSHER TRAVEL RATE (MAX.)	190.80 mm/min
FLOW RATE (MIN.)	5 pL/min (0.5 $\mu$ L syringe)
FLOW RATE (MAX.)	215.803 mL/min (140 mL syringe)
DIMENSIONS	8.89 x 25.4 x 27.94 cm (3.5 x 10 x 11 in.)
WEIGHT	4.9 kg (10.75 lb.)
CONNECTORS	RS-232 - 9 Pin D-Sub connector, S485 - IEEE-1394 6 pos, USB - Type B, I/O & TTL - 15 Pin D-Sub connector

## ORDERING INFORMATION

<b>SPLG200</b>	SPL Syringe Pump, Infuse Only
<b>504576</b>	Small Syringe Multi Rack (for six 30-60 mL syringes or ten 0.5 $\mu$ L-20 mL syringes)
<b>504577</b>	Large Syringe Multi Rack (for up to four 60-140 mL plastic syringes)
<b>504578</b>	Software Adagio/USB Key

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
Brazil: (013) 406-29703 • info@brazil.wpiinc.com • www.wpiinc.com

Germany: +49 (0)6031 1602171 • wpide@wpi-europe.com • www.wpi-europe.com  
China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net



# Affordable Syringe Pumps

Syringe pumps for high metering precision at low, pulse-free rates

## Features

- Sturdy and reliable
- Extremely simple to set up and use
- Surprisingly affordable.

## Benefits

- Automatic dispensing of small volumes
- Very precise flow rate control
- Hands-free operation with foot pedal
- Better flow performance with accuracy  $\pm 0.35\%$

## Applications

- Drug administration
- Chemical applications with slow incorporation of fixed volumes of fluids

On the SP series syringe pumps, the liquid crystal displays (LCDs) prompt you through setup:


1. Select syringe from table stored in the pump's memory and displayed on the LCD.
2. Enter the volume to be dispensed.
3. Enter the flow rate and press Start.

It's fast and simple. Your settings are permanently stored in memory — there's no need to re-enter them each day. SP pumps feature preset rate and volume control. Just set the volume you want dispensed. Volume is tracked continuously on the LCD display. Then, when the preset volume has been dispensed, the pump shuts off automatically. The easy-to-read digital display provides real-time readings using both parameters and values for clearer, mistake-free readings. The SP200 Series pumps offer TTL and RS-232C interfaces and automatic shutoff under stall conditions.

## SP100i Single-Syringe Infusion Pump

This inexpensive single-syringe infusion pump combines precision and simplicity with outstanding ease of use and durability.

- Holds any size syringe, 10  $\mu$ L to 50 mL
- Automatic volume control and shutoff
- Simple menu-driven setup: dispense volume, dispense flow rate, syringe diameter
- Last settings stored in permanent memory
- Continuous dispense volume display



SP100i

## SP120p Two-Syringe Push-Pull Pump

A second syringe mount has been added to the basic SP100i, with both syringes activated by a single pusher block for simultaneous infusion and withdrawal.

- All the features of SP100i
- Holds two syringes, from 10  $\mu$ L to 10 mL.



CE

SP120p

## SP210c Continuous Cycle Syringe Pump

The SP210c holds up to four syringes and can cycle continuously back and forth in a push-pull action. As two syringes are infusing, two other syringes are withdrawing at the same rate. At the end of the set volume the direction is automatically reversed and the next cycle begins. With the use of 2-way valves, the pump can empty and refill syringes for continuous dispensing.



SP210 c

CE

- Holds four syringes, 10 mL to 60 mL each

## SP220i Multi-Syringe Infusion Pump

Ideal for applications requiring multiple syringes, the SP220i is an adaptation of the SP200i and has been modified to hold up to 10 syringes.

- All features of SP200i
- Accommodates 10 syringes up to 10 mL or 6 syringes up to 50 mL or 4 syringes up to 140 mL.



SP220i

CE

## ORDERING INFORMATION

<b>SP100i</b>	Syringe Pump, Infusion (Single) 95-135V
<b>SP100iZ</b>	Syringe Pump, Infusion (Single) 220-240V
<b>SP101i</b>	Syringe Pump, Microdialysis (Double, Slow Speed) 95-135V
<b>SP101iZ</b>	Syringe Pump, Microdialysis (Double, Slow Speed) 220-240V
<b>SP120p</b>	Syringe Pump, Infusion-Withdrawal (Double) 95-135V
<b>SP120pZ</b>	Syringe Pump, Infusion-Withdrawal (Double) 220-240V
<b>SP200i</b>	Syringe Pump, Infusion (Double) 95-135V
<b>SP200iZ</b>	Syringe Pump, Infusion (Double) 220-240V
<b>SP210c</b>	Syringe Pump, Infusion-Withdrawal (Continuous) 95-135V
<b>SP210cZ</b>	Syringe Pump, Infusion-Withdrawal (Continuous) 220-240V
<b>SP210iw</b>	Syringe Pump, Infusion & Withdrawal (Double) 95-135V
<b>SP210iwZ</b>	Syringe Pump, Infusion & Withdrawal (Double) 220-240V
<b>SP220i</b>	Syringe Pump, Infusion (Multiple) 95-135V
<b>SP220iZ</b>	Syringe Pump, Infusion (Multiple) 220-240V
<b>SP230iw</b>	Syringe Pump, Infusion & Withdrawal (Multiple) 95-135V
<b>SP230iwZ</b>	Syringe Pump, Infusion & Withdrawal (Multiple) 220-240V
<b>SP250i</b>	Syringe Pump, Infusion (Multiple, Mixed Volumes) 95-135V
<b>SP250iZ</b>	Syringe Pump, Infusion (Multiple, Mixed Volumes) 220-240V
<b>SP260p</b>	Syringe Pump, I/W (Double) Single Cycle Action, 95-135V
<b>SP260pZ</b>	Syringe Pump, I/W (Double) Single Cycle Action, 220-240V

All 240-volt pumps are CE-approved.

- ####-A Audible Alarm (add "A" to pump part number)
- ####-P Programmable Ramp Option (SP200 Series)

## OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>15623</b>	Serial Cable, SP Pump-to-IBM 9-pin "D" Connector
<b>13685</b>	SP Pump-to-Pump "Daisy-Chain" Linking Cable, 7 ft.
<b>13962</b>	Foot Switch for SP200 Series Pumps

# Dual Rate Syringe Pump

*Two separate syringe pumps in one instrument*

## Features

- TTL, USB & RS-232 Communication
- Compatible with a wide variety of syringes
- Large touch screen display allows quick and easy setup

## Benefits

- Two separate syringe pumps in one instrument
- High or low-pressure operation
- Continuous fluid delivery

## Applications

- Continuous injection for long-term toxicology testing
- The injection of dyes, perfumes and flavoring in industrial applications
- Applications with liquids or viscous materials in micro-manufacturing
- Continuous injections of reactants into reactor vessels
- Simultaneous samplings from two sites



SPLG-G88PLUS

CE

Various applications can benefit from the. All three conditions come with the Gemini 88 Plus. No additional upgrades are required.

The intuitive Gemini 88 Plus graphical user interface controlled with a large 7" LCD color touch screen display allows quick and easy setup. The display run screen presents the user with all key dispensing parameters in real time. Syringe tables containing all major syringe manufacturers allow simple selection of any compatible syringe size. Audible Alarms, Adjustable Force and Screen Lock are all features that are available with a touch of the screen.

The SPLG-G88PLUS Dual Rate Syringe Pump opens whole new pumping possibilities. Three operating conditions named Independent, Reciprocating and Twin are available as standard features to satisfy various applications in life science.

## Connectivity & Communication

The Gemini 88 Plus comes standard with USB and RS-232 for PC communication and RS-485 for pump-to-pump communication. An entire suite of ASCII commands is available to control the pump remotely with a PC. The pump contains a foot switch input and digital input/output for each independent pumping channel.

## Independent Condition

Independent Condition allows the Gemini 88 Plus to operate as two separate syringe pumps named P1 and P2. Each syringe will operate independently with different syringe types, sizes, force, or target (volume or time, mode dependent).

## Reciprocating Condition

In Reciprocating Condition, both syringe channels move in opposite directions at the same rate using the same syringe size and type. When combined with a valve box, the Reciprocating Condition can provide the continuous fluidic delivery or a peristaltic pump with the accurate, pulseless, low flow rates provided by a syringe pump.

## Twin Condition

Twin Condition allows both syringes to operate in the same mode using the exact same syringe type, syringe size, force, target (volume or time) and flow rate settings. The pump also allows the user to combine both flows for higher speed and volume infusion applications

## SPLG-G88PLUS SPECIFICATIONS

TYPE	Microprocessor dual rate infuse/withdraw/continuous syringe pump
ACCURACY	±0.25%
REPRODUCIBILITY	±0.05%
SYRINGE TYPE	Glass, plastic and stainless steel
SYRINGE SIZE MINIMUM	0.5 µL (0.103 mm minimum inner diameter)
SYRINGE SIZE MAXIMUM	60 mL (32.573 mm maximum inner diameter)
FLOW RATE MINIMUM	1.02 pL/min (0.5 µL syringe, 0.103 mm inner diameter)
FLOW RATE MAXIMUM	106 mL/min (60 mL syringe, 32.573 mm diameter)
FLOW RATE DISPLAY	7 in. color display with touchscreen
CONNECTORS	USB - Type B RS-232 - 9-pin D-Sub connector RS-485 - IEEE-1394, 6 pos for pump-pump communication TTL Input/Output - Two 15-pin D-sub connectors, one for reach pump mechanism Foot switch - Two phonojack inputs, one for each pump mechanism
AVERAGE LINEAR FORCE	70 lb. (31.75 kg) at 100% force setting up to a flow rate of 90 mL/min using up to a 60 mL syringe with a 32.573 mm inner diameter. 50 lb. (22.6 kg) at 100% Force Setting for flow rates 90 mL/min to 106 mL/min using the same size syringe.
POWER SUPPLY	Input 100-240 VAC, 50-60 Hz, Output 30 V 1.66A 50 Watt
WEIGHT	21 lb. (9.09 kg)
DIMENSIONS (L x D x H)	11 x 15 x 8 in. (28 x 39 x 21 cm)
CLASSIFICATION	Class I
POLLUTION	Degree 1
INSTALLATION	Category II
REGULATORY CERTIFICATIONS	CE, ETL (UL & CSA), CB Scheme, EU RoHS, WEEE

## ORDERING INFORMATION

**SPLG-G88PLUS** Dual Rate Syringe Pump

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
Brazil: (013) 406-29703 • info@brazil.wpiinc.com • www.wpiinc.com

Germany: +49 (0)6031 1602171 • wpide@wpi-europe.com • www.wpi-europe.com  
China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net

# Programmable Syringe Pump

*The best value in syringe pumps on the market!*

## Features

- Cost-effective
- Program sequences without a computer
- Rich command set
- RS-232 and TTL control

## Benefits

- Automatic dispensing of small volumes
- Very precise flow rate control
- Hands-free operation with foot pedal

## Applications

- Drug administration
- Chemical applications with slow incorporation of fixed volumes of fluids



*Aladdin*

The Aladdin pump series arguably offers the best VALUE for any syringe pump on the market.

These pumps are available in single, dual, 4, 6, and 8 syringe capacities. Although it is one of the most cost-effective pumps available, the Aladdin pump series boasts a versatility that is unmatched in its price range with features typically found only on pumps costing more than twice as much.

The Aladdin series provides the capability to both inject and withdraw. These pumps can be used for manually triggered injections, but also have the capability to be programmed in multi-step pumping sequences without the requirement for a computer. Even so, this series also has both RS-232 and TTL inputs for external control and status feedback to a computer, if desired. The pumps can be daisy chained for multi pump applications involving push/pull protocols or multiple independent channels. The performance characteristics are admirable, considering the price.

For the budget minded lab looking for a versatile and reliable pump, the Aladdin series is an excellent option.

The Aladdin pump series will accept syringes from Becton Dickinson, Monoject, Terumo and Air-Tite.

## ORDERING INFORMATION

<b>AL-300</b>	Non-Programmable Syringe Pump, Infuse Only
<b>AL-1000</b>	Programmable Syringe Pump
<b>AL-1000HP</b>	Programmable Syringe Pump, High Pressure (100 lb.)
<b>AL-1000HP2</b>	Programmable Syringe Pump, High Pressure (200 lb.)
<b>AL-2000</b>	Two AL1000 Syringe Pumps <i>Includes CBL-Dual-3 Interconnecting Cable for push/pull or continuous pumping. Valves not included.</i>
<b>AL-4000</b>	Programmable Multiple (up to 4) Syringe Pump
<b>AL-6000</b>	Programmable Multiple (up to 6) Syringe Pump
<b>AL-6000H</b>	Programmable Multiple (up to 6) Syringe Pump
<b>AL-8000</b>	Programmable Multiple (up to 8) Syringe Pump

*Specify line voltage. When ordering 220V models, specify UK, Euro or Australian line cord.*

## OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>GN-PC7</b>	PC to Pump Cable, 7 ft.
<b>GN-PC25</b>	PC to Pump Cable, 25 ft.
<b>GN-NET7</b>	Pump-to-Pump Network Cable, 7 ft.
<b>GN-NET25</b>	Pump-to-Pump Network Cable, 25 ft.
<b>GN-TTL</b>	Pump-to-Pump Reciprocating Cable
<b>ADPT2</b>	Foot Switch

## ALADDIN SPECIFICATIONS

	AL-300	AL-1000	AL-1000 HP	AL-4000	AL-6000	AL-8000
SYRINGE SIZES	Plastic syringes up to 60 mL	Plastic syringes up to 60 mL and selected glass micro syringes from 0.5 to 500 µL.	Plastic syringes up to 60 mL and selected glass micro syringes from 0.5 to 500 µL.	Plastic syringes up to 60 mL	1-60 mL, or 0.5-5 µL micro syringes	Up to 100 mL
NUMBER OF SYRINGES	1	1	1	2 (may be different sizes)	6	8
MOTOR TYPE	Step Motor, 1/8 to 1/2 step modes	Step Motor, 1/8 to 1/2 step modes		Step Motor, 1/8 to 1/2 step modes	Step Motor, 1/8 to 1/2 step modes	Step Motor
STEPS PER REVOLUTIONS	400	400	200	200	400	200
STEPPING (max./min.)	0.21 µm to 0.850 µm	0.21 µm to 0.850 µm		0.425 µm to 1.7 µm, depending on motor speed	0.21 µm to 0.850 µm	0.132µm to 0.265 µm, depending on motor speed
MOTOR TO DRIVE SCREW RATIO	15/28	15/28	15/28	15/28	15/28	5/1
SPEED (max./min.)	3.774 cm/min. to 0.004 cm/hr	5.1 cm/min / 0.0042 cm/hr	18.370 cm/min / 0.008 cm/hr	18.370 cm/min / 0.008 cm/hr	5.1 cm/min / 0.0042 cm/hr	3.4917cm/min. / 0.0026cm/hr.
PUMPING RATES	1257 mL/hr with 60 mL syringe, to 0.73 µL/hr with 1 mL syringe	1699 mL/hr with 60 mL syringe, to 0.73 µL/hr with 1 mL syringe	6120 mL/hr with 60 mL syringe, to 1.459 µL/hr with 1 mL syringe	6120 mL/hr with 60 mL syringe, to 1.459 µL/hr with 1 mL syringe	1699 mL/hr with 60 mL syringe, to 0.73 µL/hr with 1 mL syringe	0.454 µL/hr. with a B-D 1cc syringe to 1163mL/hr. with a B-D 60 cc syringe
MAXIMUM FORCE	35 lb. at minimum speed, 18 lb. at maximum speed	35 lb. at min. speed, 18 lb. at maximum speed	100 lb. at minimum speed, 18 lb. at maximum speed	100 lb. at minimum speed, 15 lb. at maximum speed	35 lb. at min. speed, 18 lb. at max. speed	160 lbs. at min. speed, 30 lbs. at max. speed
PROGRAM PHASES	n/a	41	41	41	41	41
RS-232 PUMP NETWORK	n/a	100 pumps maximum	100 pumps maximum	100 pumps maximum	100 pumps maximum	100 pumps maximum
POWER SUPPLY	Wall adapter 12V DC @ 850 mA	Wall adapter 12V DC @ 850 mA	Wall adapter 12V DC @ 1000 mA	Wall adapter 12V DC @ 1000 mA	Wall adapter 12V DC @ 850 mA	Unregulated linear external wall adapter, country and power source specific
DIMENSIONS	22.9 x 14.6 x 11.4 cm (8.75 x 5.75 x 4.5 in.)	22.9 x 14.6 x 11.4 cm (8.75 x 5.75 x 4.5 in.)	22.9 x 14.6 x 11.4 cm (8.75 x 5.75 x 4.5 in.)	22.9 x 14.6 x 11.4 cm (8.75 x 5.75 x 4.5 in.)	22.9 x 14.6 x 11.4 cm (8.75 x 5.75 x 4.5 in.)	26x38.1x12.7cm (10.25 x15 x5 in.)
WEIGHT	1.6 kg (3.6 lb.)	1.6 kg (3.6 lb.)	1.6 kg (3.6 lb.)	1.6 kg (3.6 lb.)	1.6 kg (3.6 lb.)	10.125 lbs. (4.595kg)

WORLD PRECISION INSTRUMENTS

www.wpiinc.com

ORDER TOLL-FREE: (866) 606-1974 (U.S. only) • Tel: (941) 371-1003 • Fax: (941) 377-5428 • E-mail: sales@wpiinc.com

65

# Microfluidic Pumps

## Microfluidic Solenoid Pump



CX-KIMA

### Features

- Controlled by the iKima™ application for use with the iPod Touch (included) and iPhone
- Connects easily to Vena8 Endothelial+ biochip and other manufacturer's flow chambers
- Delivers pulses of fresh media to cells seeded inside the microchannels
- Flow rate: 15 - 35 mL/hr ± 4%
- Dead volume: < 300 µL

### Benefits

- Fits inside standard CO<sub>2</sub> incubators, maintaining temperature and humidity
- Recirculating long-term perfusion pump
- WiFi communication and control via iPod Touch, even with Kima pump incubated

### Applications

- Cell culture under shear stress/flow
- Biofilm studies
- Cell culture in biochips with adherent cells (HUVECs)
- Stem cells
- HepG2 cells (human liver cancer cell line)

**Kima** pump is a microfluidic pump designed to aid cell culture (e.g., epithelial and endothelial cells) under physiological conditions (shear flow) in various biochips and flow chambers, including Vena8 Endothelial+ biochips where it is possible to culture eight cell monolayers simultaneously over 24-48 hours. The solenoid pump produces a pulsatile flow to mimic the cardiovascular system, delivering a preset volume of liquid as programmed.

### ORDERING INFORMATION

<b>CX-KIMAKIT-IPOD</b>	Kima Kit for Cell Culture/Biofilm Studies <i>Includes Kima Pump, iPod Touch with Controller, 100 mL bottle &amp; cap, tubing kit, Vena Endothelial+ biochips (10-pack).</i>
<b>CX-KIMAKIT-PC</b>	Kima Kit for Cell Culture/Biofilm Studies <i>Includes Kima Pump, PC Controller, 100 mL bottle &amp; cap, tubing kit, Vena Endothelial+ biochips (10-pack).</i>
<b>CX-KIMA-I</b>	Kima Microfluidic Pump with iPod Touch & Controller
<b>CX-KIMA-P</b>	Kima Microfluidic Pump with PC Control Software and Controller
<b>CX-KIMA</b>	Kima Microfluidic Pump

## Precision Syringe Pump



CX-EXIGO

### Features

- Virtually pulse-free low volume delivery
- Fast response time
- Flow sensor feedback to pump
- Touch control

### Benefits

- Pulse-free flow
- Precise multichannel mixing
- Excellent long-term flow stability
- Rapid flow change (ms range)
- WiFi communication and control via iPad mini or USB communication with LabVIEW based PC control software
- Use standard tubing for connection to any microfluidic biochip

### Applications

- Microfluidics and nanofluidics
- Droplet generation and manipulation
- Laminar and multilaminar flow studies
- Cell-based shear flow studies
- Cell and particle manipulation studies

**ExiGo** is a precision syringe pump based on the 5-phase stepper motor drive that has more microsteps per revolution of the lead screw vs. standard syringe pumps on the market. ExiGo has 250,000 microsteps/revolution and even at low rotational speed/low flow rates, it has a very low pulsation and high accuracy. When coupled with the flow sensor and active PID feedback; this results in very fast response times for changing flow rates. A standard syringe pump typically has a smaller number of microsteps and so usually the only way a standard "microfluidic" pump can achieve pulse-free flow control is to use small syringes; e.g. 0.5 µL; 1 µL; 5 µL; etc., to achieve non-pulsatile stable flow rates in the nanoliter/minute range. By comparison, the ExiGo pump with the flow sensor can use a standard 250 µL glass syringe to produce stable non-pulsatile flow rates of 10nL/min – 1 mL/min; or a 5 mL plastic syringe to produce stable non-pulsatile flow rates of 100nL/min – 20 mL/min. ExiGo can be used in conjunction with expandable (flexible) element and fluidic resistance in order to dampen any pulsation occurring during the stepper motor operation. As it employs active feedback, the response time of the pumps still remains fast.

### ORDERING INFORMATION

<b>CX-EXIGO-IM-7</b>	Exigo with iPad, Manifold, 7µL/min FS
<b>CX-EXIGO-PM-7</b>	Exigo with LabView, Manifold, 7µL/min FS
<b>CX-EXIGO-I-7</b>	Exigo with iPad, 7µL/min FS
<b>CX-EXIGO-P-7</b>	Exigo with LabView, 7µL/min FS
<b>CX-EXIGO-IM</b>	Exigo with iPad, Manifold
<b>CX-EXIGO-PM</b>	Exigo with LabView, Manifold
<b>CX-EXIGO-I</b>	Exigo with iPad, SmartFlo app
<b>CX-EXIGO-P</b>	Exigo with LabView
<b>CX-EXIGO-M</b>	Exigo with Manifold, Tubing Kit
<b>CX-EXIGO</b>	Exigo with Tubing Kit

### WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
Brazil: (013) 406-29703 • info@brazil.wpiinc.com • www.wpiinc.com

Germany: +49 (0)6031 1602171 • wpide@wpi-europe.com • www.wpi-europe.com  
China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net

# Microfluidic Pumps

## Microfluidic Pressure Pump



CX-UNIGO

## Microfluidic syringe pump for shear flow studies



CX-MIRUS

### Features

- Precise flow control with active feedback via plug-and-play flow sensor (required add-on)
- Flow rate: 1  $\mu\text{L}/\text{min}$  – 1  $\text{mL}/\text{min}$ ; unidirectional (push)
- iPad mini or PC (LabVIEW, Matlab, Python, etc.) control which can control/program up to four pump modules independently
- External compressor (required)

### Benefits

- 2 modes of pumping: manual flow rate set or preprogrammed flow rate operation
- Preprogrammed mode includes: constant, ramp, step, sine functions
- Side port connections to dock up to four pumps (combination of UniGo and ExiGo) into one setup
- Use standard tubing for connection to any microfluidic biochip.

### Applications

- Microfluidics where accurate and stable flow rate delivery is required. The pressure pump component is based on controlled air injection.

The **UniGo™** Microfluidic Pump is a precision, microfluidic, single-channel pressure pump for a variety of microfluidic applications, where accurate and stable flow rate delivery is required. The pressure pump component is based on controlled air injection. The UniGo pump requires a plug-and-play flow sensor for active feedback and increased flow control. **SmartFlo** application executed on the iPad Mini or LabView-based interface communicates with up to four UniGo™ Microfluidic Pumps racked together, allowing simultaneous control and independent programming of each pump's flow profile. Uniquely, the UniGo™ may be docked together with the **ExiGo™** microfluidic syringe pump, combining the best features of both instruments in one microfluidic set-up. *Note: it is necessary to purchase a flow sensor with the UniGo™ pump.*

### ORDERING INFORMATION

<b>CX-UNIGO-ECI-80</b>	UniGo with iPad, External Compressor, 80 $\mu\text{L}/\text{min}$ Flow Sensor
<b>CX-UNIGO-ECP-80</b>	UniGo with LabView, External Compressor, 80 $\mu\text{L}/\text{min}$ Flow Sensor
<b>CX-UNIGO-EC-80</b>	UniGoExternal Compressor, 80 $\mu\text{L}/\text{min}$ Flow Sensor
<b>CX-UNIGO-I-80</b>	UniGo with iPad, No Compressor, 80 $\mu\text{L}/\text{min}$ Flow Sensor
<b>CX-UNIGO-P-80</b>	UniGo with LabView, No Compressor, 80 $\mu\text{L}/\text{min}$ Flow Sensor
<b>CX-UNIGO-80</b>	UniGo, No Compressor, 80 $\mu\text{L}/\text{min}$ Flow Sensor

Five flow cells are available for your UniGo system. Your choice of one flow sensor is included with your UniGo system. Contact WPI at [wpi@wpiinc.com](mailto:wpi@wpiinc.com) for more information.

### Features

- Unlimited volume
- Low per-channel cost
- Flow rate: 100nL/min - 10 mL/min  $\pm$  1%
- Dead volume:  $\sim$ 600  $\mu\text{L}$

### Benefits

- Higher throughput enabling eight assays in parallel
- Flow damper to decrease syringe pump pulses
- PC controlled via VenaFluxAssay software

### Applications

- Microfluidic applications
- Single cell analysis
- Microfluidic syringe pump for cell analysis under shear flow in biochips
- Suitable for cell samples and whole blood samples

**Mirus** is a precision syringe pump, which uses the combination of an expandable (flexible) element and fluidic resistance in order to dampen the pulsation of the syringe pump stepping. The Mirus has a very stable flow profile. Additionally, Mirus is equipped with 3-way valve allowing automatic recharging and washout of syringe. Mirus is also provided with an 8-way flow splitter, allowing multiple executions (eight parallel experiments) simultaneously.

### ORDERING INFORMATION

<b>CX-MIRUS-PRO</b>	Microfluidic Syringe Pump Includes tubing kit, VenaFlux Assay Software PC control and MultiFlow8
<b>CX-MIRUS</b>	Microfluidic Syringe Pump Includes tubing kit, VenaFlux Assay Software PC control
<b>CX-MIRUS-MULTI8</b>	MultiFlow8 Attachment for Mirus Evo Nanopump

# Microfluidic Pumps and Accessories

## 4-channel Microfluidic Pump



CX-4U

### Features

- Precise flow control with active feedback via plug-and-play flow sensor (required add-on)
- Flow rate: 1  $\mu$ L/min – 1 mL/min; unidirectional (push)
- iPad mini or PC (LabVIEW, Matlab, Python, etc.) control which can control/program up to four pump modules independently
- External compressor (required)
- Working pressure- Max 10 bars/145 psi

### Benefits

- Two modes of pumping: manual flow rate set or preprogrammed flow rate operation
- Preprogrammed mode includes: constant, ramp, step, sine functions
- Side port connections to dock up to four pumps (combination of UniGo and ExiGo) into one setup
- Use standard tubing for connection to any microfluidic biochip

### Applications

- Microfluidics where accurate and stable flow rate delivery is required. The pressure pump component is based on controlled air injection.

**4U** is a 4-channel compact, precise, economical microfluidic pressure pump for a variety of microfluidic applications, where accurate and stable flow rate delivery is required. **SmartFlo** application executed on LabView-based interface communicates with the four separate channels of the **4U** pump, allowing simultaneous control and independent programming of each pump's flow profile. The **4U** pump provides you with the advantage of using four independent flow rates at the same time.

### ORDERING INFORMATION

<b>CX-4UECIPAD-4X80FS</b>	4-channel Microfluidic Pump, External Compressor, iPad, 80 $\mu$ L/min flow sensors (4)
<b>CX-4UECPC-4X80FS</b>	4-channel Microfluidic Pump, External Compressor, PC, 80 $\mu$ L/min flow sensors (4)
<b>CX-4UNCIPAD-4X80FS</b>	4-channel Microfluidic Pump, iPad, 80 $\mu$ L/min Flow Sensors (4)
<b>CX-4UNCPC-4X80FS</b>	4-channel Microfluidic Pump, PC, 80 $\mu$ L/min Flow Sensors (4)
<b>CX-EC</b>	External Compressor

## ExiGo™ Pump Manifold

### Features

- Allows you to direct the fluid to three (3) independent ports
- Plug and play connection to the ExiGo™ pump
- Can be programmed to automatically switch between fluidic channels using SmartFlo PC software



CX-MF

### Benefits

- Refill syringes automatically

### Applications

- Microfluidics and Nanofluidics
- Droplet generation and manipulation
- Laminar and multilaminar flow studies
- Cell-based shear flow studies
- Cell and particle manipulation studies

The Manifold is a specialized microfluidic channel selector which allows the **ExiGo** pump to direct fluid to one of three microfluidic channels at a time. Accurate flow switching and low dead volume provide exceptional performance.

### ORDERING INFORMATION

<b>CX-MF</b>	ExiGo Pump Manifold
--------------	---------------------

## Microfluidic Biochips for a Variety of Applications

We offer a range of biochips for many applications. All biochips mimic human capillaries by working in tandem with our different instrumentation platforms.

Our new DropChips are either hydrophilic or hydrophobic based. DropChips comes with 5 different channel sizes from 30  $\mu$ m to 80  $\mu$ m. Please refer to the chart on page 69 for information on specific biochips.

### ORDERING INFORMATION

#### Optional Accessories/Replacement Parts

<b>CX-018</b>	Vena8 with Glass Coverslip
<b>CX-023</b>	Biochip-CONN, Biochip SNGL IN
<b>CX-024</b>	Biochip-CONN, Biochip SNGL OUT
<b>CX-027</b>	Tygon Tubing for Biochip Conn
<b>CX-032</b>	Kima-Connect Tubing Set
<b>CX-045</b>	Kima Pump Cell Culture Bottle
<b>CX-067</b>	Flow Sensor Exigo 1 mL/min
<b>CX-073</b>	Dropchip, Hydrophobic, 40 $\mu$ m x 40 $\mu$ m
<b>CX-084</b>	MF8-Connect Biochip1 Inlet CBL
<b>CX-087</b>	MULTIFLOW8 Nozzles, Standard



Vena8 Fluoro+



Vena8 Endothelial



VenaT4



Vena8 Glass Cover slip



Vena Delta Y1



Vena Delta Y2

### WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
 Brazil: (013) 406-29703 • info@brazil.wpiinc.com • www.wpiinc.com

Germany: +49 (0)6031 1602171 • wpide@wpi-europe.com • www.wpi-europe.com  
 China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net

# Microfluidic Biochips

## DROPCHIP (HYDROPHILIC/HYDROPHOBIC) BIOCHIPS ORDERING INFORMATION

Chip Size (µm)	30 x 30	40 x 40	50 x 50	60 x 60	70 x 70	80 x 80
Junction Width (cm)	0.003	0.004	0.005	0.006	0.007	0.008
Height (cm)	0.003	0.004	0.005	0.006	0.007	0.008
Length (cm)	2.8	2.8	2.8	2.8	2.8	2.8
Substrate Thickness	0.5	0.5	0.5	0.5	0.5	0.5
Hydrophilic (Pack of 2) Hydrophobic (Pack of 2)	<b>CX-026</b> <b>CX-034</b>	<b>CX-028</b> <b>CX-073</b>	<b>CX-029</b> <b>CX-035</b>	<b>CX-030</b> <b>CX-036</b>	<b>CX-031</b> <b>CX-037</b>	<b>CX-033</b> <b>CX-038</b>
Hydrophilic (Pack of 5) Hydrophobic (Pack of 5)	<b>CX-039</b> <b>CX-046</b>	<b>CX-040</b> <b>CX-047</b>	<b>CX-041</b> <b>CX-048</b>	<b>CX-042</b> <b>CX-049</b>	<b>CX-043</b> <b>CX-050</b>	<b>CX-044</b> <b>CX-051</b>
Hydrophilic (Pack of 10) Hydrophobic (Pack of 10)	<b>CX-052</b> <b>CX-058</b>	<b>CX-053</b> <b>CX-059</b>	<b>CX-054</b> <b>CX-060</b>	<b>CX-055</b> <b>CX-061</b>	<b>CX-056</b> <b>CX-062</b>	<b>CX-057</b> <b>CX-063</b>

## STANDARD BIOCHIPS ORDERING INFORMATION

	Vena8 Fluoro+	Vena8 Endothelial+	VenaT4	Vena8 Glass Coverslip		Vena Delta	
				Low Flow Rates	High Flow Rates	Y1	Y2
Channel width, b (cm)	0.04	0.08	0.08	0.16	0.08	0.008	0.008
Channel height, h (cm)	0.01	0.012	0.01	0.016	0.008	0.012	0.012
Channel length, L (cm)	2.8	2.8	2.8	2.8	2.8	2.8	2.8
Channel volume (cm <sup>3</sup> )	0.00112	0.00269	0.00224	0.00717	0.00179	0.00269	0.00269
Channel volume (µL)	1.12	2.69	2.24	7.17	1.79	2.69	2.69
Thickness of bottom substrate (mm)	0.17	0.5		0.17	0.17	0.5	0.5
# of channels / chip	8	8	4	8	8	4	4
# of assays / pack	40/80	40/80	40	80	80	40	40
Pack of 5	<b>CX-002</b>	<b>CX-004</b>					
Pack of 10	<b>CX-001</b>	<b>CX-003</b>	<b>CX-005</b>	<b>CX-009</b>	<b>CX-010</b>	<b>CX-007</b>	<b>CX-008</b>

## SPECIFICATIONS

Brightfield, phase contrast, immunostaining	✓	✓	✓	✓	✓	✓	✓
Confocal microscopy	✓			✓	✓		
<b>Cell types:</b> whole blood (human, animal); PRP; platelets; PBMC; monocytes; T-cells (primary and cell lines); eosinophils; neutrophils etc.	✓	✓	✓	✓	✓	✓	✓
<b>Cell types:</b> adherent cells e.g. endothelial – HUVECs; HMVECs etc.; HepG2; stem cells; muscle cells etc.	✓	✓	✓	✓	✓	✓	✓
<b>Protein coatings:</b> collagen, fibronectin, fibrinogen, vWF, VCAM, ICAM, selectins, MadCAM etc.	✓	✓	✓	✓	✓	✓	✓

## APPLICATIONS

Platelet adhesion, aggregation and thrombi formation; leukocyte rolling, adhesion and migration; thrombosis; immunology (inflammation); infectious diseases (e.g. malaria); sickle cell disease; respiratory (asthma and COPD)	✓			✓	✓		
Cell adhesion and culture under perfusion / shear flow; leukocyte cell-cell rolling, adhesion and migration; oncology (melanoma, breast cancer etc.); cardiovascular (atherosclerosis, drug eluting stents); immunology (inflammation); respiratory (asthma and COPD)		✓					
Biofilm assays, microbe seeding and culturing; biochips with glass coverslips (attached / not attached; treated / non-treated); biochips for the attachment of coupons for biofilm studies				✓	✓		
Chemotaxis, transmigration and invasion assays; 2D and 3D cell culture; mimicking tumour microenvironment with gels (ECM gel, hydrogel, matrigel, collagen gel)			✓				

All biochips are:

- Disposable plastic; some with glass cover slips.
- Require no assembly; unlike many standard perfusion chambers / flow chambers.

- Require no Luer lock connections which increase dead volume. Cellix's biochips have a unique plug and play connection with tubing connections which are autoclaveable and reusable.

# Ergonomic Pipetter Design

*Making repetitive procedures more efficient and comfortable*

## Features

- Lightweight and comfortable ergonomic design
- Easy calibration using provided tool
- Easy for cleaning and parts replacement
- CE and ISO13485 Certified



REG100

## Benefits

- Less stress on your hand when you are performing repetitive operations
- Save money by ordering sets of any 5, 6 or 7 pipettors
- Made from biologically inactive and chemical inert polymers

## Applications

- Routine laboratory use for accurate and affordable pipetting of liquids and solutions

Regal Pipettes use the air displacement method, where displaced air volume by the motion of the piston is equal to the liquid volume drawn into the pipette. These new pipettes are highly accurate. Within ten complete revolutions of the dial, you can set the minimum and maximum volumes. For ease of use, the dial is fixed to the plunger. Since the light plunger action reduces fatigue, results are more precise.



### ORDERING INFORMATION (SINGLES)

Model	Volume Range $\mu\text{L}$	Increment $\mu\text{L}$	Nominal Volume $\mu\text{L}$	Tolerance %	Repeatability %
REG2	0.2 ~ 2	0.01	0.2	$\pm 12.0$	$\leq 6.00$
			0.5	$\pm 5.0$	$\leq 2.50$
			2	$\pm 2.0$	$\leq 0.70$
REG10	1 ~ 10	0.1	1	$\pm 3.0$	$\leq 1.50$
			5	$\pm 1.5$	$\leq 0.60$
			10	$\pm 1.0$	$\leq 0.40$
REG20	2 ~ 20	0.1	2	$\pm 3.0$	$\leq 1.50$
			10	$\pm 1.0$	$\leq 0.50$
			20	$\pm 1.0$	$\leq 0.30$
REG50	5 ~ 50	0.5	5	$\pm 2.0$	$\leq 1.50$
			20	$\pm 1.2$	$\leq 0.40$
			50	$\pm 1.0$	$\leq 0.20$
REG100	10 ~ 100	1	10	$\pm 2.0$	$\leq 0.50$
			50	$\pm 0.8$	$\leq 0.30$
			100	$\pm 0.8$	$\leq 0.15$
REG200	20 ~ 200	1	20	$\pm 2.0$	$\leq 0.50$
			100	$\pm 0.8$	$\leq 0.30$
			200	$\pm 0.8$	$\leq 0.15$
REG1000	100 ~ 1000	5	100	$\pm 1.5$	$\leq 0.30$
			500	$\pm 0.8$	$\leq 0.30$
			1000	$\pm 0.8$	$\leq 0.15$
REG5K	1000 ~ 5000	50	1000	$\pm 1.0$	$\leq 0.50$
			2000	$\pm 0.7$	$\leq 0.25$
			5000	$\pm 0.7$	$\leq 0.15$
REG10K	1000 ~ 10000	100	1 mL	$\pm 3.0$	$\leq 0.30$
			5 mL	$\pm 0.7$	$\leq 0.20$
			10 mL	$\pm 0.7$	$\leq 0.15$

### ORDERING INFORMATION (SETS)

<b>REGS5</b>	Regal Pipettors (set of any 5) & Stand
<b>REGS6</b>	Regal Pipettors (set of any 6) & Stand
<b>REGS7</b>	Regal Pipettors (set of any 7) & Stand
<b>504591</b>	Stand for Regal Pipettors (holds 8)



# Universal Pipette Tips


Ultra-clear, certified RNase/DNase-free

# Mini Vacuum Pump

Small, reliable, durable and accurate

WPI's Universal Pipette Tips are for use with Regal and most other pipettors, including Gilson, Oxford Benchmate, Socorex, and SealPette.

\* Tips 500197 and 500198 fit Regal, Eppendorf, and BioHit pipettors.



500192



801566

Same as leading brands—  
**at about half the price!**

## ORDERING INFORMATION UNIVERSAL FILTER TIPS (STERILE)

Tip Volume	For Pipetter	Rack	Part No.
0.1 - 10 µL	REG2	960 (10 racks of 96)	<b>500199</b>
	REG10		
	REG20		
10 - 200 µL	REG20	960 (10 racks of 96)	<b>500200</b>
	REG50		
	REG100		
	REG200		

## ORDERING INFORMATION UNIVERSAL TIPS

Tip Volume	For Pipetter	Bulk	Part No.	Rack	Part No.
0.1 - 10 µL	REG2	Bag of 1000	<b>500191</b>	960	<b>500192</b>
	REG10			(10 racks of 96)	
	REG20				
5 - 200 µL	REG20	Bag of 1000	<b>500193</b>	960	<b>500194</b>
	REG50			(10 racks of 96)	
	REG100				
	REG200				
100-1000 µL	REG1K	Bag of 1000	<b>500195</b>	1000 (10 racks of 100)	<b>500196</b>
500 - 5000 µL	REG5K	Bag of 250	<b>500197 *</b>	500 (10 racks of 50)	<b>500198 *</b>

## Features

- Durable aluminum exterior
- Minimal vibration
- Low noise
- Extremely long life time

## Benefits

- Quiet operation
- Compact (18 x 7 x 7cm) unit takes little of your bench space
- Oil free, maintenance free

## Applications

- Commonly sold with Muscle Physiology setups
- Excellent accessory for use with WPI's PV830 Pneumatic PicoPump with vacuum
- Ideal for any application requiring a small, reliable pump that provides vacuum pressure up to 250 mbar

This miniature vacuum pump is durable and accurate. The industrial-strength aluminum exterior, neoprene diaphragm and neoprene/silicone valves ensure this pump will stand up to daily use.

## MINI VAC SPECIFICATIONS

POWER SOURCE	230 (50 Hz)	120 (60 Hz)
FREE FLOW	4.0 L/min.	3.0 L/min.
AT -100 MBAR	2.0 L/min.	1.5 L/min.
MOTOR TYPE	Vibrating	
POWER	4.0 W	
MAXIMUM PRESSURE	-	
MAXIMUM VACUUM	-250 mbar	
PUMP HEAD CONSTRUCTION	Aluminum	
DIAPHRAGM	CR-neoprene	
VALVES	CR-neoprene/FPM (Viton)/Silicone	
DIMENSIONS	185 x 72 x 72 mm	
WEIGHT	850 g	

## INFORMATION ORDERING

<b>801566</b>	Mini Vacuum Pump, 110 V, US power plug
<b>801963</b>	Mini Vacuum Pump, 220 V, EURO mains plug

# Luer Valve Assortment Kit

*Build your own liquid flow experiment*

## Features

- Over 300 assorted parts
- Luer fittings for quick and easy connect and disconnect

## Benefits

- Sold individually or in kits

## Applications

- Liquid flow experimental setups

A useful kit (right) for building your own liquid flow experiment. It provides the means to start, stop, add, divide and control a flow of liquid or gas. Included in the kit are **over 300 assorted parts** such as one-way and three-way stopcocks, manifolds, Y-connectors, injection sites,

male and female luer caps, check valves, syringe-activated check valves, slide clamps, roller clamps, and pinch clamps. All (except clamps) have a luer fitting for quick and easy connecting and disconnecting. Includes assorted luer fittings for use with flexible tubing.

The **14011** component parts are also sold individually. See page 220.



## ORDERING INFORMATION

**14011** Luer Valve Assortment Kit

# Luer-to-Tubing Coupler Assortment Kit

*Quick connects in nylon and polypropylene*

## Features

- Over 250 assorted parts in each kit
- Valves are polycarbonate, and the valve handles are polyethylene. Do not autoclave those parts.

## Benefits

- Polypropylene parts (**504954**) can be autoclaved repeatedly at 121°C/15PSI, 15 min. cycle
- Polypropylene fittings are chemically inert and resistant to most organic and inorganic solvents
- Nylon fittings are strong and can be bonded with adhesive.

## Applications

- Liquid flow experimental setups



504954



504955

Assemble quick-disconnect luer fittings for use with flexible tubing with internal diameters of 1/16", 3/32" and 1/8". A variety of quick-disconnect connectors can be quickly made for connecting small diameter flexible tubing; 3-way connections can be made with the use of the 3-way luer tee; luer plugs, tees, connectors, bulk-head mounts, color coding rings, locking nuts, male and female luers—are all included to enhance the versatility of this kit. The kit has over 250 assorted parts and is offered in two different types of materials. Nylon parts are not autoclavable.

## ORDERING INFORMATION

**504954** Luer-to-Tubing Coupler Assortment Kit (Polypropylene)


**504955** Luer-to-Tubing Coupler Assortment Kit (Nylon)

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
 Brazil: (013) 406-29703 • info@brazil.wpiinc.com • www.wpiinc.com

Germany: +49 (0)6031 1602171 • wvide@wpi-europe.com • www.wpi-europe.com  
 China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net

# Electrophysiology



## Over 50 years of experience with researchers

From our early beginnings in 1967 working with Yale researchers, electrophysiology has been at the heart of our business. WPI amplifiers, stimulators and isolators are designed with quality components so you get a reliable, low-noise signal every time. Our time-tested designs give you affordable solutions for electrophysiology equipment, electrodes, data acquisition and accessories for applications like:

- Intracellular/Extracellular Recording
- Voltage Clamp for Ussing
- Stimulation and Isolation
- Optogenetics
- Digital Filtering

# Bioamplifiers for Electrophysiology

*A family of very low noise battery-operated amplifiers*

## Features

- Battery powered to eliminate line noise
- High pass and low pass filtering
- Single ended or differential operation
- DC/AC amplification
- Variable output positioning
- Constructed of high quality components to ensure minimal intrinsic (shot) noise
- Portable
- Rack mountable

## Benefits

- Very low internal noise
- Ultra quiet DC power supply — no AC required
- Intrinsic low susceptibility to ground loops
- Small footprint
- Cost effective
- Electrostatic Discharge Protection!

## Applications

- Amplifying biopotentials from metal electrodes
- Brain slice field stimulation
- EAG (Electroantennogram)
- ERG (Electroretinogram)

WPI's **DAM** series amplifiers are well known as a standard of the industry for amplification of extracellular potentials. These battery powered bio-amplifiers are designed with a compact chassis profile that enables you to locate the unit closer to the preparation and thereby minimize long lead lengths which contribute to noise. Each amplifier is equipped with selectable high and low filters, and a position control to offset galvanic potentials which may develop during recording.

DAM series amplifiers can be used as stand-alone units on any tabletop or use optional clamp-mounting hardware to locate them conveniently within the work area. Alternatively, a pair of amplifiers can be mounted into a standard equipment rack with a rack mount kit (3484). A variety of hook up accessories are available to configure your application.

### DAM80 Overview

**DAM80**, an AC amplifier only, features a very low noise active headstage probe which can be mounted in micromanipulators for up-close cortical recording, for extracellular recording from high impedance glass or metal microelectrodes. The unit also provides a gated current for tissue marking. Microelectrode holder **MEH7W-XX** (sold separately) is recommended for glass microelectrodes. The DAM80 is perfect for gated or manual current generation for histological marking, iontophoresis or cell stimulation. It includes a very low noise remote active headstage that is useful for very high impedance amplification utilizing glass or metal electrodes.



CE

**DAM80**— With low-noise headstage **DAM80P**

Included with the DAM80 is a **Startup Kit** containing the following accessories needed for basic metal electrode electrophysiology research:

<b>CBL102</b>	(2) Cable, BNC-to-3.5 mm plug, 6 ft (2m)
<b>5469</b>	(2) Adapter, mini-banana to 0.031 socket
<b>13388</b>	(2) Adapter, mini-banana to 2mm socket
<b>3294</b>	Cable, ground clip to wire, 3 ft
<b>2033</b>	Mini-banana plug, black
<b>2034</b>	Mini-banana plug, red
<b>2035</b>	(2) Mini-banana plug solderable turrent
<b>EP1</b>	Ag/AgCl pellet (70 mm wire) 1 mm diam x 2.5 mm long
<b>M3301EH</b>	(2) Electrode Holder, 14cm
<b>5470</b>	0.031-inch jack on 12-inch wire (package of 4)

### FEATURE COMPARISON

	DAM50	DAM80
INPUT MODE	AC/DC	AC
INPUT CONFIGURATION	differential/single ended	differential
GAIN RANGE	100-10K (AC), 10-1K (DC)	100-10K (AC)
HIGH/LOW FILTERS	yes	yes
OFFSET POSITION CONTROL	yes	yes
CURRENT GENERATOR	No	Yes
REMOTE ACTIVE HEADSTAGE	No	Yes
OUTPUT CONNECTION	BNC	3.5 mm mini phone
STANDARD INPUT CONNECTION*	terminated wire	mini banana
POWER SUPPLY	(2) 9V alkaline batteries	(2) 9V alkaline batteries

\*See optional accessories for additional alternatives

### WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com Germany: +49 (0)6031 1602171 • wpide@wpi-europe.com • www.wpi-europe.com  
Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • www.wpiinc.com China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net

## DAM50 References

Liu, Y., Wang, Y., Zhu, G., Sun, J., Bi, X., & Baudry, M. (2016). A calpain-2 selective inhibitor enhances learning & memory by prolonging ERK activation. *Neuropharmacology*, 105, 471–477. <http://doi.org/10.1016/j.neuropharm.2016.02.022>

Ztaou, S., Maurice, N., Camon, J., Guiraudie-Capraz, G., Kerkerian-Le Goff, L., Beurrier, C., ... Amalric, M. (2016). Involvement of Striatal Cholinergic Interneurons and M1 and M4 Muscarinic Receptors in Motor Symptoms of Parkinson's Disease. *Journal of Neuroscience*, 36(35).

Kentish, S. S. J., Frisby, C. L., Kritas, S., Li, H., Hatzinikolas, G., O'Donnell, T. A., ... Ahern, G. (2015). TRPV1 Channels and Gastric Vagal Afferent Signalling in Lean and High Fat Diet Induced Obese Mice. *PLoS One*, 10(8), e0135892. <http://doi.org/10.1371/journal.pone.0135892>

## DAM80 References

Donnelly, W. T., Bartlett, D., & Leiter, J. C. (2016). Serotonin in the solitary tract nucleus shortens the laryngeal chemoreflex in anaesthetized neonatal rats. *Experimental Physiology*, 101(7), 946–961. <http://doi.org/10.1113/EP085716>

Feng, B., Joyce, S. C., & Gebhart, G. F. (2016). Optogenetic activation of mechanically insensitive afferents in mouse colorectum reveals chemosensitivity. *American Journal of Physiology - Gastrointestinal and Liver Physiology*, 310(10).


Kaldenbach, F., Bleckmann, H., & Kohl, T. (2016). Responses of infrared-sensitive tectal units of the pit viper *Crotalus atrox* to moving objects. *Journal of Comparative Physiology A*, 202(6), 389–398. <http://doi.org/10.1007/s00359-016-1076-1>

Mustafina, A. N., Koroleva, K. S., Giniatullin, R. A., & Sitdikova, G. F. (2016). Acid Sensitive Ion Channels as Target of Hydrogen Sulfide in Rat Trigeminal Neurons. *BioNanoScience*, 1–3. <http://doi.org/10.1007/s12668-016-0237-6>

Orton, L. D., Pappasavvas, C. A., & Rees, A. (2016). Commissural Gain Control Enhances the Midbrain Representation of Sound Location. *The Journal of Neuroscience*, 36(16), 4470–81. <http://doi.org/10.1523/JNEUROSCI.3012-15.2016>

## DAM SERIES SPECIFICATIONS

INPUT IMPEDANCE	10 <sup>12</sup> Ω
INPUT LEAKAGE CURRENT	50 pA (typical)
MAX. DC DIFFERENTIAL SIGNAL	± 2.5 V (DAM 50)
GAIN	AC: 100x, 1000x, 10000x DC: 10x, 100x, 1000x (DAM50)
COMMON MODE REJECTION RATIO	100 dB @ 50/60 Hz
INPUT CAPACITANCE	20 pF
AC MODE NOISE	0.4 μV RMS (2 μV p-p) 0.1-100 Hz
AC MODE NOISE	2.6 μV RMS (10 μV p-p) 1 Hz-10 kHz
DC MODE NOISE (DAM50)	7.5 μV RMS (30 μV p-p) 3-10 kHz
BANDWIDTH FILTER SETTINGS	
AC Mode	Low frequency, 0.1, 1, 10, 300 Hz
AC Mode (DAM80)	High frequency, 0.1, 1, 3, 10 kHz
DC Mode (DAM50)	High frequency, 0.1, 1, 3, 10 kHz
OUTPUT CONNECTORS	BNC on DAM50; 3.5 mm MiniPhone connector on DAM80
OUTPUT VOLTAGE SWING	±8 V
OUTPUT IMPEDANCE	470 Ω
BATTERY TEST	Audible tone
CALIBRATOR SIGNAL	10 Hz square wave
POSITION	Approximately 250 mV
CURRENT SOURCE	
DAM80: DC Generator	0 to ±50 μA, variable
EXTERNAL COMMAND	Input Voltage ±10 V commands
AC OR DC CURRENT WAVEFORM	±50 μA max. amplitude @ 200 KΩ
BATTERIES	2 x 9V alkaline (included)
DIMENSIONS	
DAM50	8 x 4 x 1.75 in. (20.3 x 10.2 x 4.4 cm)
DAM80	7 x 4 x 1.75 in. (17.8 x 10.2 x 4.4 cm)
SHIPPING WEIGHT	3.5 lb (1.6 kg)



DAM50—Basic Amplifier

(Optional 5447 Electrode Adapter)

## ORDERING INFORMATION

**SYS-DAM50** Bio-amplifier with shielded electrode cable


**SYS-DAM80** Bio-amplifier with active probe (DAM80P)

## OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>DAM80P</b>	Replacement Probe
<b>3072</b>	6 Replacement Modular Cables (DAM50)
<b>3517</b>	2 Optional Shielded Modular Cables (DAM50)
<b>CBL102</b>	3.5 mm Phone plug-to-BNC Cable
<b>2851</b>	BNC-to-BNC Cable
<b>2033</b>	Black Insulated Mini-Banana Plug
<b>2034</b>	Red Insulated Mini-Banana Plug
<b>2035</b>	Uninsulated Mini-Banana Plug
<b>3484</b>	Rack Mount Kit (for 1 or 2 DAM preamps)
<b>3485</b>	Ringstand Mounting Kit
<b>5447</b>	Electrode Adapter (DAM50)
<b>5469</b>	Metal Microelectrode Adapter for DAM80 (mini-banana plug to 0.031 in. (0.79 mm) socket)
<b>5489</b>	Adapter for Metal Microelectrode (DAM50)
<b>13388</b>	Adapter, mini-banana plug to 2mm socket
<b>5371</b>	Cable, Low Noise (2 mm pin to 2 mm pin)
<b>3578</b>	Adapter Cable for Ag/AgCl pellets (2 mm pin)
<b>300102</b>	Metal Electrode Extension, 4-in. (2 mm to 0.031 in.)
<b>300647</b>	Shielded Electrode Cable Assembly
<b>MEH7W-XX</b>	Microelectrode Holder- 1.0, 1.5 or 2.0 mm OD

See *Cables and Connectors*, page 96.

See *Metal Microelectrodes*, page 94.



Optional probe 5489 (non-active) for use with DAM50 also includes microelectrode adapter 5469.

# Which Amplifier is Right for You?

## AMPLIFIER COMPARISON CHART

Amplifier	AC/DC	Differential	Active Headstage	Stimulation	Isolated	Multi-channel	Battery Powered	Connectors
<b>Intracellular Bioamplifiers</b>								
<b>Electro 705</b>	DC		◆				◆	2 mm pin
<b>Duo773</b>	DC	◆	◆	◆		2		2 mm pin
<b>Extracellular Bioamplifiers</b>								
<b>ISODAM8A</b>	DC	◆	opt		◆	4 - 8		Mini Banana or 8-pin DIN
<b>ISO80</b>	AC	◆	◆	◆	◆		◆	Mini Banana
<b>DAM50</b>	AC/DC	◆					◆	RJ-11
<b>DAM80</b>	AC	◆	◆	◆			◆	Mini Banana
<b>Transducer Amplifiers</b>								
<b>BRIDGE8</b>	DC	◆				4 - 8		8-pin DIN WPI transducers
<b>TBM4M</b>	DC	◆				4		8-pin DIN WPI transducers

## WPI's Low-Noise Amplifiers Outperform Cheap Imitations

An amplifier is an electronic device that magnifies an input signal. However, the way an amplifier is designed to handle noise and bandwidth limitations greatly affects the quality and sustainability of the final output signal.

### Defining Terms

To knowledgeably discuss amplifiers, let's define a few common terms.

**Gain** – The gain is the multiplier defining how much the amplitude of an input signal is increased. A signal with an  $\times 1$  gain is not amplified. A  $\times 10$  gain produces an output signal ten times greater than the input signal.

**Noise** – Any unwanted signal fluctuations are called noise. While noise can also result from external sources, for the purpose of this discussion, we are primarily concerned with the noise resulting from the inner workings of the electronic device, our amplifier. This intrinsic noise is called shot (or schott) noise.

**Signal to Noise Ratio (SNR)** – The ratio of the output signal to the noise of the amplifier is called the signal to noise ratio. The smaller the shot noise signal in an amplifier in comparison with the output signal, the easier the desired signal is to discriminate. When engineering an amplifier, the SNR may be improved by boosting the first stage gain to yield a larger output signal or by using quality components to minimize the shot noise level of the amplifier.

**Output Range** – The output range determines the maximum output signal that can be generated with the amplifier. It is determined by the maximum voltage of the power supply. If the amplitude of the



output signal is too large for the output range, part of the signal is cut off (clipped).

**Rail** – The upper or lower limit of the amplifier range is called a rail. Signals that exceed the rail cannot be faithfully reproduced.

**DC Offset** – DC offsets can appear in biological preparations. This offset is the amount the output signal is displaced away from a zero reference point. It is usually a result of the potential difference at the electrode's tip.

### WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com Germany: +49 (0)6031 1602171 • wptide@wpi-europe.com • www.wpi-europe.com  
Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • www.wpiinc.com China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net

## How Amplifiers Work

### Power Supply Rails Limit the Range

In a perfect world an input signal can be infinitely multiplied by the gain factor to determine the output signal. For example:


Input Signal	Gain	Output Signal
2 mV	×1	2 mV
2 mV	×2	4 mV
2 mV	×10	20 mV
2 mV	×100	200 mV
2 mV	×10,000	20 V

In the real world, the power supply rails limit the possible output range of the amplifier. For example, a bio-amplifier could have a range of  $\pm 5.0V$ . In order for the output signal to be faithfully reproduced, the input signal times the gain factor must fall within the voltage window set by the power rails. Otherwise, the output signal goes off scale, and the input signal is not faithfully reproduced. This is called "hitting the rail."

In our example, a  $1.0 \mu V$  input signal at an  $\times 10^6$  gain would generate a  $1.0V$  output signal. Since the power supply is rated up to  $+5.0V$ , this output signal is clearly visible. If the input signal in this example is greater than  $5.0 \mu V$ , the output signal would be greater than  $+5.0V$ . Since  $5.0V$  is the top of the range that the power supply is capable of producing, the output signal hits the upper rail and gets cut off. This amplifier gives a  $+5.0V$  DC output signal for all input signals greater than or equal to  $5.0 \mu V$ . So, a smaller gain factor should be used to bring the output signal back into the dynamic output range of the amplifier.

### Noise Limits Amplifier Usability

All electronic devices produce their own internal electronic noise, an unavoidable signal that can mask the output signal. For example, if



*The higher the signal to noise ratio, the more discernible the desired signal.*

the input signal is  $2mV$  and the noise is  $1 mV$ , the signal to noise ratio is two to one (2:1), and the output signal would be undetectable. It is nearly impossible to discern which part of this output is generated by noise and which part is the desired signal. (See figure.)

Ideally, the signal to noise ratio should be at least 50 to 1 to produce a quality output signal. A good signal to noise ratio can be achieved in one of two ways:

- Boost the output signal by increasing the gain.
- Reduce the noise.

While increasing the gain is the simplest solution, too much gain can impose a limitation on the dynamic range of the amplifier. Reducing noise is a more complicated solution, but it offers a greater range and more stability.

### Two-Stage Amplifiers

Bio-amplifiers usually involve multiple stages of amplification.

**Stage One** – The unadulterated signal coming into the amplifier is unaffected by the intrinsic noise of the amplifier. Then, it runs through the critical first stage of amplification where the signal is boosted by the primary gain factor to produce an output signal with the desired signal

to noise ratio. The intrinsic noise is not amplified in the first stage. Higher gain factors used in the first stage of amplification can seriously limit the dynamic range available at output stage. Large stage one gains also limit the gain factor available in the second stage of amplification.

**Stage 2** – The stage one output signal enters the second stage of amplification where both the signal and the noise from the first stage are amplified together by the second stage gain factor so that the signal is large enough to be seen on a chart recorder or data acquisition system. The second stage amplification is the gain the user controls. It does not change the signal to noise ratio.

Instead of using high gains in the first stage of amplification, a well constructed bio-amplifier that uses high quality components, like WPI's DAM series amplifiers, minimizes the noise in the first stage of amplification so that the dynamic range is retained throughout the amplification process. Poorly designed amplifiers simply increase the gain of the first stage amplification until the desired signal to noise ratio is reached.

### Boost the power rails?

Theoretically, increasing the voltage rails powering the amplifier will increase the available dynamic output range. It would seem natural to increase the power supply rails coming into the amplifier in order to provide the capability for greater first stage gains. However, most data acquisition systems are limited to a maximum input signal ranging between  $\pm 10.0V$ . Therefore, it is not practical to increase the power rails of bio-amplifier beyond  $\pm 10.0V$ . Since the industry standard limits us to  $\pm 10.0V$  power supply rails, the only way to improve the signal to noise ratio is to minimize the shot noise in the first stage of amplification. This is why high quality amplifier components are imperative.

### Why a Signal Flatlines




Regardless of the amplifier used, biological potentials are often accompanied by a DC offset, because the electrodes polarize over time. The DC offset naturally increases over time. Since the poorly constructed amplifier that utilizes greater first stage gain has restricted its dynamic range, it has limited ability to handle this offset. As the offset continues to increase, the output signal may eventually be forced by the offset into the rail causing the flat line (clipping the signal). (See Figure.)

The amplifier that minimizes the noise in the first stage amplification offers a larger dynamic output range and handles a much greater offset value.

### WPI's Amplifiers

The purchase of a low-noise amplifier pays dividends in the end. WPI's amplifiers were engineered for the bio-medical researcher. While  $20-30 \mu V$  of noise is common in bio-amplifiers, WPI's DAM series amplifiers generate  $0.4 \mu V$  RMS (root mean squared) at  $0.1-100$  Hz. (That's equal to  $2 \mu V$  peak-to-peak.)

### Result of Amplifier Using Gain to Control Signal to Noise Ratio



*As the offset naturally increases over time, a poorly constructed amplifier will not be able to faithfully reproduce the signal. This offset can also result from gain drift as temperature rises.*

# Isolated Differential Amplifier

Excellent recording performance for extracellular nerve AP

## Features

- Battery powered, rechargeable
- High pass and low pass filtering
- Active remote headstage
- AC only amplification
- Electrode impedance test function
- Electrode current generation with polarity select
- Variable output positioning

## Benefits

- Ultra quiet DC power supply
- Intrinsic low susceptibility to ground loops
- High signal to noise ratio due to remote head stage
- Small footprint
- Stimulation/histological marking current

## Applications

- Amplifying bio-potentials using metal microelectrodes
- Brain slice field potentials
- EAG (Electroantennogram)
- ERG (Electroretinogram)
- Monitor extracellular nerve action potentials
- Use for cell marking, stimulation or electrode cleaning
- *In vivo* cortical recording

The **ISO-80** provides low noise AC coupled amplification and offers excellent recording performance for monitoring extracellular nerve action potentials *in vitro* and in living animals. The **ISO-80** is provided with a remote headstage (1 m cable) which incorporates an electrode impedance test function and a constant current stimulator. The constant current stimulator can be used for cell marking, stimulation or electrode cleaning. Typical applications include measuring EMG, EEG, extracellular and action potentials *in vitro* or *in vivo*. The ISO-80 system is DC isolated from the subject ground and employs state of the art electro-magnetic shielding for improved noise rejection. The amplifier employs both high pass and low pass filtering with gain from 100 to 10,000. The lowest low-pass setting is 5 Hz and the upper passband is 10 kHz.

Included with the **ISO-80** is a **Startup Kit** containing the following accessories needed for basic metal electrode electrophysiology research:


<b>CBL102</b>	Cable, BNC-to-3.5 mm plug, 6 ft (2m) (two)
<b>5469</b>	Adapter, mini-banana to 0.031 skt. (two)
<b>13388</b>	Adapter, mini-banana to 2mm skt. (two)
<b>3294</b>	Cable, ground clip to wire, 3 ft
<b>2033</b>	Mini-banana plug, black
<b>2034</b>	Mini-banana plug, red
<b>2035</b>	Mini-banana plug solderable current (two)
<b>EP1</b>	Ag/AgCl pellet (70 mm wire) 1 mm diam x 2.5 mm long
<b>M3301EH</b>	Electrode Holder, 14cm (two)
<b>5470</b>	0.031-inch jack on 12-inch wire (package of 4)

## References

Lynn, P. A., et al. (2003). Rectal intraganglionic laminar endings are transduction sites of extrinsic mechanoreceptors in the guinea pig rectum. *Gastroenterology*, 125(3), 786-94. Retrieved from <http://www.ncbi.nlm.nih.gov/pubmed/12949724>


Vastani, N., et al. (2013). Sensitivities of rat primary sensory afferent nerves to magnesium. *European Journal of Anaesthesiology*, 30(1), 21-28. <http://doi.org/10.1097/EJA.0b013e32835949ab>

Vastani, N., et al. (2013). Sensitivities of rat primary sensory afferent nerves to magnesium. *European Journal of Anaesthesiology*, 30(1), 21-28. <http://doi.org/10.1097/EJA.0b013e32835949ab>



Single-Ended Application

Optional Differential Application



ISO-80

## ISO-80 SPECIFICATIONS

INPUT RESISTANCE	>10 <sup>11</sup> Ω, Common mode and differential
INPUT LEAKAGE CURRENT	50 picoamperes, max.
AMPLIFICATION	×10 <sup>2</sup> , ×10 <sup>3</sup> , ×10 <sup>4</sup>
COMMON MODE REJECTION RATIO	100 dB typ. @ 50/60 Hz
EQUIVALENT NOISE SIGNAL INPUT	0.4 μV rms (0.1-100 Hz) 2.0 μV rms (1 Hz - 10 kHz)
FILTER SETTINGS	
Low frequency	5, 10, 100, 300 Hz
High frequency	100 Hz, 1, 3, 10 kHz
MAX. OUTPUT VOLTAGE SWING	±8 volts
ELECTRODE IMPEDANCE RANGE	100 kΩ - 10 MΩ @ 300 Hz
STIMULATION CURRENT	0 to ±20 μA (constant current)
MAXIMUM STIMULATION VOLTAGE	±15 V
MAXIMUM ELECTRODE VOLTAGE	±40 V
DISPLAY	3½-digit LCD
BATTERY TEST	Low battery display
POWER	(2) 9V Ni-MH batteries & charger, supplied
SHIPPING WEIGHT	4 lb (1.8 kg)

## ORDERING INFORMATION

<b>ISO-80</b>	Isolated Bioamplifier w/ active probe (ISO80P) <i>Specify line voltage</i>
---------------	---

## OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>ISO-80P</b>	Replacement ISO-80 Probe
<b>CBL102</b>	3.5 mm phone plug-to-BNC cable

WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • [wpiuk@wpi-europe.com](mailto:wpiuk@wpi-europe.com) • [www.wpi-europe.com](http://www.wpi-europe.com) Germany: +49 (0)6031 1602171 • [wptide@wpi-europe.com](mailto:wptide@wpi-europe.com) • [www.wpi-europe.com](http://www.wpi-europe.com)  
Brazil: 011 55 13 40629703 • [info@brazil.wpiinc.com](mailto:info@brazil.wpiinc.com) • [www.wpiinc.com](http://www.wpiinc.com) China: +86 21 6888 5517 • [chinasales@china.wpiinc.com](mailto:chinasales@china.wpiinc.com) • [www.wpiinc.net](http://www.wpiinc.net)



# Low Noise Modular Amplifier System

## Isolated, low noise bio-amplifier

### Features (Iso-DAM8A)

- High pass and low pass filtering
- Optional active remote head stage
- AC/DC amplification
- Variable gain adjustment
- Input is optically isolated
- 50/60 Hz notch filter
- Pre-optical isolation DC offset
- Post-optical isolation zeroing
- Independent module power switch

### Benefits (Iso-DAM8A)

- Chassis accepts combination of bioamplifiers and transducer amplifiers
- Flexible channel count (1–8) allows expandability
- Notch filter targets AC line noise sources
- Variable gain output amplitude
- Wide  $\pm 10V$  output range

### Applications (Iso-DAM8A)

- Amplifying biopotentials using metal microelectrodes
- Brain slice field potentials
- EAG (Electroantennogram)
- ERG (Electroretinogram)

### Iso-DAM8A Modules Isolated Low Noise Bio-Amplifier

The **ISO-DAM8A** is a compact modular standard rack-mountable DC amplifier system. Each channel is electrically isolated from the others and from ground. No current can flow from the input terminals and electrodes, thus, the instrument is intrinsically safe and cannot cause any electrical stimulus or shock to the preparation, in addition ground loop noise is minimized. Systems can be purchased with one, two, three or up to eight preamplifier modules or mixed with **Bridge8** transducer amplifier modules (see next page). The user can then select an appropriate low pass filter setting, gain and offset on the channel amplifier panel. A notch filter has been added to reduce line frequency interference. An optional headstage preamplifier (10x gain) allows low noise extracellular (DC) recording with **Iso-DAM8A** and adds greater signal bandwidth than a shielded cable of the same length. The **Iso-DAM8A** amplifier and headstage configuration is optimally suited for use with our metal microelectrodes and can be easily configured for many applications. Each amplifier channel has a coaxial (BNC) connector located on the rear panel.

### ISO-DAM8A SPECIFICATIONS

EACH CHANNEL	
INPUT IMPEDANCE	$> 10^{12}$ Ohms
INPUT LEAKAGE CURRENT	10 pA (typical)
INPUT DC OFFSET	$\pm 100$ mV
GAIN	$\times 10, \times 100, \times 1000, \times 10,000$
COMMON MODE REJECTION	$> 100$ dB @ 50/60 Hz
EQUIVALENT NOISE SIGNAL INPUT	$< 0.36 \mu V$ rms (1.8 $\mu V$ p-p) 0.1-10 Hz, Gain $>10$ $< 1 \mu V$ rms (5 $\mu V$ P-P) 0.1-10 kHz
BANDWIDTH FILTER SETTINGS	
High Filter (Low Pass) (kHz)	0.1, 0.5, 1, 3, 10
Low Filter (High Pass) (Hz)	0.1, 1, 10, 300
Notch Filter (Hz)	50, 60
OUTPUT VOLTAGE SWING	$\pm 7.5$ Volts
OUTPUT RESISTANCE	220 Ohms
ENCLOSURE DIMENSIONS	7 $\times$ 17 $\times$ 9.2 in. (18 $\times$ 43 $\times$ 23 cm)
SHIPPING WEIGHT	10 to 21 lb (4.5 to 9.5 kg)



74020 CE

## Low noise transducer amplifier

### Features (BRIDGE8)

- Wide range of fixed gains with independent variable gain adjustment
- Low pass filter
- Single ended or differential transducer compatibility
- Dual range output offset correction
- Independent module power switch
- Provides  $\pm$  voltage excitation to transducers

### Benefits (BRIDGE8)

- Chassis accepts combination of bioamplifiers and transducer amplifiers
- Flexible channel count (1–8) allows expandability
- Output LEDs confirm transducer output balance

### Applications (BRIDGE8)

- WPI force transducers
- Wheatstone bridge transducers
- Muscle force measurement



BRIDGE8 CE

### Bridge8 Modules Low Noise Transducer Amplifier

**Bridge8** is a modular, rack-mountable amplifier system. It is specifically designed for use as a signal conditioning amplifier with strain gages and other powered transducers. **Bridge8** includes differential amplifiers featuring high input impedance, high common mode rejection and low current leakage input terminals for low noise operation. It features a half bridge switch and channel offset. A wide variety of WPI transducers are available for force, temperature, pressure and light measurements. The **Bridge8** amplifier is a clear choice for convenience and quality.

### BRIDGE8 SPECIFICATIONS

INPUT IMPEDANCE	$> 10^{12}$ Ohms
AMPLIFICATION	1, 10, 50, 100, 500, 1000 & Adjustable
INPUT LEAKAGE CURRENT	0.1 pA at 25 °C
VOLTAGE OFFSET ADJUSTMENT	$\pm 50$ mV (low); $\pm 100$ mV (high)
AMPLIFIER OUTPUT VOLTAGE	$\pm 4.4$ V (10 mA, max.)
EXCITATION VOLTAGE	10 V ( $\pm 5.0$ V) 100 mA, max.
EQUIVALENT NOISE SIGNAL INPUT	$< 0.4 \mu V$ RMS (2 $\mu V$ p-p) 0.1-10 Hz Gain $> 10$ $< 3 \mu V$ RMS (15 $\mu V$ p-p) 0.1-100 Hz Gain $> 10$
LOW PASS FILTER BAND (KHz)	0.03, 0.1, 0.3, 1, 5, 10, "Wide Band" R-C Butterworth 6 dB /octave)

### ORDERING INFORMATION

<b>74020</b>	Iso-DAM8A Single Channel Module
<b>74030</b>	ISDB chassis and power supply
<b>74040</b>	Iso-DAM8A Active Headstage (separate)
<b>BRIDGE8</b>	Bridge8 Transducer Amplifier Module

### OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>74050</b>	ISDB Blank panels
<b>74016</b>	Replacement Cable, Input bare 5-ft wire
<b>2933</b>	ISDB Rack Mount Kit, 5/8-in. High
<b>FORT10g</b>	Force Transducer 10 g
<b>FORT25</b>	Force Transducer 25 g
<b>FORT100</b>	Force Transducer 100 g
<b>FORT250</b>	Force Transducer 250 g
<b>FORT1000</b>	Force Transducer 1000 g
<b>500184</b>	BNC to BNC 10 foot cable
<b>3161</b>	8-pin DIN plug
<b>3718</b>	Package of 4, 8-pin DIN (startup kit)
<b>3491</b>	Extension Cable (DIN male, DIN female), 5 ft (1.5 m)

# Dual Microprobe Intracellular Amplifier

2-Channel intracellular amplifier for dual and differential studies



## Features

- Two channels for differential or intracellular ISE
- Built in DC current generator with external control input
- Built in low pass filter
- Bridge balance circuit to null out electrode voltage drop
- Tickle circuit
- Built in test ports for each channel
- Dual capacitance compensations and output offset controls

## Benefits

- Dual channel, single ended recording
- Differential recording
- Bridge circuit nulls electrode voltage drop
- Assign low pass filter to either channel
- Very high impedance channel can be used with intracellular ISE

## Applications

- Intracellular electrophysiology using sharp micropipettes
- Brain slice intracellular recording
- In vivo intracellular recording from brain and spinal cord

For intracellular dual or differential studies, the **Duo773** has separate negative capacity controls and built-in active filtering that allows the precise balancing of time constants for artifact-free differential measurement. Comes complete with two probe headstages,  $10^{15}\Omega$  and  $10^{11}\Omega$  probes to monitor signals from ion-specific micro-electrodes as well as KCl-filled electrodes.

## Headstage for precise positioning

Two gold-plated, epoxy sealed miniature active probes can be positioned directly to the measurement site. Microelectrode holders containing an Ag/AgCl electrochemical half-cells plug directly into the probes. Stray capacitance can be reduced by placing the included driven guard shield over the microelectrode holder at the end of the probe.

## Capacity compensation

Channel A can compensate up to 10 pF of electrode shunt capacity and Channel B can compensate up to 50 pF.

## Tickler circuit for penetration

A Tickler Circuit assists in cell penetration. The frequency and amplitude of the oscillations may be varied for differences in membrane thickness or cell size. The duration of tickle can be controlled either by using the momentary switch, a foot switch, or by applying a signal to the remote tickler input.

## Active filters

Low pass settings on a -40 dB/decade active filter vary the cutoff from 1 to 30 kHz. Either probe or bridge outputs may be selected for filtering.

## Current injection

Channel B can eject current through the microelectrode by applying a command signal to the stimulus input connector. The resulting output from the probe will be a constant current replica of the input signal. Two ranges of current delivery are provided: 50 nA and 500 nA or by an external source. This source can be useful for delivering hyperpolarizing currents to stabilize the cell membrane potential and as a holding current for microiontophoresis.


## Bridge balance

Subtracts the excess electrode voltage associated with delivering current through the recording micropipette. Electrode resistances up to 1000 M $\Omega$  can be balanced in two ranges. The balanced signal is available from x10 or x50 front panel output connectors.

## Independent outputs


The **Duo773** has an output for each probe independent of gain filtering or balancing. In addition the **Duo773** has a 10x and a 50x output for easy integration to most data acquisition programs.

## Typical setup



See cables and connectors, page 96  
See Dri-Ref, page 158.

## Optional holders for intracellular amplifiers



See Microelectrode Holders, page 98.

## References

- Huo, Q., Chen, M., He, Q., Zhang, J., & Li, B. (2016). Prefrontal Cortical GABAergic Dysfunction Contributes to Aberrant UP-State Duration in APP Knockout Mice. *Cerebral Cortex*. Retrieved from <http://cercor.oxfordjournals.org/content/early/2016/08/23/cercor.bhw218.short>
- Spong, K. E., Rodríguez, E. C., & Robertson, R. M. (2016). Spreading depolarization in the brain of *Drosophila* is induced by inhibition of the Na<sup>+</sup>/K<sup>+</sup>-ATPase and mitigated by a decrease in activity of protein kinase G. *Journal of Neurophysiology*, *116*, 00353.2016. <http://doi.org/10.1152/jn.00353.2016>
- Zhang, J., Chen, M., Li, B., Lv, B., Jin, K., & Zheng, S. (2016). Altered striatal rhythmic activity in cylindromatosis knock-out mice due to enhanced GABAergic inhibition. *Neuroscience Letters*. Retrieved from <http://www.sciencedirect.com/science/article/pii/S002839081630274X>
- Bredeloux, P., Finday, I., & Pasqualin, C. (2016). 0194: Functional consequences of adrenergic receptors activation in the rat pulmonary veins and left atria. *Archives of Biochemistry and Biophysics*. Retrieved from <http://www.sciencedirect.com/science/article/pii/S1878648016304311>

## DUO 773 SPECIFICATIONS

HEADSTAGE (PROBE)	712P (red, port "B")	715P (blue, port "A")
ACTIVE PROBE INPUT IMPEDANCE	>10 <sup>11</sup> Ω	10 <sup>15</sup> Ω
GAIN	x1, x10	x1
OUTPUT RESISTANCE	100 Ω	100 Ω
OUTPUT VOLTAGE RANGE	±10 V	±10V
MAXIMUM INPUT VOLTAGE	±15 V	±15 V
PROBE LEAKAGE CURRENT	5 X 10 <sup>-12</sup> A	10 <sup>-14</sup> A
DC POSITION ADJUST RANGE	± 300 mV	± 300 mV
ELECTRODE RESISTANCE TEST CURRENT	1 nA	1 pA, 1 nA selectable
INPUT CAPACITY COMPENSATION	+10 to -50 pF	0 to -10 pF
NOISE		
Input shorted	<50 μV p-p 10kHz bandwidth	<50 μV p-p 10kHz bandwidth
20 MΩ carbon resistor	<200 μV p-p 10kHz bandwidth	<200 μV p-p 10kHz bandwidth
RISE TIME		
10-90% direct input small signal	1 μs, typical	
10-90% through 20 MΩ (-C "on")	25 μs, typical	
CURRENT INJECTION	(712P only)**	
Internal DC Current	± 50 nA low range, ± 500 nA high range	
Externally commanded Current	± 500 nA low range, ± 5 μA high range	
External current command factor	20 mV/nA low range, 2 mV/nA high range	
Current monitor	100 mV/nA low range, 10 mV/nA high range	
Compliance	3V low range, 10V high range	
Bridge balance	0-100 MΩ, 0-1000 MΩ	
Bridge amplifier gain	x 10, x 50	
LOW PASS FILTER	40 dB/decade, continuously variable 1-30 kHz	
METER SECTION		
Display	3.5-digit LED	
Ranges	200 mV, 2000 mV, 20 V, 200 nA, 2000 nA	
Accuracy and resolution	1 digit	
DIMENSIONS		
Instrument	17 x 5.25 x 10 in. (43 x 13 x 25 cm)	
Probe	Diameter: 12 mm Length: 34 mm	
POWER	95-135 V or 220-240 V, 50/60 Hz	
SHIPPING WEIGHT	15 lb (7 kg)	
CERTIFICATION	CE, CSA	

\* Although injected currents are "constant," the maximum current in a given situation will always be limited by the system compliance of 10 V.

\*\*The 712P headstage may be used on either A or B channels, however Current Injection specifications do not apply when used on channel A. The 715P headstage may not be used on the B channel.

## ORDERING INFORMATION

**SYS-773** Duo 773 Electromete  
Includes two probes (712P and 715P or two 712P) with driven guard shields and eight MEH15F microelectrode holders for 1.0, 1.2, 1.5 or 2.0 mm glass electrodes. Specify line voltage

### OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>712P</b>	Replacement probe (includes calibration)*
<b>715P</b>	Replacement probe (includes calibration)*
	*Instrument should be returned to WPI for free calibration with new probe.
<b>2933</b>	Rack Mount Kit, 5¼-in. high
<b>2547</b>	Driven Guard Shield for 712P & 715P Probes
<b>15790</b>	Replacement Probe Handle
<b>TW100F-4</b>	Glass capillary with filament
<b>TW150F-4</b>	Glass capillary with filament
<b>3259</b>	Foot Switch

# Dual Channel Differential Electrometer

*Electrochemical measurements with ion specific electrodes*

## Features

- Dual channel with very high input impedance
- Separate outputs for Channel A, B and A-B (Differential)
- Independent DC offset controls
- Test port
- Standby mode

## Benefits

- Measure changes in intracellular ion content electrochemically
- Stable and drift free
- Excellent amplification with low noise
- Driven guard shield for reduced noise and stray capacitance
- Set probe leakage current

## Applications

- Measure intracellular ion concentrations for K<sup>+</sup>, Ca<sup>2+</sup>, H<sup>+</sup> and other

The **FD223a** electrometer was designed specifically for use with intracellular ion selective electrodes fabricated using glass micropipettes and liquid ion exchangers.

The active head stages allow the researcher to locate the probes directly at the measurement site to minimize noise that would normally be picked up by longer cable runs. Driven guard shields cover the micropipette holders to further reduce the potential for interference from external sources of electromagnetic noise.

The **FD223a** is equipped with a test resistance port which is used to measure and adjust each probe for minimum leakage current. Each channel has a standby mode which clamps the head stage input voltage to zero, preventing extreme saturation or possible damage to the high impedance input amplifier.

## References

**E. Ermolayeva, H. Hohmeyer, E. Johannes, D. Sanders**  
 "Calciumdependent membrane depolarisation activated by phytochrome in the moss *Physcomitrella patens*" *Planta* 199. 1996: 352-358

### ORDERING INFORMATION

**FD223A** FD223a Dual Channel Differential Electrometer  
 2 probes, driven guard shields and micropipette holder MEH1SF  
 included for all glass microelectrodes O.D. 1.0, 1.2, 1.5 or 2.0 mm.

### OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>M3301L</b>	Left-hand Micromanipulator
<b>M3301R</b>	Right-hand Micromanipulator
<b>M-3</b>	80° Tilting base
<b>RC1T</b>	Reference cell (Ag/AgCl)
<b>2547</b>	Driven guard shield for FD223AP Probe
<b>MEH1SF</b>	Microelectrode holder
<b>FD223AP</b>	Replacement probe (includes calibration)


See cables and connectors, page 96.  
 See microelectrode holders, page 101.  
 See capillary glass, page 104.



### FD223A SPECIFICATIONS

INPUT IMPEDANCE	> 10 <sup>15</sup> Ω, shunted by 0.5 pF
INPUT CAPACITANCE	1 pF, nominal
LEAKAGE CURRENT	75 fA max
GAIN	1.000 ± 0.1%
OUTPUT RESISTANCE	50 Ω
INPUT SWING VOLTAGE	±10 V
RISE TIME (10 TO 90%)	5 μs, small signal
NOISE (0.1 Hz TO 10 KHz)	<100 μV p-p, input shorted
BASELINE STABILITY	±0.1 mV/day
POSITION CONTROLS RANGE	±600 mV
PHYSICAL DIMENSIONS	Case: 8.8 x 21.0 x 17.5 cm (H x W x D) Probe: 12.7 x 65 mm (D x L), 1.8 m cable
POWER	90-265 VAC, 50/60 Hz, 10 VA
PROBE HANDLE	6.5 x 65 mm (D x L)
SHIPPING WEIGHT	2.5 kg

**OPERATING CONDITIONS:** Equipment is intended to be operated in a controlled laboratory environment. Temperature: 0-40 °C; altitude: sea level to 2000 m; relative humidity: 0-95%.



**2547** Driven Guard Shield

### WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
 Germany: +49 (0)6031 1602171 • wptide@wpi-europe.com • www.wpi-europe.com  
 Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • www.wpiinc.com  
 China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.com

# Ultra Quiet Intracellular Amplifier

High quality, intracellular amplifier perfect for students

## Features

- Driven guard shield
- Test port
- Ground port
- Portable
- Remote headstage

## Benefits

- Cost effective
- Battery powered
- Capacitance compensation

## Applications

- Measure intracellular action potentials

**Electro 705**, a battery operated, low noise, wide band electrometer preamplifier, is designed for intracellular voltage measurement. Two 705's can be linked together to form a high impedance differential electrometer pair. Each instrument includes a miniature gold plated active probe to which a microelectrode can be attached using the WPI microelectrode holder supplied.

### Remote headstage

Easily mounted in any manipulator, this compact probe, containing the first stage of amplification, includes a microelectrode holder, which plugs directly into the probe input.

### Battery power

Four 9V alkaline batteries (included) power the **Electro 705** for approximately 500 hours giving a clean, low noise source of power, making the **Electro 705** the quietest amplifier available. Batteries can be easily tested by the press of a button.

### Capacitance Compensation

Corrects for loss of rise time caused by the presence of electrode capacity. Up to 50 pF of electrode shunt capacity may be neutralized.

### Driven Guard Shield

Stray capacitance can be further reduced by placing the driven guard shield (included) over the microelectrode holder at the input end of the probe.

### Test Features

A Tickler Circuit offers a momentary oscillation that helps achieve cell penetration. The **Electro 705** provides a 1 nA electrode test current. Electrode resistance is monitored at the 1X output as a voltage (1 mV/M). The Probe Test Port allows the convenience of testing the amplifier's intrinsic noise and gain without cumbersome external test hookups. Head stage leakage current can also be adjusted with minimum effort. The Baseline Position Control adds or subtracts up to 300 mV to the headstage output, allowing artifact voltages such as liquid junction potentials to be nulled prior to recording.

### Differential Output

Two **Electro 705** units can be connected in tandem to create an optional differential amplifier probe system.

## References

Wan, E., Kushner, J. S., Zakharov, S., Nui, X.-W., Chudasama, N., Kelly, C., ... Marx, S. O. (2013). Reduced vascular smooth muscle BK channel current underlies heart failure-induced vasoconstriction in mice. *FASEB Journal: Official Publication of the Federation of American Societies for Experimental Biology*, 27(5), 1859-67. <http://doi.org/10.1096/fj.12-223511>

Gokina, N. I., Bonev, A. D., Gokin, A. P., & Goloman, G. (2013). Role of impaired endothelial cell Ca<sup>2+</sup> signaling in uteroplacental vascular dysfunction during diabetic rat pregnancy. *American Journal of Physiology - Heart and Circulatory Physiology*, 304(7).



Photo shows two units.  
Manipulators sold separately.

Thomas, R. C., & Bers, D. M. (2013). How to make calcium-sensitive minielectrodes. *Cold Spring Harbor Protocols*, 2013(4), 370-3. <http://doi.org/10.1101/pdb.prot072850>

## ELECTRO 705 SPECIFICATIONS

Probe Input Resistance	10 <sup>12</sup> Ω
Probe Input Leakage Current	± 10 pA, adjustable to 0
Input Capacitance Compensation	0 to 50pF
Noise level	500 uV peak to peak
Rise time	15 uS, 10 - 90%
Output impedance	100 ohms, both outputs
Tickle	9V p-p adjustable from 200 Hz to 3000 Hz
Stepped Voltage Rise Time	<25 μs, 10-90%, through 20 MΩ*
Voltage Gain	x1 ±0.1%
Input Voltage Range	±5V
Position Range	±300mV
Electrode R Test	1mV/M ohm
Common Mode Rejection	> 104, in paired operation
Power	Four 9V alkaline cells, supplied
Dimensions	21.6 x 8.9 x 14.3cm (8.5 x 3.5 x 5.6")

\* Dependent on measuring technique.

## ORDERING INFORMATION

**SYS-705** Electro 705 Electrometer  
Probe, driven guard shield and micropipette holder MEH1SF included for glass microelectrodes O.D. 1.0 mm, 1.2 mm, 1.5 mm, or 2.0 mm.

### OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>M3301L</b>	Left-hand Micromanipulator
<b>MM3301R</b>	Right-hand Micromanipulator
<b>M-3</b>	80° Tilting base
<b>RC1T</b>	Reference cell (Ag/AgCl)
<b>2541</b>	Driven guard shield for 705PF Probe
<b>MEH1SF</b>	Microelectrode holder
<b>15804</b>	Replacement Probe Handle & Clip
<b>705PF</b>	Replacement probe (includes calibration)*

\*Instrument must be returned to WPI for free calibration with new probe.

See cables and connectors, page 96.

See microelectrode holders, page 98.

See capillary glass, page 104.

# Single Channel Pulse Generator

*The accuracy of digital electronics and convenience of analog controls*



CE

SYS-A310

3259 Optional foot switch

## Features

- Single channel pulse generator with train capability
- TTL and variable voltage output

## Benefits

- Variety of pulses: continuous run, single-shot, train/burst
- Multiple outputs available: monitor, isolator, sync and variable

## Applications

- Electrophysiology

The **A310** pulse generator/stimulator combines the reproducibility and accuracy of digital electronics with the fine resolution and continuous adjustment possible with analog circuitry. All timing parameters are entered with high resolution, ten-turn potentiometers and six-position range switches. Timing is accurate to within 1% of the set value.

### Variety of Pulses

Pulses can be created in continuous run, single-shot or train/burst modes. Duration of the train/burst is controlled using the onboard envelope generator or by using either of two external gating inputs. Used in conjunction with the **A360**, **A365**, **A385** or **A395**, constant current pulses and trains can be created easily. A foot switch allows hands-free, manual triggering.

### Multiple Outputs Available

Five separate standard BNC outputs are available on the front panel. The isolator output sends full pulse width control signals to any TTL triggered stimulus isolator, such as WPI's **A360**, **A365** or **A385** and others. The monitor output sends synchronized large scale full pulse width signals to recording or monitoring instrumentation such as a data acquisition system or oscilloscope. The sync output provides an additional synchronized 5  $\mu$ s TTL pulse for triggering external instrumentation. A variable voltage output provides two separate full pulse width signals in both positive and negative polarities in two ranges for applications that require a specific output voltage other than TTL.

## References

Cha, R., Marescaux, J., & Diana, M. (2014). Updates on gastric electrical stimulation to treat obesity: Systematic review and future perspectives. *World Journal of Gastrointestinal Endoscopy*, 6(9), 419–31. <http://doi.org/10.4253/wjge.v6.i9.419>

## SPECIFICATIONS

### TIMING PARAMETERS

EVENT INTERVAL	100 $\mu$ s to 1000 s*
EVENT DELAY	10 $\mu$ s to 100 s *
PULSE WIDTH	10 $\mu$ s to 100 s *
TRAIN DURATION (ENVELOPE)	100 $\mu$ s to 1000 s*
PULSE INTERVAL	20 $\mu$ s to 100 s*

### OUTPUTS

SYNC	5 $\mu$ s, TTL, and 5 V CMOS compatible, 20 mA max.	
MONITOR	10-15 V, 50 mA max.	
ISOLATOR	TTL & 5 V CMOS compatible, 20 mA max.	
VARIABLE (Pos or Neg)		
PULSED/DC	LOW RANGE	HIGH RANGE
Range	0 to $\pm$ 1 V	0 to $\pm$ 10 V
Resolution	1 mV	10 mV
NOISE		
Pulsed at 100 kHz bandwidth	<500 $\mu$ V	
DC Wide Band	<500 $\mu$ V	
OUTPUT IMPEDANCE	<1 $\Omega$	

### INPUTS

EXTERNAL SYNC	Accepts 1- $\mu$ s minimum pulses TTL, CMOS compatible
EXTERNAL GATE	Accepts 1- $\mu$ s pulse to continuous TTL, CMOS compatible

### POWER

95-130 V or 190-260 V, switch selectable single phase, 50/60 Hz

### DIMENSIONS

17 x 5.25 x 10 in. (43 x 13 x 25 cm)

### SHIPPING WEIGHT

14 lb (6.4 kg)

\*Continuously variable in six ranges. All accuracies better than 1% of set value. 50kHz maximum pulse frequency.

## ORDERING INFORMATION

**SYS-A310** Accupulser™ Signal Generator  
Specify line voltage

### OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>3259</b>	Foot Switch for A310
<b>2933</b>	Rack Mount Kit, 5¼ in. high

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • [wpiuk@wpi-europe.com](mailto:wpiuk@wpi-europe.com) • [www.wpi-europe.com](http://www.wpi-europe.com)  
Brazil: 011 55 13 40629703 • [info@brazil.wpiinc.com](mailto:info@brazil.wpiinc.com) • [www.wpiinc.com](http://www.wpiinc.com)

Germany: +49 (0)6031 1602171 • [wptide@wpi-europe.com](mailto:wptide@wpi-europe.com) • [www.wpi-europe.com](http://www.wpi-europe.com)  
China: +86 21 6888 5517 • [chinasales@china.wpiinc.com](mailto:chinasales@china.wpiinc.com) • [www.wpiinc.net](http://www.wpiinc.net)

# Constant Current Stimulus Isolator

Automated bipolar pulsing for zero net charge on biological preparations

## Features

- Constant current
- Unipolar and bipolar stimulation modes
- Built-in non-compliance alarm
- Input is optically isolated
- Standard TTL triggering
- DC test mode
- Powered by 9V alkaline or rechargeable batteries

## Benefits

- Compliance voltage is 100V or better
- Bipolar mode auto generates alternating positive and negative pulses from TTL input
- Test mode simplifies performance verification
- Optical isolation enhances safety of the preparation and reduces noise susceptibility

## Applications

- Electrophysiology
- Brain slice stimulation
- *In vivo* brain and CNS stimulation

Activated by conventional logic-level commands, Model **A365** can be gated by any pulse generator, stimulator, or computer output; automated bipolar pulsing for zero net charge on biological preparations.

### Dual Tone Audible Alarms

A tone sounds when an open electrode circuit is detected or when system compliance is reached. A second optional tone sounds when a signal is applied to the input. A test switch is also provided to check battery charge.

### Current Delivery up to 10 mA at More than 100V

Stimulus currents are set using a three-digit control knob and a three-position range switch. Output current tracks control settings to better than 1%. Output current is load independent; voltage sufficient to push the desired current through the load is automatically developed, subject only to compliance limits. Model **A365** produces up to 10 milliamperes current, in three ranges, at more than 100 volts compliance.

### Bipolar Output Polarity

Output polarity is determined by a push switch on the front panel. Bipolar current is toggled by the command waveform, setting alternating pulses as positive or negative.

## References

Lee, E., Hong, J., Park, Y.-G., Chae, S., Kim, Y., & Kim, D. (2015). Left brain cortical activity modulates stress effects on social behavior. *Scientific Reports*, 5, 13342. <http://doi.org/10.1038/srep13342>

Gindrat, A.-D., Quairiaux, C., Britz, J., Brunet, D., Lanz, F., Michel, C. M., & Rouiller, E. M. (2015). Whole-scalp EEG mapping of somatosensory evoked potentials in macaque monkeys. *Brain Structure & Function*, 220(4), 2121–42. <http://doi.org/10.1007/s00429-014-0776-y>

Younce, J. R., Albaugh, D. L., & Shih, Y.-Y. I. (2014). Deep Brain Stimulation with Simultaneous fMRI in Rodents. *Journal of Visualized Experiments*, (84), e51271–e51271. <http://doi.org/10.3791/51271>

Avila, I., & Lin, S.-C. (2014). Motivational Salience Signal in the Basal Forebrain Is Coupled with Faster and More Precise Decision Speed. *PLoS Biology*, 12(3), e1001811. <http://doi.org/10.1371/journal.pbio.1001811>

Herrera, C., Directores, R., Panetsos, F., Carlos, P., & Trueba, A. (2014). TESIS DOCTORAL Efectos de la estimulación artificial de un nervio periférico seccionado sobre la vía somatosensorial desaferezada de la rata.



A365RC

## A365 SPECIFICATIONS

OUTPUT WAVEFORM	DC or current pulse
OUTPUT CURRENT RANGES	0.1, 1.0, and 10 mA
CURRENT AMPLITUDE ERROR	0.5% of full scale, max.
CURRENT RESOLUTION	0.1% of full scale, typical
OUTPUT LOAD VOLTAGE	100 V (typically 144 V)
EXCURSION (COMPLIANCE)	100 V (typically 144 V)
EXTERNAL COMMAND VOLTAGE	5.0 V at 3.0 mA (TTL level), 10 V max.
TRIGGER THRESHOLD	2.0 V at 0.5 mA
OUTPUT POLARITY	Reversible, manual switch or automatic
CURRENT RISE TIME & DELAY	6 µs, typical (1 KΩ load)
CURRENT FALL TIME & DELAY	10 µs, typical (1 KΩ load)
OUTPUT TO GROUND RESISTANCE	10 <sup>12</sup> Ω
OPTOCOUPLER	2500 V, rated min. breakdown voltage
POWER	
Model <b>A365D</b> (Dry Cell)	16 alkaline 9 V batteries, included
Model <b>A365R</b> (RECHARGEABLE)	16 rechargeable NiMH 9 V batteries incl.
DIMENSIONS	8.5 x 3.5 x 5 in. (22 x 9 x 12 cm)
SHIPPING WEIGHT	4 lb. (1.8 kg)

## ORDERING INFORMATION

<b>SYS-A365D</b>	High Voltage Isolator, Bipolar, alkaline batteries
<b>A365RC</b>	A365R with charger (A362)
<b>SYS-A365R</b>	High Voltage Isolator, Bipolar, rechargeable
<b>SYS-A362</b>	Battery Charger for A320R, A365R, A395R

*Specify line voltage*

## OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>DRL</b>	Dummy Load Resistor Kit (set of 3)
<b>13347</b>	BNC-to-Double Banana Adapter

## A362 Battery Charger

Required for A320R, A365R, A395R

Recharges the high-voltage nickel-cadmium or NiMH battery stack in the A320R, A365R or A395R. LED lamp indicates charging status. Full charge overnight. Dimensions: 2.8 x 4.1 x 5 in. (7 x 10 x 13 cm). Shipping weight: 4 lb (1.8 kg).



# High Current Stimulus Isolator

Constant current stimulus isolator with 100 mA current range

## Features

- Constant current to 100 mA
- Unipolar or bipolar stimulation modes
- Built-in non-compliance alarm
- Input is optically isolated
- Standard TTL triggering
- DC test mode
- Powered by six rechargeable lead acid batteries
- 36V compliance
- Output polarity and output "on/off" switches

## Benefits

- 100 mA current capability
- Bipolar mode automatically generates alternating positive and negative pulses from TTL input
- Test mode simplifies performance verification
- Optical isolation enhances safety of the preparation and reduces noise susceptibility
- Battery charge status LEDs keep the experimenter informed of battery status
- Charger included at discounted price when system is purchased as **A385RC**

## Applications

- Muscle electrophysiology
- *In vivo/in vitro* muscle stimulation

The **A385** is an optically isolated current source, which can generate up to 100 mA of unipolar or biphasic constant current pulses or DC. Pulse duration is controlled manually or by an external 5V command. Output current amplitude is determined by a 3-digit 10-turn potentiometer. Maximum output voltage between the stimulating electrodes is +36V.

### Delivers Positive, Negative or Bipolar Currents

For bipolar delivery, the polarity of the output can be toggled to the opposite polarity state with each successive pulse presented to the input. Pulse duration is controlled by an externally applied voltage. The input connector is a standard BNC, allowing TTL signals from a data acquisition system to be used.

### Excellent Accuracy and Repeatability

The output amplitude is controlled by a 3-digit, ten-turn dial as a percentage of the range selected: for example, a setting of 45.6 in the 0-10 mA range translates to 4.56 mA at the output. Accuracy and repeatability are excellent. Designed for subcutaneous stimulation, maximum output voltage at the stimulating electrodes is 36 volts, reducing the possibility of accidental transcutaneous shocks. A compliance/output alarm sounds when the 36V limit is reached. Internal circuitry ensures electrodes are short-circuited during inactive periods ("electrode exhauster" feature). The **A385** is not appropriate for transcutaneous stimulation.

### Rechargeable Battery

The 1.2 amp-hour rating of the six heavy-duty lead-acid rechargeable batteries ensures that all day experiments will not be interrupted by dead batteries when charged daily. Indicator lights and audible alarms keep the user constantly apprised of the battery charge status. The batteries are recharged by the **A382** System Charger, which is designed especially for the **A385**, and included with the **A385RC**.

## References

- Lin, C., Disterhoft, J., & Weiss, C. (2016). Whisker-signaled Eyeblink Classical Conditioning in Head-fixed Mice. *Journal of Visualized Experiments*, (109), e53310–e53310. <http://doi.org/10.3791/53310>
- Li, T., Finch, E. A., Graham, V., Zhang, Z.-S., Ding, J.-D., Burch, J., ... Rosenberg, P. (2012). STIM1-Ca(2+) signaling is required for the hypertrophic growth of skeletal muscle in mice. *Molecular and Cellular Biology*, 32(15), 3009–17. <http://doi.org/10.1128>



A385RC CE

## A385 SPECIFICATIONS

OUTPUT WAVEFORM	DC or current pulse
OUTPUT CURRENT RANGES	1, 10, and 100 mA
CURRENT AMPLITUDE ERROR	0.5% of full scale, max.
CURRENT RESOLUTION	0.1% of full scale, typical
REPEATABILITY	0.1% of full scale, typical
OUTPUT LOAD VOLTAGE	36 V (typically 40 V)
EXCURSION (COMPLIANCE)	5 V at 3 mA minimum, 8.5 V max.
EXTERNAL COMMAND VOLTAGE	Reversible, manual switch, monophasic or electronically switched bipolar delivery
EXTERNAL COMMAND VOLTAGE: OUTPUT POLARITY	6 $\mu$ s, typical (1 K $\Omega$ load)
EXTERNAL COMMAND VOLTAGE: CURRENT RISE TIME AND DELAY	10 $\mu$ s, typical (1 K $\Omega$ load)
EXTERNAL COMMAND VOLTAGE: CURRENT FALL TIME AND DELAY	10 <sup>12</sup> $\Omega$
EXTERNAL COMMAND VOLTAGE: OUTPUT TO GROUND RESISTANCE	2500 V, rated minimum breakdown voltage
EXTERNAL COMMAND VOLTAGE: OPTOCOUPLER	Six rechargeable lead-acid batteries (Requires companion charger A382)
POWER	8.5 x 3.5 x 5 in. (22 x 9 x 12 cm)
DIMENSIONS	5 lb. (2.3 kg)
SHIPPING WEIGHT	

## ORDERING INFORMATION

<b>A385RC</b>	A385R with A382 Charger
<b>SYS-A385R</b>	High Current Isolator, rechargeable
<b>SYS-A382</b>	Battery Charger for A385 (see below)

*Specify line voltage*

## Smart Battery Charger

### Required for A385

An innovative three-step charger, the **A382** employs fast, medium, and trickle charges at a safe, low current, greatly extending battery life. After a fast initial phase, the charger automatically switches to a constant voltage mode. When charging is complete, the charger switches to the trickle-charge mode. LED lamps indicate charging status. (For use only in charging the lead acid batteries installed in the **A385**.)



WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • [wpiuk@wpi-europe.com](mailto:wpiuk@wpi-europe.com) • [www.wpi-europe.com](http://www.wpi-europe.com)  
Brazil: 011 55 13 40629703 • [info@brazil.wpiinc.com](mailto:info@brazil.wpiinc.com) • [www.wpiinc.com](http://www.wpiinc.com)

Germany: +49 (0)6031 1602171 • [wpi@wpi-europe.com](mailto:wpi@wpi-europe.com) • [www.wpi-europe.com](http://www.wpi-europe.com)  
China: +86 21 6888 5517 • [chinasales@china.wpiinc.com](mailto:chinasales@china.wpiinc.com) • [www.wpiinc.net](http://www.wpiinc.net)



# Linear Stimulus Isolator

Replicates a programmed waveform of any shape or polarity

## Features

- Creates a constant current replica of analog waveforms
- Amplitude of the output current is voltage controlled
- Input voltage from -10V to +10V
- 3 current ranges from 100  $\mu$ A to 10 mA
- Built-in test resistors
- Digital display shows current being delivered for non-varying currents of adequate duration
- Output offset adjustment
- $\pm 70$ V compliance range



A395RC

CE

## Benefits

- Amplitude of current is voltage controlled
- Built-in test resistances
- Error LEDs illuminate when current is less than commanded by control voltage

## Applications

- Neuroscience
- Muscle Physiology

All WPI stimulus isolators are designed to supply constant current, because current threshold (not voltage) is the most quantitatively reproducible parameter for stimulation of nerve and muscle. Model **A395** dispenses current reproducibly from its Output terminals; the amplitude being determined by the selected current RANGE and the input voltage. Current amplitude is "constant", that is, load resistance independent, provided that the  $I \times R$  (load) product does not exceed the available battery supply voltage. A visual indicator (the compliance LEDs) displays if  $I \times R$  reaches this limit. When the unit is out of compliance, one of the two LEDs (labeled - and +) illuminate, depending in which direction the current is flowing. Model **A395D** can generate a voltage of 70V or more across its OUTPUT terminals. You can be sure that the amplitude of the current is as dialed as long as the voltage drop across the load (stimulus electrode path) does not reach the magnitude of the supply voltage. The compliance LEDs will then be visible. Then, you would know that:

- Too much current was dialed for a given load or
- Inter-electrode resistance was too high or the electrode circuit path was open.


## User Defined Output Current of Various Forms

Model **A395** generates a user-defined output current of wave shape; DC, AC, pulse and combinations. Battery operated, photoelectrically-isolated from the input voltage drive, the instrument regenerates output currents which are linearly proportional to the analog voltage waveforms provided by your D/A converter or signal generator (see diagram below).

The **A395** is ideally suited for data acquisition and stimulator generators.

## Current Delivery for Selected Ranges

A 10 V input produces the maximum output current for the current range selected. (For example, 100  $\mu$ A, 1 mA, or 10 mA) Front panel controls allow DC current to be generated. Externally applied signals can be



superimposed simultaneously (DC offset). Warning lamps indicate open circuit or excessive current conditions.

## Digital Meter Shows DC or Average Output

The digital display meter shows the measures DC current or the average output current. Overload lamps indicate when output voltage has reached positive or negative compliance voltage limit.

## References

E.D. Zonnevillage, N.N. Somia, g.Perez Abadia, R.W. Stremel, C.J. Maldonado, P.M.N. Werker, M. Kon, J.H. Barker "Sequential Segmental Neuromuscular Stimulation Reduces Fatigue and Improves Perfusion in dynamic Gracilolasty" *Ann Plast Surg* 45. 2000: 292-297

## A395 SPECIFICATIONS

OUTPUT CURRENT, $I_{max}$	3 ranges: 100 $\mu$ A, 1 mA, and 10 mA
OUTPUT VOLTAGE RANGE	$\pm 70$ V
OUTPUT BANDWIDTH	10 kHz (measured across 1K $\Omega$ load R)
INPUT RESISTANCE	>20 m $\Omega$
INPUT VOLTAGE @ $I_{max}$	$\pm 10$ V
INPUT/OUTPUT LINEARITY ERROR	<0.5%
RISE, FALL TIME	26 $\mu$ s @ 10 K $\Omega$
POWER: Model A395D	17 alkaline 9 V batteries
POWER: Model A395R	17 rechargeable NiMH 9 V batteries
DIMENSIONS	6.5 x 4 x 3.5 in. (16 x 10 x 9 cm)
SHIPPING WEIGHT	4 lb. (1.8 kg)

## ORDERING INFORMATION

<b>A395RC</b>	A395R with Charger (A362)
<b>SYS-A395D</b>	Linear Stimulus Isolator
<b>SYS-A395R</b>	Linear Stimulus Isolator, Rechargeable
<b>SYS-A362</b>	Battery Charger

Specify line voltage

## A362 Battery Charger

### Required for A320R, A365R, A395R

Recharges the high-voltage nickel-cadmium or NiMH battery stack in the **A320R**, **A365R** or **A395R**.

LED lamp indicates charging status. Full charge overnight.

Dimensions: 2.8 x 4.1 x 5 in. (7 x 10.5 x 12.7 cm).  
Shipping weight: 4 lb. (1.8 kg).



SYS-A362

# 4-Channel Transducer Amplifier

*Amplify output voltage signals*

## Features

- Use with many different types of resistive based transducers
- WPI resistive force transducers plug in directly
- Supports full resistive bridge or single ended operation
- Output offset control
- Four gain ranges from 1–1000 ×
- Provides ±5V “excitation” voltage for resistive bridge transducers

## Benefits

- Provided with blank connectors to rewire to any resistive bridge transducer
- Bridge balance LEDs provide visual cue that unloaded transducers are at zero output state
- Low noise amplifier

## Applications

- Amplify signals from resistive strain gages and other resistive bridge configured transducers

Transbridge (TBM4M) is a four-channel analog transducer manifold, specifically designed to amplify output voltage signals from pressure, force, displacement, and temperature transducers as well as a wide variety of other signal sources. Analog output signals are available from each channel for input to a data acquisition system for digital signal processing in a computer. Each channel contains a regulated 10-volt power supply (+5 and -5 volts with respect to signal ground) to provide DC power to transducers, and a precision differential amplifier with selectable voltage amplification and variable position adjustment control.

Transducers can be connected to Transbridge via any of the 8-pin connectors on the front panel. Four spare 8-pin DIN plugs are provided with each instrument to allow you to rewire cables of other manufacturers' transducers and connect them to Transbridge. Each Transbridge channel may be used in either Full Bridge or the Half Bridge mode independently. For transducer types other than resistive bridges, such as active transistor circuits, magnetic, photocell or piezoelectric



*Transducers available separately*

devices, the instrument's differential amplifiers may still be used effectively for signal amplification in differential (full bridge) and single-ended (half bridge) modes.

## ORDERING INFORMATION

**SYS-TBM4M** Transbridge Transducer Amplifier  
*Specify line voltage*

## OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>13024</b>	Single Rack Mount Kit
<b>13025</b>	Dual Rack Mount Kit
<b>500184</b>	BNC-to-BNC cable, 10 ft
<b>3161</b>	8-pin DIN plug
<b>3718</b>	Package of 4, 8-pin DIN (startup kit)

# Battery Operated Impedance Measurement

*Measure mV and MΩ impedance of metal or glass microelectrodes*

## Features

- Battery operated

## Benefits

- Determine impedance of electrode during micropipette beveling process for pipette reproducibility

## Applications

- Measure impedance of metal or glass capillary microelectrodes

## ORDERING INFORMATION

**SYS-OMEGAZ** Omega-Tip-Z™ with Probe & Holder

## OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>711P</b>	Replacement Probe
<b>5468</b>	Adapter to connect metal microelectrodes to probe, 2 mm socket to .031 in. receptacle

Omega-Tip-Z™ was created especially for measuring impedance in etched tungsten, platinum-iridium\* and steel microelectrodes, as well as electrolyte-filled micropipettes. The meter's AC impedance-measuring circuit is unaffected by electrode offset or tip junction potentials. The gold-plated miniature probe lets you conveniently monitor microelectrode impedance in electrolytes, and an electrode tip cleaning feature lets you



remove buildup quickly. Omega-Tip-Z can also measure DC electrode tip potentials up to 2000 millivolts. The instrument operates for hundreds of hours without battery failure.

*\*See Metal Microelectrodes, page 94.*

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • www.wpiinc.com

Germany: +49 (0)6031 1602171 • wpide@wpi-europe.com • www.wpi-europe.com  
China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net

# Force Transducers

These rigid-lever force transducers transform applied force into proportional voltage. Using balanced strain gages, **FORT** transducers produce linear output voltage vs. applied force input with very little deflection.

To use, clamp the handle of the **FORT** transducer in a horizontal position and apply the forces to be measured to a rivet or hook mounted in the hole at the end of the flat sensing leaf.




## ORDERING INFORMATION

<b>FORT100</b>	Force Transducer (100 g)
<b>FORT250</b>	Force Transducer (250 g)
<b>FORT1000</b>	Force Transducer (1000 g)
<b>FORT5000</b>	Force Transducer (5000 g)

## FORT SPECIFICATIONS

	FORT100	FORT250	FORT1000	FORT5000
FORCE RANGES, FULL SCALE	100 g	250 g	1000 g	5000 g
OUTPUT SENSITIVITY ( $\pm 10\%$ )	7 $\mu\text{V/V/g}$	3 $\mu\text{V/V/g}$	0.84 $\mu\text{V/V/g}$	0.38 $\mu\text{V/V/g}$
INPUT & OUTPUT RESISTANCE	350 $\Omega$	350 $\Omega$	350 $\Omega$	350 $\Omega$
RESOLUTION	0.01% of full scale force	0.01% of full scale force	0.01% of full scale force	0.1% of full scale force
RESONANT FREQUENCY	300 Hz	300 Hz	300 Hz	60 Hz
LINEARITY ERROR	Less than 0.1% of full scale	Less than 0.1% of full scale	Less than 0.1% of full scale	Less than 0.1% of full scale
MAX. OPERATING VOLTAGE	10 V AC or DC	10 V AC or DC	10 V AC or DC	10 V AC or DC
MAXIMUM APPLIED FORCE	3 $\times$ rated full scale force	3 $\times$ rated full scale force	3 $\times$ rated full scale force	3 $\times$ rated full scale force
DRIFT	thermally compensated	thermally compensated	thermally compensated	thermally compensated
DIMENSIONS	0.3 inch diam $\times$ 4 in. (7.6 mm diam $\times$ 10.2 mm)	0.3 inch diam $\times$ 4 in. (7.6 mm diam $\times$ 10.2 mm)	0.3 inch diam $\times$ 4 in. (7.6 mm diam $\times$ 10.2 mm)	0.3 inch diam $\times$ 4 in. (7.6 mm diam $\times$ 10.2 mm)
WEIGHT (excluding cable)	0.3 oz (8 g)	0.3 oz (8 g)	0.3 oz (8 g)	0.3 oz (8 g)



## 10g & 25 Force Transducers

### FORT SPECIFICATIONS

	FORT10g	FORT25
FORCE RANGE, FULL SCALE	0-10 g	0-25 g
OUTPUT SENSITIVITY	2.75 mV/g, nominal	3 mV/g, nominal
INPUT & OUTPUT RESISTANCE	4000 $\Omega$	1500 $\Omega$
RESOLUTION	< 1 mg	< 2 mg
RESONANT FREQUENCY	480 Hz	450 Hz
LINEARITY ERROR	<0.2% of full scale	<0.2% of full scale
MAXIMUM OPERATING VOLTAGE	10 V DC (-5V ~ +5V or 0 ~ 10V)	10 V DC (-5V ~ +5V or 0 ~ 10V)
MAXIMUM APPLIED FORCE	2 $\times$ rated full scale force	3 $\times$ rated full scale force
DRIFT	<30 mg/hr	<50 mg/hr
DIMENSIONS	46 $\times$ 22 $\times$ 15 mm Handle 100 mm $\times$ 9.8 mm diam	40 $\times$ 22 $\times$ 17 mm Handle 109 mm $\times$ 9.8 mm diam
WEIGHT	135 g	200 g

These 10 g and 25 g force transducers are reliable tools for high precision force measurement. Using balanced semiconductor strain gages, both produce linear output voltage vs. applied force input with very little deflection. The rigid lever force transducer transforms the applied force into a proportional voltage. Featuring a temperature-compensated, full-bridge configuration with four high sensitivity semiconductor strain gages, these transducers have broad dynamic measuring range and very high sensitivity.

To use, clamp the handle of the **FORT10** or **FORT25** transducer in a horizontal position and apply the forces to be measured to a rivet or hook mounted in the hole at the end of the flat sensing leaf.

# Metal Microelectrodes

*Superior microelectrodes for outstanding extracellular recording— tungsten, iridium, platinum-iridium and Elgiloy®*

## Features

- Available in Tungsten, Platinum/Iridium, Elgiloy and Pure Iridium metal, and are insulated with a thin film of vapor-deposited Parylene-C
- Four different tip profiles also available (Standard, Heat-Treated, Blunt, and Fine tips)
- High corrosion resistance offers consistent long-term performance.

## Benefits

- Many standard types available (web) and custom
- Connecting pin fits Amphenol series 220-223 connectors.


## Applications

- Type C: Excellent for bipolar stimulation.
- For acute and chronic recording.

## Concentric Bipolar Electrodes

**Excellent for shielded macro recording as well as evoked potentials — especially well suited for bipolar stimulation**


The tungsten, platinum-iridium and elgiloy electrode is sharpened to a point and is 75 µm in diameter. The outer stainless steel conductor is insulated with Polyimide tubing to within 0.2 mm of the end of the stainless steel tube. Also available without the outer Polyimide insulation.



Insulated metal conductor with exposed concentric surface


## Concentric Bipolar, Rhodes Style Tips

This electrode is designed to replicate the Peter Rhodes-produced semi-micro concentric, and features the same length, thickness and tip dimensions. These probes feature a stepped tip with well-defined element lengths, allowing surface area to be easily determined for charge density calculations. Featuring a shaft diameter considerably thinner than the standard bipolar electrodes, this electrode is ideal for localized stimulation or population recording in fixed preparations and acute animal experiments where implant size is a critical concern.



## Concentric Bipolar with Extension, Rhodes Style Tips

These Concentric Bipolar Electrodes are designed for Neurological recording and stimulation and are ideal for acute applications (part numbers without the X) or chronic applications (part number ends with X). The Concentric Bipolar Electrodes, designed for neurological recording and stimulation and are ideal for acute applications. The Stainless Steel external tubing (E) provides support for accurate placement and manipulation during acute research. The extra extension may be cut off after insertion leaving only a small protrusion.



Nominal Impedance	EXPOSED TIP DIMENSIONS (nominal)			
	Tungsten	Elgiloy	Platinum Iridium	Pure Iridium
10 kΩ	250 µm	—	—	—
50 kΩ	200 µm	—	—	—
0.1 MΩ	100 µm	120 µm	60 µm	45 µm
0.5 MΩ	55 µm	66 µm	18 µm	14 µm
1.0 MΩ	30 µm	36 µm	10 µm	10 µm
2.0 MΩ	12 µm	15 µm	6 µm	5 µm
5.0 MΩ	5 µm	6 µm	3 µm	2.5 µm

**Heat Treated Tip** is ideal for penetrating tough membranes (not recommended for chronic implantation). This process is performed using a microforge in which the heating element is positioned in close proximity to the tip in order to melt the Parylene-C distal to the exposed metal. It provides a smooth transition and produces better adherence of the Parylene-C to the metal.




Kapton\* tubing, indicated by "KT" in the part number, extends from the connector to within 5 mm of the tip, providing stiffness and additional insulation to the electrode shaft. Kapton-clad electrodes are recommended when the electrode is to be inserted through a cannula for extra deep penetration.

\* Parylene is a trade mark of Union Carbide. Kapton is a trade mark of DuPont. Elgiloy is a trade mark of Elgiloy Ltd.

## Monopolar



## Bipolar



NOTE: Electrode diagrams are not shown to scale.

## CONCENTRIC ELECTRODES\* ORDERING INFORMATION

Item	Metal Core	Length	Imp	Probe Outer Diameter (total)	Tip Diam.	Core diam.	Y dim.	X dim. w/ polyimide	(pkg. of 5)
<b>TM33CCNON</b>	Tungsten	3" (76 mm)	10-15 KΩ	0.013" uninsulated (325 μm)	3-4 μm	.003" (76 μm)	0.4 mm	.005" (127 μm)	
<b>TM33CCINS</b>	Tungsten	3" (76 mm)	10-15 KΩ	0.016" insulated (400 μm)	3-4 μm	.003" (76 μm)	0.4 mm	.005" (127 μm)	
<b>TM53CCINS</b>	Tungsten	5" (127 mm)	10-15 KΩ	0.018" insulated (450 μm)	3-4 μm	.005" (127 μm)	0.4 mm	.008" (203 μm)	
<b>PTM3CC02INS</b>	Pt/Ir NS fine	3" (76 mm)	200 KΩ	0.013" insulated (325 μm)	2-4 μm	0.002" (50.8 μm)	.25 mm	.004" (114 μm)	

\*All have a stainless steel outer shaft

## METAL ELECTRODES ORDERING INFORMATION

Item	Length	Insul. Thick	Shaft Diam.	Nominal Impedance (± 20%)	Tip Diam.	Typical Use	Package of 10
<b>Tungsten — Profile A</b>							
<b>TM31A10</b>	76 mm	1 μm	0.127 mm	1.0 MΩ	1 μm	Multi unit and single unit recording	
<b>TM31A20</b>	76 mm	1 μm	0.127 mm	2.0 MΩ	1 μm	Multi unit and single unit recording	
<b>TM33A05</b>	76 mm	3 μm	0.127 mm	0.5 MΩ	1 μm	Multi unit and single unit recording	
<b>TM33A10</b>	76 mm	3 μm	0.127 mm	1.0 MΩ	1 μm	Multi unit and single unit recording	
<b>TM33A20</b>	76 mm	3 μm	0.127 mm	2.0 MΩ	1 μm	Multi unit and single unit recording	
<b>TM33B01</b>	76 mm	3 μm	0.254 mm	0.1 MΩ	1-2 μm	Single and multi unit recording	
<b>TM33B05</b>	76 mm	3 μm	0.254 mm	0.5 MΩ	1-2 μm	Single and multi unit recording	
<b>TM33B10</b>	76 mm	3 μm	0.254 mm	1.0 MΩ	1-2 μm	Single and multi unit recording	
<b>TM33B20</b>	76 mm	3 μm	0.254 mm	2.0 MΩ	1-2 μm	Single and multi unit recording	
<b>TM33C10</b>	76 mm	1 μm	0.085 mm	1.0 MΩ	1 μm	Single and multi unit recording	
<b>Tungsten — Profile C</b>							
<b>TM33A10KT</b>	76 mm	3 μm	0.216 mm	1.0 MΩ	1 μm	Multi unit and single unit recording	
<b>TM33B01KT</b>	76 mm	3 μm	0.356 mm	0.1 MΩ	1-2 μm	Single and multi unit recording	
<b>Elgiloy®/Stainless — Profile A</b>							
<b>SSM33A70</b>	76 mm	3 μm	0.229 mm	7.0 MΩ	1-2 μm	Recording and Stimulating (Prussian blue staining)	
<b>SSM33A20KT</b>	76 mm	3 μm	0.356 mm	2.0 MΩ	1-2 μm	Recording and Stimulating (Prussian blue staining)	
<b>Tungsten — Profile B</b>							
<b>TST33A001KT</b>	76 mm	3 μm	0.356 mm	10 kΩ	1 μm	Tissue slice stimulation	
<b>TST33A05KT</b>	76 mm	3 μm	0.356 mm	0.5 MΩ	1 μm	Stereotrode / Bipolar, differential measurements	
<b>TST33A10KT</b>	76 mm	3 μm	0.356 mm	1.0 MΩ	1 μm	Stereotrode / Bipolar, differential measurements	
<b>TST33A20KT</b>	76 mm	3 μm	0.356 mm	2.0 MΩ	1 μm	Stereotrode / Bipolar, differential measurements	
<b>TST33C05KT</b>	76 mm	3 μm	0.216 mm	0.5 MΩ	1 μm	Stereotrode / Bipolar, diff. meas. — extra fine (75 μm separation)	
<b>TST53A10KT</b>	127 mm	3 μm	0.356 mm	1.0 MΩ	1-2 μm	Stereotrode / Bipolar, differential measurements	
<b>Pure Iridium — Profile A</b>							
<b>IRM23E10</b>	50 mm	3 μm	0.106 mm	1.0 MΩ	1-2 μm	Single and multiunit recording and stimulation	
<b>IRM23E15</b>	50 mm	3 μm	0.106 mm	1.5 MΩ	1-2 μm	Single and multiunit recording and stimulation	
<b>IRM23E25</b>	50 mm	3 μm	0.106 mm	2.5 MΩ	1-2 μm	Greater selectivity - small cells	
<b>IRM23E30</b>	50 mm	3 μm	0.106 mm	3.0 MΩ	1-2 μm	Greater selectivity - small cells	
<b>Pure Iridium — Profile B</b>							
<b>IRM23E20KT</b>	50 mm	3 μm	0.180 mm	2.0 MΩ	1-2 μm	Greater selectivity & microstimulation	
<b>IRM23E30KT</b>	50 mm	3 μm	0.180 mm	3.0 MΩ	1-2 μm	Greater selectivity - small cells	

Elgiloy Steel \*Cobalt/chromium/nickel alloy. The KT suffix refers to Kapton™ cladding.

All Metal Microelectrodes are available in custom lengths, blunt or heat treated (extra charge).

Assorted kits may be customized.

### Ordering:

Add the **B** suffix where blunt electrodes are desired. (For example, an IRM123A10KT ordered as a blunt will be IRM123A10KT<sub>B</sub>.)


Add the **H** suffix where heat treated electrodes are desired. (For example, an IRM123A10KT ordered with heat treatment will be IRM123A10KT<sub>H</sub>.)

## ORDERING INFORMATION


### OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>300102</b>	Micromanipulator holder, 4 in., 2mm to 0.031 socket
<b>5468</b>	2 mm receptacle to 0.031-inch jack (for Omega-TipZ)
<b>5469</b>	Adapts mini banana plug (DAM80) to 0.031-inch receptacle (metal microelectrode)
<b>5470</b>	0.031-inch jack, 28 ga. wire, 12 inch (pkg. of 4)
<b>5482*</b>	Pins, 0.031-inch, gold-plated (pkg. of 50)
<b>5483*</b>	Sockets, 0.031-inch gold-plated (pkg. of 50)

\*Gold-plated pins (#5482) and sockets (#5483) may be attached to 24-, 26-, or 28-gauge wire.



Gold-plated pins (5482) and sockets (5483) may be attached to 24-, 26-, or 28-gauge wire.



# Metal Electrodes, Rhodes Style Tips with Extension


## Concentric bipolar electrodes for neurological recording and stimulation

### Features

- Rhodes-style tips are great for those that require a precisely controlled electrode surface area
- The stainless steel external tubing provides support for accurate placement and manipulation during acute research
- These are an excellent replacement for long-term Rhodes electrodes (originally made by Peter Rhodes)

### Benefits


- Concentric bipolar electrodes are ideal for bipolar stimulation paradigms
- Tools are available for shielded macro recordings, as well as evoked potential
- Electrodes may be acutely reused many times (> 20-30 insertions), if they are properly cared for:
  - Avoid bending the shaft and blunting/hooks the tip
  - Clean electrode carefully between session




### Applications

- Designed for neurological recording and stimulation and are ideal for chronic applications
- The Concentric Bipolar Electrodes, designed for Neurological recording and stimulation, are ideal for acute applications (without the X) or chronic applications (ends with X). The stainless steel external tubing (E) provides support for accurate placement and manipulation during acute research. The extra extension may be cut off after insertion leaving only a small protrusion.

### Concentric Bipolar Electrodes with Extension



### Semi-Micro Concentric Bipolar Electrodes with Extension



### ORDERING INFORMATION

Item	Type	Metal Core	Length (s)	Probe Total OD (z)	OD Inner sleeve (x)	Extension (E)
SSM04RC100	Semi micro concentric with extension	Stainless Steel	0.4" (10 mm)	Insulated (411 µm)	140 µm	50 mm
PTM04RC100	Semi micro concentric with extension	Platinum Iridium	0.4" (10 mm)	Insulated (411 µm)	140 µm	50 mm
SSM041RC225	Concentric bipolar with extension	Stainless steel	0.4" (10 mm)	insulated (625 µm)	305 µm	50 mm
PTM041RC225	Concentric bipolar with extension	Platinum iridium	0.4" (10 mm)	insulated (625 µm)	305 µm	50 mm


\*All have a stainless steel outer shaft. Package of 5.

### WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
 Germany: +49 (0)6031 1602171 • wptide@wpi-europe.com • www.wpi-europe.com  
 Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • www.wpiinc.com  
 China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net


# Metal Electrodes, Rhodes Style Tips, No Extension

Concentric bipolar electrodes for neurological recording and stimulation



## Concentric Bipolar Electrodes

Choose stainless steel or platinum iridium.




### ORDERING INFORMATION

Item	Metal Core	Length (s)	Probe Total OD (z)	OD of Inner Sleeve (x)	Outer Polyimide	Extension (E)
SSM321RC225	Stainless steel	3.125" (80 mm)	insulated (625 µm)	225 µm	Yes	30 mm
PTM321RC225	Platinum iridium	3.125" (80 mm)	insulated (625 µm)	225 µm	Yes	30 mm

## Semi-Micro Concentric Bipolar Electrodes

Choose stainless steel or platinum iridium.



### ORDERING INFORMATION

Item	Metal Core	Length (s)	Probe Total OD (z)	OD of inner sleeve (x)	Outer Polyimide	Extension (E)
SSM321RC100	Semi micro concentric	Stainless steel	3.125" (80 mm)	insulated (411 µm)	140 µm	30 mm
PTM321RC100	Semi micro concentric	Platinum iridium	3.125" (80 mm)	insulated (411 µm)	140 µm	30 mm

\*All have a stainless steel outer shaft.

# Nerve Cuff Electrodes


*For acute/chronic experiments*

## Features

- For both acute and chronic experiments
- 1-24 electrode sites
- Broad range of available inner diameters, from 5 mm down to as small as 56 µm
- Platinum, stainless steel and platinum/iridium metal electrodes

## Benefits

- Customizable electrode arrangements, including concentric and tri-bipolar
- Adaptable to interface with the acquisition system of your choice



## Applications

- Suitable for rodent, feline, bird and primate research

Nerve cuff electrodes are designed for reliable recording and/or stimulation of any peripheral nerve. They can be used acutely or for chronic implantation, with a wide selection of inner diameters available to be selected based on the nerve diameter. The nerve cuffs are designed to provide flexibility in electrode contact location for recording and stimulation protocols. Three different electrode variations are available:

- Standard, Micro Cuff and Nanocuff
- Concentric Nerve Cuff
- X-wide Nerve Cuffs

### 1. Standard, Micro Cuff and Nano Cuff

Many custom combinations of inner diameters and contact arrangements are available. Standard nerve cuffs have a variety of configurable parameters.



### 2. Concentric Nerve Cuff


Concentric electrodes have multiple contacts around a single point of the nerve, allowing recording or stimulation at different locations around the same point.



### AVAILABLE OPTIONS


Configuration	Metal Type	ID
Standard	125µm stranded Stainless Steel	1.0 – 5.0 mm
	100µm Platinum	1.0 – 2.0 mm
	250µm Platinum*	2.0 – 5.0 mm
Micro Cuff	100µm Platinum	0.5 – 0.75 mm
	50µm Platinum/Iridium	0.3 mm
Nano Cuff	25µm Platinum/Iridium	160 - 250 µm
	12.5µm Platinum/Iridium	56 - 140 µm

\* Recommended for stimulation



### Configuration Options for Standard, Micro Cuff, Nano Cuff & Concentric Nerve Cuffs


Many custom combinations of inner diameters and contact arrangements are available for standard, micro and nano nerve cuffs.



- 1 Optional Sutures
- 2 1 to 24 Contacts (Platinum, Platinum/Iridium or Stainless Steel)
- 3 Custom diameter
- 4 Custom length
- 5 Stranded stainless steel leads. The angle is selectable.



### 3. X-Wide Contact Nerve Cuffs



- d** Inner diameter  
**x** Spacing between contacts  
**n** Number of contacts  
**y** Distance from the last contacts to the edge of the cuff  
**s** Stranded stainless steel leads with Teflon insulation 300 mm long (typical)  
**w** Width of the contact, Pt/Ir 90/10

- Extra large platinum ribbon contacts maximize electrode surface area while maintaining the reliability and customizability
- Ideal for any electrical stimulation studies, especially those using KHz stimulation to block activity of peripheral nerves
- Lower thresholds and increased safety margins provide you with new freedom to develop even more effective stimulation paradigms
- The same reliability enhancement features as our standard cuffs, including reinforced bonds and specially-designed helical leads, providing exceptional implanted longevity
- Fully-customizable specifications, including contact spacing, cuff length and lead features, with a wide range of available inner diameters from over 6 mm down to 300  $\mu$ m.
- Configurable for both acute and chronic use
- Connector options for any available stimulation or recording system

The key utility for these cuffs is provided by the very large contact surface area, which is ideal for studies with high current levels and when using kilohertz electrical stimulation to achieve block of activity in peripheral nerves. The reduced impedances ensure that the output of stimulation sources is not attenuated, that no destructive charge builds up, and that true block of nerve activity is achieved.

### Types of Nerve Cuff Electrodes

Application	Electrode Type	d (mm)	x (mm)	y (mm)	W (mm)
Sciatic	Rat, 1.4 mm avg. nerve OD, 2.5 mm cuff ID	2.5	1	1.5	0.5
Tibial	Rat, 1.2 mm avg. nerve OD, 1.5 mm cuff ID	1.5	1	1.5	0.5
Cervical Vagus	Rat, 0.5 mm avg. nerve OD, 0.75 mm cuff ID	0.75	1	1.5	0.5
Subdiaphragmatic Vagus	Rat, 0.3 mm avg. nerve OD, 0.5 mm cuff ID	0.5	1	1.5	0.5
Custom Standard		1.0-5.0			
Custom Microcuff		0.5 or 0.75			

### Sterilization

Nerve cuffs are made entirely of autoclavable materials — silicone rubber, Teflon and stainless steel. They can be steam autoclaved without special precautions. If gas sterilization (ETO) is preferred, pack the nerve cuffs in a gas-permeable bag and allow adequate outgassing time (at least 48 hrs.) to be sure all toxic gases have been desorbed from the silicone rubber.

Order online at  
[www.wpiinc.com/nce](http://www.wpiinc.com/nce)

### Nerve Cuff Connectors

Omnetic connectors may be ordered separately. Other connectors must be ordered with nerve cuff electrodes. Connector options include:



Banana Terminal



Touch Proof Connector



Invivo1

- No connector
- Invivo1 (not sold separately)
- Banana Terminal (not sold separately)
- Touch Proof Connector (not sold separately)
- Omnetic Connectors

Male Implant Connector	Description
	<b>MP-A11365-001</b> Male Omnetics connector - 8 channel, 10-pin 0.050" pitch
	<b>MP-A12623-001</b> Male Omnetics connector- 8 channel, 10-pin, 0.050" pitch with latching mechanism
	<b>MP-A8393-001</b> Male Omnetics connector- 8 channel, 10-pin, 0.025" pitch, 2 guide post
	<b>MP-A79000-001 (NPS-09-DD-GS)</b> Male Omnetics connector- 7 channel, 9-pin, 0.050" pitch
	<b>MP-A79014-001 (NPD-18-DD-GS)</b> Male Omnetics connector- 16 channel, 18-pin, 0.025" pitch, 6 guides posts
	<b>MP-A70242-001</b> Male Omnetics connector- 16 channel, 18-pin, 0.025" pitch, 6 guides posts; Nickel free, MRI compatible
	<b>MP-A79038-001 (NPD-18-DDGS)</b> Male Omnetics connector- 16 channel, 18-pin, 0.025" pitch, 2 guide posts
	<b>MP-A79022-001 (NPD-36-DD-GS)</b> Male Omnetics connector -32 channel, 36-pin, 0.025" pitch, 4 guide posts
	<b>MP-A72312-001</b> Male Omnetics connector -32 channel, 36-pin, 0.025" pitch, 4 guide posts; Nickel free, MRI compatible

### References

- Guillemin, K., & Ahmed, S. H. (2017). Preference for Cocaine is Represented in the Orbitofrontal Cortex by an Increased Proportion of Cocaine Use-Coding Neurons. *Cerebral Cortex*, 28(3), 819–832. <http://doi.org/10.1093/cercor/bhw398>
- Uy, J., Yu, M., Jiang, X., Jones, C., Shen, B., Wang, J., ... Tai, C. (2017). Glutamatergic Mechanisms Involved in Bladder Overactivity and Pudendal Neuromodulation in Cats. *The Journal of Pharmacology and Experimental Therapeutics*, 362(1), 53–58. <http://doi.org/10.1124/jpet.117.240895>
- Khurram, A., Ross, S. E., Sperry, Z. J., Ouyang, A., Stephan, C., Jiman, A. A., & Bruns, T. M. (2017). Chronic monitoring of lower urinary tract activity via a sacral dorsal root ganglia interface. *Journal of Neural Engineering*, 14(3), 036027. <http://doi.org/10.1088/1741-2552/aa6801>
- Restaino, S. M., Abliz, E., Wachrathit, K., Krauthamer, V., & Shah, S. B. (2014). Biomechanical and functional variation in rat sciatic nerve following cuff electrode implantation. *Journal of NeuroEngineering and Rehabilitation*, 11(1), 73. <http://doi.org/10.1186/1743-0003-11-73>
- Lissandrello, C. A., Gillis, W. F., Shen, J., Pearre, B. W., Vitale, F., Pasquali, M., ... Gardner, T. J. (2017). A micro-scale printable nanoclip for electrical stimulation and recording in small nerves. *Journal of Neural Engineering*, 14(3), 036006. <http://doi.org/10.1088/1741-2552/aa5a5b>


# Cables & Connectors

For wiring any laboratory setup



*BNC Cables*  
 2851 (6 ft)  
 500184 (10 ft)  
 500257 (6 in.)  
 500258 (12 in.)  
 500259 (18 in.)

3492



5374



13685



13776



3517



5375



13854



3142



3578



5385



14254




3161



3670



13324



15623



3294



13347



15975



300040



15976



5371



13388




500256



3417-10



13451



13620



3491



5373



## ORDERING INFORMATION

Part #	Application/Description	Connector A	Connector B	Cable Length
<b>2026-10</b>	2 mm Socket, unwired (pkg. of 10) (Not Shown)	2 mm socket	unwired	none
<b>2851</b>	Standard BNC Cable	BNC (male)	BNC (male)	6 ft. (1.8 m)
<b>3142</b>	Mini-Banana Adapter	Screw Terminals	Dual Mini-Banana	none
<b>3161</b>	Connector for input to TBM4M and BP-1	DIN (male)	unwired	none
<b>3294</b>	Ground Wire for DAM80 probe	Clip	none	3 ft. (0.9 m)
<b>3417-10</b>	2 mm Plug, unwired (pkg. of 10)	2 mm pin	unwired	none
<b>3491</b>	Extension for any 8-pin DIN	DIN (male)	DIN (female)	5 ft. (1.5 m)
<b>3492</b>	Connector, adapts WPI transducers to non-WPI equipment	DIN (female)	unwired	none
<b>3517</b>	DAM50, DAM60, DAM70, shielded (two cables/pkg)	Modular phone plug, 4 wire	none	3 ft. (0.9 m)
<b>3578</b>	Adapter cable for Ag/AgCl pellets	2 mm pin	none	5 ft. (1.5 m)
<b>3670</b>	Double banana plug with solder turret terminals	Dual Banana (male)	Dual Banana (female)	none
<b>5371</b>	Low-noise cable for microelectrode holders	2 mm gold pin/jack	2 mm gold pin/jack	2 ft. (0.6 m)
<b>5373</b>	Low-noise cable for microelectrode holders	2 mm gold pin/jack	2 mm gold pin/jack	2 ft. (0.6 m)
<b>5374</b>	Low-noise cable for microelectrode holders	BNC (male)	2 mm gold pin	4 ft. (1.2 m)
<b>5375</b>	Low-noise cable for microelectrode holders	BNC (male)	2 mm gold jack	4 ft. (1.2 m)
<b>5385</b>	Cable, shielded transducer stock	none	none	25 ft. (7.6 m)
<b>13324</b>	Adapter	Double-banana (female)	BNC (male)	none
<b>13347</b>	ISO2 (chart recorder adapter)	Double-banana (male)	BNC (female)	none
<b>13388</b>	Electrode adapter for DAM probes	Miniature banana (male)	2 mm jack	none
<b>13451</b>	Adapter: Iso-DAM, Iso-DAM8	BNC (female)	two 2 mm pins	6 in. (15 cm)
<b>13620</b>	Low-noise cable for microelectrode holders	2 mm gold pin	2 mm gold jack	2 ft. (0.6 m)
<b>13685</b>	SP Series pump-to-pump linking cable	Modular phone plug	Modular phone plug	7 ft. (2.1 m)
<b>13776</b>	Adapts reference electrode to VF4 ground jack	Banana (male)	2 mm jack	none
<b>13854</b>	BNC T-connector, male to:	BNC (female)	BNC (female)	none
<b>14254</b>	BNC Straight Adapter	BNC (female)	BNC (female)	none
<b>15623</b>	Serial cable and adapter, SP Series pump	SP Pump	IBM 9-pin "D" connector	5 ft. (1.5 m)
<b>15975</b>	Adapter	2 mm socket	1 mm pin	none
<b>15976</b>	Adapter	1 mm socket	2 mm pin	none
<b>300040</b>	Adapter Extension	2 mm socket	2 mm socket	4 in. (10 cm)
<b>500184</b>	Standard BNC Cable	BNC (male)	BNC (male)	10 ft. (3 m)
<b>500256</b>	BNC Right Angle Adapter	BNC (male)	BNC (female)	none
<b>500257</b>	Standard BNC Cable	BNC (male)	BNC (male)	6 in. (15 cm)
<b>500258</b>	Standard BNC Cable	BNC (male)	BNC (male)	12 in. (30 cm)
<b>500259</b>	Standard BNC Cable	BNC (male)	BNC (male)	18 in. (46 cm)
<b>503301</b>	Cable, Extension	8-pin miniDIN (male)	8-pin miniDIN (female)	10 ft. (3 m)
<b>503536</b>	Cable, USB	USB (male)	USB (female)	10 ft. (3 m)
<b>504713</b>	Cable (red and black pair)	Banana (male)	Banana (male)	36 in. (91 cm)
<b>504714</b>	Cable (red and black pair)	Banana (male)	Mini-Gator	36 in. (91 cm)
<b>504715</b>	Cable (red and black pair)	Banana (male)	Mini-Clip	36 in. (91 cm)
<b>504716</b>	Cable (red and black pair)	Banana (male)	Micro-Clip	36 in. (91 cm)
<b>CBL100</b>	MiniPhone Patch Cable	3.5 mm MiniPhone plug	3.5 mm MiniPhone plug	6 ft. (1.8 m)
<b>CBL102</b>	DAM Series, PM Series	3.5 mm MiniPhone plug	BNC (male)	6 ft. (1.8 m)



504716



504713



CBL100



CBL102



503536



504715



504714

# Micropipette Holders & Half-Cells

Couple fluid filled glass micropipettes to high input impedance amplifiers

## Features

- Connector pin, jack, luer ports, Ag/AgCl half-cells, wire selections available
- Optional handles (for some varieties) allow for easy manipulation. Handles are sold separately.
- Screw cap for tight fitting of glass electrode to avoid drifting of electrode

## Benefits

- Available for 1.0, 1.2, 1.5 and 2.0 mm OD glass electrodes
- Spare gaskets available
- Optional handles available in two different sizes for some holders

## Applications

- Microinjection
- Electrophysiology recording
- Fluid handling

◆ denotes holder manufactured for you as custom orders. Call for price.

The most popular micropipette holders are stock items. Custom holders (designated by ◆) can be manufactured on demand but require an additional setup fee. Call for a quote. See all the options at [www.wpiinc.com/MEH](http://www.wpiinc.com/MEH).






WPI's microelectrode holder-half-cells couple fluid-filled glass micropipettes to high input impedance amplifiers. A Ag/AgCl pellet (or a silver wire) molded into the holder body provides stable potential. Electrical connection is made via male 2 mm pins or female 2 mm sockets. The pipette may be mounted axially or at right angles to the holder. Pipettes are held with screw-caps or rubber gaskets (without caps). Filling WPI microelectrode holders with electrolytes containing chloride results in stable electrode potential. Suitable electrolytes include KCl, NaCl and CaCl<sub>2</sub>. Holders are supplied for standard WPI single capillary tubing of 1.0, 1.2, 1.5 and 2.0 mm outside diameters. (Call WPI regarding custom designs for other glass diameters.) The holder style you select will depend on your experimental application, space, and instrumentation.

### Hints for selecting and ordering micropipette holders





















1. Determine the required electrical connection on the holder: for example, if you wish to connect the holder to a 2 mm pin you should select a holder equipped with a 2 mm jack. Most WPI probes require a holder equipped with a 2 mm jack.
2. Decide on the required alignment of the electrical connection: either in-line with the glass pipette, or at a right angle to it. Space considerations in your experimental setup and requirements imposed by other pieces of equipment typically determine which alignment is appropriate.
3. Determine if you want to hold the glass pipette by a rubber gasket (e.g., MEH1S) or a screw-cap (e.g., MEH3S). Rubber gaskets offer easier insertion and removal of glass pipettes whereas screw-caps provide more secure mounts for micropipettes.

4. Choose a holder with either a silver wire or a silver/silver chloride pellet for the metal/liquid coupling. Silver/silver chloride pellets provide a more stable low-noise baseline which is important for low-noise DC recording. Pellets require the glass pipette and holder to be free of air bubbles to achieve a good connection. Silver wire holders are durable and are easier to use when the holder is equipped with a pressure port because the fluid in the pipette does not have to be filled to the top of the pipette to achieve a good electrical connection
5. Choose a holder equipped with a pressure port only when you want to pressure inject liquid from the pipette. Two types of ports are available: 2.0 mm O.D. and standard "syringe-style" Luer. The Luer port is often recommended because it makes assembly and disassembly much easier. Quick-connect Luer fittings for four common sizes of tubing (1/16", 1/8", 5/32" I.D.) are included with each Luer-equipped holder.
6. Some non-WPI preamplifiers or headstages cannot be mounted on micromanipulators. In such cases, a holder equipped with a rod (e.g., MEH8) permits the holder to be conveniently mounted on a micromanipulator.
7. Finally, remember to specify the O.D. of the glass you will be using when you place your order.

MEH6RF/SF is designed primarily for use with the Model 900A Micropressure System; EHB1 for use in electrode beveling; and MEH3SW for microtitration of chloride with a silver wire as the electrode and a solution of silver nitrate filling the holder. MPH models do not contain Ag/AgCl half-cells and are used for pressure injection of substances through microelectrodes. PicoNozzle, used for pressure injection with PV800 Series PicoPumps, includes an MPH6S holder — which may also be used to couple a micropipette to a syringe.



















◆ denotes holder manufactured for you as custom orders. Call for price.



# Micropipette Holders & Half-Cells

Couple fluid filled glass micropipettes to high input impedance amplifiers



♦ denotes holder sizes manufactured for you as custom orders. Call for price.

**MICROELECTRODE HOLDERS ORDERING INFORMATION**

◆ = custom order (call for price)

Order Number Replace XX with glass OD *	Custom Order	Electric Connection Angle	Connector	Half-Cell	Pressure Port	Screw Cap	Designed for WPI Products
<b>EHBF</b>		Straight	Female	Wire	No Port	N/A	MBS, 48000
<b>FOIMPH</b>		Straight	Fiber Optic	None	No Port	w/Cap	MBS, 48000
<b>MEH1Sxx</b>		Straight	Male	Pellet	No Port	No Cap	
<b>MEH1SFxx</b>		Straight	Female	Pellet	No Port	No Cap	705, 773, 767, 721, FD223
<b>MEH2Rxx</b>	◆	Right	Male	Pellet	Male Luer	w/Cap	
<b>MEH2RFxx</b>	◆	Right	Female	Pellet	Male Luer	w/Cap	705, 773, 767, 721, FD223
<b>MEH2RFWxx</b>	◆	Right	Female	Wire	Male Luer	w/Cap	705, 773, 767, 721, FD223
<b>MEH2RWxx</b>	◆	Right	Male	Wire	Male Luer	w/Cap	
<b>MEH2Sxx</b>		Straight	Male	Pellet	Male Luer	w/Cap	
<b>MEH2SFxx</b>		Straight	Female	Pellet	Male Luer	w/Cap	705, 773, 767, 721, FD223
<b>MEH2SFWxx</b>	◆	Straight	Female	Wire	Male Luer	w/Cap	705, 773, 767, 721, FD223
<b>MEH2SWxx</b>		Straight	Male	Wire	Male Luer	w/Cap	
<b>MEH345xx</b>		45°	Male	Pellet	No Port	w/Cap	
<b>MEH3F45xx</b>	◆	45°	Female	Pellet	No Port	w/Cap	705, 773, 767, 721, FD223
<b>MEH3FW45xx</b>	◆	45°	Female	Wire	Port	w/Cap	
<b>MEH3Rxx</b>		Right	Male	Pellet	No Port	w/Cap	
<b>MEH3RFxx</b>		Right	Female	Pellet	No Port	w/Cap	705, 773, 767, 721, FD223
<b>MEH3RFWxx</b>	◆	Right	Female	Wire	No Port	w/Cap	705, 773, 767, 721, FD223
<b>MEH3RWxx</b>	◆	Right	Male	Wire	No Port	w/Cap	
<b>MEH3Sxx</b>		Straight	Male	Pellet	No Port	w/Cap	
<b>MEH3SBWxx</b>	◆	Straight	Banana	Wire	No Port	w/Cap	ISO-80, ISO-DAM8A
<b>MEH3SFxx</b>		Straight	Female	Pellet	No Port	w/Cap	705, 773, 767, 721, FD223
<b>MEH3SFWxx</b>		Straight	Female	Wire	No Port	w/Cap	705, 773, 767, 721, FD223
<b>MEH3SWxx</b>		Straight	Male	Wire	No Port	w/Cap	
<b>MEH3W45xx</b>	◆	45°	Male	Wire	No Port	w/Cap	705, 773, 767, 721, FD223
<b>MEH6RFxx</b>		Right	Female	Pellet	2.0-mm Port	w/Cap	705, 773, 767, 721, FD223
<b>MEH6RFWxx</b>	◆	Right	Female	Wire	2.0-mm Port	w/Cap	705, 773, 767, 721, FD223
<b>MEH6SFxx</b>		Straight	Female	Pellet	2.0-mm Port	w/Cap	705, 773, 767, 721, FD223
<b>MEH6SFWxx</b>	◆	Straight	Female	Wire	2.0-mm Port	w/Cap	705, 773, 767, 721, FD223
<b>MEH7xx</b>		Right	Male	Pellet	2.0-mm Port	w/Cap	
<b>MEH7Wxx</b>	◆	Right	Male	Wire	2.0-mm Port	w/Cap	
<b>MEH8xx</b>		Right	Male	Pellet	No Port	w/Cap	
<b>MEH900Rxx</b>		Right	Male	Pellet	2.0-mm Port	w/Cap	900A
<b>MEH900Sxx</b>		Straight	Male	Pellet	2.0-mm Port	w/Cap	900A
<b>MPH1xx</b>		—	None	None	Female Luer	w/Cap	
<b>MPH3xx</b>		—	None	None	Male Luer	w/Cap	
<b>MPH4xx</b>		—	None	None	2.0-mm Port	w/Cap	
<b>MPH6Pxx</b>	◆	Right	Male	Pellet	Female Luer	w/Cap	Piconozzle Kit (5430-XX)
<b>MPH6Rxx</b>	◆	Right	Male	Wire	Female Luer	w/Cap	Piconozzle Kit (5430-XX)
<b>MPH6Sxx</b>		—	None	None	Female Luer	w/Cap	Piconozzle Kit (5430-XX)

\* Specify O.D. of glass (1.0, 1.2, 1.5 or 2.0 mm) by replacing XX in the Order Number with 10, 12, 15 or 20.

**OPTIONAL ACCESSORIES/REPLACEMENT PARTS**

Handle **2505** is for use with WPI manipulators. The smaller diameter handle **5444** is required for use with Narishige and Zeiss manipulators.

<b>2505</b>	1/4-in (6.3 mm) diameter handle
<b>5444</b>	3/16-in (4.8 mm) diameter handle
<b>GO1-100</b>	Replacement gasket 1.0 mm, Package of 100
<b>GO2-100</b>	Replacement gasket 1.2 mm, Package of 100
<b>GO3-100</b>	Replacement gasket 1.5 mm, Package of 100
<b>GO4-100</b>	Replacement gasket 2.0 mm, Package of 100
<b>1571</b>	Clear Silicone Rubber Sealant (-4.7 oz-)

# Precious Metal and Specialty Wire

Bare and coated metal wire for most laboratory applications



Micro coaxial cables (MAXxxxx) are ideal for microelectrode fabrication and construction of similar research tools. The dual shielding eliminates electrical interference caused by radio frequencies (RF), electrostatic and microphonics (e.g., bending and vibration). Available with single or dual (twin) conductors.

Teflon-coated stainless steel (type 304) wire (SSTxxxx) is available in 25-ft and 50-ft lengths. The Teflon coating is 150 micro-in. thick (4 µm). The Teflon coating is designed to reduce surface friction, only. It is not insulation.

Carbon wire (C3005) is a single 30-micron fiber of electrochemically activated carbon. This fiber is especially useful in micro-electrochemical experiments.

Platinum/Iridium wire — uncoated (PTxxxx) and Teflon-coated (PTTxxxx) — is an alloy of 90% platinum and 10% iridium, giving excellent tensile strength and corrosion resistance. Uncoated pure platinum wire (PTPxxx) is 99.95% pure. Indium wire (IN1003) is 99.99% pure, with a melting point of 156.4°C.

Annealed silver wire (AGWxxxx), 99.99% pure, is available in five diameters; three of those sizes are also available with a Teflon coating (AGTxxxx).

Tungsten wire (TGWxxxx), available in three diameters, is 99.95% pure.

Gold wire (AUWxxxx) is 99.99% pure.

Stainless steel wire (SSxxxx) is type 316.

## Wire Cutters and Scissors



### ORDERING INFORMATION

Catalog No.	Metal	Coating	AWG*	Diameter	Precut Length
AGT0510	Silver	Teflon	36	0.005 in. (0.125 mm) <sup>1</sup>	10 ft (3 m)
AGT0525	Silver	Teflon	36	0.005 in. (0.125 mm) <sup>1</sup>	25 ft (7.6 m)
AGT05100	Silver	Teflon	36	0.005 in. (0.125 mm) <sup>1</sup>	100 ft (30 m)
AGT1010	Silver	Teflon	30	0.010 in. (0.25 mm) <sup>1</sup>	10 ft (3 m)
AGT1025	Silver	Teflon	30	0.010 in. (0.25 mm) <sup>1</sup>	25 ft (7.6 m)
AGT10100	Silver	Teflon	30	0.010 in. (0.25 mm) <sup>1</sup>	100 ft (30 m)
AGT1510	Silver	Teflon	26-27	0.015 in. (0.38 mm) <sup>1</sup>	10 ft (3 m)
AGT1530	Silver	Teflon	26-27	0.015 in. (0.38 mm) <sup>1</sup>	30 ft (9.1 m)
AGW0510	Silver	—	36	0.005 in. (0.125 mm)	10 ft (3 m)
AGW0530	Silver	—	36	0.005 in. (0.125 mm)	30 ft (9.1 m)
AGW1010	Silver	—	30	0.010 in. (0.25 mm)	10 ft (3 m)
AGW1030	Silver	—	30	0.010 in. (0.25 mm)	30 ft (9.1 m)
AGW1510	Silver	—	26-27	0.015 in. (0.38 mm)	10 ft (3 m)
AGW1530	Silver	—	26-27	0.015 in. (0.38 mm)	30 ft (9.1 m)
AGW2010	Silver	—	24	0.020 in. (0.5 mm)	10 ft (3 m)
AGW2030	Silver	—	24	0.020 in. (0.5 mm)	30 ft (9.1 m)
AGW4010	Silver	—	18	0.040 in. (1.0 mm)	10 ft (3 m)
AUW0170	Gold	—	50	0.001 in. (0.025 mm)	70 ft (21 m)
AUW201	Gold	—	24	0.020 in. (0.5 mm)	1 ft (30 cm)
C3005	Carbon	—	49	0.0012 in. (30 µm)	5 ft (1.5 m)
PT1002	Platinum / Iridium	—	30	0.010 in. (0.25 mm)	2 ft (61 cm)
PT0402	Platinum / Iridium	—	38	0.004 in. (0.102 mm)	2 ft (61 cm)
PT0203	Platinum / Iridium	—	44	0.002 in. (0.051 mm)	3 ft (91 cm)
PT0110	Platinum / Iridium	—	50	0.001 in. (0.025 mm)	10 ft (3 m)
PTP101	Platinum	—	30	0.010 in. (0.25 mm)	1 ft (30 cm)
PTP201	Platinum	—	24	0.020 in. (0.5 mm)	1 ft (30 cm)
PTP401	Platinum	—	18	0.039 in. (1.0 mm)	1 ft (30 cm)
PTP406	Platinum	—	18	0.039 in. (1.0 mm)	0.5 ft (15.2 cm)
PTT0502	Platinum / Iridium	Teflon	36	0.005 in. (0.125 mm) <sup>1</sup>	2 ft (61 cm)
PTT0203	Platinum / Iridium	Teflon	44	0.002 in. (0.051 mm) <sup>1</sup>	3 ft (91 cm)
PTT0110	Platinum / Iridium	Teflon	50	0.001 in. (0.025 mm) <sup>1</sup>	10 ft (3 m)
SS31605	Stainless Steel	—	36	0.005 in. (0.125 mm)	50 ft (15.2 m)
SS31614	Stainless Steel	—	27	0.014 in. (0.36 mm)	30 ft (9.1 m)
SST30407-25	Stainless Steel	Teflon	33	0.007 in. (0.18 mm) <sup>3</sup>	25 ft (7.6 m)
SST30407-50	Stainless Steel	Teflon	33	0.007 in. (0.18 mm) <sup>3</sup>	50 ft (15.2 m)
TGW0325	Tungsten	—	40	0.003 in. (0.075 mm)	25 ft (7.6 m)
TGW0515	Tungsten	—	36	0.005 in. (0.125 mm)	15 ft (4.6 m)
TGW1510	Tungsten	—	26-27	0.015 in. (0.38 mm)	10 ft (3 m)

### MICROCOAXIAL CABLES

MAX3820	Tinned Cu Alloy	Coaxial		0.0173 in. (0.44 mm)	20 ft (6 m) <sup>4</sup>
MAX4020	Tinned Cu Alloy	Twin Coaxial	0.0158x0.024 in. (0.4x0.61 mm)		20 ft (6 m) <sup>5</sup>

\*Brown & Sharpe

<sup>1</sup> Plus 0.002 in. for Teflon coating

<sup>3</sup> Teflon adds 0.00015 in. (4 µm) to diameter

<sup>4</sup> Impedance: 50 ohm; capacitance: 95 pF/m; resistance: 5 ohm/m

<sup>5</sup> Impedance: 100 ohm; capacitance: 54 pF/m; resistance: 1.9 ohm/m

### ORDERING INFORMATION

504749	Ergonomic Micro-Shear Flush Cutters, 12.7 cm (5 in.) Micro-shear flush cutters for delicate wires
504750	Ergonomic Mini-Scissors, 12.7 cm (5 in.) for cutting fine or delicate items with a clean, square edge. Handle design is advantageous for users with arthritic hands.
504751	Ergonomic Micro-Shear Flush Cutters, 12.7 cm (5 in.) for delicate wires. ESD safe. Extra tough high carbon steel blades. Sized for smaller hands and maximum maneuverability

### WORLD PRECISION INSTRUMENTS


UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com Germany: +49 (0)6031 1602171 • wvide@wpi-europe.com • www.wpi-europe.com  
 Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • www.wpiinc.com China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net



# Ag/AgCl Half-Cells






## Sintering pellets with low resistance and high strength

New, improved sintered pellets with lower resistance and high strength. Stable and well balanced in the presence of current, these small and inexpensive half-cells are easy to work with as bath electrodes.



### ORDERING INFORMATION

<b>RC1</b>	Reference Cell with 1.5 m lead
<b>RC1T</b>	Reference Cell, 1.5 m lead, 2 mm pin
<b>RC2</b>	Reference Cell with 2.0 mm pin
<b>RC2F</b>	Reference Cell with female connector
<b>RC3</b>	Ref. Cell with epoxy body, 4.5 mm diam x 50 mm
<b>RC6</b>	Reference Cell with glass body, 1.5 mm diam x 50 mm
<b>EP05</b>	Ag/AgCl Electrode 0.5 mm diam x 20 mm
<b>EP08</b>	Ag/AgCl Electrode 0.8 mm diam x 20 mm
<b>EP1</b>	Ag/AgCl Electrode 1.0 mm diam x 3 mm
<b>EP2</b>	Ag/AgCl Electrode 2.0 mm diam x 4 mm
<b>EP4</b>	Ag/AgCl Electrode 4.0 mm diam x 1 mm
<b>EP8</b>	Ag/AgCl Electrode 8.0 mm diam x 1 mm
<b>EP12</b>	Ag/AgCl Electrode 12.0 mm diam x 1 mm
<b>3578</b>	Adapter Cable for Ag/AgCl Pellets



# Glass Capillaries

Quality glass, superior prices for microinjection/microelectrodes

## Features

- Quality borosilicate glass capillaries
- Large variety available, including fire polished, filaments, thin wall, specialty glass and multi-barrel

## Benefits


- Superior pricing
- Most glass orders ship within 48 hours

## Applications

- Microinjection
- Electrophysiology
- Patch clamp
- Fluid Handling

## Fire Polishing

Fire-Polished glass capillaries are easier to insert into microelectrode holders without damaging the gasket. More importantly, fire-polished glass won't scratch the chloridized wire used in a recording electrode. Fire-polishing does not affect the glass's mechanical or electrical properties.



## Making Uniform, Reproducible Microelectrodes


Borosilicate glass capillaries: Close dimensional tolerances assure microelectrode uniformity and reproducibility. Capillaries are available in 1, 2, 3, 5 and 7-barrel configurations, complete range of single barrel thin-wall sizes and a variety of special configurations. Capillaries with filaments contain a solid filament fused to the inner wall, which speeds filling of electrodes. Capillaries with or without inner filaments are available for making microelectrodes in a wide range of diameters.

## Filament Glass Capillaries


Single Barrel standard wall thickness capillaries are offered either with or without inner filaments for quick filling in a variety of lengths and diameters.

## Thin Wall Glass Capillaries

Thin Wall single barrel capillaries are offered both with or without inner filaments.



Specialty glass is also available. See page 105.



## ORDERING INFORMATION

Length	OD(mm)	ID(mm)	Filament	Fire-Polished	Quantity	Item	
Single Barrel Standard Borosilicate Glass	3 in. (76 mm)	1.0	0.58	✓	500	<b>1B100F-3</b>	
	3 in. (76 mm)	1.0	0.58		500	<b>1B100-3</b>	
	3 in. (76 mm)	1.2	0.68	✓	350	<b>1B120F-3</b>	
	3 in. (76 mm)	1.2	0.68		350	<b>1B120-3</b>	
	3 in. (76 mm)	1.5	0.84	✓	225	<b>1B150F-3</b>	
	3 in. (76 mm)	1.5	0.84		✓	300	<b>1B150-3</b>
	4 in. (100 mm)	1.0	0.58	✓	✓	500	<b>1B100F-4</b>
	4 in. (100 mm)	1.0	0.58		✓	500	<b>1B100-4</b>
	4 in. (100 mm)	1.2	0.68	✓	✓	400	<b>1B120F-4</b>
	4 in. (100 mm)	1.2	0.68			350	<b>1B120-4</b>
	4 in. (100 mm)	1.5	0.84	✓	✓	300	<b>1B150F-4</b>
	4 in. (100 mm)	1.5	0.84		✓	300	<b>1B150-4</b>
	4 in. (100 mm)	2.0	1.12	✓		125	<b>1B200F-4</b>
	4 in. (100 mm)	2.0	1.12		✓	200	<b>1B200-4</b>
	6 in. (152 mm)	1.0	0.58	✓		500	<b>1B100F-6</b>
	6 in. (152 mm)	1.0	0.58			500	<b>1B100-6</b>
6 in. (152 mm)	1.2	0.68	✓		350	<b>1B120F-6</b>	
6 in. (152 mm)	1.2	0.68			350	<b>1B120-6</b>	
6 in. (152 mm)	1.5	0.84	✓		225	<b>1B150F-6</b>	
6 in. (152 mm)	1.5	0.84			225	<b>1B150-6</b>	
6 in. (152 mm)	2.0	1.12	✓		125	<b>1B200F-6</b>	
6 in. (152 mm)	2.0	1.12			125	<b>1B200-6</b>	
Thin-Wall Single-Barrel Standard	3 in. (76 mm)	1.0	0.75	✓	500	<b>TW100F-3</b>	
	3 in. (76 mm)	1.0	0.75		500	<b>TW100-3</b>	
	3 in. (76 mm)	1.2	0.90	✓	✓	400	<b>TW120F-3</b>
	3 in. (76 mm)	1.2	0.90			350	<b>TW120-3</b>
	3 in. (76 mm)	1.5	1.12	✓		225	<b>TW150F-3</b>
	3 in. (76 mm)	1.5	1.12		✓	300	<b>TW150-3</b>
	4 in. (100 mm)	1.0	0.75	✓		500	<b>TW100F-4</b>
	4 in. (100 mm)	1.0	0.75		✓	500	<b>TW100-4</b>
	4 in. (100 mm)	1.2	0.90	✓		350	<b>TW120F-4</b>
	4 in. (100 mm)	1.2	0.90			350	<b>TW120-4</b>
	4 in. (100 mm)	1.5	1.12	✓		225	<b>TW150F-4</b>
	4 in. (100 mm)	1.5	1.12		✓	300	<b>TW150-4</b>
	6 in. (152 mm)	1.0	0.75	✓		500	<b>TW100F-6</b>
	6 in. (152 mm)	1.0	0.75		✓	500	<b>TW100-6</b>
	6 in. (152 mm)	1.2	0.90	✓	✓	400	<b>TW120F-6</b>
	6 in. (152 mm)	1.2	0.90			350	<b>TW120-6</b>
6 in. (152 mm)	1.5	1.12	✓		225	<b>TW150F-6</b>	
6 in. (152 mm)	1.5	1.12		✓	300	<b>TW150-6</b>	

Single barrel glass is Kimble N51A. All thin wall glass is Schott Duran 8330. Tolerance ±10%


## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com Germany: +49 (0)6031 1602171 • wptide@wpi-europe.com • www.wpi-europe.com  
 Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • www.wpiinc.com China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net

## Multi-Barrel Glass Capillaries

Multi-barrel configurations are designed especially for micro-iontophoresis. Because the capillaries are fused together during manufacturing, you will not need to twist them while pulling to seal the tips together. An inner filament in each barrel makes filling easy and fast.

Also see PolyFil for a novel way to connect multi-barrel pipettes



### ORDERING INFORMATION


#### Multi-Barrel Borosilicate Glass Tubing with Filaments

Length	Description	OD/ID (per barrel)	Filament	Qty	Item
4" (102mm)	2-Barrel	1.5/0.84 mm	✓	100	<b>2B150F-4</b>
4" (102mm)	3-Barrel	1.2/0.68 mm	✓	100	<b>3B120F-4</b>
4" (102mm)	5-Barrel	1.2/0.68 mm	✓	65	<b>5B120F-4</b>
4" (102mm)	7-Barrel	1.0/0.58 mm	✓	60	<b>7B100F-4</b>
4" (102mm)	7-Barrel	1.2/0.68 mm	✓	75	<b>7B120F-4</b>
6" (152mm)	2-Barrel	1.5/0.84 mm	✓	100	<b>2B150F-6</b>
6" (152mm)	3-Barrel	1.2/0.68 mm	✓	100	<b>3B120F-6</b>
6" (152mm)	5-Barrel	1.2/0.68 mm	✓	65	<b>5B120F-6</b>
6" (152mm)	7-Barrel	1.0/0.58 mm	✓	60	<b>7B100F-6</b>


## Special Configuration Glass Capillaries

### Septum Theta

Septum Theta offers superior cell impalement. The natural bevel resulting from the prominent spear-like projection of the septum gives microelectrodes a sharp, spear-point tip. This style has low resistance for use as a single microelectrode, and it can be used to make superior double-tipped microelectrodes with low transtip coupling. The natural bevel of Septum Theta also significantly increases the effective tip cross-section. As supplied, the width of the septum is approximately 0.2 mm; wall thickness is approximately 0.2 mm.



### Piggyback



Piggyback glass consists of a pair of borosilicate capillaries fused together during manufacturing. One barrel is larger than the other, and both have inner filaments for quick filling. Piggyback glass makes it simple to fabricate two-barrel electrodes with a significant tip diameter differential.

### ORDERING INFORMATION

#### Special Configuration Borosilicate Glass Tubing

Description	OD/ID (mm)	Length	Qty	Item
Septum Theta	1.5/1.02	6 in. (152 mm)	100	<b>TST150-6</b>
Piggyback	1.51/0.84 0.75/0.35	4 in. (102 mm)	50	<b>PB150F-4</b>
Piggyback	1.51/0.84 0.75/0.35	6 in. (152 mm)	50	<b>PB150F-6</b>

## Borosilicate Glass Rod

1.0 mm diameter — for making tools, probes, tips

### ORDERING INFORMATION

Borosilicate Glass Rod				
Description	OD (mm)	Length	Qty	Item
Glass Rod	1.0	4 in. (102mm)	500	<b>GR100-4</b>
Glass Rod	1.0	6 in. (152mm)	500	<b>GR100-6</b>

## Multi-barrel pipette coupling kit

Securely couple multi-barrel pipettes to a pressure source



### Features

- Complete kit (all in one) for making multi-barrel pipettes
- Secure coupling of multi-barrel micropipettes to a pressure source

### Benefits

- Allows for microinjection to multi-barrel micropipette independently for one pressure source
- Pressure safe and convenient luer lock connections

### Applications

- Multi-port microinjection

**PolyFil™** allows easy and secure coupling of a multi-barrel micropipette to a pressure source. Coupling is achieved by bonding temperature-resistant and flexible **MicroFil** to the capillary tube with hot melt adhesive. The luer end of each **MicroFil** is connected to PVC tubing (200 PSI rated). Kits also include a five-port manifold that allows use of a single **PV800** Series PicoPump to drive up to six micropipette barrels independently by switching on only the barrels to be injected. All connections are locking luers — pressure safe and convenient.

### ORDERING INFORMATION

<b>5440</b>	PolyFil™ Multi-Barrel Micropipette Coupling Kit <i>Includes: 1 pipette holder/handle, plastic; 7 pieces MF28G MicroFil; 7-pieces tubing with male luer lock fittings; 1 flow-thru manifold with five luer lock ports; 1 hot melt glue gun (110V only); 3 glue sticks.</i>
<b>13316</b>	Mini Glue Gun and (3) glue sticks


# Patch Clamp Capillary Glass

*Evaluated for best softening temperature, electrical properties, sealability, leachable components*

**PG52151-4** and **PG52165-4** are prepared from Schott #8250 glass (equivalent to Corning #7052), one of the most widely used patch clamping glasses. This is a specially formulated borosilicate glass with a softening temperature that is 110°C lower than regular borosilicate glass (Corning 7740, or Pyrex). It has excellent sealing properties for most cells. Electrical properties are also very good.

**PG10150-4** and **PG10165-4** are composed of Corning #0010 glass, a high lead content (22% PbO) glass. Its thermal and electrical performance is between the Schott #8250 and Corning #8161 glasses described above. It is much more economical than Corning #8161 glass. It has been found that this glass causes much less alteration in channel behavior than Corning #8161 and Schott #8250 glass (Furman and Tanaka, *Biophys. J.* 53, p287, 1988).

*Patch clamp capillaries do not have microfilaments.*



## ORDERING INFORMATION

### Patch Clamp Capillary Glass

Catalog#	Glass Type	OD/ID (mm)	Dielectric Constant	Softening Point °C	QTY
<b>PG52151-4</b>	#8250	1.5/1.0	4.9	720°	100
<b>PG52165-4</b>	#8250	1.65/1.1	4.9	720°	100
<b>PG10150-4</b>	#0010	1.5/0.75	6.7	625°	100
<b>PG10165-4</b>	#0010	1.65/1.1	6.7	625°	100

# Glass Handling Forceps

*Firmly hold your capillary glass without risk of breaking it*

Ever had difficulty picking up a glass capillary? Special tips on these forceps solve the problem, holding glass firmly without risk of breakage. They also keep the glass clean and avoid contamination from skin oils.



## ORDERING INFORMATION

**77020** Glass Handling Forceps

### WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com Germany: +49 (0)6031 1602171 • wpide@wpi-europe.com • www.wpi-europe.com  
 Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • www.wpiinc.com China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net

# Borosilicate Glass Micropipettes

Eliminate the cost and trouble of making your own micropipettes



## Features

- Schott Duran borosilicate glass
- 0.5 µm and smaller ID micropipettes include internal glass fiber for easy filling
- Tip inner diameter tolerance ±20%
- Short taper yields high strength
- Nominal length ≈ 50 mm
- OD:ID = 1.33:1
- Standard capillary outer diameters are 1.0 mm (thin-wall) or 1.14 mm
- Every pipette individually tested and inspected
- Vacuum sealed packs of 10

## Benefits

- Plain Shank or Luer Fittings

## Applications

- Injection of dyes or proteins into cells, oocytes or other biomedical laboratory applications

WPI can quickly supply your need for consistently sized pre-pulled glass micropipettes.

Tip diameters (ID) range from 0.1 to 15 µm.

Silanized Tips (Luer Shank) are available. Silanization waterproofs the glass to retard water when inserting into cell. This will not let the outside fluid run down the pipette and get inside so easily.

## ORDERING INFORMATION

Shank	Tip I.D.	Silanize Coating Length	Glass O.D.	Filament	Fire Polished	Catalog #
PLAIN	0.1 µm	—	1.0 mm Thin-Wall	Yes	No	TIP01TW1F
	0.2 µm	—	1.0 mm Thin-Wall	Yes	No	TIP02TW1F
	0.3 µm	—	1.0 mm Thin-Wall	Yes	No	TIP03TW1F
	0.4 µm	—	1.0 mm Thin-Wall	Yes	No	TIP04TW1F
	0.5 µm	—	1.0 mm Thin-Wall	Yes	No	TIP05TW1F
	1 µm	—	1.0 mm Thin-Wall	No	Yes	TIP1TW1
	2 µm	—	1.0 mm Thin-Wall	No	Yes	TIP2TW1
	5 µm	—	1.0 mm Thin-Wall	No	Yes	TIP5TW1
	10 µm	—	1.0 mm Thin-Wall	No	Yes	TIP10TW1
	10 µm	—	1.14 mm A203W glass *	No	Yes	TIP10XV119
	—	—	1.0 mm Thin-Wall	No	Yes	TIP10LT Long Taper †
	—	—	1.0 mm Thin-Wall	Yes	Yes	TIP10FLT Long Taper †
	—	—	1.5 mm Thin-Wall	No	Yes	TIP15LT Long Taper †
	—	—	1.5 mm Thin-Wall	Yes	Yes	TIP15FLT Long Taper †
	—	30 µm	—	1.0 mm Thin-Wall	No	Yes
LUER	0.1 µm	—	1.0 mm Thin-Wall	Yes	—	TIP01TW1F-L
	0.2 µm	—	1.0 mm Thin-Wall	Yes	—	TIP02TW1F-L
	0.3 µm	—	1.0 mm Thin-Wall	Yes	—	TIP03TW1F-L
	0.5 µm	—	1.0 mm Thin-Wall	Yes	—	TIP05TW1F-L
	1 µm	—	1.0 mm Thin-Wall	No	—	TIP1TW1-L
	2 µm	—	1.0 mm Thin-Wall	No	—	TIP2TW1-L
	5 µm	—	1.0 mm Thin-Wall	No	—	TIP5TW1-L
	10 µm	—	1.0 mm Thin-Wall	No	—	TIP10TW1-L
	30 µm	—	1.0 mm Thin-Wall	No	—	TIP30TW1-L
	LUER/SILANIZED	5 µm	1 inch	1.0 mm Thin-Wall	No	—
5 µm		2 inch	1.0 mm Thin-Wall	No	—	TIP5TW1LS02
10 µm		1 inch	1.0 mm Thin-Wall	No	—	TIP10TW1LS01
10 µm		2 inch	1.0 mm Thin-Wall	No	—	TIP10TW1LS02
30 µm		1 inch	1.0 mm Thin-Wall	No	—	TIP30TW1LS01
30 µm		2 inch	1.0 mm Thin-Wall	No	—	TIP30TW1LS02

\* 10 µm (ID), 1.14 mm capillary pipettes are for use in WPI's Nanoliter 2010.

† Long Taper micropipettes are pulled with a 12-15 mm taper which the customer cuts back to obtain the desired tip diameter.

## µTIP SAMPLER ASSORTMENTS

TIPMIX01-05	Two each, 0.1, 0.2, 0.3, 0.4, 0.5 µm ID, plain shank
TIPMIX05-10	Two each, 0.5, 1, 2, 5, 10 µm ID, plain shank
TIPMIX01-05-L	Two each, 0.1, 0.2, 0.3, 0.4, 0.5 µm ID, Luer
TIPMIX05-10-L	Two each, 0.5, 1, 2, 5, 10 µm ID, Luer



## Micropipette Storage Jar

Protect your micropipettes until ready to use them

Stores up to 30 micropipettes, filled or unfilled, up to three inches in length. A gentle sliding action inserts or removes pipettes without damage to the delicate tips.

## ORDERING INFORMATION

E210	Storage Jar for 1.0 mm OD Micropipettes
E212	Storage Jar for 1.2 mm OD Micropipettes
E215	Storage Jar for 1.5 mm OD Micropipettes
E220	Storage Jar for 2.0 mm OD Micropipettes

## REPLACEMENT PARTS

1965	Foam Ring for 0.75 - 1.0 mm glass
1966	Foam Ring for 1.2 - 1.5 mm glass
1967	Foam Ring for 2.0 mm glass

# Microelectrode Puller

*A compact, versatile and reliable workhorse*



PUL-1000

## Features

- Program sequences up to four steps
- Store up to 95 programs in memory for easy recall
- Supplied with a vial of capillary glass

## Benefits

- Tempered glass cover to reduce the effects of humidity on puller reproducibility
- Switchable power supply ensures that line voltage fluctuations don't affect reproducibility

## Applications

- Pull your own microelectrodes and micropipettes

**PUL-1000** is a microprocessor controlled horizontal puller for making glass micropipettes or microelectrodes used in intracellular recording, patch clamp studies, microperfusion or microinjection. The puller was designed with tight mechanical specifications and precision electronics for complete control of the pulling process and accurate reproducibility. It offers programmable sequences of up to four steps with heating, force, movement and cooling time. This allows graduated cycles for applications like patch clamp recording.

This puller is a reasonably priced, compact, versatile and reliable workhorse. The microprocessor, combined with the LCD display, makes the **PUL-1000** easy to use.

### Tempered Glass Cover

The cover of pulling chamber is made with tempered glass to minimize the humidity effect on the reproducibility of pulled pipettes.

### Switchable Power Supply

**PUL-1000** has a high quality switching power supply for use anywhere in the world without worry about the line voltage differences. Pulling

reproducibility is unaffected by line voltage fluctuation. Heating voltage can be controlled to within 0.1% accuracy even when line voltage fluctuates from 90 to 240 VAC.

## References

**Plautz, C. Z., Williams, H. C., & Grainger, R. M.** (2016). Functional Cloning Using a Xenopus Oocyte Expression System. *Journal of Visualized Experiments*, (107), e53518–e53518. <http://doi.org/10.3791/53518>

**Komarova, Y., Peloquin, J., & Borisy, G.** (2011). Components of a microinjection system. *Cold Spring Harbor Protocols*, 2011(8), 935–9. <http://doi.org/10.1101/pdb.ip27>

## PUL-1000 SPECIFICATIONS

HEATER ELEMENT	Platinum/Iridium
PULLING FORCE	Solenoid, adjustable
CAPILLARY OD RANGE	1.0–2.0 mm
MAXIMUM CAPILLARY LENGTH	170 mm
MINIMUM CAPILLARY LENGTH	55 mm
MEMORY SETS	95
POWER	90–240 VAC, 50/60 Hz
DIMENSIONS	34 x 24 x 12 cm
SHIPPING WEIGHT	15 lb.
REPLACEMENT FILAMENTS	2.5 mm Square Box Filament, 2.5 mm wide (13834)

## ORDERING INFORMATION

<b>PUL-1000</b>	Micropipette Puller
<b>13834</b>	Filaments, Platinum/Iridium, 2.5mm wide 2.5mm Square Box
<b>504951</b>	Filaments, Platinum/Iridium, 4.5mm wide Trough

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • [wpiuk@wpi-europe.com](mailto:wpiuk@wpi-europe.com) • [www.wpi-europe.com](http://www.wpi-europe.com)  
Brazil: 011 55 13 40629703 • [info@brazil.wpiinc.com](mailto:info@brazil.wpiinc.com) • [www.wpiinc.com](http://www.wpiinc.com)

Germany: +49 (0)6031 1602171 • [wptide@wpi-europe.com](mailto:wptide@wpi-europe.com) • [www.wpi-europe.com](http://www.wpi-europe.com)  
China: +86 21 6888 5517 • [chinasales@china.wpiinc.com](mailto:chinasales@china.wpiinc.com) • [www.wpiinc.net](http://www.wpiinc.net)

# Micropipette Puller for Quartz and Glass

*CO<sub>2</sub> laser-based heat source for pulling quartz pipettes*

new

## Features

- Capable of pulling quartz, borosilicate and aluminosilicate glass
- Fully programmable, including heating filament characteristics
- Pulls electrodes with tip diameters smaller than 0.03  $\mu\text{m}$
- Preprogrammed sample programs for intracellular and patch pipettes

## Benefits

- The laser has no melting point limit as with conventional metal filaments and cannot be burned out
- Optimized velocity sensing circuit for maximized sensitivity and reproducibility
- Operating life of the CO<sub>2</sub> laser is expected to be in excess of ten years with normal use, after which the laser can be refurbished
- Individual programs can be write-protected in order to secure them from inadvertent changes
- The total time that the heat is on during the pull is displayed for improved program development and troubleshooting
- A date and time stamp is displayed to show the last time that a program has been changed



## Applications

- Patch clamp, single isolated and whole cell
- Intracellular recording
- Microinjection
- Nano probe research
- SECM
- P-2000/F is ideal for applications such as nanospray and NSOM

The **SU-P2000** micropipette puller represents a significant advance in the technology of fabrication of micropipettes, optical fiber probes, and nanospray tips. The **SU-P2000** integrates a CO<sub>2</sub> laser-based heat source with the technology derived from our extensive experience with conventional pullers. This system offers capabilities unmatched by other pullers.

While the **SU-P2000** is suitable for working with most conventional glasses, its primary advantage is the ability to work with quartz glass (fused silica). Quartz offers superior material properties for a variety of research applications. Quartz is stronger than other glasses and can facilitate penetration through tough tissues which would normally break conventional pipettes. For applications requiring a low noise glass, users will find that quartz is the lowest noise glass available. Quartz contains none of the metals used in conventional glasses. Optically, quartz is virtually free from fluorescence when illuminated.

A CO<sub>2</sub> laser was selected as the heat source for the **SU-P2000** for several reasons:

- The nominal emission wavelength of the laser approximates the resonant frequency of the SiO<sub>2</sub> lattice in glass. Thus, quartz and other conventional glasses can be melted when the appropriate laser power is supplied.
- Laser heat is clean and leaves no metal residue on the pipette as do conventional heating filaments.
- Laser heat can be turned off instantly, leaving no residual filament heat.
- The user can program the amount and distribution of heat supplied to the glass.
- Laser heat source means there are no filaments to burn out or replace.

The **SU-P2000** can store up to 100 separate programs, with each program consisting of up to 8 command lines. Programmable parameters include: laser power level, scan width, trip velocity, delay/laser on time, and hard pull strength.

One important consideration for the use of the **SU-P2000** is the diameter of the glass used. The **SU-P2000** is designed to produce even heating on glass up to 1.2 mm in outside diameter. Larger diameter glass can be used with the **SU-P2000** (up to 1.5 mm quartz and 1.8 mm conventional glass), but the performance is best with glass that is 1.2 mm diameter or less.

## SPECIFICATIONS

DIMENSIONS	30 x 14.25 x 13.25 in. (76 x 36 x 33.5 cm)
WEIGHT	90 lbs. (41 kg)
ELECTRICAL	115/230 V, 50/60 Hz
*Patent No.4,600,424	
CLASS 1 LASER PRODUCT	

## ORDERING INFORMATION

**SU-P2000** Micropipette Puller

## References

- Munoz, J.L. and Coles, J.** Quartz micropipettes for intracellular voltage microelectrodes and ion selective microelectrodes. *Journal of Neuroscience Methods*: 22:57-64, 1987.
- Rae, J.L. and Levis, R. A.** A method for exceptionally low noise single channel recordings. *European Journal of Physiology - Pflügers Archiv*: 420:618-620, 1992.
- Zuazaga, C. and Steinacker, A.** Patch-clamp recording of ion channels: Interfering effects of patch pipette glass. *News in Physiological Science*: 5:155-159, 1990.
- Levis, R.A. and Rae, J. L.** The use of quartz patch pipettes for low noise single channel recording. *Biophysical Journal*: 65:1666-1677, 1993.

# Next Generation Micropipette Puller

*Flaming/Brown style puller with color touch screen control*

## Features

- Color touch-screen display
- Help topics and error detection
- Copy and paste functions for writing new programs
- Record of last two pull results
- Two symmetrical pipettes with each pull
- Self-contained air supply with filtration system and humidity control chamber
- Memory storage for up to 100 programs
- Two cooling modes: time and delay

## Benefits

- Safe heat mode to protect and extend filament life
- Pre-heat mode improves stability
- Line repeat mode simplifies multi-line programming
- Pipette Cookbook program directory
- Glossary with micropipette and puller terminology
- Jaw temperature sensor helps define ideal pulling conditions
- Ramp test more easily accessed and can now be stored and referenced within each program. Helps to establish program heat settings and protect filament.

## Applications

- Patch pipettes
- Sharp electrodes
- Pronuclear injection
- Zebrafish injection
- Insect egg injection
- Aspiration pipettes

The **SU-P1000** micropipette puller was developed through years of experience with the Flaming/Brown style micropipette pullers and infused with leading edge technology. The most obvious new feature is the color touchscreen display that provides an intuitive and full-featured interface.

The extensive library of programs found in the popular Sutter Cookbook has been incorporated into the puller and is available to the user via the touchscreen display. You need only specify the glass, filament, and type of pipette you require and a suitable program will be identified and available for installation. This takes the guesswork out of pipette pulling and simplifies programming.

The Pre-heat mode actively heats and controls the jaw temperature and assures that the jaws have reached a specific temperature before the glass is pulled. This can increase the stability of the program from pull to pull. Copy and paste functions assure that programs can be easily written and the line repeat mode simplifies multi-line programming. A safe heat mode is an additional feature that helps the user avoid using heat settings that might damage or burn out the filament. When the safe heat mode is turned on, the puller will "check" the installed heat and alert the user when a given heat value is too low or too high in relation to the ramp value. User notes can be added to each program for annotating important information.

## New features

- Diagnostic testing of all puller components
- Built-in error detection of air pressure loss or filament burnout
- Easy access to ramp test



- Measurement of jaw temperature
- Access to previous pull results with the heat on times for each cycle of the program
- Help topics are pre-loaded to assist with on-site troubleshooting and the built-in glossary includes text, pictures, and diagrams explaining the terms used in micropipette fabrication
- Rotary dial is offered as an alternative to the keypad for numerical entry

## SU-P1000 SPECIFICATIONS

DIMENSIONS	21 x 14 x 12" (53 x 36 x 30 cm)
WEIGHT	41 lb. (18.59 kg)
POWER	115/230 V, 50/60 Hz

## ORDERING INFORMATION

**SU-P1000** Pipette Puller

### OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>505094</b>	2 mm Box Filament, 1.5 mm wide
<b>505095</b>	2 mm Box Filament, 2 mm wide
<b>505096</b>	2 mm Box Filament, 3 mm wide
<b>505097</b>	2.5 mm Box Filament, 4.5 mm wide
<b>13834</b>	2 mm Box Filament, 2.5 mm wide
<b>505098</b>	3 mm Box Filament, 1.5 mm wide
<b>14074</b>	3 mm Box Filament, 2 mm wide
<b>506122</b>	3 mm Box Filament, 3 mm wide
<b>505099</b>	Custom Filament
<b>505101</b>	3 mm Loop Filament, 1.5 mm wide
<b>505103</b>	3 mm Loop Filament, 2 mm wide
<b>505105</b>	3 mm Loop Filament, 2.5 mm wide
<b>505106</b>	Fire Polishing Spacer
<b>505108</b>	Trough Filament, 1.5 mm wide
<b>13835</b>	Trough Filament, 2 mm wide
<b>13836</b>	Trough Filament, 3 mm wide
<b>505109</b>	Trough Filament, 4.5 mm wide

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com Germany: +49 (0)6031 1602171 • wpide@wpi-europe.com • www.wpi-europe.com  
Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • www.wpiinc.com China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net



# Flaming/Brown Pipette Puller

Create custom micropipettes, patch pipettes and microinjection needles

## Features

- Environmental chamber for humidity control
- Programmable air pressure
- Memory storage for up to 100 programs
- Write protection and date stamp for each program
- Two symmetrical pipettes with each pull
- Two cooling modes: time and delay
- Vacuum fluorescent display

## Benefits

- Preprogrammed sample programs for intracellular and patch pipettes. Special programming on request
- Ramp test — self test for establishing program heat settings when a new filament or glass is introduced
- Constant current power supply for filament and pull solenoid
- Looping pull cycle for fabrication of patch type micropipettes
- Self-contained air supply with filtration system and humidity control
- Consistent and reliable electrodes with tip diameters less than 0.1  $\mu\text{m}$
- Control over the time and pressure at which the air is delivered
- Optimized velocity sensing circuit for maximized sensitivity and reproducibility
- Quality control, SEM photograph of a tip pulled with each puller. Criterion is tip measurement less than 0.1  $\mu\text{m}$  and typically is  $\sim 0.06 \mu\text{m}$ .

## Applications

- Patch pipettes
- Sharp electrodes
- Pronuclear injection
- Zebrafish injection
- Insect egg injection
- Aspiration pipettes

The **SU-P97** Flaming/Brown type micropipette puller is ideal for fabricating micropipettes, patch pipettes and microinjection needles. While retaining many of the features of earlier models, the **SU-P97** offers improvements in mechanical, electronic and software design. The result is better control of the pulling process and a higher degree of reproducibility.

The **SU-P97** combines a proven mechanical system with a sophisticated, programmable microprocessor controller. This programmable control of the pulling parameters allows the investigator to design application specific pipettes from a wide range of glass compositions and sizes.

A number of other features have been incorporated in the design of the **SU-P97**. Most apparent is the environmental chamber which surrounds the heating filament. This environmental chamber is designed to minimize the effect of changing humidity on the reproducibility of pulled pipettes. A 25% increase in power over the P-87 allows for the use of larger heating filaments, larger diameter glass and multi-barreled glass. The metal jaws that clamp the heating filament have also been redesigned to minimize heat retention. A gas delivery mode switch provides for extended cooling for large diameter and multi-barreled glass. A spring-loaded clamping mechanism has been added for easier loading of glass. A vacuum fluorescent display has been added that allows easy viewing.

Software improvements on the **SU-P97** include a display of the total heat-on time to assist in program development and troubleshooting. Up to 100 programs can now be written and stored in memory, which makes the **SU-P97** suitable for multiple users. These programs can now be write-protected, adding security to prevent programs from being changed or altered inadvertently. The display shows the last date and time the



program was written or edited. In addition, the air pressure is included as a programmable parameter.

## SU-P97 SPECIFICATIONS

DIMENSIONS	21 x 16 x 12 (53 x 40.6 x 30 cm)
WEIGHT	50 lb. (23 kg)
ELECTRICAL	115/230 V, 50/60 Hz

## ORDERING INFORMATION

**SU-P97** Pipette Puller

### OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>505094</b>	2 mm Box Filament, 1.5 mm wide
<b>505095</b>	2 mm Box Filament, 2 mm wide
<b>505096</b>	2 mm Box Filament, 3 mm wide
<b>505097</b>	2.5 mm Box Filament, 4.5 mm wide
<b>13834</b>	2 mm Box Filament, 2.5 mm wide
<b>505098</b>	3 mm Box Filament, 1.5 mm wide
<b>14074</b>	3 mm Box Filament, 2 mm wide
<b>506122</b>	3 mm Box Filament, 3 mm wide
<b>505099</b>	Custom Filament
<b>505101</b>	3 mm Loop Filament, 1.5 mm wide
<b>505103</b>	3 mm Loop Filament, 2 mm wide
<b>505105</b>	3 mm Loop Filament, 2.5 mm wide
<b>505106</b>	Fire Polishing Spacer
<b>505108</b>	Trough Filament, 1.5 mm wide
<b>13835</b>	Trough Filament, 2 mm wide
<b>13836</b>	Trough Filament, 3 mm wide
<b>505109</b>	Trough Filament, 4.5 mm wide

# Digital Microforge

*Microforging, Micropipette Calibration and Microinjection — in a single device!*

## Features

- Microprocessor-controlled microforge
- Digital signal processor technology precisely controls the polish heating time
- Unique digital pneumatic pressure feature polishes the tip without changing the size
- **W30S-LED** Microscope (Optional)

## Benefits

- Digital Signal Processor (DSP) technology
- Complete system package available
- Kohler illuminator and Abbe condenser for less glare and sharper images.
- Pneumatic pressure polishing that allows the preparation of blunt tips without change of tip ID
- Heating filament is attached to the microscope objective so they move together
- Pipette holder sits on the microscope stage to simplify the locating and polishing of the pipette

## Applications

- Polishing patch pipettes
- Microforging holding pipettes
- Microforging beveled injection pipettes
- Pipette tip calibration and microinjection


The **DMF1000** is a microprocessor-controlled microforge offering unmatched performance. Designed for fabrication of both small patch clamp glass pipettes and larger injection pipettes, the **DMF1000** should find many uses in the laboratory. The **DMF1000** is based on a design similar to that first used in WPI's extremely popular microforge model, the **MF200**. The extensive improvements incorporated into the **DMF1000** greatly increase its versatility and performance, making it one of the most powerful microforges on the market.

### Digital Signal Processor (DSP) Technology

The **DMF1000** is powered by the latest digital signal processor (DSP) technology. A digital timer is used to precisely control the polish heating time. Ten memories can be used to store settings of the heating power and heating duration. All of the settings are controlled and displayed digitally for better accuracy and reproducibility. Two different operating modes are provided: Manual and Auto. In the Manual mode, the DSP will memorize the duration of the time that is used to achieve a desired polishing. In Auto mode, the heat will be applied for the duration of the timer setting.

### Complete System Available


The **DMF1000** system includes a specially configured WPI model **W30S-LED** research grade compound microscope (optional) equipped with a high quality metallurgic 40x long-working distance objective and a pair of 10x eyepieces. The long working distance objective reduces the danger of damage to the objective lens during the heating process.




Fire Polishing




Large Tip Sharpening (contact stretching)




Tip Sealing




Tip Reducing (holding pipettes)



Tip Bending



Carbon Fiber Sealing in Plastic Sensor



DMF1000-M1

*Filament Holder mounts directly to objective to provide precise control of heating element position.*

### Kohler Illuminator and Abbe Condenser

Other benefits of the **DMF1000** design include the use of a Kohler illuminator and Abbe condenser, which provide the reduced glare and sharper image contrast necessary when polishing pipettes as small as half a micron (0.5  $\mu\text{m}$ ) in diameter.

### Pressure Polishing

The **DMF1000** incorporates a unique digital pneumatic pressure feature that enables pressurized air to be delivered through the pipette during fire polishing. In the fabrication of patch pipettes, the pressurized air can be used to blunt the taper at the pipette tip without changing the size of the tip opening. This reduces electrical resistance of the tip, leading to lower noise during patch-clamp recordings (Goodman & Lockery, 2000).

### The Heating Filament

With a conventional microforge often the most difficult and time-consuming part of using a high magnification objective is being able to move both the heating filament and the pipette into the same viewing area. Finding and moving both the heating filament and the pipette without collision can be a challenge. However, this difficulty is eliminated with the **DMF1000** because the heating filament is directly attached to the

microscope's objective. Hence it can be easily adjusted to any position within the viewing area.

The low heat capacity and low thermal coefficient of linear expansion of the filaments are key design features. The low heat capacity of the filament allows it to reach fire-polishing temperatures without excessive heat. This permits you to bring the pipette tip close to the filament during polishing without fear of collapsing the pipette tip. Low heat capacity eliminates the need for an auxiliary air-cooling system. The low coefficient of expansion characteristic of the filament ensures minimal displacement of the filament during heating. This feature eliminates much of the guesswork out of tip placement in relation to the filament.

Two different heating filaments are provided to accommodate various applications. The H5 filament is large gauge and can be reformed into a "U" for fabrication of pipettes, air forming of patch pipettes and other applications. The H4 is a smaller gauge filament and is ideal for polishing patch clamp pipettes.

### Pipette Holder Sits on the Microscope Stage

The pipette rests on a specially designed holder that sits on top of the microscope stage. The position of the pipette, relative to the heating filament, is controlled by the (X, Y, Z) adjustment of the stage. This unique design makes locating and polishing the pipette extremely easy. The stage of the microscope has a high quality rail that gives precise, smooth and stable control of the pipette's movement. This configuration also eliminates the need and expense of an additional micromanipulator to control pipette movement.

## Typical applications of the DMF1000

### Polishing the Patch Pipettes

Proper fire polishing of patch pipettes is the single most important factor for forming a stable giga-seal in patch clamp recording. This is even more important than the type of glass capillary used. Difficulties often arise in forming giga-seals because the polishing of patch pipettes using a conventional low magnification microforge is inadequate. However, the DMF1000 uses a 40X long-working distance objective. Pipette polishing is much more accurately controlled. Both whole cell patch pipettes and single channel patch pipettes can be conveniently polished to the highest quality and reproducibility achievable with any microforge.

### Microforging Holding Pipettes

A holding pipette with a large blunt tip and a small opening is used to hold a floating cell in place prior to microinjection by applying suction to the rear of the pipette. The procedure for making holding pipettes involves three steps: squaring off, large bore flame polishing and tip reducing. These steps are accomplished with a larger heating filament.

### Microforging Beveled Injection Pipettes

Occasionally, a beveled large bore pipette is not sharp enough to penetrate a cell without damaging the area around the pipette. With the DMF1000 and the large heating filament, a sharp point can be formed on the beveled tip to assist the penetration of the cell. This process is referred to as contact stretching.

### Pipette Tip Calibration & Microinjection

The integrated digital pneumatic pressure system can be used to calibrate the precise diameter (I.D.) of a micropipette tip, based on a technique described previously (Hagag & Randolph 1990, Bowman & Ruknudin 1999). The pressure system can also be used separately as a simple but highly accurate controller for microinjection applications.

## References

Wu, Z.-Z., Chen, S.-R., & Pan, H.-L. (n.d.). Differential Sensitivity of N-and P/Q-Type Ca<sup>2+</sup> Channel Currents to an Opioid in Isolectin B 4 -Positive and -Negative Dorsal Root Ganglion Neurons. <http://doi.org/10.1124/jpet.104.073429>



CE

## SPECIFICATIONS

AC POWER MODULE	100-240 VAC 50/60 Hz
TIMER RANGE (for heater & timer)	0.01 to 360 sec
NUMBER OF MEMORYS	10
PRESSURE ADJUSTING RANGE	0.5 – 60 PSI (3.5 – 414 kPa)
PRESSURE RESOLUTION	0.1 PSI (0.7 kPa)
FILAMENTS: H4	Small filament for working with 40x long working distance objective.
FILAMENTS: H5	Large filament for working with 10x objective. Filament adjustment assembly provided for both objectives.
HEATER AND TIMER CONTROL	Auto or Manual via Pushbutton, TTL, or Optional Foot switch.
DIMENSIONS: Control Unit	4 x 7 x 17 in. (10.2 x 17.8 x 4.8 cm)
SHIPPING WEIGHT	4 lb. (1.8 kg)
MICROSCOPE	See W30S, page 187
MICROSCOPE: SHIPPING WEIGHT	16 lb. (7.3 kg)

## ORDERING INFORMATION

<b>DMF1000-1</b>	Microforge, w W30S-LED Microscope (110 V)
<b>DMF1000-2</b>	Microforge, w W30S-LED Microscope (220 V)
<b>DMF1000-M1</b>	Microforge without microscope (110 V)
<b>DMF1000-M2</b>	Microforge without microscope (220 V)

\*Above DMF1000 microforges include 40X long working distance objective

## OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>500329</b>	25x Long Working Distance Objective, 5 mm 0.5 NA
<b>13142</b>	Optional foot switch
<b>800292</b>	40x Long Working Distance Objective, 3mm 0.25 NA
<b>503513</b>	21 mm 10X Eyepiece with 100/10 reticle
<b>DMF1000-H5</b>	Replacement heating filament (large gauge)
<b>MF200-H4</b>	Replacement Heating filament (small gauge)
<b>75050</b>	Replacement Micropipette Slide
<b>75040</b>	Replacement Filament Cable

# Analog Microforge

*Sometimes the simplest designs work best*

## Features

- Simple, reliable and economically priced
- Analog temperature controller
- W30S-LED microscope (optional)

## Benefits

- Includes 40x long-working distance objective and 10x eyepiece
- Kohler illuminator and Abbe condenser for less glare and sharper images.

## Applications

- Patch pipette tip polishing
- Micropipette tip size reduction
- Contact stretching in in vitro fertilization pipette production

The **MF200** Microforge is a versatile instrument designed specifically for the fabrication of glass micropipettes and other related tools. The system was developed in collaboration with Dr. Ming Li of the Department of Pharmacology, University of South Alabama. The **MF200** is simple, reliable and priced economically.

### 40x LWD Objective Included

The **MF200** system includes: An easy to use analog temperature controller, a specially configured WPI model **W30S-LED** research grade compound microscope (optional), 40x long-working distance objective and 10x eyepiece. 40x magnification is essential when polishing pipettes as small as half a micron (0.5  $\mu\text{m}$ ) in diameter. Compared to a conventional 40x objective, the long working distance objective reduces the danger of damage to the pipette and/or objective lens during the polishing process.

### Kohler Illuminator

It is the only commercial microforge using the Kohler illuminator and Abbe condenser for illumination.

This provides less glare and sharper image of the pipette than the frosted glass illuminator, which is used on other commercial microforges.

## References

Guillou, L., Babataheri, A., Puech, P.-H., Barakat, A. I., & Husson, J. (2016). Dynamic monitoring of cell mechanical properties using profile microindentation. *Scientific Reports*, 6, 21529. <http://doi.org/10.1038/srep21529>

Vasaskas, A. A., Chen, H., Wu, S., & Cioffi, D. L. (2014). The serine-threonine phosphatase calcineurin is a regulator of endothelial store-operated calcium entry. *Pulmonary Circulation*, 4(1), 116–27. <http://doi.org/10.1086/675641>



## SPECIFICATIONS

AC POWER MODULE	100-240 VAC 50/60 Hz
FILAMENTS (3)	H2, H3, H4
FILAMENT ON	Pushbutton Controlled or Optional Foot Switch Controlled
FILAMENT ADJUSTMENT ASSEMBLY	For 40x and 25x Long-Working Distance Objectives: mounts on objective
FILAMENT ADJUSTMENT ASSEMBLY: OBJECTIVE	40x Long-Working Distance (3 mm)
FILAMENT ADJUSTMENT ASSEMBLY: OPTIONAL	25x Long-Working Distance (5 mm)
EYEPIECE	10x (pair)
EYEPIECE: RETICLE (10x eyepiece only) (OPTIONAL)	1.25 $\mu\text{m}$ /division (at 40x), 0-90° Angle at 5°/division
EYEPIECE: OPTIONAL EYEPIECE	15x (pair)
GLASS HOLDER	Mounts on Microscope Stage
DIMENSIONS: Control Unit	4 x 7 x 1.875 in. (10.2 x 17.8 x 4.8 cm)
SHIPPING WEIGHT	3 lb. (1.4 kg)
MICROSCOPE: Note	See W30S-LED
MICROSCOPE: SHIPPING WEIGHT	16 lb. (7.3 kg)

## ORDERING INFORMATION

<b>MF200-1</b>	Microforge System w. W30S-LED Microscope (110 V)
<b>MF200-2</b>	Microforge System w W30S-LED Microscope (220 V)
<b>MF200-M1</b>	MF 200 without microscope (110 V)
<b>MF200-M2</b>	MF 200 without microscope (220 V)

\*Above MF200 microforges include 40X long working distance objective

## OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>500329</b>	25x Long-Working Distance Objective (fits most microscopes with a 160 mm Focal Length)
<b>13142</b>	Optional foot switch
<b>MF200-H3</b>	Replacement heating filament (medium gauge)
<b>MF200-H4</b>	Replacement heating filament (small gauge)
<b>75090</b>	Filament Adjustment Assembly, 22mm OD Objectives
<b>75050</b>	Replacement Micropipette Slide
<b>75040</b>	Replacement Filament Cable

## FEATURE COMPARISON

	MF200	DMF1000
W30S-LED Microscope	✓	✓
40x Long Working Distance Objective	✓	✓
Analog Controller	✓	✓
Digital Controller		✓
Pressurized Air Control		✓
Microinjection Capability		✓
Optional Foot Switch	✓	✓
Memory		✓
Auto-sense of Filament Type		✓
Digital Temperature Control		✓

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • [wpiuk@wpi-europe.com](mailto:wpiuk@wpi-europe.com) • [www.wpi-europe.com](http://www.wpi-europe.com)  
 Brazil: 011 55 13 40629703 • [info@brazil.wpiinc.com](mailto:info@brazil.wpiinc.com) • [www.wpiinc.com](http://www.wpiinc.com)

Germany: +49 (0)6031 1602171 • [wptide@wpi-europe.com](mailto:wptide@wpi-europe.com) • [www.wpi-europe.com](http://www.wpi-europe.com)  
 China: +86 21 6888 5517 • [chinasales@china.wpiinc.com](mailto:chinasales@china.wpiinc.com) • [www.wpiinc.net](http://www.wpiinc.net)

# Micropressure System

## Measure hydrostatic pressure in small vessels and oocytes

### Features

- Utilizes a liquid filled micropipette (2–5  $\mu\text{m}$  tip opening) for sensing pressure
- Pressure range from 1–100 mmHg (pressure range to 350 mmHg is available)
- Lower limit 1 mmHg (133 Pa)
- Includes calibration/test chamber
- Tubing and fittings for interconnecting system sub-components are provided
- Pressure in the pipette can be manually set to positive or negative relative to the outside
- Probe holder for mounting on micromanipulator is included
- 10 pre-pulled pipettes are included

### Benefits

- Measures biological pressures in very small liquid (aqueous) filled spaces
- Pre-configured pipettes are available for convenience

### Applications

- Pressure in kidney tubules
- Embryonic blood pressure
- Intracellular pressure
- Mouse intraocular pressure

The **900A** system is designed to measure liquid pressures dynamically in aqueous biological micro-environments, such as in kidney tubules or intracellular pressures. A liquid filled micropipette is used as a pressure probe, and pressure external to the pipette is measured at the tip. The outside diameter of the micropipette tip typically measures between 2–7  $\mu\text{m}$ . Pressure measurement is achieved by monitoring the pipette's electrical resistance. The resistance changes according to changes in the pressure outside the pipette tip via displacement of an electrolyte concentration gradient. As the position of the concentration gradient changes, the resistance of the pipette changes. The resistance signal from the pipette is used as feedback to control a pressure source that applies pressure to the inside of the pipette to counterbalance pressure from the outside. The feedback loop forces the gradient to a neutral balance point, which is user-defined at atmospheric pressure beforehand. The internal pressure required to equally balance the external pressure to the neutral point is readily measured, and it is converted into an analog voltage available at the pressure output BNC and displayed numerically on the LED meter.

### System requirements

The **900A** requires stable sources of both pressure and vacuum, which are essential for the system to rapidly counteract changing pressures encountered at the pipette tip. Pressure and vacuum sources are not provided with the **900A** instrument because some labs are already equipped with suitable sources of pressure and vacuum. For researchers who do not possess pressure or vacuum sources, a cylinder of compressed air or inert gas with a dual stage regulator serves very well as a pressure source. Vacuum must be very stable. It is often best supplied by a quality vacuum pump. WPI offers a very quiet continuous duty vacuum pump well suited for use with the system. In addition, a vacuum regulation kit is recommended to fine tune the vacuum source to the ideal level (**900A-VAC**).

A manometer or meter for independent pressure measurement is necessary to calibrate the pressure and vacuum sources, as well as for validation of the performance of the **900A** system prior to experimental use. A pressure measurement device capable of measuring within a range of +300 mmHg and –150 mmHg is recommended (**PM015D** or **PM015R**). For system performance validation at pressures well below 100 mmHg, the **PM01D** or **PM01R** is recommended, because it provides higher resolution at low pressure.

For transient response performance evaluation, a rapid burst of air or water pressure is required. WPI's **PV830** or **PV820** series PicoPumps provide this capability. Rapidly occurring transient pressure measurements are typically captured on a data acquisition system. For details about testing and measurement of rapidly occurring pressure phenomena, contact a WPI sales representative for additional information.



### References

- Inamoto, R., et al.** (2017). The difference in endolymphatic hydrostatic pressure elevation induced by isoproterenol between the ampulla and the cochlea. *Auris Nasus Larynx*, 44(3), 282–287. <http://doi.org/10.1016/j.anl.2016.07.018>
- Wei, J., et al.** (2017). Role of intratubular pressure during the ischemic phase in acute kidney injury. *American Journal of Physiology - Renal Physiology*, 312(6), F1158–F1165. <http://doi.org/10.1152/ajprenal.00527.2016>
- Inamoto, R., et al.** (2017). The difference in endolymphatic hydrostatic pressure elevation induced by isoproterenol between the ampulla and the cochlea. *Auris Nasus Larynx*, 44(3), 282–287. <http://doi.org/10.1016/j.anl.2016.07.018>
- Wei, J., et al.** (2017). Role of intratubular pressure during the ischemic phase in acute kidney injury. *American Journal of Physiology-Renal Physiology*, 312(6), F1158–F1165. <http://doi.org/10.1152/ajprenal.00527.2016>
- Petrie, R. J., et al.** (2017). Activating the nuclear piston mechanism of 3D migration in tumor cells. *The Journal of Cell Biology*, 216(1), 93–100. <http://doi.org/10.1083/jcb.201605097>
- Warmerdam, T., et al.** (n.d.). Perilymphatic and endolymphatic pressures during endolymphatic hydrops. *European Archives of Oto-Rhino-Laryngology*, 260(1), 9–11. <http://doi.org/10.1007/s00405-002-0518-2>
- Petrie, R. J., & Koo, H.** (2014). Direct measurement of intracellular pressure. *Current Protocols in Cell Biology / Editorial Board, Juan S. Bonifacino ... [et al.]*, 63, 12.9.1–9. <http://doi.org/10.1002/0471143030.cb1209s63>

### SPECIFICATIONS

PRESSURE RANGE	$\pm 0$ –100 mmHg (900APP100)
LINEARITY	< $\pm 0.5\%$ from a straight line
STABILITY	$\pm 0.1$ mmHg up to 1 hour or more
ACCURACY	$\pm 0.5\%$ of full scale
RISE TIME	>10 ms (10–90%), based on residual volume
OUTPUT ("Pressure Signal")	10 mV/mmHg
AMPLIFIER PROBE	Input Resistance >10 <sup>10</sup> $\Omega$ , Voltage Gain 1.0
DIMENSIONS (MAIN FRAME)	43.2 x 13.3 x 25.4 cm (17 x 5.25 x 10 in.)
DIMENSIONS (PRESSURE POD)	9.4 x 2.5 x 5.7 cm (3.7 x 1 x 2.25 in.)
POWER	110 VAC/220 VAC

### ORDERING INFORMATION

**SYS-900A** Micropressure System  
Specify line voltage. See training videos at [www.wpiinc.com/sys-900a](http://www.wpiinc.com/sys-900a).

### OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>900AP</b>	Replacement Probe
<b>CAL900A</b>	Pressure Calibration Chamber
<b>3491</b>	Probe Extension Cable
<b>2933</b>	Rack Mount Kit
<b>5332</b>	Replacement Liquid Trap
<b>MEH6RF</b>	Micropipette Holder (1.0, 1.2, 1.5 or 2.0 mm – Specify OD)
<b>TIPTW900A</b>	Pre-pulled Micropipette for 900A (1 mm thin-wall, 2 $\mu$ Tip) (pkg. of 10)
<b>900APP100</b>	Replacement Pressure Pod
<b>SYS-PM015D</b>	Pressure Manometer (15 psi)
<b>SYS-PM01D</b>	Pressure Manometer (1 psi)
<b>801566</b>	Vacuum Pump, 110V
<b>801963</b>	Vacuum Pump, 220V
<b>900A-VAC</b>	Vacuum Pump Regulation Kit

# Pressure Manometer

For measuring hydrostatic pressures

## Features

- Measure vacuum and pressure in gases
- Ranges available:  $\pm 1$  PSI,  $\pm 15$  PSI,  $\pm 100$  PSI
- Measure in PSI or kPa on the 100 PSI unit or PSI and mmHg on the 15PSI unit
- Battery powered
- Includes tubing and mini-phone to BNC cable

## Benefits

- Easy and accurate measurements

## Applications

- Measure pressure of non-corrosive gases

Hand-held and battery operated, PM Series pressure manometers monitor vacuum and pressure in non-corrosive gases. An integral transducer and digital display allow easy and accurate pressure readings. Three versions measure pressures in the range of  $\pm 1$  PSI,  $\pm 15$  PSI or  $\pm 100$  PSI. A range switch allows measurement in units of PSI or kPa for the 100 PSI version, and PSI or mmHg for the 15 psi version. Pressure can be read on the built-in LCD display or relayed to a chart recorder, oscilloscope, or computer.

PM Series pressure manometers come with 4 feet of  $\frac{1}{8}$ -in. ID soft vinyl tubing. A mini-phone-to-BNC cable for the recorder output is also available (**CBL102**). Standard versions are equipped with a nine-volt alkaline battery.



PM100

CE

## PRESSURE MANOMETER SPECIFICATIONS

	PM01	PM015	PM100
PRESSURE RANGES	$\pm 1$ psi ( $\pm 52$ mm Hg)	$\pm 15$ psi ( $\pm 775$ mm Hg)	$\pm 100$ psi ( $\pm 690$ kPa)
MAX. PRESSURE	20 psi (1035 mm Hg)	30 psi (1550 mm Hg)	150 psi (1035 kPa)
RESOLUTION	0.001 psi (0.1 mm Hg)	0.01 psi (1 mm Hg)	0.1 psi (1 kPa)
OUTPUT	1 V/psi	100 mV/psi	10 mV/psi
OUTPUT RANGE	$\pm 1.0$ V	$\pm 1.5$ V	$\pm 1.0$ V
LINEARITY		0.5% full-scale	
TEMPERATURE EFFECT		1.0% full-scale (0-70°C)	
ZERO		Screwdriver-adjust	
RESPONSE TIME		30 ms	
POWER		Nine-volt battery	
BATTERY LIFE		Alkaline, 200 hours; rechargeable, 25 hours	
RECORDER OUTPUT		Mini-phone jack, 0.141 inch (3.5 mm)	
OUTPUT IMPEDANCE		1 k $\Omega$	
PNEUMATIC CONNECTORS		Barbed, for $\frac{1}{8}$ -inch or $\frac{3}{16}$ -inch ID soft tubing	
DIMENSIONS		3 x 6 x 1 inches (8 x 15 x 4 cm)	
SHIPPING WEIGHT		3 lb (1.4 kg)	

## ORDERING INFORMATION

<b>SYS-PM01D</b>	Pressure Manometer (1 psi)
<b>SYS-PM01R</b>	Pressure Manometer (1 psi), Rechargeable*
<b>SYS-PM015D</b>	Pressure Manometer (15 psi)
<b>SYS-PM015R</b>	Pressure Manometer (15 psi), Rechargeable*
<b>SYS-PM100D</b>	Pressure Manometer (100 psi)
<b>SYS-PM100R</b>	Pressure Manometer (100 psi), Rechargeable*
<b>CBL102</b>	Mini-Phone-to-BNC Cable

Specify line voltage

\*Rechargeable versions come with nickel/cadmium battery and charger

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
 Germany: +49 (0)6031 1602171 • wptide@wpi-europe.com • www.wpi-europe.com  
 Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • www.wpiinc.com  
 China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net

# Animal Physiology



## Specialized Equipment for Small Animal Research

WPI offers a broad range of instruments and laboratory supplies that are designed for small animal physiology. Let us help you get the laboratory equipment that's right for your research application.

# Motorized Stereotaxic Frame

*Automatically calculates coordinates and accurately places the probe*

## Features

- Accurate microstepping motor drive for high resolution placement
- Digital controller
- No computer required
- Touch screen controller for easy, intuitive control with better than 10  $\mu\text{m}$  precision
- Graphic controller display for instant operational feedback
- Hand controller for complete manual control

## Benefits

- Reliable positioning reduces errors
- Increased precision and repeatability of motion
- Feature rich for your convenience of use
- Multiple modes of operation for precise probe positioning

## Applications

- Stereotaxic surgery for small animal research

When precision and repeatability of motion are critical, the **MTM-3** Motorized Stereotaxic Frame and digital display outperforms manual models, and it greatly reduces human error when performing routine stereotaxic surgery. The motorized axis of the **MTM-3** provides precise, controlled, 3-dimensional placement of any probe or accessory within the working space of a stereotaxic frame.

No computer is required. The **MTM-3** 3-axis manipulator supplied with WPI stereotaxic frames is compatible with standard stereotaxic frames and may be adapted to existing frames of other manufacturers.

Single and dual manipulator arm motorized systems are available. This allows you to mount a stereotaxic drill and a probe simultaneously. Both 3-Axis Wheel Control and Touchscreen Controller are included.

### Reliable positioning for reduced errors


The controller automatically calculates the brain map coordinates for precise probe placement during a stereotaxic surgery. This eliminates the errors resulting from reading Vernier scales.

### Increased precision and repeatability of motion

The MTM-3 offers better than 10  $\mu\text{m}$  accuracy with the microstepping motor drive for high resolution placement of your probe, drill or electrode. Coordinate distances are automatically calculated.

### Convenient to use

Brain atlas coordinates may be entered into the controller, and no computer required. The simple to use touch screen controller offers ease of control. The graphic controller display provides instant operational feedback during a stereotaxic surgery. For example, the "final approach" speed may be set between 2 mm/sec and 0.02 mm/sec.



The **MTM-3** is available in both left- and right-hand versions, or get the **MTM-6**, which includes both manipulator arms and dual controllers.

## Multiple modes for positioning

The arms of the motorized stereotaxic frame may be positioned manually using the digital interface or hand controller. Use the Coordinates mode to control the arm positions based on specific mapped coordinates. Set your retracted position, and store up to three reference (origin) positions. Then, display the probe position with respect to any of your references or the absolute coordinates. If you prefer to write short scripts, computer control is also available which uses a simple, text based command set. For example, if your brain research requires that you routinely drill in the same location of the skull and implant a probe to the same prescribed depth, a script could dramatically improve your process.

## SPECIFICATIONS

TRAVEL	80 mm on each axis
PRECISION	10 $\mu\text{m}$

## ORDERING INFORMATION

<b>MTM-3</b>	3-Axis Stereotaxic Frame (Left or Right)
<b>MTM-6</b>	6-Axis Stereotaxic Frame (Left and Right)
<b>MTM</b>	3-Axis Stereotaxic Arm/Controller only (Left), Upgrade
<b>MTM-R</b>	3-Axis Stereotaxic Arm/Controller only (Right), Upgrade

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
 Germany: +49 (0)6031 1602171 • wptide@wpi-europe.com • www.wpi-europe.com  
 Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • www.wpiinc.com  
 China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net



# Digital Stereotaxic Frame with LCD Display

Battery powered, low noise for electrically quiet operation



## Features

- Large LCD display from sealed sensors on each axis
- Adaptors available for use with rats, mice, birds, cats, guinea pigs and other species
- 80 mm of travel in three planes
- 180° rotation and lock at any vertical angle
- 360° rotation and lock at any level angle
- Syringe pump and drill can be attached directly
- Stable and accurate movement
- Targets specific coordinates by zeroing
- 10 µm precision in all directions
- Extended base plate (400 mm x 255 mm) suitable for various animals
- Vertical lock and fixing knob are separated to ensure accurate position at any angle

## Benefits

- Quiet, battery-powered operation
- Battery-powered sensors, without electronic noise, are suitable for electrophysiology experiments

## Applications

- Stereotaxic surgery for small animal research
- Electrophysiology applications

This new Digital Stereotaxic Frame features sealed electronic sensors and an easy-to-read LCD display with 10 µm resolution in all three axes. A zeroing function aids in targeting specific coordinates.

- The precisely designed rotary knob and U frame allow sufficient space for the anterior-posterior operation
- Laser engraved scales and a darkened rod for easy to read scales
- Curved design of nose clamp fixes the head of the animal more securely

### ORDERING INFORMATION

<b>502300</b>	Digital Stereotaxic Frame with 18° Ear Bars
<b>502350</b>	Digital Stereotaxic Frame with 45° Ear Bars
<b>502303</b>	Dual Manipulator Digital Stereotaxic Frame, 18° Ear Bars
<b>502353</b>	Dual Manipulator Digital Stereotaxic Frame, 45° Ear Bars
<b>TAXIC-300</b>	Digital Stereotaxic Frame with 18° Ear Bars & UMP3T-1
<b>TAXIC-303</b>	Dual Manipulator Digital Stereotaxic Frame with 18° Ear Bars with UMP3T-1 syringe pump
<b>TAXIC-350</b>	Digital Stereotaxic Frame with 45° Ear Bars & UMP3T-1
<b>TAXIC-353</b>	Dual Manipulator Digital Stereotaxic Frame with 45° Ear Bars with UMP3T-1 syringe pump

### UPGRADE MANUAL STEREOTAXIC TO DIGITAL

<b>502360</b>	Manipulator, 3 Axes, Right Hand, Digital Display
<b>502361</b>	Manipulator, 3 Axes, Left Hand, Digital Display

# Precision Stereotaxic Instruments

## For small animal research

### Features

- 5 mm linear movement per revolution on each axis
- Absolute lock at 90° (vertical)
- Entire frame is electrically continuous, ideal for electrophysiology
- Accessories available for wide variety of small animals

### Benefits

- Versatile positioning
- Easy to read scales
- Convenient for electrophysiology because the frame may be grounded

### Applications


- Stereotaxic surgery for small animal research
- Electrophysiology research

WPI's Precision Stereotaxic Instrument is built around the time-proven U-frame design concept, providing stability and adaptability to most species. Precision alignment ensures accurate placement of electrodes, micropipettes and other devices. It is ideal for researchers in need of a versatile, reliable instrument for stereotaxic procedures with small animals.

The manipulator arm controls mediolateral and vertical positioning via lead screws with 80 mm of travel. This allows the fastest positioning possible, consistent with lining up the scales easily at a given coordinate. The anteroposterior movement is controlled via a dovetail slide movement, with 80 mm of travel possible in each direction. A universal joint allows the investigator to change the angle of the probe up to 90° in either the anteroposterior or mediolateral planes. The locking mechanism will hold any angle position without drift or creep. It also provides an absolute lock at 90° vertical.

All scales are oriented to be read easily from the open end of the "U." This is the position from which most scientists prefer to work. The numerals on the scales are clear and easy to read. Precise alignment with facing Vernier scales gives accurate resolution to 0.1 mm.

The entire stereotaxic frame, including the dovetails, manipulator arms and base are electrically continuous. Grounding of the entire frame including the base plate can be accomplished by connecting the provided grounding stud to earth. This is ideal for electrophysiological studies where the animal and surrounding structures need to be grounded to reduce electrical noise.




**TAXIC-600**  
Stereotaxic Frame with  
UltraMicroPump Injection  
System



502600

502603



502612

*Stereotaxic  
Frame for Rat,  
Portable, Single  
Manipulator, 18°  
Ear Bars*

502614

*Stereotaxic Frame  
for Rat, Rail-  
mounted, Single  
Manipulator, 18°  
Ear Bars*

### ORDERING INFORMATION


<b>502600</b>	Precision Stereotaxic Frame with 18° Ear Bars
<b>502650</b>	Precision Stereotaxic Frame with 45° Ear Bars
<b>502603</b>	Dual Manipulator Stereotaxic Frame with 18° Ear Bars
<b>502653</b>	Dual Manipulator Stereotaxic Frame with 45° Ear Bars
<b>502612</b>	Stereotaxic Frame, Single Manipulator, for Rat, 18° Ear Bars
<b>502614</b>	Stereotaxic Frame, Single Manipulator, Rail Mounted, for Rat, 18° Ear Bars
<b>TAXIC-600</b>	Stereotaxic Frame with 18° Ear Bars plus UMP3T-1
<b>TAXIC-650</b>	Stereotaxic Frame with 45° Ear Bars plus UMP3T-1
<b>TAXIC-603</b>	Dual Manipulator Stereotaxic Frame with 18° Ear Bars plus UMP3T-1 System
<b>TAXIC-653</b>	Dual Manipulator Stereotaxic Frame with 45° Ear Bars plus UMP3T-1 System

### WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
 Germany: +49 (0)6031 1602171 • wpide@wpi-europe.com • www.wpi-europe.com  
 Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • www.wpiinc.com  
 China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.com

# Mice and Neonatal Rats

*Specialized equipment for the smallest subjects*



## Features

- Base plate dimensions: 25.5 cm × 25.25 cm
- Unique and light ear bars composed of polyoxymethylene resin
- Ear bar and palate clip height regulation range is 20 mm (0 mm~+20 mm)

## Benefits

- Laser engraved scales shown on the ear bars make the unit easy to operate (1 mm accuracy)
- Two types of ear bar tips provide a secure hold without compressing the sides of the skull inward
- Dual arm is optional
- Digital LED display is optional

## Applications

- Small animal surgery
- Electrophysiology applications

This stereotaxic instrument uses an adapter block to replace the traditional U-frame assembly, and it holds the mouse or neonatal rat firmly and accurately. This instrument provides easy access to the animal and offers an open working area for extra surgical instruments. This unit is an upgrade of the digital stereotaxic frame shown on page 119.

**502610**  
*Stereotaxic Frame for Mouse, Portable, Single Manipulator*

**504926**  
*Mouse/Neonatal frame with the optional digital controller*

## ORDERING INFORMATION

<b>502610</b>	Stereotaxic Frame, Single Manipulator, for Mouse
<b>504926</b>	Digital Mouse Frame, Single
<b>504927</b>	Digital Mouse Frame, Dual

# Parallel Rail Stereotaxic Instruments

*For large animal research*

## Features

- 5 mm linear movement per revolution on each axis
- Entire frame is electrically continuous, ideal for electrophysiology
- Includes the U-frame base plate, manipulator(s), cat/monkey adaptor, 18° ear bars and swivel mount

## Benefits

- Versatile positioning
- Easy to read scales
- Convenient for electrophysiology because the frame may be grounded
- Ability to expand to accommodate a variety of large animals

## Applications

- Stereotaxic surgery for large animals
- Electrophysiology applications



502232

WPI's Parallel Rail Stereotaxic Frame systems are heavy-duty research instruments for large laboratory animals such as cats, monkeys and dogs. The solid, large frame and superior rigidity ensure the precise alignment of animals for stereotaxic surgery, injection and recording. The system can accommodate up to four manipulators with 100 µm resolution on each axis. Each manipulator can smoothly move to and lock at any location on both parallel rails in a range of 20 cm.

Parallel Rail Stereotaxic Frame system for large animals includes the Parallel Rail Frame, base plate, manipulator(s), cat/monkey or dog adaptor and swivel mount.



502616

*MRI Stereotaxic Instrument for Dogs and Monkeys*

## ORDERING INFORMATION

<b>502227</b>	Stereotaxic Frame, 1 Manipulator for Cat/Monkey
<b>502228</b>	Stereotaxic Frame, 2 Manipulators for Cat/Monkey
<b>502229</b>	Stereotaxic Frame, 3 Manipulators for Cat/Monkey
<b>502230</b>	Stereotaxic Frame, 4 Manipulators for Cat/Monkey
<b>502231</b>	Stereotaxic Frame, 1 Manipulator for Dog
<b>502232</b>	Stereotaxic Frame, 2 Manipulators for Dog
<b>502233</b>	Stereotaxic Frame, 3 Manipulators for Dog
<b>502234</b>	Stereotaxic Frame, 4 Manipulators for Dog

## ORDERING INFORMATION


<b>502616</b>	MRI Stereotaxic Frame, for Dogs and Monkeys
---------------	---

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • [www.wpi-europe.com](http://www.wpi-europe.com) Germany: +49 (0)6031 1602171 • wptide@wpi-europe.com • [www.wpi-europe.com](http://www.wpi-europe.com)  
 Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • [www.wpiinc.com](http://www.wpiinc.com) China: +86 21 6888 5517 • chinasales@china.wpiinc.com • [www.wpiinc.net](http://www.wpiinc.net)

# Stereotaxic Accessories


Instruments for small animal research



502063

The WPI Mouse and Neonatal Rat Adaptor (502063) employs light, Delrin® adjustable ear bars with tapered points on one end and thumbscrew on the other to facilitate surgery on mice and rat pups. Adjustable ear bars may be independently adjusted in height to level the skull. Laser engraved scales show the vertical positions of the ear bars. A tooth bar and nose clamp secures the nose. A well in the thick aluminum body may be filled with dry ice and alcohol for hypothermic anesthesia of neonatal animals. The adaptor clamps securely on the right side of the "U" frame of the stereotaxic instrument. When gas anesthesia is needed, a gas mask assembly (502213) may be mounted on the adaptor.

502241



502265




502259—Rat Spinal Adaptor



502266

## ORDERING INFORMATION (ADAPTORS)

502063	Mouse and Neonatal Rat Adaptor
502213	Platform, Gas Anesthesia, with Mouse Mask (use with 502063)
502265	Mouse Adaptor with Nose Clamp
502266	Rat Adaptor
502226	Cat/Monkey Adaptor for 502600 series
502259	Spinal Adaptor for Rat
502060	Guinea Pig Adaptor for 502600 series
502241	Dog/Monkey Adaptor for Parallel Rail Instruments, with Ear Bars, 18°
502259	Rat Spinal Adaptor, for Microscope Stage



**ORDERING INFORMATION (PROBE HOLDERS)**

<b>502210</b>	Probe Holder with Corner Clamp, 0.3-1.5 mm opening
<b>502067</b>	Probe Holder with Side Clamp, 0.3-3.5 mm opening
<b>502070</b>	Cannula Holder, Metal Tip, opens to 3.4 mm (See cannulas, page 134)
<b>502256</b>	Cannula Holder, plastic tip, opens to 3.4 mm (See cannulas, page 134)
<b>502068</b>	Large Probe Holder, 6-12 mm Opening
<b>502237</b>	Extra Large Holder for OmniDrill35 Microdrill
<b>502236</b>	Microdialysis Probe Holder, 1.5 mm Hole
<b>502244</b>	Micrometer Adjust Electrode Holder, 10 µm Resolution, 25 mm Travel, 0.3-1.5 mm Opening
<b>502245</b>	Manual Microsyringe Injection Holder, 10 µm Resolution, 25 mm Travel

502235      502056      502225      502242



**ORDERING INFORMATION (EAR BARS)**

<b>502055</b>	Ear Bars, Rat, 18° (pair)
<b>502056</b>	Ear Bars, Rat, 45° (pair)
<b>502224</b>	Ear Bars, Cat, 18° (pair)
<b>502225</b>	Ear Bars, Cat, 45° (pair)
<b>502235</b>	Ear Bars, Mouse, 60°, Non-rupture (pair)
<b>502242</b>	Ear Bars, Rat, Hollow, 1.5 mm Hole for Auditory Stimulation
<b>502249</b>	Non-Rupture, Rubber Ear Bars for Mice
<b>502250</b>	Non-Rupture, Cuff Ear Bars for Mice



**502257**  
MRI Rat Adaptor

**ORDERING INFORMATION (OTHER ACCESSORIES)**

<b>502257</b>	MRI Rat Adaptor, 36 x 12.5 x 9.6 cm
<b>502201</b>	V-Clamp, with 10/32 Screw for Mounting UMP3
<b>504608</b>	V-Clamp, with 10/32 Thumb Screw
<b>502213</b>	Platform, Gas Anesthesia with Mouse Mask (use with 502063)
<b>502243</b>	Adjustable Stage Platform for 502600 series (5 cm high)
<b>503598</b>	Micro-Drill, 35K rpm, with a set of bits, 110/220VAC
<b>503599</b>	Micro-Drill, 35K rpm, with a set of bits, 240VAC, UK plug
<b>61840</b>	Heating Plate for 502063, 4x15 cm (use with ATC2000)

# Stereotaxic Anesthesia Facemask

*Ensures a snug fit for mice and rats with zero leakage and zero dead space*

## Features

- Facemask custom fitted to a toothbar that fits a Cunningham Stereotaxic Frame
- Ports swivel 360°
- Optically clear
- Includes OC-tubing and OC-adaptor (requires EZ-1130)

## Benefits

- Anatomically correct design for a tight fit for mice and rats
- Zero leakage
- Zero dead space
- Provides access to eyes and head of the rodent
- Mask provides stability of the animal for examinations and surgery
- Reduces animal stress by limiting handling
- Reduces technician's exposure to gas

## Applications

- Small animal stereotaxic surgery

Our custom toothbar for the Cunningham adaptor is designed to accommodate the use of a soft gas mask. The included mask securely fits the toothbar to eliminate gas leakage.

Choose either small or medium size. The **502063-SGM** and **502063-MGM** are customized versions of our standard rodent gas masks.

The OC gas masks allow you access to the eyes and head of the rodent while sealing the nose and mouth. The gas masks provide stability of the animal for examinations and surgery, and they can be easily incorporated in existing examination procedures. The waste-line connects to gravity or vacuum filters, and swiveling ports allow for greater flexibility in positions. They are easily adapted for stereotaxic applications.

Thanks to the anatomical face mask design, you have head-to-tail access. The face mask reduces the technician's exposure to gas by sequestering and filtering waste, minimizing dead space and leakage.

You can reduce rodent stress by limiting handling required when using injectable anesthesia. By choosing these gas masks with isoflurane and thereby reducing the use of Schedule II drugs, you can simplify your lab management. Setting up and changing masks is simple, and the swivel adaptors allow you to make adjustments quickly and even modify the mask for stereotaxic uses. Gas anesthesia reduces recovery time, making for a faster and more efficient lab.

*The 502063-MGM gas mask fits onto the 502063 Cunningham Stereotaxic Frame.*



### RODENT FACEMASKS

\*Requires the use of OC-TUBING, OC-ADAPTER and EZ-1130.

<b>OC-SFM</b>	Small Mask	Small Mouse 0.5 – 2 months 8 mm Aperture x 25 mm overall length
<b>OC-MFM</b>	Medium Mask	Mouse 2+ months and Small Rat 10 mm Aperture x 30 mm overall length
<b>OC-LFM</b>	Large Mask	Adult Rat up to 250 g 13 mm Aperture x 35 mm overall length
<b>OC-XLFM</b>	Extra Large Mask	Adult Rat from 250 – 500 g 18 mm Aperture x 42 mm overall length
<b>502063-SGM</b>	Small Mask for Cunningham Stereotaxic Frame *Includes OC-TUBING and OC-ADAPTER	Small Mouse 0.5 – 2 months 8 mm Aperture x 25 mm overall length
<b>502063-MGM</b>	Medium Mask for Cunningham Stereotaxic Frame *Includes OC-TUBING and OC-ADAPTER	Mouse 2+ months and Small Rat 10 mm Aperture x 30 mm overall length



### ORDERING INFORMATION

**502063-SGM** Small Mask for Cunningham Stereotaxic Frame  
*Includes OC-TUBING and OC-ADAPTER*

**502063-MGM** Medium Mask for Cunningham Stereotaxic Frame  
*Includes OC-TUBING and OC-ADAPTER*

**PLEASE NOTE:** This part is made to fit the WPI neonatal frame only. If you have a neonatal frame of another make, please contact WPI.

# Stereotaxic Frame Anesthesia Masks

*Active and passive anesthesia masks for use with stereotaxic equipment*

Two kinds of stereotaxic frame nosecone masks are available:

- Passive masks have collinear inlets and outlets.
- Active masks have separate (non-collinear) inlets and outlets on both sides

The passive mask is used with a charcoal canister only, and the active mask can be used with an evacuation apparatus.


## Active Anesthesia



502046 (above) — 20 mm ID  
502047 (not shown) — 22 mm ID



502042 (above) — 10 mm ID  
502043 (not shown) — 11.5 mm ID



502262

### ORDERING INFORMATION (ACTIVE SCAVENGING)

<b>502042</b>	Stereotaxic Anesthesia/Evacuation Mask, Mice or Neonatal Rats, < 30 g
<b>502043</b>	Stereotaxic Anesthesia/Evacuation Mask, Mice or Neonatal Rats, 30-70 g
<b>502046</b>	Stereotaxic Anesthesia/Evacuation Mask, Rats, < 300 g
<b>502047</b>	Stereotaxic Anesthesia/Evacuation Mask, Rats, > 300 g
<b>502262</b>	Mouse Gas Anesthesia Head Holder with Rubber Ear Bars

## Passive Anesthesia



502213



502054

502053

### ORDERING INFORMATION (PASSIVE SCAVENGING)

<b>502213</b>	Platform, Gas Anesthesia, with Mouse Mask (use with 502063)
<b>502053</b>	Stereotaxic Mask, Gas Anesthesia, Mouse or Neonatal Rat (< 70g)
<b>502054</b>	Stereotaxic Mask, Gas Anesthesia, Rat (< 300 g)



# Animal Temperature Controller

Maintain optimal animal temperature during research procedures

## Features

- Animal Temperature Controller with Adaptive mode—auto adjust PID regardless of animal size
- Plate's internal temperature sensor prevents localized hot spots under animal with maximum temperature stability
- Compatible with RTD (resistive temperature device) and thermocouple probes
- Extremely quiet DC heater to facilitate electrophysiological recordings
- Three temperature sensor inputs
- Auto shutdown if the plate reaches 45°C

## Benefits

- Low noise temperature control system
- Adapts to different size animals automatically
- Accepts RTD and thermocouple probes

## Applications

- Small animal surgery
- Stereotaxic surgery

The **ATC2000** Animal Temperature Controller is a low-noise heating system for maintaining animal body temperature during experimental procedures. The DC heater is extremely quiet in terms of electromagnetic radiation. This is essential in electrophysiological recordings which are very sensitive to electromagnetic interference. The controller uses proportional, integral, derivative (PID) technology to provide precise and stable control of a subject's temperature. Compared with switched on/off type controllers, PID controllers provide a much more precise and stable control of temperature. The PID approach is also more immune to the variation of the experimental conditions such as change in animal size and unexpected disturbances. Our unique adaptive mode technology automatically senses and adjusts the PID values based on the size of the animal.

## Multiple temperature sensing inputs

The **ATC2000** has three temperature sensing inputs.

- The resistive temperature device (RTD) probe input can be used to monitor an RTD rectal probe to control the animal temperature or to monitor ambient temperature, induction chamber temperature or any other temperature.
- When using a thermocouple probe, the thermocouple (TC) probe input can be used just like the RTD input. (A T type thermocouple must be used.)
- The heater plate also has an internal temperature sensor. The **ATC2000** monitors this sensor to prevent the localized hot spots under animal.

## Operational modes

The controller has three operational modes:

- Normal mode uses the configured sensor (RTD or TC) or the plate sensor to drive the control loop.
- Adaptive mode uses the temperature of the heated plate and the temperature of the subject to control. This approach is less prone to overshoot, but somewhat slower than the normal mode, depending on the sampling rate used.
- Shutdown is a fail safe mode used if the plate temperature ever exceeds 45°C.

## Required Accessories

The **ATC2000** is tuned at the factory. However, the PID parameters may also be set manually. The temperature resolution of the controller is 0.1°C. The rectal temperature probe has a 5-ft. ultra-flexible, shielded cable and an RTD sensor.



ATC2000

The metal heating plates (available separately) have built-in temperature sensors. The smaller heating plate (WPI# **61840**) is compatible with stereotaxic systems; the rigid, flat surface of the warming pad fits under the U-frame. Our homeothermic warming pads are washable with water and detergent.

A heating plate and a probe are required for use with this unit.

## References

- Nguyen Chi, V., Müller, C., Wolfenstetter, T., Yanovsky, Y., Draguhn, A., Tort, A. B. L., & Brankač, J. (2016). Hippocampal Respiration-Driven Rhythm Distinct from Theta Oscillations in Awake Mice. *Journal of Neuroscience*, 36(1).
- Okun, M., Lak, A., Carandini, M., Harris, K. D., Buzsaki, G., Stevenson, I., ... Kaufman, M. (2016). Long Term Recordings with Immobile Silicon Probes in the Mouse Cortex. *PLOS ONE*, 11(3), e0151180. <http://doi.org/10.1371/journal.pone.0151180>
- Gaylo, A., Overstreet, M. G., & Fowell, D. J. (2016). Imaging CD4 T Cell Interstitial Migration in the Inflamed Dermis. *Journal of Visualized Experiments*, (109), e53585–e53585. <http://doi.org/10.3791/53585>



**61840**—Small heating plate (15x4 cm) for use with stereotaxic frames.

## ORDERING INFORMATION

**ATC2000** Animal Temperature Controller

### RECOMMENDED ACCESSORIES (select one sensor and one plate)

**61800** Heating Plate with Built-in RTD Sensor, 15 x 25 cm  
**61830** Heating Plate with Built-in RTD Sensor, 15 x 10 cm  
**61840** Heating Plate with Built-in RTD Sensor\*, 15 x 4 cm

\*for stereotaxic frame

### OPTIONAL ACCESSORIES/REPLACEMENTS

**61824** RTD Rectal Temp Probe, 1.25 mm Shaft 2.5 mm Ball Tip  
**RET-2** TC Rat Rectal Temp Probe, 1 mm Shaft, 3.2 mm Ball Tip  
**RET-3** TC Mouse Rectal Temp Probe, 1 mm shaft, 1.6 mm Ball Tip  
**503573** Silicone Pad for ATC2000 (10 x 15 cm)  
**MT-29/1HT** 29 ga 1 cm Needle Microprobe, 5-ft Cable  
**MT-23/3HT** 23 ga 3 cm Needle Microprobe, 5-ft Cable  
**MT-D** Needle Microprobe, 5-ft Cable  
**IT-18** Flexible Implantable Probe, 0.025-inch Diam., 3-ft Cable  
**IT-23** Flexible Implantable Probe, 0.009-inch Diam., 3-ft Cable  
**IT-1E** Flexible Implantable Probe, 0.025-inch Diam., 3-ft Cable

# Small Animal Anesthesia

Complete systems and all the accessories

EZ-B800 Basic System



Components may not be as pictured.

## Features

- Safe for surgical personnel, 90% below OSHA isoflurane limit
- Time efficient and cost-effective
- Virtually stress free for the animals
- Speedy recovery time

## Benefits

- Compact and portable
- Easy to set up and use, simplifying the training of new staff and reducing the risk of human error
- Turnkey (all-in-one) system or plug-and-play system

## Applications

- Small animal surgery

**EZ-Anesthesia** is the system of choice for anesthetizing small animals, and it comes with a variety of choices. Animals to be anesthetized are placed in the acrylic induction chamber, and the system delivers a precisely blended mixture of oxygen and isoflurane. An activated charcoal air filter canister at the top of the chamber releases safe, filtered air back into the room. A water-heated cage warmer or warming plate (**ATC2000**) is used to retain the animal body temperature while in the induction chamber. After the initial anesthetizing, the animal may be moved to the heated surgical water bed and positioned properly in the snugly fitted nose cone. A highly sensitive valve regulated by the animal's breathing works with the nose cone to ensure non-rebreathing efficiency. It allows safe anesthesia for up to several hours.

The breathing device also includes an air filter that releases safe, filtered air back into the room.

Oxygen and liquid isoflurane are not supplied. These are required to operate the system. Each system comes with all necessary components and connections for immediate use, including:

- Oxygen regulator
- Vaporizer unit
- Induction chamber (standard **EZ-178**)
- Breathing device (standard **EZ-103**)
- Case of charcoal filters
- Connecting tubing

Other available options include:

- Heated water beds
- Additional breathing devices
- Chambers in a variety of sizes
- Key fill isoflurane vaporizer is available separately or may be substituted in a complete system (call for details)

## EZ-FF9000 Fixed Flow System

The Fixed Flow System is our most advanced system and provides preset, fixed flow rates with no need for adjustment. The system offers five gas outlets, each with an individual ON/OFF switch. The system ensures consistent, precise gas flow when connected to any pressurized gas source higher than 7 psi.

- Induction chamber output is fixed at 1 LPM
- Four breathing circuit outputs are fixed at 0.5 LPM

This unit features an oxygen flush button that purges the induction chamber with pure oxygen, thereby protecting personnel from exposure to anesthesia gas when opening the chamber. A handle on the vaporizer makes the unit easy to safely transport.

## EZ-7000 Classic System

The Classic System has a user-friendly flow meter. The system offers five gas outlets that can supply a single induction chamber and four breather circuits simultaneously.

The system also features an oxygen flush button that purges the induction chamber with pure oxygen, thereby protecting personnel from exposure to anesthesia gas when opening the chamber. A handle on the vaporizer makes the unit easy to safely transport.

## EZ-B800 Basic System

The Basic System is designed for use with a single animal. It utilizes one output directly from the vaporizer into a Y-splitter which creates a dual feed that can be directed to the induction chamber or the breather circuit.

This unit incorporates an oxygen flush system that purges the induction chamber with pure oxygen, thereby eliminating personnel exposure to anesthesia gas when opening the chamber.

### ORDERING INFORMATION

<b>EZ-FF9000</b>	Fixed Flow System
<b>EZ-7000</b>	Classic System
<b>EZ-B800</b>	Basic System

### WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • www.wpiinc.com

Germany: +49 (0)6031 1602171 • wpide@wpi-europe.com • www.wpi-europe.com  
China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net

# Anesthesia Accessories

## Induction Chambers

These chambers incorporate a positive seal O-ring gasket for containment during use.

### ORDERING INFORMATION

<b>EZ-177</b>	Sure-Seal Mouse Chamber (5"L x 4.75"W x 4.38"H)
<b>EZ-178</b>	Sure Seal Mouse/Rat Chamber (9.75"L x 4.75"W x 4.38"H)
<b>EZ-1785</b>	Large Mouse/Rat Chamber (7"W x 11"D x 6"H inside)
<b>EZ-179</b>	Rabbit/Guinea Pig Chamber (18.75" x 12.75" x 12.75")

## Ventilator Connection Kit

This kit is used to connect the **SAR-830** ventilators with EZ-Anesthesia systems. It includes all required components, pre-assembled for simple connection between the ventilator and the anesthesia system.

### ORDERING INFORMATION

<b>EZ-830</b>	Ventilator Connection Kit
---------------	---------------------------

## Connection Kit

This kit is used for connecting additional components to the EZ-Anesthesia Systems. It includes 6 ft of clear PVC tubing with a quick-disconnect fitting.

### ORDERING INFORMATION

<b>EZ-1130</b>	Connection Kit
----------------	----------------

## EZ Oxygen Regulators

Regulators are preset to 50 psi.

- **EZ-320** utilizes a CGA-540 connection for large "H" tanks.
- **EZ-330** utilizes a CGA-870 connection for small "E" tanks.

### ORDERING INFORMATION

<b>EZ-320</b>	Oxygen Regulator for large tanks
<b>EZ-330</b>	Oxygen Regulator for small tanks

## Custom Hose Assemblies Available

Custom built hose assemblies for wall or ceiling outputs specific to facility needs: Chemetron, Ohmeda, Schrader or DISS.

## Other Accessories



### ORDERING INFORMATION

<b>EZ-104A</b>	Versaflex Non-Rebreathing Unit
<b>EZ-103A</b>	Microflex Non-Rebreathing Unit
<b>EZ-107A</b>	Rat Stereotaxic Non-Rebreathing Unit
<b>EZ-109</b>	Multi-Animal Non-Rebreathing Unit
<b>EZ-211</b>	Mouse/Rat Thin-Line Heated Waterbed

## Mobile Workstations

Two mobile workstations, constructed of heavy-duty stainless steel with locking casters, integrate all your EZ-Anesthesia components into one portable unit. Open side shelves accommodate 20 lb. cylinders, and convenient 2" port holes allow for easy rigging of gas and electrical lines. Below the work surface of each mobile workstation is an open shelf and a locking cabinet. **EZ-E25000** provides a 42"x24" work surface and holds up to four cylinders. **EZ-E27000** has a 22"x21" work surface and holds up to two cylinders. These systems are easy to set up and provide maximum flexibility and mobility.

### ORDERING INFORMATION

<b>EZ-E25000</b>	Mobile Workstation, 42" x 24" Top
<b>EZ-E27000</b>	Mobile Workstation, 22" x 21" Top

## Euthanex Lids

The Euthanex Lid system has become an industry standard. You no longer need to transfer animals. The lids are available in five sizes to accommodate virtually all plastic cage sizes designed for small lab animals. These heavy duty stainless steel lids have a stem fitting for connection to the quick-disconnect fitting on the hose from the gas source. A foam lid gasket ensures a good seal on the cage. Multiple lids may be used to treat several cages at once.



### ORDERING INFORMATION


<b>EZ-20027</b>	Small Lid (13" x 9") Fits old-style mouse cages
<b>EZ-20028</b>	Small Lid (16" x 10") Fits new-style mouse cages
<b>EZ-20030</b>	Square Lid (13" x 13") Fits Thorn cages
<b>EZ-20032</b>	Medium Lid (20.5" x 11") Fits rat cages
<b>EZ-20034</b>	Large Lid (23" x 16.5") Fits guinea pig cages
<b>EZ-20029</b>	Lid Storage Bracket (wall-mounted, holds up to four lids)
<b>EZ-20027G</b>	Small Lid Gasket (13" x 9")
<b>EZ-20028G</b>	Small Lid Gasket (16" x 9")
<b>EZ-20030G</b>	Square Lid Gasket (13" x 13")
<b>EZ-20032G</b>	Medium Lid Gasket (20.5" x 11")
<b>EZ-20034G</b>	Large Lid Gasket (23" x 16.5")

## Vaporizer Pole Mount

The **EZ-E28000** is a mounting option for the EZ-Anesthesia system, combining system portability with a small footprint. Constructed of aluminum and stainless steel, it features a five-leg base for maximum stability and an oxygen "E" tank mounting bracket. Anesthesia system must be purchased separately.

### ORDERING INFORMATION

<b>EZ-E28000</b>	Vaporizer Pole Mount
------------------	----------------------



# Blood Pressure Monitor and Transducer

*Audible monitor with variable pitch provides feedback*

## Features

- Variable amplification
- Displays systolic, diastolic or average blood pressure
- Audio monitor
- BNC output for data logging

## Benefits

- **BLPR2** can be used for the direct arterial and venous pressure measurement in animal blood vessels

## Applications

- Invasive monitoring of animal arterial or venous blood pressure

**BP1** accepts WPI's **BLPR2** blood pressure transducer (below) as well as other blood pressure transducers. An audio monitor provides a signal with variable pitch and amplitude, allowing you to hear changes in blood pressure. Digital LCD display provides average or peak signal values from 0 to 1999 mV. With an optional pressure gauge (not provided — see **PM015D**, page 116), you may calibrate the display to read mmHg. Recorder output connector allows direct connection to a pen recorder, oscilloscope or computer via a data acquisition system.

Supplied sterile, **BLPR2** is accurate, linear and stable with temperature. May be sterilized cold with Rapicide OPA or a similar bactericide.

**BLPR2** is equipped with a 12-ft cable and connector compatible with WPI's four-channel signal conditioning unit, **TBM4M** Transbridge, and the single-channel **BP1** blood pressure monitor. The cable has a moisture-resistant locking connector. A continuous, uniform lumen eliminates places for bubbles to form and lodge. The clear fluid path is easy to inspect. Easy to mount — slotted transducer body accepts Velcro strap.

To facilitate setup and operation, a four-way stopcock that allows easy filling, flushing and zeroing of the transducer is included. Typically, the stopcock is located between the transducer and the animal catheter where it can be used to quickly zero, flush or de-bubble the transducer.



Stopcock **14036** included with **BLPR2**

## Micro Cannula KZ1101



- 0.4 mm O.D., 0.2 mm I.D. tubing
- Autoclavable
- Biocompatible perfluorocarbon tubing material

See description on page 221.



**SYS-BP1**

## BP1 SPECIFICATIONS

AMPLIFICATION	x1, x10, x100, variable (x5 to x1000)
OUTPUT VOLTAGE SWING	± 5 V
MAXIMUM OUTPUT CURRENT	2 mA
INPUT IMPEDANCE, EACH INPUT	100 kΩ    0.01 μF
TRANSDUCER APPLIED VOLTS	10 V nominal, varies with load. 25 mA, maximum
POWER	95-135 V or 220-240 V, 50/60 Hz
DIMENSIONS	8.5 x 5.12 x 10 in. (21.6 x 13 x 25.44 cm)
SHIPPING WEIGHT	11 lb. (5 kg)

## BLPR2 SPECIFICATIONS

WORKING PRESSURE	-50 to +300 mm Hg
OVERPRESSURE	-400 to +4000 mm Hg
EXCITATION VOLTAGE	1-10 VDC or RMS to 5 kHz
SENSITIVITY	5 μV/mm Hg
DYNAMIC RESPONSE	100 Hz
EIGHT HOUR DRIFT	1 mm Hg after 5 minute warm-up
MAXIMUM ERROR	Total combined effects of Sensitivity, Linearity, Hysteresis (at 25°C and 5 μV/mm Hg) do not exceed ±2% or 1 mm Hg, whichever is greater.
SHIPPING WEIGHT	1 lb.

## ORDERING INFORMATION

**SYS-BP1** Pressure Monitor (transducer & cable not included)  
*Specify line voltage*

### OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>BLPR2</b>	Blood Pressure Transducer & Cable
<b>BPCABLE2</b>	Cable (12 ft) with DIN connector for BLPR2
<b>503067</b>	BLPR2 Transducer without cable
<b>13024</b>	Single Rack Mount Kit
<b>13025</b>	Dual Rack Mount Kit
<b>3491</b>	Extension Cable, 5 ft
<b>500184</b>	BNC-to-BNC Cable, 10 ft
<b>14036-15</b>	4-Way Luer Stopcock, Blue Tint (package of 15)
<b>KZ1101</b>	Micro Cannula, 3 inch

*NOTE: BLPR2 is intended for animal research only and may not be used for human blood pressure measurement.*

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • www.wpiinc.com

Germany: +49 (0)6031 1602171 • wptide@wpi-europe.com • www.wpi-europe.com  
China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net

# Blood Pressure Measurement for Rodents

Unique environmental control system and single tail cuff allows tail movement, minimizing stress in the animal subjects



## Features

- Test up to 200 animals at a time
- Sensor is MRI compatible
- Quick and accurate blood pressure measurement at temperatures as low as 32°C
- Highly sensitive photoelectric sensor for blood pressure detection
- Monitor, record, store or export real time systolic, diastolic, mean and heart rate

## Benefits

- Automatic evaluation of data
- Touch screen or software control
- Data can be collected and stored on and displayed from a USB flash drive, which is included

## Applications

- Non-invasive blood pressure monitoring of rodents

This revolutionary design brings non-invasive blood pressure testing to a new level — a true turn-key system for accurate, consistent blood pressure measurement on mice, rats or any other laboratory animal test subject.

It is a compact, simple yet versatile system that can test from one to 200 animals at a time with independent control of each channel. Simple daisy-chaining allows expansion of up to 200 independently controlled systems. All components are built into one small unit — controls inflation of tail cuff, warming environment with whisper-quiet fans — providing an ideal system for teaching facilities and for the pharmaceutical industry when high throughput is a must.

Single animal systems are controlled from the touch screen, which allows keying in all necessary test setups. Touch screen control allows ease of operation, supplying automatic evaluation of test results — systolic, diastolic, mean and heart rate.

Data is collected, stored, displayed and can be transferred to the supplied memory stick. The USB interface allows for software control of multichannel systems. For single animal systems, built-in software lets you view and export data. Reports are in an Excel-style format and may be easily exported.

No computer is required. However, the analog output may be interfaced with your own data acquisition software.

System is easily cleaned. Removable trays are included with each system.

**Tail Cuff Sensors** have the ability to detect mouse blood pressure or rat blood pressure via the tail photoelectrically. The photoelectric sensor, light source and sensor are built into a unit which attaches to a restrainer with nylon hand screws. Sensors are nonmagnetic. The ambient temperature that a reading can be taken starts at 30°C. The maximum temperature needed for a test is usually 32°C for rats and 34°C for mice. The systolic readings are within five mmHg of Cannulation and Telemetry.

In addition to the standard one-year warranty on the system, tail cuff sensors have a lifetime warranty.



## ORDERING INFORMATION

**II-MRBP-M** Mouse Blood Pressure System

**II-MRBP-R** Rat Blood Pressure System

*Includes one tail cuff sensor and one restrainer (select when ordering).*

*Call for pricing on multiple channels.*

## OPTIONAL ACCESSORIES

<b>II-B64-1</b>	Tail Cuff Sensor, 650-800 grams, ID 1 in.
<b>II-B63-075</b>	Tail Cuff Sensor, 500-650 grams, ID ¾ in.
<b>II-B60-044</b>	Tail Cuff Sensor, 150-380 grams, ID ⅞ in.
<b>II-B60-038</b>	Tail Cuff Sensor, 70-160 grams, ID ⅝ in.
<b>II-B60-025</b>	Tail Cuff Sensor, 20-80 grams, ID ¼ in.
<b>II-78</b>	Rodent Restrainer, 750 to 1000 grams, ID 4" (105 mm), 13" Long
<b>II-79</b>	Rodent Restrainer, 600 to 800 grams, ID 3½" (88 mm), 13" Long
<b>II-80</b>	Rodent Restrainer, 400 to 600 grams, ID 3" (75 mm), 12" Long
<b>II-81</b>	Rodent Restrainer, 250 to 400 grams, ID 2½" (63 mm), 9" Long
<b>II-82</b>	Rodent Restrainer, 180 to 270 grams, ID 2" (50 mm), 8" Long
<b>II-83</b>	Rodent Restrainer, 70 to 170 grams, ID 1½" (38 mm), 6" Long
<b>II-84</b>	Rodent Restrainer, adult mice and neonate rats, 20 to 45 grams, ID 1" (25 mm), 4" Long
<b>II-84XL</b>	Rodent Restrainer, 45 to 70 grams, ID 1¼" (30 mm), 5" Long

# Microprobe Thermometers

*Instrument of choice for biological and laboratory use*



## Features

- Super accuracy
- Fast response
- Analog output signal
- Multiple inputs
- Differential temperature measurement

## Benefits

- Portable
- A variety of probes available

## Applications

- Temperature monitoring
- Cryogenic measure for blood banking
- Skin temperature measurement in exercise experiments

A Microprobe Thermometer is the instrument of choice for biological and laboratory temperature measurements. These thermometers are very versatile, providing fast response, high accuracy and stability with digital display and analog signal for connection to a computer or recorder. With the wide selection of probes, the instruments can be used in almost any application.

**BAT-12:** This thermometer has a sealed construction making it water, dust and fume resistant. The **BAT-12** has a single microprobe input and a single range with the same high accuracy as the **BAT-10**. Comes complete with carrying case.

The thermometers can be used with any "Type T" thermocouple. Select a temperature microprobe on the following page for your specific application.

**BAT-10:** This is the most versatile thermometer available. The instrument has a wide temperature range and fast response with most microprobes. The **BAT-10** accuracy is NIST traceable and in each of the two temperature ranges, the accuracy is the same as the resolution. There are three microprobe inputs, 1 and 2 can be selected as separate inputs while 2 and 3 will read the differential temperature measurement between the two. The instrument has automatic warnings for low battery or faulty probes on the digital display. The linearized analog output (LOP) signal allows ease of connection to a data acquisition system or recorder.

## MICROPROBE THERMOMETER SPECIFICATIONS

	BAT-10	BAT-12
TEMPERATURE RANGE & RESOLUTION	-200°C to +400°C, 1°C resolution -100°C to +199.9°C, 0.1°C resolution	-100°C to +199.9°C, 0.1°C resolution
DIFFERENTIAL TEMP. RANGE	-19.99°C to +19.99°C Linearization centered at 40°C 0.01°C resolution	N/A
ACCURACY	1° Range 0.1° Range Diff. Range	0.1°C ± 1 digit between 0-50°C 0.1% ± 1 digit over full range
REPEATABILITY	± 1 least significant digit	
CALIBRATION CONFORMITY	Conforms to NIST tables	Follows NIST thermocouple tables within 1 digit
DISPLAY	3½ Digit LCD	3½ Digit LCD
INPUT SOCKET	Miniature, quick disconnect, copper-constantan	Miniature, quick disconnect, copper-constantan
ANALOG OUTPUT	Non-linearized set at 1.6 V, corresponding to temperature of 401°C	≈ 10 mV per degree C
POWER SUPPLY/BATTERIES	BAT-10: 4 alkaline "C" cells (life: 1000 hr) BAT-10R: 4 Ni-Cad "C" cells (rechargeable unit)	BAT-12: 9V cell BAT-12R: 9V Ni-Cad with charger
SENSORS	Three Type T thermocouple inputs	One Type T thermocouple input
AMBIENT OPERATING RANGE	15-45°C	Auto-compensated to 0.1°C from 0°C to 50°C
DIMENSIONS	21.6 x 22.9 x 8.9 cm (8.5 x 9 x 3.5 in.)	12.7 x 6.4 x 15.2 cm (5 x 2.5 x 6 in.)
WEIGHT	1.6 kg (3.5 lb), including carrying case	1 kg (2 lb), including carrying case

## ORDERING INFORMATION

**BAT-10R/LOP** Multiple Input Type T Thermocouple Thermometer, rechargeable NiCad batteries and 110 VAC adapter (microprobes ordered separately)

**BAT-10R/LOP220** Multiple Input Type T Thermocouple Thermometer, rechargeable NiCad batteries and 220 VAC adapter (microprobes ordered separately)

**BAT-12R** Single Input Type T Thermocouple Thermometer, rechargeable NiCad batteries and 110 VAC adapter (microprobes ordered separately)

**BAT-12R-220** Single Input Type T Thermocouple Thermometer, rechargeable NiCad batteries and 220 VAC adapter (microprobes ordered separately)

## OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>EXT-6</b>	Probe Extension Lead, 180 cm long
<b>501608</b>	Tripod Stand for BAT-12

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • www.wpiinc.com

Germany: +49 (0)6031 1602171 • wptide@wpi-europe.com • www.wpi-europe.com  
China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net

# Temperature Probes

A variety of probes for any application

## Features

- Flexible Teflon microprobes are used for implantation in tissue, in spectrophotometer cuvettes, rectally in neonatal mice, in water baths, PCR thermal cyclers, etc.
- Monitors animal rectal temperatures during surgical procedures and pyrogen testing


## Benefits

- Thermocouple wire used meets NIST standards (certification is available for an additional cost)

## Applications

- Rectal temperature monitoring
- Instant read in tissue, semi-solids, small specimen and other materials

When precise temperature measurements are required, WPI can provide you with a very accurate monitor and thermocouple microprobes. WPI monitors have both resolution and accuracy of 0.1°C in the 0-50°C range and are traceable to NIST standards, whereas, other competitive electronic thermometers have an accuracy that is usually to 0.5°C or worse. Furthermore, all our type T clinical probes are guaranteed accurate to 0.1°C, due to our stringent wire standards. These are five times more accurate than competitive probes made with regular "Special Limits" wire.



### NEEDLE MICROPROBES ORDERING INFORMATION

Fast-response needle probes for instant readings in tissue, semi-solids, liquids, very small specimens, powders and materials. Needle tip is sealed to ensure only stainless steel contacts specimen. Maximum autoclavable temperature is 135°C.

Probe Type	Size	Style	Time Constant	Isolated	Max. Temp.	Lead Length	Description
MT-29/1HT	29 ga / 1 cm	A	0.125 sec	No	200°C	5 ft	29 gauge approximately 0.013 in.
MT-29/2HT	29 ga / 2 cm	A	0.125 sec	No	200°C	5 ft	-
MT-29/3HT	29 ga / 3 cm	A	0.125 sec	No	200°C	5 ft	-
MT-29/5HT	29 ga / 5 cm	A	0.125 sec	No	200°C	5 ft	-
MT-26/2HT	26 ga / 2 cm	A	0.1 sec	No	200°C	5 ft	26 gauge approximately 0.018 in.
MT-26/4HT	26 ga / 4 cm	A	0.1 sec	No	200°C	5 ft	-
MT-26/6HT	26 ga / 6 cm	A	0.1 sec	No	200°C	5 ft	-
MT-23/3HT	23 ga / 3 cm	A	0.15 sec	No	200°C	5 ft	23 gauge approximately 0.125 in.
MT-23/5HT	23 ga / 5 cm	A	0.15 sec	No	200°C	5 ft	-
MT-23/8HT	23 ga / 8 cm	A	0.15 sec	No	200°C	5 ft	-
MT-4	29 ga / 1 cm	A	0.025 sec	No	200°C	5 ft	Similar to MT-29/1 but has a blunt tip. Good for instant skin and surface temperatures, liquids
MT-D	—	C	0.025 sec	No	200°C	5 ft	Fast response surface probe (stainless steel for locating inflammation, arteries, etc. Also for dental use.

### FLEXIBLE IMPLANTABLE PROBES ORDERING INFORMATION

Designed for high accuracy on extremely small specimens such as insects, seeds, etc. Maximum insertion depth 1/8". Totally sheathed in chemical resistant Teflon.

Probe Type	Sensor Lead Diameter	Style	Time Constant	Isolated	Max. Temp.	Lead Length	Description
IT-14	0.050" dia	D	0.3 sec	Yes	150°C	3 ft	-
IT-18	0.025" dia	D	0.1 sec	Yes	150°C	3 ft	-
IT-18EXLONG	0.025" dia.	D	-	Yes	150°C	5 ft	-
IT-21	0.016" dia	D	0.08 sec	Yes	150°C	1 ft	-
IT-23	0.009" dia	E	0.005 sec	Yes	150°C	3 ft	For ultra fast measurements and for use on micro-size specimens. Tissue implantable with 23ga. Needle (supplied). Rather fragile. Teflon coated.
IT-1E	0.025" dia	F	0.005 sec	Yes	150°C	3 ft	As IT-18 sensor except bead exposed. Combines ultra-fast response of IT-23 with sheath strength of IT-18.


### RECTAL PROBES ORDERING INFORMATION

Probe Type	Sensor Lead Diameter	Style	Time Constant	Isolated	Max. Temp.	Lead Length	Description
RET-2	-	G	0.8 sec	No	125°C	5 ft	Rectal probe for rats typically for fast intermittent measurements. Smooth ball tip (0.125 in. dia.) with 1" long (0.59 in. dia) stainless steel shaft.
RET-3	-	G	0.5 sec	No	125°C	5 ft	Rectal probe for mice similar to RET-2 except tip diam. 1.6mm (0.063") and shaft 3/4" long (0.028" diam.)

### GENERAL PURPOSE ORDERING INFORMATION

Probe Type	Sensor Lead Diameter	Style	Time Constant	Isolated	Max. Temp.	Lead Length	Description
HT-1		H	0.5 sec	No	400°C	5 ft	"Workhorse" probe for liquids, gases, semi-solids. Plastic handle with straight stainless steel shaft. Not good for surface temperatures.
HT-2		H	0.5 sec	No	400°C	5 ft	Like HT-1 except shaft length is 9 in.

Due to the fragility of these probes, warranty is limited:  
Defective probes may be returned within 14 days of receipt.



These probes may also be used with ATC2000 Animal Temperature Controller.  
See page 127.

# Neuroscience Cannulas

For in vivo investigation of rodents

## Features

- Exceptional quality
- Best prices
- Rapid order response

## Benefits

- Beveled inside and out and polished to remove burrs and ensure that the inside diameter is perfectly cylindrical. This limits trauma to tissue and ensures smooth operation.
- Exceptional quality at the best prices

## Applications

- Drug administration
- Optogenetics



Cannula assembly: fixing screw, internal cannula with attached tubing, and guide cannula.



Internal cannula secured with fixing screw.

These cannulas for neuroscience study and pre-clinical research include an entire range of cannula options, including the Guide Cannula, the Internal Cannula and the Dummy Cannula (cap).

### Guide Cannulas



Gauge	OD	ID
22	0.64	0.46
24	0.56	0.40
26	0.48	0.32

The Guide Cannula is a surgical grade, stainless steel tube that is implanted into a rodent's skull and cemented in place using dental cement and screws. It guides the Internal Cannula to the specific injection site.

### Internal Cannula



Gauge	OD	ID
22	0.41	0.26
24	0.36	0.21
26	0.30	0.16

Insert the Internal Cannula into the Guide Cannula to sample or inject fluid.

### Dummy Cannula



The Dummy Cannula has a stainless steel wire core, and it is placed in the Guide Cannula when the Internal Cannula is removed. It seals the opening and prevents tissue from entering the bottom of the Guide Cannula. The Dummy Cannula is threaded to securely tighten it so that the animal will not unscrew

it while grooming. Note that the DUM26 cap is used for the 24 ga cannulas also.


### Ordering

Order the Guide Cannulas based on the gauge and length from the base and the Internal or Dummy Cannulas based on the length of the Guide Cannula and the projection from the Guide Cannula tip.

### Understanding Part Numbers

#### Cannula Type

- GC—Guide Cannula
- INC—Internal Cannula
- DUMC—Dummy Cannula



**Gauge**—Choose the gauge of the Guide Cannula to be used, even if you are ordering an Internal or Dummy Cannula. Choices include 22, 24 and 26 gauge.


**Length**—Guide Cannulas can be ordered in lengths from 1.0–9.9 mm. Length can be specified to 0.1 mm, with a tolerance of ±0.07mm. For example, a **GC24-60** is a 24 gauge cannula that is 6.0 mm long.

For other cannulas, the length is determined by the desired projection length and the guide cannula length. The projection can extend beyond the tip of the Guide Cannula up to 1.0 mm.

$$DUMC\ Length = GC\ Length + Projection$$

$$INC\ Length = GC\ Length + Projection$$

For example, for a Dummy Cannula flush with the end of a 26 ga Guide Cannula that is 6.0 mm long, order a **DUMC26-60**. For an internal cannula that projects 0.5 mm beyond the 6.0 mm Guide Cannula (6.5 mm), order an **INC26-65**.



\* DUMC Length=GC Length + Projection. If the cap is screwed on too tightly, the projection will be longer than expected.  
\*\* INC Length=GC Length + Projection. Internal cannula mounts flush and does not screw into place.

## ORDERING INFORMATION

### Cannulas

<b>GC22-X*</b>	Guide Cannula, 22 ga, X.0 mm
<b>INC22-X*</b>	Internal Cannula, 22 ga, X.0 mm
<b>DUMC22-X*</b>	Dummy Cannula, 22 ga, X.0 mm
<b>GC24-X*</b>	Guide Cannula, 24 ga, X.0 mm
<b>INC24-X*</b>	Internal Cannula, 24 ga, X.0 mm
<b>GC26-X*</b>	Guide Cannula, 26 ga, X.0 mm
<b>INC26-X*</b>	Internal Cannula, 26 ga, X.0 mm
<b>DUMC26-X*</b>	Dummy Cannula, 26 ga, X.0 mm

\* X indicates the length of the cannula or projection. See Understanding Part Numbers.

### Flexible PE Tubing

<b>504278</b>	0.25 mm ID, 0.5 mm OD, 1 m long
<b>504279</b>	0.42 mm ID, 0.85 mm OD, 1 m long
<b>504280</b>	0.6 mm ID, 1.1 mm OD, 1 m long

### Miscellaneous Accessories

<b>504281</b>	Fixing Screw-connect INC & GC
---------------	-------------------------------

Flexible PE tubing is recommended to connect to Internal Cannula (INC) only, NOT Guide Cannula (GC) or Dummy Cap (DUMC).

**504278** matches **INC26-XX**

**504279** matches **INC24-XX**

**504280** matches **INC22-XX**

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • www.wpiinc.com

Germany: +49 (0)6031 1602171 • wpi@wpi-europe.com • www.wpi-europe.com  
China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net



# Electronic von Frey Anesthesiometer

## Assessing mechanical allodynia with 15 SuperTips

### Features

- LCD readout (floating or last maximum/minimum)
- Rigid tips up to 800 g
- Supertips™ up to 65 g
- 1,000 g probe available
- Optional analog output cable for chart recorder
- Pipette tips can be customized to any specification
- Microprocessor electronics 0.1 g plug-in probes

### Benefits

- Plug up to three probes into a single unit
- Independent from temperature

### Applications

- Small animal analgesia testing

To assess mechanical allodynia, which is a painful response to a light touch or pressure from a stimulus that is not normally painful, the Electronic von Frey Anesthesiometer was developed. The Electronic von Frey meter uses one of 15 different flexible von Frey hairs called "SuperTips™" (or rigid tips up to 800 grams). Each hair, regardless of model chosen, is exactly 0.8 mm in diameter. This uniformity of design eliminates false readings and allows for comparison of test results. The Electronic von Frey can be used with chart recorders and analog/digital converters, and it never needs calibration. This system includes either a 90, 800 or 1,000 gram probe. Mesh stands are available in a variety of sizes for large group studies.



II-2391

### ORDERING INFORMATION

<b>II-2390</b>	Electronic von Frey Anesthesiometer, Rigid Tips, 90 g Range
<b>II-2391</b>	Electronic von Frey Anesthesiometer, Rigid Tips, 800 g Range
<b>II-2392</b>	Electronic von Frey Anesthesiometer, Rigid & 15 Supertips, 90 g Range
<b>II-2393</b>	Electronic von Frey Anesthesiometer, Rigid & 15 Supertips, 800 g Range
<b>II-23931</b>	Electronic von Frey Anesthesiometer, Custom Rigid Tips, 1000 g Range
<b>II-2394</b>	von Frey Probe, 90 g Range
<b>II-2395</b>	von Frey Probe, 800 g Range
<b>II-2396</b>	von Frey Probe, 1000 g Range
<b>II-2397</b>	MRI Probe Option ( <i>add to price of probe above</i> )
<b>II-2400</b>	Analog Output Cable

## Quattro

### Four test systems



II-2889

This special package offers four tests, including Electronic von Frey, Plethysmometer, Randall Selitto and the Grip Strength Meter. You get all four test modules and the electronic controller that is interchangeable with all four systems. The electronic controller has up to three inputs, so you can perform up to four unique tests with only one electronic system. If you prefer, you may build your system as you grow. Because of the modular design of these four systems, you need to order only one complete system. Then, the modules for the other three tests, which integrate into the system, can be purchased separately, as needed. The stand and sling for the Randall Selitto test are sold separately.

### ORDERING INFORMATION

**II-2889** Quattro 4-in-1 System

## Trio

### Three test systems



II-2888

Get three test systems in one package with the Trio, featuring the Electronic von Frey, Plethysmometer and Randall Selitto Meters. Just like the Quattro package, the modular design allows these three test systems to communicate with the same electronic controller. The stand and sling for the Randall Selitto test are sold separately.

### ORDERING INFORMATION

**II-2888** Trio 3-in-1 System

# Plethysmometer (Paw Volume Meter)

## Test effectiveness of anti-inflammatory agents

### Features

- No wetting solution needed
- One calibration/year
- Battery-operated or line powered controller
- One-year warranty

### Benefits

- A sensor in the water notes a pressure change when the paw is immersed
- Pressure is calibrated in 0.1 milliliters and displays on the battery-powered monitor

### Applications

- Measure the effectiveness of anti-inflammatory agents and agents to reduce edemic conditions

Test the effectiveness of anti-inflammatory agents and agents to reduce edemic conditions with the Plethysmometer. Simply insert the rat or mouse paw into water. The acrylic stand offers full visibility of the subject throughout the testing.



### ORDERING INFORMATION

**II-520MR** Paw Volume Meter for Mouse & Rat

# Randall Selitto Paw Pressure Meter

## Get instantaneous, live readings

### Features

- Visible force limit indicator
- Portable electronic display
- Battery-operated or line powered controller
- No calibration required

### Benefits

- Hands-free operation with foot switch
- 3-oz. hand held probe with accuracy of 0.5%

### Applications

- Analgesic drug studies

The Randall Selitto Paw Pressure Meter for analgesia testing is digitally controlled. Use the handheld instrument to apply force to an animal's extremity and get instantaneous, live readings. You can even view the last maximum force



800 g pressure applicator

applied during the test. The new limit indicator lets you select the maximum force limit, and then indicates with a warning light when the system reaches that limit. This unit comes with an acrylic stand to allow for easy viewing of the display. Stand and sling are sold separately.

### ORDERING INFORMATION

**II-2500** Randall Selitto Paw Pressure Meter

# Grip Strength Meter

## Measure muscle hyperalgesia in rodents

### Features


- Mouse and rat wire mesh grid plates included
- Maximum force range is 2,000 g (10% over range allowed). Higher ranges are available by special order.
- 1 gram increment readings
- Suction feet on the heavy, anodized base plate resist even large pulling forces

### Benefits

- Battery-operated or line powered controller
- One-year warranty

### Applications

- Measure muscle hyperalgesia in rats and mice



Measure muscle hyperalgesia in rats and mice with the Grip Strength Meter, which gauges the forelimb grip force using a digital force transducer. Simply hold the animal by the tail and gently dangle it over the wire mesh plate until the animal grasps the plate with its forepaws. The force transducer, connected with the wire mesh plate, measures the strength of the animal at the time of the test. The battery-operated, electronic control device calculates the average of three measurements, and it holds the last maximum force in a "peak and hold" type readout until you reset it.

### ORDERING INFORMATION

**II-2200** Grip Strength Meter for Mouse & Rat

### WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • www.wpiinc.com

Germany: +49 (0)6031 1602171 • wptide@wpi-europe.com • www.wpi-europe.com  
China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net

# Incremental Hot Cold Plate Analgesia Meter

## *Latency and threshold-based nociception*

### Features

- Heat or cool, from 0-70 °C
- Ramping temperatures for threshold & latency results
- Rapid increase or decrease in temperature
- Printout of data
- Includes clear animal enclosure
- Plate size 4" x 8"
- Two-year warranty

### Benefits

- Precise programmable digital control
- Temperature stability is 0.1 °C

### Applications

- Analgesia research
- Latency and threshold-based nociception



II-PE34

This safe, humane device for rats and mice is used for latency and threshold-based nociception, ramping temperatures for 0-70 °C. Because this hot cold plate is incremental, it measures latencies of much more than just the strong narcotic agents, broadening dramatically the range of analgesia research with devices of this type. Microprocessor-controlled, the Incremental Hot Cold Plate can heat or cool in increments of 0.1 °C, at a rate of 1-10 °C per minute. With

uniform heating and cooling and upper/lower cut-off limits, this device is predictable and safe. It can also function as a constant temperature plate with great stability (0.1 °C). As soon as a reaction is observed from the chosen paw, the unit reverses to the standby temperature.

#### ORDERING INFORMATION

**II-PE34** Incremental Hot Cold Plate Analgesia Meter

# Incapacitance Meter for Mouse & Rat

## *Test pain and inflammation in small animal hind limbs*

### Features

- Precise programmable digital control
- Start, stop and reset test from controller's front panel
- 180-270 gram holder included (other sizes available)
- Reaction detected automatically
- Manual override of all timer functions
- Alphanumeric readout
- 5 to 999 seconds test period
- Alphanumeric readout
- All functions and parameters entered via keypad
- Two-year warranty

### Benefits

- Overcomes the drawback of placing unnecessary stress on the animals

### Applications

- Test pain and inflammation in the hind limbs of mice, rats or birds

The Incapacitance Meter uses a technique called dual channel weight averaging, which tests both hind limbs. This gives you a clean, stress-free correlation of the paw pressure test. Conduct control and testing of the animal at the same time. Place the animal in



II-600 mR

the holder with its hind limbs resting on the two weight-averaging platform pads. The controller records the average weight (grams) over the test period as the animal shifts its weight from each pad.

#### ORDERING INFORMATION

**II-600MR** Incapacitance Meter for Mouse & Rat

# Hot Plate Analgesia Meter

## Latency testing in rodents

### Features

- Includes both mouse and rat enclosure
- Temperature indicated in 0.1 °C increments
- Holding accuracy is  $\pm 0.1$  °C
- Digitally controlled
- Two-year warranty

### Benefits

- Plate has consistent temperature throughout a procedure, ensuring accurate testing

### Applications

- Latency testing in rats and mice

To use the Hot Plate Analgesia meters, simply place the animals on the black anodized aluminum plate (11" x 10.5" x 0.75", 275 x 263 x 15 mm) and set the plate's surface temperature



II-39

to the desired set point (up to 75 °C). The plate maintains a consistent temperature throughout the test.

### ORDERING INFORMATION

II-39 Hot Plate Analgesia Meter for Mouse & Rat

# Plantar Test Apparatus/Tail Flick Test Analgesia

## Meter to test narcotics and strong non-narcotic drugs

### Features

- Includes three acrylic animal enclosures that each hold two rats or four mice
- Precise programmable digital control
- User-defined humane cutoff feature
- Adjustable beam intensity in 1% increments up to 250 °C
- Reaction is detected automatically
- Alphanumeric readout
- Manual override of all timer functions
- All functions and parameters entered via keypad
- Heated glass option
- Tail temperature monitor option (for use with the Tail Flick Meter)

### Benefits

- Plantar test and tail flick test applications are combined in one platform
- Stimulate other body parts by adjusting the height of the glass

### Applications

- Test the properties of narcotics on unrestrained mice/rats in Plantar mode

This unit, which is designed for testing narcotics and strong non-narcotic drugs, offers both Plantar (Hargreaves Method) and Tail Flick testing with a single unit. Either testing system is also available individually. In plantar mode, the visible light/heat source is directed at the paw or other desired body part, and in tail flick mode it is directed at the subjects' tails. Test up to 12 mice or 6 rats simultaneously. If desired, other animals like cats and rabbits may also be used. Tests are simple to setup. The focused, radiant heat/light source creates a 4 x 6 mm intense



II-336TG

spot. Because the light is visible, you know when the test starts and ends. The equipment is silent (no whining or clicking sounds) to avoid triggering an automatic response in conditioned animals. You can set a humane cutoff timer that automatically shuts off the heat if no response is observed during the designated time frame.

When an animal is placed on cold glass, its reaction time may be slower. This unique system offers a heated glass option that prevents the glass enclosure from acting as a heat sink, giving a more accurate reading. An optional tail temperature monitor can also be selected for use with the Tail Flick Meter. This option actually preheats the tail before experimentation. Once the preset tail temperature is reached, the test and timer automatically begin. A glass stand is also available in two sizes for large group studies.

### ORDERING INFORMATION

II-336T	Combination Plantar/Tail Flick Meter (non-heated glass and tail temperature for mouse and rat)
II-336TG	Combination Plantar/Tail Flick Meter (tail temperature and heated glass for mouse and rat)
II-390	Plantar Test Analgesia Meter (non-heated glass for mouse and rat)
II-390G	Plantar Test Analgesia Meter (heated glass for mouse and rat)

### WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
 Germany: +49 (0)6031 1602171 • wpide@wpi-europe.com • www.wpi-europe.com  
 Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • www.wpiinc.com  
 China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net

## ADHESIVES APPLICATION GUIDE

WPI Part #	Description	Curing Time	Useful Applications and Characteristics
<b>Epoxies</b> Form strong bonding. Used in wire bonding applications.			
<b>4898</b>	Silver Filled Conductive Epoxy	12 hr @ 50° C – 5 min @ 150° C	Connecting conductors that can't be soldered. Constructing or connecting Ag/AgCl pellets.
<b>7335</b>	Carbon Filled Conductive Epoxy	48 hr @ 25° C – 5 min @ 150° C	Constructing carbon electrodes
<b>4886</b>	High Performance Structural Epoxy	12 hr @ 25° C	Forms a strong and slightly flexible bond on plastic, metal, & glass. Bonds some low surface.
<b>Hot melt (EVA)</b> Easy to use for bonding, needs large gap filling			
<b>13316</b>	Mini Glue Gun with glue sticks	As soon as it cools down	Bonds wood, glass, metals, and many plastics
<b>Silicone Adhesives/Sealants/Primers</b> Good moisture resistant and elastic. Low toxic.			
<b>1571</b>	Room Temperature Vulcanizing (RTV) adhesive. Acyloxy/moisture cure system. Acetic acid is cure by-product.	24 hr @ 25° C	Has the best adhesion property in this silicone family. Will bond to many materials.
<b>7128</b>	RTV sealant. Alkoxy/moisture cure system. Methanol is cure by-product.	72 hr @ 25° C	Good for bonding or sealing electronics circuits (metal)
<b>SYLG184</b>	Sylgard. Two-part vinyl/platinum cure sealant. Hydrogen is cure by-product. Very low toxicity.	24 hr @ 25° C – 15 min @ 150° C	Coating patch clamp electrodes, cell culture dish, making dissection pads
<b>KWIK-SIL</b>	Two-part adhesive. Vinyl/platinum system. Hydrogen is cure by-product. Very low toxicity.	< 5 min @ 25° C	Live tissue and nerve studies. Medium strength adhesion.
<b>KWIK-CAST</b>	Two-part sealant. Vinyl/platinum cure system. Hydrogen is cure by-product. Very low toxicity.	< 5 min @ 25° C	Sealant for live tissues. Embedding peripheral nerves with electrodes.
<b>6820</b>	Primer for silicone	N/A	Enhances adhesion of silicone adhesives for difficult to bond plastic surfaces
<b>Cyanoacrylate</b> Forms an instantaneous bonding.			
<b>7341</b>	Ethyl cyanoacrylate, low viscosity 90-120 cps	<10 seconds	Mounting rat/mouse brain slices. Ideal for relatively small gaps and smooth surfaces. Bonds plastic, metals and rubber. Package of 10 vials, each approximately 1.5 mL.
<b>7342</b>	Ethyl cyanoacrylate, high viscosity 1100-1600 cps	<30 seconds	Use in brain slice experiments. Ideal for larger gaps, allows slightly longer bonding time. Bonds plastic, metals and rubber. Package of 10 vials, each approximately 1.5 mL.
<b>VETBOND</b>	Butyl cyanoacrylate, low toxicity	<10 seconds	Bonds tissues, alternative to suture, helps small wound healing. Antimicrobial effect. Used in forensic science.
<b>503763</b>	Octyl cyanoacrylate, low toxicity	<15 seconds	Suitable for surface wound bonding, protection, holding a sensor or other device on the tissue.

## Kwik-Gard™

Kwik-Gard™ is specially packaged **Sylgard 184** silicone for quicker and easier application, eliminating the messy procedure of preparing the mixture before application. Its special cartridge controls the precise mixing ratio to ensure proper curing. The disposable tip mixes resin and hardener as they are dispensed. Since no air is introduced during mixing, the resin does not need degassing for most applications. The mixed silicone is applied directly to the site, reducing preparation time and material waste.

Each Kwik-Gard cartridge contains 37 mL of resin and hardener. The dispensing tip has a dead volume of 0.75 mL.



KWIKGARD

### ORDERING INFORMATION

<b>KWIKGARD</b>	Kwik-Gard™ Start-up Kit (incl. dispenser, 1 cartridge, 5 tips)
<b>KWIKGLUE</b>	Kwik-Gard™ Refill (2 cartridges, 10 dispensing tips)
<b>KWIKMIX</b>	Dispensing Tips (pack of 10)
<b>KWIKGUN</b>	Kwik-Gard™ Dispenser

## Scotch-Weld 2216 Structural Epoxy

Scotch-Weld 2216 remains the best epoxy for bonding plastic, often used as the benchmark for testing the binding strength of other adhesives. The slightly rubbery texture also makes it less easy to break off. *It is the only epoxy known that can bond PEEK.*

Color: gray  
Cures at room temperature.  
Shipping weight: 1 lb. (0.5 kg)



4886

### ORDERING INFORMATION

<b>4886</b>	Scotch-Weld 2216 (2 oz.)
-------------	--------------------------

# Low Toxicity Adhesive

Ideal for neuroscience applications, nerve studies and more

## Features

- Bio-compatible adhesive for live tissue and nerve studies
- Pre-mixing tips simplify use
- Medium strength adhesion
- Low toxicity
- Rapid curing silicone adhesive, cure on contact
- Cures without producing heat
- Includes 10 mixing tips
- Volume discounts - save up to 15%

## Benefits

- Low toxicity
- Rapid cure time

## Applications

- Neuroscience and nerve studies
- Biomedical applications



**Kwik-Sil** and **Kwik-Cast** silicones have very low toxicity before, during and after curing. The by-product of curing is a small amount of hydrogen gas, which is much less toxic to cells than acetic acid or alcohol from traditional RTV silicone systems.

**Kwik-Sil** and **Kwik-Cast** curing speed is hundreds of times faster than traditional RTV silicones. A curing time of a few minutes at room temperature is especially useful for encapsulation of live tissue or implanting into a live animal.

Unlike many vinyl-based silicones in which the platinum complex catalysts are easily poisoned by contamination from amines and animal tissue, **Kwik-Sil** and **Kwik-Cast** are not sensitive to contamination from animal tissue.

**Kwik-Sil™** is a translucent, medium-viscosity silicone adhesive, developed for chronic peripheral nerve studies such as anterograde tracing with fluorescent indicators or electrode recording. Good adhesion and mechanical properties (tear strength and elongation) allow days of study without breaking of the bonding. Curing speed is very reproducible.

**Kwik-Cast™** is a very low viscosity silicone sealant developed to embed peripheral nerves with electrodes for acute multi-fiber recordings. It flows easily, filling the small spaces around the nerve and leaving no channels through which peritoneal fluid can travel and thus short the nerve/electrode contact. Equally important is the ability of the material to flow into itself and create one continuous mass from underneath the nerve all the way to the top of the nerve/electrode contact to ensure long-term recording stability. **Kwik-Cast** is color-coded to make the mixing foolproof. The catalyst is yellow and the base is blue. When uniformly mixed, it is green. **Kwik-Cast** can be applied and cured underneath mineral oil. After recording, electrodes are easily recovered due to the low tear strength.



KWIK-CAST



KWIK-SIL

## References

Maura Casciola, Shu Xiao & Andrei G. Pakhomov. Damage-free peripheral nerve stimulation by 12-ns pulsed electric field. *Scientific Reports* 7, Article number: 10453. <https://www.nature.com/articles/s41598-017-10282-5>

Tiran, E., Ferrier, J., Deffieux, T., Gennisson, J.-L., Pezet, S., Lenkei, Z., & Tanter, M. (2017). Transcranial Functional Ultrasound Imaging in Freely Moving Awake Mice and Anesthetized Young Rats without Contrast Agent. *Ultrasound in Medicine & Biology*, 43(8), 1679–1689. <http://doi.org/10.1016/j.ultrasmedbio.2017.03.011>

Lee, D., & Lee, A. K. (2017). Efficient Method for Whole-Cell Recording in Freely Moving Rodents Using Ultraviolet-Cured Collar-Based Pipette Stabilization. *Cold Spring Harbor Protocols*, 2017(4), pdb.prot095810. <http://doi.org/10.1101/pdb.prot095810>

Welle, T., Hoekstra, A. T., Daemen, I. A. J. M., Berkers, C. R., & Costa, M. O. (2017). Metabolic response of porcine colon explants to in vitro infection by *Brachyspira hyodysenteriae*: a leap into disease pathophysiology. *Metabolomics*, 13(7), 83. <http://doi.org/10.1007/s11306-017-1219-6>

Lee, D., & Lee, A. K. (2017). In Vivo Patch-Clamp Recording in Awake Head-Fixed Rodents. *Cold Spring Harbor Protocols*, 2017(4), pdb.prot095802. <http://doi.org/10.1101/pdb.prot095802>

Suwabe, T., Fukami, H., & Bradley, R. M. (n.d.). Synaptic Responses of Neurons Controlling the Parotid and von Ebner Salivary Glands in Rats to Stimulation of the Solitary Nucleus and Tract. <http://doi.org/10.1152/jn.01115.2007>

## KWIK-CAST & KWIK-SIL SPECIFICATIONS

	Kwik-Sil	Kwik-Cast
MIX RATIO	1 to 1	1 to 1
WORKING TIME	< 5 minutes*	4 minutes
SETTING TIME (ROOM TEMP., 1:1 RATIO)	5–10 minutes**	<10 minutes
CURE TIME	~15 minutes	
VISCOSITY, CPS	15,000	10,000
GUARANTEED SHELF LIFE AT 23 °C	6 months	6 months
VOLUME	5 mL	5 mL
NUMBER OF MIXING TIP	10	10
DEAD VOLUME OF THE MIXING TIP	<0.12 mL	<0.12 mL
<b>AFTER CURING 24 HOURS:</b>		
TEAR STRENGTH, PPI	90	44
ELONGATION %	650	60
DUROMETER (SHORE A-2)	30	36
COLOR	translucent	green
VOLUME RESISTIVITY, W/CM	1x10 <sup>15</sup>	1x10 <sup>15</sup>

\* 3 minutes average with about 90 seconds of liquidity

\*\* no longer mixable at this point

## ORDERING INFORMATION

<b>KWIK-SIL</b>	Silicone Adhesive Compound (two 5-mL syringes)
<b>KWIK-CAST</b>	Silicone Casting Compound (two 5-mL syringes)
<b>600022</b>	Replacement KWIK Mixing Tips (pack of 10)

Quantity discounts available

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • [wpiuk@wpi-europe.com](mailto:wpiuk@wpi-europe.com) • [www.wpi-europe.com](http://www.wpi-europe.com)  
Brazil: 011 55 13 40629703 • [info@brazil.wpiinc.com](mailto:info@brazil.wpiinc.com) • [www.wpiinc.com](http://www.wpiinc.com)

Germany: +49 (0)6031 1602171 • [wptide@wpi-europe.com](mailto:wptide@wpi-europe.com) • [www.wpi-europe.com](http://www.wpi-europe.com)  
China: +86 21 6888 5517 • [chinasales@china.wpiinc.com](mailto:chinasales@china.wpiinc.com) • [www.wpiinc.net](http://www.wpiinc.net)

# "Super" Adhesives for Life Science Research

*Four times stronger than butyl cyanoacrylate and less toxic*

Cyanoacrylate adhesives have been on the market since 1958. Most industrial or household grade cyanoacrylate is made of shorter alkyl chain derivatives such as methyl or ethyl cyanoacrylate (WPI's **7341** and **7342**). They are very useful for temporarily holding tissues such as mounting specimens for microtome sectioning. However, they are not suitable for bonding wounds on live animals. The difficulties of using cyanoacrylate for bonding live animals are: (1) a strong, irritating odor; (2) quick loss of bonding strength due to breakdown of the bonding by hydration; (3) the breakdown products, cyanoacetate and formaldehyde, are toxic and can cause inflammatory reactions; and (4) they have low flexibility and tend to be brittle.

To overcome these problems, several longer alkyl chain cyanoacrylates have been developed especially for veterinary and human use. The first longer alkyl chain product is butyl cyanoacrylate. This product has been used for animal and human applications outside the USA since 1970. It is much less toxic and has a lower odor than the methyl or ethyl cyanoacrylate. The butyl cyanoacrylate offered by WPI is **Vetbond™**.

A family of adhesives containing octyl cyanoacrylate, a plasticizer and stabilizer, was developed in the 1990's (one of them approved by FDA for human use). When bonding to tissue, these new adhesives are four times stronger and less toxic than butyl cyanoacrylate. Compared with the traditional suture, the new super adhesive has several advantages. On average, it takes only one-tenth of the time to close an incision. The bonding strength is equal to 5-0 monofilament suture. It also has a



**VETBOND**

mysterious antimicrobial effect that can decrease infection rates in contaminated wounds. Bonding will slough off naturally in 5 to 7 days. Cosmetic appearance of the healed incision is also better.

Glutire Topical Tissue Adhesive **503763** forms a strong and flexible film and is thus more suitable for surface wound bonding, protection, and holding a sensor or other device on the tissue.

Setting time is about 10 seconds, which gives ample time for application. It can also be used for temporarily holding a live tissue. For example, there is a report of using it to hold nematodes on a glass slide for patch-clamp neurons recording.

All of the products offered by WPI are veterinary grade (not suitable for human application). Though very similar to the grade for human use, they are not sterile and do not have FDA approval.



## ORDERING INFORMATION

<b>503763</b>	Glutire Topical Tissue Adhesive (10 tips), 1.5 mL
<b>7341</b>	Cyanoacrylate Adhesive, low viscosity—90-120 cps (pack of 10 vials, each approx. 1.5 mL)
<b>7342</b>	Cyanoacrylate Adhesive, high viscosity—1100-1600 cps (pack of 10 vials, each approx. 1.5 mL)
<b>VETBOND</b>	3M Vetbond™ Adhesive (3 mL)

## Sylgard



**SYLG184**

A two-part silicone elastomer, ideal for potting and encapsulating applications. Very low dielectric constant sealing compound used in patch clamping and many other lab applications. After cure, will withstand -55° to 200 °C.

## ORDERING INFORMATION

**SYLG184** Sylgard, 0.5 kg (1.1 lb)

## Silicone Dissecting Pad Kit



**501986**

Make your own silicone dissection pads easily and quickly. Mix the 2-part silicone right in the plastic Petri dishes and allow to cure 24 hours at room temperature. Kit includes enough two-part Sylgard silicone elastomer to prepare 20 dishes; pins; and 20 plastic Petri dishes with lids, 65 mm.

## ORDERING INFORMATION

**501986** Silicone Dissecting Pad Kit

## Electrically Conductive Silver Epoxy



Two-component silver-filled epoxy for electrical connections which cannot be soldered, such as Ag/AgCl pellets. This widely used silver-filled epoxy features low viscosity and smooth flowing character. Pure silver is dispersed in both resin and hardener. Cures in 15 minutes at 120 °C. Mix ratio 1:1. May be premixed and frozen for later use.

### ORDERING INFORMATION

**4898** Silver Epoxy, 28 g (1 oz)

## Electrically Conductive Carbon Epoxy



Two-component carbon-epoxy, curable at room and elevated temperatures. Ideal for electrostatic discharge protection and EMI/RFI shielding. 1:1 mix ratio. May be premixed and frozen for later use.

### ORDERING INFORMATION

**7335** Carbon Epoxy, 56 g (2 oz)

## Silicone RTV Adhesive (non-acidic)



Because it is non-corrosive, this material is ideal for use on metal, for encapsulating small circuits on connectors. After cure, will withstand -55° to 200 °C. No mixing required.

### ORDERING INFORMATION

**7128** RTV Coating, 90 mL (3 fl oz)

## Silicone RTV Adhesive



Clear silicone sealant provides good bonding to plastic. After cure, will withstand -55 to 200 °C. No mixing required. A handy, general purpose laboratory sealant. (Releases acetic acid during curing.)

### ORDERING INFORMATION

**1571** RTV Sealant, 139 mL (4.7 fl oz)



## RTV Prime Coat

Enhances adhesion of silicone adhesives to many difficult-to-bond plastic surfaces.

### ORDERING INFORMATION

**6820** RTV Prime Coat, 400 mL (13.5 fl oz)



## Mini Glue Gun

Comes with three sticks of special formula hot melt glue. UL approved. 110V 60 Hz only.

### ORDERING INFORMATION

**13316** Mini Glue Gun

### WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • www.wpiinc.com

Germany: +49 (0)6031 1602171 • wptide@wpi-europe.com • www.wpi-europe.com  
China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net



# Biosensing



## Selection includes high selectivity and low detection limit sensors

WPI's biosensors are unique, because they offer a high selectivity and low detection limit (down to nM concentration) with a broad dynamic range; covering physiological concentration of species with different sizes from nM to mM. The majority of our sensors are the only commercially available sensors in the world. Scientists across a variety of disciplines have relied on our sensors for over 25 years. These scientists use WPI's sensors for research performed in universities, hospitals, biomedical research labs, pharmaceutical companies, food/ nutrition research labs, environmental monitoring centers and military labs. Our popular biosensors are listed in thousands of publications.

# Leading Scientist Heads WPI Biosensing Division




## About Dr. Xueji Zhang

The notable expert in the area of electrochemical sensors, Xueji Zhang, PhD, Sr. Vice President and Chief Scientist at WPI, is leading the development for microsensors (including the **ISO-H2S-100-Cxx**) for many applications. His research includes bio-analysis, electrochemical sensors and biosensors, microelectrodes, ultramicroelectrodes and nanoelectrodes, nanosensors, free radical sensors, nitric oxide sensors, cancer diagnostics, design and application of biomedical instrumentation.

Dr. Zhang has written hundreds of published research papers. Listed below is a select list of Dr. Zhang's papers covering the detection of Hydrogen Sulfide. We are very proud to share the research findings with you. Please email [wpi@wpiinc.com](mailto:wpi@wpiinc.com) for a copy of the paper(s).

- **Zhang XJ, et al.**, Electrochemical Sensors Biosensors and their Applications, Elsevier 2008.
- **Zhang XJ, et al.**, Electrochemical hydrogen sulfide biosensors. *Analyst*. 2016;141:1185-95.
- **Zhang XJ, et al.**, A novel enzymatic method for determination of homocysteine using electrochemical hydrogen sulfide sensor, *Frontiers in Bioscience*, 2007, 12, 3774

Dr. Zhang is a national chair professor and director at the University of Science & Technology, Beijing, China. He is also an adjunct professor of University of Tokyo, University of South Florida, Chinese Academy of Sciences, PLA 301 hospital. He is president of the Chinese Alliance of Biodetection & Biomonitoring. He is the honorary president of MITBJ. He has published over 400 papers in peer-reviewed journals, authored 8 books and 80 patents. He has developed numerous biomedical sensors, instruments and devices for commercialization. Dr. Zhang also serves as the editor-in-chief of *Frontiers in Bioscience* ([www.bioscience.org](http://www.bioscience.org)) and *The American Journal of Biomedical Sciences* ([www.nwpii.com/ajbms](http://www.nwpii.com/ajbms)) and has been a member of the advisory editor board of 19 other international journals. He has received numerous awards and honors, including academican of Russian Academy of Engineering (2013), Fellow of American Institute for Medical & Bioengineering (2016), Fellow of Royal Chemistry Society (2014), Fellow of Chinese Chemical Society (2014), Scientist of the Year in China (2015), Engineer of the Year in China (2016), National Innovation Award (2017) and W. Simon Fellow of ICSC-World lab, United Nations (1996).



### WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • [wpiuk@wpi-europe.com](mailto:wpiuk@wpi-europe.com) • [www.wpi-europe.com](http://www.wpi-europe.com)  
Brazil: (013) 406-29703 • [info@brazil.wpiinc.com](mailto:info@brazil.wpiinc.com) • [www.wpiinc.com](http://www.wpiinc.com)

Germany: +49 (0)6031 1602171 • [wpide@wpi-europe.com](mailto:wpide@wpi-europe.com) • [www.wpi-europe.com](http://www.wpi-europe.com)  
China: +86 21 6888 5517 • [chinasales@china.wpiinc.com](mailto:chinasales@china.wpiinc.com) • [www.wpiinc.net](http://www.wpiinc.net)

# Biosensors Designed for Selective Detection

## Selecting microsensors for your application

Choosing the right sensor for your application is critical for successful research. The manufacturer's specifications provide valuable information for selecting your sensor. Consider the following five performance factors in reference to your application:

**Response**—Electrochemical electrodes produce changes in current in response to changes in concentration. "Response" is most often specified in terms of the amount of current per concentration unit: nA/μM or pA/nM, etc. The larger the current per concentration the higher the sensitivity of the sensor.

**Detection Limit**—Detection limit is the minimum change in concentration that can be reliably seen. This specification is directly related to the noise of the sensor. A sensor with a 100nA/μM response but a 3 μM detection limit is not as good as a 10nA/μM response sensor with a 1 μM detection limit.

**Free Radical Detection**—The best sensors have low detection limits and high sensitivity. A sensor can have a low detection limit and a good response, however, to be useful in long term studies it must be stable

when temperature and concentration are constant. A drifting baseline, if monotonic, can be corrected, but wandering baselines limit the utility of sensors to short experiments.

**Selectivity**—It is a rare instance that the ion species of interest is the only ion in the medium to be measured. In a perfect world your sensor would respond ONLY to the ion of interest. In reality there is always some contribution from competing species. The lower the contribution the better.

**Linearity**—For an electrode to be useful and easy to calibrate the response must be "linear" with changes in concentration over the range of interest. Non-linear behavior requires special curve fit software to calibrate the sensors. This approach is more time consuming and can be unreliable. "Good" linearity is expressed by a R<sup>2</sup> of 0.980 or higher. (1.00 is perfect.) All of the electrochemical sensors made by WPI are 0.98 or better.

The table below presents the specifications of WPI's macro sensors (2 mm sensor) and microsensors.

### MACRO SENSORS

SPECIES	Carbon Monoxide	Nitric Oxide	Hydrogen Peroxide	Oxygen	Hydrogen Sulfide
Order Number	<b>ISO-COP-2</b>	<b>ISO-NOP</b>	<b>ISO-HPO-2</b>	<b>ISO-OXY-2</b>	<b>ISO-H2S-2</b>
Available Diameters	2 mm	2 mm	2 mm	2 mm	2 mm
Response Time	< 10 sec	< 5 sec	< 5 sec	< 10 sec	< 5 sec
Detection Limit/Range	10 nM to 10 μM	1 nM to 40 μM*	< 100 nM to 100 μM	0.1%-100%	< 5 nM to 100 μM
Sensitivity	~0.5 pA/nM	≤ 2 pA/nM	8 pA/μM	0.3-0.6nA/%	2 pA/nM
Drift	< 1 pA/min	< 1 pA/min	< 1 pA/min	≥ 1%/min	< 1 pA/min
Temperature Dependent	Yes	Yes	Yes	Yes	Yes
Physiological Interference	nitric oxide	NaNO <sub>2</sub> (10 <sup>-6</sup> or better)	None	None	None
Replacement Sleeves (pkg. of 4)	<b>95620</b>	<b>5436</b>	<b>600012</b>	<b>5378</b>	<b>600016</b>
Filling Solution	<b>95611</b>	<b>7325</b>	<b>100042</b>	<b>7326</b>	<b>100084</b>
Start-up Kit	<b>95699</b>	<b>5435</b>	<b>600011</b>	<b>5377</b>	<b>600015</b>

\*Higher detection limit available on request — call for custom pricing.

### MICRO SENSORS

Species	Nitric Oxide										H <sub>2</sub> O <sub>2</sub>		H <sub>2</sub> S
	ISO-NOPF200	ISO-NOPF200-Lxx <sup>3</sup>	ISO-NOPF100	ISO-NOPF100-Lxx <sup>3</sup>	ISO-NOPF500-Cxx	ISO-NOP3005	ISO-NOP3020	ISO-NOP007	ISO-NOPNM	ISO-HPO-100	ISO-HPO-100-L	ISO-H2S-100-Cxx	
Package Quantity	(pkg. of 2)	(pkg. of 2)	(pkg. of 2)	(pkg. of 2)	(pkg. of 2)	(pkg. of 3)	(pkg. of 3)	(pkg. of 3)	(pkg. of 3)	(pkg. of 3)	(pkg. of 3)	(pkg. of 2)	(pkg. of 2)
Fiber Diameter (μm)	200	200	100	100	500	30	30	7	7 Conical tip: 100nm	100	100	100	
Tip Length <sup>2</sup> (mm)	1-5 <sup>1</sup>	1-10 <sup>1</sup>	1-5 <sup>1</sup>	1-5 <sup>1</sup>	5-10	0.5	2	2	2	1-5 <sup>1</sup>	1-5 <sup>1</sup>	2-5 <sup>1</sup>	
Response Time (sec.)	< 5	< 5	< 5	< 5	< 10	< 3	< 3	< 3	< 3	< 5	< 5	< 5	
Lowest Detection Limit/Range (nM)	0.2	0.2	0.2	0.2	0.2	1	1	0.5	0.5	1	1	< 5	
Nominal Sensitivity-New Sensor <sup>2</sup> (pA/nM)	≥ 10	≥ 20	≥ 6	≥ 20	≥ 20	≥ 0.2	≥ 1	≥ 0.5	≥ 0.5	≥ 1	≥ 1	1-4	
Baseline Drift (pA/min.)	none	none	none	none	none	none	none	none	none	< 2.0	< 2.0	< 2	

<sup>1</sup>Sensor available in 1 mm length increments (for example, 1 mm, 2mm, 3mm...)

<sup>2</sup>Sensor sensitivity varies with length and diameter.

Any 100 μm sensor can be purchased with a hypodermic sheath. Add **-H** to the end of the part number (for example, **ISO-HPO-100-H**).

<sup>3</sup>L-shaped sensor for use with a tissue bath.

Some nitric oxide sensors are available in custom lengths. When ordering custom lengths, use the part numbers **ISO-NOPF100-Cxx** or **ISO-NOPF200-Cxx** and replace the **xx** with the desired length. For example, for a 1 mm flexible sensor tip, the part number should be **ISO-NOPF200-C01**. Sensors can be ordered in the following custom lengths: 1 mm, 2mm, 3mm, 4mm or 5 mm (**ISO-H2S-100**: 2-5mm only).

# Nitric Oxide Sensors

For routine detection of NO at ultra low concentrations

## Features

- Excellent selectivity to NO
- Rapid response time
- Highly sensitive

## Benefits

- NO sensors in different sizes from 100 nm to a few mm can be used for many NO detection applications, such as single cell measurement, *in vivo* measurement in tissues (even in animals), and NO release from drugs

## Applications

- Cell culture, cell suspensions
- Measurements at the cellular level
- Arteries, microvessels, *in vivo* applications
- Tissue bath applications

WPI offers the most extensive range of nitric oxide (NO) sensors available. Developed over a decade of extensive research in the field of nitric oxide, the result is a superior range of NO sensors that enable routine detection of nitric oxide at ultra low concentrations.

The ideal NO sensor should be insensitive to other reactive species likely to be present within the measurement environment. The conventional Nafion coated carbon fiber NO micro sensor exhibits a large response to such species. WPI's unique NO sensor technology utilizes a novel surface membrane which amplifies the response to NO while eliminating responses to a vast range of reactive species, including nitrite, ascorbic acid, hydrogen peroxide, catecholamines, and much more.

## ISO-NOP


The original nitric oxide probe – ideal for cell cultures, cell



### suspensions and many other applications

The ISO-NOP is a popular, robust and high performance sensor encased within a 2 mm diameter disposable stainless steel protective sleeve. The tip of the sleeve is covered with a NO-selective membrane. Replacement membrane sleeves can be purchased separately (5436) and require an internal electrolyte (7325).

A simple change in experimental protocol enables the ISO-NOP to be conveniently used for indirect rapid accurate determination of nitrite (NO<sub>2</sub>) and nitrate (NO<sub>3</sub>) concentration in samples. Using this method a detection limit for NO<sub>2</sub> or NO<sub>3</sub> as low as 1 nM is routinely possible.



Abdominal X-ray showing the apparatus consisting of two customized ISO-NOP nitric oxide probes, 4-channel pH catheter, and Teflon nasogastric tube. (Courtesy Prof. K.E.L. McColl, University Department of Medicine and Therapeutics, Western Infirmary, Glasgow, Scotland.) Iijima, K., et al. *Gastroenterology* 2002; **122**: 1248-1257.

## References


- Zhang X.J. Real time and *in vivo* monitoring of nitric oxide by electrochemical sensors—From dream to reality, *Frontiers in Bioscience*, 9,3434-3446(2004)
- Zhang X.J., Ju H.X., Wang J. Electrochemical Sensors, *Biosensors and Their Biomedical Applications*, Elsevier, 2008
- Zhang X.J., et al. Nanometer size Electro for nitric oxide and s-nitrosothiols measurement, *Electrochem. Commun.* 4, 11-16 (2002)
- Zhang X.J. Construction and Characterization of a new *in vivo* nitric oxide microsensor, *Electroanal.* 9, 640-643, (2004)
- Zhang X.J. An integrated nitric oxide sensor based on carbon fiber electrode coated with selective membranes, 12. 1113-1117(2000)

## ISO-NOPNM


The world's smallest nitric oxide NanoSensor, designed for measurement of NO at the cellular level.

The ISO-NOPNM NanoSensor has a tip diameter of just 100 nm (0.1 μm) and a detection limit for NO of less than 0.5 nM — making it indisputably the smallest and most sensitive NO sensor in the world!

The ISO-NOPNM is based on a novel design in which an electrochemically



“activated” composite graphite nanofiber is used as the NO-sensing element. The surface of the NanoSensor is then modified using a




Amperometric response of the NO nanosensor (ISO-NOPNM) to the successive additions of 2 nM, 4 nM, 8 nM NO into 0.1 M CuCl<sub>2</sub>.

unique multi-layered NO-selective membrane. Figure at right illustrates the response of the ISO-NOPNM following successive additions of nanomolar concentrations of NO. The ultra-low noise of the ISO-NOPNM (0.5 pA) enables a detection limit of just 0.5 nM NO. The response time of ISO-NOPNM is less than 3 s.


## ISO-NOPF

Unique flexible NO sensor! Designed for arteries, microvessels, *in vivo* applications and similar applications.

ISO-NOPF electrodes are available in 100 μm and 200 μm diameters. Utilizing the latest advances in nano-technology and material science, scientists at WPI's Sensor Laboratory have created these completely flexible and virtually unbreakable NO sensors. The new sensors are based on a composite graphite NO-sensing element combined with a



reference electrode. The surface of the sensor is then coated with a unique multi-layered NO-selective membrane.




## NITRIC OXIDE SENSOR SPECIFICATIONS

	ISO-NOPF	ISO-NOP	ISO-NOPF100-L10	ISO-NOPF200-L10	ISO-NOP30	ISO-NOP007	ISO-NOPNM
APPLICATION	In Vivo	Cell Cultures, NO <sub>2</sub> , NO <sub>3</sub>	Tissue Bath	Tissue Bath	Microvessels	Microvessels	Single Cell
SENSOR DIAMETER	100, 200, or 500 μm	2 mm	30 μm	70 μm	30 μm	7 μm	100 nm
RESPONSE TIME	< 5 sec	< 5 sec	< 3 sec	< 3 sec	< 3 sec	< 3 sec	< 3 sec
LOWEST DETECTION LIMIT	0.2 nM	1 nM	1 nM	1 nM	1 nM	0.5 nM	0.5 nM
TEMPERATURE SENSITIVITY	yes	yes	yes	yes	yes	yes	yes
DRIFT	none	none	none	none	none	none	none
SENSITIVITY	8 pA/nM	1–2 pA/nM	≥1 pA/nM	≥1 pA/nM	1.4 pA/nM	0.5 pA/nM	0.5 pA/nM
PHYSIOLOGICAL INTERFERENCE	none	none	none	none	none	none	none


All WPI NO sensors are 100% compatible with the TBR4100 and TBR1025 Free Radical Analyzers. They are also compatible with the discontinued ISO-NO Mark II (NOMK2), APOLLO 4000 and APOLLO 1000.

### ISO-NOP007, ISO-NOP30xx


7 and 30 micron sensors with exceptional performance — ideal for tissues and microvessels



The ISO-NOP007 has a tip diameter of just 7 microns and a length of 2 mm. The ISO-NOP30 has a tip diameter of 30 microns and is available in two different tip lengths (i.e., ISO-NOP3020 has tip length of 2 mm, ISO-NOP3005 has tip length of 0.5 mm). The response of the ISO-NOP007 and ISO-NOP30 is linear over a wide dynamic concentration range of NO. The design of both electrodes is based on a single carbon fiber coated with WPI's NO-selective membrane. A detection limit of approximately 1 nM NO makes these electrodes ideal for use in tissues and microvessels.



The response of a 7 μm NO sensor (ISO-NOP007) to successive additions of NO (100 nM). Inset shows the linearity of the resulting calibration plot.




Simultaneous measurement of force (top trace) and changes of NO concentration (lower trace) in (A) the rat superior mesenteric artery relaxed with ACh and (B) a small human artery relaxed with ACh and SNAP using ISO-NOP30. In this artery oxyhaemoglobin (oxyHb) partly reversed the increase in NO concentration, with only a small change in force. [U. Simonsen, et al., J. Physiol., 1999, 516: 271-282.]

### ISO-NOPF100-L10, ISO-NOPF200-L10, ISO-HPO-100-L

L-shaped sensors for tissue bath & cell culture studies

These unique L-shaped nitric oxide sensors are designed specifically for use in tissue bath studies and similar applications. The shape of the sensor has been engineered to facilitate placement of the electrode within the lumen of the tissue vessel under study. The ISO-NOPF200-L10 is designed specifically for cell culture studies.



### ORDERING INFORMATION

ISO-NOP	Replacement 2 mm shielded sensor and cable
ISO-NOPNM*	100 nm NanoSensor, pkg. of 3 (requires cable 91580)
ISO-NOPF100*	100 μm Flexible NO Sensor, pkg. of 2
ISO-NOPF200*	200 μm Flexible NO Sensor, pkg. of 2
ISO-NOPF500*	500 μm Flexible NO Sensor, pkg. of 2
ISO-NOP007*	7 μm Nitric Oxide Sensor (pkg. of 3)
ISO-NOP3020*	30 μm Sensor Tips (2 mm length), pkg. of 3 (requires 91580)
ISO-NOP3005*	30 μm Sensor Tips (0.5 mm fiber), pkg. of 3 (requires 91580)
ISO-NOPF100-L10*	NO Sensor, 100 μm Flexible L-shaped (pkg. of 2)
ISO-NOPF200-L10*	NO Sensor, 200 μm Flexible L-shaped (pkg. of 2)
ISO-HPO-100-L*	HPO Sensor, L-Shaped 100 μm (pkg. of 2)
5435	ISONOP Startup Kit (recommended with first purchase)
5436	Replacement Sleeve Kit for 2mm sensor, pkg. of 4
7325	ISO-NO Electrolyte (10 mL)
7521	ISO-NO Electrolyte, CO2-insensitive (10 mL)
5399	T-Adapter Kit (pkg. of 3) for ISO-NOP
7357	Nitrite Standard Solution, 1 gram/liter (100 mL)
91580	Microsensor Adapter Cable
SNAP50	SNAP, 50 mg vial

\*Tip diameters given do not include the coatings. Tips are bare wire.

# Hydrogen Sulfide Sensors

Low detection limit sensors to record *in vitro*

## Features

- High selectivity and sensitivity to H<sub>2</sub>S with fast response time
- Broad linearity


## Benefits

- A dry sensor that works like a 2 mm wet sensor
- Use like an ISO-H2S-2 without sleeves and filling solution


## Applications

- Measure biological H<sub>2</sub>S *in vivo* or *in vitro* in tissues or biological media
- Measure H<sub>2</sub>S released from drugs
- Cell culture, cell suspensions, arteries, *in vivo* applications

Although hydrogen sulfide (H<sub>2</sub>S) is generally thought of as a poisonous gas, it is endogenously produced in many mammalian tissues. It has



Stepped response to increasing concentrations of H<sub>2</sub>S are linear (R=0.9963).



The sensor is insensitive to competing species such as ascorbic acid, dopamine, nitrate, nitrite and glutathione.

been detected in micromolar amounts in blood and brain tissue. Hydrogen sulfide is reported as having a broad range of biological functions and although its potential to participate in cell signaling is clear, this biological role is not well understood. H<sub>2</sub>S is strongly analogous to nitric oxide (NO) because they share several physical and metabolic properties.

Like NO, H<sub>2</sub>S is a potent vascular signal that can mediate vasoconstriction or vasorelaxation depending on the O<sub>2</sub> level and tissue. In the rat aorta, H<sub>2</sub>S concentrations that mediate rapid constriction at one O<sub>2</sub> level will cause rapid relaxation at lower O<sub>2</sub> levels.

The ISO-H2S sensor is a low detection limit sensor to record H<sub>2</sub>S *in vitro*. This is the only sensor available that measures H<sub>2</sub>S. The ISO-H2S-100 is a hydrogen sulfide sensor with a 100 μm diameter tip. It is designed like the dry microsensors, however, it works like a traditional 2 mm sensor. The sensor can be ordered in a variety of lengths from 2–5 mm. It incorporates WPI's proprietary combination electrode technology in which the hydrogen sulfide-sensing element and separate reference electrode are encased within a single shielded sensor design. The ISO-H2S-100 offers several advantages:


- Requires no sleeves or filling solutions like 2mm macrosensors
- Durable for long-term use, because of its platinum wire construction
- Rapid response time
- Broad linear range (Range is based on the length of the sensor tip)

## ORDERING INFORMATION

<b>ISO-H2S-100-Cxx</b>	Hydrogen Sulfide Micro Sensor (pkg. of 2)
<b>ISO-H2S-2</b>	2 mm Shielded Hydrogen Sulfide Macro Sensor

# Hydrogen Peroxide Sensors

Direct quantitative measurements in biological samples



## Features

- Gas permeable polymer membrane sleeve blocks liquids, ions and particulate matter
- Incorporated reference electrode

## Benefits

- Non breakable integrated hydrogen peroxide sensor with tip dimension of 100 μm, detection limit down to nM range. Tip size can be custom made as small as a few micrometers.

## Applications

- Cell culture, cell suspensions
- Cell tissue measurements

Despite the recognized importance of this oxidant in biology, real-time measurements at low concentration have been difficult. The hydrogen peroxide sensors developed by WPI are designed to complement existing high sensitivity fluorescent approaches with direct quantitative measurement in biological samples in the low nM range.

The ISO-HPO-2 is a 2.0 mm stainless steel sensor, with replaceable membrane sleeves (600012) and an internal refillable electrolyte (100042). It is designed for use in cell cultures and similar applications.

The ISO-HPO-100 is a 100 μm tip diameter hydrogen peroxide micro

sensor designed for use in tissues and similar applications. The design is based on a platinum wire sensing electrode coated with a proprietary membrane to enhance H<sub>2</sub>O<sub>2</sub> detection.

These sensors incorporate WPI's proprietary combination electrode technology whereby the hydrogen peroxide sensing element and separate reference electrode are encased within a single Faraday-shielded probe design. This design has been shown to enhance performance during measurements and minimizes overall sensor size. Our hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) sensors work with the TBR4100 and TBR1025 free radical analyzers.

## H<sub>2</sub>O<sub>2</sub> SENSOR SPECIFICATIONS

	ISO-HPO-2	ISO-HPO-100	ISO-HPO-100 H	ISO-HPO-100-L
APPLICATION	Cell Cultures, Cell Suspensions	Tissue/Microvessels	Hypodermic Sheath	Tissue Bath
SENSOR DIAMETER	2.0 mm	100 μm	100 μm	100 μm
RESPONSE TIME	< 5 SEC (90%)	< 5 SEC (90%)	< 5 SEC (90%)	< 5 SEC (90%)
DETECTION LIMIT	< 100 nM to 100 μM	1 nM to 1 mM	< 10 nM to 1 mM	1 nM to 100 μM
DRIFT	< 0.2 pA/min	< 2.0 pA/min	< 2.0 pA/min	< 2.0 pA/min
SENSITIVITY	8 pA/μM	1 pA/nM	1 pA/nM	1 pA/nM
PHYSIOLOGICAL INTERFERENCE	none	Contact WPI	Contact WPI	Contact WPI

## ORDERING INFORMATION

<b>600011</b>	ISO-HPO Startup Kit (recommended with first purchase)
<b>ISO-HPO-2</b>	2 mm Shielded HPO Sensor & Cable
<b>ISO-HPO-100</b>	100 μm HPO Sensor*, pkg. of 3
<b>ISO-HPO-100-L</b>	100 μm HPO Sensor, L-shaped*, pkg. of 2
<b>ISO-HPO-100H</b>	100 μm HPO Sensor in hypodermic sheath*, pkg. of 3
<b>600012</b>	Replacement Sleeve Kit for ISO-HPO-2, pkg. of 4
<b>100042</b>	ISO-HPO-2 Electrolyte (10 mL)
<b>91580</b>	Microsensor Adapter Cable

\* Requires 91580 Microsensor Adapter Cable

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
 Brazil: (013) 406-29703 • info@brazil.wpiinc.com • www.wpiinc.com

Germany: +49 (0)6031 1602171 • wpide@wpi-europe.com • www.wpi-europe.com  
 China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net

# Oxygen Sensors

Make direct quantitative measurements in biological samples



ISO-OXY-2

## Features

- Gas permeable polymer membrane sleeve blocks liquids, ions and particulate matter
- Incorporated reference electrode

## Benefits

- Clark type of oxygen sensor with tip size of 2 mm, can be used for very small volume O<sub>2</sub> measurement


## Applications

- Cell culture, cell suspensions
- Cell tissue measurements

This sensor incorporates WPI's proprietary combination electrode technology whereby the oxygen-sensing element and separate reference electrode are encased within a single shielded sensor design. A gas-permeable polymer membrane is fitted over the end of the sleeve, which allows oxygen to pass while blocking liquids, ions and particulate matter.

Oxygen diffuses through the membrane. The voltage applied to the sensor is held at -0.7V when the monitoring device is on and the sensor is properly connected. The magnitude of the generated electrical current is determined by the rate of diffusion through the membrane. The rate is proportional to the partial pressure of oxygen outside the membrane. The current serves as a measure of the partial pressure of O<sub>2</sub>.

The **ISO-OXY-2** is a 2.0 mm stainless steel sensor, with replaceable membrane sleeves (5378) and an internal refillable electrolyte (**7326**).



## O<sub>2</sub> SENSOR SPECIFICATIONS

APPLICATION	Cell cultures, cell suspensions.
SENSOR DIAMETER	2.0 mm
RESPONSE TIME	<10 SEC (90%)
DETECTION LIMIT	0.1 % to 100%
DRIFT	<1%/min
SENSITIVITY	N/A
PHYSIOLOGICAL INTERFERENCE	None

## ORDERING INFORMATION

<b>ISO-OXY-2</b>	2 mm Shielded Oxygen Sensor & Cable
<b>5377</b>	ISO-OXY Startup Kit (recommended with first purchase)
<b>5378</b>	Replacement Electrode Sleeve Kit, pkg. of 4
<b>7326</b>	ISO2 Filling Solution (electrolyte)

# Implantable Glucose Sensor

Measure glucose in vivo or in vitro over long term

## Features

- Implantable sensor for long term studies
- Incorporated reference electrode

## Benefits

- Implantable microsensor,
- Biocompatible
- Long term monitoring



IGS100

## Applications

- *In vivo* long term measurement of glucose in animals or tissues

Measuring glucose *in vivo* over the long term is challenging and difficult. Previous measurement systems were limited to acute studies or a few days at best. WPI introduces a new kind of implantable glucose sensor based on a patented technology. This sensor provides a tool for researchers to directly detect glucose in chronic studies *in vitro* or *in vivo*. The sensor is fully compatible with WPI's TBR systems.

## GLUCOSE SENSOR SPECIFICATIONS

<i>IN VITRO</i> PRECISION	Coefficient of Variation (CV) $\leq 5\%$
GLUCOSE RANGE	36-450 mg/dl (or 2-25 mM/L)
RESPONSE TIME (SEC)	100-300 s
<i>IN VIVO</i> CALIBRATION	<i>In vivo</i> calibration
INTERFERENCE SPECIES	Acetaminophen, ascorbic acid, uric acid
LENGTH	5 cm
SENSOR SIZE	0.6 × 0.7 mm
REFERENCE ELECTRODE	Ag/AgCl
POLARIZATION VOLTAGE (V)	0.65 -0.7 V vs. Ag/AgCl
SENSOR LIFE	3-4 months in solutions at room temperature under continuous polarization; 15-30 days <i>in vivo</i>
SHELF LIFE	6 months
OPERATION CONDITIONS	20° to 40° C (68° to 104° F)
STORAGE CONDITIONS	10° to 25° C (50° to 77° F)

## ORDERING INFORMATION

<b>IGS100</b>	Implantable Glucose Sensor (pkg. of 2)
<b>91580</b>	Microsensor Adapter Cable

# Four-Channel Free Radical Analyzer

*Fast, reliable, real-time detection— Measure redox-reactive species*



TBR4100

## Features

- Real-time detection using electrochemical sensors
- Integrated system includes one temperature sensor, your choice of two additional sensors and a start-up kit
- Current measurement range from 300 fA to 10  $\mu$ A (four ranges) permits wide dynamic range for detection
- Wide bandwidth allows recording of fast events
- Measure carbon monoxide from 10 nM to 10  $\mu$ M
- Measure nitric oxide from < 0.3 nM to 100  $\mu$ M
- Measure hydrogen peroxide < 10 nM to 100 mM
- Measure hydrogen sulfide
- Measure glucose
- Measure oxygen from 0.1% to 100%
- Isolated architecture allows Lab-Trax interface to simultaneously measure free radical and independent analog data (for example, ECG, BP, etc.) on any channel

## Benefits

- Measure up to four different species and temperature in the same preparation or simultaneous measurement in four different preparations
- Lab-Trax data acquisition system is flexible

## Applications

- Free radical detection (NO, H<sub>2</sub>O<sub>2</sub>, H<sub>2</sub>S, CO, O<sub>2</sub> and glucose)

Real-time detection and measurement of a variety of redox-reactive species is fast and easy using the electrochemical (amperometric) detection principle employed in the **TBR4100**. This optically isolated four-channel free radical analyzer has ultra low noise and independently operated channels.

### Measure four species simultaneously

For use with WPI's wide range of nitric oxide, hydrogen peroxide, hydrogen sulfide and oxygen sensors, the **TBR4100** can measure four different species simultaneously in the same preparation. Simply plug a sensor into any one of the input channels on the front panel and select the current range. Poise voltage can be selected from a range of values tuned for optimal response from WPI sensors. An independent output for real-time monitoring of temperature is also included.


## Lab-Trax data acquisition system is flexible

The **TBR4100** analyzer utilizes PC-based data acquisition via our **Lab-Trax** interface. Data traces are displayed and recorded in real-time. The LabScribe software comes pre-configured for single or multiple electrode recording; filters, gains and smoothing are all set for optimal results. Data can be viewed making adjustments to smoothing and filter settings without affecting the original stored raw data. Electrode calibration from multiple concentration readings can be input into the software's Multipoint Calibration utility quickly, providing a plot and slope calculation for electrode sensitivity determination. Alternately, the **Lab-Trax** data interface can be used for providing simultaneous acquisition of Free Radical data along with other physiological data (ECG, HR, BP, etc.) as each of the four input channels has its own independent input, filters and 24-bit converter. See [www.wpiinc.com/TBR4100](http://www.wpiinc.com/TBR4100) for more information on Lab-Trax data acquisition.

## References

- Zhang X.J. Real time and *in vivo* monitoring of nitric oxide by electrochemical sensors—From dream to reality, *Frontiers in Bioscience*, 9,3434-3446(2004)
- Zhang X.J. Ju HX. Wang J. Electrochemical Sensors, *Biosensors and Their Biomedical Applications*, Elsevier, 2008
- Murine strain differences in inflammatory angiogenesis of internal wound in diabetes. (2017). *Biomedicine & Pharmacotherapy*, 86, 715–724. <https://doi.org/10.1016/j.biopha.2016.11.146>
- Xie, L., Feng, H., Li, S., Meng, G., Liu, S., Tang, X., ... Ji, Y. (2016). SIRT3 Mediates the Antioxidant Effect of Hydrogen Sulfide in Endothelial Cells. *Antioxidants & Redox Signaling*, 24(6), 329–343. <http://doi.org/10.1089/ars.2015.6331>
- Nguyen, T.-K., Selvanayagam, R., Ho, K. K. K., Chen, R., Kutty, S. K., Rice, S. A., ... Boyer, C. (2016). Co-delivery of nitric oxide and antibiotic using polymeric nanoparticles. *Chem. Sci.*, 7(2), 1016–1027. <http://doi.org/10.1039/C5SC02769A>
- Chen, G., Yang, L., Zhong, L., Kutty, S., Wang, Y., Cui, K., ... Bin, J. (2016). Delivery of Hydrogen Sulfide by Ultrasound Targeted Microbubble Destruction Attenuates Myocardial Ischemia-reperfusion Injury. *Scientific Reports*, 6, 30643. <http://doi.org/10.1038/srep30643>
- Xu, T., Scafa, N., Xu, L.-P., Zhou, S., Abdullah Al-Ghanem, K., Mahboob, S., ... Zhang, X. (2016). Electrochemical hydrogen sulfide biosensors. *The Analyst*, 141(4), 1185–1195. <http://doi.org/10.1039/C5AN02208H>
- Zong, Y., Huang, Y., Chen, S., Zhu, M., Chen, Q., Feng, S., ... Jin, H. (2015). Downregulation of Endogenous Hydrogen Sulfide Pathway Is Involved in Mitochondrion-Related Endothelial Cell Apoptosis Induced by High Salt. *Oxidative Medicine and Cellular Longevity*, 2015, 1–11. <http://doi.org/10.1155/2015/754670>
- Cho, Y., Park, Y. M., Barate, A. K., Park, S.-Y., Park, H. J., Lee, M. R., ... Holden, D. (2015). The role of rpoS, hmp, and ssrAB in Salmonella enterica Gallinarum and evaluation of a triple-deletion mutant as a live vaccine candidate in Lohmann layer chickens. *Journal of Veterinary Science*, 16(2), 187. <http://doi.org/10.4142/jvs.2015.16.2.187>





Multipoint electrode calibration and slope determination can be quickly derived from recorded calibration data using Lab-Trax data acquisition software.

### TBR4100 SPECIFICATIONS

POWER	100 – 240 VAC, 50-60 Hz, <15 W		
OPERATING TEMPERATURE (ambient)	0–50°C (32–122°F)		
OPERATING HUMIDITY (ambient)	15 – 70% RH non-condensing		
WARM UP TIME	<5 minutes		
DIMENSIONS	135 X 419 X 217 mm (5.25" X 16.5" X 8.16")		
WEIGHT	1.35 kg (3 lb.)		
DISPLAY FUNCTIONS	18 mm (0.7") LCD readout, 4.5 digit Polarization Voltage (mV) Current input (nA, $\mu$ A)		
CONTROLS	Power (on/off) Current Input Range Polarization Voltage		
ANALOG OUTPUT RANGE	$\pm$ 10 V (continuous)		
ANALOG OUTPUT IMPEDANCE	10 K $\Omega$		
CHANNEL TO CHANNEL ISOLATION	>10 G $\Omega$		
CHANNEL TO OUTPUT ISOLATION	>10 G $\Omega$		
POWER SUPPLY TO AC LINE ISOLATION	>100 M $\Omega$		
ANALOG OUTPUT DRIFT	<10 pA/h		
<b>TEMPERATURE INPUT</b>			
NUMBER OF CHANNELS	1		
SENSING ELEMENT	Platinum RTD, 1000 $\Omega$		
RANGE	0–100°C		
ACCURACY	$\pm$ 1°C		
RESOLUTION	0.1°C		
ANALOG OUTPUT	31.25 mV/°C (continuous)		
<b>AMPEROMETRIC INPUT</b>			
NUMBER OF AMPEROMETRIC CHANNELS	4		
SIGNAL BANDWIDTH	0-3 Hz		
POLARIZATION VOLTAGE (SELECTABLE VIA ROTARY SWITCH)			
Nitric Oxide	865 mV		
Hydrogen Sulfide	150 mV		
Hydrogen Peroxide	450 mV		
Glucose	600 mV		
Oxygen	700 mV		
ADJ (user adjustable)	$\pm$ 2500 mV		
POLARIZATION VOLTAGE ACCURACY	$\pm$ 5 mV		
POLARIZATION VOLTAGE DISPLAY RESOLUTION	$\pm$ 1 mV		
<b>CURRENT MEASUREMENT PERFORMANCE</b>			
Range	Analog Output	Noise @ 3Hz	Noise @ 0.3 Hz*
$\pm$ 10 nA	1 mV / 1 pA	< 1 pA	< 0.3 pA
$\pm$ 100 nA	1 mV / 10pA	< 7 pA	< 3 pA
$\pm$ 1 $\mu$ A	1 mV / 100pA	< 70 pA	< 30 pA
$\pm$ 10 $\mu$ A	1 mV / 1 $\mu$ A	< 700 pA	< 300 pA

\*Instrument performance is measured as the (max-min) over 20-second period with open input. Typical values are given at 3 Hz and 0.3 Hz bandwidth.

Typical sensor performance with TBR4100: ISO-NOPF100 NOISE . . . . . 0.2 nM NO (<2 pA)\*\*

\*\*Sensor noise is measured as the (max-min) over a 20-second period with the sensor immersed in 0.1 M CuCl<sub>2</sub> solution.

### Temperature Sensor



The temperature sensor (ISO-TEMP-2) is based on a 2.0 mm tip diameter high quality miniature platinum RTD (Resistance Temperature Detector) electrode. This design has been shown to provide greater accuracy, stability and interchangeability during temperature measurements than traditional thermistor and thermocouple sensors. The ISO-TEMP-2 is included with the purchase of a system.



TBR1025

*Don't need four channels? The single-channel TBR1025 packs the power of the full-size 4-channel unit in a small, economical package.*

### ORDERING INFORMATION

**TBR4100-416** Four-Channel Free Radical Analyzer with Lab-Trax 4/16 Data Acquisition System

*Includes TBR4100 analyzer & power cord, Lab-Trax-4/16 data acquisition system & USB cable, 4 BNC cables, 1 electrode adapter cable, 1 temperature probe, 2 sensors of your choice, and sensor start-up kit(s), if applicable.*

**TBR1025** Single-Channel Free Radical Analyzer

*Includes 1 sensor of your choice & 1 start-up kit*

### RECOMMENDED ACCESSORIES / REPLACEMENT PARTS

<b>SNAP50</b>	SNAP S-Nitroso-N-acetyl-D-penicillamine, 50 mg vial
<b>ISO-TEMP-2</b>	2 mm Platinum RTD Temperature Sensor (requires 91580)
<b>91580</b>	Microsensor Adapter Cable

# Optical Oxygen Sensor

## Oxygen detection in small samples

### Features

- Oxygen detection using the phase shift between a reference signal and measured signal
- Common sterilization procedure via gamma radiation, ethanol, para-acetic acid, ethylene oxide
- Easy cleaning procedure with 3% H<sub>2</sub>O<sub>2</sub>, ethanol or common soap solution
- Compatibility with ethanol, methanol and isopropanol
- No cross-sensitivity in static aqueous solutions with pH 1-14, CO<sub>2</sub>, H<sub>2</sub>S, SO<sub>2</sub> or any other ionic species
- Cross-sensitivity for organic solvents such as acetone, toluene, chloroform or methylene chloride and any chlorine gas
- Available with SMA or ST connector
- Small sensor tip diameter ≤ 50 μm

### Benefits

- Designed for use in gaseous and aqueous solutions
- Easy calibration via a two-point calibration with an oxygen-free environment (nitrogen, sodium sulfite) and an air-saturated environment
- Fast response time  $t_{90} \leq 3$  seconds (in gas phase)
- Retractable up to 10 mm
- No oxygen consumption during measurements

### Applications

- Oxygen measurement in small samples
- Bioprocess control

**BioOxy** is a new and innovative technology for measuring oxygen in gaseous and aqueous phase. **BioOxy** is an optical oxygen sensor with important advantages over using common Clark type electrodes.

- **BioOxy** does not consume oxygen while making measurements.
- Similar performance overall at a lower cost. This helps lower research or operational expenses compared to competitive consumable sensors.
- **BioOxy** may be gamma sterilized. Sensors can be irradiated, as well as ETOH and PAA, which may be ideal for biological applications.
- Very good resistance to long term photobleaching allows for extended use, a periodic re-calibration is recommended for optimal results.
- Sensor-to-sensor repeatability because of the very small unit-to-unit variation allows much more confident data from test to test using different sensors.

The **BioOxy** is ideally adapted for examination of small sample volumes thanks to its miniature sensor tip diameter of less than 50 μm. The **BioOxy** is well suited for long-term measurements and biotechnological applications. The signal does not depend on the flow rate of the sample and the measurement principle requires no contact.

The physical principle is based on quenching of luminescence, which is caused by the collision between molecular oxygen and a luminescent dye molecule in the excited state. There is a direct relationship between the measured phase angle and oxygen concentration. The phase angle decreases when oxygen is present in your sample. The transfer of the emitted light signal is done optically using high quality fused silica fibers, ensuring a constant performance from lot-to-lot.



### BIOOXY SPECIFICATIONS

OXYGEN MEASUREMENT	Gaseous & Dissolved O <sub>2</sub>
PROBE SIZE	240 μm with active area ~50 μm Retractable 10 mm
MEASUREMENT RANGE	0 – 100 % O <sub>2</sub> 0 – 1000 hPa
LIMIT OF DETECTION	0.05 % Oxygen
RESOLUTION	± 0.09 % O <sub>2</sub> or better
ACCURACY @ 20 °C	± 0.4 % O <sub>2</sub> @ 20.9 % O <sub>2</sub> ± 0.05 % O <sub>2</sub> @ 0.2 % O <sub>2</sub>
MEASUREMENT TEMPERATURE RANGE	0 to + 80 °C
RESPONSE TIME (T <sub>90</sub> )	< 3 sec. (in gas phase)

### ORDERING INFORMATION

- 02-T0050L250-12012** BioOxy Sensor, 50 μm Tip, ST Connector, pack of 2  
**02-T0050L250-12022** BioOxy Sensor, 50 μm Tip, SMA Connector, pack of 2

### COMPATIBLE WITH:

- OXY-MICRO-AOT** Fiber Optic Oxygen Meter with ST Connection  
**FLOX-PROBE** In Situ Oxygen Monitoring Kit with SMA Connection

### WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
 Brazil: (013) 406-29703 • info@brazil.wpiinc.com • www.wpiinc.com

Germany: +49 (0)6031 1602171 • wptide@wpi-europe.com • www.wpi-europe.com  
 China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net

# Oxygen Monitoring Kit

Detect O<sub>2</sub> for real-time analysis

## Features

- Fluorescence-sensing detector for optical sensors, a viable alternative to traditional chemical sensing devices
- Measure fluorescence lifetime, phase and intensity
- LED based photometry makes this unit affordable—Half the cost of comparable phase measurement systems
- Excellent stability, extremely low drift and phase noise
- Simple calibration, setup and control

## Benefits

- Self-contained, benchtop system—invariant to fiber bending or stray light
- Easy setup and control, combined with its stability and sensitivity make this unit perfect for long term studies
- Three coatings for probes and patches give you many options.

## Applications

- Oxygen sensing applications requiring stability and sensitivity to drift which must be undisturbed for lengths of time
- Monitor O<sub>2</sub> partial pressure in gas and aqueous solutions and in non-aqueous vapors and solutions
- Monitor traces of oxygen in gas and liquids

**FLOX** is a device for measuring fluorescence lifetime, phase and intensity. It uses LED excitation and photodiode detection with filter-based wavelength selection for easy experimental set-up and control.

### Self contained unit

The compact, self-contained unit makes it invariant to fiber bending and stray light. It also has a wide dynamic range of optical intensity, as well as low optical and electronic crosstalk, and low drift and phase noise.

**FLOX** is especially useful for oxygen sensing applications where stability and sensitivity to drift is important and where sample set-ups must be left undisturbed for long periods of time.

### Easy setup, great for long term studies

The new oxygen sensing system measures fluorescence lifetime, phase and intensity, using LED excitation and photodiode detection with filter-based wavelength selection. The system is simple to set-up and control. When stability and sensitivity to drift are important in your oxygen sensing experiment, this unit is ideal. It is perfect for applications where sample set-ups must be left undisturbed for long periods of time.

### Three coatings available for probes and patches

The optical sensors consist of transducer materials, applied to the tips of optical fibers or to substrates such as patches or cuvettes, which change optical properties in response to specific analytes in their immediate environment.

**OXY**—The standard oxygen sensor designed for monitoring oxygen partial pressure in gas and aqueous solutions is a fiber optic fluorescence probe with a proprietary oxygen sensing coated tip.

**HIOXY**—Designed for monitoring oxygen partial pressure in non-aqueous vapors and solutions. The sensor coating chemistry is compatible with oils, alcohols, and hydrocarbon-based vapors and liquids.

**FOSPOR**—A new generation of highly sensitive sensor coating can be used for monitoring traces of oxygen in gas and liquids.

The oxygen sensor probes are low-power and offer high sensitivity, reversibility and stability, ideal for remote monitoring. The thin coating on the probe tips consumes no oxygen, allowing for continuous contact with the sample. They are ideal for viscous samples and are immune to interference caused by pH, ionic strength or salinity fluctuations or biofouling.



FLOX-PROBE



### OR125 g and OR125 gT

These probes are 1/8" OD optical fiber probes used as direct replacement for 1/8" OD O<sub>2</sub> electrodes. • 1000 µm optical fiber, stainless steel ferrule; 3.175 mm OD, 63.5 mm length

### R Sensor Probe

1000 µm fiber in a stainless steel 1/16" ferrule; 1.587 mm OD, 152.4 mm length

### AL300

A 500 µm OD (300 µm core diameter) aluminum-jacketed optical fiber probe for applications that require fine spatial resolution. • 300 µm aluminum-jacketed fiber assembly; 500 µm OD, 1 m length

### OR125

The OR125 is a 1/8" OD optical fiber probe used as a direct replacement for 1/8" OD O<sub>2</sub> electrodes. • 100 µm optical fiber, stainless steel ferrule; 3.175 mm OD, 63.5 mm length

### P1600


The P1600 is a silicone-jacketed, polyimide-coated optical fiber probe used in environments where a non-metallic probe is required. • 200 or 600 µm optical fiber with silicone jacketing; 710 µm OD, 2 m length

## ORDERING INFORMATION

<b>FLOX-PATCH</b>	Non-Invasive Oxygen Monitoring Kit, including phase measurement system, temperature probe; select sensor patches when ordering
<b>FLOX-PROBE</b>	In Situ Oxygen Monitoring Kit, including phase measurement system, temperature probe; select sensor when ordering

# Carbon Fiber Microelectrodes

## Electrochemical detection of oxidizable compounds



### Features

- Sensitive, renewable/durable and economical carbon fiber electrodes for electrochemical detection of oxidizable compounds
- Excellent linearity to the oxidizable compounds

### Benefits

- Precision tip size and length of the CF electrodes
- Renewable

### Applications

- Detection of neurotransmitters and oxidizable compounds *in vivo* or *in vitro*

Carbon fiber microelectrodes have been used in both the detection of oxidizable compounds (Gonon, *et al.*, 1978; Cahill and Wightman, 1995) and extracellular single-unit recording (Armstrong-James and Millar, 1979). WPI's ultra-sensitive and low-noise carbon fiber (CF) electrodes can be applied, with our Micro-C Potentiostat or similar instruments, in the electrochemical detection of catecholamines (epinephrine, norepinephrine and dopamine), indolamines (serotonin, 5-HT or melatonin), ascorbic acid, Fe (II) and other oxidizable compounds.

CF electrodes (diameter of 10 or 30 μm) respond with an excellent linearity to the oxidizable compounds (see figure below) and can detect the compounds as low as 0.2 nM. While the shorter (25-100 μm) CF electrodes are suitable for *in vivo* amperometric and voltammetric measurements, the longer CF electrodes provide higher sensitivity and are especially useful for the *in vitro* studies (amperometric or differential pulse voltammetry). When used with the Micro-C Potentiostat, these CF electrodes can be activated and renewed

Dopamine Concentration/Response Curve




Fig. 1 — Excellent linearity in the response of carbon fiber electrode (CF30-500) to dopamine recorded on Micro-C. Courtesy: Drs. D. Yeomans and X.-T. Wang, University of Illinois at Chicago.

in sensitivity for multiple use. The selective detection of catecholamines can be achieved with our Nafion-coated CF electrodes. For selective detection of 5-HT and ascorbic acid, please contact WPI for more information.

### References

- P. S. Cahill, R. M. Wightman, *Anal. Chem.*, **67**, 2599-2605 (1995).  
 F. Gonon, *et al.*, *Hebd Seances Acad. Sci. Ser.* **286**, 1203 (1978).  
 M. Armstrong-James, J. Millar, *J. Neurosci. Methods*, **1**, 279 (1979).




Fig. 2 — Extracellular recording using a carbon electrode in CA1 region of the hippocampus in an anesthetized rat shows ultra-low noise (<5 μV). Courtesy: Dr. Carolyn Harley of Memorial University, Newfoundland, Canada.

### ORDERING INFORMATION

#### CARBON FIBER MICROELECTRODES, UNCOATED

	Diameter	Length	(pack of 5)
CF10-100	10 μm	100 μm	
CF10-250	10 μm	250 μm	
CF10-500 *	10 μm	500 μm	
CF30-50 *	30 μm	50 μm	
CF30-100	30 μm	100 μm	
CF30-500 *	30 μm	500 μm	
CF30-1000 *	30 μm	1000 μm	

#### CARBON FIBER MICROELECTRODES, NAFION-COATED

	Diameter	Length	(pack of 5)
CFN10-50 *	10 μm	50 μm	
CFN10-100 *	10 μm	100 μm	
CFN10-250 *	10 μm	250 μm	
CFN30-50 *	30 μm	50 μm	
CFN30-100 *	30 μm	100 μm	
CFN30-250 *	30 μm	250 μm	
CFN30-1000 *	30 μm	1000 μm	

\* Built to order — allow up to 4 weeks manufacturing time.

# Fiber Optic Oxygen Meter and Sensors

*A generation of sensors based on luminescence lifetime*

## Features

- Oxygen is not consumed during the experiment
- Immune to electrical and magnetic interference
- Excellent long-term stability
- No lengthy polarization necessary (like Clark-type O<sub>2</sub> electrodes)
- Fast response time < 0.5 s for MicroTip sensors
- Probe size of MicroTip sensors as small as 50 μm
- Measurement is feasible in dry gas
- Optical isolation of sensor tip available for fluorescent or photosynthetically active samples



OXY-MICRO-AOT

## Benefits

- Unaffected by light source stability and intensity fluctuations, because it is based on luminescence lifetime detections.
- Compact, portable meters may be used inside or outside

## Applications

- Process control like bottling plant in breweries and quality control of packages (OxyMini)
- Biotechnology like control of cell culture media and non-invasive control of bioreactors (OxyMini)
- Implantation of oxygen sensors into soil and trees (OxyMini)
- Oxygen profiles of marine sediment, soils, or tissue (OxyMicro)
- Implantation in living tissue like heart or muscle tissue (OxyMicro)
- Control of cell culture media in Biotechnology (OxyMicro)

The measurement principle of the sensor system is based on the detection of oxygen concentration as a function of luminescence lifetime, either in dissolved or gaseous phase environments.

### Unaffected by light source stability/intensity

Conventional fiberoptic oxygen sensor systems based on intensity measurements are limited in their accuracy by light source stability and ambient light fluctuations. Using a luminescence lifetime detection, measurements are not affected by light source stability, intensity fluctuations caused by fiber bending or changes of the optical properties of the sample (turbidity, refractive index, coloration, etc.).

Calibration: The sensors can be calibrated by a simple two point calibration, 100% air-saturation and 0% air saturation.


### Compact, portable meters, inside and outside

The **OxyMini** and **OxyMicro** fiber optic oxygen meters are compact and easy to transport. Designed for in/outdoor use, they can be connected to a PC via a RS232 interface. Data can be visualized, analyzed and stored with supplied software. A full range of sensors are available.


## OxyMini systems

The OxyMini is a one channel fiber optic O<sub>2</sub> meter for fiber optic O<sub>2</sub> minisensors. These sensors are based on 2 mm polymer optical fibers and have a length of 2.5 m.

**MiniTip** – This dipping probe (501641) has a tip diameter of 4 mm and consists of a polymer optical fiber with an O<sub>2</sub> sensitive coating. Its range is 0–100%. It has a response time (t<sub>90</sub>) of ~40 s.



**MiniFlow** – The MiniFlow probe (501642) is a miniaturized fiber optic chemical sensor integrated in a standard T-shape flow through cell which can be easily connected via Luer-Lock adapters to external tubings. Liquids can be pumped through the cell. It has a response time (t<sub>90</sub>) of approximately 40 s and an excellent long-term stability.



**MiniFoil** – Sensor material on a 1 cm<sup>2</sup> support disk made of polyester can be glued inside glass vials. Measure oxygen concentration non-invasively and non-destructively from



outside through the wall of the vial. Illuminate the sensor foil with a plastic fiber optic cable (501644, 01645). The wall of the flask must be transparent/non-fluorescent (response time (t<sub>90</sub>) of ~50s). The material can be implanted into animal tissues or custom-made housings.


## OxyMicro systems

**OxyMicro** is a one channel O<sub>2</sub> meter for fiber optic O<sub>2</sub> microsensors.


**MicroTip** – The MicroTip (501656) is a needle-type (27 ga.) oxygen micro sensor designed for applications where a small tip size (50 μm) and fast response time (t<sub>90</sub>) of 1 s are required. The oxygen sensitive sensor tip consists of 140 μm fiber tapered to a 50 μm tip. Housed inside a stainless steel needle 22 mm long and 0.4 mm diameter, it can penetrate a septum rubber or similar material. They are ideal for oxygen profiling in sediments and biofilms.



**MicroFlow** – The MicroFlow oxygen sensor (501657) is a miniaturized fiber optic chemical sensor optimized for fast response time (t<sub>90</sub> < 1 sec in gases, < 5 sec in liquids). The tiny probe has a tip size of 50 μm and is integrated in a T-shape flow cell for easy connection via Luer-Lock adapters to external tubings. Liquids (like water, blood, etc.) can be pumped through the cell.



**Microlmplant** – The Microlmplant oxygen sensor (501658) is an implantable probe (IMP) with a probe tip size 50 μm, an exposed fiber length of 5 mm and a jacket diameter of 900 μm. It has been successfully implanted in crabs, fishes and soil.



## SPECIFICATIONS

	MiniTip	MiniFlow	MiniSpot	MicroTip	MicroFlow	Microlmplant
MEASURE RANGE DISSOLVED/ GASEOUS	0-45 ppm, 0-100% 0-760 mmHg	0-45 ppm, 0-100% 0-760 mmHg	0-45 ppm, 0-100% 0-760 mmHg	0-45 ppm, 0-100% 0-760 mmHg	0-45 ppm, 0-100% 0-760 mmHg	0-45 ppm, 0-100% 0-760 mmHg
RESPONSE TIME [T <sub>90</sub> ] DISSOLVED/ GASEOUS	40 s 10 s	40 s 10 s	40 s 10 s	< 2 s < 0.5 s	< 2 s < 0.5 s	< 2 s < 0.5 s
STERILIZATION ETOH, H <sub>2</sub> O <sub>2</sub> AUTOCLAVABLE*	Y N	Y Y	Y Y	Y N	Y Y	Y Y
DRIFT**	< 0.1%	< 0.1%	< 0.1%	< 0.3%	< 0.3%	< 0.3%
ACCURACY**	0.2%					
RESOLUTION***	2.75 ±0.01 ppm, 9.00 ±0.05 ppm, 220 ±0.15 ppm, 45.0 ±0.25 mmHg, 150 ±0.75 mmHg, 375 ±2.6 mmHg					
TEMP RANGE	-10°C to 50°C					
PROBE ASSY LENGTH	2.5 m					

\*130°C, 1.5 atm \*\*100,000 data points, 20°C \*\*\*20°C

## ORDERING INFORMATION

### MINISENSOR SYSTEM

**OXY-MINI-AOT** Fiber-optic Oxygen Meter for Minisensors \*

### MINISENSORS (not interchangeable with Microsensors)

- 501641** MiniTip, fiber-optic oxygen sensor
- 501642** MiniFlow, flow-through cell with integrated planar oxygen sensor
- 503090** MiniSpot, planar oxygen-sensitive spot, 5 mm diam. (includes 10) Requires **501644**
- 501644** Polymer optical fiber with 1 SMA connector

### MICROSENSOR SYSTEM

**OXY-MICRO-AOT** Fiber-optic Oxygen Meter for Microsensors \*

### MICROSENSORS (not interchangeable with Microsensors)

- 501656** MicroTip, needle-type housing fiber-optic oxygen sensor, 50 μm tip
- 501656-C** MicroTip, needle-type housing, 50 μm tip, optical isolation
- 501656-F** MicroTip, needle-type housing, 140 μm flat tip
- 501657** MicroFlow, flow-through housed oxygen microsensor
- 501658** Microlmplant, implantable oxygen microsensor, 50 μm tip
- 501658-F** Microlmplant, 140 μm flat tip

\*Meter contains two analog outputs and one trigger input

# 4-Channel Data Acquisition System with Software

Low noise, high resolution system with 8 analog input and 3 analog output channels

## Features

- Powerful low-noise (<1 mV RMS) and high-resolution (16 bits) data acquisition system for sampling up-to 8 analog input channels and 3 analog output channels simultaneously, using standard BNC connections
- MDAC software provides easy to use interface controlling, with extensible standard and customized Data Processing and Analysis Tools



## Benefits

- Online Channel Math operations, general purpose Fast Fourier analysis (FFT) and digital filtering of Analog In channels
- Numerous basic signal forms can be combined to design experimental protocols, for most physiological applications
  - Factory designed standard or customized protocols
  - Semi-automated data analysis toolbox
- Protocol repeat function to avoid time consuming protocol programming of extended experiments

## Applications

- Muscle physiology (Can be used with **SI-MTM** Muscle Testing Platform, **SI-CTS200** Cell Tester System, **SI-HTB2** Horizontal Tissue Bath and **SI-BF-100** Biofluorometer)
- Stand alone general data recorder for Spectroscopy, Neuroscience and Electrophysiology (Can be used with **TBR4100** Free Radical Analyzer, Extracellular Bioamplifiers like **SYS-DAM50**, **SYS-DAM80**, **SYS-900A**, **ISO-80**, **EVOM2™** Volt Ohm Meter, **ATC2000** Animal Temperature Controller, **BP-1** Blood Pressure Monitor or the **BAT-12** Microprobe Thermometer)
- Instrument control for software triggered devices like **A365/A385/A395** Constant Current Stimulators, **MPS-2** Perfusion System, **SYS-PV820/SYS-PV830** Pneumatic PicoPumps, **Duo 773** Intracellular amplifiers and the **SYS-TBM4M** Transbridge Transducer Amplifier (e.g. for FORT force transducers)

Knowledge of the physiological characteristics of muscle tissue can be useful to quantify beneficial or adverse effects of drug supply on muscle function in pre-clinical and toxicological studies, evaluating muscle dystrophies, training effects in sports and rehabilitation (disuse vs. overuse) and advanced physiology and biomedical research.

This is usually achieved by quantifying the contractile and/or the elastic properties of muscle tissue. This needs the programming of different and specific experimental protocols (isometric, concentric and eccentric, isokinetic or isotonic), so that the physiological structure of interest can be quantified. **LabTrax-MDAC** data acquisition software was designed for use with WPI's Muscle Physiology line to test physiological characteristics of muscle tissues in various conditions, using factory designed standard or customized protocols. The semi-automated Data Analysis Toolbox of standard protocols gives quick access to user-friendly, readable and interpretable results of the experiments.

### Variety of muscle physiology applications

The physiological response of muscle tissue to training, disuse, nutrition, drug supply and others factors may be studied by adding accessories to the system, like:

- Study of the muscle's force production capacities in combination with the  $Ca^{2+}$  release from the sarcoplasmic reticulum (SR) and ATPase consumption. The perfect instrument for this is WPI's Biofluorometer (**SI-BF-100**) in combination with any system of WPI's Muscle Physiology line, controlled via **LABTRAX-MDAC**.

- Study of the muscle's force production capacities from direct muscle or peripheral nerve stimulation. For this experiment, use WPI's programmable isolated current stimulators (**A365**, **A380** or **A395**), controlled via **LABTRAX-MDAC**.

**LABTRAX-MDAC** provides easily used continuous stimulation protocols, so that especially cardiac cells/tissue remain intact during experimental resting periods.

**LABTRAX-MDAC** is also well suited for other software triggered instruments or as a stand-alone general data recorder for selected WPI Instruments.



The back panel of the Lab-Trax-8/16 has four analog outputs, digital inputs or outputs, a USB port, power socket and power switch.

## LAB-TRAX-8/16 SPECIFICATIONS

ANALOG INPUTS	8 BNC connections
INPUT RANGE	± 10V
SYSTEM NOISE	< 1 mV RMS
ISOLATION	1,500V
OPERATING CURRENT	800 mA maximum
ANALOG OUTPUTS	4 BNC connections
OUTPUT RANGE	± 10V
IMPLEMENTED FILTER	5 <sup>th</sup> order low-pass Bessel filter with 3dB cut-off frequency
OUTPUT IMPEDANCE	100Ω
OUTPUT CURRENT	15 mA
DIGITAL I/O	16/16 TTL (BNC or DB-9 Connector)
LOGIC HIGH VOLTAGE	3.3 V minimum
LOGIC LOW VOLTAGE	1.0 V maximum
ANALOG & DIGITAL INPUTS	Operating voltage protected to ±30V
PC INTERFACE	USB 2.0
RESOLUTION	16 bits
POWER SOURCE	12V DC

## ORDERING INFORMATION

<b>LABTRAX-MDAC</b>	Lab-Trax-8/16 with MDAC software
<b>LAB-TRAX-4</b>	4-Channel General Data Acquisition System
<b>2851</b>	BNC to BNC Cable

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
Brazil: (013) 406-29703 • info@brazil.wpiinc.com • www.wpiinc.com

Germany: +49 (0)6031 1602171 • wvide@wpi-europe.com • www.wpi-europe.com  
China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net

# Novel Fiber Optic pH System

Referenced measurements with single excitation

## Features

### pHOptica Meter

- Single-channel, compact, easy to transport fiber-optic meter for pH measurements with miniature sensors
- Two 12-bit, programmable analog outputs, with electrical isolation
- One external trigger input, with electrical isolation
- Computer with RS232 interface required for operation

### pH Optical Sensors

- Immune to electrical interferences and magnetic fields
- Low drift
- High spatial resolution due to small tip size
- Measurement in very small sample volumes
- Additional optical isolation of the sensor tip is available for measurements in colored or photosynthetically active samples

## Benefits

- User-friendly software saves and visualizes measured values
- Several pHOptica meters can be connected to one computer
- Temperature variation is recorded using a temperature sensor
- No reference electrode is needed


## Applications

### pHOptica Micro System

The pH Optica micro system is a single channel pH system for use with fiber optic micro sensors. The applications include:

- Penetration or implantation into living tissue (heart, muscle or animal blood vessels).
- Soil implantation for pH measurement.
- Implantation into customer-made housing.

TIP: To protect the small glass fiber tip against breaking, suitable housings and tubings around it, depending on the respective application, were designed.



### pHOptica Mini System

The pHOptica mini system is a one channel pH system for use with fiber optic mini sensors, foil and spot surface sensors for applications like:

- Non-invasive and non-distractive pH measurements from outside through flask walls (cell culture)
- Online pH monitoring by flow through cells
- Dipping probe pH measurements

### pH Micro Sensors


- Tip size 140 micrometer
- Drift of 0.1 pH units for 2000 measurements (16 hours measurement in the 30 sec data update mode)



**Implantable sensor**—without any housings implantation into animal blood circuits; soil implantation; implantation in custom-made housings



**Needle-Type Housing Sensor**—the glass-fiber with its pH-sensitive tip is protected inside a stainless steel needle (18 ga.); fiber has to be extended during measurement; penetration through septum.



PH-OPTICA-MINI

### pH Mini Sensors

- OD of the dipping sensor is 4 mm
- Sterilization of the pH sensor spots via gamma radiation
- pH mini sensor meter is based on 2 mm PMMA waveguides
- Drift of 0.1 pH units for 10,000 measurements (4 days measurement in the 30 sec data update mode).

pHOptica™ is a pH measuring system which uses fiber optic sensors and patented DLR technology. This method allows referenced measurements with single excitation to be implemented.

Two different housings and sensor spots (sensorfoils) are available.



**POF Coated with a pH-Sensitive Foil**—Small and robust pH dipping sensor; no reference electrode needed.



**Flow-Through Cell with Integrated pH Sensor**—On-line monitoring; can be easily connected via Luer-Lock adapters.

## SPECIFICATIONS

DATA INTERFACE	RS232
SAMPLE RATE	1 sample per sec
MEASURING PH RANGE	5 - 9
RESOLUTION (20°C)	± 0.03 (microsensors) ± 0.01 (minisensors)
RESPONSE TIME	<1 min.
DIMENSIONS	185 x 110 x 45 mm
WEIGHT	630 g
POWER SUPPLY	100-220 VAC

## ORDERING INFORMATION

### MINISENSOR SYSTEM (cannot be used with microsensors)

<b>PH-OPTICA-MINI</b>	Fiber Optic pH Meter for Minisensors, Foils and Spots
<b>503538</b>	pH MiniTip, fiber optic pH sensor dipping probe, disposable (4 mm OD), pkg. of 3. <i>Requires cable 503110.</i>
<b>501644</b>	Polymer Optical Fiber with 1 SMA connector
<b>503110</b>	Fiber Optic Cable with 1 SMA connector

### MICROSENSOR SYSTEM (cannot be used with minisensors)

<b>PH-OPTICA-MICRO</b>	Fiber Optic pH Meter for Microsensors
<b>502124</b>	pH Microimplant, fiber optic pH implantable sensor (140 µm OD), pkg. of 3

# Reference Electrodes

Low electrolyte leakage, stable potential with low resistance

## Features

- Extremely low electrolyte leakage
- Stable, reproducible potential with low resistance
- Chemically resistant to strong acids and bases

## Benefits

- May be used with ion selective electrodes without contamination from the reference electrode

## Applications

- Small volume, low salt concentration measurement (SDR)

**Dri-Ref™** reference electrodes were developed by WPI to have extremely low electrolyte leakage properties, hence the name "Dri-Ref." In addition to this key feature, these electrodes exhibit stable and reproducible potential and low resistance. Stored in KCl when not in use, they have a long life expectancy.

### May be used with ion selective electrodes


Although the internal filling solution contains KCl, the low fluid leakage means Dri-Ref may be used in combination with ion selective electrodes, including those for  $K^+$  and  $Cl^-$ , without significant contamination from the reference electrode.

The Dri-Ref electrodes are chemically resistant to strong acids and alkalines. Dri-Ref electrodes are not suitable for use in organic solvents. In addition, the long, thin FLEXREF may be easily manipulated to accommodate a difficult experimental setup.

**SUPER-Dri-Ref** – With a diameter of 2mm, SUPER-Dri-Ref *does not leak electrolyte at all*. Exhibiting the electrical stability of a classic flowing junction reference cell, this electrode exhibits low resistance and a stable half-cell potential essentially independent of sample electrolyte concentration. SUPER-Dri-Ref is ideal for small volume and low salt concentration measurements.

**Micro-Reference Electrode** – Only 450  $\mu m$  in diameter and one inch long, the new **DRIREF-450** reference electrode can be used along with other sensors in space-restricted areas and very small sample volumes.

**Luer-Tip Reference** – The male luer fitting at the front of the **DRIREF-L** allows it to be easily connected to a female luer port (see WPI's luer fittings kit, page 219 to form a tight seal — a very convenient installation for a flow-through system.



## DRI-REF ELECTRODE SPECIFICATIONS

	DRIREF-450	DRIREF-5	DRIREF-2	FLEXREF	SDR	DRIREF-L	DRIREF-5SH	DRIREF-2SH
LENGTH	2.54 cm	9 cm	13 cm	13 cm	9 cm	7.5 cm	3.5 cm	2.54 cm
DIAMETER	450 $\mu m$	4.7 mm	2 mm	1.5 mm	2 mm	Standard Luer	4.7 mm	2 mm
CONSTRUCTION	Coated Glass	Epoxy	Isoplast™	Teflon™	PVC	Polypropylene	Epoxy	Isoplast™
RESISTANCE (TYPICAL)	<5 K $\Omega$	-500 $\Omega$	-2.7 K $\Omega$	-2.7 K $\Omega$	<5 K $\Omega$	-500 $\Omega$	-500 $\Omega$	-2.7 K $\Omega$
ELECTROLYTE LEAKAGE (ML/HR)	—	-7.4x10 <sup>-7</sup>	-5.7x10 <sup>-8</sup>	-5.7x10 <sup>-8</sup>	—	7.4x10 <sup>-7</sup>	-7.4x10 <sup>-7</sup>	-5.7x10 <sup>-8</sup>
LEAD LENGTH	30 in. (76 cm)							
CONNECTOR	2 mm pin							
Filling Solution	KCl							

## ORDERING INFORMATION

<b>FLEXREF</b>	Flexible Dri-Ref, 1.5 mm diam.
<b>DRIREF-2</b>	Dri-Ref, 2 mm diam.
<b>DRIREF-2SH</b>	Dri-Ref, 2 mm diam. (Short)
<b>DRIREF-5</b>	Dri-Ref, 4.7 mm diam.
<b>DRIREF-5SH</b>	Dri-Ref, 4.7 mm diam. (Short)
<b>SDR2</b>	SUPER-Dri-Ref, 2 mm diam.
<b>DRIREF-450</b>	Micro-Dri-Ref, 450 $\mu m$ diam.
<b>DRIREF-L</b>	Reference Electrode with Luer Tip

# Calcium Calibration Solutions

## CALBUF-1

### For use with calcium electrodes

A set of eight calcium buffers covering the range of concentration from  $10^{-1}$  to  $10^{-8}$  M  $Ca^{++}$ . Each buffer contains 20 mL of solution and enough potassium chloride to set the ionic strength to 0.1 M. Limited shelf life; use within 30 days.

Concentration:  $1 \times 10^{-1}$ ,  $1 \times 10^{-2}$ ,  $1 \times 10^{-3}$ ,  $1 \times 10^{-4}$ ,  $1 \times 10^{-5}$ ,  $1 \times 10^{-6}$ ,  $1 \times 10^{-7}$ ,  $1 \times 10^{-8}$  M at 20°C. Limited shelf life; use within 30 days.

## ORDERING INFORMATION

**CALBUF-1** Kit of 8 Calcium Buffer Solutions

## CALBUF-2

### Use with calcium fluorescent indicators

CALBUF-2 is especially suitable for calibrating fluorescent  $Ca^{++}$  indicators. It provides eleven buffer standards in the  $10^{-4}$  to  $10^{-8}$  M  $Ca^{++}$  range, whereas other commonly used fluorescent  $Ca^{++}$  indicators have the apparent  $K_d$  in the range of 100 to 300 nM. As with any ionic sensitive indicator, the sensitivity range of these indicators is about 1.0 log unit above and below the  $K_d$ . CALBUF-2 provides seven calibration points in this sensitivity range. It has an osmolarity of 0.305, which is isotonic with most mammalian cells.

Concentration:  $1 \times 10^{-8}$ ,  $4 \times 10^{-8}$ ,  $1 \times 10^{-7}$ ,  $2.5 \times 10^{-7}$ ,  $5 \times 10^{-7}$ ,  $7.5 \times 10^{-7}$ ,  $1 \times 10^{-6}$ ,  $4 \times 10^{-6}$ ,  $1 \times 10^{-5}$ ,  $4 \times 10^{-5}$ , and  $1 \times 10^{-4}$  M at 20°C. Ionic strength: 0.150 M. 11 bottles, 20 mL each. Limited shelf life; use within 30 days.

## ORDERING INFORMATION

**CALBUF-2** Kit of 11 Calcium Buffer Solutions

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
Brazil: (013) 406-29703 • info@brazil.wpiinc.com • www.wpiinc.com

Germany: +49 (0)6031 1602171 • wptide@wpi-europe.com • www.wpi-europe.com  
China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net



# Ion Selective Electrodes

Accurately measure calcium, potassium, hydrogen or TPP

## Features

- Fast, accurate, economical
- Superior, stable PVC membrane
- Fast response
- 2mm diameter tips
- Interchangeable tip holder
- Each kit includes 3 electrode tips and MicroFil filling syringe

## Benefits

- Inexpensive
- Use to measure various ions in biological media

## Applications

- Detection of ions *in vivo* or *in vitro* for biological applications

These highly stable electrodes accurately measure calcium, potassium, hydrogen and TPP (Tetraphenylphosphonium) ion activity. Tips consist of 2 mm diameter plastic tubes sealed at one end with an ion-sensitive membrane. After filling with electrolyte solution, you can insert the tube into the holder and connect it to a pH meter. Tips and holders are interchangeable, so one tip may be replaced with another sensitive to a different ion. Replacing a tip takes less than a minute. Electrode tips normally last several months, when stored properly in saline solution. When replacement is necessary, only the tip needs to be replaced.

Kwik-Tip electrodes are available separately and as kits. Each "KWIK" Electrode Holder kit includes a reusable holder and three removable tips. In addition to a 4-foot BNC cable and an electrolyte filling syringe; "TIP" Electrode Kits contain three electrode tips for a specific ion. **A separate reference electrode, such as WPI's Dri-Ref™, is also required.**



## ORDERING INFORMATION

### COMPLETE KITS

<b>KWIKCAL-2</b>	Holder & 3 Calcium Electrodes
<b>KWIKH-2</b>	Holder & 3 Hydrogen Electrodes
<b>KWIKPOT-2</b>	Holder & 3 Potassium Electrodes
<b>KWIKTPP-2</b>	Holder & 3 TPP (Tetraphenylphosphonium) Electrodes

### HOLDERS AND REPLACEMENT TIPS

<b>KWIK-2</b>	Electrode Holder with BNC cable
<b>TIPCA</b>	Calcium Electrode Tips (3)
<b>TIPH</b>	Hydrogen Electrode Tips (3)
<b>TIPK</b>	Potassium Electrode Tips (3)
<b>TIPTPP</b>	TPP+ (Tetraphenylphosphonium) Electrode Tips (3)

## KWIK-TIP ELECTRODE SPECIFICATIONS

Part #	Electrode	Color Code	Recommended Filling Solution	Min. Slope/Decade	Concentration Range	Selectivity Coefficients (-log)
<b>TIPCA</b>	Calcium	Green	0.1 M CaCl <sub>2</sub>	28 mV	0.1 M - 10 <sup>-6.75</sup> M	Na <sup>+</sup> 5.5, K <sup>+</sup> 5.4, Mg <sup>++</sup> 4.9
<b>TIPH</b>	Hydrogen	Red	1 M Citric Acid, 0.01 M NaCl, pH 5.6	54 mV	pH 5.0 - 12	Na <sup>+</sup> 10.4, K <sup>+</sup> 9.8, Ca <sup>++</sup> 11.1
<b>TIPK</b>	Potassium	Yellow	0.1 M KCl	54 mV	0.1 M - 10 <sup>-4.5</sup> M	Na <sup>+</sup> 4.0, Ca <sup>++</sup> 3.9, Mg <sup>++</sup> 3.0
<b>TIPTPP</b>	TPP <sup>+</sup>	Purple	10 mM TPP <sup>+</sup>	54 mV	0.001 M - 10 <sup>-4</sup> M	K <sup>+</sup> 6.0

# Liquid Ion Exchangers

Make micropipettes to record cellular concentrations

WPI's Liquid Ion Exchangers (LIX), for use with the **FD223A** Electrometer, allow intracellular measurements to be made for cations (hydrogen, potassium and calcium).

## SPECIFICATIONS

ION	H <sup>+</sup>	K <sup>+</sup>	Ca <sup>++</sup>
CATALOG NO.	<b>IE 010</b>	<b>IE 190</b>	<b>IE 200</b>
<b>SELECTIVITY COEFFICIENTS*</b>			
Na <sup>+</sup>	12.7	1.97	5.5
Mg <sup>++</sup>	—	2.95	4.9
K <sup>+</sup>	—	—	5.4
Ca <sup>++</sup>	—	2.7	—
USEFUL pH RANGE	2-10	4-10	4-10
SLOPE	56 mV	58 mV	28 mV
LINEAR RANGE	pH 4-12	pK 0-3	pCa 1-7
APPROX. EQUIV.	—	Corning 477317	ETH1001

\*Selectivity Coefficients are expressed here as -log K<sub>ij</sub> or pK<sub>ij</sub>.



When used in micropipettes to record cellular ion concentrations, consider using WPI's Duo 773 electrometer (channel A).

## ORDERING INFORMATION

<b>IE010</b>	Hydrogen Ion Exchanger (0.1 mL)
<b>IE190</b>	Potassium Ion Exchanger (1.0 mL)
<b>IE200</b>	Calcium Neutral Ion Exchanger (0.1 mL)

# Multi-Port Measurement Chamber

4-port closed chamber for measurements of NO, O<sub>2</sub>, H<sub>2</sub>O<sub>2</sub> & other species in cell culture, temperature stabilized

## Features

- Four port (**NOCHM-4**) chamber accommodates WPI's 2 mm sensors for nitric oxide (**ISO-NOP**), oxygen (**ISO-OXY-2**), hydrogen peroxide (**ISO-HPO-2**) and WPI's KWIK-TIP ion selective electrodes in combination with WPI's 2 mm Dri-Ref™ reference electrodes
- Two additional top ports for injection of reagents using WPI's MicroFil™ syringe needles
- Temperature control through an external circulating bath
- The chamber can be used for nitric oxide and other species calibration at temperatures from 4-40 °C

## Benefits

- Closed chamber design greatly reduces the surface area of the solution exposed to air
- One top port and up to three side ports configuration provides adequate space for convenient sample and electrode manipulation

## Applications

- Simultaneously measurement of free radicals such as NO, H<sub>2</sub>O<sub>2</sub>, H<sub>2</sub>S, O<sub>2</sub> and other ions at controlled conditions for cultured cell, cell suspensions or biological media

## References

**Olson, K. R., et. al.** (2017). Catalase as a sulfide-sulfur oxido-reductase: An ancient (and modern?) regulator of reactive sulfur species (RSS). *Redox Biology*, 12, 325–339. <http://doi.org/10.1016/j.redox.2017.02.021>

**Zhou, D., et. al.** (2017). Oxygen binding and nitric oxide dioxygenase activity of cytoglobin are altered to different extents by cysteine modification. *FEBS Open Bio*, 7(6), 845–853. <http://doi.org/10.1002/2211-5463.12230>

**Liu, X., et. al.** (2017). Cytoglobin regulates blood pressure and vascular tone through nitric oxide metabolism in the vascular wall. *Nature Communications*, 8, 14807. <http://doi.org/10.1038/ncomms14807>



NOCHM

**Santos, S. S., et. al.** (2017). NO production and potassium channels activation induced by *Crotalus durissus cascavella* underlie mesenteric artery relaxation. *Toxicol*, 133, 10–17. <http://doi.org/10.1016/j.toxicol.2017.04.010>

**DeLeon, E. R., et. al.** (2016). A case of mistaken identity: are reactive oxygen species actually reactive sulfide species? *American Journal of Physiology-Regulatory, Integrative and Comparative Physiology*, 310(7), R549–R560. <http://doi.org/10.1152/ajpregu.00455.2015>

## ORDERING INFORMATION

<b>NOCHM-4</b>	Four-Port Closed Chamber, for use with WPI's 2.0 mm electrodes (e.g., <b>ISO-NOP</b> and <b>ISO-OXY-2</b> , etc.)
<b>NOCHM-P</b>	Spare Plug-adapter for ISO-NOP nitric oxide electrode
<b>800100-5</b>	Spare Center Chamber Gasket (package of 5)

# NSA Pre-polarizer

Keep extra NO sensors ready to use

Achieve a stable background current quickly. This small battery-powered device applies a potential to the NO electrode equivalent to the potential applied by the ISO-NO meter. A sensor, which has been connected to the activator, may be transferred to the meter for immediate use. **For use with all WPI NO electrodes.**



## ORDERING INFORMATION

<b>NSA-3</b>	ISO-NO Activator
--------------	------------------

# SNAP

S-Nitroso-N-acetyl-D-penicillamine

SNAP is a stable green crystalline S-nitrosothiol compound that mimics the action of nitric oxide *in vivo*. It has vasodilatory properties and has been shown to relax isolated bovine coronary artery rings by activating soluble granulate cyclase. This reagent also actuates apoptosis in mouse thymocytes and has been accounted for reversible inactivation of protein Kinase C. SNAP can be used for calibration of all WPI NO sensors.

M.W. 220.2 • Purity > 98% by NMR or TLC

## ORDERING INFORMATION

<b>SNAP25</b>	SNAP, 25 mg vial
<b>SNAP50</b>	SNAP, 50 mg vial
<b>SNAP100</b>	SNAP, 100 mg vial

# GSNO

S-nitrosoglutathione

GSNO has been identified *in vivo* as a potential storage and transport vehicle for NO in the body. GSNO has been used in clinical trials to treat a form of preeclampsia and to prevent platelet aggregation. It also has considerable potential as an NO donor in medicine.

M.W.336.3 • C<sub>10</sub>H<sub>16</sub>N<sub>4</sub>O<sub>7</sub>S • Purity > 98% • Soluble in water or DMSO • Storage: -20°C

## ORDERING INFORMATION


<b>GSNO-50</b>	GSNO 50 mg vial
<b>GSNO-100</b>	GSNO 100 mg vial

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • [wpiuk@wpi-europe.com](mailto:wpiuk@wpi-europe.com) • [www.wpi-europe.com](http://www.wpi-europe.com)  
Brazil: (013) 406-29703 • [info@brazil.wpiinc.com](mailto:info@brazil.wpiinc.com) • [www.wpiinc.com](http://www.wpiinc.com)

Germany: +49 (0)6031 1602171 • [wptide@wpi-europe.com](mailto:wptide@wpi-europe.com) • [www.wpi-europe.com](http://www.wpi-europe.com)  
China: +86 21 6888 5517 • [chinasales@china.wpiinc.com](mailto:chinasales@china.wpiinc.com) • [www.wpiinc.net](http://www.wpiinc.net)

# Optical Measurement Systems



## All you need for optical measurements

WPI now provides a full range of components and solutions, from proprietary qualified fibers to optical assemblies and fully integrated spectroscopy measurement systems. WPI knows the importance of making measurements where the science happens, and WPI has a full range of reflectance and transmission dipping probes, flowcells and long pathlength liquid waveguide capillary cells (LWCCs) with up to 500 cm pathlength that eliminate the need to take a sample to an instrument. Our Photonics Center of Excellence in Friedberg, Germany also allows us to build proven custom configurations quickly. Whether that is a new optical design for a probe, an optical assembly or a customized measurement system.

# Solutions for All Your Optical Needs



Optical measurement is important to researchers. Whether it is fluorescence studies or measuring absorbance of trace compounds, WPI is your solution partner who has the tools you need, and the technical understanding to help you solve your application challenges.

Starting in 1993 with our first patented liquid waveguide capillary cell (**LWCC**) for precise absorbance measurement of very low concentrations, WPI now provides a full range of components and solutions, including proprietary qualified fibers, optical assemblies and fully integrated spectroscopy measurement systems.

WPI's Photonics Center of Excellence in Friedberg, Germany, allows WPI to build proven custom configurations quickly. The close collaboration with key universities such as the THM University of Applied Sciences gives WPI a leading edge.

## Features

- High quality silica/silica fibers for DUV/UV/VIS/NIR
- High quality silica/polymer clad fiber for high NA and broad wavelength range
- High flexible PMMA fibers for VIS wavelength range, large core fiber and high NA
- Industry leading solarization resistance at DUV-UV wavelengths: 180-1200 nm
- UV-enhanced broad range: 260-2200 nm
- High laser damage resistance
- Broad temperature range
- Easy customization by: core size, connector, jacketing, and configuration
- Assemblies manufactured with biocompatible materials (connectors, epoxies and furcation tubing)
- Radiation resistance – dose dependent
- Sterilizable by common methods

## Benefits


- Robust fiber and assemblies for UV use. WPI's fiber resists degradation due to UV exposure better than virtually any fiber on the market. This means less changeouts and better data comparison run to run.
- Unique extended range fiber that covers the high UV to long NIR wavelength range
- Qualified fiber allows true "plug and play" replacement
- Qualified **QVLUX** fiber is ideal for OEM and applications requiring field service or where multiple instruments must be used for a battery of tests
- Wide selection of fiber materials and sizes to manufacture the optimum fiber assembly for your application
- WPI's fiber optic assemblies are uniquely designed for biology and physiology applications with small size and excellent photon efficiency
- Custom configurations for your specific application are easy to specify and quick to deliver
- Customer support with in-house engineering and manufacturing for optical assemblies
- For critical biological work, many of our assemblies can be sterilized with gamma radiation, H<sub>2</sub>O<sub>2</sub>, ETOH or autoclave. Ask for details.

The **QVLUX** proprietary fiber series is WPI's most advanced fiber product line. **QVLUX-DUV** silica fibers are very resistant to UV light down to 185 nm but also have the uniformity and performance that OEM applications require. **QVLUX-NIR** fibers perform exceptionally well in the 700 to 2200 nm wavelength region. WPI is the only company offering DIN 58145:2018-01 certified low solarization UV fibers and assemblies for 185 m to 340 nm. For less demanding applications, we offer our standard **WVLUX** fiber in both silica and polymer cores. While WPI offers **QVLUX** and **WVLUX** assemblies in many standard configurations, custom configurations can be quickly manufactured in house.

## Basics of using optical fibers in science

Different than the fiber used for telecommunications, scientific optical fibers efficiently transmit light at many wavelengths. The amount of light that can enter a fiber and be transmitted is determined mostly by three factors:

- Core diameter size
- Material of the core
- Material of the cladding.



BIFURCATED FIBER

The difference in the refractive index between core and cladding sets the numerical aperture, i.e. the maximum angle that light can enter the fiber and transmits through the fiber.

In addition, fibers are specified by wavelength region: typically, UV-VIS (190-800 nm) or VIS-NIR (400-2200 nm). For sensing and light delivery applications we offer fused silica, polymer clad fused silica and solid polymer fibers, whichever are most suitable for your application. Silica fiber has a lower numerical aperture, but more efficiently transports light, especially at UV wavelengths. Polymer core fiber has a higher numerical aperture allowing more light to enter the fiber and may be a better choice for normal temperature conditions in a laboratory, especially in fluorescence detection.

## Qualified fibers for UV applications

Our **QVLUX** fiber has been engineered to not only provide superior transmission in the UV but have excellent consistency over length. This makes **QVLUX** a perfect choice for critical measurements or where "plug-and-play" capability is required without having to recalibrate after a fiber assembly change. Standard silica fiber exhibits a "solarization impact" when transmitting light in the UV below ~280 nm. Over time, depending on the optically transmitted power, the fiber's ability to transmit UV light drops off. Our **QVLUX** fiber offers superior resistance to solarization, lowest basic attenuation at 200 nm (< 1 dB/m) and the lowest solarization at 214 nm on the market. WPI is the only company offering DIN 58145:2018-01 certified low solarization UV fibers and assemblies for 185 nm to 340 nm. Our standard **WVLUX** fused silica optical fibers are designed for value, but keep in mind that they lose their transmission over time and at high power due to the solarization effects of UV light below 280 nm.

## Qualified fibers for NIR applications

We also offer **QVLUX** fiber for the 700-2200 nm range that works well in demanding petro-chemical, food and pharmaceutical applications. This region is very difficult for standard fibers from traditional suppliers since they cannot control their processes or glass chemistry sufficiently to have the purity required for high-standard of uniformity and superior optical throughput that these applications demand.

### WORLD PRECISION INSTRUMENTS

## Match your Fiber to your Needs

With WPI you can match the fiber, the assembly and even the instrument exactly to your application needs:

Fiber	QVLUXDUV	QVLUXNIR	WVLUXDUV
FIBER MATERIAL	Qualified Silica	Qualified Silica	Specialty Silica
WAVELENGTH RANGE (nm)	180-1200	260-2200	200-1200
AVAILABLE CORE DIAMETERS (µm)	400, 600	200, 400, 600	100, 200, 300, 400, 600
KEY SPECIFICATIONS	<ul style="list-style-type: none"> <li>Lowest attenuation @ 180 nm</li> <li>Lowest solarization @ 214 nm</li> <li>NA = 0.22 @ 633 nm</li> </ul>	<ul style="list-style-type: none"> <li>Low solarization @ 248 nm</li> <li>NA = 0.22 @ 633 nm</li> </ul>	<ul style="list-style-type: none"> <li>Low solarization @ 214 nm</li> <li>NA = 0.22 @ 633 nm</li> </ul>



OPTICAL FIBER WITH SILICONE MONOCOIL TUBING

Fiber	WVLUXPV	WVLUXUVIS	WVLUXNXR	WVLUXHNA
FIBER MATERIAL	Silica/PolyClad	Silica	Silica	PMMA
WAVELENGTH RANGE (nm)	300-1200	260-1200	320-2200	365-720
AVAILABLE CORE DIAMETERS (µm)	200, 300, 400, 600, 800, 1000	100, 200, 300, 400, 600, 800, 1000	200, 400, 600	250, 500, 750, 1000, 1500, 2000, 3000
KEY SPECIFICATIONS	<ul style="list-style-type: none"> <li>365 nm for low power UV applications</li> <li>NA = 0.37 @ 633 nm</li> </ul>	<ul style="list-style-type: none"> <li>Great general purpose fiber for UV-VIS spectroscopy</li> <li>NA = 0.22 @ 633 nm</li> </ul>	<ul style="list-style-type: none"> <li>Extended NIR range</li> <li>NA=0.22 @ 633 nm</li> <li>Custom core diameters available</li> </ul>	<ul style="list-style-type: none"> <li>Solid polymer fiber (PMMA)</li> <li>Highly flexible</li> <li>Largest range of core diameters</li> </ul>

## Fiber Selection Guide

WPI offers standard, bifurcated and cross fiber assemblies. For Y assemblies and X assemblies, the split point is approximately at 50 cm. The standard jacketing is silicone monocoil.

Ordering assemblies with our fiber is simple! Standard lengths are 100 cm and terminated with SMA connectors.

**1** Choose fiber type.

**2** Choose style.

**3** Choose diameter (in microns).

The available diameters are limited by the fiber type chosen in step 1. Refer to the tables above for the options.

**4** Select length (in centimeters)

If your assembly uses all the same fibers, start with the standard part number defined above. The standard fiber is 100 cm.

**5** Choose connector type.

For example, the part number for a standard QVLUX-UV fiber with a Y assembly and 600 µm core diameter would be **QVLUXDUV-Y-0600**. For The part number for a custom cable of 250 cm with FC connectors using 750 µm high NA polymer fiber would be **WVLUXHNA-S-0750-250-F125**.


For a mixed fiber, mixed connectors, X-cross assemblies or a special configuration, contact WPI (wpi@wpiinc.com).



1x1 Standard – S  
Optical fibers with 2 connectors and furcation tubing

1x2 Bifurcated – Y  
Split or combines similar intensities by mixing different or same fiber type.

2x2 Cross – X  
Cross fiber assembly with mixed fiber types (call for details)



Standard Option  
Core diameters available depend on fiber chosen.

# High Performance Spectrophotometer

*A unique multiple long pathlength sample cell for absorbance spectroscopy*



## Features

- Process Control & Oceanography
- Rugged system for laboratory and onboard measuring
- Portable & easy to use
- User-selected optical path lengths: 2, 10, 50 & 200 cm
- Highly sensitive and stable

## Benefits

- Designed with NASA for colored dissolved organic matter in seawater and fresh water
- Highly sensitive with extended dynamic range for UV and VIS absorbance measurements, with fewer complications associated with standard long pathlength systems
- Portable system for field operations and mobile lab environments

## Application

- CDOM - Colored dissolved organic matter
- QFT - Quantitative Filter Technique

**UltraPath™** is a unique high-performance spectrophotometer system offering selectable optical path lengths of 2, 10, 50 and 200 cm. The instrument operates in the wavelength range of 250 to 730 (**UPUV**) or 380 to 730 nm (**UPVIS**) and has an exceptional dynamic range. Designed for the detection of low absorbing species in aqueous solutions, **UltraPath** is an ideal tool for any study requiring precise and highly sensitive spectroscopic determination of analytes, either in the lab or in the field.

### Designed with NASA for CDOM

**UltraPath** was developed by WPI under a collaborative agreement with NASA (Stennis Space Center) for the spectroscopic determination of colored dissolved organic matter (CDOM) in seawater and fresh water environments. It can be used in the laboratory and in the field (i.e., at sea). CDOM concentrations vary significantly between open ocean samples with low CDOM (e.g., 0.007 m<sup>-1</sup> at 380 nm), and high CDOM freshwater environments (e.g., 10-20 m<sup>-1</sup> at 380 nm). To address these problems, the design requirements of **UltraPath** mandated the development of a rugged portable system capable of high sensitivity measurements across a wide dynamic range. The **UltraPath** system meets these stringent design criteria and enables reliable measurement of CDOM in the range of 0.002 m<sup>-1</sup> to 200 m<sup>-1</sup> (250 to 730 nm).

## Highly Sensitive, Extended Dynamic Range

**UltraPath** has four optical pathlengths contained within a single sample cell (i.e., 2 cm, 10 cm, 50 cm and 200 cm). The pathlengths are selectable, offering a very high sensitivity and an extended dynamic range for UV and VIS absorbance measurements. The fluid path of the sample cell is optimized to produce a laminar flow that is virtually free of interference from trapped air bubbles and adherence of dissolved substances to the cell wall. In particular, the design greatly minimizes the problems commonly found with flow cells of long optical pathlengths.

- Reduces the risk of trapping dust particles.
- Reduces contamination of fibers or particulate matter inside the cell.

The UV/VIS **UltraPath** system (**UPUV**) includes a low noise photodiode array-based **TIDAS S-300 UV/VIS** spectrophotometer with a deuterium/halogen light source. The VIS **UltraPath** system includes a low noise photodiode array-based **TIDAS E-Base UV/VIS** spectrometer module and a FO-6000 visible light source. Light is coupled from the light source to the sample cell and from the sample cell to the detector using two fused silica fibers. A peristaltic pump (**PeriPro-4LS**) draws the sample into the **UltraPath** sample cell.

A standard PC or laptop (not included) is connected to the **TIDAS** units using an RJ-45 Ethernet interface.

### Portable system for field environments

The system is designed for mobility. The components of the **UltraPath** system are designed to function over a broad range of laboratory and field environments.



*The Mini Star is a compact, lightweight peristaltic pump that fits just about anywhere. It can be mounted directly on the bench, in a regular rack or to a post. See page 53.*

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com Germany: +49 (0)6031 1602171 • wpide@wpi-europe.com • www.wpi-europe.com  
Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • www.wpiinc.com China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net

## References

M. G. Novak, A. Mannino, N. Nelson, E. D'Sa, R. Miller, J. Werdell, R. Del Vecchio, C. Del Castillo, J. Chaves, J. F. Berthon, E. Boss, M. Tzortziou, A. Neeley, S. Freeman, A. Bricaud, R. Röttgers, A. Matsuoka, M. Belz, N. Blough, Measurements of CDOM absorption spectra using different Instruments and techniques: A round robin exercise and extensive field data set., *Ocean Optics XXII*, (2014)

N. B. Nelson, D. A. Siegel, C. A. Carlson, C. Swan, W. M. Smethie Jr. and S. Khatiwala. 2007. Hydrography of chromophoric dissolved organic matter in the North Atlantic. *Deep-Sea Res. I*. 54: 710 – 731.

V. Kitidis, A. P. Stubbins, G. Uher, R. C. Upstill Goddard, C. S. Law, E. M. S. Woodward, "Variability of chromophoric organic matter in surface waters of the Atlantic Ocean", *Deep Sea Research Part II: Topical Studies*, Vol. 53, Issue 14-16, 2006, p. 1666-1684.

R. L. Miller, M. Belz, C. Del Castillo, R. Trzaska, "Determining CDOM Absorption Spectra in Diverse Coastal Environments Using a Multiple Pathlength, Liquid Core Waveguide System", *Continental Shelf Research*, July 2002, 22:9, p 1301-1310.

"System Analyzes Water Samples at Sea", NASA Aerospace Technology Innovation, 2001, 9 (5). <http://nctn.hq.nasa.gov/innovation/innovation95/3-tchtrans2.html>

R. L. Miller and E. D'Sa. "Evaluating the influence of CDOM on the remote sensing signal in the Mississippi River Bight". In *Eos Transactions AGU Ocean Sciences*, 2002. Honolulu, HI, p. 171.

E. D'Sa, R.L. Miller and R. Trzaska. "Aparent Optical Properties in Waters Influenced by the Mississippi River", *Proceedings of the Seventh Thematic Conference, Remote Sensing for Marine and Coastal Environments*, 2002, 6 pg, Miami, FL.

R. L. Miller, C. Hall, C. Del Castillo, B. McKee and M. Dagg. "Bio-optical Properties of the Mississippi River Plume and Adjacent Shelf." *ASLO Aquatic Sciences*, Albuquerque, NM, 2001.

R. L. Miller, M. Belz and S. Y. Liu, "Measuring the absorption of CDOM in the field using a multiple pathlength liquid waveguide system", *Ocean Optics XV*, paper 1308, Monaco, October 2000.

## ULTRAPATH SPECIFICATIONS

DYNAMIC RANGE	5 $\mu\text{AU}/\text{cm}$ to 1 $\text{AU}/\text{cm}$ 0.002 $\text{m}^{-1}$ to 200 $\text{m}^{-1}$
WAVELENGTH RANGE	250 nm – 730 nm (UPUV) 380 nm – 730 nm (UPVIS)
WAVELENGTH RESOLUTION (FWHM)	5 nm
NOISE (PEAK TO PEAK)	< 0.2 $\text{mAU}$
DRIFT	< 1 $\text{mAU}/\text{h}$
OPTICAL PATHLENGTH	2, 10, 50 & 200 cm (user selectable)
SAMPLE CELL INNER DIAMETER	2 mm
CELL VOLUME	10 mL (at 200 cm pathlength)
SAMPLE INLET / OUTLET	1/8"
FIBER INPUT/OUTPUT	600 $\mu\text{m}$
SOLVENT RESISTANCE	Most organic and inorganic solvents
SHIPPING WEIGHT	UPUV: 44 lb. (20 kg) UPVIS: 33 lb. (15 kg)

## ORDERING INFORMATION

<b>UPVIS</b>	Ultrapath System, Visible Light
<b>UPUV</b>	Ultrapath System, Ultraviolet & Visible Light
<i>The UltraPath system includes Multiple pathlength cell, Tidas S300 (UPUV) with integrated light source or Tidas E Base (UPVIS) with TidasDAQ3 software, two WVLXUVIS-S-600 optical fibers, PeriStar Pro peristaltic pump, silicone tubing, sample injector and Waveguide Cleaning Kit.</i>	
<b>501609</b>	Waveguide Cleaning Kit
<b>89575</b>	QFT1, Fiber Optic Holder for Glass Fiber Filters
<b>WVLXUVIS-S-1000</b>	Fiber Optic Cable 1 m, SMA, 1000 $\mu\text{m}$ core, UV enhanced

## Applications

### Simple measurements for particulate absorption

WPI's filter holder for particulate absorption measurements is specially designed for field use. It is rugged and portable. It performs as well with a laboratory-based spectrophotometer. It can be directly connected to WPI's line of fiber optic spectrometers and light sources. Improve your particulate absorption measurements on site, instead of loss of spectral information by transporting your sample to a laboratory (Sosik, *Limnol Oceanogr.* Vol 44, 1139-1141, 1999).

### How does it work?

Particulate absorption of fresh and seawater can be determined by filtering a known amount of sample through a Glass Fiber Filter (GF/F) and measuring the particulate absorption coefficient concentrated on the filter. This technique is called quantitative filter technique (QFT) and corrects for pathlength amplification, an effect of scattering (Mitchell, *SPIE Vol 1302*, 137-148, 1990).

### Advantages of QFT


WPI's filter holder for glass fiber filters is designed for particulate absorption or cuvette measurements in a single instrument. A further advantage of the filter holder is its large beam diameter of 5 mm, resulting in "averaging out" of larger non-organic particles frequently found on the filter pad when using natural samples. The removable filter fixture allows simple filter alternation and cleaning.

### QFT for particulate absorption

Particulate absorption can be measured by the well established Quantitative Filter Technique (QFT). WPI now offers a fiber optic filter holder for Glass Fiber Filters (QFT1, page 168) which can be used with the spectrometer (Tidas E Base) and light source (D4H or FO6000) supplied with the UltraPath. With this accessory, particulate absorption can be measured on site, avoiding loss of spectral information due to freezing and shipping particulate samples to a laboratory.



### Ponceau S measured with 200 cm pathlength



*Ponceau S absorption measured with UltraPath (200 cm cell). Ponceau S was dissolved in Millipore water.*

Special attention should be drawn to the exceptional sensitivity of UltraPath enabling detection of CDOM absorption below  $0.03 \text{ m}^{-1}$ . To exemplify the performance of the UltraPath in laboratory chemistry and process control, Ponceau S absorbance was measured with the 200 cm pathlength of an UltraPath. Normalizing the Ponceau absorbance graph to  $\text{AU}/\text{cm}$ , the range of this measurement is  $150 \mu\text{AU}$  with a noise level below  $2 \mu\text{AU}$  peak to peak. Sub-nanomolar concentration of this dye can clearly and reliably be detected, which is a novelty in absorbance based spectroscopy.

# Liquid Waveguide Capillary Cell

## Long pathlengths for small sample volumes

### Features

- 10–500 cm pathlength
- 10–500 fold sensitivity improvement in comparison to 1 cm cuvette
- 0.55 mm ID for low sample volume sampling
- 2 mm ID for unfiltered liquid samples
- SMA 905 fiber optic connections
- 250 nm – 720 nm wavelength range with MilliPore water (see datasheet for details)

### Benefits

- Adapts to most fiber optic detection systems
- 20 years of manufacturing experience
- Low UV drift

### Applications

- Trace detection of nutrients (nitrite, nitrate, phosphate, iron) in seawater
- Environmental and oceanographic monitoring
- Drinking water analysis
- Colored dissolved organic matter (CDOM)
- Process control

UV/VIS/NIR absorbance spectroscopy is governed by Beer's Law, where the absorbance signal is proportional to chemical concentration, light path length and the compound's specific molar absorption coefficient. Typical optical pathlengths of cuvettes and flow cells are between 0.1 cm and 10 cm. Longer pathlengths are difficult to achieve due to mechanical constraints. Liquid Waveguide Capillary Cells (LWCCs) fill this gap. LWCCs are fiber optic flow cells that combine an increased optical pathlength (10–500 cm) with small sample volumes ranging from 2.4  $\mu$ L to about 3 mL. Compared with a standard 1 cm cell, a 1 mAU signal is enhanced one hundred fold with a 100 cm flowcell to 100 mAU, using WPI's patented aqueous waveguide technology.\* They can be connected via optical fibers to a spectrophotometer with fiber optic capabilities. Ultra-sensitive absorbance measurements can be performed in the ultraviolet (UV), visible (VIS) and near-infrared (NIR) to detect low sample concentrations in a laboratory or process control environment.

### Your sample is the core of a light guide

WPI's Liquid Waveguide Capillary Cells are made of fused silica tubing with an outer coating of a low refractive index polymer. Your liquid sample is guided through the capillary and represents the core of the waveguide. The hydrophilic character of the fused silica capillary inner wall results in high signal stability and easy removal of air bubbles trapped in the flow cell. However, the transmission of the LWCC is mainly dependent on the intrinsic attenuation of the sample liquid. In case of water, a usable wavelength range from 250 nm to 720 nm wavelength can be observed in a 100 cm pathlength LWCC. Using a 500 cm pathlength LWCC will reduce that transmission range from 300 nm to about 700 nm. However, when switching from water to methanol as a solvent, transmission into the NIR are possible with suitable light sources and detectors.

### Connections

The LWCC-3xxx series of flow cells uses traditional HPLC type 10-32 coned port fittings with 1/32 inch tubing for liquid connection and 500  $\mu$ m SMA fiber optic adapters for light input and output. The LWCC-4xxx series of flow cells uses 1/4-28 flangless flat bottom fittings with 0.125" tubing and 600  $\mu$ m SMA fiber optic adapters. Liquid can be pumped into the flow cells using (in the simplest case) a sample injector (58006) and a ministar peristaltic pump (MINISTAR). The LWCC may be connected directly to a fluid injection analysis (FIA) system or to a gas segmented fluid injection analysis (GFIA) system via a debubbler. Finally, for routing discrete measurements, WPI's LWCC Injection system (89372) may be used where the sample is injected into a constant flow via an injection loop of 3–4 times the internal flow cell volume to ensure a stable baseline and avoid the introduction of micro air bubbles into the flow cell.

### Applications

When the entire spectral shape of an absorbance curve is required for analysis, WPI's TIDAS E Base spectrometer with a D4H or a FO-6000, or the TIDAS S300 spectrophotometer can be used. LWCCs have



MINISTAR,  
see page 53.

LWCC-3050

been used in a variety of applications, such as liquid chromatography, stopped-flow and colorimetric detection, drinking water analysis, as well as environmental and oceanographic monitoring systems.

### Accessory: LWCC Injection System

For flow analysis, including simple fluid injection analysis (FIA) setups, add WPI's LWCC injection system (89372). A selection valve provides baseline or cleaning solutions to the sample stream. The injection valve injects a sample into the stream, avoiding the introduction of air bubbles or changes of flow rate.



89372 LWCC Injector System (pump and LWCC not included)

### Related Patents

Micro Chemical Analysis Employing Flow Through Detectors, 1995, U.S. Patent No. 5,444,807.

Aqueous Fluid Core Waveguide, 1996, U.S. Patent No. 5,507,447.

Long Capillary Waveguide Raman Cell, 1997, U.S. Patent No. 5,604,587.

Chemical Sensing Techniques Employing Liquid-Core Optical Fibers, U.S. Patent No. 6,016,372

### References

Yuhan Liu, Keding Lu, Huabin Dong, Xin Li, Peng Cheng, Qi Zou, Yusheng Wu, Xingang Liu, Yuanhang Zhang, (2016) In situ monitoring of atmospheric nitrous acid based on multi-pumping flow system and liquid waveguide capillary cell. *Journal of Environmental Sciences Volume 43*, 2016, Pages 273–284 <http://dx.doi.org/10.1016/j.jes.2015.11.034>

Ye, C., Zhou, X., Pu, D., Stutz, J., Festa, J., Spolaor, M., ... Knotte, C. (2016). Rapid cycling of reactive nitrogen in the marine boundary layer. *Nature*, 532(7600), 489–491. <http://doi.org/10.1038/nature17195>

Catelani, T. A., Tóth, I. V., Lima, J. L. F. C., Pezza, L., & Pezza, H. R. (2014). A simple and rapid screening method for sulfonamides in honey using a flow injection system coupled to a liquid waveguide capillary cell. *Talanta*, 121, 281–7. <http://doi.org/10.1016/j.talanta.2013.12.034>

### WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • [wpiuk@wpi-europe.com](mailto:wpiuk@wpi-europe.com) • [www.wpi-europe.com](http://www.wpi-europe.com)  
Brazil: 011 55 13 40629703 • [info@brazil.wpiinc.com](mailto:info@brazil.wpiinc.com) • [www.wpiinc.com](http://www.wpiinc.com)

Germany: +49 (0)6031 1602171 • [wptide@wpi-europe.com](mailto:wptide@wpi-europe.com) • [www.wpi-europe.com](http://www.wpi-europe.com)  
China: +86 21 6888 5517 • [chinasales@china.wpiinc.com](mailto:chinasales@china.wpiinc.com) • [www.wpiinc.net](http://www.wpiinc.net)



## LWCC SPECIFICATIONS

	LWCC-3050	LWCC-3100	LWCC-3250	LWCC-3500	LWCC-4010	LWCC-4050	LWCC-4100
OPTICAL PATHLENGTH	50 cm	100 cm	250 cm	500 cm	10 cm	50 cm	100 cm
INTERNAL VOLUME	125 $\mu$ L	250 $\mu$ L	625 $\mu$ L	1250 $\mu$ L	0.31 mL	1.57 mL	3.1 mL
FIBER CONNECTION	500 $\mu$ m SMA	500 $\mu$ m SMA	500 $\mu$ m SMA	500 $\mu$ m SMA	600 $\mu$ m SMA	600 $\mu$ m SMA	600 $\mu$ m SMA
TRANSMISSION @254nm*	20	10	5	-	5	3	2
TRANSMISSION @540 nm*	35	30	25	20	5	4	4
NOISE [mAU]**	<0.1	<0.2	<0.5	<1.0	<0.1	<0.2	<0.2
MAXIMUM PRESSURE	100 PSI	100 PSI	100 PSI	100 PSI			
WETTED MATERIAL	PEEK, Fused Silica, PTFE						
LIQUID INPUT	LWCC-3xxx series: 10-32 coned port fittings for 1/16" tubing LWCC-4xxx series: 1/4-28 flangeless flat bottom fitting for 1/8" tubing						

\* Referenced using coupled 500  $\mu$ m fibers or LWCC-3xxx series and 600  $\mu$ m fibers for LWCC-4xxx series

\*\* Measured using ASTM E685-93

\*\*\* A one-meter waveguide of 550  $\mu$ m internal diameter requires approximately 1.5PSI for water flow of 1.0 mL/min.



501609

**Waveguide Cleaning Kit (#501609)**, above, includes the most commonly needed cleaning solutions for the LWCC waveguides. The **LWCC Start-up Kit (#KITLWCC)**, at right, includes two fiber optic cables (#WVLXUVIS-S-600), Sample Injector Assembly (#58006), MiniStar™ Peristaltic Pump, and Waveguide Cleaning Kit (#501609).

## ORDERING INFORMATION

<b>LWCC-3050</b>	Liquid Waveguide Capillary Cell, 50 cm pathlength, 0.55 mm ID
<b>LWCC-3100</b>	Liquid Waveguide Capillary Cell, 100 cm pathlength, 0.55 mm ID
<b>LWCC-3250</b>	Liquid Waveguide Capillary Cell, 250 cm pathlength, 0.55 mm ID
<b>LWCC-3500</b>	Liquid Waveguide Capillary Cell, 500 cm pathlength, 0.55 mm ID
<b>LWCC-4010</b>	Liquid Waveguide Capillary Cell, 10 cm pathlength, 2 mm ID
<b>LWCC-4050</b>	Liquid Waveguide Capillary Cell, 50 cm pathlength, 2 mm ID
<b>LWCC-4100</b>	Liquid Waveguide Capillary Cell, 100 cm pathlength, 2 mm ID

### Accessories

A sample injector assembly can be used to conveniently fill an LWCC with sample solution using a peristaltic pump. The LWCC requires two optical fibers to connect to spectrophotometer system. Choose between anti-solarized 400  $\mu$ m core or UV-enhanced cables (ordered in 1 or 3 meter lengths).

<b>89372</b>	LWCC Injection System
<b>58006</b>	Sample Injector Attachment
<b>PERIPRO-4LS</b>	Peri-Star™ Pro Peristaltic Pump (See page 52.)
<b>MINISTAR</b>	Miniature Peristaltic Pump, 1-channel (See page 53.)
<b>WVLXUVIS-S-600</b>	Fiber Optic Cable, 1 m, SMA, 600 $\mu$ m core, UV-enhanced
<b>501609</b>	Waveguide Cleaning Kit (available only in USA)
<b>KITLWCC</b>	LWCC Start-up Kit*
<b>58450</b>	Kit, Adapter Syringe, LWCC

\*includes WVLXUVIS-S-600 (two), 58006, MINISTAR, 501609

## Low Volume Flow Cell

### For FIA, HPLC and process analysis

#### Features

- UV/VIS flow cell for absorbance
- Low internal volume
- Fits 500 and 600  $\mu$ m fibers
- High UV transmission
- 0–10 mL/min flow rate

#### Benefits

- High efficiency coupling
- Low refractive index offset
- Fits WPI TIDAS systems

#### Applications

- FIA, GFIA, HPLC, Optofluidics
- Process control

**MicroLWCC** is a new fiber optic low volume flow cell for UV/VIS/ NIR absorbance analysis. Based on WPI's established liquid core waveguide technology, the analyte solution functions as the core of a fluid filled light waveguide. Wetted parts in the sample cell light path are PEEK, fused silica and PTFE. Optical fibers are used



to transport light to and from the sample cell. The cell can be used in biochemistry for DNA, RNA & protein quantification, colorimetric nutrient and trace metal analysis, drug discovery and dissolution testing, process control, and HPLC analysis.

## LWCC-M SPECIFICATIONS

	LWCC-M-10	LWCC-M-50	LWCC-M-100
OPTICAL PATHLENGTH	10 mm	50 mm	100 mm
INTERNAL VOLUME	2.4 $\mu$ L	12 $\mu$ L	24 $\mu$ L
REFRACTIVE INDEX @ 280 nm**	< 15 mAU	< 20 mAU	< 30 mAU
TRANSMISSION @ 254 nm *	25%	20%	15%
TRANSMISSION @ 500 nm	40%	35%	30%
WAVELENGTH RANGE	200 – 1000 nm		
FIBER CONNECTION [ $\mu$ m]	500 (SMA)		
MAXIMUM PRESSURE	50 Bar		
WETTED MATERIALS	PEEK, Fused Silica, PTFE		

\* Reference: 2 \* 600  $\mu$ m Fiber, butt-coupled

\*\* Measured using ASTM E 685 - 93


WPI U.S. Patents: 5,444,807; 5,570,447; 5,604,587; 6,603,556; 6,385,380.

## ORDERING INFORMATION

<b>LWCC-M-10</b>	Low Volume Flow Cell, 10 mm pathlength
<b>LWCC-M-50</b>	Low Volume Flow Cell, 50 mm pathlength
<b>LWCC-M-100</b>	Low Volume Flow Cell, 100 mm pathlength

# Measure Colored Dissolved Organic Matter

Detect lower concentrations of solutes & conduct broader range of measurement



Colored Dissolved Organic Matter (CDOM) is organic matter whose optical properties are measurable using WPI's Liquid Waveguide Capillary Cell (LWCC). CDOM occurs naturally in water systems and is derived from organic tannins. CDOM concentration depends on the location where samples are taken, with coastal waters showing higher CDOM concentrations compared to open-ocean waters. In addition, CDOM absorption depends on open-ocean water depth.

LWCCs are fiber optic flow cells that combine an increased optical pathlength (1–500 cm) with small sample volumes ranging from 2.4  $\mu\text{L}$  to about 3 mL.

## Features

### Select models:

- 10–250 cm pathlength (depending on LWCC model)
- 0.31–3.1 mL internal volume (depending on LWCC model)
- Up to 230–730 nm wavelength range (depending on LWCC model)

### All models:

- 10 to 500 fold sensitivity improvement in comparison to 1 cm cuvettes
- 0.55 mm ID for low sample volume sampling
- 2 mm ID for filtered and unfiltered liquid samples
- SMA 905 fiber optic connections

## Benefits

- Improved dynamic range for broader range of absorbance measurements
- Improved sensitivity of measurement
- Detect low concentrations of solutes
- Measurements can be made using small sample volumes
- Compact, portable system for real time measurement aboard ship

## Determining LWCC pathlength

Selecting the proper equipment is imperative when setting up your system. Here are a few considerations:

- Select the LWCC pathlength based on the desired absorption range. Here are some reference values for different water types:
  - Fresh water with absorption range  $> 4.0 \text{ m}^{-1}$
  - Coastal-ocean waters with absorption range 1.0–4.0  $\text{m}^{-1}$
  - Open ocean water with absorption range  $< 1.0 \text{ m}^{-1}$
- Select the usable internal volume.

Next, you can select components to complete your CDOM analysis system, depending on the selected LWCC. The effective pathlength of

WPI's LWCC is defined as the equivalent pathlength of the cell, if it is assumed that the LWCC strictly follows the Beer-Lambert law.

Typically, the longer LWCC pathlength is used to increase the sensitivity when the maximum absorbance values are supposed to be  $< 0.1$  AU (Absorbance Unit). Inversely, when absorbance measurements are above 1.4 AU, the LWCC pathlength should be decreased to ensure that measurements still remain within the linear range of the LWCC detection system.

LWCC Type	Pathlength (cm)	Noise (mAU)	Absorbance Range (mAU)	Absorption Range ( $\text{m}^{-1}$ )	Internal Volume (mL)
LWCC-4010	10	$< 0.1$	0.5 - 1400	0.012 - 32.0	0.31
LWCC-4050	50	$< 0.2$	1.0 - 1400	0.005 - 6.4	1.57
LWCC-4100	100	$< 0.5$	2.5 - 1400	0.006 - 3.2	3.1
LWCC-3050	50	$< 0.1$	0.5 - 1400	0.002 - 6.4	0.125
LWCC-3100	100	$< 0.2$	1.0 - 1400	0.002 - 3.2	0.250
LWCC-3250	250	$< 0.5$	2.5 - 1400	0.002 - 1.2	0.625

*The useful absorption range calculation is based on the absorbance detection limits of the LWCC, considering a wavelength range of 300–700 nm.*

Absorbance measurements obtained with WPI's LWCC and Tidas S300 UV/VIS spectrophotometer are linear up to 1.4 AU. The measured absorbance can be converted to the spectral absorption coefficient  $a(\lambda)$ , commonly used in oceanography for CDOM measurements. Absorbance and spectral absorption are related by the formula:  $a(\lambda) = 2.303 A(\lambda) / L$ , where 2.303 is the conversion factor from decimal to natural logarithmic,  $A(\lambda)$  is the absorbance at wavelength  $\lambda$  and  $L$  the LWCC pathlength.

## Common applications of CDOM detection



- Determine the biogeochemical cycles, e.g., the organic carbon-based cycle in the ocean.
- Monitor and map surface-water masses.
- Measure the UV light penetration into the ocean to determine:
  - Photosynthesis reaction with effects on phytoplankton population
  - Effect on oceanic food chains
  - Atmospheric oxygen concentration
- Monitor the light absorption of CDOM as it relates to heat storage and the decline of sea ice.

Two typical absorption spectra recorded with an UltraPath (UPUV) of a seawater and a fresh water sample collected in November 2007 are shown in the figures. Due to their high absorbance, both samples were analyzed in the 10 cm pathlength. The CDOM sample labeled Mayagüez Bay is from oligotrophic, low productive waters with high salinity collected off the west coast of Puerto Rico in the Mayagüez Bay and was measured with a 200 cm pathlength.

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • www.wpiinc.com

Germany: +49 (0)6031 1602171 • wptide@wpi-europe.com • www.wpi-europe.com  
China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net



Two typical absorption spectra measured using UltraPath. The sample labeled "Sarasota Bay" is a CDOM sample with 34 PSU salinity collected from Sarasota Bay (Nov. 2007), and the sample labeled "Pond" is a highly concentrated CDOM sample collected from a local pond in Sarasota, Florida (Nov. 2007).

CDOM Sample "Mayagüez Bay" was collected from the high salinity oligotrophic waters of Mayagüez Bay on the west coast of Puerto Rico (2001). Data courtesy of NASA Stennis Space Center.

## System Configurations

Refer to the table on the previous page for selecting the LWCC pathlength.

### The CDOM-FRESH System > 4 m<sup>-1</sup>

LOW VOLUME		HIGH VOLUME	
Product Description	Item #	Product Description	Item #
<b>CDOM-FRESH-LV System Includes:</b>	CDOM-FRESH-LV	<b>CDOM-FRESH-HV System Includes:</b>	CDOM-FRESH-HV
LWCC, 50 cm pathlength	LWCC-3050	LWCC, 10 cm pathlength	LWCC-4010
Photo Diode Array (PDA) Spectrophotometer System, UV/VIS (190-720 nm), with integrated deuterium/halogen light source	505067	PDA Spectrophotometer System, UV/VIS (190-720 nm), with integrated deuterium/halogen light source	505067
(2) UV-Enhanced Fiber Optic Cables, 1 m, 600 µm Core	WWLUXUVIS-S-600	(2) UV-Enhanced Fiber Optic Cables, 1 m, 600 µm Core.	WWLUXUVIS-S-600
Ministar Peristaltic Pump	504011	PeriStar Pro Pump	PeriPro-4LS
LWCC Injection System	89372	Injector Kit	72100
TTL Control Module for Ministar and/or Peristar	503120	TTL Control Module for Ministar and/or Peristar	503120
UV-Enhanced Fiber Optic Cable	WWLUXUVIS-S-200	UV-Enhanced Fiber Optic Cable	WWLUXUVIS-S-200
(2) SMA Bulkhead feed through Connector/Coupler	13395	(2) SMA Bulkhead feed through Connector/Coupler	13395

### The CDOM-COAST System 1-4 m<sup>-1</sup>

LOW VOLUME		HIGH VOLUME	
Product Description	Item #	Product Description	Item #
<b>CDOM-COAST-LV System Includes:</b>	CDOM-COAST-LV	<b>CDOM-COAST-HV System Includes:</b>	CDOM-COAST-HV
LWCC, 100 cm pathlength	LWCC-3100	LWCC, 50 cm pathlength	LWCC-4050
PDA Spectrophotometer System, UV/VIS (190-720 nm) with integrated deuterium/halogen light source	505067	PDA Spectrophotometer System, UV/VIS (190-720 nm) with integrated deuterium/halogen light source	505067
(2) UV-Enhanced Fiber Optic Cables. 1 m, 600 µm Core.	WWLUXUVIS-S-600	(2) UV-Enhanced Fiber Optic Cables. 1 m, 600 µm Core.	WWLUXUVIS-S-600
Ministar Peristaltic Pump	504011	PeriStar Pro Pump	PeriPro-4LS
LWCC Injection System	89372	Injector Kit	72100
TTL Control Module for Ministar and/or Peristar	503120	TTL Control Module for Ministar and/or Peristar	503120
UV-Enhanced Fiber Optic Cable	WWLUXUVIS-S-200	UV-Enhanced Fiber Optic Cable	WWLUXUVIS-S-200
(2) SMA Bulkhead feed through Connector/Coupler	13395	(2) SMA Bulkhead feed through Connector/Coupler	13395

### The CDOM-Ocean System < 1 m<sup>-1</sup>

LOW VOLUME		HIGH VOLUME	
Product Description	Item #	Product Description	Item #
<b>CDOM-OCEAN-LV System Includes:</b>	CDOM-OCEAN-LV	<b>CDOM-OCEAN-HV System Includes:</b>	CDOM-OCEAN-HV
LWCC, 250 cm pathlength	LWCC-3250	LWCC, 100 cm pathlength	LWCC-4100
PDA Spectrophotometer System, UV/VIS (190-720 nm) with integrated deuterium/halogen light source	505067	PDA Spectrophotometer System, UV/VIS (190-720 nm) with integrated deuterium/halogen light source	505067
(2) UV-Enhanced Fiber Optic Cables. 1 m, 600 µm Core	WWLUXUVIS-S-600	(2) UV-Enhanced Fiber Optic Cables. 1 m, 600 µm Core.	WWLUXUVIS-S-600
Ministar Peristaltic Pump	504011	PeriStar Pro Pump	PeriPro-4LS
LWCC Injection System	89372	Injector Kit	72100
TTL Control Module for Ministar and/or Peristar	503120	TTL Control Module for Ministar and/or Peristar	503120
UV-Enhanced Fiber Optic Cable	WWLUXUVIS-S-200	UV-Enhanced Fiber Optic Cable	WWLUXUVIS-S-200
(2) SMA Bulkhead feed through Connector/Coupler	13395	(2) SMA Bulkhead feed through Connector/Coupler	13395

# Filter Holder for Glass Fiber Filters

*Designed for field use*

## Features

- Simple measurements for particulate absorption
- Rugged and portable
- Performs as well as a laboratory based spectrophotometer

## Benefits

- QFT corrects for the effect of scattering
- Averages out larger, non-organic particles

## Applications

- Measurement of UV/VIS absorption of particulate matter in seawater

WPI's filter holder for particulate absorption measurements is specially designed for field use. It is rugged and portable. It performs as well as a laboratory based spectrophotometer. It can be directly connected to WPI's line of fiber optic spectrometers and light sources. Instead of collecting your samples, transporting them to a laboratory, and accepting the loss of spectral information associated with it (Sosik, 1999), particulate absorption can now be measured on site.

### QFT corrects for the effect of scattering

Particulate absorption of fresh water and seawater can be determined by filtering a known amount of sample through a Glass Fiber Filter (GF/F) and measuring the particulate absorption coefficient  $a_p(\lambda)$  concentrated on the filter. This technique is called quantitative filter technique (QFT) and corrects for the pathlength amplification, an effect of scattering. The correction of the pathlength amplification and the correction of the non-linear relationship between the optical density of samples on a Whatman GF/F filter and in suspension are discussed in Mitchell (1990).

### Averages out larger, non-organic particles

A significant advantage of the filter holder is its large beam diameter of 5 mm, resulting in "averaging out" of larger non-organic particles frequently found on the filter pad when using natural samples. The removable filter fixture allows simple filter alternation and cleaning.

### System Requirements

The optical throughput of QFT1 equipped with a classical GF/F filter is very low and requires a matched light source/spectrometer system. TIDAS E Base UV/VIS (504718) in combination with a D4H UV/VIS light source



89575

or a **FO-6000** VIS light source, as it is offered in the Ultrathin system, are ideally suited, as well as TIDAS S300 UV/VIS (505067). Light should be coupled into the **QFT1** with a 1000  $\mu\text{m}$  fiber (**WVLUXUVIS-S-1000**) and from the **QFT1** to a spectrometer with a 600  $\mu\text{m}$  fiber (**WVLUXUVIS-S-600**).

## SPECIFICATIONS

GF/F FILTER DIAMETER	25 mm
WAVELENGTH RANGE	280-730 nm *
FIBER OPTIC CONNECTION	Input: 1000 $\mu\text{m}$ Output: 600 $\mu\text{m}$
MATERIAL IN CONTACT WITH FILTER PAD	Delrin
WEIGHT	0.5 kg (1 lb)

\* Using a Tidas E Base spectrometer and D4H UV/VIS light source.

## References

M. Belz, K. Larsen, K.-F. Klein, "Fiber optic sample cells for polychromatic detection of dissolved and particulate matter in natural waters", *Proc. SPIE*, Vol. 6377, Oct 2006, 63770X

## ORDERING INFORMATION

**89575** QFT1, Fiber Optic Holder for Glass Fiber Filters

# In-Line Fiber Optic Filter Holder

*Insert optical filters into a fiber optic pathway*

## Features

- In-line filter for SMA terminated fibers
- Collimators for UV/VIS/NIR
- Filter diameter 8–25.4 mm
- Filter thickness 2–10 mm

## Benefits

- No stray light
- Solid design

## Applications

- Stray light filtering
- Removal of excitation light in fluorescence detection



56200

This In-Line Fiber Optic Filter Holder allows the insertion of optical filters within a fiber optic pathway. The connectors of the filter holder assembly are compatible with WPI's range of fiber optic jumper cables and can be coupled using SMA connectors.

Filters with outer diameters from 8 to 25.4 mm and thicknesses from 2 to 10 mm can be accommodated. The design limits lateral and axial movement of the filter when secured in the holder.

Two fiber optic collimators are internally mounted in the holder to pass collimated light through the filter and then refocus the filtered light into the aperture of the output fiber. Spectral range will be

largely limited by the bandpass of the optical fibers (from UV to near IR using WPI UV-enhanced cables).

## ORDERING INFORMATION

**56200** In-Line Fiber Optic Filter Holder (SMA)

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • www.wpiinc.com

Germany: +49 (0)6031 1602171 • wptide@wpi-europe.com • www.wpi-europe.com  
China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net

# Smallest Fiber Optic Dipping Probe

Perfect for UV/Vis spectroscopy



## Features

- Compact and efficient dipping probe
- Tip diameter 1.8 mm

## Benefits

- Connects with most standard spectrometers using a 600  $\mu\text{m}$  connection
- Perfect for mobile applications

## Applications

- Protein and DNA sample measurements
- Dissolution system

**Mini DipTip™** is a miniature transmission probe for microliter spectroscopic sampling. **Mini DipTip's** tip diameter is only 1.5 mm—the size of a 17-gauge needle. It will fit into all micro centrifuge tubes on the market. Microliter samples can be analyzed cost effectively when you combine the **Mini DipTip™** with one of the following:

- TIDAS E Base with **FO-6000** or **D4H**
- TIDAS S300 series
- Compatible with most fiber coupled spectrometers
- Ideal for multi channel applications with LEDSpec

**DIP-UV-MINI-10** 10 mm light pathlength

**DIP-UV-MINI-5** 5 mm light pathlength

**DIP-UV-MINI-2** 2 mm light pathlength

## DIPTIP SPECIFICATIONS

TIP DIAMETER	1.8 mm
LIGHT PATHLENGTH	2, 5, 10 mm
WAVELENGTH RANGE (nm)	200-1000
SAMPLE VOLUME REQUIRED	20-50 $\mu\text{L}$
DISTANCE FROM TIP TO UPPER EDGE OF SAMPLE WINDOW	7 mm
FIBER LENGTH	1.0 m
FIBER OPTIC CONNECTION	SMA 905
LAUNCH FIBER BUNDLE (7 x 200 $\mu\text{m}$ )	680 $\mu\text{m}^*$
RETURN FIBER BUNDLE (7 x 200 $\mu\text{m}$ )	680 $\mu\text{m}^*$


\*Circular packaging of the fiber bundle results in an active area equivalent to a fiber with a core diameter of 680  $\mu\text{m}$ . Using a 600  $\mu\text{m}$  connection is recommended and will result in negligible light loss.

## ORDERING INFORMATION

<b>DIP-UV-MINI-2</b>	Mini DipTip™ for UV/VIS/NIR (2 mm path)
<b>DIP-UV-MINI-5</b>	Mini DipTip™ for UV/VIS/NIR (5 mm path)
<b>DIP-UV-MINI-10</b>	Mini DipTip™ for UV/VIS/NIR (10 mm path)

# Fiber Optic Collimator

Collimate an emitted light beam or couple light into an optical fiber



## Features

- Maximum coupling efficiency with fused silica fibers
- Easily adjust the distance between the lens and the optical fiber for focusing
- Not suitable for single mode laser applications

## Benefits

- UV/VIS/NIR light collimation

## Applications

- Generates plan-parallel beam with a 5mm beam diameter from UV/ VIS light guided in a fiber

WPI's Fiber Optic Collimator can be used for both collimating a light beam emitted by an optical fiber or coupling light from a collimated light beam into an optical fiber. The numerical aperture of the collimator is optimized for maximum coupling efficiency into typical fused silica fibers. The collimator can, for example, be used to guide a parallel light beam through a sample cuvette or an optical filter with virtually no optical losses. In this application, one collimator collimates the light into a parallel beam 5 mm in diameter, enabling it to pass a long distance without losing the energy. After the light passes the sample media, a second collimator can be used to collect the beam into the receiving fiber. A unique design feature of this collimator is that the distance between the lens and the optical fiber can be easily adjusted. This permits it to be used as a focusing device or for fine-tuning the color balance when coupling light from a light source into multimode fibers.

## COLLIMATOR SPECIFICATIONS

LENS DIAMETER	5 mm
LENS FOCAL DISTANCE	10 mm
LENS MATERIAL	Ultraviolet grade synthetic fused silica (KU-1)
WAVELENGTH RANGE	220 nm-2 $\mu\text{m}$
MOUNTING THREADS	3/8-24 UNF
DIVERGENCE	< 0.1 rad for 1 mm core fiber
FIBER CONNECTOR INTERFACE	SMA or ST

## ORDERING INFORMATION

<b>300051</b>	Fiber Optic Collimator (SMA)
<b>300052</b>	Fiber Optic Collimator (ST)

## OPTIONAL ACCESSORIES/REPLACEMENT PARTS


<b>13395</b>	SMA Bulkhead feed through connector/coupler, D-hole
<b>13370</b>	SMA half-length bulkhead coupler/connector
<b>CC-3-UV</b>	Cosine Corrector



**CC-3-UV**




**13395**



**13370**

# Z-Dimensions Are Not Created Equal

Cuvettes come in a variety of shapes and sizes, but one of the most important specifications of a cuvette is its Z-dimension. The Z-dimension of an instrument (cuvette holder or spectrometer) is the distance from the bottom of the cuvette chamber floor to the center of its light beam (see image). A cuvette's Z-dimension must match the Z-dimension of the instrument with which it will be used.



Each manufacturer designs its instruments

with a specific Z-dimension. Common Z-dimensions include 8.5 and 15 mm, and sometimes 20 mm. When purchasing small volume cuvettes, the correct Z-dimension becomes critical. Matching the Z-dimension of the cuvette to the Z-dimension of the instrument ensures that the light beam passes through the center of small samples. The table below shows the standard Z-dimension of the spectrometer sample compartments for many manufacturers.

Manufacturer	Z-Dimension
Agilent®	15 mm
Avantes®	15 mm
Beckman®	8.5 mm
Bio-Rad®	8.5 mm
Cecil®	15 mm
Eppendorf®	8.5 mm
Hewlett – Packard®	15 mm
Hitachi®	8.5 mm
Jasco®	11 mm
J & M®	8.5 mm
Ocean Optics®	15 mm
Perkin – Elmer®	15 mm
Pharmacia®	15 mm
Shimadzu®	15 mm
Spectronics®	8.5 mm
Stellarnet®	15 mm
Turner®	8.5 mm
Varian®	20 mm
WPI	15 mm

To determine the Z-dimension of a cuvette holder:

- Use strips of heavy paper that will fit neatly into a cuvette (for example, 12 mm x 50 mm) and not allow light to pass through the cuvette.
- Poke a tiny hole in each paper “sample.” For example, one paper sample could have a hole at 8.5 mm, one at 15 mm, one at 20 mm.
- One at a time, insert the paper samples into the cuvette and place the cuvette into the cuvette holder. The paper sample with the pin hole at the instrument’s Z-dimension will allow light to pass. The other paper samples will not allow light to pass.

If you have an instrument that is not on the list and need to know its Z-dimension, please contact WPI at (941) 371-1003 or wpi@wpiinc.com.

# Optical Glass and Quartz Cuvettes

*For spectrophotometry and fluorometry*

## Features

- Absorbance and fluorescence cuvettes
- 2–100 mm path length
- Fused silica and glass cuvettes
- Flow cuvettes
- Microliter cuvettes

## Benefits


- Eight different styles
- Inexpensive single use glass cuvettes

## Applications


- Absorbance spectrophotometers
- Fluorescence fluorometers
- Chemistry, Biochemistry
- FIA, GFIA, Quality control

WPI's glass and synthetic quartz cuvettes are ideal for UV/VIS/NIR absorbance or fluorescence experiments.

Synthetic quartz can be used in deep UV applications and is recommended for fluorescent applications, as it does not exhibit background fluorescence. Quartz cuvettes (absorbance, fluorescence and flow) are shipped individually packaged, glass cuvettes are shipped in packages of 10 cuvettes. These economic quartz and glass cuvettes are ideal for precision measurements because of their high quality materials used and their low manufacturing tolerances. Typical transmission curves of glass and synthetic quartz cuvettes are shown below.



*Transmission curves of Glass and Synthetic Quartz Cuvettes. The cuvettes were empty, thickness 1.25 mm x 2, including surface reflections, measured with a TIDAS II against air as reference.*



*A complete transmission spectrum of glass and synthetic quartz cuvettes from 190 nm to 4 mm is shown. Cuvettes were empty, thickness 1.25 mm x 2, including surface reflections, measured with a TIDAS II.*

## SPECIFICATIONS

Cuvette Material	Spectral Range (>80%)	Transmission Difference Between Different Cuvettes
OPTICAL GLASS	350 – 2500 nm	Less than 1%
SYNTHETIC QUARTZ	200 – 2500 nm	Less than 1%

Style G



ORDERING INFORMATION

WPI PN	Style	material	Polished windows	path [mm]	Dimensions [mm]	volume [mL]	Beam width [mm]
<b>STANDARD RECTANGULAR CUVETTES</b>							
CUV2101-1*	B	Quartz	2	1	3.5x12.5x45	0.35	10
CUV2102-1*	B	Quartz	2	2	4.5x12.5x45	0.7	10
CUV2011-1*	B	Quartz	2	5	7.5x12.5x45	1.7	10
CUV1022-10	C	Optical Glass	2	10	12.5x12.5x45	3.5	10 <b>pack of 10</b>
CUV2012-1	C	Quartz	2	10	12.5x12.5x45	3.5	10
CUV2105-1	C	Quartz	2	20	22.5x12.5x45	7	10
CUV2106-1	C	Quartz	2	30	32.5x12.5x45	10.5	10
CUV2107-1	C	Quartz	2	40	42.5x12.5x45	14	10
CUV2108-1	C	Quartz	2	50	52.5x12.5x45	17.5	10
*89341	Cuvette spacer for 1-mm cuvettes (part CUV2101-1)						
*89342	Cuvette spacer for 2-mm cuvettes (part CUV2102-1)						
*89337	Cuvette spacer for 5-mm cuvettes (part CUV2011-1, CUV2023-1, CUV2063-1)						
<b>SELF MASKING SEMI MICRO CELL CUVETTE</b>							
CUV2023-1*	D	Quartz	2	5	7.5x12.5x45	0.7	4
CUV2031-1	D	Quartz	2	10	12.5x12.5x45	1.4	4
CUV2025-1	D	Quartz	2	20	22.5x12.5x45	2.8	4
CUV2032-1	D	Quartz	2	10	12.5x12.5x45	1	3
CUV2033-1	D	Quartz	2	10	12.5x12.5x45	0.7	2
CUV2034-1	D	Quartz	2	10	12.5x12.5x45	0.35	1
<b>SELF MASKING CONTINUOUS FLOWTHROUGH CELL</b>							
CUV2063-1*	E	Quartz	2	5	7.4x12.5x45	0.035	∅ 3
CUV2061-1	E	Quartz	2	10	12.5x12.5x45	0.07	∅ 3
CUV2065-1	E	Quartz	2	20	22.6x12.5x45	0.14	∅ 3
CUV2066-1	E	Quartz	2	30	32.6x12.4x45	0.21	∅ 3
CUV2062-1	F	Quartz	2	10	12.5x12.5x45	0.48	4x12
<b>SELF MASKING CONTINUOUS FLOW THROUGH CELL, SMALL INPUT, LARGE OUTPUT Z=8.5 MM</b>							
CUV2614-1	H	Quartz	2	10	12.4x12.4x35.6	0.03	∅ 2
<b>MICRO CELL WITH BLACK WALLS</b>							
CUV2674-1	J	Quartz	2	10	12.5x12.5x45	0.05	2
<b>FLUORESCENCE</b>							
CUV2051-1	A	Quartz	4	10	12.5x12.5x45	3.5	10
CUV2052-1	A	Quartz	4	10	12.5x12.5x45	1.4	4
<b>LONG PATH CUVETTE</b>							
CUV2071-1	G	Quartz	2	100	102.5 x 22 ∅	28	19



# Photodiode Array Based Spectrophotometers

*High performance fiber optic spectrometer systems*

CE



504719



## Features

**Plug and Play** – with help of standard accessories the TIDAS® E sits ready to be immediately operational

**Increased Flexibility** – the patented cuvette holder allows a variety of measurement configurations that are seamlessly integrated within one instrument

**Diode Array Technology** – enables quick and precise measuring with a low level of scattered light

**Software** – modular and user friendly, permits the data acquisition and evaluation of individual tasks

- Fiber-optic diode-array spectrometer
- Measures in transmission and absorption with external light source
- Supports external immersion probes, temperature-controlled cuvette holders and flow cells
- Supports external reflection measurement heads
- Performs film thickness measurement using white light interference
- Suitable for color measurement and analysis
- Chemometric modelling
- Certified instrument quality (wavelength and photometric accuracy)
- Software for instrument control, data acquisition and data evaluation included

## Benefits

- Compact, rugged and transportable instruments
- Multimodal and combined methods
- High dynamic range
- Fast and reliable results
- Powerful, yet easy to use software

## Applications

- Standard chemistry and biochemistry
- Trace level nutrient analysis
- Detection in flow injection analysis (FIA)
- HPLC analysis with WPI's **LWCC-M** flow cells
- CDOM detection with **UltraPath**

The **TIDAS** spectrometer and spectrophotometer series is based on a high performance monolithic spectrometer module manufactured by Zeiss for UV, UV/VIS and VIS/NIR applications. The **Tidas E Base** spectrometers are ideal for modular spectrometer systems with separate detectors, light sources and sample cells. The **TIDAS S300** spectrophotometers already include a light source.


The **TIDAS** series instruments are specialized for research and instructional concepts. This Diode Array based technology allows for fast and precise measurements, and the instruments are universally applicable. The **TIDAS E Base** and the **TIDAS S300** are suitable for routine measurements in the laboratory and outperform conventional bench-based spectrophotometers and CCD-based spectrometer modules, when it comes to high precision fiber optic sampling. They rely on a monolithic optical bench made by Zeiss, which is optimized for fiber optic applications. Most cuvette-based standard spectrometers lose more than 90% of light through expensive prism decoupling. The Tidas systems are designed for fiber optic sampling cells. Using suitable light sources and sample cells, spectral detection in the wavelength range of 190 – 720 nm can be performed.

## Applications

The Tidas systems are ideally suited for WPI's fiber optic sampling equipment. High sensitivity detection systems for flow analysis can be assembled using WPI's Liquid Waveguide Capillary Cells (**LWCC**) with effective pathlengths ranging from 50 to 500 cm. These setups are frequently used in fluid injection analysis systems for nutrient analysis (nitrite, nitrate, phosphate, iron) in oceanographic applications. Microliter sampling systems for UV/VIS applications can be assembled using WPI's **DipTip™** dipping probes.

## Software

**TidasDAQ 3** software is included with each instrument for data collection and data analysis. TidasDAQ is used to run the spectrometer module, collect spectra in either single or continuous mode, control the digital I/Os, save the experimental data to disk, and analyze the data. Further, TidasDAQ can export data directly into GRAMS/AI, a feature very useful for advanced data analysis for pharmaceutical applications and requirements.




TIDASDAQ acquisition window, showing an absorbance baseline.

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • www.wpiinc.com

Germany: +49 (0)6031 1602171 • wpide@wpi-europe.com • www.wpi-europe.com  
China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net





Spectra may be displayed and analyzed in 2D and 3D format. This allows the user to conveniently interpret "time acquisition" data typically done with a TIDAS-E-BASE-LWCC flow system.

### TidasDAQ: Data Collection & Instrument Control


With TidasDAQ, high precision intensity, absorbance, transmittance or normalized spectra can be obtained in less than a second. Only a few parameters need to be adjusted to obtain spectral data. Sampling of single scans, continuous full spectra scans or triggered scans is possible. Chromatograms can be displayed and logged to disk at up to four wavelengths. Data Export of 2D and 3D Spectrograms, as well as Chromatograms is supported in ASCII, Spectralys/SpectraView, Excel and Grams/AI formats. Light sources and other sampling instrumentation can be controlled via the TTL level digital outputs, and data collection can be triggered by TTL leveled external inputs of the Tidas E Base.

Spectra can be recorded in 2D and 3D view. Mathematical computation, Derivation, Smoothing, Quantification and other functions are available to work with your data. The Quantification module allows single point and multiple point analysis, multiple linear regression, partial least square and principle component analysis. Data can be exported out of a 3D analysis file into separate scans. Furthermore, chromatograms as well as spectrograms can be copied directly into Excel for further data analysis.

### Typical Usage


These high performance spectrometer modules and spectrophotometers combine modularity, flexibility, ease of use and high quality in one unit. They are ideally suited for WPI's fiber optic sampling equipment. High sensitivity detection system for flow injection analysis (FIA) and gas segmented flow injection analysis (GFIA) systems can be assembled with WPI's proprietary Liquid Waveguide Capillary (LWCC) flow cells with effective pathlengths ranging from 50 to 500 cm. Typical applications include colorimetric trace detection of nitrite, nitrate, phosphate and iron in oceanographic research.

Small flow volume samples, like those that may be found in Optofluidics and High Performance Liquid Chromatography (HPLC), can be measured conveniently and accurately with WPI's LWCC-M series of flowcells. LWCC-M flow cells feature pathlengths of 10, 50 and 100 mm with corresponding sample cell volumes of 2, 4, 12 and 24  $\mu\text{L}$ . Finally, small discrete samples of 2-40  $\mu\text{L}$  volume often found in biochemistry (like, Protein or DNA) may be analyzed using WPI's DIP-UV-MINI miniature dipping probes with 2, 5 and 10 mm pathlength. Other applications for



Mean measured absorption spectra for marine samples, measured with TIDAS1 and a 100-cm pathlength LWCC (adapted from Beine et al., J Geophys Res, Vol. 117, D00R15).

the TIDAS series spectrometers and spectrophotometers include material science (forensics, semiconductor technology), quality control (automotive and food), as well as new research in the field of nanotechnology and rapid kinetics.



Absorbance spectrum displayed using TIDASDAQ software

### SPECIFICATIONS

	TIDAS S300/E BASE UV	TIDAS S300/E BASE UV/VIS	TIDAS S300/E BASE VIS/NIR
WAVELENGTH RANGE	190-390 nm	190-720 nm	300-1100 nm
SPECTRAL RESOLUTION	<3 nm	<7 nm	103 nm
WAVELENGTH ACCURACY	$\pm 1$ nm	$\pm 1$ nm	<3 nm
PHOTOMETRIC ACCURACY UV	$\pm 10$ mAU	$\pm 10$ mAU	-
PHOTOMETRIC ACCURACY VIS	-	$\pm 10$ mAU	-
WAVELENGTH REPRODUCIBILITY	<0.1 nm	<0.1 nm	<0.1 nm
BASELINE DRIFT @250NM	$5 \times 10^{-4}$ AU/h according to ASTM E685	$5 \times 10^{-4}$ AU/h according to ASTM E685	$5 \times 10^{-4}$ AU/h according to ASTM E685
SIGNAL-TO-NOISE RATIO	$<3 \times 10^{-5}$ AU according to ASTM E685	$<3 \times 10^{-5}$ AU according to ASTM E685	$<3 \times 10^{-5}$ AU according to ASTM E685
NUMBER OF DIODES	256	256	256
LIGHT SOURCE	Deuterium	Deuterium & Halogen	Halogen
	<i>Not included for Tidas E Base</i>		
OPTICAL BENCH	Monolithic spectrometer module with concave aberration corrected holographic grating;		
DETECTOR ARRAY	Hamamatsu photodiode array, 256 pixel		
A/D RESOLUTION	16 Bit		
FIBER OPTIC CONNECTION	SMA905, 600 $\mu\text{m}$		
SOFTWARE	Tidas DAQ3 (included)		
SCRIPT LANGUAGE FOR METHOD DEVELOPMENT	Yes		
SYSTEM REQUIREMENT	Windows 7, 8, 10		

### ORDERING INFORMATION

<b>504717</b>	TIDAS E Base, UV 190-390 nm
<b>504718</b>	TIDAS E Base, UV/VIS 190-720 nm
<b>504719</b>	TIDAS E Base, VIS/NIR 300-1100 nm
<b>505066</b>	TIDAS S300, UV 190-390 nm with deuterium lamp
<b>505067</b>	TIDAS S300, UV/VIS 190-720 nm with deuterium/halogen lamp
<b>505068</b>	TIDAS S300, VIS/NIR 300-1100 nm with halogen lamp
<b>505069</b>	Tidas S300 VIS/NIR 300-1100 nm (without halogen lamp)

*Systems includes power supply, TIDAS DAQ software, RJ-45 cable and manual.*

# Deuterium Halogen Fiber Light Source

For a continuous spectrum in the UV, VIS and NIR range

## Features

- Continuous spectrum from 200 nm–1700 nm
- Integrated shutter with switch and TTL control
- SMA fiber optic connection
- Separate UV and VIS bulb control

## Benefits

- Matched deuterium and halogen bulbs
- Matched optical output for WPI's flow cells and fiber optic probes
- Low drift < 1 mAU/h @ 254 nm
- Simple exchange of light bulbs

## Applications

- UV/VIS/NIR absorbance spectroscopy
- Excitation light source for fluorescence applications
- Colored Dissolved Organic Matter (CDOM) detection in seawater
- Nutrient analysis in fresh water, drinking water and seawater using WPI's **LWCCS**
- Protein detection using WPI's **DIP UV MINI** fiber optic probes

The **D4H** is a combined deuterium and halogen light source for UV/VIS and NIR applications. This light source is ideally suited to work with WPI's spectrometer modules and sample cells. It supplies a continuous spectrum in the UV, VIS and NIR range from 200 nm to 1100 nm. The **D4H** is equipped with an integrated electrical shutter, which can be controlled by a switch or a TTL signal.



D4H




Replacement Deuterium Lamp  
**503847**



Replacement Halogen Lamp  
**503848**

### D4H Spectrum



## LIGHT SOURCE SPECIFICATIONS

	D4H	FO-6000
APPLICATION	UV/VIS/NIR	VIS/NIR
SPECTRAL RANGE	200–1100 nm	380–1700 nm
DEUTERIUM LAMP LIFE	2000 hr	NA
TUNGSTEN/HALOGEN LAMP LIFE	2000 hr	3000* hr
STABILITY	1-2 mAU/h	<0.5 mAU/h
POWER CONSUMPTION	140 W	6 W
POWER REQUIREMENTS	110/240V, 50-60 Hz, 1A	12VDC/1A
SHUTTER/TTL TRIGGER	Yes	Yes
MAX. FIBER OUTPUT	1000 μm	1000 μm
CONNECTIONS	SMA	SMA
SHIPPING WEIGHT	13.2 lb (6 kg)	1.3 lb (0.6 kg)
DIMENSIONS (W/H/L)	7 x 6.2 x 9.8 in. (17.8 x 15.7 x 25 cm)	4.8 x 2.8 x 7.5 in. (12 x 7 x 19 cm)

\*Lamp life is dependent upon internal power settings.

## ORDERING INFORMATION

<b>D4H</b>	Deuterium Halogen Light Source (200 nm-1100 nm)
<b>503848</b>	Halogen Replacement Lamp for D4H
<b>503847</b>	Deuterium Replacement Lamp for D4H (> 215 nm)

## OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>TIDAS-D2</b>	Replacement Deuterium Lamp, for Tidas II
<b>TIDAS-H</b>	Replacement Halogen Lamp (Type 1), for Tidas II
<b>TIDAS-H2</b>	Replacement Halogen Lamp (Type 2), for Tidas II
<b>D2H-DB</b>	Replacement Deuterium Lamp, for D2H
<b>D2H-HB</b>	Replacement Halogen Lamp, for D2H
<b>D2H-HBER</b>	Replacement Deuterium Lamp, Extended Range, for D2H

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • www.wpiinc.com

Germany: +49 (0)6031 1602171 • wpide@wpi-europe.com • www.wpi-europe.com  
China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net

# Tungsten Fiber Light Source

High color temperature



## Features

- Visible light source
- 380 nm – 1700 nm
- Low drift < 0.5 mAU/h
- SMA fiber connection
- Electrical shutter control (switch & TTL)


## Benefits

- Bulb with 10,000-hour lifetime
- Fits all of WPI's fiber optic probes and flow cells
- Design matched to **LOW DRIFT** applications
- Temperature controlled optical bench

## Applications

- Low noise VIS absorbance measurements
- Trace analysis of nutrients in seawater and freshwater
- CDOM analysis in seawater
- Analytical chemistry, environmental science and life science

The **FO-6000** is a continuous fiber optic light source featuring an extended visible part of the light source (380 nm - 1700 nm). It has a SMA fiber optic connector. The shutter and lamp can both be controlled via a switch or external TTL triggering. This light source offers a wide assortment of applications. A special feature of the **FO-6000** is its color balancing optics, which shifts the usable range of the light source from traditional 420 nm down to 380 nm wavelength. Due to its thermally controlled optical bench, it is particularly suitable for low noise and low drift applications.



## ORDERING INFORMATION

<b>FO-6000</b>	Fiber Optic Light Source
<b>800120</b>	Replacement Lamp for FO-6000

# Absorbance Detection

Detection of organic compounds in water analysis

Absorption of light correlates to the energy of a photon that is taken up by electrons of the substance atom. The electromagnetic energy is transformed into internal energy of the absorbent substance. The absorbance of a substance quantifies how much of the incident light is absorbed by it (instead of being reflected or refracted). Precise measurements of the absorbance at many wavelengths allow the identification of a substance via absorption spectroscopy, where a sample is illuminated from one side, and the intensity of the light that exits from the sample in every direction is measured (see Fig. 1). A few examples of absorption are ultraviolet-visible (UV-Vis) spectroscopy or infrared (IR) spectroscopy.




Fig 1. Concept of absorbance spectroscopy using white light and optical components to filter out light of a specific wavelength that interacts with molecules in the solution. Absorbance at this specific wavelength by the molecules in the solutions is detected as a decrease in light intensity (Spectrophotometer-Source: <http://chemwiki.ucdavis.edu/>).

Absorption is the amount of light that a substance takes in and does not allow to pass through it. Spectrophotometers actually measure transmission, the amount of light that passes through a sample, but this is converted into absorption by comparing the bulb output to the light that has passed through the sample. Light sources that can be used for absorbance spectroscopy depend strongly on the used substance to label a specific molecule and can span the entire electromagnetic spectrum of light.

Protein detection uses the UV- spectrum (typically, 260 nm and 280 nm, while further information is obtained at 230 nm and 320 nm, but compensated by a selection of ratios and background corrections) and is most commonly used to estimate DNA or RNA concentration and to analyze the purity of the preparation. Further application spans the measurement of the light scatter at 600 nm to monitor the growth rate of a cultured bacterial population and to identify the peak concentration.


Industrial applications that also cover the UV spectrum for protein detection is linked to food analysis (for example, characterizing the grading of olive oil as extra virgin, virgin oil, etc., as set out by European regulations) or quality control in the pharmaceutical industry. In addition, industrial applications of absorbance spectrophotometry cover the characterization of water purity or waste water analysis, in addition to the determination of specific organic molecules like nitrate, nitrite or phosphor. This last application is interesting when using the visible light (Vis) spectrum.

# Biofluorometer

## Reliable, 2-channel, LED-based fluorometer

### Features

- Two photomultiplier inputs
- Light excitation with high power LEDs
- Modes: single excitation & single emission, single excitation & dual emission, dual excitation & single emission
- 2 channel mode (two single excitation & single emission)
- Optical connections: Liquid Light Guides (LLGs) and SMA terminated fibers
- Sampling rates up to 1 kHz (1000 ratios/second)
- Automatic LED light drift correction for long term measurements
- Automatic room light correction
- Optional fiber optic probes for horizontal tissue bath applications
- Optional imaging probes for Langendorff systems
- Optional attachments for direct connection to fluorescence microscopes via epifluorescence port (excitation) and C-Mount (Emission) via liquid light guides



The Biofluorometer can be connected to the epifluorescence port of an inverted microscope and its high intensity LED light source is used for excitation of the corresponding dye.

### Benefits

- Versatile instrument for horizontal tissue bath, Langendorff and microscope applications
- Warmup time < 1 minute
- Low bleaching mode (5% LED On time)
- Can be combined with imaging based Sarcomere Detection System (OptiSarc)
- SMA and Liquid Light Guide (LLG) may be exchanged (optional)
- Customized analysis techniques in WPI's MDAC data acquisition software

### Applications

- Auto-fluorescence imaging

The **SI-BF-100** is an LED-based fluorometer for life science applications. It is ideally suited for ratiometric calcium detection (FURA-8) and ATPase detection (via NADH fluorescence). With up to three LED modules (wavelengths), the **SI-BF-100** covers many fluorometric applications in neuroscience and cell biology.


The **SI-BF-100** enables the detection and analysis of fluorescence signals in four different modes:

- Single excitation/single emission–In this classical mode, a fluorophore is excited at one wavelength and the fluorescence signal is detected at a single higher wavelength using one photomultiplier. The concentration of the analyte is directly proportional to the intensity of the detected signal.
- Dual excitation/single emission–A fluorophore is excited at two wavelengths and the fluorescence signal is detected at one wavelength using one photomultiplier. The concentration of the analyte is proportional to the ratio of the two detected fluorescence signals. This ratiometric concept minimizes the effect of indicator dye bleaching and motion artifact in experiments. A typical example is the detection of free calcium in muscle tissue using the indicator dye Fura-8™.
- Single excitation/dual emission–A fluorophore is excited at one wavelength and the fluorescence signal is detected at two wavelength using two photomultipliers.
- Dual excitation/dual emission–Two separate fluorophores are excited at different wavelengths and the fluorescence signal of each fluorophore is detected at two separate wavelengths using two photomultipliers.

### Fluorescence Probe

The principle of a fluorescence probe is to bring light to a sample and to detect reflected or fluorescent light from the sample surface. These probes may be used for the detection of the transient response of free ion concentrations, like calcium, potassium, sodium, magnesium. Fluorescence probes can also be designed to detect pH and membrane potential. Auto-fluorescence like the detection of ATPase activity via NADH or FAD is another application. If a fluorescent dye is used to label a molecule, a custom fluorescence probe can be designed to detect the emission.

To select the appropriate fluorescent probe, consider the detection principle and the dye to be used. All probes use fibers with 300 μm core diameter. Excitation fibers have 1000 μm SMA connectors for excitation and double emission probes, and 1500 μm SMA connectors for excitation and single emission probes. The standard probes come with a rectangular head (3 mm × 0.7 mm) or round head (2 mm diameter). The fibers that are used for such probes are highly flexible plastic fibers or fused silica fibers. Custom fibers for your application may be designed with different core fiber diameters.



(Left) Customized single emission round probe with a 5 mm diameter for one excitation and one emission wavelength.


(Right) Customized rectangular probe.

Application	Probe	Excitation	Emission
Free Calcium via Indo-1	Double Emission	365nm	405nm/486nm
Free Calcium via Fura-8	Single Emission	365nm/410nm	525nm
ATPase via NADH/TAMRA	Double Emission	365nm	472/572nm
FAD	Single Emission	470nm	525nm

## Auto-fluorescence imaging


### Genetically modified frogs

Fluorescence signaling with single emission probe on living frogs of wildtype tissue (wt) and genetically modified tissue from frog expressing a protein of "auto-fluorescence" called EGFP. Data below showing an increase in autofluorescence signal for EGFP frogs.



**A** Wildtype (wt, left) and EGFP (right) frog strain, showing the difference in color between both strains. Photo courtesy of Dr. Xu-Feng, Qi Jinan University.

**B** Experimental setup for frog dorsum *in vivo* fluorescence measurement with single emission probe. Observe the LED beam/spot as measurement area of frog dorsum. Photo courtesy of Dr. Xu-Feng, Qi Jinan University.




Fluorescence response of frog dorsum *in vivo* with single emission probe from whole frog dorsum measurements for EGFP (red) vs. wildtype (wt, green). Note the distinguishable increase in auto-fluorescence signal for EGFP frog dorsum but confounded with tissue color backreflection

### Genetically modified mice isolated hearts

Fluorescence signaling with LLG Imaging probe of wildtype tissue (wt) and genetically modified tissue from mice isolated heart expressing a protein of "auto-fluorescence" called EGFP. Data below showing an increase in autofluorescence signal for EGFP heart at two LED current intensity.



Experimental setup for isolated mouse heart fluorescence measurement with LLG Imaging Probe.




Fluorescence signal from isolated mouse heart with LLG Imaging Probe for EGFP (red) vs. wildtype (wt, green). Note that auto-fluorescence signals can be distinguished between EGFP vs. wildtype, although the signals and difference are weak. Note that back reflection can be considered as constant, as heart and dish support colors are the same. Photo courtesy of Dr. Xu-Feng, Qi Jinan University.

## Dye loading concept for measuring fluorescence signal on a leaf

Often calcium fluorescence is measured in leaves via specific dye loading techniques, like microinjection. A simplified concept for dye loading can also be the use of the root fluid transport system. In this example Fluorescein dye was used to exemplify this technique. Measurements of Fluorescein in the leaf was done with the BF-100 and an SMA single emission probe for direct sensing under a dark condition. Backreflection is considered as constant as booth leaves had the same coloring.

For the setup, the plant with root was incubated in Fluorescein solution (using physiological transportation) for ~1hr or 2~hrs in Fluorescein concentration of 16 mg/l and 100 mg/l. Measurements were taken on the side without vein. Clear distinguishable effects of the different Fluorescein concentration can be observed, notably regarding the control measurement with Milli-Q water.

Example of plant dye loading via the plant root for two different Fluorescein concentration. Photo courtesy of Professor Yao Youli of Yangzhou University, China.



Plant dye loading of fluorescein at two different concentrations (green, red) via the plant root. For comparison see only Milli-Q water (blue). Note that dye loading time was 1 hr (green) and 2 hrs (red). Fluorescence measurement was conducted by shining directly on the leaf under a dark condition.

### SI-BF-100 SPECIFICATIONS

MEASUREMENT PRINCIPLE	Fiber optic fluorometer with 2 inputs and 1 output
DETECTOR INPUTS	2 PMTs
EXCITATION	High Power LED Modules: 365 nm, 420 nm, 470 nm, 530 nm (select any 3 modules, when ordering)
ANALOG OUTPUT RANGE	0-10V
OPTICAL CONNECTIONS	Choice of Liquid Light Guide (LLG) or SMA connections
POWER	12 V/2 A (includes external 100 – 240 V / 50 – 60 Hz power supply)
DIMENSIONS (h x w x d)	3.5 x 17 x 13 in. (88 x 431 x 330 mm)

### ORDERING INFORMATION

<b>SI-BF-100LLG</b>	Biofluorometer with LLG Optical Connections
<b>SI-BF-100SMA</b>	Biofluorometer with SMA Fiber Optic Connections
<b>99261</b>	C-Mount Microscope Attachment for 2x PMTs
<b>99259</b>	C-Mount microscope attachment for 1x camera & 2x PMTs (includes 1x camera C-Mount adapter with adjustable aperture,

### OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>SI-BF-SMA-UPGRADE</b>	Biofluorometer Upgrade Kit for SMA optical probes
<b>SI-BF-LLG-UPGRADE</b>	Biofluorometer Upgrade Kit for LLG connections
<b>802407</b>	Liquid Light Guide (LLG), 3mm diameter, 6' long
<b>M3301</b>	Manual Manipulator for securing the probe
<b>M10</b>	Magnetic Base
<b>94650</b>	Single Emission, Small Tissue Probe
<b>94689</b>	Dual Emission, Small Tissue Probe

Check our website for new LED modules, emission filters and dichroic mirrors for specific applications.

# Precision Stereo Zoom Microscope

*Modular, two parallel beam path design with high quality optics*

## Features

- Planachromatic objectives, no optical distortion
- Large zoom ratio: 8:1
- Large zoom range: 0.62x-5.0x
- Step and continuous zoom

## Benefits

- True trinocular version available
- High-contrast imaging, ideal for observing transparent, low-contrast objects
- Long working distance option available

## Applications

- Integrated optical and biological research

The fourth generation of WPI's precision stereo zoom microscope uses modular, two parallel beam path design and high quality optical system. The advanced optical design with planachromatic objectives provides a sharp and distortion-free contrast image throughout the entire zoom range and comes with an impressive 5-year warranty.

The **PZMIV** is available in a binocular or trinocular version. In addition, an extensive list of optional accessories is available that makes the PZMIV suitable for integrated optical and biological research.

The microscope comes with a track stand, standard 10x eyepieces (wide-field, distortion-free and high eye point) and a 1x planachromatic distortion-free objective.

All **PZMIV** and **PZMTIV** microscopes come with 10x eyepieces, built-in 1x auxiliary lens and light ring adapter.

### Long working distance option

A common application for the PZMIV and the **PZMIII** is to add a 0.5X objective to allow for more working distance under the lens. This increases the PZMIV working distance from 80 mm to 189 mm. The PZMIII working distance increases from 100 mm to 177mm.


The use of this auxiliary lens drops the magnification range to half, so we also recommend acquiring the 20X eyepieces to return the microscope to its standard magnifications.

See the specification tables marked areas for more details. On trinocular versions of the microscope with the LWD option, the camera view is nearly the same as the eyepiece visual field of view.

**USBCAM50** digital camera included with version shown, also available separately

**Z-LITE** fiber optic illuminator now included with systems below differs from unit in photo — See page 191.

Computer and display not included



**PZMTIV-BS-LWD-DIG50**

*The trinocular version is a true trinocular, with continuous operation of both eyepieces and photo tube simultaneously. There is no need to block the right eyepiece to use the photoport.*

**504134 LED Ringlight**



See page 191

## PZMIV & PZMTIV Eyepieces and Objectives

Objective	10x Eyepiece		16x Eyepiece		20x Eyepiece		25x Eyepiece		Working Distance
	Mag	Field (mm) (Video Field)	Mag	Field (mm) (Video Field)	Mag	Field (mm) (Video Field)	Mag	Field (mm) (Video Field)	
0.32x	1.9x - 16x	106 - 13.1 (49.8 - 6.1)	3.2x 25.6x	70.6 - 8.8 (49.8 - 6.1)	3.9x - 32x	55.4 - 6.9 (49.8 - 6.1)	5x - 40x	45.4 - 5.6 (49.8 - 6.1)	296 mm
0.5x	3.1x - 25x	67.7 - 8.4 (31.8 - 3.95)	5x - 40x	45.2 - 5.6 (31.8 - 3.95)	6.2x - 50x*	35.5 - 4.4 (31.8 - 3.95)	7.8x - 62.5x	29 - 3.6 (31.8 - 3.95)	189 mm*
0.63x	3.9x - 31.5x	53.8 - 6.7 (25.3 - 3.15)	6.2x - 50.4x	35.8 - 4.4 (25.3 - 3.15)	7.8x - 63x	28.2 - 3.5 (25.3 - 3.15)	9.8x - 78.8x	23 - 2.9 (25.3 - 3.15)	149 mm
1.0x (inc)	6.2x - 50x	33.9 - 4.2 (15.9 - 1.97)	9.9x - 90x	22.6 - 2.8 (15.9 - 1.97)	12.4x - 100x	17.7 - 2.2 (15.9 - 1.97)	15.5x - 125x	14.5 - 1.8 (15.9 - 1.97)	80 mm

*The Video Field is based on a 1/2-inch CCD (8 mm diagonal) and a 0.5x camera adapter.*

*\*Long working distance (LWD) configuration*

### WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
 Germany: +49 (0)6031 1602171 • wvide@wpi-europe.com • www.wpi-europe.com  
 Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • www.wpiinc.com  
 China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net



**ORDERING INFORMATION**

<b>PZMIV</b>	Precision Stereo Zoom Binocular Microscope (Model IV), on Track Stand
<b>PZMIV-BS</b>	PZMIV Microscope on Boom Stand
<b>PZMTIV</b>	Precision Stereo Zoom Trinocular Microscope (Model IV), on Track Stand
<b>PZMTIV-BS</b>	PZMTIV Microscope on Boom Stand
<b>* PZMTIV-DIG50</b>	PZMTIV Microscope System, including PZMTIV, USBCAM50 USB Computer Camera, 0.5× CCD Camera Coupler, Z-LITE Optical Illuminator, Bifurcated Optical Fiber Light Guide
<b>* PZMTIV-BS-DIG50</b>	PZMTIV Microscope System (80 mm WD), including PZMTIV, USBCAM50 USB Computer Camera, 0.5× CCD Camera Coupler, Z-LITE Optical Illuminator, Bifurcated Optical Fiber Light Guide, Boom Stand
<b>* PZMTIV-BS-LWD-DIG50</b>	PZMTIV Microscope System (189 mm WD), including PZMTIV, USBCAM50 USB Computer Camera, 0.5× CCD Camera Coupler, Z-LITE Optical Illuminator, Bifurcated Optical Fiber Light Guide, Boom Stand, 0.5× Objective, 20× Eyepieces for Long Working Distance Viewing

See website for complete configurations.  
 \* Z-LITE-Z may be substituted in EU countries.

**OPTIONAL ACCESSORIES/REPLACEMENT PARTS**

<b>Optional Bodies</b>	
<b>502000</b>	PZMIV Binocular Body, 10× Eyepieces, 1× Objective, Eye Guards
<b>502001</b>	PZMTIV Trinocular Body, 10× Eyepieces, 1× Objective, Eye Guards
<b>Optional Stands, Mounts, Bases</b>	
<b>502002</b>	76mm Rectangular Track Stand
<b>502004</b>	Boom Stand (Heavy) without Focus Mount (requires 502009 Focus Mount for PZMIV)
<b>502005</b>	Ball Bearing Boom Stand (Heavy) without Focus Mount (requires 502009 Focus Mount for PZMIV)
<b>502006</b>	Boom Clamp Stand (Heavy) (requires 502009 Focus Mount for PZMIV)
<b>502009</b>	Universal Focus Mount for 76 mm PZMIV (Required for BS, AAC, BBS, and BCS) (5/8" pin)
<b>504123</b>	Extension for Heavy Clamp Stand
<b>504929</b>	LED Microscope Stand, 12.5-in. Post
<b>Optional Eyepieces</b>	
<b>502010</b>	10× Wide Field Eyepiece for PZMIV (pair)
<b>502011</b>	16× Wide Field Eyepiece for PZMIV (pair)
<b>502012*</b>	20× Wide Field Eyepiece for PZMIV (pair)
<b>502013</b>	25× Wide Field Eyepiece for PZMIV (pair)
<b>500264</b>	10× Eyepiece with Reticle (matches 10× eyepiece #502010)
<b>500266</b>	20× Eyepiece with Reticle (matches 20× eyepiece #502012)
<b>Optional Objectives</b>	
<b>502016</b>	0.32×, Planachromatic Objective (Distortion-free) (278 mm WD)
<b>502017*</b>	0.50×, Planachromatic Objective (Distortion-free) (174 mm WD)
<b>502018</b>	0.63×, Planachromatic Objective (Distortion-free) (138 mm WD)
<b>502019</b>	1.0×, Planachromatic Objective (Distortion-free) (73mm WD)
<b>Other Accessories</b>	
<b>500028</b>	1× CCD Camera Coupler
<b>500261</b>	0.35× CCD Camera Coupler, C-Mount (Use with USBCAM33)
<b>500262</b>	0.5× CCD Camera Coupler, C-Mount (Use with COLCAM, USBCAM50)
<b>Z-LITE</b>	Z-Lite Fiber Optic Illuminator
<b>Z-LITE-186</b>	Z-Lite Fiber Optic Illuminator with Bifurcated Light Guide and Lenses
<b>500186</b>	Bifurcated Light Guide with Lenses
<b>R-8-8-WPI01</b>	Ring Light Guide
<b>502015</b>	Ring Light Adapter for PZMIV (For R-8-8-WPI01 Ring Light Guide)
<b>505131</b>	Replacement Trinocular Tube
<b>505132</b>	Replacement Eye Guards for 30 mm Eyepieces (pair)

\*Long working distance (LWD) option with the 502012 and 502017 on a trinocular microscope shows nearly the same scene on the screen as the viewer sees in the eyepieces.

**PZMIV SPECIFICATIONS**

EYEPIECES	WFH 10×
AUXILIARY LENSES	1×
ZOOM RANGE	0.62× - 5×
TOTAL MAGNIFICATION	6.2× - 50×
ZOOM RATIO	8 : 1
FIELD OF VIEW	∅33.9-∅4.2 mm
WORKING DISTANCE	80 mm
BINOCULAR TUBE	Inclined 45°
INTERPUPILARY DISTANCE	50 - 75 mm
DIOPTER ADJUSTMENT	± 5 Diopter
MICROSCOPE BODY	Rotatable 360°
<b>OPTIONAL ACCESSORIES</b>	
Eyepieces	16×, 20×, 25×
Auxiliary lenses	0.32×, 0.5×, 0.63×
Total Magnification	1.9× - 125×
Field of view	∅106 - ∅1.8 mm
Working Distance	80 - 296 mm
SHIPPING WEIGHT	23 lb.

# Precision Stereo Zoom Microscope

*Quality and precision to improve your vision*

## Features

- Stereo viewing with ample working distance when used with the long working distance (LWD) option
- Affordably priced

## Benefits

- Advanced optics
- True trinocular version available
- Perfect for a boom stand or articulating arm

## Applications

- Tissue Dissection
- Cell injection
- Specimen manipulation
- Electrode inspection

WPI's third-generation stereo microscope, **PZMIII**, is an ideal tool for tissue dissection, cell injection, specimen manipulation, electrode inspection, and many other applications that require a magnified, stereo viewing and ample working distance. It offers the leading brand's quality and performance at an affordable price. Advanced optics provide the sharpest image that can only be found among the best of this class. It is superior to many stereo microscopes costing almost twice as much. Zooming is achieved by a spring-loaded knob that is smooth and effortless. The compact size and light weight make it more stable and easily manipulated on the boom stand. A specially designed photo/video module is used in the trinocular version of the microscope (**PZMTIII**) for photo, video, or digital imaging. In addition, an extensive list of optional accessories is available that can make the PZMIII suitable for almost any bio-research applications requiring a stereo microscope. See next page for options.


The PZMTIII trinocular version of this microscope offers a true trinocular view. Both eyepieces and the trinocular ports are active all the time. There is no need for the right eye piece to switch off. The left eye piece dims slightly, and the right eye piece is 50/50.

*Z-LITE fiber optic illuminator now included with systems below differs from unit in photo — see page 191.*



*Computer and display not included*

**PZMTIII-BS-DIG50** includes Boom Stand and USB Digital Camera



**503051 Manual Stage** — Mounts in the circular opening in the PZMIII base. XY travel distance: 75 x 56 mm. Glass size: 116 x 96 mm. Active diameter: 37.6 mm. Dimensions: 180 x 155 x 27 mm. **Fits 503102 base only.**



**503102 Stand**

## PZMIII Eyepieces and Objectives

Objective	10x Eyepiece		15x Eyepiece		20x Eyepiece		25x Eyepiece		Working Distance (mm)
	Mag.	Field (mm) (Video Field)†	Mag.	Field (mm) (Video Field)†	Mag.	Field (mm) (Video Field)†	Mag.	Field (mm) (Video Field)†	
0.3x	2x - 13.5x	114 - 17 (53.6 - 8)	3x - 20.3x	84 - 13 (53.6 - 8)	4x - 27x	69 - 10.3 (53.6 - 8)	5x - 33.8x	44.8 - 6.7 (53.6 - 8)	287 mm
0.5x	3.4x - 22.5x	69 - 10 (32.4 - 4.7)	5x - 33.8x	51 - 7 (32.4 - 4.7)	6.7x - 45x*	42 - 6.2 (32.4 - 4.7)	8.4x - 56.3x	26.9 - 4.0 (32.4 - 4.7)	177 mm*
0.75x	5x - 33.8x	45 - 7 (21.1 - 3.3)	7.5x - 50.6x	34 - 5 (21.1 - 3.3)	10x - 67.5x	28 - 4.2 (21.1 - 3.3)	12.6x - 84.4x	17.9 - 2.7 (21.1 - 3.3)	117 mm
1.0x	6.7x - 45x	34 - 5 (16 - 4.7)	10x - 67.5x	25 - 3.7 (16 - 4.7)	13.4x - 90x	21 - 3.1 (16 - 4.7)	16.8x - 112.5x	13.4 - 2.0 (16 - 4.7)	100 mm
1.5x	10x - 67.5x	23 - 3.4 (10.8 - 1.6)	15x - 101.3x	17 - 2.5 (10.8 - 1.6)	20.1x - 135x	14 - 2.1 (10.8 - 1.6)	25.1x - 168.8x	9.0 - 1.3 (10.8 - 1.6)	47 mm
2.0x	13.4x - 90x	12 - 2.5 (5.6 - 1.17)	20.1x - 135x	13 - 1.8 (5.6 - 1.17)	26.8x - 180x	10 - 1.5 (5.6 - 1.17)	33.5x - 225x	6.7 - 1.0 (5.6 - 1.17)	26 mm

† The video field of view is based on a 1/2-inch (8 mm diagonal) CCD camera and a 0.5x camera adapter.

\*Long working distance (LWD) configuration

### WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com Germany: +49 (0)6031 1602171 • wpide@wpi-europe.com • www.wpi-europe.com  
 Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • www.wpiinc.com China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.com





**502007**  
Articulated Arm

**502163** Wall-Mount Plate  
— Mount the microscope  
on the wall for convenient  
storage when space is tight.

Max. clamp  
opening: 40 mm

**504134** LED Ringlight

See page 191.



**502006** Boom Clamp Stand —  
Combines the stability of a boom  
stand with the versatility of an  
articulated arm.



PZMIII-BS

### PZMIII SPECIFICATIONS

EYEPIECES	WFH 10×
ZOOM RANGE	0.67× - 4.5×
TOTAL MAGNIFICATION	6.77× - 45×
FIELD OF VIEW	∅ 34 mm - ∅ 5 mm
WORKING DISTANCE	100 mm
BINOCULAR TUBE	Inclined 45°
INTERPUPILLARY DISTANCE	Adjustable 47-70 mm
DIOPTER ADJUSTMENT	±5 Diopter (both eyepieces)
MICROSCOPE BODY	Rotatable 360°
AUXILIARY LENSES	Total magnification 2× - 225× Biggest Field of View ∅ 110 mm Working Distance 26-287 mm
SHIPPING WEIGHT	23 lb.

### ORDERING INFORMATION

<b>PZMIII</b>	Precision Stereo Zoom Microscope (Model III), on Post Stand
<b>PZMIII-BS</b>	PZMIII Microscope on Boom Stand
<b>PZMIII-AAC</b>	PZMIII Microscope on Articulated Arm with Table Clamp
<b>PZMIII</b>	Precision Stereo Zoom Trinocular Microscope (Model III)
<b>*PZMIII-DIG50</b>	PZMIII Microscope System (100 mm WD) Including PZMIII, USBCAM50 Computer Camera, 0.5× CCD Camera Coupler, Z-LITE Optical Illuminator, Bifurcated Optical Fiber Light Guide
<b>*PZMIII-LWD-DIG50</b>	PZMIII Microscope System (177 mm WD) Including PZMIII, Pair of 10× Eyepieces and Eye Guards, Ring Light Adapter, USBCAM50 Camera, 0.5× CCD Camera Coupler, Z-Lite Optical Illuminator With Bifurcated Optical Fiber Light Guide, 0.5× Objective, 20× Eyepieces
<b>*PZMIII-BS-DIG50</b>	PZMIII Microscope System (100 mm WD) Including PZMIII, USBCAM50 Computer Camera, 0.5× CCD Camera Coupler, Z-LITE Optical Illuminator, Bifurcated Optical Fiber Light Guide, Boom Stand
<b>*PZMIII-BS-LWD-DIG50</b>	PZMIII Microscope System (177 mm WD) Including PZMIII, USBCAM50 Computer Camera, 0.5× CCD Camera Coupler, Z-LITE Optical Illuminator, Bifurcated Optical Fiber Light Guide, 0.5× Objective, 20× Eyepieces for Long Working Distance Viewing
<b>PZMIII-BS</b>	PZMIII Microscope on Boom Stand
<b>PZMIII-AAC</b>	PZMIII Microscope on Articulated Arm with Table Clamp All PZMIII and PZMIII microscopes come with 10× eyepieces and built-in 1× auxiliary lens.
<b>PZMIII-AAC-DIG50</b>	PZMIII Microscope System Including PZMIII, 10× Eyepieces and Eyeguards, USBCAM50 Camera, 0.5× CCD Camera Coupler, Ring Light Adapter, and Articulated Arm With Table Clamp

\* Z-LITE-Z may be substituted in EU countries.

### OPTIONAL ACCESSORIES/REPLACEMENT PARTS

#### Optional Bodies

<b>501352</b>	PZMIII Binocular Body, 10× eyepieces (pair), eye guards 13338 Ring Light Adapter NOT included
<b>501379</b>	PZMIII Trinocular Body, 10× eyepieces & eye guards True trinocular view — 13338 Ring Light Adapter NOT included

#### Optional Stands, Mounts, Bases

<b>501353</b>	Fan Post Stand with 76mm Focus Mount
<b>502004</b>	Boom Stand (Heavy) without Focus Mount
<b>502005</b>	Ball Bearing Boom Stand (Heavy) without Focus Mount
<b>502006</b>	Boom Clamp Stand (Heavy) (requires 502009 Focus Mount)
<b>502007</b>	Articulated Arm, 40 cm Table Clamp, without Focus Mount
<b>502009</b>	Universal Focus Mount, 76 mm ID for PZMIII, PZMIV body
<b>502163</b>	Wall-Mount Plate, 6" × 6" (or 15.24 cm × 15.24 cm)
<b>503051</b>	Manual Stage for PZMIII
<b>503102</b>	76 mm Rectangular Base Post Stand for PZMIII
<b>504123</b>	Post Extension for Heavy Boom Stand
<b>504929</b>	LED Microscope Stand, 12.5-in. Post

#### Optional Eyepieces

<b>501369</b>	Wide Field 10× Eyepieces (pair)
<b>501370</b>	Wide Field 15× Eyepieces (pair)
<b>501371*</b>	Wide Field 20× Eyepieces (pair)
<b>501372</b>	Wide Field 25× Eyepieces (pair)
<b>504128</b>	10× Eyepiece with Reticle (matches 10× eyepiece on PZMIII)
<b>504129</b>	20× Eyepiece with Reticle (matches 20× eyepiece 501371)

#### Optional Objectives

<b>501373</b>	0.3× Long Working Distance Objective Lens
<b>501375*</b>	0.5× Long Working Distance Objective Lens
<b>501376</b>	0.75× Long Working Distance Objective Lens
<b>501377</b>	1.5× Long Working Distance Objective Lens
<b>501378</b>	2.0× Long Working Distance Objective Lens

#### Other Accessories

<b>13338</b>	Ring Light Adapter for PZMIII Series (included with all microscope configurations on previous page)
<b>501381</b>	0.5× C-Mount CCD Camera Coupler
<b>504596</b>	76 mm Halogen-Halogen Dual Illuminated Track Stand
<b>504597</b>	Replacement Lamp for 504596

\*Long working distance (LWD) option with the 501371 and 501375 on a trinocular microscope shows nearly the same scene on the screen as the viewer sees in the eyepieces.

# Stereo Microscope with LED Illuminated Base

*Articulating mirror ideal for Brightfield and Darkfield applications*

## Features

- Illuminated base
- Articulating mirror
- Includes base, microscope head and focus mount
- Trinocular version available
- Includes opaque, black/white contrast stage plate

## Benefits

- Oblique, transmitted light for low magnification view of tissues, cells or embryo transfer
- Adjustable brightfield and darkfield illumination
- Pseudo DIC illumination

## Applications

- Microinjection
- Nematode research

**PZMIII-MI** Microscope with Illuminated Base and Articulating Mirror is perfect for microinjection and transfection. It includes our standard stereo microscope head mounted on a research grade Brightfield/Darkfield pole type stand. It has a large stable work surface and a rotatable lens/mirror system which provides transmitted LED intensity illumination.

The sliding mirror is gimbaled, allowing for a full range of movement front to back, as well as rotation. The mirror rotates 360° on one axis and can slide for further lighting effect directionally, front to back. A knob on the right of the base adjusts the mirror and a locking ring holds the desired mirror position.

The articulating mirror lets you vary the microscope illumination from Brightfield LED to Darkfield LED at an appropriate angle. It is an effective tool for viewing live bacteria. At low magnifications, view tissues, cells or embryo transfer where oblique, transmitted illumination is critical.

## Microscope options

This unit is sold with our standard PZMIII binocular stereo microscope head. Other options are available. Speak with a specialist today to configure a microscope for your needs



**PZMIII-MI**

## Systems include

The scope mounting pole diameter is 32 mm.

This platform includes the PZMIII microscope, a focus mount, two stage clips, 94.5 mm glass stage plate, opaque black/white contrast stage plate and an automatic voltage sensing power supply and a 25.6 cm vertical post. The base may be fitted to other microscope heads.

## PZMTIII-MI SPECIFICATIONS

DIMENSIONS	13 x 11 x 3.37 in. (33 x 28 x 8.5 cm)
PILLAR HEIGHT	10 in. (25.6 cm)
GLASS STAGE PLATE	3.75 in. (9.5 cm) diameter
STAGE CLIPS	75 mm clip length, 4 mm pin, paired
ILLUMINATION PORT	Built in incident and transmitted variable intensity LED illumination
ROTATABLE MIRROR	Dual-reflection lens/mirror system provides transmitted brightfield/ pseudo-darkfield illumination
STAND MANUFACTURER	Made in Japan
WARRANTY	Limited lifetime warranty
SHIPPING WEIGHT	10 lb. (4.5 kg)

## ORDERING INFORMATION

**PZMIII-MI** Stereo Microscope with LED Illuminated Base Stand with Articulating Mirror  
Includes PZMIII microscope, a focus mount, 2 stage clips, 94.5 mm glass stage plate, black/white stage plate, automatic voltage sensing power supply

**PZMTIII-MI** Trinocular Microscope with LED Illuminated Base Stand with Articulating Mirror  
Includes PZMTIII microscope, a focus mount, 2 stage clips, 94.5 mm glass stage plate, black/white stage plate, automatic voltage sensing power supply



**Left:** If you want to add a camera or connect with a video monitor, select a PZMTIII trinocular head for your lighted microscope base.



**Right:** This trinocular microscope head has a camera attached. It also shows an M3301 mounted to the base using a M4C coupler for microinjection studies. A setup like this facilitates classroom or collaborative environments so everyone can see on the remote screen what the researcher is viewing through the microscope.

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • [www.wpi-europe.com](http://www.wpi-europe.com)  
Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • [www.wpiinc.com](http://www.wpiinc.com)

Germany: +49 (0)6031 1602171 • wpide@wpi-europe.com • [www.wpi-europe.com](http://www.wpi-europe.com)  
China: +86 21 6888 5517 • chinasales@china.wpiinc.com • [www.wpiinc.net](http://www.wpiinc.net)

# Precision Surgical Microscope

*Ideal for small animal surgery*

## Features

- Motorized focusing system, allows hands-free operation
- Optional video adapter
- Improved optics 119 lp/mm
- Convenient handles
- New head tilting mechanism

## Benefits

- Light weight, compact and easy to maneuver—weighs only 70 lb.
- Dual bulbs prevent illumination failure during surgery
- Five magnification steps

## Applications

- Small animal surgery
- Veterinary surgery

WPI's popular precision SurgioScope with five magnification steps is a portable, high quality surgical microscope offering outstanding image quality and value. Incorporating an agile extension arm and excellent working distance objectives, the SurgioScope provides convenient movement and maneuverability necessary for accurate positioning. These important features, together with a high quality optical system, provide sharp image contrast and enhanced large field of vision. The SurgioScope comes fully equipped with a foot-controlled motorized focusing system, normally only found in more expensive surgical microscopes. A unique dual lamp housing enables safe and rapid changing of the lamp during an operation, without the need to power down. The optional video port on the trinocular version permits operational procedures to be monitored or recorded simultaneously using a video recorder and a COLCAM-HD video camera or digital stills with USBCAM50.



## SURGIOSCOPE SPECIFICATIONS

TOTAL MAGNIFICATION (F200)	3.2x — 25x
ADJUSTABLE DIOPTRER	± 6 Diopter
ADJUSTABLE INTERPUPILLARY DISTANCE	min. 50 mm — max. 70 mm
EYEPiece	12.5x
FINE FOCUS ADJUSTMENT RANGE	30 mm
WORKING HEIGHT (Arm Movement Range Above Floor)	
89 cm Post.....Focus from 34.5" (88 cm) to 51" (130 cm) above floor *	
103 cm Post.....Focus from 40.5" (103 cm) to 57" (146 cm) above floor *	
* Subtract Working Distance for height above specimen, 103 cm post recommended for F350 objective.	
RANGE OF MOTION	
Maximum Stretch Radius of Arm ..	870 mm
Vertical Movement Range .....	700-1100 mm
ILLUMINATION	
Spot =	42 mm
Dual lamp housing with quick-change spare and internal coaxial fiber optic cable.	
HALOGEN-TUNGSTEN LAMP	12 V, 100 W, with cold reflection
OPTIONAL CAMERA	COLCAM-HD, COLCAM-HD1080P, USBCAM50 (1/2" CCD) USBCAM33 (1/3" CCD)
POWER	110 V, 50-60 Hz, or 220 V, 50-60 Hz
SHIPPING WEIGHT	94 lb. (43 kg)

## OBJECTIVES

Objective	Working Distance	Magnification step	Visual Field of view (mm)	Camera field 1/2" CCD (mm)	Camera field 1/3" CCD (mm)
<b>F100</b> #504284	90 mm	6.4, 10, 16, 26, 40x	25, 15.5, 10, 6, 4	25, 15.5, 10, 6, 4.5	17.5, 11.5, 7, 4.6, 2.8
<b>F200</b> (included)	190 mm	3.2, 5, 8, 13, 20x	50, 31, 20, 12, 8	50, 31, 20, 12, 8	35, 23, 14, 9, 5.5
<b>F250</b> #504285	240 mm	2.6, 4, 6.4, 10.4, 16x	65, 40, 25, 16, 10	63, 40, 25, 16, 10	45, 28, 18, 11, 7
<b>F300</b> #504286	290 mm	2.1, 3.3, 5.3, 8.7, 13x	75, 46.5, 30, 18, 12	75, 46.5, 30, 18, 12	52.5, 34.5, 21, 13.5, 8.3
<b>F350</b> #504287	340 mm	1.8, 2.9, 4.6, 7.4, 11x	91, 57, 36, 22, 14	88, 55, 35, 21, 13	60, 38, 24, 15, 9.5

For additional objectives and specifications, see [www.wpiinc.com](http://www.wpiinc.com)

## ORDERING INFORMATION

<b>PSMB5N</b>	Binocular SurgioScope, F200 objective (Specify post height)
<b>PSMT5N</b>	Trinocular SurgioScope, beam splitter, standard video adapter, F200 objective (Specify post height)

**Specify 89 cm or 103 cm post**  
**Specify line voltage**

## OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>501636</b>	1/2" CS-mount Adapter (requires Beam Splitter 501637)
<b>501637</b>	Beam splitter
<b>504284</b>	F100 Objective
<b>504285</b>	F250 Objective
<b>504286</b>	F300 Objective
<b>504287</b>	F350 Objective
<b>500162</b>	Replacement lamp, 12V, 100 W

# Inverted Trinocular Microscope

*Excellent for video recording and photography*

## Features

- Phase contrast at 10X and 20X
- 160 x 260 mm stage
- Infinity Optical System
- Trinocular head for attaching camera
- Abbe condenser
- Coarse and fine focusing
- Multiple objectives included

## Benefits

- Stage includes inserts for holding standard culture dishes and slides
- Halogen illumination, offering an image with true to life coloring
- Fixed Stage—optics move during focusing—excellent for patch clamp and brain slice recording

## Applications

- Video recording/photography
- Patch clamping
- Brain slice recording

The **INV-101** is an affordable inverted microscope for pathologists, biologists and medical researchers. Able to perform intricate and varied applications for vital tissue cultures, it is perfect for observation and education or for professional research in a clinical lab.

Accessories included: Green, blue and neutral filters, dust cover, immersion oil, manual and warranty card.



INV-101

## TRINOCULAR HEAD

EYEPIECES	EW 10× /22 extra wide field eyepieces
CONDENSER	Abbe condenser with 10× / 20× phase
FOCUS	Coarse adjustment range +8 to -3 mm Fine adjustment: 0.002 mm
STAGE	160×250 mm
INSERTS	35 mm round, 50 mm round, 87×46 rectangular
STAGE DRIVE	Coaxial drive controls
MECHANICAL STAGE	X-Y coaxial control; 120×78mm range of traverse Nosepiece
NOSEPIECE	Quintuple nosepiece
OBJECTIVES	Infinity optical system Plan 4× and 40×, Plan Phase 10× and 20×
ILLUMINATION	6V/30 W halogen bulb Centering eyepiece for phase objective
FIXED STAGE	Optics move during focusing—excellent for patch clamp and brain slice recording

## ORDERING INFORMATION

<b>INV-101</b>	Trinocular Inverted Microscope
<b>503510</b>	30 mm 10× Eyepiece with 100/10 reticle
<b>503520</b>	Replacement lamp

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • [www.wpi-europe.com](http://www.wpi-europe.com) Germany: +49 (0)6031 1602171 • wptide@wpi-europe.com • [www.wpi-europe.com](http://www.wpi-europe.com)  
Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • [www.wpiinc.com](http://www.wpiinc.com) China: +86 21 6888 5517 • chinasaes@china.wpiinc.com • [www.wpiinc.net](http://www.wpiinc.net)

# LED Illuminated Microscopes

*Superb optics, durable and high performance microscopes*

## Features

- Binocular or trinocular
- Infinity plan achromat optics
- 20 mm FOV
- 3 W LED illuminator

## Benefits

- Super bright LED illuminator

## Applications

- Student laboratories
- Life science research

The ACCU-SCOPE® 3000-LED Microscope Series delivers outstanding optical performance, value and resolution to meet the exacting standards of life science professionals and students. With a newly designed infinity plan optical system, a best-in-class 20 mm field of view and a super-bright 3-watt LED illuminator, the 3000-LED Series provides high contrast images with outstanding resolution, precision design and enhanced illumination.



504221

## 3000-LED SERIES SPECIFICATIONS

OPTICAL SYSTEM	Infinity Optical System, f=180 mm, Anti-Mold
VIEWING HEAD	Siedentopf type, inclined 30°; interpupillary distance adjustment 48-75 mm
EYEPIECES	HWF Plan 10× eyepiece, 20 mm field of view with built-in diopter adjustment; a pointer is standard in one eyepiece
NOSEPIECE	Rear facing quadruple
OBJECTIVES	Infinity Plan achromat 4× (N.A. 0.10), 10× (N.A. 0.25), 40×R* (N.A. 0.65), 100×R oil* (N.A. 1.25) are standard *spring-loaded
STAGE	Stage size 140 mm × 132mm with X-Y movement range of 76mm × 50 mm
PHASE CONTRAST (OPTION)	Phase sliders for 10×/BF/40×
CONDENSER	N.A. 1.25 Abbe condenser
ILLUMINATION	3 watt LED with variable intensity control
ACCESSORIES	Dust cover, immersion oil and instruction manual, universal power supply 110v-240v
STAND	Cast alloy aluminum; coaxial coarse and fine focusing controls.
DIMENSIONS	15.25 × 7.75 × 15.4 in. (387 × 196 × 391 mm)
WEIGHT	16 lb. (7.26 kg)

## ORDERING INFORMATION

<b>504221</b>	Binocular, Infinity Plan Achromat 4×, 10×, 40×R and 100×R oil objectives
<b>504443</b>	Trinocular, Infinity Plan Achromat 4×, 10×, 40×R and 100×R oil objectives
<b>504445</b>	Binocular, Plan Phase Contrast 10× and 40×R objectives
<b>504444</b>	Trinocular, Plan Phase Contrast 10× and 40×R objectives
<b>504416</b>	0.50× C-mount Adapter for 1/2" Sensors, adjustable focus
<b>504417</b>	0.35× C-mount Adapter for 1/3" Sensors, adjustable focus

# Professional-Grade Microscope

*Best-seller in universities, medical schools, research laboratories*

## Features

- Titanium finished DIN or semi-plan optic
- 30-year antifungal coating

## Benefits

- Affordable research grade microscope

## Applications

- Universities, medical schools, and research laboratories

The **W30S** professional-grade microscope is a best-seller in universities, medical schools, and research laboratories. Equipped for performance, its features include titanium-finished DIN or Semi-Plan optics and a 30-year antifungal coating. The **W30S** is the choice for superior performance at a great price.



W30S-LED

## W30S SPECIFICATIONS

HEAD	Binocular (Seidentopf), True Trinocular Inclined 30°, rotates 360° Dual diopter adjustment, Interpupillary distance range 55-75 mm 10×/18 wide field eyepieces
NOSEPIECE	Quadruple forward-facing nosepiece
OBJECTIVES	DIN Plan, antifungal 4×, 10×, 40×, 100×R (oil) Parfocal, parcentric, color-coded
STAGE	Mechanical stage (140 mm × 140 mm) Coaxial drive controls XY Movement: 73 mm × 43 mm
FOCUS	Coarse adjustment: range of 30 mm Fine adjustment: graduation of 2 μm Tension control knob
ILLUMINATION	Moveable Abbe condenser, NA 1.25, Iris diaphragm Variable LED light source (3W bulb) 110 V/220 V switchable electronics
ACCESSORIES INCLUDED	Replacement 0.5 amp fuses, mirror attachment (for field use), blue and green filters, dust cover, immersion oil
DIMENSIONS AND WEIGHT	15" (38 cm) × 9" (23 cm) × 7" (17.8 cm) 14 lb. (6.4 kg)

## ORDERING INFORMATION

<b>W30S-LED</b>	Binocular Microscope
<b>W30ST-LED</b>	Trinocular Microscope
<b>503513</b>	21 mm 10X Eyepiece with 100/10 reticle
<b>500828</b>	Stage Micrometer, 1 mm scale, 200 div. at 10 μm
<b>504606</b>	Stage Micrometer, 50-0.5 mm scale and 10 μm scale

# High Definition Camera & Monitoring System

*HD imaging for scientific and industrial applications*

## Features

- High definition camera and monitoring system
- 30 fps (HDMI) and 15 fps (USB 2.0), Aptina sensor and an 4 GB SD card
- Includes HD camera, 1/2.5" (7.182 diagonal), 11.6" HD display screen, 4 GB SD card, USB 2.0 mouse, CaptaVision PC Imaging Software, HDMI Cable, USB 2.0 Cable, 12 V Power Adapter, HDMI Adapter, Y Power Splitter, Mounting Brackets and Hardware
- Specify line voltage

## Benefits

- Connect your way: use directly, USB to PC or HDMI to a projector
- Superior performance and color
- Built-in mouse control software for HDMI viewing and recording without a computer
- Fluorescence imaging with 3D noise reduction
- Many built-in camera functions

## Applications

- Scientific and industrial imaging

The new **PRO-300HDS** sets a new standard for excellence in high definition imaging for scientific and industrial applications. This full featured HD camera offers super fast frame rates in video preview, with unrivaled color fidelity and on-board image capturing capability.

**PRO-300HDS** lets you view and capture images and video directly to the supplied SD card without the need for a computer or separate monitor. The 11.6-in. HD display offers beautiful, crystal-clear image quality, and vibrant, true-to-life color with exceptional viewing from all angles.

### Connect your way

**PRO-300HDS** can be used as a stand-alone system or connected to a PC via USB cable so images can be displayed simultaneously on a PC. For even more flexibility—and ideal for teaching environments—the HDMI architecture allows the device to be connected to an HDMI-enabled projector as well. Each camera is supplied with the PC/Mac-compatible on-board image capture software, as well as full-featured CaptaVision PC image analysis software.

### Superior performance and color

**PRO-300HDS** incorporates dual FPGA processors and unique algorithms that produce perfect color reproduction. It offers full 1080p video preview for HD output with no lag time or compression.

### Built-in mouse control software for HDMI viewing and recording without a computer

An external HDMI port allows you to connect directly to an HD monitor for live/real time viewing, capturing and saving of images to an SD card without being connected to a computer. On-board software lets you control the camera with a click of the mouse rather than searching for buttons on the camera, making **PRO-300HDS** the ideal choice for teaching, group presentations or when a computer set-up is not an option.

### Fluorescence imaging

The ultra-high signal-to-noise ratio sensor lets you set exposure time from 1 ms to up to 10 seconds and adjust 20 scales of gain value. The incredibly efficient 3D noise reduction performance delivers detailed low light images, making **PRO-300HDS** a great choice for fluorescence applications.



**PRO-300HDS**  
camera & screen  
mounted on  
**PZMTIII-MI**  
(sold separately).

### Built-in camera functions

**Auto or Manual Exposure Time Capability**—You are in full control of exposure and gain. Use the Auto Exposure function or set the exposure time from 1 ms up to 10s and adjust 20 scales of gain value.

**3D Noise Reduction**—Longer exposure times increase image noise. The integrated 3D noise reduction keeps the image clean and sharp.

**1080P Video Recording**—Just click on the video record icon to start recording 1080P videos at 30fps. The recorded video files are saved directly to the high speed SD card. You can playback videos directly from the SD card.

**Magnify (ROI), Rotate and Flip Images**—Image operation buttons on the right side of the screen allow you to select a ROI (Region of Interest), as well as flip or rotate an image.

**Side by Side Image Comparison**—The image comparison function allows you to choose one image, move the image position or select the ROI area to compare with the live image.

**Browse Captured Images and Video**—Easily browse images on the SD card, zoom in on images or delete them. You can even playback video files saved to the SD card from the Browse feature.

**On Board Imaging Tools**—No driver is needed when you connect to a PC or Mac via the USB 2.0 camera port. You'll be up and running in no time with the on-board software that features basic imaging tools:

- Flip Horizontally
- Flip Vertically
- Zoom in
- Zoom Out
- Crop
- Divison
- Cancel
- Insert a Crossline
- Compare Images
- Browse Images



## ORDERING INFORMATION

**PRO-300HDS** High Definition Camera with Screen

*Specify line voltage*

*Includes HD camera, 11.6" HD display screen, 4GB SD card, USB 2.0 mouse, CaptaVision PC Imaging Software, HDMI Cable, USB 2.0 Cable, 12V Power Adapter, HDMI Adapter, Y Power Splitter, Mounting Brackets and Hardware*

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • [www.wpi-europe.com](http://www.wpi-europe.com)  
Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • [www.wpiinc.com](http://www.wpiinc.com)

Germany: +49 (0)6031 1602171 • wpide@wpi-europe.com • [www.wpi-europe.com](http://www.wpi-europe.com)  
China: +86 21 6888 5517 • chinasales@china.wpiinc.com • [www.wpiinc.net](http://www.wpiinc.net)

# Digital Microscope Cameras

Photograph or share your research images

## USBCAM133/USBCAM152/USBCAM202

### Features

- Ultra-Compact USB cameras with color CCD

### Benefits

- Rear mini-USB connector
- Hardware & Software Trigger

### Applications

- Image capturing from a microscope



USBCAM133

### Options

- **USBCAM133** with a 1/3-inch color CCD, 6mm diagonal
- **USBCAM152** with a 1/2-inch color CCD, 8mm diagonal
- **USBCAM202** with a 1/1.8-inch color CCD, 8.9 mm diagonal

All three cameras are cased models with a rear mini-USB connector, hardware and software triggering image capture, digital zoom and a feature-rich, user based menu setup and control. Sentech USB cameras include a SDK, DirectX, Twain and Linux driver, as well as the Sentech Viewing Software.

SPECIFICATIONS			
	USBCAM133	USBCAM152	USBCAM202
IMAGE SENSOR	1/3" Interline SXVGA color progressive CCD	1/2" Interline SXGA color progressive CCD	1/1.8" Interline UXGA color progressive CCD
CELL SIZE	3.75 (H) x 3.75 (V) $\mu$ m, 6mm diagonal	4.65 (H) x 4.65 (V) $\mu$ m, 8mm diagonal	4.40 (H) x 4.40 (V) $\mu$ m, 8.9 mm diagonal
SCANNING SYSTEM	Progressive	Progressive	Progressive
RESOLUTION	1280 (H) x 960 (V)	1360 (H) x 1024 (V)	1600 (H) x 1200 (V)
MIN. SCENE ILLUMINATION	11 Lux at F1.2	18 Lux at F1.2	7.7 Lux at F1.2
SPEED	22.4 Frames per Second	19.26 Frames per Second	15.3 Frames per Second
ELECTRONIC SHUTTER	Auto/Manual (software selectable)	Auto/Manual (software selectable)	Auto/Manual (software selectable)
GAIN	Auto/Manual (software selectable)	Auto/Manual (software selectable)	Auto/Manual (software selectable)
GAMMA	Manual (software selectable)	Manual (software selectable)	Manual (software selectable)
WHITE BALANCE	Auto/Manual/One shot (software selectable)	Auto/Manual/One shot (software selectable)	Auto/Manual/One shot (software selectable)
INPUT/OUTPUT	USB 2.0 High Speed	USB 2.0 High Speed	USB 2.0 High Speed
POWER	+5 VDC through USB connector, < 300 mA	+5 VDC through USB connector, < 420 mA	+5 VDC through USB connector, < 450 mA
DIMENSIONS	28 (W) x 28 (H) x 37 (D) mm (excluding connector)	28 (W) x 28 (H) x 42 (D) mm (excluding connector)	28 (W) x 28 (H) x 42 (D) mm (excluding connector)
LENS MOUNT	CS mount	C mount	C mount
WEIGHT	Approximately 45 g	Approximately 45 g	Approximately 45 g
INTERFACE CONNECTOR	USB: mini-B USB connector IO signal: 6 pin connector (HR10A-7R-6PB or equivalent)	USB: mini-B USB connector IO signal: 6 pin connector (HR10A-7R-6PB or equivalent)	USB: mini-B USB connector IO signal: 6 pin connector (HR10A-7R-6PB or equivalent)
RoHS	RoHS Compliant	RoHS Compliant	RoHS Compliant

### ORDERING INFORMATION

<b>USBCAM133</b>	Digital Microscope Camera
<b>USBCAM152</b>	Digital Microscope Camera
<b>USBCAM202</b>	Digital Microscope Camera
<b>504570</b>	Replacement USB Cable, 3 m (10 ft.)

## USBCAM33/USBCAM50

### Features

- Automatic and manual controls
- C/CS mount
- Software included

### Benefits

- Connect via USB

### Applications

- Image capturing from a microscope



USBCAM33

Record images directly to your computer. These digital microscope cameras offer flexibility, with a range of configurations for image capture, a choice of mount option (C or CS) and file output alternatives. Since both cameras connect via the USB port, installing the image capture software is simple. Either camera can be used on WPI's stereo microscopes **PZMTIV**, **PZMTIII**, compound microscopes **W30ST** and **GPL-T** and also the **PSMT5** Surgical Microscope. Choose from the one third-inch CCD with 1024x768 resolution and 30 frames per second (**USBCAM33**) or one half-inch CCD with 1280x960 resolution and 15 frames per second (**USBCAM50**).

These cameras include **IC Imaging Control** software that has:

- Real-time video preview
- Text and graphics can be drawn on a live video stream
- Scroll and Zoom
- Acquisition of single frames
- Capture pause, for intermittent image capture
- Timestamps

### SPECIFICATIONS

	USBCAM33	USBCAM50
IMAGE SENSOR	1/3" Sony CCD, progressive scan	1/2" Sony CCD, progressive scan
MAX RESOLUTION	1024 x 768	1280 x 960
SIZE	4.65 $\mu$ m x 4.65 $\mu$ m, 6mm diameter	4.65 $\mu$ m x 4.65 $\mu$ m, 8mm diameter
SPEED (PC DEPENDENT)	30 fps, 15 fps, 7.5 fps or 3.75 fps	15 fps, 7.5 fps or 3.75 fps
SENSITIVITY	0.5 lux @ 1/15 s	0.5 lux @ 1/7.5 s
EXPOSURE, SHUTTER CONTROL, WHITE BALANCE	Automatic/Manual	Automatic/Manual
INTERFACE	USB 2.0 cable	USB 2.0 cable
SYSTEM REQUIREMENT	Windows Vista (32 & 64 bit), Windows 7 (32 & 64 bit), Windows 8 (32 & 64 bit) or Windows 10 (32 & 64 bit)	Windows Vista (32 & 64 bit), Windows 7 (32 & 64 bit), Windows 8 (32 & 64 bit) or Windows 10 (32 & 64 bit)
SOFTWARE	IC Imaging Control Software	IC Imaging Control Software
LENS MOUNT	C/CS-Mount	C/CS-Mount
CAMERA BODY	50.6 x 50.6 x 50 mm	50.6 x 50.6 x 50 mm
WEIGHT	265 g (9.5 oz.)	265 g (9.5 oz.)

### ORDERING INFORMATION

<b>USBCAM33</b>	Digital Microscope Camera, 1/3-in. CCD
<b>USBCAM50</b>	Digital Microscope Camera, 1/2-in. CCD
<b>503536</b>	Cable, USB Extension (male-female)

# Color CCD Video Camera for Microscopy

Watch live video images

## COLCAM-HD1080P

### Features

- 16:9 Aspect Ratio, 1:1, no scaling
- DVI signal output via HDMI cable
- Improved design of COLCAM-HD with better low light sensitivity
- Low cost alternative to existing HD and 3 CCD cameras, with outstanding image quality



COLCAM-HD1080P

### Benefits

- Directly connect to HDTV, no PC required and no software to load

### Applications

- Microscopy video output to HDTV

The **COLCAM-HD1080P** is a CMOS-based camera that outputs a true HD 1080P or 720P image at 60fps in the 16x9 format, with the capability to program individual DSP profiles accessed via remote hand held controller.

If still image or HDTV recording is required, contact WPI.

### COLCAM-HD1080P SPECIFICATIONS

IMAGER	1/2.33" 14 MP CMOS Progressive
HD ACTIVE PICTURE ELEMENT	1920 (H) x 1080 (V)
CHIP SIZE	7.8 mm diagonal
MINIMUM SCENE ILLUMINATION	650 mV/Lux-sec
SYNC SYSTEM	Internal
VIDEO OUTPUT	DVI 1.0 conformity 1080P RGB, 1920 H x 1080 V 60/50/30 Hz
GAIN	AGC *
SHUTTER SPEED	Auto *
GAIN	AGC or Fixed gain *
GAMMA	presets or manual gamma *
WHITE BALANCE	Auto/Manual/Push-to-set *
POWER	12V power jack; 5.5x2.1 mm
DIMENSIONS	40 (W) x 45 (H) x 41.1 (D) mm
OPTICAL FILTER	IR cut filter included
LENS MOUNT	C/CS mount
VIDEO OUTPUT	HDMI connector
WEIGHT	88 g
ROHS	RoHS compliant

\* Selectable via the UART communication

### ORDERING INFORMATION

<b>COLCAM-HD1080P</b>	Color Video Camera, 1920 x 1080 HD resolution Includes 3m DV Cable, Power Supply, C/CS Mount
<b>504136</b>	3 meter HDMI Cable

## C-Mount Eyepiece Adapters



503097

Camera not included.

For 1/3-inch (6 mm diagonal) and 1/2-inch (8 mm diagonal) video cameras and eyepiece camera conversion, this lens and its accessories make it possible to connect a typical video or C-mount camera to almost any microscope on the market. The lens fits right into the ocular socket of standard 23.2mm microscopes and the 30 mm adapter allows for use on the typical stereo zoom microscope. If you already have a trinocular microscope you can add the included C-Adapter to the top of an existing 1X C-mount (no lens) adapter.

### ORDERING INFORMATION

<b>503097</b>	C-Mount to Eyepiece Adapter Kit, 0.45X for 1/3-in. and 1/2-in. Video Cameras, 30 mm Stereo Adapter, 1X C-Mount Adapter
---------------	--

## SLR Digital Camera-to-Microscope Eyepiece Adapter



503099

This adapter connects T-mount SLR digital cameras to almost any microscope on the market. The adapter is built to 23.2 mm ocular tubes that are found on most high magnification (upright, inverted, standard) microscopes. The 30 mm adapter allows mounting on most stereo zoom microscopes that use 30 mm oculars. If you already have a trinocular microscope, you can add this adapter to the top of an existing 1X C-mount adapter.

The 2X magnification of this microscope adapter yields an approximate 65% field of view from the visual field as measured on a Canon 80D Digital camera. (CCD Sensor size = 22.7 x 15.1 mm). 35 mm film reference size is 24 x 36mm.

Please contact your camera dealer for a suitable T-mount to bayonet adapter for your camera.

### ORDERING INFORMATION

<b>503099</b>	C-Mount to Eyepiece Adapter Kit, 2.0x for SLR digital Cameras Includes 30 mm Stereo adapter, 1X C-mount adapter
---------------	--

### WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • www.wpiinc.com

Germany: +49 (0)6031 1602171 • wptide@wpi-europe.com • www.wpi-europe.com  
China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net



# High Intensity Fiber Optic Illumination Source

*Uninterrupted, directed light for microscopes and other applications*

## Features

- Reliable, uninterrupted high-intensity light
- Use with microscopes
- Intensity controlled by rotary knob
- Use with ring light, or single/bifurcated guides, sold separately

## Benefits

- Contains full spectrum of visible light

## Applications

- General laboratory use
- Microscopy illumination

The **Z-LITE** Fiber Optic Illuminator provides reliable, uninterrupted high-intensity light for microscopes. **Z-LITE** allows a continuous range of subdued or concentrated lighting controlled by a rotary dimmer on the front panel. **Z-LITE** may be used with a ring light and single or bifurcated flexible fiber bundles, enabling the light beam to be placed exactly where it is needed. Forced air cooling prolongs the lamp life. Lamp color temperature is 3350°K. An interlock switch automatically cuts off power when the front panel is opened to replace the bulb.

## Light guides available separately



Ring Light **R-8-8-WPI01** can be used with PZM Stereo Microscope for shadow-free illumination.  
18-in. (46 cm) flexible cable

ID:  
58 mm

This single pipe light guide with focusing lens has a fiber diameter of 6mm. It is 18" (460 mm) long with a 5/8" standard adapter for attaching to a fiber optic illuminator like the Z-Lite series.



504931



Shown with Bifurcated Light Guide **500186**

Z-LITE

## LED Ring Light

*For PZMIII & PZMIV stereo microscopes*



504134

## Features

- "White" light illumination — 72 LED bulbs
- Maximum opening 61 mm
- Ring light is divided into four areas and each area is turned on and off separately
- Brightness adjustable
- ESD safe
- Power supply AC 90-264V, 50/60 Hz, US plug only

## Benefits

- Direct top illumination plus four zones for shadows
- Low cost
- Less bulky than Halogen
- Cold light with no heat

## Applications

- Stereo microscope use



## Z-LITE SPECIFICATIONS

LAMP	150 W quartz halogen (EKE lamp)
SIZE	30.5×25×25 cm (12×10×10 in.)
POWER	115 VAC, 50/60 Hz, 3 A
WEIGHT	5.9 kg (13 lb.)

## ORDERING INFORMATION

<b>Z-LITE-186</b>	Z-Lite & Bifurcated Light Guide (115 V, 60 Hz, beige case)
<b>Z-LITE-Z186</b>	Z-Lite & Bifurcated Light Guide (230 V, 50 Hz, black case) *
<b>Z-LITE</b>	Z-Lite Fiber Optic Illuminator (115 V, 60 Hz, beige case)
<b>Z-LITE-Z</b>	Z-Lite Fiber Optic Illuminator (230 V, 50 Hz, black case) *

## OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>500186</b>	Bifurcated Light Guide (with lenses)
<b>504930</b>	Flexible Light Guide with Focusing Lens, 18" (46cm)
<b>504931</b>	Flexible Light Guide with Focusing Lens, 24" (61cm)
<b>R-8-8-WPI01</b>	Ring Light Guide for PZM and PZMIII Series <i>Ring Light Guide requires adapter #13338 for use with PZM, PZMII and PZMIII, included with each PZMIII and PZMIV microscope system.</i>
<b>13338</b>	Ring Light Adapter (48 mm Ø) for PZM, PZMII, PZMIII
<b>502015</b>	Ring Light Adapter for PZMIV
<b>EJA</b>	Replacement Halogen Lamp, 150W, 3350°K, 40-hour
<b>EKE</b>	Replacement Lamp, 150W, 3250°K, 200-hour


\* Not available in EU countries

## ORDERING INFORMATION

**504134** LED Ringlight

# Microdissection Instruments

MICRODISSECTION



## Quality you can count on

WPI surgical instruments are manufactured with the highest quality materials and craftsmanship to provide you with reliable instruments at a cost-effective price. We offer a wide variety of instruments, including surgical kits. Whether you are looking for a pair of quality European surgical instruments, marine grade stainless steel forceps or precision scissors, the quality of our surgical instruments is backed by our 100% satisfaction guarantee.



WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
Germany: +49 (0)6031 1602171 • wpide@wpi-europe.com • www.wpi-europe.com  
Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • www.wpiinc.com  
China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net

# Swiss Tweezers

Our Swiss surgical tweezers are manufactured from fine Swiss steel. We offer a variety of tweezers and forceps.

## WPI SWISS TWEEZERS

- Used in microsurgical procedures to grasp fine tissue
- Sharp tips may puncture tissue
- 10.5 cm (4.1 in.) Long
- Extra fine 90° angled tips

<b>504508</b>	Stainless Steel
<b>504509</b>	Acid-Resistant/Anti-magnetic



Tips: 0.06 × 0.02 mm

**504508**



**504508**  
**504509**

## WPI SWISS TWEEZERS

- 11 cm (4.3 in.) Long
- Handcrafted to perfection in Switzerland
- For use in general laboratory and medical applications, microscopy, precision manufacturing, electronic industry, watch making
- 0.1 × 0.06 mm Straight tips

<b>504506</b>	Stainless Steel
<b>504507</b>	Acid-Resistant/Anti-magnetic



**504506**



**504506**  
**504507**

## WPI SWISS TWEEZERS

- 11 cm (4.3 in.) Long
- Diamond Like Carbon coated tips to improve hardness, toughness and wear characteristics

<b>504524</b>	Ni-Cr-Mo Superalloy
<b>504516</b>	Ni-Cr-Mo Superalloy, Diamond Coated Tips
<b>504515</b>	Anti-magnetic, Diamond Coated Tips



**504516**



**504516**



**504515**



**504524**

## WPI SWISS TWEEZERS

- 11.5 cm (4.5 in.) Long
- Stainless Steel

<b>504505</b>	Acid-Resistant/Anti-magnetic, Curved Tips
---------------	---



**504505**



**504505**

## WPI SWISS CUTTING TWEEZERS

- 11.5 cm (4.5 in.) Long
- Angled blades
- Stainless Steel

<b>504744</b>	Stainless Steel
<b>504745</b>	High Precision, Tungsten Carbide Tips



**504745**



**504744**



**504745**

# Dumont Forceps

High quality Dumont tweezers are available in a variety of metals.

**Dumostar** – Anti-magnetic. Dumostar is more elastic and more corrosion-resistant than the best stainless steel.

**Dumoxel** – Extremely flexible. Dumoxel is 95% anti-magnetic and stain resistant.

**Stainless Steel** – This is the standard metal for most surgical instruments.

## MINI DUMONT #M5S

- 8.2 cm (3.25 in.) Long
- 0.09 x 0.05 mm Tips
- Stainless Steel

**501764** Mini Dumont #M5S



501764

## DUMONT #3C

- 11 cm (4.3 in.) Long
- 0.08 x 0.04 mm Tips
- Non-magnetic
- Non-corrosive
- Dumostar Steel

**500064** Dumont #3C



500064

## DUMONT #4

- 11 cm (4.3 in.) Long
- 0.06 x 0.02 mm Tips
- Non-magnetic
- Non-corrosive
- Dumostar Steel

**500339** Dumont #4



500339

## DUMONT #4

- 11 cm (4.3 in.) Long
- Stainless Steel

**500340** 0.13 x 0.08 mm Tips  
**500231** 0.02 x 0.06 mm Tips



500340

500231

500340

## DUMONT #5

- 11 cm (4.3 in.) Long
- Non-magnetic
- Non-corrosive
- Dumostar Steel

**500085** 0.025 x 0.015 mm Tips  
**500233** 0.1 x 0.06 mm Tips



500085

# Dressing Forceps


Dressing forceps are used when dressing wounds. They hold gauze and other dressings. They may also be used during wound debridement to remove infected or necrotic tissue or debris from the wound. They may also be used for suturing.

## ADSON FORCEPS

- Used to grasp tissue during stapling of the skin
- 12 cm (4.75 in.) Long
- Serrated
- Stainless Steel

<b>14226</b>	1.5 mm Tips
<b>14226-G</b>	1.5 mm Tips, <b>German</b>
<b>500222</b>	Tungsten Carbide Jaws
<b>500222-G</b>	Tungsten Carbide Jaws, <b>German</b>

Standard and German versions may differ slightly in appearance.



500222

## IRIS FORCEPS

- Used for grasping tissue and dressing applications in microscopy
- 10 cm (4 in.) Long
- 0.8 mm Tips
- 1x2 Teeth
- Stainless Steel

<b>15916</b>	Straight
<b>15916-G</b>	Straight, <b>German</b>
<b>15917</b>	Curved
<b>15917-G</b>	Curved, <b>German</b>

Standard and German versions may differ slightly in appearance.



15916

# Hemostatic Forceps

Ring forceps (also called hemostatic forceps) are hinged and look like ring scissors. Frequently, hemostatic forceps have a locking mechanism, called a ratchet, which is used for clamping. The jaws of the locking forceps gradually come together as each increment of the ratchet is employed.

- Used to compress blood vessels or other tubular structures to obstruct the flow of blood or fluids



- Used for small, delicate vessels
- Full serrated jaws
- Used in ophthalmology, biopsy, skin grafting or general surgery and lab procedures

## HALSTED MOSQUITO HEMOSTATIC FORCEPS

- Serrated
- Stainless Steel

<b>15920</b>	12.5 cm (5 in.), Straight
<b>15920-G</b>	12.5 cm (5 in.), Straight, <b>German</b>
<b>15921</b>	12.5 cm (5 in.), Curved
<b>15921-G</b>	12.5 cm (5 in.), Curved, <b>German</b>
<b>501705</b>	9 cm (3.5 in.), Straight
<b>501705-G</b>	9 cm (3.5 in.), Straight, <b>German</b>
<b>501291</b>	9 cm (3.5 in.), Curved
<b>501291-G</b>	9 cm (3.5 in.), Curved, <b>German</b>

Standard and German versions may differ slightly in appearance.



15920

# Titanium Forceps & Needle Holders


100% anti-magnetic, corrosion-resistant, lightweight and strong, titanium alloy is ideal for biological and medical applications. Titanium has the tensile strength of carbon steel and is completely resistant to corrosion from nitric acid, chloride, salt water, and industrial and organic chemicals. Titanium is more flexible and 40% lighter than Innox. When heated or cooled, the dimensions of titanium alloy change less than half of what stainless steel alloys will, making titanium surgical instruments much more durable. Titanium is stain-free and temperature resistant up to 430°C. Titanium tools are the premium choice for corrosive environments or MRI applications.

- 100% non-corrosive (great for sea water procedures)
- 40% lighter than stainless steel (reduces hand fatigue)
- 100% non-magnetic (MRI compatible)
- Anodized, non-glare blue finish

## FORCEPS

- 8.5 cm (3.3 in.) Long
- 12 mm Tips
- 4.5 mm Tying platform
- Straight
- **Titanium**

**WP5000** Forceps





WP5000

## FORCEPS WITH TYING PLATFORM

- 8.5 cm (3.3 in.) Long
- 12 mm Tips
- 5 mm Tying platform
- Straight with 0.12 mm teeth (1x2)
- **Titanium**

**WP2000** Forceps with Tying Platform



WP2000

## TROUTMAN NEEDLE HOLDER

- 11.5 cm (4.5 in.) Long
- Curved
- 10 mm Jaws
- Blunt
- **Titanium**

**WP3620** Troutman Needle Holder




WP3620

## NEEDLE HOLDER, SERRATED JAWS

- 16 cm (6.3 in.) Long
- Straight
- Serrated Jaws
- **Titanium**

**503391** Needle Holder



503391

# Spring Scissors

*Vannas scissors – small, sharp blades designed for researchers working in tight spaces.*

- May be used with either left or right hand
- Designed for neurosurgical, vascular and ophthalmic procedures
- Cuts by shearing, reduces crushing of tissue
- Fine tips are ideal for use in very restricted spaces
- Curved tips are useful to avoid cutting of underlying tissue
- Heavier blades are helpful when a thick tissue or vessel is to be cut
- Length of the scissor tips can help with the depth of the incision:  
Shorter tips for near-to-surface incision  
Longer tips for deep incision


## VANNAS SCISSORS

- Small, long blades for uses in ophthalmology, neurosurgery, microsurgery
- Super fine
- 8 cm Long
- 3 mm Blades
- 0.015 mm Tips
- Stainless Steel

<b>501778</b>	Straight
<b>501839</b>	Curved



501778



## VANNAS SCISSORS


- Diamond shaped tips
- 9 cm (3.5 in.) Long
- 5 mm Blades
- 0.1 mm Tips
- Stainless Steel

<b>14003</b>	Straight
<b>14003-G</b>	Straight, German

Standard and German versions may differ slightly in appearance.



14003-G



## NOYES SCISSORS

- 12 cm Long
- 15 mm Tips
- Originally designed to remove the iris (iridectomy)
- Fine for multi-purpose use
- Best for trimming tissue from nerves

<b>500228</b>	Sharp/Sharp Tips, Straight
<b>500228-G</b>	Sharp/Sharp Tips, Straight, German
<b>501236</b>	Sharp/Sharp Tips, Curved
<b>501236-G</b>	Sharp/Sharp Tips, Curved, German
<b>503305</b>	Sharp/Blunt Tips, Straight
<b>503306</b>	Blunt/Blunt Tips, Straight
<b>501237</b>	14 cm, 20 mm Tips, Straight
<b>501238</b>	14 cm, 20 mm Tips, Curved

Standard and German versions may differ slightly in appearance.



500228



500228



501236



503305



503306



501237




501238

## KATENA-VANNAS SCISSORS

- 11 cm (4.3 in.) Long
- Stainless Steel

<b>504024</b>	7 mm Thin Blades, Straight
<b>504025</b>	7 mm Thin Blades, Curved
<b>504026</b>	10 mm Thin Blades, Straight
<b>504027</b>	10 mm Thin Blades, Curved




## MICRO SCISSORS – DIAMOND LIKE CARBON COATED BLADES

Diamond like carbon coated blades increase the longevity of the cutting edge and improve the quality of the cut. Used to reduce wear on razor blades and industrial metal cutting tools, this coating is now being used in biomedical applications. Some studies have shown that even after double the normal usage of the cutting implement, no blade wear was detected.

- Hardness of coated surface resists wear and increases life of cutting edge
- Anti-reflective coating resists glare and assists in better view of surgical field
- Coating is corrosion resistant
- Reduced friction at cutting site decreases necrosis of surrounding tissue and bone
- 12 cm (4.75 in.)
- Sharp coated tips
- 8 mm Blades
- Stainless Steel,


<b>503365</b>	Straight
<b>503364</b>	Curved



## SPRING SCISSORS

- 12 cm (4.75 in.) Long
- 12 mm Extra fine, long blades
- Stainless Steel

<b>14125</b>	Straight
<b>14125-G</b>	Straight, German
<b>14126</b>	Curved



Standard and German versions may differ slightly in appearance.



# Titanium Scissors

100% anti-magnetic, corrosion-resistant, lightweight and strong, titanium alloy is ideal for biological and medical applications. Titanium has the tensile strength of carbon steel and is completely resistant to corrosion from nitric acid, chloride, salt water, and industrial and organic chemicals. Titanium is more flexible and 40% lighter than Inox. This reduces hand fatigue during long procedures. When heated or cooled, the dimensions of titanium alloy change less than half of what stainless steel alloys will, making titanium surgical instruments much more durable. Titanium is stain-free and temperature resistant up to 430°C. Titanium instruments are the premium choice for corrosive environments or MRI applications.

- 100% non-corrosive (great for seawater procedures)
- 40% lighter than stainless steel (reduces hand fatigue)
- 100% non-magnetic (MRI compatible)
- Anodized, non-glare blue finish

## VANNAS SCISSORS

- 8.2 cm (3.2 in.) Long
- Sharp 9.5 mm tips
- **Titanium**

<b>WP5050</b>	Straight
<b>WP5070</b>	Angled Up



**WP5070**




**WP5050**

**WP5070**

## IRIS SCISSORS

- 9 cm (3.5 in.) Long
- Curved
- Sharp 11 mm tips
- **Titanium**

<b>WP4410</b>	Iris Scissors
---------------	---------------




**WP4410**

## VANNAS CAPSULOTOMY SCISSORS

- 9.5 cm (3.7 in.) Long
- Angled
- Sharp
- 10.5 mm Tips
- **Titanium**

<b>WP5070L</b>	Vannas Capsulotomy Scissors
----------------	-----------------------------




**WP5070L**

## VANNAS CAPSULOTOMY SCISSORS

- 9.5 cm (3.7 in.) Long
- Curved
- Sharp
- 12.5 mm Tips
- **Titanium**

<b>WP5060</b>	Vannas Capsulotomy Scissors
---------------	-----------------------------



**WP5060**

# Ring Scissors

Our standard scissors all have ring handles. You can typically recognize the construction style of scissors based on the following color code:

- Black handled rings indicate SuperCut scissors (with one serrated blade).
- Gold handles indicate that the scissors have tungsten carbide inserts.
- One gold handle and one black handle indicates that the scissors have serrations in one blade (black handle) and tungsten carbide on the other (gold handle).
- Blue instruments are typically titanium.

## IRIS DISSECTING SCISSORS WITH PROBE TIPS

- Multi-purpose
- Ideal for detailed dissection of delicate tissues and can be used in suture removal
- 10 cm (4 in.) Long
- Blunt probe tips
- Stainless Steel

<b>500366</b>	Straight
<b>500366-G</b>	Straight, <b>German</b>
<b>500367</b>	Curved
<b>500367-G</b>	Curved, <b>German</b>

Standard and German versions may differ slightly in appearance.




## OPERATING SCISSORS

- Recommended for cutting and dissecting tissue
- Also used for cutting surgical drapes
- 16 cm (6.25 in.) Long
- Stainless Steel

<b>501223</b>	Sharp/Blunt, Straight
<b>501223-G</b>	Sharp/Blunt, Straight, <b>German</b>
<b>501224</b>	Blunt/Blunt, Straight
<b>501224-G</b>	Blunt/Blunt, Straight, <b>German</b>
<b>501225</b>	Sharp/Sharp, Straight
<b>501225-G</b>	Sharp/Sharp, Straight, <b>German</b>
<b>501226</b>	Sharp/Blunt, Curved
<b>501226-G</b>	Sharp/Blunt, Curved, <b>German</b>
<b>501227</b>	Blunt/Blunt, Curved
<b>501227-G</b>	Blunt/Blunt, Curved, <b>German</b>
<b>501228</b>	Sharp/Sharp, Curved
<b>501228-G</b>	Sharp/Sharp, Curved, <b>German</b>

Standard and German versions may differ slightly in appearance.



# SuperCut Scissors

Our black handled surgical scissors designate our SuperCut scissors. These scissors have one on razor sharp blade and one micro-serrate blade. This construction offers a couple of advantages for surgeons.


- Sharp edge gives a clean cut with minimal tissue damage
- Serrated edge actually holds the tissue to prevent it from slipping while you are making an incision

## IRIS SUPERCUT SCISSORS

- One edge micro serrated, one edge honed to the sharpness of a knife edge—prevents tissue slippage
- Stainless Steel

<b>14218</b>	10 cm (4 in.), Straight
<b>14218-G</b>	10 cm (4 in.), Straight, <b>German</b>
<b>14219</b>	10 cm (4 in.), Curved
<b>14219-G</b>	10 cm (4 in.), Curved, <b>German</b>
<b>503259</b>	11.5 cm (4.5 in.), Straight
<b>503260</b>	11.5 cm (4.5 in.), Curved
<b>14225</b>	12.5 cm (5 in.), Straight
<b>14225-G</b>	12.5 cm (5 in.), Straight, <b>German</b>
<b>14224</b>	12.5 cm (5 in.), Curved
<b>14224-G</b>	12.5 cm (5 in.), Curved, <b>German</b>

Standard and German versions may differ slightly in appearance.



# Tungsten Carbide SuperCut Scissors

Our black handled surgical scissors designate our SuperCut scissors. These scissors have one on razor sharp blade and one micro-serrate blade. Surgical scissors with one black handle and one gold handle are SuperCut scissors with tungsten carbide inserts. Tungsten carbide is harder than steel, which means the instruments last longer. The SuperCut edges offer a couple of advantages for surgeons.

- Sharp edge gives a clean cut with minimal tissue damage
- Serrated edge actually holds the tissue to prevent it from slipping while you are making an incision

## IRIS SUPERCUT SCISSORS, TUNGSTEN CARBIDE

- Gold blade has tungsten carbide insert — black blade is serrated
- Produces clean cut without tearing the tissue
- Used for cutting delicate or heavy tissue
- 11 cm (4.3 in.) Long
- Sharp tips
- **Tungsten Carbide** blades
- Stainless Steel

<b>501263</b>	Straight
<b>501263-G</b>	Straight, <b>German</b>
<b>501264</b>	Curved
<b>501264-G</b>	Curved, <b>German</b>

Standard and German versions may differ slightly in appearance.




# Ring-Handle Needle Holders

## OLSEN-HEGAR NEEDLE HOLDER WITH SUTURE SCISSORS

- Straight
- Serrated
- **Tungsten Carbide** jaws
- Stainless Steel

<b>501312</b>	12 cm (4.75 in.)
<b>500227</b>	14 cm (5.5 in.)
<b>500227-G</b>	14 cm (5.5 in.), <b>German</b>
<b>501725</b>	16.5 cm (6.5 in.)
<b>501725-G</b>	16.5 cm (6.5 in.), <b>German</b>

Standard and German versions may differ slightly in appearance.



500227-G

## OLSEN-HEGAR NEEDLE HOLDER

- 14 cm (5.5 in.) Long
- Straight
- Serrated jaws with suture scissors
- Stainless Steel

<b>500023</b>	Olsen-Hegar Needle Holder
<b>500023-G</b>	Olsen-Hegar Needle Holder, <b>German</b>

Standard and German versions may differ slightly in appearance.





500023

## CRILE-WOOD NEEDLE HOLDER

- Straight
- Serrated
- **Tungsten Carbide** jaws
- Stainless Steel

<b>500224</b>	14.5 cm (5.75 in.)
<b>500224-G</b>	14.5 cm (5.75 in.), <b>German</b>
<b>503741</b>	15 cm (5.9 in.), Left-Handed

Standard and German versions may differ slightly in appearance.




500224

## MAYO-HEGAR NEEDLE HOLDER

- 15 cm (6 in.) Long
- Straight
- Stainless Steel

<b>V503382</b>	Mayo-Hegar Needle Holder
----------------	--------------------------



V503382

### WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com Germany: +49 (0)6031 1602171 • wptide@wpi-europe.com • www.wpi-europe.com  
Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • www.wpiinc.com China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net

# Black Instruments

## JAFFE FORCEPS

- 10 cm (3.9 in.)
- Curved 0.52 mm Tip
- Stainless Steel with black coating

**WPB505210** Jaffe Forceps, Black



WPB505210

## IRIS SUPERCUT SCISSORS, TITANIUM POWDER COATED

- More resistance to rust and body fluids
- Anti slippage serrated blade
- Used in microsurgery
- 11 cm (4.3 in.) Long
- Anti-reflective
- One serrated blade
- Stainless Steel

**504487** Straight  
**504500** Curved



504487

## BLACK EPOXY COATED TWEEZERS WITH PATENTED EZ LEVER

- 11 cm (4.25 in.) Long
- Strong blades
- Stainless Steel with black coating

**WPB315AB** 4.5 mm from tip to bend  
**WPB315AS** 6.5 mm from tip to bend



WPB315AB



WPB315AB



WPB315AS

## Skin Stapler


The **WPSS01** skin stapler is an affordable, reliable solution for ideal veterinary wound closures.

- Made in the U.S.A.
- Designed specifically for veterinarians
- Clear handle makes it easy to spot contaminants
- Uses reloadable cartridges containing 35 stainless steel staples
- Autoclavable handle
- Improves patient safety by considerably reducing surgery and anesthesia time
- Angled tip design helps you maintain a clear view of the incision line at all times
- Lightweight but durable with a comfortable ergonomic grip design
- Simple to use — simple to reload
- Designed to reduce problems with rotated staples

**WPSS01** Skin Stapler (one staple cartridge included)  
**WPSS02** Staple Cartridge (35 staples)



WPSS01



Stainless steel skin staples are 3.6 mm high by 6.9 mm wide after closure.

# Reusable Biopsy Punches with Plunger

## Features

- Individual Sterile Packs, Limited Reuse
- Tip material: 304 stainless steel
- Sterilization: ETO (ethylene oxide)

## Applications

- **Electrophysiology**—Patch recording (various regions of the specimen can be punched prior to separation).
- **Biomedical Pharmacology**—Collect samples from the brain neurotransmitters for analyzing the metabolic changes in response to active substances, or skin samples to analyze the response of anti aging medications.
- **Anatomy**—Sample brain regions for analysis of neurotransmitter concentrations of mRNA levels.
- **Forensic Sampling**—WPI punches can easily cut, retrieve, store and eject material like skin, gel, films, paint chips and paper. (Better



504529

samples are collected using one of our self-healing cutting mats **504620** or **504621**.)

- **Microfluidics**—Punch flow inlets in microfluidics chambers molded of PDMS (polydimethylsiloxane).

Surgeons and laboratory technicians use biopsy punches and curettes to remove tissue samples. These punches are for use in research laboratories.

<b>504638</b>	ID 0.35 mm, OD 0.63 mm, Wall 0.15 mm
<b>504528</b>	ID 0.50 mm, OD 0.80 mm, Wall 0.15 mm
<b>504529</b>	ID 0.75 mm, OD 1.07 mm, Wall 0.15 mm
<b>504646</b>	ID 1.0 mm, OD 1.26 mm, Wall 0.15 mm
<b>504530</b>	ID 1.2 mm, OD 1.5 mm, Wall 0.15 mm
<b>504647</b>	ID 1.5 mm, OD 1.9 mm, Wall 0.20 mm
<b>504531</b>	ID 2.0 mm, OD 2.4 mm, Wall 0.20 mm
<b>504648</b>	ID 2.5 mm, OD 2.9 mm, Wall 0.20 mm

<b>504649</b>	ID 3.0 mm, OD 3.4 mm, Wall 0.20 mm
<b>504650</b>	ID 3.5 mm, OD 3.9 mm, Wall 0.20 mm
<b>504651</b>	ID 4.0 mm, OD 4.4 mm, Wall 0.20 mm
<b>504652</b>	ID 5.0 mm, OD 5.5 mm, Wall 0.25 mm
<b>504653</b>	ID 6.0 mm, OD 6.5 mm, Wall 0.25 mm
<b>504654</b>	ID 7.0 mm, OD 7.5 mm, Wall 0.25 mm
<b>504655</b>	ID 8.0 mm, OD 8.5 mm, Wall 0.25 mm

## Reusable Rapid-Punch Biopsy Kit

The Rapid-Punch consists of a razor sharp cutting tip designed to cut, retrieve and store cored samples from source materials such as tissue, gels, paper, cloth, leaves, paint chips, films or other thin substrates. The Rapid-Punch is recommended as the special DNA sampling tool. It is ideal for forensic and most laboratory sampling purposes. Designed for heavy, numerous and precise sampling operations. It is durable, lasting up to 2000 punches.

<b>504639</b>	Reusable Rapid Punch Kit with 0.5 mm Tip
<b>WP1010</b>	Reusable Rapid Punch Kit with 1.0 mm Tip
<b>WP1212</b>	Reusable Rapid Punch Kit with 1.2 mm Tip
<b>WP2020F</b>	Reusable Rapid Punch Kit with 2.0 mm Tip
<b>WP3030</b>	Reusable Rapid Punch Kit with 3.0 mm Tip
<b>504640</b>	Replacement Tip for Rapid Punch size 0.5 mm
<b>504641</b>	Replacement Tip for Rapid Punch size 1.0 mm
<b>504642</b>	Replacement Tip for Rapid Punch size 1.2 mm
<b>504643</b>	Replacement Tip for Rapid Punch size 2.0 mm
<b>504644</b>	Replacement Tip for Rapid Punch size 3.0 mm
<b>504735</b>	Replacement Plunger for Rapid Punch, 0.5 mm
<b>504736</b>	Replacement Plunger for Rapid Punch, 1.0 mm
<b>504737</b>	Replacement Plunger for Rapid Punch, 1.2 mm
<b>504738</b>	Replacement Plunger for Rapid Punch, 2.0 mm
<b>504739</b>	Replacement Plunger for Rapid Punch, 3.0 mm



WP1010



Rapid-Punch with 0.5 mm tip

### WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com Germany: +49 (0)6031 1602171 • wpide@wpi-europe.com • www.wpi-europe.com  
Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • www.wpiinc.com China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net

# Clips and Clamps

Save time in surgery using our reflex clips. They are a fast, effective suture alternative for wound closure. Choose the 9 mm clips for use in rats or 7 mm clips for mice and young rats. The spring-loaded clip applicator is easy to use. You get minimal tissue trauma and maximum holding power. Use the clip remover to extract the clips when the wound heals. The clips, applicator and remover are all autoclavable.

Vessel clips are used to clamp arteries or veins during surgical procedures. These vascular clamps come in a variety of shapes and sizes.

## REFLEX CLIP APPLIER

- 11 cm (4.3 in.) Long

<b>500343</b>	For 7 mm Clips
<b>500345</b>	For 9 mm Clips



500343

## REFLEX CLIP

- 100 per box
- Non-sterile
- Stainless Steel

<b>500344</b>	7 mm
<b>500346</b>	9 mm



500344

## VESSEL CLIPS

- 5 per pack
- Stainless Steel

<b>15911</b>	0.8 x 5 mm Jaws, 10 g Pressure
<b>15913</b>	1 x 6 mm Jaws, 15 g Pressure
<b>14120</b>	1.5 x 5 mm Jaws, 30 g Pressure
<b>14121</b>	1.5 x 8 mm Jaws, 60 g Pressure



15911



14120



15913



14121

# Disposable Scalpels

- Carbon steel blades
- Non-slip handles
- Used for incisions and underlying tissue
- Packed 10 per box
- Sterile
- Individually wrapped

<b>500348</b>	#10
<b>500349</b>	#11
<b>500350</b>	#12
<b>500351</b>	#15
<b>500352</b>	#20
<b>500353</b>	#21
<b>500354</b>	#22
<b>500355</b>	#23



500348



500351



500355

# Rodent Accessory

## Features


- Large, stable base
- Hardened blades for long service
- Ambidextrous configuration
- 3 sizes
- Cuts through bone

## Benefits

- Blades are drawn together by magnetic force to ensure a clean and precise cut through very strong bones and skin
- Spring action, locking device prevents blades from accidentally falling
- Fluoropolymer coated base for easy cleanup

## Applications

- Decapitator




The small animal guillotine is easy to use with extra added safety features. There is a large base for stability, long handle for extra leverage, spring action so the blades can not fall down unexpectedly, hardened stainless blades for endurance, simplified construction for easy maintenance. The fluoropolymer coated surface on the base makes cleaning easy. The guillotine is considered one of the most humane methods to sacrifice a subject.

<b>DCAP</b>	For Rodents/Other Small Animals, 1.5 x 1.5 in. Opening
<b>DCAP-M</b>	For Large Rodents/Other Medium Animals, 2.5 x 2.5 in. Opening
<b>DCAP-L</b>	For Larger Animals, 4.0 x 4.0 in. Opening

# Rodent Brain Matrices

WPI offers a large selection of acrylic and stainless steel brain matrices. These matrices are sturdy and can be heated, chilled, autoclaved, scrubbed and stand up to rigorous daily use. Coronal matrices have the

additional feature of a mid-line sagittal cut to facilitate splitting of the left and right hemispheres. Sections can be as fine as 1 mm. The olfactory/spinal/notch is cut into each matrix.



## SPECIFICATIONS

Order #	Subject	Material	Section	A	B	C	D	E	Cavity Depth	Weight
<b>RBMA-200C</b>	Adult Mouse, 40-75 g	Acrylic	Coronal	3.18	11.1	8.73	19.1	12.2	7.4	0.5 lb.
<b>RBMA-200S</b>	Adult Mouse, 40-75 g	Acrylic	Sagittal	3.18	11.1	8.73	19.1	12.2	7.4	0.5 lb.
<b>RBMA-300C</b>	Rat, 175-300g	Acrylic	Coronal	4.76	15.9	12.7	36.6	23.8	7.61	0.5 lb.
<b>RBMA-300S</b>	Rat, 175-300g	Acrylic	Sagittal	4.76	15.9	12.7	36.6	23.8	7.61	0.5 lb.
<b>RBMA-600C</b>	Rat, 300g-600g	Acrylic	Coronal	4.76	19.8	14.7	36.6	24.7	10.91	0.5 lb.
<b>RBMA-600S</b>	Rat, 300g-600g	Acrylic	Sagittal	4.76	19.8	14.7	36.6	24.7	10.91	0.5 lb.
<b>RBMS-200C</b>	Adult Mouse	Stainless Steel	Coronal	3.18	11.1	8.73	19.1	12.2	7.4	1.0 lb.
<b>RBMS-200S</b>	Adult Mouse	Stainless Steel	Sagittal	3.18	11.1	8.73	19.1	12.2	7.4	1.0 lb.
<b>RBMS-300C</b>	Rat, 175-300g	Stainless Steel	Coronal	4.76	15.9	12.7	36.6	23.8	7.61	1.0 lb.
<b>RBMS-300S</b>	Rat, 175-300g	Stainless Steel	Sagittal	4.76	15.9	12.7	36.6	23.8	7.61	1.0 lb.
<b>RBMS-600C</b>	Rat, 300g-600g	Stainless Steel	Coronal	4.76	19.8	14.7	36.6	24.7	10.91	1.0 lb.
<b>RBMS-600S</b>	Rat, 300g-600g	Stainless Steel	Sagittal	4.76	19.8	14.7	36.6	24.7	10.91	1.0 lb.

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com  
 Germany: +49 (0)6031 1602171 • wpide@wpi-europe.com • www.wpi-europe.com  
 Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • www.wpiinc.com  
 China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net




# Surgical Kits

## PHYSIOLOGY KIT I

### Kit Includes:

- SuperCut Tenotomy Scissors, Curved (14396)
- Dumont Forceps (500342)
- Rongeur, 3 mm jaw (14091)
- Filter Forceps, 11 cm, Straight, Flat Jaw (500456)
- Utility Scissors (501322)
- Probe, 1.0 mm Diameter, Blunt (501313)
- Operating Scissors, Straight, Sharp/Blunt (14192)

**KIT-PHYSIO-I** Physiology Kit I (7 instruments)








## PHYSIOLOGY KIT II

### Kit Includes:

- Vannas Scissors (14003)
- Stevenson Retractor (14131)
- SuperCut Iris Scissors, Straight (14218)
- Iris Forceps, Curved, Serrated (15915)
- Rongeur, 1.3 mm Jaw (14292)
- Adson Forceps, 1x2 Teeth (500092)
- Utility Scissors, Straight, Sharp/Blunt (501322)
- Olsen-Hegar Needle Holder (500227)
- Probe, 1.0 mm Diameter, Blunt (501313)
- Dumont Forceps (500342)
- Operating Scissors, Straight, Sharp/Blunt (14192)
- Portfolio (501319)

**PHYSIO-II** Physiology Kit II (11 instruments plus zipper case)




## MOUSE DISSECTING KIT

### Kit Includes:

- Dumont Forceps (500342)
- Wire Retractor (14130)
- Vannas Scissors (14003)
- Needle Holder (14109)
- Iris Forceps, Curved, Serrated (15915)
- Probe, 1.0 mm Diameter, Blunt (501313)
- Dissecting Scissors, 10 cm, Straight (14393)
- Portfolio (503294)

**MOUSEKIT** Mouse Kit (7 instruments plus zipper case)  
**504741** Mouse Kit for Telemetry Training (includes **500453**)



14109




14393




500342



15915



14130



501313



Vessel Cannulation Forceps (500453)  
 Included only with Kit **504741**



503294



14003

## RAT DISSECTING KIT


### Kit Includes:

- Dumont Forceps 500342
- Alm Retractor 14240
- Vannas Scissors 14124
- Needle Holder 14110
- Iris Forceps, Curved, Serrated 15915
- Probe, 1.0 mm Diameter, Blunt 501313
- SuperCut Iris Scissors, Straight 14218

**RATKIT** Rat Kit (7 instruments)  
**504740** Rat Kit for Telemetry Training (includes **500453**)



15915



14110



14240



500342



501313



Vessel Cannulation Forceps (500453)  
 Included only with Kit **504740**



14218



14124

## STORAGE PORTFOLIO

**501319** Large Surgical Instrument Storage Portfolio, Holds up to 40 Instruments

**503294** Small Surgical Instrument Storage Portfolio, Holds up to 10 Instruments



503294



501319

### WORLD PRECISION INSTRUMENTS

# OmniDrill35 Micro Drill

*Drill, grind or finish bone or teeth*

## Features

- Excellent tool for grinding, finishing, and drilling bone and other material
- High-torque motor (35,000 rpm) is quiet
- Includes foot switch

## Benefits

- Removable nose cone that can be cleaned and sterilized
- Consistent power for the duration of its use
- Variable speeds so you control the amount of heat generated

## Applications

- Grinding, finishing, cutting and drilling bone, teeth and other material
- This line-powered micro drill will make easy work of grinding, finishing, cutting and drilling bone, teeth and other material. The high-torque 35,000 rpm (maximum) motor is quiet and has minimal vibration which reduces wear on the motor and provides greater comfort for you. It also features a forward and reverse switch, "E Type" handpiece and handpiece holder. The handpiece has a removable nose cone that can be cleaned and sterilized. It accepts 3/32" and 2.33mm bur shanks. Unlike battery-powered drills, this unit maintains consistent power for the duration of use. The wide range of speeds allows you to control the amount of heat generation.



503598

## ORDERING INFORMATION

**503598** OmniDrill35 Micro Drill System, 110 V

**503599** OmniDrill35 Micro Drill System, 220 V

*Stand, tips, cutoff disk, mandrels and foot switch are all included.  
Replacements available below.*

## OPTIONAL ACCESSORIES/REPLACEMENT PARTS

<b>501851</b>	Abrading Tip, Stone, pkg. of 5
<b>501852</b>	Accessory Stand
<b>501853</b>	Ball Mill, Carbide, #1, .031" Diameter, pkg. of 5
<b>501854</b>	Ball Mill, Carbide, #2, .039" Diameter, pkg. of 5
<b>501855</b>	Ball Mill, Carbide, #3, .047" Diameter, pkg. of 5
<b>501856</b>	Ball Mill, Carbide, #4, .055" Diameter, pkg. of 5
<b>501857</b>	Ball Mill, Carbide, #5, .063" Diameter, pkg. of 5
<b>501858</b>	Ball Mill, Carbide, #6, .071" Diameter, pkg. of 5
<b>501842</b>	Ball Mill, Carbide, #7, .083" Diameter, pkg. of 5
<b>501860</b>	Ball Mill, Carbide, #1/4, .019" Diameter, pkg. of 5
<b>501861</b>	Ball Mill, Carbide, #1/2, .027" Diameter, pkg. of 5
<b>501862</b>	Cutoff Disk, pkg. of 20
<b>501863</b>	Mandrel, Screw, pkg. of 5
<b>501864</b>	Mandrel, Threaded, pkg. of 5
<b>504459</b>	Foot switch
<b>502237</b>	Stereotaxic Holder for OmniDrill35 Microdrill

## OMNIDRILL35 SPECIFICATIONS

INPUT	110 V, 50/60 Hz ( <b>503598</b> ) 240 V, 50/60 Hz ( <b>503599</b> )
OUTPUT	0-32 VDC
FUSE	1 A
OPERATING SPEED	0-35,000 RPM
DIMENSIONS	178 x 114 x 89 mm (7 x 4.5 x 3.5in.)
WEIGHT	1.7 kg (3.75 lbs.)

# Economy Electrosurgical Unit

*Cut and coagulate with a press of the foot switch*

## Features

- 10 levels of output intensity
- Three operational modes (cut, coagulate, cut/coagulate)
- Choice of electrodes

## Benefits

- Ready to use out of the box.
- Includes a variety of electrodes, which are also sold individually

## Applications

- Cutting and coagulation of tissue

Electrosurgery utilizes alternating current at radio frequencies to cut and coagulate. Using this method, the current enters the subject's body and the subject becomes part of the electrical circuit. This requires the use of a return, or ground, plate. This economically priced electrosurgical



501274

unit has 10 levels of output intensity, three operational modes (cut, coagulate, and cut/coagulate) and various choices of electrodes. The unit comes complete and ready-to-use with a handpiece, ground plate, foot switch and one of each electrodes. All accessories can also be ordered separately.

## ELECTROSURGICAL UNIT SPECIFICATIONS

OPERATION FREQUENCY	1.5 MHz
STABLE & FINE POWER SETTING	10 Steps
POWER SUPPLY	115 V ± 10% - 50/60 Hz 1.8 A, 210 VA 230 V ± 10% - 50/60 Hz 0.9 A, 210 VA
OUTPUT POWER	70 W ± 5%
WORKING FREQUENCY	1.5-1.7 MHz ± 5%
DIMENSIONS	24 cm x 22 cm x 8.5 cm (lwxh)
SHIPPING WEIGHT	10 lb. (4.5 kg)

## ORDERING INFORMATION

**501274** Electrosurgical Unit, 110 V

**501284** Ball Electrode, Ø 1.6 mm shaft

# Dry Sterilizer

*Heat sterilize instruments in seconds*

## Features

- Large LED control panel
- Adjustable temperature control up to 300°C
- Sterilizes forceps, needles and other instruments
- Choose Bead capacity: 300 g (**ST5193**), 150 g (**ST5191**)

## Benefits

- No chemicals
- No flames
- No risk of burns
- No disinfectant fluids

## Applications

- Sterilize your microdissecting and tissue culture instruments, thoroughly and conveniently, in seconds

WPI's **Steri-Lite Micro Bead Sterilizer** is ideal for the sterilization of small research instruments (forceps, scissors, etc.). When the chamber is filled with the included glass beads, high temperatures (up to 300°C) can be used to eliminate bacteria, spores and other microorganisms.

Chamber temperature is displayed on the large LED control panel and may be adjusted with the control knob. The outer casing of the unit stays cool to the touch, even when it is used continually. The supplied safety cover can be placed over the chamber to conserve energy and reduce the heating/ramping rate. The Steri-Lite Sterilizer does not use gas, flame or chemicals and is safe for use in hoods.

With adjustable temperature, the sterilizer can be used for incubations and the decontamination of plastic objects/tools.

This product is designed for RESEARCH and DEVELOPMENT USE ONLY.




### STERI-LITE SPECIFICATIONS

TEMPERATURE RANGE:	100 to 300°C
PREHEAT TIME:	Approximately 25 min. to 300°C
CHAMBER DIAMETER:	40 mm
CHAMBER HEIGHT:	ST5191: 80 mm ST5193: 140 mm
DIMENSIONS:	ST5191: 14x14.5x15.5 cm (5.5x5.7x6.1 in.) ST5193: 14X14.5X21.5 cm (5.5X5.7X8.7 in.)
BEAD CAPACITY:	ST5191: 150 g ST5193: 300 g
WEIGHT:	ST5191: 2 Kg ST5193: 2.5 Kg
ELECTRICAL:	ST5191: 115 V or 230 V, 50-60 Hz, 120 W ST5193: 115 V or 230 V, 50-60 Hz, 250 W
WARRANTY:	2 Years

### ORDERING INFORMATION

<b>ST5193</b>	Steri-Lite Micro Bead Sterilizer, large, 300 g capacity
<b>ST5191</b>	Steri-Lite Micro Bead Sterilizer, small, 150 g capacity



*Request the surgical instrument catalog to see the full line of instruments:*

[www.wpiinc.com/catalog](http://www.wpiinc.com/catalog)

#### WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • [wpiuk@wpi-europe.com](mailto:wpiuk@wpi-europe.com) • [www.wpi-europe.com](http://www.wpi-europe.com) Germany: +49 (0)6031 1602171 • [wpide@wpi-europe.com](mailto:wpide@wpi-europe.com) • [www.wpi-europe.com](http://www.wpi-europe.com)  
Brazil: 011 55 13 40629703 • [info@brazil.wpiinc.com](mailto:info@brazil.wpiinc.com) • [www.wpiinc.com](http://www.wpiinc.com) China: +86 21 6888 5517 • [chinasales@china.wpiinc.com](mailto:chinasales@china.wpiinc.com) • [www.wpiinc.net](http://www.wpiinc.net)

# Ultrasonic Cleaning Systems

## Quantrex® bath with a timer

### Features

- 60 Minute Timer
- Vinyl-clad steel and stainless steel
- 14 quality inspection steps for strength and durability

### Benefits

- Provides super-strength cleaning every time
- When used with L&R's specialty formulated solutions, the self-contained, compact unit offers efficient trouble-free cleaning
- Each Quantrex machine comes standard with increased power—strength you can see as soon as you turn the unit on
- Stainless steel drain with multi-positional outlet for easy removal of solution

### Applications

- Versatile enough for a variety of cleaning applications



504216

### SPECIFICATIONS

DIMENSIONS	26.0 x 16.5 x 21.0 cm (10.25 x 6.5 x 8.25 in.)
SHIPPING WEIGHT	4.5 kg (10 lb.)
TANK CAPACITY	3.2 L (0.85 gal.)
INTERNAL DIMENSIONS	23.8 x 13.7 x 10.2 cm (9.38 x 5.38 x 4.0 in.)

### ORDERING INFORMATION

504216	Quantrex Ultrasonic Cleaning System <i>Includes Timer, Drain and Cover. Heater Optional.</i>
--------	---

## Economy ultrasonic cleaner



UBATH-Y

### Features

- Half liter stainless steel tank
- Durable and compact

### Benefits

- Robust all-metal construction allows for continuous duty

### Applications

- Microdissection and veterinary instrument cleaning

### UBATH SPECIFICATIONS

INPUT POWER	22W
PEAK OUTPUT	70W, 55 kHz
TANK CAPACITY	0.53L (18 oz.)
TANK I.D.	12.1 x 8.6 x 6.6 cm, (4.75 x 3.375 x 2.625 in.)
TANK O.D.	13.7 x 10.5 x 12.1 cm, (5.375 x 4.125 x 4.75 in.)
SHIPPING WEIGHT	1.8 kg (4 lb.)



504766

504767

504768

### ORDERING INFORMATION

UBATH-Y	Ultrasonic Cleaner (110 V)
13740	Ultrasonic Detergent (4 lb)
504766	Sterilization mesh cassette 40x40x20 mm
504767	Sterilization mesh cassette 80x80x34mm
504768	Sterilization mesh cassette 105x70x25 mm



*Request a complete catalog of sterilization baskets and trays:*

[www.wpiinc.com/catalog](http://www.wpiinc.com/catalog)

# Lab Supplies



## Affordable laboratory equipment

Looking for a quality source for affordable laboratory equipment? From capillary glass and adhesives up to micromanipulators, pumps and microscopes, we offer the full range for the well stocked scientific laboratory.

### WORLD PRECISION INSTRUMENTS

**UK:** +44 (0)1462 424700 • [WPIUK@WPI-EUROPE.COM](mailto:WPIUK@WPI-EUROPE.COM) • [WWW.WPI-EUROPE.COM](http://WWW.WPI-EUROPE.COM)  
**BRAZIL:** (013) 406-29703 • [INFO@BRAZIL.WPIINC.COM](mailto:INFO@BRAZIL.WPIINC.COM) • [WWW.WPIINC.COM](http://WWW.WPIINC.COM)

**GERMANY:** +49 (0)6031 1602171 • [WPIDE@WPI-EUROPE.COM](mailto:WPIDE@WPI-EUROPE.COM) • [WWW.WPI-EUROPE.COM](http://WWW.WPI-EUROPE.COM)  
**CHINA:** +86 21 6888 5517 • [CHINASALES@CHINA.WPIINC.COM](mailto:CHINASALES@CHINA.WPIINC.COM) • [WWW.WPIINC.NET](http://WWW.WPIINC.NET)

# Which Adhesive is Right for Me?

## ADHESIVES APPLICATION GUIDE

WPI Part #	Description	Curing time	Useful Applications and Characteristics
<b>Epoxies</b> Form strong bonding. Used in wire bonding applications.			
<b>4898</b>	Silver filled conductive Epoxy	12 hr @ 50°C; 5 min @ 150°C	Connecting conductors that can't be soldered. Constructing or connecting Ag/AgCl pellets.
<b>7335</b>	Carbon filled conductive Epoxy	48 hr @ 25°C; 5 min @ 150°C	Constructing carbon electrode.
<b>4886</b>	High performance Structural Epoxy	12 hr @ 25°C.	Forms a strong and slightly flexible bond on plastic, metal, & glass. Bonds some low surface.
<b>Hot Melt (EVA)</b> Easy to use for bonding. For large gap filling.			
<b>13316</b>	Mini Glue Gun with glue sticks	As soon as it cools down	Bonds wood, glass, metals, and many plastics.
<b>Silicone Adhesives/Sealants/Primers</b> Good moisture resistance and elasticity. Low toxicity.			
<b>1571</b>	Room temperature vulcanizing (RTV) adhesive. Acyloxy/moisture cure system. Acetic acid is cure by-product.	24 hr @ 25°C	Has the best adhesion property in this silicone family. Will bond to many materials.
<b>7128</b>	RTV sealant. Alkoxy/Moisture cure system. Methanol as cure by-product.	72 hr @ 25°C	Good for bonding or sealing electronics circuits (metal).
<b>SYLG184</b>	Sylgard, two parts, vinyl/platinum cure sealant. Hydrogen as cure by-products. Very low toxic	24 hr @ 25°C, 15 min. @ 150°C	Coating Patch Clamp electrodes, Cell culture dish, making dissection pads.
<b>KWIK-SIL</b>	Two part, adhesive. Vinyl/platinum system, Hydrogen as cure by-products. Very low toxic.	< 5 min. @ 25°C	Live tissue and nerve studies. Medium strength adhesion.
<b>KWIK-CAST</b>	Two part sealant. Vinyl/platinum cure system. Hydrogen as cure by-products. Very low toxic.	< 5 min. @ 25°C	Sealant for live tissues. Embedding peripheral nerves with electrodes.
<b>6820</b>	Primer for silicone	N/A	Enhances adhesion of silicone adhesives for difficult to bond plastic surfaces
<b>Cyanoacrylate</b> Forms an instantaneous bond.			
<b>7341</b>	Ethyl Cyanoacrylate, low viscosity 90-120 cps	<10 seconds	Mounting rat/mouse brain slices. Ideal for relatively small gaps and smooth surfaces. Bonds plastic, metals and rubber. Package of 10 vials, each approximately 1.5 mL.
<b>7342</b>	Ethyl Cyanoacrylate, high viscosity 1100-1600 cps	<30 seconds	Use on brain slice exp. Ideal for larger gaps, allows slightly longer bonding time. Bonds plastic, metals and rubber. Package of 10 vials, each approximately 1.5 mL.
<b>VETBOND</b>	Butyl Cyanoacrylate, low toxic	<10 seconds	Bonds tissues, alternative to suture, helps small wound healing. Antimicrobial effect. Used in forensic science.
<b>503763</b>	Octyl Cyanoacrylate, low toxic	<15 seconds	Suitable for surface wound bonding, protection, holding a sensor or other device on the tissue.

## Easy Application Silicone Structural Epoxy

**Kwik-Gard™** is specially packaged **Sylgard 184** silicone for quicker and easier application, eliminating the messy procedure of preparing the mixture before application. Its special cartridge controls the precise mixing ratio to ensure proper curing. The disposable tip mixes resin and hardener as they are dispensed. Since no air is introduced during mixing, the resin does not need degassing for most applications. The mixed silicone is applied directly to the site, reducing preparation time and material waste.



Each Kwik-Gard cartridge contains 37 mL of resin and hardener. The dispensing tip has a dead volume of 0.75 mL. Perfect for PDMS.

Biopsy punches (0.5–8mm) are available online at [www.wpiinc.com/punch](http://www.wpiinc.com/punch).

### ORDERING INFORMATION

<b>KWIKGARD</b>	Kwik-Gard Start-up Kit (incl. dispenser, 1 cartridge, 5 tips)
<b>KWIKGLUE</b>	Kwik-Gard Refill (2 cartridges, 10 dispensing tips)
<b>KWIKMIX</b>	Dispensing Tips (pkg. of 10)
<b>KWIKGUN</b>	Kwik-Gard Dispenser

## Easy Application Silicone Structural Epoxy

Scotch-Weld 2216 remains the best epoxy for bonding plastic, often used as the benchmark for testing the binding strength of other adhesives. The slightly rubbery texture also makes it less easy to break off. *It is the only epoxy known that can bond PEEK.*

Color: gray  
Cures at room temperature.  
Shipping weight: 1 lb. (0.5 kg)



### ORDERING INFORMATION

<b>4886</b>	Scotch-Weld 2216 (2 oz.)
-------------	--------------------------

# Low Toxicity Adhesive

Ideal for neuroscience applications, nerve studies and more

## Features

- Bio-compatible adhesive for live tissue and nerve studies
- Pre-mixing tips simplify use
- Medium strength adhesion
- Low toxicity
- Rapid curing silicone adhesive, cure on contact
- Cures without producing heat
- Includes 10 Mixing Tips
- Volume discounts - Save up to 15%!

## Benefits

- Low toxicity
- Rapid cure time

## Applications

- Neuroscience and nerve studies
- Biomedical applications



600022

**Kwik-Sil** and **Kwik-Cast** silicones have very low toxicity before, during and after curing. The by-product of curing is a small amount of hydrogen gas, which is much less toxic to cells than acetic acid or alcohol from traditional RTV silicone systems.

**Kwik-Sil** and **Kwik-Cast** curing speed is hundreds of times faster than traditional RTV silicones. A curing time of a few minutes at room temperature is especially useful for encapsulation of live tissue or implanting into a live animal.

Unlike many vinyl-based silicones in which the platinum complex catalysts are easily poisoned by contamination from amines and animal tissue, **Kwik-Sil** and **Kwik-Cast** are not sensitive to contamination from animal tissue.

**Kwik-Sil™** is a translucent, medium-viscosity silicone adhesive, developed for chronic peripheral nerve studies such as anterograde tracing with fluorescent indicators or electrode recording. Good adhesion and mechanical properties (tear strength and elongation) allow days of study without breaking of the bonding. Curing speed is very reproducible.

**Kwik-Cast™** is a very low viscosity silicone sealant developed to embed peripheral nerves with electrodes for acute multi-fiber recordings. It



KWIK-SIL



KWIK-CAST

flows easily, filling the small spaces around the nerve and leaving no channels through which peritoneal fluid can travel and thus short the nerve/electrode contact. Equally important is the ability of the material to flow into itself and create one continuous mass from underneath the nerve all the way to the top of the nerve/electrode contact to ensure long-term recording stability. **Kwik-Cast** is color-coded to make the mixing foolproof. The catalyst is yellow and the base is blue. When uniformly mixed, it is green. **Kwik-Cast** can be applied and cured underneath mineral oil. After recording, electrodes are easily recovered due to the low tear strength.

## KWIK-CAST & KWIK-SIL SPECIFICATIONS

	Kwik-Sil	Kwik-Cast
MIX RATIO	1 to 1	1 to 1
WORKING TIME	< 5 minutes*	4 minutes
SETTING TIME (ROOM TEMP., 1:1 RATIO)	5-10 minutes**	<10 minutes
CURE TIME	~15 minutes	
VISCOSITY, CPS	15,000	10,000
GUARANTEED SHELF LIFE AT 23 °C	6 months	6 months
VOLUME	5 mL	5 mL
NUMBER OF MIXING TIP	10	10
DEAD VOLUME OF THE MIXING TIP	<0.12 mL	<0.12 mL
AFTER CURING 24 HOURS:		
TEAR STRENGTH, PPI	90	44
ELONGATION %	650	60
DUROMETER (SHORE A-2)	30	36
COLOR	translucent	green
VOLUME RESISTIVITY, W/CM	1x10 <sup>15</sup>	1x10 <sup>15</sup>

\* 3 minutes average with about 90 seconds of liquidity

\*\* no longer mixable at this point

## ORDERING INFORMATION

<b>KWIK-SIL</b>	Silicone Adhesive Compound (two 5 mL syringes)
<b>KWIK-CAST</b>	Silicone Casting Compound (two 5 mL syringes)
<b>600022</b>	Replacement KWIK Mixing Tips (pkg. of 10)

Quantity discounts available

# MicroFil™ Non-metallic Syringe Needle

Perfect for filling micropipettes

WPI's **MicroFil™** fills micropipettes easily and reliably. Its long and fine tip allows you to start the filling very close to the pipette tip, eliminating both air bubble formation and clogging due to the washing down of dust particles. The transparent amber **MicroFil™** needle is constructed from a combination of plastic and fused silica — no metal components are used. The **MicroFil™** needle can be stored for days with the filling solution inside without clogging. The **MicroFil's** tip elasticity is sturdy and very flexible though not unbreakable. Since it is more flexible than stainless steel needles, moderate bending will not block or damage the **MicroFil™** needle. The combination of plastic and fused silica in the **MicroFil™** tip is sturdier than plastic tips, allowing easy and repeated insertions into micropipettes. **MicroFil's** luer fitting allows easy coupling to syringes and syringe filters.



U.S. PATENT  
5,294,325

## ORDERING INFORMATION

<b>MF34G-5</b>	MicroFil™, 34 ga., 67 mm long (pkg. of 5)
<b>MF28G-5</b>	MicroFil™, 28 ga., 97 mm long (pkg. of 5)
<b>MF28G67-5</b>	MicroFil™, 28 ga., 67 mm long (pkg. of 5)

## Custom MicroFil™

All **MicroFil™** products, including custom orders, can be shipped immediately. Custom orders for special needs can be made using nine sizes of **MicroFil™** tubing in lengths up to 50 cm — except for **CMF90U** which has a maximum length of 10 cm because of its high resistance to flow. Quantity discounts available. Specify length when ordering.

<b>CMF20G</b>	MicroFil™, 20 ga, 700 µm ID, 850 µm OD (pkg. of 4)
<b>CMF22G</b>	MicroFil™, 22 ga, 530 µm ID, 700 µm OD (pkg. of 4)
<b>CMF23G</b>	MicroFil™, 23 ga, 530 µm ID, 665 µm OD (pkg. of 4)
<b>CMF26G</b>	MicroFil™, 26 ga, 320 µm ID, 430 µm OD (pkg. of 4)
<b>CMF28G</b>	MicroFil™, 28 ga, 250 µm ID, 350 µm OD (pkg. of 4)
<b>CMF31G</b>	MicroFil™, 31 ga, 100 µm ID, 238 µm OD (pkg. of 4)
<b>CMF34G</b>	MicroFil™, 34 ga, 100 µm ID, 164 µm OD (pkg. of 4)
<b>CMF35G</b>	MicroFil™, 35 ga, 75 µm ID, 144 µm OD (pkg. of 4)
<b>CMF90U</b>	MicroFil™, ~36 ga, 20 µm ID, 90 µm OD (pkg. of 4)

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • WPIUK@WPI-EUROPE.COM • WWW.WPI-EUROPE.COM  
BRAZIL: (013) 406-29703 • INFO@BRAZIL.WPIINC.COM • WWW.WPIINC.COM

GERMANY: +49 (0)6031 1602171 • WPIDE@WPI-EUROPE.COM • WWW.WPI-EUROPE.COM  
CHINA: +86 21 6888 5517 • CHINASALES@CHINA.WPIINC.COM • WWW.WPIINC.NET



# "Super" Adhesives for Life Science Research

*Four times stronger than butyl cyanoacrylate and less toxic*

Cyanoacrylate adhesives have been on the market since 1958. Most industrial or household grade cyanoacrylate is made of shorter alkyl chain derivatives such as methyl or ethyl cyanoacrylate (**7341** and **7342**). They are very useful for temporarily holding tissues such as mounting specimens for microtome sectioning. However, they are not suitable for bonding wounds on live animals. The difficulties of using cyanoacrylate for bonding live animals include:

- A strong, irritating odor
- Quick loss of bonding strength due to breakdown of the bonding by hydration
- Breakdown products (cyanoacetate and formaldehyde) are toxic and can cause inflammatory reactions
- Low flexibility and tend to be brittle.

To overcome these problems, several longer alkyl chain cyanoacrylates have been developed especially for veterinary and human use. The first longer alkyl chain product is butyl cyanoacrylate. This product has been used for animal and human applications outside the USA since 1970. It is much less toxic and has a lower odor than the methyl or ethyl cyanoacrylate. The butyl cyanoacrylate offered by WPI is **Vetbond™**.

A family of adhesives containing octyl cyanoacrylate, a plasticizer and stabilizer, was developed in the 1990's (one of them approved by FDA for human use). When bonding to tissue, these new adhesives are four times stronger and less toxic than butyl cyanoacrylate. Compared with the traditional suture, the new super adhesive has several advantages. On average, it takes only one-tenth of the time to close an incision.



**VETBOND**

The bonding strength is equal to 5-0 monofilament suture. It also has a mysterious antimicrobial effect that can decrease infection rates in contaminated wounds. Bonding will slough off naturally in 5 to 7 days. Cosmetic appearance of the healed incision is also better.

Glutire Topical Tissue Adhesive **503763** forms a strong and flexible film and is thus more suitable for surface wound bonding, protection, and holding a sensor or other device on the tissue. Setting time is about 10 seconds, which gives ample time for application. It can also be used for temporarily holding a live tissue. For example, there is a report of using it to hold nematodes on a glass slide for patch-clamp neurons recording.

All of the products offered by WPI are veterinary grade (not suitable for human application). Though very similar to the grade for human use, they are not sterile and do not have FDA approval.



**503763**

## ORDERING INFORMATION

<b>503763</b>	Glutire Topical Tissue Adhesive (10 tips), 1.5 mL
<b>7341</b>	Cyanoacrylate Adhesive, Low Viscosity—90-120 cps (package of 10 vials, each approximately 1.5 mL)
<b>7342</b>	Cyanoacrylate Adhesive, High Viscosity—1100-1600 cps (pkg. of 10 vials, each approximately 1.5 mL)
<b>VETBOND</b>	3M Vetbond™ Adhesive (3 mL)

## Sylgard



**SYLG184**

A two-part silicone elastomer, ideal for potting and encapsulating applications. Very low dielectric constant sealing compound used in patch clamping and many other lab applications. After cure, will withstand -55° to 200 °C.

## ORDERING INFORMATION

**SYLG184** Sylgard, 0.5 kg (1.1 lb)

## Silicone Dissecting Pad Kit



**501986**

Make your own silicone dissection pads easily and quickly. Mix the 2-part silicone right in the plastic petri dishes and allow to cure 24 hours at room temperature. Kit includes enough 2-Part Sylgard silicone elastomer to prepare 20 dishes, pins and 20 plastic petri dishes with lids, 65 mm.

## ORDERING INFORMATION

**501986** Silicone Dissecting Pad Kit

## Electrically Conductive Silver Epoxy



Two-component silver-filled epoxy for electrical connections which cannot be soldered, such as Ag/AgCl pellets. This widely used silver-filled epoxy features low viscosity and smooth flowing character. Pure silver is dispersed in both resin and hardener. Cures in 15 minutes at 120 °C. Mix ratio 1:1. May be premixed and frozen for later use.

### ORDERING INFORMATION

**4898** Silver Epoxy, 28 g (1 oz)

## Electrically Conductive Carbon Epoxy



Two-component carbon-epoxy, curable at room and elevated temperatures. Ideal for electrostatic discharge protection and EMI/RFI shielding. 1:1 mix ratio. May be premixed and frozen for later use.

### ORDERING INFORMATION

**7335** Carbon Epoxy, 56 g (2 oz)

## Silicone RTV Adhesive (non-acidic)



Because it is non-corrosive, this material is ideal for use on metal, for encapsulating small circuits on connectors. After cure, will withstand -55° to 200 °C. No mixing required.

### ORDERING INFORMATION

**7128** RTV Coating, 90 mL (3 fl oz)

## Silicone RTV Adhesive



Clear silicone sealant provides good bonding to plastic. After cure, will withstand -55 to 200 °C. No mixing required. A handy, general purpose laboratory sealant. (Releases acetic acid during curing.)

### ORDERING INFORMATION

**1571** RTV Sealant, 139 mL (4.7 fl oz)



## RTV Prime Coat

Enhances adhesion of silicone adhesives to many difficult-to-bond plastic surfaces.

### ORDERING INFORMATION

**6820** RTV Prime Coat, 400 mL (13.5 fl oz)



## Mini Glue Gun

Comes with three sticks of special formula hot melt glue. UL approved. 110V 60 Hz only.

### ORDERING INFORMATION

**13316** Mini Glue Gun

### WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • WPIUK@WPI-EUROPE.COM • WWW.WPI-EUROPE.COM  
BRAZIL: (013) 406-29703 • INFO@BRAZIL.WPIINC.COM • WWW.WPIINC.COM

GERMANY: +49 (0)6031 1602171 • WPIDE@WPI-EUROPE.COM • WWW.WPI-EUROPE.COM  
CHINA: +86 21 6888 5517 • CHINASALES@CHINA.WPIINC.COM • WWW.WPIINC.NET

# Glass Capillaries

Quality glass, superior prices for microinjection/microelectrodes

## Features

- Quality borosilicate glass capillaries
- Large variety available, including fire polished, filaments, thin wall, specialty glass and multi-barrel

## Benefits

- Superior pricing
- Most glass orders ship within 48 hours

## Applications

- Microinjection
- Electrophysiology
- Patch clamp
- Fluid Handling

## Fire Polishing

Fire-Polished glass capillaries are easier to insert into microelectrode holders without damaging the gasket. More importantly, fire-polished glass won't scratch the chloridized wire used in a recording electrode. Fire-polishing does not affect the glass's mechanical or electrical properties.



## Making Uniform, Reproducible Microelectrodes



Borosilicate glass capillaries: Close dimensional tolerances assure microelectrode uniformity and reproducibility. Capillaries are available in 1, 2, 3, 5 and 7-barrel configurations, complete range of single barrel thin-wall sizes and a variety of special configurations. Capillaries with filaments contain a solid filament fused to the inner wall, which speeds filling of electrodes. Capillaries with or without inner filaments are available for making microelectrodes in a wide range of diameters.

## Filament Glass Capillaries


Single Barrel standard wall thickness capillaries are offered either with or without inner filaments for quick filling in a variety of lengths and diameters.

## Thin Wall Glass Capillaries

Thin Wall single barrel capillaries are offered both with or without inner filaments.



Specialty glass is also available. See page 105.



ORDERING INFORMATION							
	Length	OD(mm)	ID(mm)	Filament	Fire-Polished	Quantity	Item
Single Barrel Standard Borosilicate Glass	3 in. (76 mm)	1.0	0.58	✓		500	<b>1B100F-3</b>
	3 in. (76 mm)	1.0	0.58			500	<b>1B100-3</b>
	3 in. (76 mm)	1.2	0.68	✓		350	<b>1B120F-3</b>
	3 in. (76 mm)	1.2	0.68			350	<b>1B120-3</b>
	3 in. (76 mm)	1.5	0.84	✓		225	<b>1B150F-3</b>
	3 in. (76 mm)	1.5	0.84		✓	300	<b>1B150-3</b>
	4 in. (100 mm)	1.0	0.58	✓	✓	500	<b>1B100F-4</b>
	4 in. (100 mm)	1.0	0.58		✓	500	<b>1B100-4</b>
	4 in. (100 mm)	1.2	0.68	✓	✓	400	<b>1B120F-4</b>
	4 in. (100 mm)	1.2	0.68			350	<b>1B120-4</b>
	4 in. (100 mm)	1.5	0.84	✓	✓	300	<b>1B150F-4</b>
	4 in. (100 mm)	1.5	0.84		✓	300	<b>1B150-4</b>
	4 in. (100 mm)	2.0	1.12	✓		125	<b>1B200F-4</b>
	4 in. (100 mm)	2.0	1.12		✓	200	<b>1B200-4</b>
	6 in. (152 mm)	1.0	0.58	✓		500	<b>1B100F-6</b>
6 in. (152 mm)	1.0	0.58			500	<b>1B100-6</b>	
6 in. (152 mm)	1.2	0.68	✓		350	<b>1B120F-6</b>	
6 in. (152 mm)	1.2	0.68			350	<b>1B120-6</b>	
6 in. (152 mm)	1.5	0.84	✓		225	<b>1B150F-6</b>	
6 in. (152 mm)	1.5	0.84			225	<b>1B150-6</b>	
6 in. (152 mm)	2.0	1.12	✓		125	<b>1B200F-6</b>	
6 in. (152 mm)	2.0	1.12			125	<b>1B200-6</b>	
Thin-Wall Single-Barrel Standard	3 in. (76 mm)	1.0	0.75	✓		500	<b>TW100F-3</b>
	3 in. (76 mm)	1.0	0.75			500	<b>TW100-3</b>
	3 in. (76 mm)	1.2	0.90	✓	✓	400	<b>TW120F-3</b>
	3 in. (76 mm)	1.2	0.90			350	<b>TW120-3</b>
	3 in. (76 mm)	1.5	1.12	✓		225	<b>TW150F-3</b>
	3 in. (76 mm)	1.5	1.12		✓	300	<b>TW150-3</b>
	4 in. (100 mm)	1.0	0.75	✓		500	<b>TW100F-4</b>
	4 in. (100 mm)	1.0	0.75		✓	500	<b>TW100-4</b>
	4 in. (100 mm)	1.2	0.90	✓		350	<b>TW120F-4</b>
	4 in. (100 mm)	1.2	0.90			350	<b>TW120-4</b>
	4 in. (100 mm)	1.5	1.12	✓		225	<b>TW150F-4</b>
	4 in. (100 mm)	1.5	1.12		✓	300	<b>TW150-4</b>
	6 in. (152 mm)	1.0	0.75	✓		500	<b>TW100F-6</b>
	6 in. (152 mm)	1.0	0.75		✓	500	<b>TW100-6</b>
	6 in. (152 mm)	1.2	0.90	✓	✓	400	<b>TW120F-6</b>
6 in. (152 mm)	1.2	0.90			350	<b>TW120-6</b>	
6 in. (152 mm)	1.5	1.12	✓		225	<b>TW150F-6</b>	
6 in. (152 mm)	1.5	1.12		✓	300	<b>TW150-6</b>	

Single barrel glass is Kimble N51A. All thin wall glass is Schott Duran 8330.  
Tolerance ±10%

# Borosilicate Glass Micropipettes

Eliminate the cost and trouble of making your own micropipettes

*µTip*™



## Features

- Schott Duran borosilicate glass
- 0.5 µm and smaller ID micropipettes include internal glass fiber for easy filling
- Tip inner diameter tolerance ±20%
- Short taper yields high strength
- Nominal length ≈ 50 mm
- OD:ID = 1.33:1
- Standard capillary outer diameters are 1.0 mm (thin-wall) or 1.14 mm
- Every pipette individually tested and inspected
- Vacuum sealed packs of 10

## Benefits

- Plain Shank or Luer Fittings

## Applications

- Injection of dyes or proteins into cells, oocytes or other biomedical laboratory applications

WPI can quickly supply your need for consistently sized pre-pulled glass micropipettes.

Tip diameters (ID) range from 0.1 to 10 µm.

Silanized Tips (Luer Shank) are available. Silanization waterproofs the glass to retard water when inserting into cell. This will not let the outside fluid run down the pipette and get inside so easily.

## ORDERING INFORMATION

Shank	Tip I.D.	Silanize Coating Length	Glass O.D.	Filament	Fire Polished	Catalog #
PLAIN	0.1 µm	—	1.0 mm Thin-Wall	Yes	No	TIP01TW1F
	0.2 µm	—	1.0 mm Thin-Wall	Yes	No	TIP02TW1F
	0.3 µm	—	1.0 mm Thin-Wall	Yes	No	TIP03TW1F
	0.4 µm	—	1.0 mm Thin-Wall	Yes	No	TIP04TW1F
	0.5 µm	—	1.0 mm Thin-Wall	Yes	No	TIP05TW1F
	1 µm	—	1.0 mm Thin-Wall	No	Yes	TIP1TW1
	2 µm	—	1.0 mm Thin-Wall	No	Yes	TIP2TW1
	5 µm	—	1.0 mm Thin-Wall	No	Yes	TIP5TW1
	10 µm	—	1.0 mm Thin-Wall	No	Yes	TIP10TW1
	10 µm	—	1.14 mm A203W glass *	No	Yes	TIP10XV119
—	—	—	1.0 mm Thin-Wall	No	Yes	TIP10LT Long Taper †
	—	—	1.0 mm Thin-Wall	Yes	Yes	TIP10FLT Long Taper †
	—	—	1.5 mm Thin-Wall	No	Yes	TIP15LT Long Taper †
	—	—	1.5 mm Thin-Wall	Yes	Yes	TIP15FLT Long Taper †
	—	—	1.5 mm Thin-Wall	Yes	Yes	TIP15FLT Long Taper †
	—	—	1.0 mm Thin-Wall	No	Yes	TIP30TW1
LUER	0.1 µm	—	1.0 mm Thin-Wall	Yes	—	TIP01TW1F-L
	0.2 µm	—	1.0 mm Thin-Wall	Yes	—	TIP02TW1F-L
	0.3 µm	—	1.0 mm Thin-Wall	Yes	—	TIP03TW1F-L
	0.5 µm	—	1.0 mm Thin-Wall	Yes	—	TIP05TW1F-L
	1 µm	—	1.0 mm Thin-Wall	No	—	TIP1TW1-L
	2 µm	—	1.0 mm Thin-Wall	No	—	TIP2TW1-L
	5 µm	—	1.0 mm Thin-Wall	No	—	TIP5TW1-L
	10 µm	—	1.0 mm Thin-Wall	No	—	TIP10TW1-L
	30 µm	—	1.0 mm Thin-Wall	No	—	TIP30TW1-L
	LUER/ SILANIZED	5 µm	1 inch	1.0 mm Thin-Wall	No	—
5 µm		2 inch	1.0 mm Thin-Wall	No	—	TIP5TW1LS02
10 µm		1 inch	1.0 mm Thin-Wall	No	—	TIP10TW1LS01
10 µm		2 inch	1.0 mm Thin-Wall	No	—	TIP10TW1LS02
30 µm		1 inch	1.0 mm Thin-Wall	No	—	TIP30TW1LS01
30 µm		2 inch	1.0 mm Thin-Wall	No	—	TIP30TW1LS02

\* 10 µm (ID), 1.14 mm capillary pipettes are for use in WPI's Nanoliter 2010.


† Long Taper micropipettes are pulled with a 12-15 mm taper which the customer cuts back to obtain the desired tip diameter.

## µTip Sampler Assortments

TIPMIX01-05	Two each, 0.1, 0.2, 0.3, 0.4, 0.5 µm ID, Plain Shank
TIPMIX05-10	Two each, 0.5, 1, 2, 5, 10 µm ID, Plain Shank
TIPMIX01-05-L	Two each, 0.1, 0.2, 0.3, 0.4, 0.5 µm ID, Luer
TIPMIX05-10-L	Two each, 0.5, 1, 2, 5, 10 µm ID, Luer

## Slides

These clean glass microscope slides are 25 x 75 mm, 1.0-1.2 mm thick with 90° ground edges. They are available as frosted and red ended. The frosted end slides feature a fine 20 mm frosted area on both sides of one end for easy marking. The red frosted slides feature a 20 mm colored end useful for identifying hazardous materials.



## ORDERING INFORMATION

503506	Frosted Glass Microscope Slides, Box of 144
503507	Red Frosted Glass Microscope Slides, Box of 144

## Cover Slips

These cover slips (made of German glass) can be used for growing and culturing cells that normally have poor adhesion to plastic surfaces. They are small enough to be placed in the micro plate or other cell culture devices. The 5 mm size will fit inside the 96-well culture plate and leave enough room to pick it up from the bottom of the well with forceps. The 8mm size fits inside the 24-well plates.



## ORDERING INFORMATION

Order #	Diameter	Thickness	Quantity
502040	5 mm	#1.5 (0.16 - 0.19 mm)	100
502041	8 mm	#1.5 (0.16 - 0.19 mm)	100
503508	25 mm	#1.5 (0.16 - 0.19 mm)	100

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • WPIUK@WPI-EUROPE.COM • WWW.WPI-EUROPE.COM  
BRAZIL: (013) 406-29703 • INFO@BRAZIL.WPIINC.COM • WWW.WPIINC.COM

GERMANY: +49 (0)6031 1602171 • WPIDE@WPI-EUROPE.COM • WWW.WPI-EUROPE.COM  
CHINA: +86 21 6888 5517 • CHINASALES@CHINA.WPIINC.COM • WWW.WPIINC.NET

# Luer Valve Assortment Kit

*Build your own liquid flow experiment*

## Features

- Over 300 assorted parts
- Luer fittings for quick and easy connect and disconnect

## Benefits

- Sold individually or in kits

## Applications

- Liquid flow experimental setups

A useful kit (right) for building your own liquid flow experiment. It provides the means to start, stop, add, divide and control a flow of liquid or gas. Included in the kit are **over 200 assorted parts** such as one-way and three-way stopcocks, manifolds, Y-connectors, injection sites, male and female luer caps, check valves, syringe-activated check valves, slide clamps, roller clamps, and pinch clamps. All (except clamps) have a luer fitting for quick and easy connecting and disconnecting. Includes assorted luer fittings for use with flexible tubing.



14011

## ORDERING INFORMATION

**14011** Luer Valve Assortment Kit

# Luer-to-Tubing Coupler Assortment Kit

*Quick connects in nylon and polypropylene*

## Features

- Over 250 assorted parts in each kit
- Valves are polycarbonate, and the valve handles are polyethylene. Do not autoclave those parts.

## Benefits

- Polypropylene parts (**504954**) can be autoclaved repeatedly at 121°C/15PSI, 15 min. cycle
- Polypropylene fittings are chemically inert and resistant to most organic and inorganic solvents
- Nylon fittings are strong and can be bonded with adhesive.

## Applications

- Liquid flow experimental setups



504954



504955

## ORDERING INFORMATION

**504954** Luer-to-Tubing Coupler Assortment Kit (Polypropylene)  
**504955** Luer-to-Tubing Coupler Assortment Kit (Nylon)

Assemble quick-disconnect luer fittings for use with flexible tubing with internal diameters of 1/16", 3/32" and 1/8". A variety of quick-disconnect connectors can be quickly made for connecting small diameter flexible tubing; 3-way connections can be made with the use of the 3-way luer tee; luer plugs, tees, connectors, bulk-head mounts, color coding rings, locking nuts, male and female luers—are all included to enhance the versatility of this kit. The kit has over 250 assorted parts and is offered in two different types of materials. Nylon parts are not autoclavable.

# Luer Valve Assortment Kit 14011 Parts

*Kit parts are also available individually*



**14034-40**  
Injection Site  
Male Luer lock  
Pack of 40



**14039-10**  
Check Valve  
Pack of 10



**14044-5**  
Syringe  
Activated Dual  
Check Valve  
Pack of 5



**14045-20**  
Syringe  
Slip Luer Valve  
Activated Check  
Pack of 20



**13822-10**  
0.135"/"3.4 mm  
OD Tubing  
Pack of 10



**14041-60**  
Roller Clamp  
3/16" Tubing  
Pack of 60



**7465-20**  
Pinch Clamp  
Large Bore  
Pack of 20



**14040-50**  
Pinch Clamp for  
7mm Tubing  
Pack of 50



**3742-20**  
Female T Luer  
Pack of 20



**14047-10**  
4-Port Infusion Y  
Swivel Thread  
Pack of 10



**14048-20**  
3-Port Infusion Y  
Swivel Thread  
Pack of 20



**14057-10**  
4-Way Stopcock,  
Luer Lock  
Pack of 10



**14036-15**  
4-Way Luer  
Stopcock  
Pack of 15



**14058-10**  
4-Way Stopcock,  
Luer Lock  
Pack of 10



**14035-10**  
4-Way Stopcock,  
Luer lock,  
Pack of 10



**14051-100**  
Pinch Clamp for 5  
mm Tubing  
Pack of 100



**14038-10**  
1-Way Stopcock  
Luer Lock,  
Pack of 10



**14054-10**  
1-Way Stopcock,  
Luer Slip  
Pack of 10



**14055-2**  
4-Port Manifold (6 Female Ports)  
Pack of 2



**14059-2**  
3-Port Manifold (5 Female Ports)  
Pack of 2



**13156-100**  
Female Luer  
Fitting for  
1/16" ID Tubing  
Pack of 100



**13157-100**  
Female Luer  
Fitting for  
3/32" ID Tubing  
Pack of 100



**13158-100**  
Female Luer Fitting  
for  
1/8" ID Tubing  
Pack of 100



**13159-100**  
Female Luer  
Fitting for  
5/32" ID Tubing  
Pack of 100



**13160-100**  
Male Luer  
Fitting for  
1/16" ID Tubing  
Pack of 100



**13161-100**  
Male Luer  
Fitting for  
3/32" ID Tubing  
Pack of 100



**13162-100**  
Male Luer  
Fitting for  
1/8" ID Tubing  
Pack of 100



**13163-100**  
Male Luer  
Fitting for  
5/32" ID Tubing  
Pack of 100



**14061-60**  
Male/Female  
Luer Plug  
Pack of 60



**14042-100**  
Slide Clamp  
for 2.5 mm  
O.D. Tubing  
Pack of 100

*Parts in kit may differ slightly in appearance from those pictured.*

**WORLD PRECISION INSTRUMENTS**

**UK:** +44 (0)1462 424700 • [WPIUK@WPI-EUROPE.COM](mailto:WPIUK@WPI-EUROPE.COM) • [WWW.WPI-EUROPE.COM](http://WWW.WPI-EUROPE.COM)  
**BRAZIL:** (013) 406-29703 • [INFO@BRAZIL.WPIINC.COM](mailto:INFO@BRAZIL.WPIINC.COM) • [WWW.WPIINC.COM](http://WWW.WPIINC.COM)

**GERMANY:** +49 (0)6031 1602171 • [WPIDE@WPI-EUROPE.COM](mailto:WPIDE@WPI-EUROPE.COM) • [WWW.WPI-EUROPE.COM](http://WWW.WPI-EUROPE.COM)  
**CHINA:** +86 21 6888 5517 • [CHINASALES@CHINA.WPIINC.COM](mailto:CHINASALES@CHINA.WPIINC.COM) • [WWW.WPIINC.NET](http://WWW.WPIINC.NET)

# BNC Cables & Connectors

For wiring any laboratory setup



## ORDERING INFORMATION

PART #	APPLICATION/DESCRIPTION	CONNECTOR A	CONNECTOR B	CABLE LENGTH
2851	Standard BNC cable	BNC (male)	BNC (male)	6 ft (1.8 m)
5374	Low-noise cable for microelectrode holders	BNC (male)	2 mm gold pin	4 ft (1.2 m)
5375	Low-noise cable for microelectrode holders	BNC (male)	2 mm gold jack	4 ft (1.2 m)
13324	Adapter	Double-banana (female)	BNC (male)	none
13347	ISO2 (chart recorder adapter)	Double-banana (male)	BNC (female)	none
13451	Adapter: Iso-DAM, Iso-DAM8	BNC (female)	two 2 mm pins	6 in. (15 cm)
13854	BNC T-connector, male to:	BNC (female)	BNC (female)	none
14254	BNC Straight Adapter	BNC (female)	BNC (female)	none
500184	Standard BNC Cable	BNC (male)	BNC (male)	10 ft (3 m)
500256	BNC Right Angle Adapter	BNC (male)	BNC (female)	none
500257	Standard BNC Cable	BNC (male)	BNC (male)	6 in. (15 cm)
500258	Standard BNC Cable	BNC (male)	BNC (male)	12 in. (30 cm)
500259	Standard BNC Cable	BNC (male)	BNC (male)	18 in. (46 cm)
CBL102	DAM Series, PM Series	3.5 mm MiniPhone plug	BNC (male)	6 ft (1.8 m)

## Micro Cannula

Quickly connect to existing experimental plumbing

### Features

- 0.4 mm O.D., 0.2 mm I.D. tubing
- Autoclavable
- Biocompatible Perfluorocarbon tubing material

### Benefits

- May be used with a pressure transducer (BLPR2) or a micro syringe injection system (UMP3 or MMP)

### Applications

- Cannula for carotid or femoral arteries of rodents and small animal blood vessels

This micro cannula is ideal for placement in the carotid or femoral artery of mice, rats, and other small animal blood vessels. It can be used with a pressure transducer (WPI's **BLPR2**) for blood pressure measurement, or in conjunction with a micro-syringe injection system (like WPI's **UMPIII** or

**MMP** pumps). The incorporated standard female luer fitting makes connecting to existing experimental plumbing quick and easy. The cannula is provided with a contoured-tip stainless steel stylet (trocar) to facilitate placement using established techniques. A movable "shoulder" ring provides a tie-in point to prevent accidental removal. The cannula may be left in place for 2 hours or more, and with proper care and cleaning, may be re-used multiple times. Instructions for use included.

KZ1101

## ORDERING INFORMATION

**KZ1101** Micro Cannula, 3-inch

# Precious Metal and Specialty Wire

Bare and coated metal wire for most laboratory applications



Micro coaxial cables (**MAXxxxx**) are ideal for microelectrode fabrication and construction of similar research tools. The dual shielding eliminates electrical interference caused by radio frequencies (RF), electrostatic and microphonics (e.g., bending and vibration). Available with single or dual (twin) conductors.

Teflon-coated stainless steel (type 304) wire (**SSTxxxx**) is available in 25-ft and 50-ft lengths. The Teflon coating is 150 micro-in. thick (4  $\mu$ m). The Teflon coating is designed to reduce surface friction, only. It is not insulation.

Carbon wire (**C3005**) is a single 30-micron fiber of electrochemically activated carbon. This fiber is especially useful in micro-electrochemical experiments.

Platinum/iridium wire — uncoated (**PTxxxx**) and Teflon-coated (**PTTxxxx**) — is an alloy of 90% platinum and 10% iridium, giving excellent tensile strength and corrosion resistance. Uncoated pure platinum wire (**PTPxxx**) is 99.95% pure. Indium wire (IN1003) is 99.99% pure, with a melting point of 156.4°C.

Annealed silver wire (**AGWxxxx**), 99.99% pure, is available in five diameters; three of those sizes are also available with a Teflon coating (**AGTxxxx**).

Tungsten wire (**TGWxxxx**), available in three diameters, is 99.95% pure.

Gold wire (**AUWxxxx**) is 99.99% pure. Stainless steel wire (**SSxxxxx**) is type 316.

## Wire Cutters and Scissors



504749

### ORDERING INFORMATION

Catalog No.	Metal	Coating	AWG*	Diameter	Precut Length
AGT0510	Silver	Teflon	36	0.005 in. (0.125 mm) <sup>1</sup>	10 ft (3 m)
AGT0525	Silver	Teflon	36	0.005 in. (0.125 mm) <sup>1</sup>	25 ft (7.6 m)
AGT05100	Silver	Teflon	36	0.005 in. (0.125 mm) <sup>1</sup>	100 ft (30 m)
AGT1010	Silver	Teflon	30	0.010 in. (0.25 mm) <sup>1</sup>	10 ft (3 m)
AGT1025	Silver	Teflon	30	0.010 in. (0.25 mm) <sup>1</sup>	25 ft (7.6 m)
AGT10100	Silver	Teflon	30	0.010 in. (0.25 mm) <sup>1</sup>	100 ft (30 m)
AGT1510	Silver	Teflon	26-27	0.015 in. (0.38 mm) <sup>1</sup>	10 ft (3 m)
AGT1530	Silver	Teflon	26-27	0.015 in. (0.38 mm) <sup>1</sup>	30 ft (9.1 m)
AGW0510	Silver	—	36	0.005 in. (0.125 mm)	10 ft (3 m)
AGW0530	Silver	—	36	0.005 in. (0.125 mm)	30 ft (9.1 m)
AGW1010	Silver	—	30	0.010 in. (0.25 mm)	10 ft (3 m)
AGW1030	Silver	—	30	0.010 in. (0.25 mm)	30 ft (9.1 m)
AGW1510	Silver	—	26-27	0.015 in. (0.38 mm)	10 ft (3 m)
AGW1530	Silver	—	26-27	0.015 in. (0.38 mm)	30 ft (9.1 m)
AGW2010	Silver	—	24	0.020 in. (0.5 mm)	10 ft (3 m)
AGW2030	Silver	—	24	0.020 in. (0.5 mm)	30 ft (9.1 m)
AGW4010	Silver	—	18	0.040 in. (1.0 mm)	10 ft (3 m)
AUW0170	Gold	—	50	0.001 in. (0.025 mm)	70 ft (21 m)
AUW201	Gold	—	24	0.020 in. (0.5 mm)	1 ft (30 cm)
C3005	Carbon	—	49	0.0012 in. (30 $\mu$ m)	5 ft (1.5 m)
PT1002	Platinum / Iridium	—	30	0.010 in. (0.25 mm)	2 ft (61 cm)
PT0402	Platinum / Iridium	—	38	0.004 in. (0.102 mm)	2 ft (61 cm)
PT0203	Platinum / Iridium	—	44	0.002 in. (0.051 mm)	3 ft (91 cm)
PT0110	Platinum / Iridium	—	50	0.001 in. (0.025 mm)	10 ft (3 m)
PTP101	Platinum	—	30	0.010 in. (0.25 mm)	1 ft (30 cm)
PTP201	Platinum	—	24	0.020 in. (0.5 mm)	1 ft (30 cm)
PTP401	Platinum	—	18	0.039 in. (1.0 mm)	1 ft (30 cm)
PTP406	Platinum	—	18	0.039 in. (1.0 mm)	0.5 ft (15.2 cm)
PTT0502	Platinum / Iridium	Teflon	36	0.005 in. (0.125 mm) <sup>1</sup>	2 ft (61 cm)
PTT0203	Platinum / Iridium	Teflon	44	0.002 in. (0.051 mm) <sup>1</sup>	3 ft (91 cm)
PTT0110	Platinum / Iridium	Teflon	50	0.001 in. (0.025 mm) <sup>1</sup>	10 ft (3 m)
SS31605	Stainless Steel	—	36	0.005 in. (0.125 mm)	50 ft (15.2 m)
SS31614	Stainless Steel	—	27	0.014 in. (0.36 mm)	30 ft (9.1 m)
SST30407-25	Stainless Steel	Teflon	33	0.007 in. (0.18 mm) <sup>3</sup>	25 ft (7.6 m)
SST30407-50	Stainless Steel	Teflon	33	0.007 in. (0.18 mm) <sup>3</sup>	50 ft (15.2 m)
TGW0325	Tungsten	—	40	0.003 in. (0.075 mm)	25 ft (7.6 m)
TGW0515	Tungsten	—	36	0.005 in. (0.125 mm)	15 ft (4.6 m)
TGW1510	Tungsten	—	26-27	0.015 in. (0.38 mm)	10 ft (3 m)

### MICROCOAXIAL CABLES

MAX3820	Tinned Cu Alloy	Coaxial		0.0173 in. (0.44 mm)	20 ft (6 m) <sup>4</sup>
MAX4020	Tinned Cu Alloy	Twin Coaxial		0.0158x0.024 in. (0.4x0.61 mm)	20 ft (6 m) <sup>5</sup>

\*Brown & Sharpe

<sup>1</sup> Plus 0.002 in. for Teflon coating

<sup>3</sup> Teflon adds 0.00015 in. (4  $\mu$ m) to diameter

<sup>4</sup> Impedance: 50 ohm; capacitance: 95 pF/m; resistance: 5 ohm/m

<sup>5</sup> Impedance: 100 ohm; capacitance: 54 pF/m; resistance: 1.9 ohm/m

### ORDERING INFORMATION

504749	Ergonomic Micro-Shear Flush Cutters, 12.7 cm (5 in.) Micro-shear flush cutters for delicate wires
504750	Ergonomic Mini-Scissors, 12.7 cm (5 in.) For cutting fine or delicate items with a clean, square edge. Handle design is advantageous for users with arthritic hands.
504751	Ergonomic Micro-Shear Flush Cutters, 12.7 cm (5 in.) For delicate wires. ESD safe. Extra tough high carbon steel blades. Sized for smaller hands and maximum maneuverability

### WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • WPIUK@WPI-EUROPE.COM • WWW.WPI-EUROPE.COM  
BRAZIL: (013) 406-29703 • INFO@BRAZIL.WPIINC.COM • WWW.WPIINC.COM



GERMANY: +49 (0)6031 1602171 • WPIDE@WPI-EUROPE.COM • WWW.WPI-EUROPE.COM  
CHINA: +86 21 6888 5517 • CHINASALES@CHINA.WPIINC.COM • WWW.WPIINC.NET



# FrameWorks

## Non-Magnetic Bases, Stainless Steel Rods, & Clamps

These high quality components are made of stainless steel and polymer that resist organic solvents and corrosion. They can be easily assembled to make a stand-alone setup for student labs or to make a complicated frame for research labs.




### ORDERING INFORMATION

<b>503041</b>	Large Clamp with Rod (157 mm), opens up to 85 mm
<b>503042</b>	Medium Clamp with Rod (157 mm), opens up to 45 mm
<b>503086</b>	Small Clamp with Rod (157 mm), opens up to 16 mm
<b>14073-4</b>	Open-sided Frame Clamp (pkg. of 4)
<b>503082-4</b>	Board Frame Clamp, opens to 8.5 mm (pkg. of 4)
<b>503080-4</b>	Frame Clamp with Parallel Surface Mount (includes mounting screws) (pkg. of 4)
<b>502193-4</b>	Parallel Frame Clamp (pkg. of 4)
<b>503078-4</b>	T-joint Frame Clamp (pkg. of 4)
<b>503079-4</b>	In-line Frame Clamp (pkg. of 4)
<b>502190</b>	Heavy Rectangular Base (with M8 Thread mount and thumbscrew mount), 23×15.6 cm, 4 lb.
<b>503083</b>	Light Rectangular Base (with M8 Thread mount and thumbscrew mount), 23×15.6 cm, 0.5 lb.
<b>503085</b>	Large 10-in. V-base with M8 Thread Mount
<b>503084</b>	Small 6-in. V-base with M8 Thread Mount
<b>503081-4</b>	Vertical Surface Mount, M8 Threaded (pkg. of 4)
<b>503070</b>	Polished Stainless Steel Post, 12 mm OD, 25 cm long, No Thread
<b>503071</b>	Polished Stainless Steel Post, 12 mm OD, 50 cm long, No Thread
<b>503072</b>	Polished Stainless Steel Post, 12 mm OD, 75 cm long, No Thread
<b>503073</b>	Polished Stainless Steel Post, 12 mm OD, 25 cm long, M8 Thread
<b>502191</b>	Polished Stainless Steel Post, 12 mm OD, 50 cm long, M8 Thread
<b>503075</b>	Polished Stainless Steel Post, 12 mm OD, 60 cm long, M8 Thread
<b>503076</b>	Polished Stainless Steel Post, 12 mm OD, 75 cm long, M8 Thread
<b>503077</b>	Polished Stainless Steel Post, 12 mm OD, 80 cm long, M8 Thread

# Digital Caliper

For easy and accurate measurements



The high quality electronic digital caliper is a useful tool — no laboratory should be without one because it is more accurate and easier to use than the traditional analog devices. Measure in either inches or millimeters at the touch of a button. The floating zero feature allows

you to read the increment without calculation. The caliper is made from hardened stainless. The caliper measures up to 150 mm (6 in.) with 0.01 mm (or 0.0005") resolution. The caliper is designed to be water-resistant to IP54 as defined in the IEC529 standard. The IP54 code's first digit "5" means dust-protected. The second digit "4" means that the caliper is protected against splashing water — it can withstand a shower from every direction for 10 minutes.

## ORDERING INFORMATION

<b>501601</b>	Digital Caliper
<b>502157</b>	Replacement Battery (package of 10)

# Non-rotating Spindle Digital Micrometer Head

Build your own precision micro-positioning device



The new non-rotating spindle digital micrometer head allows you to create your own micro-positioning instrument. With micron-level accuracy, it gives higher precision than a normal micromanipulator. Since the spindle does not rotate as it advances, instruments can be directly attached without the need for a complicated decoupling device. The digital display eliminates the need to squint at the notational scale. Readings can be clearly seen in either inches or millimeters. You can read both absolute position and the increment relative to a previously chosen point.

## 502102 SPECIFICATIONS

TOTAL TRAVEL DISTANCE	25 mm
RESOLUTION	0.001 mm
ACCURACY	± 0.003 mm
SPINDLE	Ø 8 mm
MOUNTING	Ø 12 mm x 10 mm
TOTAL LENGTH	166 mm
MEASUREMENT MODE	Absolute and incremental
DIGITAL READOUT	mm or inch
ANALOG READOUT	mm
DATA OUTPUT	RS232
ENVIRONMENTAL PROTECTION	IP54
SHIPPING WEIGHT	0.51 kg (1.12 lb.)

## ORDERING INFORMATION


<b>502102</b>	Non-Rotating Spindle Micrometer Head
---------------	--------------------------------------

# Powerful Ball Joint Rare Earth Magnet

Construct holding devices for small parts/equipment

- Small but very powerful: holds 2 kilograms (~5 pounds)!
- Steel ball rotates freely 360° on a 180° axis
- M3 mounting screw on ball for attachment to equipment
- Magnet base threaded (M3) for mounting onto a base or equipment

This novel magnetic ball joint has phenomenal holding power for up to 2kg of attached weight while permitting the ball a full 360° rotation on a 180° axis. You can freely orient your equipment to an infinite number of positions within this rotation. This is made possible by the combination of a steel ball (10 mm diameter) and a powerful rare earth magnet contained in the magnet cylinder (ø 10 x



20 mm). Convenient M3 attachment sites are provided on both the ball (male) and the magnet base (female). For use with micromanipulators for the positioning and holding of optical instruments including various lighting sources and lasers, pipettes and any small parts that would benefit from the flexibility offered by this new magnetic ball joint.



## ORDERING INFORMATION

<b>500871</b>	Magnetic Ball Joint
---------------	---------------------

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • WPIUK@WPI-EUROPE.COM • WWW.WPI-EUROPE.COM  
BRAZIL: (013) 406-29703 • INFO@BRAZIL.WPIINC.COM • WWW.WPIINC.COM

GERMANY: +49 (0)6031 1602171 • WPIDE@WPI-EUROPE.COM • WWW.WPI-EUROPE.COM  
CHINA: +86 21 6888 5517 • CHINASALES@CHINA.WPIINC.COM • WWW.WPIINC.NET

# Culture Dishes with Optical Glass Bottom

Cover-glass bottom for observing & growing cells for imaging

Back by popular demand

Black Wall FluoroDish™

## Features

- Optical quality glass bottom for better imaging quality (RI=1.525)
- Low sample volume for expensive chemicals
- Lowest access angle for micropipette
- Black wall available for low background fluorescent measurement

## Benefits

- Multiple sizes and designs to suit your application
- Optional Poly-D-Lysine coating for neurons
- Dishes designed for low volumes or large growth areas

## Applications

- High resolution image analysis
- Microinjection
- Electrophysical recording of fluorescent-tagged cells
- Black wall available for use with confocal microscopes



WPI's FluoroDish™ tissue culture dishes provide exceptional imaging quality for many applications requiring the use of inverted microscopes such as high resolution image analysis, microinjection and electrophysical recording of fluorescent-tagged cells. Taking advantage of WPI's extensive experience with low toxicity adhesives, FluoroDish™ uses a specially formulated adhesive that is optically clear, durable and with extremely low toxicity. Tests by an independent laboratory have shown that the 96-hour surviving rate of embryos is 100% when kept in FluoroDish™, substantially better than other brands. The bottom glass has superior UV transmission (30% transmission at 300 nm, compared to less than 7% for the most popular German glass). Stringent quality control ensures that glass thickness is 0.17 ±0.01 mm.

Conventional plastic dishes and chambers limit the use of the inverted scope, because the thick plastic bottom requires a long working distance objective available only in lower magnifications. Each WPI dish has a flat (0.17mm thick), optical quality glass bottom, allowing the use of a much shorter working distance, larger numerical aperture (NA) and higher magnification (up to 100X). The larger NA and higher magnification provide superior quality imaging for both classical and fluorescence microscopy. Higher effective NA yields brighter images for fluorescence and higher resolution in image analysis. The glass bottom permits the use of immersion objectives with media such as water, glycerin or oil for the highest magnification possible. To optimize heat-exchange, WPI's glass-bottom dish is designed to be flush (flat) with the microscope stage or heating unit, eliminating the air gap that exists with modified plastic dishes where a cover slip was inserted.

## Multiple sizes

We have a 50 mm diameter dish and two types of 35 mm diameter dishes. An inner well is created within the dish by the glass bottom and the tissue culture grade polystyrene which forms the sides of the dish. They are individually packed and gamma sterilized.

The 35 mm dish has outside dimensions similar to that of a Corning 35 mm dish and has ø23.5 mm glass window (FD35) or ø10 mm glass window (FD3510). Most heaters and perfusion adapters designed for the Corning 35 mm dish also fit this dish.

## Poly-D-Lysine coating or uncoated

The 23.5 mm glass dish is available uncoated or poly-D-lysine-coated. Some cell lines (e.g., PC3 and HEK) adhere well to the uncoated glass bottom dish. The poly-D-lysine coating reportedly improves the adhesion of neuron cells. You may also apply any special coating that is best for your cell line to uncoated dishes.

## Low volume dishes or large growth areas

The 10 mm glass dish (FD3510) has low sidewalls for easy microelectrode access and low solution volume. The low microelectrode access angle is the lowest among all of 35 mm dishes on the market (very close to that of a 50 mm dish). The dish needs only about 115 µL to cover the bottom well, an important feature when using expensive drugs and chemicals.

The 50 mm dish (FD50) has a large growth area (35 mm well diameter), a low access angle for microelectrodes, and grips for easy handling.

## References

Petros Batsios, Xiang Ren, Otto Baumann, Denis A. Larochelle, Ralph Gräf, *Cells*. 2016 Mar; 5(1): 13. Published online 2016 Mar 18. doi: 10.3390/cells5010013 PMID: PMC4810098


P. Calcraft, A. Arredouani, M. Ruas, Z. Pan, X. Cheng, X. Hao, J. Tang, K. Rietdorf, L. Teboul, K. Chuang, P. Lin, R. Xiao, C. Wang, Y. Zhu, Y. Lin, C. Wyatt, J. Parrington, J. Ma, A. M. Evans, A. Galione, M. Zhu. *Nature*. 2009 May 28; 459(7246): 596–600. Published online 2009 April 22. doi: 10.1038/nature08030 PMID: PMC2761823

T Mitra-Ganguli, I. Vitko, E. Perez-Reyes, A. Rittenhouse. *J Gen Physiol*. 2009 November; 134(5): 385–396. doi: 10.1085/jgp.200910204 PMID: PMC2768804 ArticlePubReaderPDF-3.0 m

S. Khuon, L. Liang, R. Dettman, P. Sporn, R. Wysolmerski, T. Chew. *J Cell Sci*. 2010 February 1; 123(3): 431–440. Published online 2010 January 12. doi: 10.1242/jcs.053793 PMID: PMC2816188

## SPECIFICATIONS

Style	ID (mm)	OD (mm)	Glass Ø (mm)	height (inside)	Height (outside)	Access Angle
FD35	33	35.5	23.5	7.8	9	29°
FD5040	47.5	49.82	35	7.25	7.4	17°



## ORDERING INFORMATION

FD35-100	FluoroDish™ Sterile Culture Dish, clear wall, 35 mm, 23 mm well, box of 100
FD35B-100	FluoroDish™ Sterile Culture Dish, black wall, 35 mm, 23 mm well, box of 100
FD35PDL-100	FluoroDish™ Sterile Culture Dish, Poly-D-Lysine Coated, clear wall, 35 mm, 23 mm well, box of 100
FD3510-100	FluoroDish™ Sterile Culture Dish, clear wall, 35 mm, 10 mm well, low sidewall, box of 100
FD3510B-100	FluoroDish™ Sterile Culture Dish, black wall, 35 mm, 10 mm well, low sidewall, box of 100
FD5040-100	FluoroDish™ Sterile Culture Dish, clear wall, 50 mm, 35 mm well, box of 100
FD5040B-100	FluoroDish™ Sterile Culture Dish, black wall, 50 mm, 35 mm well, box of 100

# Product Index

Item	Page	Item	Page	Item	Page	Item	Page	Item	Page
1B100-3	42, 47, 104, 217	3578	97, 100, 103	13620	97	15920-G	195	500199	71
1B100-4	47, 104, 217	3579-20	9	13685	63, 97	15921	195	500200	71
1B100-6	47, 104, 217	3669	9	13740	211	15921-G	195	500222	195
1B100F-3	47, 104, 217	3670	97	13776	97	15975	97	500222-G	195
1B100F-4	47, 104, 217	3718	79, 88	13822-10	220	15976	97	500224	202
1B100F-6	47, 104, 217	3742-20	220	13834	48, 108, 110, 111	40239	59	500224-G	202
1B120	47, 104, 217	3955	9	13835	110, 111	53330-01	5	500227	202, 207
1B120-4	47, 104, 217	3960	9	13836	110, 111	56200	170	500227-G	202
1B120-6	47, 104, 217	3993	4	13854	97, 221	58006	167	500228	197
1B120F-3	47, 104, 217	4886	139, 213	13962	63	58450	167	500228-G	197
1B120F-4	47, 104, 217	4898	139, 142, 213, 216	14003	197, 207, 208	61800	127	500231	194
1B120F-6	47, 104, 217	5052	32	14003-G	197	61824	127	500233	194
1B150-3	47, 104, 217	5153	9	14011	72, 219	61830	127	500256	97, 221
1B150-4	47, 104, 217	5210	9	14034-40	220	61840	124, 127	500257	97, 221
1B150-6	47, 104, 217	5233	9	14035-10	220	74016	79	500258	97, 221
1B150F-3	47, 104, 217	5332	115	14036-15	130, 220	74020	79	500259	97, 221
1B150F-4	47, 104, 217	5361	9	14038-10	220	74030	79	500261	181
1B150F-6	47, 104, 217	5362	9	14039-10	220	74040	79	500262	181
1B200-4	47, 104, 217	5371	97, 100	14040-50	220	74050	79	500264	181
1B200-6	47, 104, 217	5373	97	14041-60	220	75040	113, 114	500266	181
1B200F-4	47, 104, 217	5374	97, 221	14042-100	220	75050	113, 114	500292	113
1B200F-6	47, 104, 217	5375	97, 221	14044-5	220	75122-210	45	500329	114
2B150F-4	105	5377	145, 149	14045-20	220	75122-310	45	500339	194
2B150F-6	105	5378	145, 149	14047-10	220	75122-410	45	500340	194
2 mm Plug	97	5385	145, 147	14048-20	220	77020	106	500342	207, 208
3-Axis Micropositioner	26	5389	37	14051-100	220	89337	173	500343	205
3-Axis Stereotaxic Frame	118	5399	147	14054-10	220	89341	173	500344	205
3B120F-4	105	5430-10	45	14055-2	220	89342	173	500345	205
3B120F-6	105	5430-12	45	14057-10	220	89575	170, 173	500346	205
4-Axis Motorized Micromanip	28	5430-15	45	14058-10	220	91580	147, 148, 149, 151	500348	205
5B120F-4	105	5430-20	45	14059-2	220	91736	4	500349	205
5B120F-6	105	5430-ALL	44, 45	14061-60	220	91750	4	500350	205
5 lb Weight	25	5435	145, 147	14073-4	223	91799	4	500351	205
6-Axis Stereotaxic Frame	118	5436	145, 147	14074	110, 111	94650	179, 21	500352	205
7B100F-4	105	5440	105	14091	207	94689	179, 21	500353	205
7B100F-6	105	5444	101	14109	208	95611	145	500354	205
7B120F-4	105	5447	75	14110	208	95620	145	500355	205
8-pin DIN Plug	79, 88	5464	25	14120	205	95699	145	500366	200
705	83	5468	88, 91	14121	205	96824	37	500366-G	200
705PF	83	5469	74, 75, 91, 92	14124	208	97204	11, 17	500367	200
711P	88	5470	74, 91, 92	14125	198	98727	37	500367-G	200
712P	81	5479	32	14125-G	198	98728	37	500453	208
715P	81	5482	91	14126	198	98991	42	500456	207
773	81	5483	91	14130	208	99259	179, 21	500475	25, 29, 33
900A	115	5489	75	14131	207	99261	179, 21	500476	25, 29, 33
900AP	115	6820	139, 142, 213, 216	14192	207	100042	145, 148	500778	59
900APP100	115	7128	139, 142, 213, 216	14218	201, 207, 208	100084	145	500828	187
900A-VAC	115	7325	145, 147	14218-G	201	300033	59	500871	30, 224
1571	101, 139, 142, 213, 216	7326	145, 149	14219	201	300040	97	501223	200
1965	107	7335	139, 142, 213, 216	14219-G	201	300051	171	501223-G	200
1966	107	7341	139, 141, 213, 215	14224	201	300052	171	501224	200
1967	107	7342	139, 141, 213, 215	14224-G	201	300102	75, 91	501224-G	200
2026-10	97	7357	147	14225	201	300411	37	501225	200
2033	74, 75, 78	7465-20	220	14225-G	201	300521	59	501225-G	200
2034	74, 75, 78	7521	147	14226	195	300523	4	501226	200
2035	74, 78	13024	88, 130	14226-G	195	300626	35, 37	501226-G	200
2505	101	13025	88, 130	14240	208	300647	75	501227	200
2541	83	13142	55, 59, 113, 114	14254	97, 221	500023	202	501227-G	200
2547	80, 81, 82	13156-100	220	14292	207	500023-G	202	501228	200
2851	15, 17, 156, 97, 151, 23	13157-100	220	14393	208	500028	181	501228-G	200
2932	45	13158-100	220	14396	207	500064	194	501236	197
2933	45, 81, 84, 115,	13159-100	220	15590	37	500085	194	501236-G	197
3072	75	13160-100	220	15623	63, 97	500092	207	501237	197
3161	79, 88, 97	13161-100	220	15790	81	500162	185	501238	197
3259	81, 84	13162-100	220	15804	83	500184	79, 88, 97, 130, 221	501263	201
3260	45	13163-100	220	15873	25, 29	500186	181, 191	501263-G	201
3294	74, 97, 100	13316	105, 139, 142, 213, 216	15911	205	500191	71	501264	201
3316	45	13324	97, 221	15913	205	500192	71	501264-G	201
3417-10	97	13338	183, 191	15915	207, 208	500193	71	501274	209
3484	75	13347	97, 100, 221	15916	195	500194	71	501284	209
3485	9, 75	13370	171	15916-G	195	500195	71	501291	195
3491	97, 100, 115, 130	13388	74, 75, 97, 100	15917	195	500196	71	501291-G	195
3492	97	13395	171	15917-G	195	500197	71	501312	202
3517	97, 100	13451	97, 221	15920	195	500198	71	501313	207, 208

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com Germany: +49 (0)6031 1602171 • wpide@wpi-europe.com • www.wpi-europe.com  
 Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • www.wpiinc.com China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.net

Item	Page	Item	Page	Item	Page	Item	Page	Item	Page
501319.....	207, 208	502043.....	126	503070.....	223	504444.....	187	505062.....	9
501322.....	207	502046.....	126	503071.....	223	504445.....	187	505063.....	9
501352.....	183	502047.....	126	503072.....	223	504459.....	209	505066.....	175
501353.....	183	502053.....	126	503073.....	223	504487.....	203	505067.....	175
501369.....	183	502054.....	126	503075.....	223	504500.....	203	505068.....	175
501370.....	183	502055.....	124	503076.....	223	504505.....	193	505069.....	27, 28, 175
501371.....	183	502056.....	124	503077.....	223	504506.....	193	505070.....	28
501372.....	183	502060.....	123	503078-4.....	223	504507.....	193	505071.....	27, 28
501373.....	183	502063.....	123	503079-4.....	223	504508.....	193	505073.....	27, 28
501375.....	183	502063-MGM.....	125	503080-4.....	223	504509.....	193	505074.....	27
501376.....	183	502063-SGM.....	125	503081-4.....	223	504515.....	193	505085.....	27, 28
501377.....	183	502067.....	124	503082-4.....	223	504516.....	193	505094.....	110, 111
501378.....	183	502068.....	124	503083.....	223	504524.....	193	505095.....	110, 111
501379.....	183	502070.....	124	503084.....	223	504528.....	204	505096.....	110, 111
501381.....	183	502102.....	224	503085.....	223	504529.....	204	505097.....	110, 111
501601.....	224	502105.....	25	503086.....	223	504530.....	204	505098.....	110, 111
501608.....	132	502109.....	51	503090.....	155	504531.....	204	505099.....	110, 111
501609.....	167, 171	502110.....	51	503097.....	190	504570.....	189	505101.....	110, 111
501622.....	38	502124.....	157	503099.....	190	504576.....	61, 62	505103.....	110, 111
501623.....	38	502125.....	51	503102.....	183	504577.....	61, 62	505105.....	110, 111
501624.....	38	502157.....	224	503110.....	157	504578.....	61, 62	505106.....	110, 111
501636.....	185	502163.....	183	503120.....	52, 53	504591.....	70	505108.....	110, 111
501637.....	185	502190.....	223	503121.....	53	504596.....	183	505109.....	110, 111
501641.....	155	502191.....	223	503122.....	53	504597.....	183	505131.....	181
501642.....	155	502193-4.....	223	503207.....	55	504606.....	187	505132.....	181
501644.....	155, 157	502201.....	55, 124	503259.....	201	504608.....	124	506122.....	110, 111
501651.....	32	502210.....	124	503260.....	201	504638.....	204	600011.....	145, 148
501652.....	32	502213.....	123, 124, 126	503294.....	208	504639.....	204	600012.....	145, 148
501653.....	32	502224.....	124	503301.....	55	504640.....	204	600015.....	145
501656.....	155	502225.....	124	503305.....	197	504641.....	204	600016.....	145
501656-C.....	155	502226.....	123	503306.....	197	504642.....	204	600022.....	140, 214
501656-F.....	155	502227.....	122	503364.....	198	504643.....	204	800100.....	160
501657.....	155	502228.....	122	503365.....	198	504644.....	204	800120.....	177
501658.....	155	502229.....	122	503391.....	196	504646.....	204	800292.....	113
501658-F.....	155	502230.....	122	503506.....	218	504647.....	204	801566.....	11, 41, 71, 115
501705.....	195	502231.....	122	503507.....	218	504648.....	204	801963.....	11, 41, 71, 115
501705-G.....	195	502232.....	122	503508.....	218	504649.....	204	802407.....	179, 21
501725.....	202	502233.....	122	503510.....	186	504650.....	204	µManager.....	18
501725-G.....	202	502234.....	122	503513.....	113, 187	504651.....	204	µTip.....	107, 218
501764.....	194	502235.....	124	503520.....	186	504652.....	204		
501778.....	197	502236.....	124	503536.....	97, 189	504653.....	204		
501839.....	197	502237.....	124, 209	503538.....	157	504654.....	204	<b>A</b>	
501842.....	209	502241.....	123	503568.....	32	504655.....	204	A310 Accupulser.....	84
501851.....	209	502242.....	124	503569.....	32	504713.....	97	A362.....	87, 88
501852.....	209	502243.....	124	503570.....	32	504714.....	97	A365D.....	85
501853.....	209	502244.....	124	503571.....	32	504715.....	97	A365R.....	85
501854.....	209	502245.....	124	503573.....	127	504716.....	97	A365RC.....	85
501855.....	209	502249.....	124	503573.....	127	504717.....	175	A385R.....	86
501856.....	209	502250.....	124	503598.....	124, 209	504717.....	175	A385RC.....	86
501857.....	209	502256.....	124	503599.....	124, 209	504718.....	175	A395D.....	87
501858.....	209	502257.....	124	503741.....	202	504718.....	175	A395R.....	87
501860.....	209	502259.....	123	503763.....	139, 141, 213, 215	504719.....	175	A395RC.....	87
501861.....	209	502262.....	126	503847.....	176	504735.....	204	Abrading Tip.....	209
501862.....	209	502265.....	123	503848.....	176	504736.....	204	Absorbance Detection.....	177
501863.....	209	502266.....	123	504011.....	53	504737.....	204	Accessory Stand.....	209
501864.....	209	502266.....	123	504024.....	198	504738.....	204	Accupulser Signal Generator.....	84
501986.....	141, 215	502300.....	119	504025.....	198	504739.....	204	Active Headstage.....	79
502000.....	181	502303.....	119	504026.....	198	504740.....	208	Adagio Software.....	61
502001.....	181	502350.....	119	504027.....	198	504741.....	208	Adapter.....	74, 78, 221
502002.....	181	502353.....	119	504123.....	181, 183	504744.....	193	Adapter Cable.....	97, 103
502004.....	181, 183	502360.....	119	504127.....	59	504745.....	193	Adapter, Spinal.....	123
502005.....	181, 183	502361.....	119	504128.....	183	504749.....	102, 222	Adhesive.....	213, 214
502006.....	181, 183	502600.....	120	504129.....	183	504750.....	102, 222	Adhesives.....	139, 140
502007.....	183	502603.....	120	504134.....	191	504751.....	102, 222	Adjustable Stage Platform.....	124
502009.....	181, 183	502610.....	121	504136.....	190	504766.....	211	ADPT2 Foot Switch.....	65
502010.....	181	502612.....	120	504216.....	211	504767.....	211	Adson Forceps.....	195
502011.....	181	502614.....	120	504221.....	187	504768.....	211	Ag/AgCl Electrode.....	103
502012.....	181	502616.....	120	504278.....	134	504926.....	121	Ag/AgCl Half-Cells.....	103
502013.....	181	502650.....	120	504279.....	134	504927.....	121	Ag/AgCl Pellet.....	74, 78
502015.....	181, 191	502653.....	120	504280.....	134	504929.....	181, 183	AGT0525.....	102, 222
502016.....	181	503022.....	52	504281.....	134	504930.....	191	AGT1010.....	102, 222
502017.....	181	503023.....	52	504284.....	185	504931.....	191	AGT1025.....	102, 222
502018.....	181	503041.....	223	504285.....	185	504949.....	59	AGT1510.....	102, 222
502019.....	181	503042.....	223	504286.....	185	504950.....	59	AGT1530.....	102, 222
502040.....	218	503049.....	52	504287.....	185	504951.....	48, 108	AGT05100.....	102, 222
502041.....	218	503050.....	52	504416.....	187	504954.....	72, 219	AGT10100.....	102, 222
502042.....	126	503051.....	183	504417.....	187	504955.....	72, 219	AGTxxxx.....	102, 222
		503067.....	130	504443.....	187	505060.....	9	AGW0530.....	102, 222
						505061.....	9	AGW1010.....	102, 222

Item	Page	Item	Page	Item	Page	Item	Page	Item	Page	
AGW1030.....	102, 222	Bipolar Stimulation.....	90	CHM1.....	9	CUV2108-1.....	173	CX-KIMAKIT-PC.....	66	
AGW1510.....	102, 222	Black Epoxy Coated Tweezers ...	203	CHM2.....	9	CUV2614-1.....	173	CX-KIMA-P.....	66	
AGW1530.....	102, 222	Black Instruments.....	203	CHM3.....	9	CUV2674-1.....	173	CX-MIRUS.....	67	
AGW2010.....	102, 222	Blood Pressure Monitor.....	130	CHM4.....	9	Cuvettes.....	172	CX-MIRUS-MULTI8.....	67	
AGW2030.....	102, 222	Blood Pressure, Rodent.....	131	CHM5.....	9	CUVxxx-1.....	173	CX-MIRUS-PRO.....	67	
AGW4010.....	102, 222	Blood Pressure System.....	131	CHM6.....	9	CX-001.....	69	CX-UNIGO-80.....	67	
AGWxxxx.....	102, 222	Blood Pressure Transducer.....	130	CHM7.....	9	CX-002.....	69	CX-UNIGO-EC-80.....	67	
AirTherm-SAT.....	35	BLPR2.....	130	CHM8.....	9	CX-003.....	69	CX-UNIGO-ECI-80.....	67	
AIRTHERM-SAT-1W.....	36	BLPR2 Transducer.....	130	Chopstick Electrode Set.....		4	CX-004.....	69	CX-UNIGO-ECP-80.....	67
AIRTHERM-SAT-2W.....	36, 37	BNC Cable.....	15, 17, 97, 221	Circulating Bath.....		9	CX-4UECIPAD-4X80FS.....	68	CX-UNIGO-I-80.....	67
AirTherm Satellite.....	37	BNC Right Angle Adapter ....	97, 221	Circulation Reservoirs.....		8, 9	CX-4UECPC-4X80FS.....	68	CX-UNIGO-P-80.....	67
AirTherm Satellite Heater.....	36	BNC Straight Adapter.....	221	Clamps.....		33, 223	CX-4UNCIPAD-4X80FS.....	68	Cyanoacrylate Adhesive ....	141, 215
AirTherm-SMT.....	35, 37	BNC T-Connector.....	97, 221	Clamp with Rod.....		223	CX-4UNCPC-4X80FS.....	68		
AIRTHERM-SMT-2W.....	37	BNC to BNC Cable ....	23, 75, 79, 88, 130, 156	Classic System, Anesthesia.....		128	CX-005.....	69	<b>D</b>	
AL-300.....	65	BNC-to-Double Banana Adapter .	85	Clip Applier.....		205	CX-007.....	69	D2H-DB.....	176
AL-1000.....	65	Board Frame Clamp.....	223	Clips and Clamps.....		205	CX-008.....	69	D2H-HB.....	176
AL-1000HP.....	65	Boom Clamp Stand.....	181	Closed Chamber.....		160	CX-009.....	69	D2H-HBER.....	176
AL-1000HP2.....	65	Boom Stand.....	181, 183	Closed Loop.....		37	CX-010.....	69	D4H.....	176
AL-2000.....	65	Borosilicate ..	43, 47, 104, 105, 107, 217, 218	CMF20G.....		214	CX-018.....	68	DAM50.....	75, 76
AL-4000.....	65	Box Filament.....	110, 111	CMF22G.....		214	CX-023.....	68	DAM80.....	74, 76
AL-6000.....	65	BP1.....	130	CMF23G.....		214	CX-024.....	68	DAM80P.....	75
AL-6000H.....	65	BPCABLE2.....	130	CMF26G.....		214	CX-026.....	69	DAM Series Bioamplifiers.....	74
AL-8000.....	65	Brain Matrices.....	206	CMF28G.....		214	CX-027.....	68	Data Acquisition.....	11, 23
Aladdin Syringe Pump.....	65	BRIDGE8.....	76, 79	CMF31G.....		214	CX-028.....	69	DC Amplifier System.....	79
Allodynia.....	135	Buffer Solutions.....	158	CMF34G.....		214	CX-029.....	69	DCAP.....	206
ALM Retractor.....	208	Butyl Cyanoacrylate.....	139	CMF35G.....		214	CX-030.....	69	DCAP-L.....	206
Amplifier.....	76			CMF90U.....		214	CX-031.....	69	DCAP-M.....	206
Amplifier System.....	13			CMOS Camera with OptiSarc.....		18	CX-032.....	68	Deuterium Halogen Light Source	176
Analgesia Meter.....	137, 138			C-Mount Adapter.....		187	CX-033.....	69	Deuterium Lamp.....	176
Anesthesia.....	128	<b>C</b>		C-Mount CCD Camera Coupler .		183	CX-034.....	69	Deuterium Replacement Lamp.....	176
Anesthesia Accessories.....	129	C3005.....	102, 222	C-Mount Eyepiece Adapters.....		190	CX-035.....	69	Diamond Coated Micro Scissors	198
Anesthesia Facemask.....	125	Cables.....	74, 78, 96	C-Mount Microscope Attachment.....		21, 179	CX-036.....	69	Digital Caliper.....	224
Anesthesia Masks.....	126	CAL900A.....	115	C-Mount to Eyepiece Adapter...		190	CX-037.....	69	Digital Filtering.....	73
Anesthesiometer.....	135	CALBUF-1.....	158	Coaxial Cables, Micro.....		102, 222	CX-038.....	69	Digital Microforge.....	112
Angled Electrode Holder.....	25, 29	CALBUF-2.....	158	COLCAM-HD1080P.....		190	CX-039.....	69	Digital Micrometer Head.....	224
Animal Physiology.....	117	Calcium Calibration Solutions...	158	Collimator.....		171	CX-040.....	69	Digital Microscope Camera.....	187, 189
Animal Temperature Controller	127	Calcium Detection.....	19	Colored Dissolved Organic Matter ...		168	CX-041.....	69	Digital Mouse Stereotaxic Frame	121
Anti-Oscillation Module.....	13, 17	Calcium Electrode.....	159	Compact Magnetic Stand.....		32	CX-042.....	69	Digital Stereotaxic Frame.....	119
Articulated Arm.....	181, 183	Calcium Neutral Ion Exchanger.	159	Concentric Bipolar Electrodes.....		90, 92, 93	CX-043.....	69	Dipping Probe.....	171
ATC2000.....	127	Caliper, Digital.....	224	Concentric Nerve Cuff.....		94	CX-044.....	69	DipTip.....	171
Auto-Fluorescence Imaging.....	179	Camera, Digital.....	189	Condenser.....		9	CX-045.....	68	DIP-UV-MINI-2.....	171
Automated Tissue Resistance		Camera, High Definition.....	188	Connectors.....		96	CX-046.....	69	DIP-UV-MINI-5.....	171
Measuring System.....	7	Cannula, Micro.....	130	Constant Load.....		13, 16	CX-047.....	69	DIP-UV-MINI-10.....	171
Auxiliary Knob Controller.....	42	Cannulas.....	134	Controller.....		55	CX-048.....	69	Dispensing Tips.....	213
AUW201.....	102, 222	Carbon Epoxy.....	142, 216	Cosine Corrector.....		171	CX-049.....	69	Disposable Scalpels.....	205
AUWxxxx.....	102, 222	Carbon Fiber Microelectrodes..	154	Cover Slips, Glass.....		218	CX-050.....	69	Dissecting Pad Kit.....	141, 215
Axis Adjustment Tool.....	25	Carbon Filled Conductive Epoxy	139	Crile-Wood Needle Holder.....		202	CX-051.....	69	Dissecting Scissors.....	208
		Carbon Monoxide.....	145	Cryose Fiber Assemblies.....		163	CX-052.....	69	DMF1000.....	112
		Carbon Wire.....	102, 222	CS-mount Adapter.....		185	CX-053.....	69	DMF1000-1.....	113
<b>B</b>		Cat/Monkey Adaptor.....	123	Culture Dishes.....		225	CX-054.....	69	DMF1000-2.....	113
Ball Bearing Boom Stand..	181, 183	CBL100.....	97	Current Sources.....		85	CX-055.....	69	DMF1000-H5.....	113
Ball Electrode.....	209	CBL102.....	74, 75, 78, 97, 116, 221	Cutoff Disk.....		209	CX-056.....	69	DMF1000-M1.....	113
Ball Joint.....	25, 29, 30, 33	CC-3-UV.....	171	CUV1022-10.....		173	CX-057.....	69	DMF1000-M2.....	113
Ball-Joint Holder Attachment.....	33	CD Camera Coupler.....	181	CUV2011-1.....		173	CX-058.....	69	DMP.....	46
Ball Joint, Magnetic.....	224	CDOM.....	168	CUV2012-1.....		173	CX-059.....	69	Dog/Monkey Adaptor.....	123
Ball Mill.....	209	CDOM-COAST-HV.....	169	CUV2023-1.....		173	CX-060.....	69	Double Banana Plug.....	97
Base.....	223	CDOM-COAST-LV.....	169	CUV2025-1.....		173	CX-061.....	69	Double-Holder Micromanip.....	29
Basic System, Anesthesia.....	128	CDOM-FRESH-HV.....	169	CUV2031-1.....		173	CX-062.....	69	Dovetail Extension.....	27, 28
BAT-10R/LOP.....	132	CDOM-OCEAN-HV.....	169	CUV2032-1.....		173	CX-063.....	69	Dressing Forceps.....	195
BAT-10R/LOP220.....	132	CDOM-OCEAN-LV.....	169	CUV2033-1.....		173	CX-064.....	68	Drill, Micro.....	209
BAT-12R.....	132	Cell Tester.....	12, 18	CUV2034-1.....		173	CX-065.....	68	Dri-Ref.....	158
BAT-12R-220.....	132	Cell & Tissue.....	3	CUV2051-1.....		173	CX-066.....	68, 69	DRIFEF-2.....	158
Battery.....	4, 224	Center Chamber Gasket.....	160	CUV2052-1.....		173	CX-067.....	68	DRIFEF-25H.....	158
Battery Charger.....	85, 86, 87	CF10-100.....	154	CUV2061-1.....		173	CX-068.....	68	DRIFEF-5.....	158
Beam Splitter.....	185	CF10-250.....	154	CUV2062-1.....		173	CX-069.....	68	DRIFEF-55H.....	158
Bifurcated Fibers.....	163	CF10-500.....	154	CUV2063-1.....		173	CX-070.....	66	DRIFEF-450.....	158
Bifurcated Light Guide.....	179, 181	CF30-50.....	154	CUV2065-1.....		173	CX-071.....	66	DRIFEF-L.....	158
Binocular Body.....	183	CF30-100.....	154	CUV2066-1.....		173	CX-072.....	66	Driven Guard Shield... 80, 81, 82, 83	
Binocular Microscope.....	187	CF30-500.....	154	CUV2071-1.....		173	CX-073.....	66	DRL.....	85
Binocular SurgioScope.....	185	CF30-1000.....	154	CUV2102-1.....		173	CX-074.....	66	Dropchip.....	68
Bio-amplifier.....	74, 75	CFN10-50.....	154	CUV2105-1.....		173	CX-075.....	66	Dry Sterilizer.....	210
Biochip.....	68	CFN10-100.....	154	CUV2106-1.....		173	CX-076.....	66	Dual Channel Differential	
Biofluorescence.....	21	CFN10-250.....	154	CUV2107-1.....		173	CX-077.....	66	Electrometer.....	82
Biofluorometer 14, 18-20, 179,		CFN30-50.....	154				CX-078.....	66	Dual Emission.....	179
BioOxy Sensor.....	152	CFN30-100.....	154				CX-079.....	66	Dual Manipulator Digital Stereotaxic	
Biopsy Punch.....	204	CFN30-250.....	154				CX-080.....	66	Frame.....	119
Biosensing.....	143	CFN30-1000.....	154				CX-081.....	66	Dual Manipulator Stereotaxic	
Biosensors.....	145	Chamber, Measurement.....	160				CX-082.....	66	Frame.....	120
Bipolar Electrodes.....	90						CX-083.....	66		

Item	Page	Item	Page	Item	Page	Item	Page	Item	Page
Dual Microprobe System.....	80	EZ-830.....	129	FORT25.....	79, 89	Hydrogen Sulfide.....	145	Iris Forceps.....	195, 208
Dual Rate Syringe Pump.....	64	EZ-1130.....	129	FORT100.....	79, 89	Hydrogen Sulfide Macro Sensor	148	Iris Scissors.....	199
Dual Tool-Holder Micromanip....	29	EZ-1785.....	129	FORT250.....	79, 89	Hydrogen Sulfide Micro Sensor	148	Iris SuperCut Scissors.....	201, 203
DUMC22.....	134	EZ-7000.....	128	FORT1000.....	79, 89	Hydrogen Sulfide Sensor.....	148	IRM23E10.....	91
DUMC26.....	134	EZ-20027.....	129	FORT5000.....	89	Hydrostatic Pressure.....	116	IRM23E15.....	91
Dummy Load Resistor Kit.....	85	EZ-20027G.....	129	FORT Transducers.....	89			IRM23E20KT.....	91
Dumont.....	194	EZ-20028.....	129	Frame Clamp.....	223	<b>I</b>		IRM23E25.....	91
Dumont #3C.....	194	EZ-20028G.....	129	Frame, Stereotaxic.....	122	IE010.....	159	IRM23E30.....	91
Dumont #4.....	194	EZ-20029.....	129	FrameWorks.....	223	IE190.....	159	IRM23E30KT.....	91
Dumont #5.....	194	EZ-20030.....	129	Free Radical Analyzer.....	150, 151	IE200.....	159	ISDB.....	79
Dumont Forceps.....	194, 208	EZ-20030G.....	129			IGS100.....	149	ISDB Blank Panels.....	79
Dumostar.....	194	EZ-20032.....	129	<b>G</b>		II-39.....	138	ISO2 Adapter.....	221
Dumostar Forceps.....	194	EZ-20032G.....	129	Gas Anesthesia.....	123, 124, 126	II-78.....	131	ISO2 Filling Solution.....	149
Duo773.....	76, 80	EZ-20034.....	129	Gasket.....	101	II-79.....	131	ISO-80.....	76, 78
Duo 773 Electrometer.....	81	EZ-20034G.....	129	Gaskets for PicoNozzle.....	45	II-80.....	131	ISO-80P.....	78
Dye Loading.....	179	EZ-Anesthesia.....	128	GC22.....	134	II-81.....	131	ISO-COP-2.....	145
		EZ-B800.....	128	GC24.....	134	II-82.....	131	Iso-DAM8A.....	76, 79
<b>E</b>		EZ-E25000.....	129	GC26.....	134	II-83.....	131	ISO-H2S-2.....	145, 148
E210.....	107	EZ-E27000.....	129	Glass Capillaries....	43, 47, 104, 217	II-84.....	131	ISO-H2S-100-Cxx.....	145, 148
E212.....	107	EZ-E28000.....	129	Glass Fiber Filters.....	170	II-84XL.....	131	ISO-HPO-2.....	145, 148
E215.....	107	EZ-FF900.....	128	Glass Handling Forceps.....	106	II-336T.....	138	ISO-HPO-2 Electrolyte.....	148
E220.....	107	EZ Lever.....	203	Glass Micropipettes.....	107, 218	II-336TG.....	138	ISO-HPO-100.....	145, 148
Ear Bars.....	124			Glass Rod.....	105	II-390.....	138	ISO-HPO-100H.....	148
Economy Manual Micromanip... ..	25	<b>F</b>		Glucose Sensor.....	149	II-390G.....	138	ISO-HPO-100-L.....	145, 147
ECU-H5.....	35	Facemask.....	125	Glue Gun.....	105, 142, 216	II-520MR.....	136	ISO-HPO Startup Kit.....	148
ECU-HC.....	35	FC-H5.....	35	GN-NET7.....	65	II-600MR.....	137	Isolated Bioamplifier.....	78
ECU-HOC.....	35	FC-HC.....	35	GN-NET25.....	65	II-2200.....	136	Isolated Differential Amplifier.....	78
EHB1.....	98	FC-HOC.....	35	GN-PC7.....	65	II-2390.....	135	ISO-NO Activator.....	160
EHBF.....	99, 101	FC-SMT.....	35	GN-PC25.....	65	II-2391.....	135	ISO-NOP.....	145, 146, 147
EJA.....	191	FD35-100.....	225	GN-TTL.....	65	II-2392.....	135	ISO-NOP007.....	145, 147
EK1.....	9	FD35B-100.....	225	GO1-100.....	101	II-2393.....	135	ISO-NOP3005.....	145, 147
EKC.....	9	FD35PDL-100.....	225	GO2-100.....	101	II-2394.....	135	ISO-NOP3020.....	145, 147
EKC Cartridges.....	9	FD223A.....	82	GO3-100.....	101	II-2395.....	135	ISO-NOPF.....	146
EKE.....	191	FD223AP.....	82	GO4-100.....	101	II-2396.....	135	ISO-NOPF100.....	145, 147
EKV.....	9	FD3510-100.....	225	Gold Plated Pins, Sockets.....	91	II-2397.....	135	ISO-NOPF100-L10.....	147
EKV Cartridges.....	9	FD3510B-100.....	225	Gold Wire.....	102, 222	II-2400.....	135	ISO-NOPF200.....	145, 147
Electro 705.....	76, 83	FD5040-100.....	225	GR100-4.....	105	II-2500.....	136	ISO-NOPF200-L10.....	147
Electrode Adapter.....	75	FD5040B-100.....	225	GR100-6.....	105	II-2888.....	135	ISO-NOPF500.....	147
Electrode Holder....	33, 74, 78, 159	Fiber Optic Assemblies.....	162	Grip Strength Meter.....	136	II-23931.....	135	ISO-NOPF500-Cxx.....	145
Electrolyte, ISO-NO.....	147	Fiber Optic Cables.....	157, 165	Ground Wire.....	97	II-B60-025.....	131	ISO-NOPNM.....	145, 146, 147
Electrometer.....	83	Fiber Optic Collimator.....	171	GSNO.....	160	II-B60-038.....	131	ISONOP Startup Kit.....	147
Electrophysiology.....	73	Fiber Optic Filter Holder.....	170	GSNO-50.....	160	II-B60-044.....	131	ISO-OXY-2.....	145, 149
Electrosurgical Unit.....	209	Fiber Optic Holder.....	168, 170	GSNO-100.....	160	II-B63-075.....	131	ISO-TEMP-2.....	151
Elgiloy Microelectrodes.....	90	Fiber Optic Illuminator.....	179, 181	Guide Cannula.....	134	II-B64-1.....	131	IT-1E.....	127, 133
Endohm.....	5	Fiber Optic Light Source.....	177	Guillotine.....	206	II-MRBP-M.....	131	IT-14.....	133
ENDOHM-6G.....	5, 7, 9	Fiber Optic Oxygen Meter.....	152, 155	Guinea Pig Adaptor.....	123	II-MRBP-R.....	131	IT-18.....	127, 133
ENDOHM-12G.....	5	Fiber Optic pH Meter.....	157			II-PE34.....	137	IT-18EXLONG.....	133
ENDOHM-24GSNAP.....	5	Fiber Optic pH System.....	157	<b>H</b>		ILS005LT.....	55	IT-21.....	133
Environmental Control.....	34	Filament.....	110, 111	Half-Cells.....	98, 100	ILS010LT.....	55	IT-23.....	127, 133
EP1.....	74, 97, 103	Filament Adjustment Assembly.....	114	Halogen-Halogen Dual Illuminated		ILS025LT.....	55	<b>J</b>	
EP2.....	103	Filament Cable.....	114	Track Stand.....	183	ILS Gas-Tight Syringes.....	55	Jaffe Forceps.....	203
EP4.....	103	Filament for Puller.....	48	Halogen Lamp.....	176, 191	Impedance Measurement.....	88	Joystick Control.....	29
EP05.....	103	Filaments, Platinum/Iridium.....	108	Halogen Replacement Lamp.....	176	Implantable Probes..	127, 133, 155	Julabo Circulating Bath.....	9
EP08.....	103	Filling Solution.....	145	Halsted Mosquito Forceps.....	195	INC2.....	134		
EP12.....	103	Filter Holder.....	170	HD imaging.....	188	INC24.....	134		
Epithelial Tissue Voltohmmeter....	4	Fire-Polished.....	47, 104, 217	HDMI Cable.....	190	INC26.....	134	<b>K</b>	
Epoxy.....	139, 142, 213, 216	Fixed Flow System.....	128	Heated Waterbed.....	129	Incapacitance Meter.....	137	Katena-Vannas Scissors.....	198
Ethyl Cyanoacrylate.....	139	Flaming/Brown Pipette Puller ..	111	Heater Controller.....	37	Incremental Hot Cold Plate Analgesia		KG Transducers.....	22
Euthanex Lids.....	129	Flexible Dri-Ref.....	158	Heating Filament.....	113	Meter.....	137	Kima Pump.....	66, 68
EVOM2.....	4, 5	Flexible Implantable Probe.....	127	Heating Plate.....	124, 127	Induction Chambers.....	129	KITE-L.....	25
Exigo.....	66	Flexible Light Guide.....	191	Homostatic Forceps.....	195	Infinity Plan Achromat.....	187	KITE-M3-L.....	25
EXT-6.....	132	Flexible NO Sensor.....	147	High Current Isolator.....	86	Infusion Pump.....	63	KITE-M3-R.....	25
Extension.....	97	FLEXREF.....	158	High Performance Structural Epoxy ..		Injection Assembly Parts Kit.....	46	Kite Manual Manipulator	24, 25, 41
Extension Cable.....	79, 130	Flow Cell, Low Volume.....	167	.....	139	injector.....	40	KITE-R.....	25
Eye Guards.....	181	FLOX-PATCH.....	153	High Voltage Isolator.....	85	Inlet and Outlet Hoses.....	37	KITLWCC.....	167
Eyepiece Adapter, SLR.....	190	FLOX-PROBE.....	152, 153	Horizontal Tissue Bath.....	11	In-Line Fiber Optic Filter Holder	170	KIT-PHYSIO-I.....	207
Eyeieces.....	181, 182, 186, 187	Fluorescence Spectroscopy.....	19	Hose Assembly.....	37	In-line Frame Clamp.....	223	KWIK-2.....	159
Eyeiece with Reticle.....	181, 183	FluoroDish.....	41, 225	HOT Plate Analgesia Meter.....	138	Input Kit.....	45	KWIKCAL.....	159
EZ-103A.....	129	Flush Cutters.....	102	HPO Sensor.....	147, 148	In Situ Oxygen Monitoring	152, 153	KWIK-CAST.....	139, 140, 213, 214
EZ-104A.....	129	FO-6000.....	177	HT-1.....	133	Internal Cannula.....	134	KWIKGARD.....	139, 213
EZ-107A.....	129	Foam Ring.....	107	HT-2.....	133	Intracellular Amplifier.....	80, 83	Kwik-Gard Dispenser.....	213
EZ-109.....	129	FOIMPH.....	99, 101	HTS Electrodes.....	6	INV-101.....	18, 186	Kwik-Gard Refill.....	213
EZ-177.....	129	Foot Switch.....	45, 55, 81, 84, 113, 114, 209	Humidifier.....	35, 37	Inverted Microscope.....	18	Kwik-Gard Start-up Kit.....	213
EZ-178.....	129	Forceps.....	196	Humidifier Kit.....	37	Inverted Trinocular Microscope	186	KWIKGLUE.....	139, 213
EZ-179.....	129	Forceps with Tying Platform.....	196	Humidity Probe.....	37	IO-KIT.....	58	KWIKGUN.....	139, 213
EZ-211.....	129	Force Transducers.....	15, 22, 79, 89	Hydrogen Electrodes.....	159	Ion Selective Electrodes.....	159	KWIKH-2.....	159
EZ-320.....	129	FORT10g.....	79, 89	Hydrogen Ion Exchanger.....	159	Iridium Microelectrodes.....	90	KWIKMIX.....	139, 213
EZ-330.....	129			Hydrogen Peroxide Sensors	145, 148	Iris Dissecting Scissors.....	200	KWIK Mixing Tips.....	140





Item	Page	Item	Page	Item	Page	Item	Page	Item	Page	
PERIPRO-4LS .....	11, 52, 167	Pre-Pulled Pipettes.....	41	Rat Stereotaxic Non-Rebreathing Unit .....	129	RTV Coating .....	142, 216	SI-TCM2 .....	16	
PERIPRO-8LS .....	11, 52	Pressure Calibration Chamber..	115	RBMA-200C.....	206	RTV Prime Coat.....	142, 216	SI-TCM2B .....	13, 16	
Peristaltic Pump .....	52, 53, 167	Pressure Manometer.....	115, 116	RBMA-200S.....	206	RTV Sealant .....	139, 142, 216	SI-TM1.....	22	
Peri-Star Pro.....	52	Pressure Monitor .....	130	RBMA-300C.....	206	Rubber Ear Bars .....	124	SI-TM2.....	22	
PE Tubing, Flexible .....	134	Pressure Pod.....	115	RBMA-300S.....	206	<b>S</b>				
PG10150-4 .....	106	PRO-300 HDS .....	41	RBMA-600C.....	206	Sample Cell .....	164	SI-TM3.....	22	
PG10165-4 .....	106	PRO-300HDS .....	188	RBMA-600S.....	206	Sample Injector Attachment.....	167	SI-TM4.....	22	
PG52151-4 .....	106	Probe.....	75, 115, 208	RBMS-200C.....	206	Scalpels, Sterile .....	205	SI-TM5.....	22	
PG52165-4 .....	106	Probe Extension Cable.....	115	RBMS-200S.....	206	Scissors .....	197	SI-TM6.....	22	
pH Implantable Sensor.....	157	Probe Extension Lead.....	132	RBMS-300C.....	206	Scotch-Weld .....	139, 213	SI-TM7.....	22	
pH Meter.....	157	Probe Handle.....	81, 83	RBMS-300S.....	206	SDR2 .....	158	SI-TM8.....	22	
pH Meter, Fiber Optic.....	157	Programmable Filter Module .....	13, 17	RBMS-600C.....	206	Semi-Micro Concentric Bipolar Electrodes.....	92, 93	SI-TM9.....	22	
pH Mini Sensors.....	157	Programmable Syringe Pump.....	65	RBMS-600S.....	206	Septum Theta.....	47, 104, 105, 217	SI-TM10.....	22	
Phone Plug-to-BNC Cable.....	75	PSMB5N.....	185	RC1.....	103	Serial Cable .....	97	SI-TM11.....	22	
PH-OPTICA-MICRO.....	157	PSMT5N.....	185	RC1T.....	83, 97, 103	SGE001RN.....	55	Skin Staplers .....	203	
PH-OPTICA-MINI.....	157	PT0110.....	102, 222	RC2.....	103	SGE005RN.....	55	Slides, Glass.....	218	
Photodiode Array Based Spectrophotometers.....	174	PT0203.....	102, 222	RC2F.....	103	SGE010RNS.....	55	SMA Bulkhead.....	171	
Photonics Center of Excellence .....	162	PT0402.....	102, 222	RC3.....	103	SGE025RN.....	55	Small Tissue Probe .....	179	
pH System.....	157	PTM3C02INS .....	91	RC6.....	103	SGE050RN.....	55	SMARTouch.....	55	
PHYSIO-II.....	207	PTM04RC100.....	92	Rectal Probes .....	133	SGE050TLL.....	147, 151, 160	SNAP.....	151	
Physiology Kit I.....	207	PTM041RC225.....	92	Rectangular Base .....	223	SGE100RN.....	55	SNAP25.....	160	
Physiology Kit II.....	207	PTM321RC100.....	93	Rectangular Base Post Stand....	183	SGE250TLL.....	55	SNAP50.....	160	
PicoNozzle Kit .....	41, 45	PTM321RC225.....	93	Reference Cell .....	82, 83, 103	SGE Syringes .....	55	SNAP100.....	160	
PicoPump.....	40, 43, 44, 45	PTP201.....	102, 222	Reference Electrodes.....	158	Shielded Electrode Cable .....	75	S-nitroso-L-glutathione .....	160	
Piggyback Glass....	47, 104, 105, 217	PTP401.....	102, 222	Reflex Clip .....	205	Shielded Modular Cables.....	75	S-Nitroso-N-acetyl-D,L-penicillamine .....	160	
Pinger.....	17	PTP406.....	102, 222	Reflex Clip Applier.....	205	SI-AOSUB.....	11, 13, 17	Socket .....	97	
Pinpoint Injection.....	42	PTPxxx.....	102, 222	REG2 .....	70	SI-BAM21-LC .....	15	SO-OXY Startup Kit.....	149	
Pinpoint Injector.....	40, 42	PTT0110.....	102, 222	REG5K.....	70	SI-BAM21-LCB .....	13, 15	SP100i.....	63	
Pipetters.....	70	PTT0203.....	102, 222	REG10.....	70	SI-BF-100.....	18, 178	SP100iZ.....	63	
Pipettor Stand .....	70	PTTxxx.....	102, 222	REG10K.....	70	SI-BF-100LLG.....	21, 179	SP101i.....	63	
Pipette Tips .....	71	PTxxxx.....	102, 222	REG20.....	70	SI-BF-100SMA .....	21, 179	SP101iZ.....	63	
Planachromatic Objective .....	181	PUL-1000.....	41, 48, 108	REG50.....	70	SI-BF-LLG-UPGRADE .....	21, 179	SP120p.....	63	
Plan Phase Contrast.....	187	Puller ....	43, 48, 108, 109, 110, 111	REG100.....	70	SI-BF-SMA-UPGRADE.....	21, 179	SP120pZ.....	63	
Plantar/Tail Flick Meter .....	138	PV820.....	40, 44, 45	REG200.....	70	SI-BMFA.....	13	SP200i.....	63	
Plantar Test .....	138	PV830.....	40, 43, 44	REG1000.....	70	SI-COLUB.....	11, 13, 16	SP200iZ.....	63	
Plantar Test Analgesia Meter ....	138	PZMIII.....	41, 183	Regal Pipetters .....	70	SI-CTS200.....	12, 18	SP210c.....	63	
Platinum/Iridium Microelectrodes .....	90	PZMIII-AAC.....	183	REG55.....	70	SI-CTS200A.....	12	SP210cZ.....	63	
Platinum/Iridium Teflon Wire.....	102, 222	PZMIII-BS.....	183	REG56.....	70	SI-CTS200B.....	12	SP210iw.....	63	
Platinum/Iridium Wire .....	102, 222	PZMIII-MI.....	41, 184	REG57.....	70	Signal Conditioning Amplifiers....	13	SP210iwZ.....	63	
Platinum Wire .....	102, 222	PZMIV.....	181	REMS-24.....	7	Signal Conditioning Amplifier System.....	13, 14, 16	SP220i.....	63	
Plethysmometer .....	136	PZMIV-BS.....	181	REMS-24M.....	7	Signal Generator .....	84	SP220iZ.....	63	
PM01.....	116	PZMTIII.....	41, 183	REMS-96.....	7	SI-HTB2.....	11	SP230iw.....	63	
PM01D.....	116	PZMTIII-AAC.....	183	REMS-96C.....	7	SI-HTB4.....	11	SP230iwZ.....	63	
PM01R.....	116	PZMTIII-AAC-DIG50.....	183	REMS AutoSampler .....	7	SI-KG2.....	15, 22	SP250i.....	63	
PM01S.....	116	PZMTIII-BS.....	183	Replacement Electrode Holder.....	25, 29	SI-KG2B.....	15, 22	SP250iZ.....	63	
PM015.....	116	PZMTIII-BS-DIG50.....	183	Replacement Electrode Sleeve..	149	SI-KG4.....	15, 22	SP260p.....	63	
PM015D.....	116	PZMTIII-BS-LWD-DIG50.....	183	Replacement Hoses.....	37	SI-KG7.....	22	SP260pZ.....	63	
PM015R.....	116	PZMTIII-DIG50.....	184	Replacement Lamp.....	177	SI-KG7a.....	15, 22	Spare Plug-Adapter .....	160	
PM100.....	116	PZMTIII-MI.....	184	Replacement Needles.....	55	SI-KG7b.....	15, 22	Specialty Micromanipulators.....	29	
PM100D.....	116	PZMTIV.....	181	Replacement Needle Kit....	147, 148	SI-KGxx.....	22	Spectrophotometers.....	174	
PM100R.....	116	PZMTIV-BS.....	181	Replacement Sleeves.....	145	SILFLEX-2.....	57	SPLG100.....	61	
PM Series.....	116	PZMTIV-BS-DIG50.....	181	Replacement Tubing.....	52	SilFlex Tubing.....	57	SPLG101.....	61	
Pneumatic PicoPump.....	42, 43, 44	PZMTIV-BS-LWD-DIG50.....	181	Restrainer.....	131	Silicone Adhesive.....	140, 214	SPLG110.....	61	
Polishing Patch Pipettes.....	113	QFT1.....	165, 170	RET-2.....	127, 133	Silicone Casting Compound.....	140, 214	SPLG210.....	61	
PolyFil.....	105	QUAD.....	28	RET-3.....	127, 133	Silicone Dissecting Pad Kit .....	215	SPLG212.....	61	
Polymer Core Fiber.....	162	Qualified Optical Fiber.....	162	Rhodes Style Tips .....	92, 93	Silicone Gasket.....	57	SPLG270.....	61	
Polymer Optical Fiber .....	155, 157	Quantrex.....	211	Ring-Handle Needle Holders....	202	Silicone Pad for ATC2000.....	127	SPLG272.....	61	
Portfolio .....	208	Quattro.....	135	Ring Light .....	191	Silicone Rubber Sealant .....	101	SPLG-G88PLUS.....	64	
Post.....	223	Quattro Tubing.....	57	Ring Light Adapter ....	183, 189, 191	Silicone Tubing.....	52, 53	SPL Syringe Pump .....	61	
Post Extension.....	183	Quattro 4-in-1 System .....	135	Ring Light Guide.....	179, 181	Silver Epoxy.....	142, 216	Spring Scissors.....	198	
Post Mounting Kit.....	9	QVLUX-DUV.....	163	Ring Scissors.....	200	SI-MOTDB.....	13, 14	SS316xx.....	102, 222	
Post Stand.....	183	QVLUX-NIR.....	163	Ringstand Mounting Kit.....	75	SI-MOT-MT.....	13, 14	SS31614.....	102, 222	
Potassium Electrodes .....	159	<b>R</b>				RN001.....	SI-MT-L.....	10	SSM04RC100.....	92
Potassium Ion Exchanger.....	159	R-8-8-WPI01.....	181, 191	Rabbit/Guinea Pig Chamber.....	129	SI-MTM-L.....	10	SSM33A20KT.....	91	
Power Frame Enclosure.....	13	Rack Mount Kit... 45, 75, 79, 81, 84, 88, 115, 130		Rodent Brain Matrices.....	206	Single-Barrel Standard Borosilicate Glass Tubing.....	42, 47, 104, 217	SSM33A70.....	91	
Precious Metal.....	222	Randall Selitto Paw Pressure Meter .....	136	Rodent Restrainer .....	131	Single Emission .....	179	SSM041RC225.....	92	
Precision Current Sources.....	85	Rapid-Punch Biopsy Kit.....	204	Rod Holder.....	27, 28	Sintering Pellets.....	103	SSM321RC100.....	93	
Precision Stereotaxic Frame.....	120	Rat Adaptor .....	123	Room Temperature Vulcanizing .....	139	SI-OSARC.....	18	SSM321RC225.....	93	
Precision Stereo Zoom Binocular Microscope.....	181	Rat Dissecting Kit.....	208	Rotating Base .....	27, 28	SI-PF100.....	11, 17	SST30407-50.....	102, 222	
Precision Stereo Zoom Microscope... ..	183	Rat KIT.....	208	RPE-KIT.....	58	SI-PF100B.....	13, 17	SSTxxxx.....	102, 222	
Precision SurgioScope .....	185	Rat Spinal Adaptor.....	123	RTD Rectal Temp Probe .....	127			ST5191.....	210	
Pre-Polarizer.....	160			RTD Temperature Sensor.....	151			ST5193.....	210	
Pre-Pulled Glass Pipettes.....	42			RTV.....	139			Stage Micrometer.....	185, 187	
Pre-Pulled Micropipette .....	115							Stagetop Environmental Control .....	36	

Item	Page	Item	Page	Item	Page	Item	Page	Item	Page
Stand Alone Optical Force Transducer Amplifier.....	15	SYS-A365R.....	85	TIPSTW1LS02.....	107, 218	TTL Control Module.....	52, 53	<b>V</b>	
Standard Cuff.....	94	SYS-A382.....	86	TIP10FLT.....	107, 218	Tubing Cartridge.....	52	V503382.....	202
Steel Base Plate.....	32	SYS-A385R.....	86	TIP10LT.....	107, 218	Tubing Kit.....	9	Vacuum Pump.....	41, 115
Stereo Microscope.....	184	SYS-A385SRC.....	86	TIP10TW1.....	107, 218	Tungsten Carbide SuperCut Scissors.....	201	Vacuum Pump Regulation Kit ...	115
Stereo Microscope with LED Illuminated Base.....	184	SYS-BP1.....	130	TIP10TW1-L.....	107, 218	Tungsten Light Source.....	177	Vannas Capsulotomy Scissors, Titanium.....	199
Stereotaxic.....	120	SYS-DAM50.....	75	TIP10TW1LS01.....	107, 218	Tungsten Microelectrodes.....	90	Vannas Scissors.....	197, 199, 208
Stereotaxic Accessories.....	123	SYS-DAM80.....	75	TIP10TW1LS02.....	107, 218	Tungsten Wire.....	102, 222	Vaporizer Pole Mount.....	129
Stereotaxic Anesthesia/Evacuation Mask, Mice or Neonatal Rats	126	SYS-OMEGAZ.....	88	TIP10XV119.....	59, 107, 218	TW100-3.....	47, 104, 217	Vascular Hooks.....	22
Stereotaxic Anesthesia/Evacuation Mask, Rats.....	126	SYS-PM01D.....	115, 116	TIP12FLT.....	41	TW100-4.....	47, 104, 217	V-Base.....	223
Stereotaxic Anesthesia/Facemask	125	SYS-PM01R.....	116	TIP15FLT.....	41, 107, 218	TW100-6.....	47, 104, 217	V-Clamp.....	124
Stereotaxic Frame.....	120, 121	SYS-PM015D.....	115, 116	TIP15LT.....	107, 218	TW100F-3.....	47, 104, 217	Ventilator Connection Kit.....	129
Stereotaxic Holder.....	209	SYS-PM015R.....	116	TIP30TW1.....	107, 218	TW100F-4.....	47, 104, 105, 217	Versaflex Non-Rebreathing Unit	129
Stereotaxic, Mice and Neonatal Rats.....	121	SYS-PM100D.....	116	TIP30TW1-L.....	107, 218	TW100F-6.....	47, 104, 217	Vertical Surface Mount.....	223
Stereo Zoom Microscope... 180, 182		SYS-PM100R.....	116	TIP30TW1LS01.....	107, 218	TW120-3.....	47, 104, 217	Vessel Cannulation Forceps.....	208
Steri-Lite Micro Bead Sterilizer .	210	SYS-PV820.....	42, 45	TIP30TW1LS02.....	107, 218	TW120-4.....	47, 104, 217	Vessel Clips.....	205
Sterilization Mesh Cassette.....	211	SYS-PV830.....	42, 45	TIPCA.....	159	TW120-6.....	47, 104, 217	VEVBOND.....	139, 141, 213, 215
Sterilizers.....	210	SYS-REMS.....	7	TIPH.....	159	TW120F-3.....	47, 104, 217	Vibration-Free Tables.....	38
ST-H5.....	36	SYS-TBM4M.....	88	TIPK.....	159	TW120F-6.....	47, 104, 217	Vibration-Free Workstation.....	38
ST-HC.....	36	<b>T</b>		TIPMIX01-05.....	107, 218	TW150-3.....	47, 104, 217	Vibration Isolation Platform.....	38
ST-HOC.....	36	T-Adapter Kit.....	147	TIPMIX01-05-L.....	107, 218	TW150-4.....	47, 104, 217	Vibration Platform.....	38
Stimulator.....	84	Tail Cuff Sensor.....	131	TIPMIX05-10.....	107, 218	TW150-6.....	47, 104, 217	Video Camera.....	190
Stimulus Isolator.....	86, 87, 88	Tail Flick Meter.....	138	TIPMIX05-10-L.....	107, 218	TW150F-3.....	47, 104, 217	VIS-NIR.....	162
Stopcock.....	220	Tail Flick Test.....	138	TIPTPP.....	159	TW150F-4.....	47, 104, 105, 217	Von Frey Anesthesiometer.....	135
Storage Jar.....	107	TAXIC-300.....	119	TIPTW900A.....	115	TW150F-6.....	47, 104, 217	Von Frey Probe.....	135
Storage Portfolio.....	208	TAXIC-303.....	119	Tissue Bath, Horizontal.....	11	Tygon Tubing.....	68		
Structural Epoxy.....	213	TAXIC-350.....	119	Tissue Mounting Hooks.....	22			<b>W</b>	
STX2.....	4, 5	TAXIC-353.....	119	Tissue Probe.....	21			W30S-LED.....	187
STX3.....	4	TAXIC-600.....	120	Titanium Forceps.....	196			W30ST-LED.....	187
STX100.....	6	TAXIC-603.....	120	Titanium Scissors.....	199	UBATH-Y.....	211	Wall-Mount Plate.....	183
STX100C.....	6	TAXIC-650.....	120	T-joint Frame Clamp.....	223	UltraMicroPump.....	55	Waterbed.....	129
STX100C96.....	6	TAXIC-653.....	120	TM31A10.....	91	UltraPath.....	164	Waveguide Cleaning Kit....	167, 171
STX100F.....	6	TBM4M.....	88, 89	TM31A20.....	91	Ultraphath System.....	165	Wide Field.....	183
STX100M.....	6	TBR1025.....	151	TM33A05.....	91	Ultrasonic Cleaner.....	211	Wire.....	102, 222
Sub-Microliter Injection.....	56	TBR4100.....	150, 151	TM33A10.....	91	Ultrasonic Cleaning System.....	211	Wire Cutters, Ergonomic.....	102
SU-P97.....	41, 48, 111	TBS.....	33	TM33A10KT.....	91	UMP3.....	55	Wire Cutters.....	102, 222
SU-P1000.....	41, 48, 110	TC Mouse Rectal Temp Probe... 127		TM33A20.....	91	UMP3-3.....	55	Wire Retractor.....	208
SU-P2000.....	109	TC Rat Rectal Temp Probe..... 127		TM33B01.....	91	UMP3-4.....	55	WP1010.....	204
Super Adhesive.....	141, 215	TEER Measurement..... 4, 6, 7, 8		TM33B01KT.....	91	UMP3T-1.....	55	WP1212.....	204
SuperCut Iris Scissors.....	208	Telemetry Training..... 208		TM33B05.....	91	UMP3T-2.....	55	WP2000.....	196
SuperCut Scissors.....	201	Temperature Control Module 13, 16		TM33B10.....	91	UMS.....	38	WP2020F.....	204
SUPER-Dri-Ref.....	158	Temperature Probes..... 37, 133		TM33B20.....	91	UniGo.....	67	WP3030.....	204
Support Stand.....	9	Temperature Sensor..... 151		TM33C10.....	91	Universal Focus Mount.... 179, 181		WP3620.....	196
SU-QUAD.....	24, 41	Test Resistor..... 4		TM33CC10.....	91	Universal Manipulator Stand.... 38		WP4410.....	199
SU-QUAD-L.....	28	TGW0515.....	102, 222	TM33CCINS.....	91	Universal Pipette Tips..... 71		WP5000.....	196
SU-QUAD-R.....	28	TGW1510.....	102, 222	TM33CCNON.....	91	Univette.....	10	WP5050.....	199
Sure-Seal Chambers.....	129	TGWxxxx.....	102, 222	TM53CCINS.....	91	UPUV.....	165	WP5060.....	199
Surgical Instruments.....	192	Thermocouple Thermometer... 132		TPP+ (Tetraphenylphosphonium) Electrodes.....	159	UPVIS.....	165	WP5070.....	199
Surgical Instrument Storage Portfolio.....	208	Thin-Line Heated Waterbed.... 129		TPP Electrode.....	159	USBCAM33.....	189	WP5070L.....	199
Surgical Kits.....	207	Thin Wall.....	47, 104, 217	Transbridge.....	88	USBCAM50.....	189	WPB315AB.....	203
SurgioScope.....	185	TIDAS-D2.....	176	Transbridge Transducer Amplifier 88		USBCAM133.....	189	WPB315AS.....	203
SU-TRIO.....	41	TIDAS-DAQ.....	174, 175	Transducer Amplifier..... 79, 88, 89		USBCAM152.....	189	WPB505210.....	203
SU-TRIO235.....	24, 27	TIDAS E Base.....	175	Transducer Amplifier Module.... 79		USBCAM202.....	189	WPSS01.....	203
SU-TRIO245.....	24	TIDAS-H.....	176	Transmission Probes..... 171		USBS1L.....	9	WPSS02.....	203
SU-TRIO245-L.....	27	TIDAS-H2.....	176	Trinocular.....	187	USS1S.....	9	WVLUX-DUV.....	163
SU-TRIO245-R.....	27	TIDAS S300.....	175	Trinocular Body.....	183	USS2L.....	9	WVLUXHNA.....	163
Swiss Tweezers.....	193	Tilting Base.....	25, 33, 82, 83	Trinocular Microscope..... 181, 187		USS3L.....	9	WVLUXNHR.....	163
Sylgard.....	139, 141, 213, 215	Tinned Cu Alloy.....	102	Trinocular Microscope with LED Illuminated Base.....	184	USS3S.....	9	WVLUXXPV.....	163
Syringe Multi Rack.....	61	TIP1TW1.....	107, 108, 42	Trinocular SurgioScope.....	185	USS4L.....	9	WVLUXUVIS.....	163
Syringe Pumps.....	60, 62, 63, 65	TIP01TW1F.....	107, 218	Trinocular Tube.....	181	USS4S.....	9		
Syringes with Luers.....	55	TIP01TW1F-L.....	107, 218	Trio Analgesia Meters.....	135	USS5L.....	9	<b>X</b>	
SYS-705.....	83	TIP1TW1-L.....	107, 218	Trio 3-Axis Manipulator..... 27		USS5S.....	9	X-Wide Contact Nerve Cuffs.....	95
SYS-773.....	81	TIP2TW1.....	107, 218	Trio 3-in-1 System.....	135	USS6L.....	9		
SYS-900A.....	115	TIP02TW1F.....	107, 218	Tripod Stand.....	132	USS6S.....	9	<b>Z</b>	
SYS-A310.....	84	TIP02TW1F-L.....	107, 218	Trough Filament.....	110, 111	USS7L.....	9	Z-Axis Vertical Extension..... 27, 28	
SYS-A362.....	87, 88	TIP03TW1F.....	107, 218	Troutman Needle Holder.....	196	USS7S.....	9	Z-Dimensions.....	172
SYS-A365D.....	85	TIP03TW1F-L.....	107, 218	TST33A001KT.....	91	USS8L.....	9	Zebrafish.....	40, 46
		TIP04TW1F.....	107, 218	TST33A05KT.....	91	USS8S.....	9	Zhang, Xueji.....	144
		TIP5TW1.....	107, 218	TST33A10KT.....	91	Ussing Current Electrode.....	9	Z-LITE.....	88, 181, 191
		TIP05TW1F.....	107, 218	TST33A20KT.....	91	Ussing Electrode Kit.....	9	Z-LITE-186.....	181, 191
		TIP05TW1F-L.....	107, 218	TST33C05KT.....	91	Ussing Electrode Kit.....	9	Z-LITE-Z.....	191
		TIP5TW1-L.....	107, 218	TST33C05KT.....	91	Ussing System.....	8, 9	Z-LITE-Z186.....	191
		TIP5TW1LS01.....	107, 218	TST53A10KT.....	91	Ussing Voltage Electrode.....	9	Z-MOLDS.....	41, 46
				TST150-6.....	105	UV-VIS.....	162	Zoom Microscope.....	180, 182

# TERMS & CONDITIONS

WPI values its customers and will make every effort to maintain the highest professional and ethical standards. However, WPI regards as invalid any special terms and conditions set forth by any other party, including purchasing, quotations, warranties or other statements. WPI's response to a quotation, purchase order or inquiry shall not be construed as acceptance of any conditions stated therein.

## Credit Terms

Net 30 days from date of invoice to customers with satisfactory credit references.

## Standard Equipment & Accessories

To insure against misunderstandings or delays, always include model number, name, and significant specifications of systems and accessories when ordering. Also specify any options desired. WPI retains the right to change specifications or discontinue products without prior notice and without responsibility for incorporating such changes into equipment previously sold. Please state the desired operating line voltage and frequency.

## Warranty

WPI strives to maintain the highest quality standards in all products. We warranty these products against defects in workmanship or materials, and all WPI products retain a warranty that is stated in the product's manual. (Typically, one year for parts and labor except for consumable items such as electrodes, glass, etc., for which warranty varies from 30 days to six months, or as otherwise specified.) However, WPI hereby specifies that no use or application of any product may be suitable or successful for a specific application and therefore waives all liability in that regard. It is the sole responsibility of the buyer/user to apply the product in a manner consistent with its intended use.

WPI products are not approved for human use unless specifically stated in writing with accompanying FDA (or applicable) documents. WPI assumes no liability or legal, moral, ethical, or fiduciary responsibility for a use of any product in the WPI catalog in treating or studying humans.

## Prices

WPI maintains competitive world wide pricing. Prices listed for standard items as described are net and subject to change without notice. Prices (in U.S. dollars) do not include shipping charges. Formal quotations are valid for a period of 30 days unless otherwise specified. Purchases may be charged to your VISA, MasterCard, or American Express account.

## Taxes

We reserve the right to add any taxes required by local law or ordinance. Institutions operating under tax-free conditions must state applicable license or contract number and be prepared to furnish necessary exemption certificate.

## Packaging and Shipping

Equipment is shipped in adequate commercial packing. All packages should be checked thoroughly before being discarded. Shipping and handling expenses, including adequate insurance coverage, are paid by the customer.

## Claims and Returns

Inspect all shipments upon receipt. Missing cartons or obvious damage to cartons should be noted on the delivery receipt before signing. Concealed loss or damage should be reported at once to the carrier and an inspection requested. All claims for shortage or damage must be made within ten (10) days after receipt of shipment. Claims for lost shipments must be made within thirty (30) days of receipt of invoice or other notification of shipment. Please save damaged or pilfered cartons until claim is settled. In some instances, photographic documentation may be required. Some items are time-sensitive; WPI assumes no extended warranty or any liability for use beyond the date specified on the container.

Do not return any goods to us without obtaining prior approval and instructions from our Returns Department. Goods returned unauthorized or by collect freight may be refused. Goods accepted for restocking will be exchanged or credited to your WPI account and are subject to a restock fee of up to 25%. Goods returned which were ordered by customers in error are subject to a 25% restock fee. Equipment which was built as a special order cannot be returned.

## Repairs

Contact our Returns Department for assistance in the repair of apparatus. Do not return goods until instructions have been received. Returned items must be securely packed to prevent further damage in transit. The Customer is responsible for paying shipping expenses, including adequate insurance on all items returned for repairs. Identification of the item(s) by model number, name, serial number and proof of purchase (packing slip number or invoice number) as well as complete description of the difficulties experienced should be written on the RMA request form and a copy of the form must be included with the item.

# ORDERING INFORMATION

Orders from **North America** should be directed to WPI's main office:

**World Precision Instruments, Inc.**

175 Sarasota Center Boulevard  
Sarasota, Florida 34240-9258

Telephone: (941) 371-1003 (collect calls not accepted). Facsimile copies of purchase orders, inquiries or other correspondence may be transmitted to WPI's Fax number: (941) 377-5428. When sending written purchase orders to confirm a telephone order, please clearly mark "CONFIRMING" on the order to avoid it being shipped a second time. Our Sales Department's e-mail address is: sales@wpiinc.com

North American customers may also order many products online through WPI's Web site — www.wpiinc.com. All customers may use the online Quote Request Form and online Purchase Order Form.

Orders from **Belgium, Denmark, Eire, Finland, Luxembourg, Norway, Portugal, Spain, Sweden, Switzerland (FR) and the United Kingdom** should be directed to

**WPI UK:** 1 Hunting Gate

Hitchin, Hertfordshire SG4 0TJ England

Tel: +44 (0)1462 424700 Fax: +44 (0)1462 424701

E-mail: wpiuk@wpi-europe.com

Customers in **France** may call +33 (0)970 44 9000 or

E-mail: wpifr@wpi-europe.com

Orders from **Austria, Germany, Greece, Italy, Malta, Netherlands, Russian States, Switzerland (DE) and Turkey** should be directed to

**WPI Germany:** Saarstrasse 23

D-61169 Friedberg, Germany

Tel: +49 (0)6031 1602171 Fax: +49 (0)6031 1602180

E-mail: wpi@wpi-europe.com

Orders from **China, Hong Kong and Taiwan** should be directed to

**WPI China:** WPI Shanghai Trading Co., Ltd.

Rm 25e, No8 Dong Fang Rd., Lu Jia Zui Financial District  
Shanghai PRC

Tel: + 86 688 85517

Email: chinasales@china.wpiinc.com

Orders from **Brazil** should be directed to

**WPI Brazil:** Av. Conselheiro Nebias, 756 sala 2611

Santos - Sao Paulo - Brasil

CEP: 11045-002

Email: info@brazil.wpiinc.com

Orders from the countries listed below should be directed to the appropriate WPI distributors. All other international orders should be directed to WPI's main office and factory in Sarasota, Florida.

To expedite shipment of orders outside North America (other than to areas served by the international offices above), payment must be made in advance via wire transfer, check payable in U.S. dollars, or charged to a credit card. Proforma Invoices will gladly be furnished upon request. Please specify proper electrical current for your order.

Contact our distributor for Poland, Baltics,  
Czechia, Slovakia, Romania, Croatia and Slovenia

**ANIMA LAB**  
animal facility and laboratory equipment • animal research models

info@animalab.eu

www.animalab.eu



WORLD  
PRECISION  
INSTRUMENTS

[www.wpiinc.com](http://www.wpiinc.com)