


Conventional cages



Type I L	Type II	Type II L	Type III	Type III H	Type IV	Type IVS	-----
Floor Area:							
360cm ²	375cm ²	553cm ²	840cm ²	820cm ²	1875cm ²	1400cm ²	-----
Polycarbonate:							
•	•	•	•	•	•	•	---
Polysulfone:							
•	•	•	•	•	•	•	----

Stainless Steel Wire Covers



Type I	Type I L	Type II	Type II L	Type III	Type III H	Type IV	Type IVS
Stainless steel:							
•	•	•	•	•	•	•	•
Wire lid with integrated divider and drinking hole protection plate:							
X	X	•	X	•	X	X	X

Bottles




200ml	250ml	300ml	400ml	450ml	500ml	750ml
Polycarbonate:						
•	•	•	•	X	•	•
Polysulphone:						
•	•	X	X	•	•	•
Suitable for Neoprenstopper						
•	•	•	•	•	•	•
Suitable for Stainless Steel Drinking cap						
•	•	X	•	•	•	•

Racks for cages




Standard Rack	Light Rack	Rack with wire shelves	Rack with full shelves
Stainless steel:			
AISI 316	AISI 316	AISI 304	AISI 304
Dimension:		According to the request	
Racks available for all type of cages			

NexGen IVC for Rats



NexGen IVC for Mice



IVC for Rats	IVC for Mice
From 30 up to 70 cages	From 48 up to 162 cages
Single or double sided	
Factory certified HEPA blower assembly	
Polysulfone plastic	
Autoclavable	

- Advantages:**
- Smart-engineered to be stronger while needing less structural support
 - Rivet technology coupled with high-strength, high-grade autoclavable plastic components
 - Streamlined assembly process improves cost and delivery times
 - Up to 25% lighter than other IVCs
 - No loose parts, easy prep for washing and autoclaving
 - Easy assembled caging systems with increased visibility