


RoVent® Jr.

Small animal automated ventilator

ANIMA LAB
animal facility and laboratory equipment • animal research models



Features & Benefits:

- **Free Offer:**
RightTemp® temperature monitor and homeothermic control module
Warm and auto-regulate the animal's body temperature
- **Fully automatic with touchscreen control**
Safely ventilate using auto-calculated settings based on animal weight
- **Ventilate animals from 3g to 1,250g**
Ideal for mice, rats and small animals
- **Volume control and pressure control**
Reduce footprint with one, dual control ventilator
- **Multi-animal capability***
Greater efficiency

**Our clinical department will be pleased to review your study protocol to confirm that our product will meet your specific study requirements.*

1-Year 100% Parts & Labor Warranty (3-Year extended warranty available)

The RoVent Jr. small animal automated ventilator is ideal for mice, rats and small animals up to 1,250 grams. RoVent Jr. is fully automatic – simply enter the animal's weight and press run.

kentscientific.com

Order your RoVent Jr.:
www.animalab.eu | info@animalab.eu

 **Kent Scientific**
CORPORATION

Choose the ventilator that fits your research needs:

Comparison Chart

RoVent® Jr. & PhysioSuite® with RoVent Module



	RoVent Jr.	PhysioSuite with RoVent Module
Volume ventilation	●	●
Pressure ventilation	●	●
Pulse ventilation mode	○	●
PEEP	○	●
Assist mode	○	●
Sigh breath	○	●
I:E ratio	●	●
Tidal volume — measured	●	●
Auto set-up	Limited	●
User profiles	○	●
Alarms (parameter)	○	●
Trigger in	○	●
Analog out	○	●
Digital out	○	●
Sync out	○	●
ADI integration	○	●
Data output	○	●
Internal data logging	○	●
Realtime clock	○	●
Temperature monitor and homeothermic control module	FREE	FREE