



Sychem is a leading expert, provider of disinfection, sterilization & decontamination solutions for animal facilities. With over 40 years of experience, manufacture wide range of workstations dedicated for animal facilities.

## SychemCHANGE - cage changing station



- ISO class 3 airflow
- H14 HEPA filters
- User-friendly touch screen interface
- Chemical and impact resistant worktop made of TRESPA®

| Models                 | Working area dimensions [mm] |       |
|------------------------|------------------------------|-------|
|                        | Width                        | Depth |
| SychemCHANGE level     | 1208                         | 666   |
| SychemCHANGE dipped    | 1208                         | 712   |
| SychemCHANGE XL level  | 1508                         | 666   |
| SychemCHANGE XL dipped | 1508                         | 712   |

- Ergonomic, electrically driven legs

## SychemSAFETY - class 2 microbiological safety cabinets




| Models           | Width [mm] |
|------------------|------------|
| SychemSAFETY-3FT | 900        |
| SychemSAFETY-4FT | 1200       |
| SychemSAFETY-5FT | 1500       |
| SychemSAFETY-6FT | 1800       |

- Active Airflow Technology
- Airflow sensors
- Maximum 47 dB(A) noise level
- Efficient HEPA filters
- High-quality customised options to suit individual requirements

## SychemFILL - bedding dispensers

- Unparalleled dispensing accuracy with a tolerance of  $\pm 4\%$
- Compatible with all free-flowing rodent bedding, with sizes ranging from 0,9–10 mm, as well as all commonly used wood and paper-based bedding (including corn cob).
- Custom-built, automatic bedding filling station – SychemFILL XL



## SychemWASTE - bedding disposal station

- Triple stage filtration allows quiet operating
- Double chute & sliding worksurface
- Available as SychemWASTE XL with macerator and conveyor system, to reduce bin change frequency and increases cage throughput.



## SychemAIR - air showers

- Active Airflow Technology with an air velocity of 25 m/s
- Particle Removal Excellence with a total of 48 air nozzles, creating an enveloping airflow
- Three Levels of filtration with a HEPA H14 filter
- Easy Installation Options - SychemAIR can be ordered as pre-built, or in a flatpack configuration, for easy installation in limited access facilities