

Detergents for pharmaceutical industry



RBS detergents offer adequate solutions for the cleaning in place CIP by immersion, circulation or by spraying of internal or external surfaces of tanks and other production equipment of medicines, paramedical products or cosmetics.

RBS cleaning agents allow the cleaning of mixers, reactors, pipes, tablets presses and all types of surfaces in contact with intermediates or end products. They assure the removal of insoluble residues, ointments, oils, greases, chemical solutions, pigments, proteins, petroleum, fragrances, creams, silicone, metal oxides, mineral salts, etc.



RBS products information

Product	Priority Product	MSDS	TDS	Composition Certificate	Conductivity Curve	TOC Test report	HPLC Test Report	ICP Test Report	MAS Test Report	Declaration of Non-Animal and Non-Human Origin	Declaration of non-us of Products of Concerns (POC)
RBS IND 458		✓	✓	✓	✓					✓	✓
RBS IND 705	✓	✓	✓	✓	✓	✓	✓			✓	✓
RBS IND 721		✓	✓	✓	✓					✓	✓
RBS IND 740		✓	✓	✓	✓		✓			✓	✓
RBS IND 500		✓	✓	✓	✓			✓	✓	✓	✓
RBS IND 750		✓	✓	✓	✓					✓	✓
RBS IND 756	✓	✓	✓	✓	✓		✓			✓	✓
RBS IND 826		✓	✓	✓	✓					✓	✓
RBS IND 854		✓	✓	✓	✓					✓	✓
RBS IND 945		✓	✓	✓	✓					✓	✓
RBS IND 950	✓	✓	✓	✓	✓					✓	✓

RBS products compatibility

Surface (**)	Inox	High-density polyethylene	Polypropylene	Polytetrafluoroethylene	Silicone	Glass
RBS IND 458	++	++	++	++	Limited	++
RBS IND 705	++	++	++	++	Limited	++
RBS IND 721	++	++	++	++	Limited	++
RBS IND 740	++	++	++	++	Limited	++
RBS IND 500	++	++	++	++	Limited	++
RBS IND 750	++	++	++	++	Limited	++
RBS IND 756	++	++	++	++	Limited	++
RBS IND 826	++	++	++	++	+	++
RBS IND 854	++	++	++	++	+	++
RBS IND 945	++	++	++	++	Poor	++
RBS IND 950	++	++	++	++	Poor	++

RBS detergents for other industries

- Production of cosmetics
- Biotechnology
- Food & beverage processing
- Electronic industry, metal working
- Glass, high-tech glass
- Optics

RBS detergents packing



4.5 kg



100 g



200 ml



5 L



20 kg



220 kg



1000 L