



WORLD  
PRECISION  
INSTRUMENTS  
*Instrumenting scientific ideas*

**ANIMA LAB**  
animal facility and laboratory equipment • animal research models

# Sterilization Baskets & Trays



# Introduction

WPI offers a comprehensive range of surgical instrument baskets and trays. Our instrument trays and baskets provide a high degree of protection for costly surgical instruments during storage, transportation and sterilization. Trays are available in a variety of sizes and designs that accommodate a wide range of surgical instruments.

We are backed by ISO 9001-2008, ISO 13485, 2003 Certification, CE & FDA. Our baskets are made from high quality stainless steel, providing both strength and durability at an affordable price. They are designed with rounded corners and edges to avoid two common problems:


- Burr-free rounded corners prevent wraps from being torn by sharp edges
- Smooth edges reduce user injury

## Exceptionally robust and resistant to accidental damage

We present a comprehensive range of surgical instrument sterilizing baskets including:

- Side Perforated Baskets
- Wire Mesh Baskets
- Fine Mesh Baskets
- Endoscopy Baskets
- Customized Baskets (for holding scopes)

We guarantee all products against poor craftsmanship and material. Any item found defective will be repaired or replaced free of charge. This guarantee does not cover any product which becomes defective through misuse, negligence or mishandling.



### WORLD PRECISION INSTRUMENTS

# Side Perforated Baskets with Wire Base


The wire base is designed for quickly drying instruments after washing and sterilizing. They are perfect for handling and storage of instruments.

Our side perforated baskets are:


- Manufactured with stainless steel 304, non-magnetic
- Electropolished
- Burr free to protect wraps and prevent injury
- Smooth rounded corners
- Single frame, double frame and detention frame lids available
- Side walls are perforated to protect sharp instrument edges and tips from damage. The hole size is 2.0 mm.
- Meets all the requirements of automatic washing and sterilization for all types of instruments
- Separate lids are available with single frame, double frame and detention shapes
- Available with drop in and folding handles
- Top and bottom frame diameter: 5.0 mm
- Base available in wire mesh with opening squares of 4.2 mm and wire diameter of 1.0 mm
- The images show the smooth round corners (left) and the inner mounted drop in handles (right)



Rounded corners



Drop in Handles mounted on the inside



| Item Code | Sizes in mm |
|-----------|-------------|
| WP-4610   | 240X120X30  |
| WP-4611   | 240X120X50  |
| WP-4612   | 240X120X60  |
| WP-4608   | 240X120X70  |
| WP-4609   | 240X120X100 |
| WP-4607   | 240X250X30  |
| WP-4613   | 240X250X50  |
| WP-4615   | 240X250X60  |
| WP-4616   | 240X250X70  |
| WP-4617   | 240X250X100 |
| WP-4656   | 405X250X30  |
| WP-4657   | 405X250X50  |
| WP-4658   | 405X250X60  |
| WP-4659   | 405X250X70  |
| WP-4660   | 405X250X100 |
| WP-4606   | 480X250X30  |
| WP-4618   | 480X250X50  |

| Item Code | Sizes in mm |
|-----------|-------------|
| WP-4619   | 480X250X60  |
| WP-4620   | 480X250X70  |
| WP-4621   | 480X250X100 |
| WP-4661   | 540X250X30  |
| WP-4662   | 540X250X50  |
| WP-4663   | 540X250X60  |
| WP-4664   | 540X250X70  |
| WP-4665   | 540X250X100 |
| WP-4622   | 300x300x50  |
| WP-4623   | 300x300x60  |
| WP-4624   | 470x400x60  |
| WP-4625   | 480x480x50  |
| WP-4626   | 480x480x60  |
| WP-4628   | 500x300x50  |
| WP-4632   | 550x430x50  |
| WP-4634   | 600x300x60  |


Many more sizes available on request. Minimum buy on special orders - 10 pieces.

# Lids for Side Perforated Baskets with Wire Base

## Lids for perforated baskets (Single Frame)

| Item Code | Sizes in mm |
|-----------|-------------|
| WP-4635   | 240X120     |
| WP-4636   | 240X250     |
| WP-4642   | 405x250     |
| WP-4637   | 480X250     |
| WP-4643   | 540x250     |
| WP-4638   | 300X300     |
| WP-4639   | 480X480     |
| WP-4640   | 500X300     |
| WP-4641   | 600X300     |


Many more sizes available on request. Minimum buy on special orders - 10 pieces.



## Lids for perforated baskets (Double Frame)

| Item Code | Sizes in mm |
|-----------|-------------|
| WP-4635D  | 240X120     |
| WP-4636D  | 240X250     |
| WP-4642D  | 405x250     |
| WP-4637D  | 480X250     |
| WP-4643D  | 540x250     |
| WP-4638D  | 300X300     |
| WP-4639D  | 480X480     |
| WP-4640D  | 500X300     |
| WP-4641D  | 600X300     |


Many more sizes available on request. Minimum buy on special orders - 10 pieces.



## Lids for perforated baskets (Detention Frame)



| Item Code | Sizes in mm |
|-----------|-------------|
| WP-4635DT | 240X120     |
| WP-4636DT | 240X250     |
| WP-4642DT | 405x250     |
| WP-4637DT | 480X250     |
| WP-4643DT | 540x250     |
| WP-4638DT | 300X300     |
| WP-4639DT | 480X480     |
| WP-4640DT | 500X300     |
| WP-4641DT | 600X300     |

Many more sizes available on request. Minimum buy on special orders - 10 pieces.



# Side Perforated Baskets with Flat Base

The flat base is designed for quick drying of the instruments after washing and sterilizing. These baskets are also used for handling and storage of instruments.



| Item Code | Sizes in mm |
|-----------|-------------|
| WP-4545   | 120X250x50  |
| WP-4546   | 120X250x60  |
| WP-4547   | 120X250x70  |
| WP-4548   | 120X250x100 |
| WP-4549   | 240x250x50  |
| WP-4550   | 240x250x60  |
| WP-4551   | 240x250x70  |
| WP-4552   | 240x250x100 |
| WP-4553   | 405X250x50  |
| WP-4554   | 405X250x60  |
| WP-4555   | 405X250x70  |
| WP-4556   | 405X250x100 |
| WP-4557   | 480x250x50  |
| WP-4558   | 480x250x60  |

| Item Code | Sizes in mm |
|-----------|-------------|
| WP-4559   | 480x250x70  |
| WP-4560   | 480x250x100 |
| WP-4561   | 540X250x50  |
| WP-4562   | 540X250x60  |
| WP-4563   | 540X250x70  |
| WP-4564   | 540X250x100 |
| WP-4565   | 480X360x50  |
| WP-4566   | 480X360x60  |
| WP-4567   | 480X360x70  |
| WP-4568   | 480X360x100 |
| WP-4569   | 480x480x50  |
| WP-4570   | 480x480x60  |
| WP-4571   | 480x480x70  |
| WP-4572   | 480x480x100 |

Many more sizes available on request. Minimum buy on special orders - 10 pieces.

# Side Perforated Baskets with Flat Base



- Manufactured with stainless steel 304
- Electropolished
- Burr free to protect wraps and prevent injury
- Smooth rounded corners
- Double frame lids available
- Side walls are perforated to protect sharp instrument edges and tips from damage. The hole size is 2.0 mm.
- Meets all the requirements of automatic washing and sterilization for all types of instruments
- Separate lids are available with single frame, double frame and detention shapes
- Available with drop in and folding handles
- Top and bottom frame diameter: 5.0 mm
- Flat base basket available in welded mesh with opening square 6.0 mm, wire diameter 1.2 mm.
- Stacking frame can also be be available on request.
- The images show the flat base from the inside (upper left), the inner mounted drop handles (upper right), the flat base from the outside (lower left) and the smooth rounded corners (lower right)..



## Lids for perforated baskets (Double Frame)

| Item Code | Sizes in mm |
|-----------|-------------|
| WP-4582   | 120X250     |
| WP-4583   | 240X250     |
| WP-4584   | 405X250     |
| WP-4585   | 405X250     |
| WP-4586   | 540X250     |
| WP-4587   | 480X360     |
| WP-4588   | 480X480     |






Many more sizes available on request. Minimum buy on special orders - 10 pieces.



# Wire Mesh Baskets with Drop In Handles

Wire Mesh Baskets are used for washing, disinfecting, sterilizing, storing and handling surgical instruments.

- Manufactured with stainless steel 304
- Electropolished
- Burr free to protect wraps and prevent injury
- Smooth rounded corners
- Handles are mounted on the inside and drop inside the basket for stacking
- Single frame, double frame and detention frame lids available
- Meets all the requirements of automatic washing and sterilization for all types of instruments
- Top and bottom frame diameter: 5.0 mm
- Wire mesh with mesh openings of 4.5 mm, wire diameter of 1.0 mm
- Baskets are stackable



| Item Code | Sizes in mm |
|-----------|-------------|
| WP-4401   | 105X240X30  |
| WP-4402   | 105X240X50  |
| WP-4403   | 105X240X60  |
| WP-4404   | 105X240X70  |
| WP-4405   | 105X240X100 |
| WP-4406   | 240X240X30  |
| WP-4407   | 240X240X50  |
| WP-4408   | 240X240X60  |
| WP-4409   | 240X240X70  |
| WP-4410   | 240X240X100 |
| WP-4411   | 240X250X30  |
| WP-4412   | 240X250X50  |
| WP-4413   | 240X250X60  |
| WP-4414   | 240X250X70  |
| WP-4415   | 240X250X100 |

| Item Code | Sizes in mm |
|-----------|-------------|
| WP-4416   | 250X250X30  |
| WP-4417   | 250X250X50  |
| WP-4418   | 250X250X60  |
| WP-4419   | 250X250X70  |
| WP-4420   | 250X250X100 |
| WP-4421   | 405X250X30  |
| WP-4422   | 405X250X50  |
| WP-4423   | 405X250X60  |
| WP-4424   | 405X250X70  |
| WP-4425   | 405X250X100 |
| WP-4426   | 480X240X30  |
| WP-4427   | 480X240X50  |
| WP-4428   | 480X240X60  |
| WP-4429   | 480X240X70  |
| WP-4430   | 480X240X100 |

| Item Code | Sizes in mm |
|-----------|-------------|
| WP-4431   | 480X250X30  |
| WP-4432   | 480X250X50  |
| WP-4433   | 480X250X60  |
| WP-4434   | 480X250X70  |
| WP-4435   | 480X250X100 |
| WP-4436   | 480X480X30  |
| WP-4437   | 480X480X50  |
| WP-4438   | 480X480X60  |
| WP-4439   | 480X480X70  |
| WP-4440   | 480X480X100 |
| WP-4441   | 540X250X30  |
| WP-4442   | 540X250X50  |
| WP-4443   | 540X250X60  |
| WP-4444   | 540X250X70  |
| WP-4445   | 540X250X100 |

Many more sizes available on request. Minimum buy on special orders - 10 pieces.

# Lids for Wire Mesh Baskets

## Lids for Wire Mesh Baskets (Single Frame)

| Item Code | Sizes in mm |
|-----------|-------------|
| WP-4455   | 105X240     |
| WP-4456   | 240X240     |
| WP-4457   | 240X250     |
| WP-4458   | 250X250     |
| WP-4459   | 405X250     |
| WP-4460   | 480X240     |
| WP-4461   | 480X250     |
| WP-4462   | 480X480     |
| WP-4463   | 540X250     |


Many more sizes available on request. Minimum buy on special orders - 10 pieces.



## Lids for Wire Mesh Baskets (Double Frame)

| Item Code | Sizes in mm |
|-----------|-------------|
| WP-4455D  | 105X240     |
| WP-4456D  | 240X240     |
| WP-4457D  | 240X250     |
| WP-4458D  | 250X250     |
| WP-4459D  | 405X250     |
| WP-4460D  | 480X240     |
| WP-4461D  | 480X250     |
| WP-4462D  | 480X480     |
| WP-4463D  | 540X250     |


Many more sizes available on request. Minimum buy on special orders - 10 pieces.



## Lids for Wire Mesh Baskets (Detention Frame)

| Item Code | Sizes in mm |
|-----------|-------------|
| WP-4455DT | 105X240     |
| WP-4456DT | 240X240     |
| WP-4457DT | 240X250     |
| WP-4458DT | 250X250     |
| WP-4459DT | 405X250     |
| WP-4460DT | 480X240     |
| WP-4461DT | 480X250     |
| WP-4462DT | 480X480     |
| WP-4463DT | 540X250     |

Many more sizes available on request. Minimum buy on special orders - 10 pieces.







# Crimped Wire Mesh Baskets with Tilt Handles

Crimped Wire Mesh Baskets are used for washing, disinfecting, sterilizing, storing and handling surgical instruments.


- Manufactured with stainless steel 304
- Electropolished
- Burr free to protect wraps and prevent injury
- Smooth rounded corners
- Tilt handles are mounted on the inside and drop inside the basket for stacking
- Single frame, double frame and detention frame lids available
- Meets all the requirements of automatic washing and sterilization for all types of instruments
- Top and bottom frame diameter: 5.0 mm
- Crimped mesh with mesh openings of 5 mm, wire diameter of 0.9 mm
- Baskets are stackable.



Stackable



Bottom stacking frame



Tilted Handles

| Item Code | Sizes in mm |
|-----------|-------------|
| WP-4401CW | 105X240X30  |
| WP-4402CW | 105X240X50  |
| WP-4403CW | 105X240X60  |
| WP-4404CW | 105X240X70  |
| WP-4405CW | 105X240X100 |
| WP-4406CW | 240X240X30  |
| WP-4407CW | 240X240X50  |
| WP-4408CW | 240X240X60  |
| WP-4409CW | 240X240X70  |
| WP-4410CW | 240X240X100 |
| WP-4411CW | 240X250X30  |
| WP-4412CW | 240X250X50  |
| WP-4413CW | 240X250X60  |
| WP-4414CW | 240X250X70  |
| WP-4415CW | 240X250X100 |

| Item Code | Sizes in mm |
|-----------|-------------|
| WP-4416CW | 250X250X30  |
| WP-4417CW | 250X250X50  |
| WP-4418CW | 250X250X60  |
| WP-4419CW | 250X250X70  |
| WP-4420CW | 250X250X100 |
| WP-4421CW | 405X250X30  |
| WP-4422CW | 405X250X50  |
| WP-4423CW | 405X250X60  |
| WP-4424CW | 405X250X70  |
| WP-4425CW | 405X250X100 |
| WP-4426CW | 480X240X30  |
| WP-4427CW | 480X240X50  |
| WP-4428CW | 480X240X60  |
| WP-4429CW | 480X240X70  |
| WP-4430CW | 480X240X100 |

| Item Code | Sizes in mm |
|-----------|-------------|
| WP-4431CW | 480X250X30  |
| WP-4432CW | 480X250X50  |
| WP-4433CW | 480X250X60  |
| WP-4434CW | 480X250X70  |
| WP-4435CW | 480X250X100 |
| WP-4436CW | 480X480X30  |
| WP-4437CW | 480X480X50  |
| WP-4438CW | 480X480X60  |
| WP-4439CW | 480X480X70  |
| WP-4440CW | 480X480X100 |
| WP-4441CW | 540X250X30  |
| WP-4442CW | 540X250X50  |
| WP-4443CW | 540X250X60  |
| WP-4444CW | 540X250X70  |
| WP-4445CW | 540X250X100 |

Many more sizes available on request. Minimum buy on special orders - 10 pieces.

# Lids for Crimped Wire Mesh Baskets

## Lids for Crimped Wire Mesh Baskets (Single Frame)

| Item Code | Sizes in mm |
|-----------|-------------|
| WP-4455CW | 105X240     |
| WP-4456CW | 240X240     |
| WP-4457CW | 240X250     |
| WP-4458CW | 250X250     |
| WP-4459CW | 405X250     |
| WP-4460CW | 480X240     |
| WP-4461CW | 480X250     |
| WP-4462CW | 480X480     |
| WP-4463CW | 540X250     |

Many more sizes available on request. Minimum buy on special orders - 10 pieces.



## Lid for Crimped Wire Mesh Baskets (Double Frame)

| Item Code  | Sizes in mm |
|------------|-------------|
| WP-4455DCW | 105X240     |
| WP-4456DCW | 240X240     |
| WP-4457DCW | 240X250     |
| WP-4458DCW | 250X250     |
| WP-4459DCW | 405X250     |
| WP-4460DCW | 480X240     |
| WP-4461DCW | 480X250     |
| WP-4462DCW | 480X480     |
| WP-4463DCW | 540X250     |


Many more sizes available on request. Minimum buy on special orders - 10 pieces.



## Lids for Crimped Wire Mesh Baskets (Detention Frame)

| Item Code   | Sizes in mm |
|-------------|-------------|
| WP-4455DTCW | 105X240     |
| WP-4456DTCW | 240X240     |
| WP-4457DTCW | 240X250     |
| WP-4458DTCW | 250X250     |
| WP-4459DTCW | 405X250     |
| WP-4460DTCW | 480X240     |
| WP-4461DTCW | 480X250     |
| WP-4462DTCW | 480X480     |
| WP-4463DTCW | 540X250     |




Many more sizes available on request. Minimum buy on special orders - 10 pieces.



# Classic Crimped Wire Mesh Baskets, Tilted Handles

Classic Crimped Wire Mesh Baskets are used for washing, disinfecting, sterilizing, storing and handling surgical instruments.

- Manufactured with stainless steel 304
- Electropolished
- Burr free to protect wraps and prevent injury
- Smooth rounded corners
- Tilt handles
- Single frame, double frame and detention frame lids available
- Meets all the requirements of automatic washing and sterilization for all types of instruments
- Top and bottom frame diameter: 5.0 mm
- Crimped mesh with mesh openings of 6 mm, wire diameter of 1.0 mm
- Smooth corners and tilt handles



| Item Code | Sizes in mm |
|-----------|-------------|
| WP-4401ER | 105X240X30  |
| WP-4402ER | 105X240X50  |
| WP-4403ER | 105X240X60  |
| WP-4404ER | 105X240X70  |
| WP-4405ER | 105X240X100 |
| WP-4406ER | 240X240X30  |
| WP-4407ER | 240X240X50  |
| WP-4408ER | 240X240X60  |
| WP-4409ER | 240X240X70  |
| WP-4410ER | 240X240X100 |
| WP-4411ER | 240X250X30  |
| WP-4412ER | 240X250X50  |
| WP-4413ER | 240X250X60  |
| WP-4414ER | 240X250X70  |
| WP-4415ER | 240X250X100 |

| Item Code | Sizes in mm |
|-----------|-------------|
| WP-4416ER | 250X250X30  |
| WP-4417ER | 250X250X50  |
| WP-4418ER | 250X250X60  |
| WP-4419ER | 250X250X70  |
| WP-4420ER | 250X250X100 |
| WP-4421ER | 405X250X30  |
| WP-4422ER | 405X250X50  |
| WP-4423ER | 405X250X60  |
| WP-4424ER | 405X250X70  |
| WP-4425ER | 405X250X100 |
| WP-4426ER | 480X240X30  |
| WP-4427ER | 480X240X50  |
| WP-4428ER | 480X240X60  |
| WP-4429ER | 480X240X70  |
| WP-4430ER | 480X240X100 |


| Item Code | Sizes in mm |
|-----------|-------------|
| WP-4431ER | 480X250X30  |
| WP-4432ER | 480X250X50  |
| WP-4433ER | 480X250X60  |
| WP-4434ER | 480X250X70  |
| WP-4435ER | 480X250X100 |
| WP-4436ER | 480X480X30  |
| WP-4437ER | 480X480X50  |
| WP-4438ER | 480X480X60  |
| WP-4439ER | 480X480X70  |
| WP-4440ER | 480X480X100 |
| WP-4441ER | 540X250X30  |
| WP-4442ER | 540X250X50  |
| WP-4443ER | 540X250X60  |
| WP-4444ER | 540X250X70  |
| WP-4445ER | 540X250X100 |

Many more sizes available on request. Minimum buy on special orders - 10 pieces.

# Fine Mesh Baskets With Lids

Fine Mesh Baskets With Lids are used for washing, sterilizing and storing micro surgical instruments.

- Manufactured with stainless steel 304
- Electropolished
- Burr free to protect wraps and prevent injury
- Smooth rounded corners
- Top and lid frame wire diameter: 3.5 mm and 4.0 mm (depending on the basket)
- Fine mesh with opening squares of 1.6 mm and wire diameter of 0.7 mm
- Wire locking clip



Many more sizes available on request.  
Minimum buy on special orders - 10 pieces.

| Item Code | Sizes in mm |
|-----------|-------------|
| WP-4343   | 40X40x20    |
| WP-4344   | 70x70x30    |
| WP-4345   | 80x40x20    |
| WP-4346   | 80x40x30    |
| WP-4347   | 80x80x34    |
| WP-4348   | 80x80x40    |
| WP-4349   | 90x90x50    |
| WP-4350   | 150x110x20  |
| WP-4351   | 160x40x20   |
| WP-4352   | 105x70x25   |
| WP-4353   | 141x50x40   |
| WP-4354   | 150x50x35   |
| WP-4355   | 180x130x20  |
| WP-4356   | 190x85x100  |
| WP-4357   | 200x100x35  |


| Item Code | Sizes in mm |
|-----------|-------------|
| WP-4358   | 220x140x50  |
| WP-4359   | 210x170x20  |
| WP-4360   | 235x155x40  |
| WP-4361   | 240x90x40   |
| WP-4362   | 240x120x35  |
| WP-4363   | 275x175x35  |
| WP-4364   | 280x170x60  |
| WP-4365   | 250x200x20  |
| WP-4366   | 300x180x20  |
| WP-4367   | 200x110x40  |
| WP-4368   | 225x110x50  |
| WP-4369   | 235x120x60  |
| WP-4370   | 270x135x40  |
| WP-4371   | 280x145x50  |
| WP-4372   | 400x250x50  |

# Micro Fine Mesh Baskets with Lids

Micro Fine Baskets are used for washing, sterilizing and storing delicate micro surgical instruments and parts.

- Manufactured with stainless steel 304
- Electrolyte polished
- Top and lid frame wire diameter: 3.5 mm
- Micro fine mesh with opening squares of 1.0 mm and wire diameter of 0.5 mm
- Wire locking clip

| Item Code  | Sizes in mm |
|------------|-------------|
| WP-4343MM  | 40X40x20    |
| WP-4344MM  | 70x70x30    |
| WP-4345MM  | 80x40x20    |
| WP-4346MM  | 80x40x30    |
| WP-4347MM  | 80x80x34    |
| WP-4348MM  | 80x80x40    |
| WP-4349MM  | 90x90x50    |
| WP-4350MM  | 150x110x20  |
| WP-4351 MM | 160x40x20   |
| WP-4352MM  | 105x70x25   |
| WP-4353MM  | 141x50x40   |
| WP-4354MM  | 150x50x35   |
| WP-4355MM  | 180x130x20  |




**WORLD PRECISION INSTRUMENTS**

# Full Perforated Trays with Lids

Micro Fine Baskets are used for washing, sterilizing and storing delicate micro surgical instruments and parts.

- Manufactured with stainless steel 304
- Electropolished
- Smooth rounded corners
- Perforated hole size of 2.0 mm
- Hinged lids
- Top and lid frame rod diameter: 3.5 mm and 4.0 mm (depending on the basket)
- Full wall perforation to protect sharp instruments' edges and tips from damage
- Fine mesh with opening squares of 1.6 mm and wire diameter of 0.7 mm
- Sheet locking clip




| Item Code | Sizes in mm |
|-----------|-------------|
| WP-2343   | 40X40X20    |
| WP-2344   | 70x70x30    |
| WP-2345   | 80x40x20    |
| WP-2346   | 80x40x30    |
| WP-2347   | 80x80x34    |
| WP-2348   | 80x80x40    |
| WP-2349   | 90x90x50    |
| WP-2350   | 150x110x20  |
| WP-2351   | 160x40x20   |
| WP-2352   | 105x40x20   |
| WP-2353   | 141x50x40   |
| WP-2354   | 150x50x35   |
| WP-2355   | 190x85x35   |
| WP-2356   | 200x100x35  |
| WP-2357   | 240x90x40   |
| WP-2358   | 200x110x40  |
| WP-2359   | 225x110x50  |

# Endoscope Sterilization Tray with Lens Basket

Endoscope Baskets provide security for costly endoscopes during washing and sterilization processes. The silicone padding also provide protection for the endoscope during handling and storage.



| Item Code | Sizes in mm |
|-----------|-------------|
| WP-4046   | 640x150x77  |

Small Lens Baskets    Locking Clip





## Arthroscope Baskets

| Item Code | Sizes in mm |
|-----------|-------------|
| WP-4040   | 180X80X52   |
| WP-4041   | 290X80X52   |
| WP-4042   | 290X80X52   |



## Endoscope Baskets




| Item Code | Sizes in mm |
|-----------|-------------|
| WP-4044   | 460x80x52   |
| WP-4043   | 460x80x52   |

# Endoscope Baskets


Endoscope Baskets provide security for costly endoscopes during washing and sterilization processes. The silicone padding also provide protection for the endoscope during handling and storage.

- Safely sterilize and store expensive endoscopes
- Silicone pads to protect endoscope during handling and storage
- Electropolished
- Burr free to prevent injury
- Locking clip
- Variety of sizes
- Top lid and frame wire diameter of 5.0 mm
- Loop wire diameter of 2.0 mm
- Autoclavable silicone inserts




| Item Code | Sizes in mm |
|-----------|-------------|
|-----------|-------------|

|         |           |
|---------|-----------|
| WP-4045 | 670x80x52 |
| WP-4047 | 600x80x52 |
| WP-4048 | 700x80x52 |



| Item Code | Sizes in mm |
|-----------|-------------|
|-----------|-------------|

|         |            |
|---------|------------|
| WP-4049 | 290x160x52 |
|---------|------------|




| Item Code | Sizes in mm |
|-----------|-------------|
|-----------|-------------|

|         |            |
|---------|------------|
| WP-4050 | 290x160x52 |
|---------|------------|



| Item Code | Sizes in mm |
|-----------|-------------|
|-----------|-------------|

|         |            |
|---------|------------|
| WP-4051 | 460x160x52 |
|---------|------------|




| Item Code | Sizes in mm |
|-----------|-------------|
|-----------|-------------|

|         |            |
|---------|------------|
| WP-4052 | 460x160x52 |
|---------|------------|




| Item Code | Sizes in mm |
|-----------|-------------|
|-----------|-------------|

|         |            |
|---------|------------|
| WP-4053 | 460x160x52 |
|---------|------------|




# Endoscope Baskets

Endoscope Baskets provide security for costly endoscopes during washing and sterilization processes. The silicone padding also provide protection for the endoscope during handling and storage.




| Item Code | Sizes in mm |
|-----------|-------------|
| WP-4054   | 310x150x77  |



| Item Code | Sizes in mm |
|-----------|-------------|
| WP-4066   | 460x225x52  |

| Item Code | Sizes in mm |
|-----------|-------------|
| WP-4067   | 460x225x52  |



| Item Code | Sizes in mm |
|-----------|-------------|
| WP-4068   | 580x160x52  |

| Item Code | Sizes in mm |
|-----------|-------------|
| WP-4069   | 580x160x52  |

## WORLD PRECISION INSTRUMENTS

UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com Germany: +49 (0)6031 1602171 • wptide@wpi-europe.com • www.wpi-europe.com  
 Brazil: (013) 406-29703 • info@brazil.wpiinc.com • https://www.wpi brasil.com.br China: +86 21 6888 5517 • chinasales@china.wpiinc.com • www.wpiinc.cn



# Micro Ophthalmic Baskets with Lids

Micro ophthalmic baskets are specially designed for delicate micro surgical instruments.

- Manufactured with stainless steel 304
- Electrolyte polished
- Silicone inserts on the bottom of the baskets accommodates micro instruments
- Lids have silicone strips to grip the instruments
- Side walls are constructed of perforated sheets
- Easy to open locking handle
- Designed for delicate micro instruments.



| Item Code | Sizes in mm |
|-----------|-------------|
| WP-22000  | 240x250x38  |
| WP-22001  | 480x250x38  |
| WP-22002  | 400x250x38  |



# Instrument Cassettes with Spring Holding Lids

These cassettes are specially designed to minimize the contact area when instruments are inside the cassette. Each instrument is held by a separate spring. The cassette holds thick and thin instruments securely in a single cassette arrangement.

| Item Code | Sizes in mm |
|-----------|-------------|
| WP-4242   | 178X64x32   |
| WP-4243   | 178X114x32  |
| WP-4244   | 203X140x38  |
| WP-4245   | 203x280x38  |
| WP-4246   | 220x125x40  |
| WP-4247   | 250x220x40  |



# Fine Mesh Cassettes for Micro Instruments


These cassettes are specially designed to micro instruments. They are perfect for cleaning instruments, especially when the water spray comes from any angle.

| Item Code | Sizes in mm |
|-----------|-------------|
| WP-4148   | 240X90x40   |
| WP-4149   | 235x80x30   |
| WP-4150   | 235x155x40  |



# Instruments Cassettes

| Item Code | Sizes in mm |
|-----------|-------------|
| WP-3010   | 180X60X22   |
| WP-3011   | 180X125X22  |



| Item Code | Sizes in mm |
|-----------|-------------|
| WP-3012   | 184x145x34  |
| WP-3013   | 184x74x34   |
| WP-3014   | 184x74x34   |




| Item Code | Sizes in mm |
|-----------|-------------|
| WP-3015   | 205x143x34  |




| Item Code | Sizes in mm |
|-----------|-------------|
| WP-3016   | 280x183x32  |

For 16 Instruments with partition plate



| Item Code | Sizes in mm |
|-----------|-------------|
| WP-3017   | 280x183x32  |

For 20 Instruments




| Item Code | Sizes in mm |
|-----------|-------------|
| WP-3018   | 280x205x32  |



## WORLD PRECISION INSTRUMENTS


| Item Code | Sizes in mm |
|-----------|-------------|
| WP-3019   | 280x183x32  |



## Micro Cassettes

Micro cassettes are great for sterilization of micro surgical instruments and small clamps.

| Item Code | Sizes in mm |
|-----------|-------------|
| WP-3040   | 40x55x20    |
| WP-3041   | 74x45x20    |
| WP-3042   | 74x59x24    |



## Clip-In

| Item Code | Sizes in mm |
|-----------|-------------|
| WP-3051   | 185x82x18   |
| WP-3052   | 185x140x18  |



## Fixing Clamp


- Adjustable clamp length
- Universally adaptable for all sizes of surgical instruments
- Easily dismantled for cleaning purposes.

| Item Code | Sizes in mm   |
|-----------|---------------|
| WP-4141   | 60/220x220x39 |
| WP-4142   | 60x220x39     |





WORLD  
PRECISION  
INSTRUMENTS  
*Instrumenting scientific ideas*



Contact our distributor

**ANIMA LAB**  
animal facility and laboratory equipment • animal research models

[info@animalab.eu](mailto:info@animalab.eu)